

# **Bio Signal Total Solution**

Fetal Monitor
Fetal & Vascular Doppler
Patient Monitor
Infant Incubator
Infant Warmer
Phototherapy
Breast Pump
Thermometer
Head Lamp



# **Fetal Monitor**



# **BT-380 Fetal & Maternal Monitor I**

- 1MHz Frequency with Twin Doppler
- 10.1" Color Touch Screen
- Fetal Monitoring : FHR, UC, DECG Maternal Monitoring : ECG, NIBP, SpO2, TEMP
- Option: DECG, Rechargeable battery, Acoustic Stimulator, Wifi, Cart

# **BT-350 LCD Fetal Monitor**

- 1MHz Frequency with Twin Doppler
- 7" color TFT LCD Display
- CTG Analysis
- Option: Rechargeable battery, Acoustic Stimulator, Carrying bag, Cart





# **BT-350 LED Fetal Monitor I**

- 1MHz Frequency with Twin Doppler
- Large 7 Segment LED Display
- Waterproof Probe (IPX 8)
- Option : Rechargeable battery, Acoustic Stimulator, Carrying bag, Cart

# I BT-300 LED Fetal Monitor

- 1MHz Frequency with Twin Doppler
- Waterproof Probe (IPX 8)
- Option : Rechargeable battery, Acoustic Stimulator, Carrying bag, Cart



# **Fetal & Vascular Doppler**

# HBT-200 Fetal Doppler

- 2/3 MHz Frequency
- Various Display Type (Sound only, LCD, Color)
- Waterproof Probe
- FHR and Sound management





# **BT-220 Fetal Doppler**

- 2/3 MHz Frequency
- 2.4" Color LCD Display with FHR Trend Display
- Acoustic Stimulator Function
- BMI and Body fat Analysis



- 2/3 MHz Frequency
- 3.2" color TFT LCD Display with FHR Trend Display
- Rechargeable Battery





# **BT-200 Vascular Doppler 1**

- 2, 4, 5, 8 MHz Frequency
- LCD Display
- Interchangeable Waterproof Probe
- Low Battery Indicator

# **Patient Monitor**

# **BT-780 Patient Monitor**

- ECG, Resp., SpO2, NIBP, Temp., CO<sub>2</sub>, IBP, Multi-gas, C.O.
- 15.6" Color Touch Screen
- Ultra slim design
- 5 hours operation on battery
- Smart Hook & stand
- HL 7 support
- Option : CO<sub>2</sub>, IBP, Multi-gas, C.O., Printer, Masimo SpO2, WIFI



# **BT-770 Patient Monitor**

- ECG, Resp., SpO2, NIBP, Temp., CO<sub>2</sub>, IBP, Multi-gas, C.O.
- 12.1" Color Touch Screen
- Ultra slim design
- 5 hours operation on battery
- Smart Hook & stand
- ■HL 7 support
- Option : CO<sub>2</sub>, IBP, Multi-gas, C.O., Printer, Masimo SpO2, WIFI





# **BT-740 Patient Monitor**

- ECG, Resp., SpO2, NIBP, Temp., CO<sub>2</sub>, IBP, Multi-gas, C.O.
- 8.4" Color Touch Screen
- Ultra slim design
- 6 hours operation on battery
- Smart Hook & stand
- HL 7 support
- Option : CO<sub>2</sub>, IBP, Multi-gas, CO, Printer, Masimo SpO2, WIFI

# **Patient Monitor**

# **BT-720 Vital Sign Monitor**

- SpO2, Pulse
- 4.3" Color Touch Screen
- 8 hours operation on battery
- Fan-less design
- Rechargeable by USB charger
- Option: NIBP, Temp





# **BT-710 Pulse Oximeter H**

- SpO2, Pulse
- Hand-held design
- 4.3" Color Touch Screen
- ■8 hours operation on battery
- Fan-less design
- Rechargeable by USB charger
- Option: CO<sub>2</sub>



# **BT-750 Patient Monitor**

- ECG, Resp., SpO2, NIBP, Temp., CO<sub>2</sub>, IBP
- 10.4" Color LCD Display
- ■4 hours operation on battery
- Multi language support
- Option : CO<sub>2</sub>, IBP, Printer, LAN, Cart, Wall mount

# **Neonate Care**



# **BT-500 Infant Incubator** 4

- Air / Skin Temperature servo control
- Humidity Servo Control
- Tilting bed
- Option : O2 Servo Control, O2 Monitoring, Masimo SpO2 & CCD Camera & External monitor, Lifting Stand, Weighing scale, Basket, IV Plate, IV Pole, Shelf

# l BT-550 Infant Warmer

- Infrared Heater
- Pre-warm / Baby / Manual Mode
- APGAR Timer
- Option : Tilting, Masimo SpO2, Lifting Stand, Weighing scale, Basket, IV Plate, IV Pole, Oxygen delivery, Suction





# **BT-400 Phototherapy**

- Blue LED Lamp
- Intensity Control (High/Low)
- 100,000 hours LED lifetime
- Timer function

# **BT-450 Phototherapy**

- Blue LED (455 nm to 465 nm)
- Intensity Control (High/Low)
- Blanket Type
- 10 hours operation on battery
- Alarm for High Temperature





# **Others**

# + BT-100 Breast Pump

- Portable Type with Dual Pumping
- Massage / Expression / Memory Mode
- Backflow Prevention
- BPA / Latex Free





# **BT-150 Breast Pump**

- Desktop Type with Dual Pumping
- Massage / Expression / Program Mode
- Various Vacuum Level and Cycle Control
- Option: Carrying bag, Car charging cable,
   Nursing Night Lamp, Additional bottle set,
   Rechargeable battery, Bluetooth for mobile app

# **4 Non-Contact Thermometer**

- Fever Alert Function (Blue/Red Color)
- Measurement of Body and Object
- Protective Cap for Sensor (BT-35)
- Faster Measurement Time (BT-36)





# **BT-410 Head Lamp I**

- Power LED Light with 6,000K Color Temperature
- Adjustable illumination (BT-410F)
- Adjustable Light Spot (BT-410A)
- Option: Loupe, Astral Lamp, Additional battery

# **BIO SIGNAL TOTAL SOLUTION**

### Bistos Co., Ltd. (Headquarter)

7th Fl., A Bldg., Woolim Lions Valley 5-cha, 302, Galmachi-ro, Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea

**Tel.** +82 (31) 750 0340 Fax. +82 (31) 750 0344



bistos@bistos.co.kr www.bistos.co.kr

## Bistos America Inc. (US Branch)

Fax. +1-949-614-8745 Site, www.bistosamerica.com **E-mail.** info@bistosamericainc.com

# Specifications: BT-770 12.1" Multi-parameter Patient Monitor

Functional Characteristics			
Display			
Туре	Color TFT touch screen LCD		
Size and resolution	12.1", 800 x 600 pixels		
LED	12.17, 000 X 000 pixels		
Alarm indicator	Yellow & red		
Adaptor power indicator	1 green		
Battery status indicator	1 green		
Audio	, g.co		
	Alarm sound (45 ~ 85dB), key pressing sound		
Speaker	QRS sound, PR sound		
Speaker .	Alarm sound meets the IEC60601-1-8		
Data Storage	Additi sound freets the recood 1 1 0		
Trend	168hours, resolution : 1min		
Alarm event	200 physiological and 100 technical alarm events		
NiBp measurement result	1,000 groups		
Function	1,000 groups		
	English, France, Spanish, Turkey		
Multi-language Trend	. ,		
Alarm	Graphic/tabular		
-	Viewel andible information and the Control		
Mode	Visual, audible, information, parameter flashing		
Alarm delay	Off, 1s, 2s, 3s, 4s, 5s, 6s, 7s, 8s		
Pause duration	1, 2, 3, 4, 5, 10, 15min or permanent		
System	Low battery		
Interface			
Auxiliary	Nurse call		
RJ45 (LAN)	CMS		
USB	S/W upgrade		
ECG			
Standard compliance	IEC60601-2-27		
Lead type	3Lead : I, II, III		
Lead type	5Lead : I, II, III, aVR, aVL, aVF, V		
Display sensitivity (gain)	Auto, 1.25, 2.5, 5, 10, 20mm/mV		
Wave sweep speed	12.5, 25, 50mm/s		
	Diagnostic mode : 0.05 ~ 130Hz		
Band width	Monitoring mode : 0.5 ~ 40Hz		
band width	Surgery mode : 1 ~ 25Hz		
	Surgery mode . 1 ZSMZ		
	Strong filter mode : 5 ~ 20Hz		
CMRR	_ ·		
CMRR Notch	Strong filter mode : 5 ~ 20Hz		
	Strong filter mode : 5 ~ 20Hz > 100dB		
Notch Differential input	Strong filter mode : 5 ~ 20Hz > 100dB 50/60Hz (can be set on or off)		
Notch Differential input Electrode polarization voltage range	Strong filter mode : 5 $\sim$ 20Hz > 100dB 50/60Hz (can be set on or off) > 5M $\Omega$ ±400mV		
Notch Differential input Electrode polarization voltage range Baseline recovery time	Strong filter mode : 5 ~ 20Hz > 100dB 50/60Hz (can be set on or off) > 5MΩ ±400mV < 5s after defibrillation (monitor and surgery mode)		
Notch Differential input Electrode polarization voltage range Baseline recovery time Calibration signal	Strong filter mode : 5 ~ 20Hz > 100dB 50/60Hz (can be set on or off) > 5MΩ ±400mV < 5s after defibrillation (monitor and surgery mode) 1mV (peak-peak), accuracy ±3%		
Notch Differential input Electrode polarization voltage range Baseline recovery time	Strong filter mode : 5 ~ 20Hz > 100dB 50/60Hz (can be set on or off) > 5MΩ ±400mV < 5s after defibrillation (monitor and surgery mode) 1mV (peak-peak), accuracy ±3% Measuring electrode : < 0.1μA		
Notch Differential input Electrode polarization voltage range Baseline recovery time Calibration signal Lead-off detection current	Strong filter mode : 5 ~ 20Hz  > 100dB  50/60Hz (can be set on or off)  > 5MΩ  ±400mV  < 5s after defibrillation (monitor and surgery mode)  1mV (peak-peak), accuracy ±3%  Measuring electrode : < 0.1μA  Drive electrode : < 1μA		
Notch Differential input Electrode polarization voltage range Baseline recovery time Calibration signal	Strong filter mode : 5 ~ 20Hz  > 100dB  50/60Hz (can be set on or off)  > 5MΩ  ±400mV  < 5s after defibrillation (monitor and surgery mode)  1mV (peak-peak), accuracy ±3%  Measuring electrode : < 0.1μA  Drive electrode : < 1μA  Adult : 15 ~ 300bpm		
Notch Differential input Electrode polarization voltage range Baseline recovery time Calibration signal Lead-off detection current HR measuring range	Strong filter mode : 5 ~ 20Hz > 100dB 50/60Hz (can be set on or off) > 5MΩ ±400mV < 5s after defibrillation (monitor and surgery mode) 1mV (peak-peak), accuracy ±3% Measuring electrode : < 0.1μA Drive electrode : < 1μA Adult : 15 ~ 300bpm Pediatric/Neonate : 15 ~ 350bpm		
Notch Differential input Electrode polarization voltage range Baseline recovery time Calibration signal Lead-off detection current HR measuring range HR measuring resolution	Strong filter mode : 5 ~ 20Hz  > 100dB  50/60Hz (can be set on or off)  > 5MΩ  ±400mV  < 5s after defibrillation (monitor and surgery mode)  1mV (peak-peak), accuracy ±3%  Measuring electrode : < 0.1μA  Drive electrode : < 1μA  Adult : 15 ~ 300bpm  Pediatric/Neonate : 15 ~ 350bpm  1bpm		
Notch Differential input Electrode polarization voltage range Baseline recovery time Calibration signal Lead-off detection current HR measuring range	Strong filter mode : 5 ~ 20Hz  > 100dB  50/60Hz (can be set on or off)  > 5MΩ  ±400mV  < 5s after defibrillation (monitor and surgery mode)  1mV (peak-peak), accuracy ±3%  Measuring electrode : < 0.1μA  Drive electrode : < 1μA  Adult : 15 ~ 300bpm  Pediatric/Neonate : 15 ~ 350bpm  1bpm  ±1bpm or ±1%, whichever is greater		
Notch Differential input Electrode polarization voltage range Baseline recovery time Calibration signal Lead-off detection current HR measuring range HR measuring resolution	Strong filter mode : 5 ~ 20Hz  > 100dB  50/60Hz (can be set on or off)  > 5MΩ  ±400mV  < 5s after defibrillation (monitor and surgery mode)  1mV (peak-peak), accuracy ±3%  Measuring electrode : < 0.1μA  Drive electrode : < 1μA  Adult : 15 ~ 300bpm  Pediatric/Neonate : 15 ~ 350bpm  1bpm  ±1bpm or ±1%, whichever is greater  Ventricular bigeminy : 80±1bpm		
Notch Differential input Electrode polarization voltage range Baseline recovery time Calibration signal Lead-off detection current HR measuring range HR measuring resolution	Strong filter mode : 5 ~ 20Hz  > 100dB  50/60Hz (can be set on or off)  > 5MΩ  ±400mV  < 5s after defibrillation (monitor and surgery mode)  1mV (peak-peak), accuracy ±3%  Measuring electrode : < 0.1μA  Drive electrode : < 1μA  Adult : 15 ~ 300bpm  Pediatric/Neonate : 15 ~ 350bpm  1bpm  ±1bpm or ±1%, whichever is greater  Ventricular bigeminy : 80±1bpm  Slow alternating ventricular bigeminy : 60±1bpm		
Notch Differential input Electrode polarization voltage range Baseline recovery time Calibration signal Lead-off detection current HR measuring range HR measuring resolution HR measurement accuracy	Strong filter mode : 5 ~ 20Hz  > 100dB  50/60Hz (can be set on or off)  > 5MΩ  ±400mV  < 5s after defibrillation (monitor and surgery mode)  1mV (peak-peak), accuracy ±3%  Measuring electrode : < 0.1μA  Drive electrode : < 1μA  Adult : 15 ~ 300bpm  Pediatric/Neonate : 15 ~ 350bpm  1bpm  ±1bpm or ±1%, whichever is greater  Ventricular bigeminy : 80±1bpm  Slow alternating ventricular bigeminy : 60±1bpm  Rapid alternating ventricular bigeminy : 120±1bpm		
Notch Differential input Electrode polarization voltage range Baseline recovery time Calibration signal Lead-off detection current HR measuring range HR measuring resolution HR measurement accuracy	Strong filter mode : 5 ~ 20Hz  > 100dB  50/60Hz (can be set on or off)  > 5MΩ  ±400mV  < 5s after defibrillation (monitor and surgery mode)  1mV (peak-peak), accuracy ±3%  Measuring electrode : < 0.1μA  Drive electrode : < 1μA  Adult : 15 ~ 300bpm  Pediatric/Neonate : 15 ~ 350bpm  1bpm  ±1bpm or ±1%, whichever is greater  Ventricular bigeminy : 80±1bpm  Slow alternating ventricular bigeminy : 60±1bpm  Rapid alternating ventricular bigeminy : 120±1bpm  Bidirectional systoles : 90±2bpm		
Notch Differential input Electrode polarization voltage range Baseline recovery time Calibration signal Lead-off detection current HR measuring range HR measuring resolution HR measurement accuracy	Strong filter mode : 5 ~ 20Hz  > 100dB  50/60Hz (can be set on or off)  > 5MΩ  ±400mV  < 5s after defibrillation (monitor and surgery mode)  1mV (peak-peak), accuracy ±3%  Measuring electrode : < 0.1μA  Drive electrode : < 1μA  Adult : 15 ~ 300bpm  Pediatric/Neonate : 15 ~ 350bpm  1bpm  ±1bpm or ±1%, whichever is greater  Ventricular bigeminy : 80±1bpm  Slow alternating ventricular bigeminy : 60±1bpm  Rapid alternating ventricular bigeminy : 120±1bpm		

HR alarm upper limit (bpm)	Adult : 16 ~ 300, 1bpm step  Pediatric/Neonate : 16 ~ 350, 1bpm step		
	Adult : 15 ~ 299, 1bpm step		
HR alarm lower limit (bpm)	Pediatric/Neonate : 15 ~ 349, 1bpm step		
	Detection range: ±2mV ~ ±700mV		
Pacing pulse identification	Pulse width: 0.2ms ~ 2.0ms		
Pacing pulse average HR	15s data		
Pacing pulse interval of HR Refreshing	Every second		
Pacing pulse HR change response time	≤ 10sec		
Pacing pulse tall T-wave suppression	2mV		
3 h	Communication, configuration, selfcheck error		
	Lead off		
	HR high/low, PVCS high		
Alarm	Asystole, VF/VTA, R on T, Tachycardia/bradicardia, PVC		
	frequent/couplet/singlr/bigeminy/trigeminy, Miss Beat		
	Pacemaker not capture/work		
	Signal weak, ST-I, II, II high/low		
Respiration			
Measurement method	Thoracic electrical bio impedance method		
Measuring lead	Lead RA-LA, RA-LL		
Wave gain	X0.5, x1, x2		
Respiratory impedance range	0.2 ~ 3 Ω		
Base line impedance	500 ~ 2,000Ω		
Gain	10 grades		
Wave sweep speed	6.25mm/s, 12.5mm/s, 25mm/s		
Measurement accuracy	±2rpm		
Measurement range	0 ~ 120rpm		
	RR high/low		
Alarm	Apnea		
	Respiration artifact		
Temperature			
Standard compliance	ISO80601-2-56		
Measurement method	Thermistor		
Measuring range	0°C ~ 50.0°C (32°F ~ 122.0°F)		
Resolution	0.1°C		
Measurement accuracy	±0.3℃		
Number of channel	2		
T1/T2 alarm upper limit	0.1°C ~ 50.0°C, 0.1°C/°F step		
T1/T2 alarm lower limit	0°C ~ 49.9°C, 0.1°C/°F step		
Temperature difference alarm upper limit	0°C ~ 50.0°C, 0.1°C/°F step		
Alarm	T1, T2 Sensor off		
NiPa	T1/T2 high/low, TD high		
NiBp Standard compliance	IEC80601-2-30		
Measurement method	Automatic oscillometric method		
Operating mode			
Useful life	Manual, automatic, continuous(STAT) 100,000times		
Measurement interval in automatic mode			
Typical measurement time	1/2/3/4/5/10/15/30/60/90/120/180/240/480min 20~40s		
Typical measurement time	Systolic : Adult(40~270), Pediatric(40~200), Neonate(40~130)		
Normal mode measuring range (mmHg)	Mean : Adult(20~230), Pediatric(20~175), Neonate(20~100)		
Troinial mode measuring range (mining)	Diastolic : Adult(10~210), Pediatric(20~173), Neonate(20~100)		
	Maximum average error: ±5mmHg		
Measurement accuracy	Maximum average error: ±5mmng  Maximum standard deviation: 8mmHg		
Resolution	1mmHg		
Nesolution	Adult: 150(default), 80~240(pressure setting range)		
Initial inflation pressure (mmHg)	Pediatric: 100(default), 80~240(pressure setting range)		
mass imader pressure (imility)	Neonate: 100(default), 60~120(pressure setting range)		
	resonate. Toolaciaary, oo Tzo(pressure setting range)		

	Adult: 200mmHa		
Overpressure protection point (software)	Adult: 300mmHg		
Overpressure protection point (software)	Pediatric: 240mmHg		
	Neonate: 150mmHg		
Overpressure protection point (hardware)	Adult: 320~330mmHg		
Overpressure protection point (hardware)	Pediatric: 265~275mmHg		
Ct-ti- Du	Neonate: 160~165mmHg		
Static Pressure accuracy	±3mmHg		
Supply voltage	10V~14VDC		
Maximum power consumption	3.6W		
Quiescent current	50mA		
Maximum current during measurement	180mA		
Maximum current during inflation	300mA		
	Communication, selfcheck, CFG error		
	System error, measurement timeout		
	Cuff loose, no, leak, type error		
Alarm	Air pressure error		
, warm	Over range, signal weak/unstable/saturated		
	Over pressure		
	Module reset failed		
	Systolic, mean, diastolic high/low		
SpO2			
Standard compliance	ISO80601-2-61		
Display range	0% ~ 100%		
SpO2 display resolution	1%		
	Adult/Pediatric: 70 ~ 100% ±2%		
SpO2 accuracy	Neonate: 70 ~ 100% ±3%		
	0 ~ 69% : Unspecified		
Wave sweep speed	12.5mm/s, 25mm/s		
Wave mode	Scan, fill		
Pulse volume	0, 1, 2, 3, 4, 5, 6, 7, 8, 9 level		
i dise volume	Upper Alarm Limit : 86% ~ 100%		
SpO2 alarm preset limits	Lower Alarm Limit: 85% ~ 99%		
CnO2 plarm procet accuracy	±1%		
SpO2 alarm preset accuracy			
SpO2 alerting signal generates delay	No Delay		
SpO2 value refresh period	1s/time		
SpO2 value refresh delay	< 10s		
	Low Sensitivity: 7 ~ 8s		
Average period	Intermediate Sensitivity : 4 ~ 6s		
	Advanced Sensitivity : 2 ~ 3s		
	Low Sensitivity: < 8s		
Alarm condition delay period	Intermediate Sensitivity: < 6s		
	Advanced Sensitivity : < 3s		
Alarm sign generates delay period	Os		
Perfusion index	0.05 ~ 20%		
PR Measurement Range	25 ~ 254bpm		
PR Resolution	±1bpm		
PR Measurement accuracy	±2% or ±2bpm, whichever is greater		
	Communication stop/error		
	No sensor/ sensor off		
Alarm	Search timeout		
	Search pulse(weak)		
	SpO2, RR high/low		
IBP (Option)	-1-1		
Standards compliant	IEC60601-2-34		
Pressure measurement range	-50 ~ 400 mmHg		
Pressure measurement accuracy	±3 mmHg or±2%, whichever is greater		
Pressure resolution	1 mmHg		
PR measurement range	35 ~ 250 bpm		
	23 L30 0pm		

PR measurement accuracy	±3bpm		
PR resolution	1bpm		
Transducer sensitivity	5μV/V/mmHg		
Transducer resistance range	300-5,000Ω		
Supply voltage	+12VDC		
Maximum power consumption	≤5W		
Scan speed	12.5mm/s, 25mm/s		
	IBP1, 2 communication stop/error		
	IBP1, 2 sensor off		
	Art-sys, PA-sys, P1-sys, P2-sys high		
Alarm	Art-dia, PA-dia, P1-dia, P2-dia high		
	Art-mean, PA-mean, CVP-mean, LAP-mean, RAP-mean, ICP-		
	mean, P1-mean, P2-mean high		
EtCO2 Mainstream & Sidestream (Opti			
Measurement parameters	EtCO2、FiCO2、AwRR		
Measuring range	EtCO2: 0~150mmHg, AwRR: 0~150rpm		
Resolution	EtCO2/FiCO2 : 1mmHg, AwRR : 1rpm		
Apnea delay	20s, 25s, 30s, 35s, 40s, 45s, 50s, 55s, 60s		
Operating mode	Standby, measure		
O2 compensation	Low, mid, high		
N2O compensation	On, off		
It - 2005.	EtCO2 lower limit : 0~149mmHg		
	EtCO2/FiCO2 upper limit: 1~150mmHg		
Alarm limit	AWRR lower limit : 0~119rpm		
	AWRR upper limit: 1~120rpm		
	Communication stop/error		
	CO2 sensor off/error		
	O2 sensor error/replace		
	adaptor/sampling line no/check		
	Parameter accuracy error		
	O2, Air calibration error		
	S/W, H/W error		
	Motor accuracy error		
	CO2 factory calibration error		
Alarm	Adaptor, sampling line replace		
	O2 port error		
	CO2, O2, N2O out of accuracy		
	CO2 temp., pressure out of accuracy		
	CO2 zero required		
	CO2 zeroing/sleeping		
	CO2 module calibrating/calibration error		
	EtCO2, FiCO2, AWRR high/low		
	Apnea		
C.O. (Cardiac Output : Option)	Αμιτέα		
( amaz z arpate z pasi,	C.O.: 0.2 ~ 20 L/min		
Measurement range	BT : 23 ~ 45°C±0.5 °C		
- 3 ·	IT : 0 ~ 20°C±0.5 °C		
Resolution factor	C.O. : 0.1L/min		
	BT, IT : 0.1°C		
	C.O.: $\pm 5\%$ or $\pm 0.1$ L/min, subject to the bigger one		
Accuracy	BT, IT: ±0.1°C (sensor exclusive)		
	BT high limit: (Low limit +0.1) ~ 43°C		
Scope of alarm limit	BT low limit: 23.0 ~ (high limit -0.1) °C		
ocepe of diam. In the	Step size : 0.1°C		
	BT sensor off		
Alarm	BT high/low		
/ MOTH	C.O. high		
	c.o. mgn		

Printer (Option)			
Type	Thermal dot array		
Print speed	12.5, 25, 50mm/s		
Paper size	58mm(W) x 42m		
Power			
Adaptor	Input : AC 100 ~ 240V (50/60Hz)		
Adaptor	Output: DC 15V/2.4A		
Consumption	13.5W		
	11.1V Li-ion 4,400mA		
Rechargeable battery	Operating Time : 5hrs		
	Charging Time : 4hrs		
Standard Configurations			
ECG cables and lead wire	1ea(5lead)		
ECG electrode for adult	1pack(25pcs)		
SpO2 adult reusable sensor	1ea		
SpO2 extension cable	1ea		
NiBp adult cuff	1ea		
NiBp extension tube	1ea		
Temperature sensor	1ea		
Power adaptor	1ea		
Bracket	1ea		
Operation manual	1ea		
Options (Function)			
IBP	Sensor cable & package		
EtCO2 Mainstream (Bistos)	Airway adaptor & module		
EtCO2 Sidestream (Bistos)	Sampling tube		
EtCO2 IRMA Mainstream (Masimo)	Airway adaptor & module		
EtCO2 ISA Sidestream (Masimo)	Sampling tube		
C.O.	Sensor cable		
Printer	Printer & paper		
Cart			
Options (Accessory)			
ECG cables and lead wire	5/3 lead		
ECG electrode	adult/neonate		
SpO2 reusable sensor	adult/pediatric/neonate		
SpO2 disposable sensor	adult/pediatric/neonate		
Skin & rectal temperature sensor	adult/pediatric/neonate		
NiBp cuff	adult(27~35cm)/pediatric(14~21.5cm)/neonate(4*9cm)		
Physical Characteristics			
Dimension			
Main unit	320(W) x 65(D) x 250(H)mm		
Packing	400(W) x 350(D) x 290(H)mm		
Weight			
Main unit	< 2.8Kg		
Packing	4.6Kg		
Environmental Conditions			
Operating temperature	10 ~ 40°C (50 ~ 104°F)		
Operating humidity	5 ~ 85% non-condensing		
Storage temperature	−20 ~ 60°C (−4 ~ 140°F)		
Storage humidity	0 ~ 95% non-condensing		
Warranty			
Main unit	2 years		
Optional sensor & accessory	1 year		
Certificates			
KFDA, CE			

# **Specifications: BT-710 Handheld Pulse Oximeter**

Functional Characteristics			
SpO2			
Display range	0% ~ 100%		
SpO2 display resolution	1%		
op o z display reservice.	Adult/Pediatric : 70 ~ 100% ±2%		
SpO2 accuracy	Neonate : 70 ~ 100% ±3%		
op or accaracy	0 ~ 69% : unspecified		
	Upper alarm limit : 86% ~ 100%		
SpO2 alarm preset limits	Lower alarm limit: 85% ~ 99%		
SpO2 alarm preset accuracy	±1%		
SpO2 alerting signal generates delay	No delay		
SpO2 value refresh period	1s/time		
SpO2 value refresh delay	< 10s		
Spoz value refresir delay	Low sensitivity: 7 ~ 8s		
Average period	Intermediate sensitivity: 4 ~ 6s		
Average period	,		
	Advanced sensitivity: 2 ~ 3s		
Alarm condition dolar pariod	Low sensitivity: < 8s		
Alarm condition delay period	Intermediate sensitivity: < 6s		
	Advanced sensitivity: < 3s		
Alarm sign generates delay period	0s		
Pulse Rate	25 250		
Measuring range	25 ~ 250bpm		
Resolution	±1bpm		
Accuracy	±2% or ±2bpm, whichever is greater		
Display			
Type	Color TFT touch screen LCD		
Size	4.3"		
Function			
Sleep mode			
Perfusion index			
Multi-language			
Trend	168hours		
Alarm			
	SpO2 high/low		
SpO2 alarms	PR high/low		
	PI high/low		
	Spo2 sensor no/off/error		
	SpO2 search timeout/pulse		
	SpO2 signal unstable		
System alarms	SpO2 board failure		
System diamis	Low perfusion		
	Too much light		
	Battery low		
	System will shutdown		
PC Interface			
SD card interface	S/W upgrade		
Others			
Liquid Inlet Protection Grade	IPX2		
Power			
	Input : AC 100 ~ 240V (50/60Hz)		
Adaptor	Output : DC 5V/2A		
	3.7V Li-ion 3,800mA		
Rechargeable battery	Operating time: 8hrs		
	Charging time : 4hrs		

Standard Configurations		
Adult SpO2 probe	1ea	
Carrying pouch	1ea	
Operation manual	1ea	
Power adaptor	1ea	
Physical Characteristics		
Dimension		
Main unit	84(W) x 34.5(D) x 158.5(H)mm	
Packing (one unit)	201(W) x 106(D) x 69(H)mm	
Carton box (10ea)	355(W) x 230(D) x 220H)mm	
Weight		
Main unit	300g	
Packing (one unit)	620g	
Carton box (10ea)	6.6Kg	
<b>Environmental Conditions</b>		
Operating temperature	10 ~ 40°C (50 ~ 104°F)	
Operating humidity	5 ~ 85% non-condensing	
Storage temperature	−20 ~ 60°C (−4 ~ 140°F)	
Storage humidity	0 ~ 95% non-condensing	
Warranty		
Main unit 2 years		
Optional sensor & accessory	1 year	
Certificates		
KFDA, CE		

ndavd Canfigurati



# BT-710 Hand held pulse oximeter Operation Manual



Keep this manual for future reference

P/N: 710-ENG-OPM-EUR-R02

### **Proprietary Material**

Information and descriptions contained in this manual are the property of Bistos Co., Ltd. and may not be copied, reproduced, disseminated, or distributed without express written permission from Bistos Co., Ltd.

Information furnished by Bistos Co., Ltd is believed to be accurate and reliable. However, no responsibility is assumed by Bistos for its use, or any infringements of patents or other rights of third parties that may result from its use. No license is granted by implication or otherwise under any patent or patent rights of Bistos Co., Ltd.

The information contained herein is subjects to change without notice.

### Prepared by:

Bistos Co., Ltd. 7<sup>th</sup> FL., A Bldg., Woolim Lions Valley 5-cha, 302, Galmachi-ro, Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea

Telephone: +82 31 750 0340 Fax: +82 31 750 0344

Revision R02 March, 2020

Printed in Korea

Copyright © Bistos Corporation 2018. All rights reserved.

### **Contents**

0. Safety information	4
0.1 General precautions, warnings and cautions	5
0.2 Shock hazards	8
0.3 Battery warnings	8
0.4 General precautions on environment	9
1. System basics	9
1.1 Intended use	9
1.2 Operating principle	11
1.3 System configurations	11
1.4 Product outlook	
1.5 Description of pulse oximeter	11
1.6 Understanding the display	
1.7 Essential performance	13
2. Preparing for operation	13
2.1 Installation	13
2.2 Connecting to power	15
3. Basic operations	15
3.1 Turn on	15
3.2 Turn off	16
4. Setup the pulse oximeter	16
5. SpO <sub>2</sub>	17
5.1 Overview	17
5.2 Safety information	18
5.3 Monitoring steps	19
5.4 Setting SpO <sub>2</sub>	
5.5 Measuring influencing factors	
5.6 Technical Description	21
6. Review	21
7. Alarm	22
7.1 Alarm types	
7.2 Alarm condition priorities	23
7.3 Alarm mode	24

BT-710 Operation manual	3
7.4 Alarm state	
7.5 Alarm information	25
8. Battery	26
8.1 Overview	26
8.2 Battery usage guide	27
8.3 Checking battery performance	28
8.4 Battery recycling	28
9. Caring and cleaning	28
9.1 Overview	28
9.2 Cleaning	28
9.3 Disinfection	29
10. Maintenance	30
10.1 Checking	30
10.2 Trouble shootings	30
10.3 Maintenance plan	31
11. Accessories	31
12. Specifications	32
12.1 Safety specification	32
12.2 Hardware specifications	32
12.3 Functional specification	32
13. Manufacturer's declaration on EMC	34
13.1 Electromagnetic emissions	
13.2 Recommended separation distances between porta	ble and mobile RF communications
equipment and BT-710	35
13.3 Electromagnetic immunity	35

### 0 Safety information

Before using BT-710 Pulse oximeter, read this entire manual and be fully understood the following safety information to prevent injury of patient and user.

### Symbols Used

The following symbols identify all instructions that are important to safety. Failure to follow these instructions can lead to injury or damage to the pulse oximeter. When used in conjunction with the following words, the symbols indicate:



Can lead to serious injury or death.



Can lead to minor injury or product/property damage

The following symbols are placed on product, label, packaging and this manual in order to stand for the information about:

<u> </u>	Used to identify safety information.
	Be well-known this information thoroughly before using BT-710.
$\triangle$	Used to identify safety information.
	Be well-known this information thoroughly before using BT-710
	Indicates the protection level against the ingress of liquid.
181/2	IPX2 is protection from some water drops when the device is tilted
IPX2	up to and including 15°.
	It correspond the device and the accessory for SpO2.
Refer to operation manual. Read manual before placing the	
~~	Indicates the production date.
***	Indicates the manufacturer.
SN	Indicates the serial number of the device.
EC REP	Indicates the authorized representative in the European Community
of manufacturer.	
<u> </u>	Indicates a BF applied part.
	Indicates CLASS II equipment.(Adapter)
24	Indicates the date after which the medical device is not to be used.
<del></del>	Indicates to keep the device dry.
Indicates the medical device that can be broken or damage	
	handled carefully.
Ш	Indicates to keep upright
V	Indicates the temperature limitation for operation, transport and
-4	storage.
<b></b>	Indicates the humidity limitation for operation, transport and
الفيز	storage.

₽	Indicates the packing material is recyclable.
<u> </u>	Indicates to not dispose the device together with unsorted munici- pal waste (for EU only). The solid bar symbol indicates that mains adapter is put on the market after 13 August 2005.

### 0.1 General precautions, warnings and cautions

- Examine the pulse oximeter and any accessories periodically to ensure that the cables
  including adapter cords and instruments do not have visible evidence of damage that
  may affect patient safety or performance. The recommended inspection interval is
  once per week or less. Do not use the pulse oximeter if there is any visible sign of
  damage.
- Only the DC power adapter supplied with the BT-710 is approved for use with the device
- Do not attempt to service the BT-710 pulse oximeter. Only qualified service personnel by Bistos Co. Ltd. should attempt any needed internal servicing. There is no user serviceable part
- Perform periodic safety testing to insure proper patient safety. This should include leakage current measurement and insulation testing. The recommended testing interval is once per year.
- If the hospital or healthcare institutions using this device fail to implement a satisfactory ry maintenance schedule, it will result in device failure and may endanger the patient's safety.
- Use the pulse oximeter under the conditions specified in this operation manual. Beyond the conditions, the pulse oximeter may not function properly and the measurement results may not accurate and may result in device failure or endangering the patient's safety.
- · During the operation, do not disconnect any cable.
- Do not operate the BT-710 pulse oximeter if it fails to pass the power on self-test procedure.
- The BT-710 pulse oximeter is intended to be used by clinical professionals or trained doctors, nurses or laboratory assistant.
- Do not service and maintain or clean the device including accessories while in use with a patient.
- · Using the device to one patient at a time.

### ▲ WARNING

- Thoroughly read and understand the manual prior to use of the BT-710. Failure to do so could result in personal injury or equipment damage.
- The device is intended for measure the clinical blood oxygen saturation via

pulse oximeter, and only trained and qualified doctors and nurses should use the device.

- The alarm volume, upper and lower alarm limits should be set according to the
  actual situation of the using environment. Do not just rely on audio alarm system while monitoring the patient, because too low alarm volume or muted
  alarm may result in notice failure of alarm situation and endanger the patient's
  safety. Please pay close attention to the actual clinical status of the patient.
- Use only the power adapter supplied with pulse oximeter.
- Do not open the enclosure to avoid an electric shock. Any repair and upgrade
  of pulse oximeter should be done by service personnel trained and authorized
  by Bistos. Co., Ltd.
- When handling packaging materials and the device, abide by local laws and regulations or hospital waste disposal regulations. Keep the away from children
- Do not use in the presence of flammable anesthetics to prevent explosion or fire.
- Install the power lines and cables of accessories carefully to avoid patient entanglement or suffocation, cables tangled or electrical interference.
- When the pulse oximeter is used together with electrosurgical devices, the user (a doctor or a nurse) should ensure the safety of the patient and instrument.
- The physiological wave, physiological parameters and alarm information displayed on the pulse oximeter are only for the doctor's reference and should not be directly used as the basis for clinical treatment.
- · This is not a therapeutic device.
- Use of accessories other than approved for use with this product may result in increased emissions or decreased immunity.
- Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this manual. In addition, portable and mobile RF communications equipment can affect medical electrical equipment.
- The equipment shall not be used adjacent to other devices unless verification of normal operation in the configuration in which it is to be used can be achieved.
- Keep matches, and all other sources of ignition, out of the room in which the pulse oximeter is located. Textiles, oils, and other combustibles are easily ignited and burn with great intensity in oxygen rich environment. Personal injury or equipment damage could occur.
- A fire and explosion hazard exists when performing cleaning or maintenance procedures in an oxygen rich environment.
- The pulse oximeter has been validated with the accessories listed in this manual and found to comply with all relevant safety and performance require-

ments applicable to the device. It is therefore the responsibility of the person or organization who makes an unauthorized modification, or incorporates an unapproved attachment to the device.

- An operator may only perform maintenance procedures specifically described in this manual.
- Do not remove the covers of a BT-710 yourself to avoid damage to the equipment and unexpected electrical shock. Only qualified Bistos service engineer must repair or replace components.

### **∆**CAUTION

- Please install or carry the pulse oximeter properly to prevent damage due to falling, collision, strong vibration or other mechanical force.
- · Avoid instrument splashed by water.
- Avoid high temperatures, the instrument should be used within a temperature range of 5  $^{\circ}$  C  $^{\circ}$  40  $^{\circ}$ C  $_{\circ}$
- Avoid using pulse oximeter in the environment such as pressure is too high, poor ventilation, dusty, or contain salt, sulfur gas and chemical.
- Before using the pulse oximeter, check the pulse oximeter and accessories if there is damage that may affect patient safety. If there is obvious damage or aging, replace the parts before use. The replacement should be made with same parts of original parts.
- Before powering on the device, make sure that the power used by the device complies with the supply voltage and frequency requirements on the equipment label or in this manual.
- Equipment should be tested at least once a year, the test should be done and recorded by trained, have safety testing knowledge and experienced personnel. If there are any problems in the tests, they must be repaired.
- When the instrument and accessories are about to exceed the useful life (expected service life: 5 years), it must be treated in accordance with relevant local laws and regulations or the hospital's rules and regulations.
- Do not connect to other equipment or network which not specified in the instruction for use, in risk of external high voltage.
- Do not connect any equipment or accessories that are not approved by the
  manufacturer or according to IEC 60601-1 to the pulse oximeter. The operation or use of non-approved equipment or accessories with the pulse oximeter
  is not tested or supported, and pulse oximeter operation and safety are not
  guaranteed in such a case.
- · Any non-medical equipment (such as the external printer) is not allowed to be

used within the patient vicinity (1.5m/6ft.).

- Parts and accessories used must meet the requirements of the applicable safety standards, and/or the system configuration must meet the requirements of the medical electrical systems standard.
- Ensure that the conductive parts of electrodes and associated connectors, including neutral electrodes, do not come in contact with earth or any other conducting objects.

### 0.2 Shock hazards

### MARNING

- Unplug the pulse oximeter from its power source prior to cleaning or maintenance to prevent personal injury or equipment damage.
- Some chemical cleaning agents may be conductive and leave a residue that
  may permit a build-up of conductive dust or dirt. Do not allow cleaning agents
  to contact electrical components and do not spray cleaning solutions onto any
  of these surfaces. Personal injury or equipment damage could occur.
- Do not expose the pulse oximeter to excessive moisture that would allow for liquid pooling. Personal injury or equipment damage could occur.
- Do not touch the patient and signal input/output parts simultaneously
- Due to the risk of electrical shock hazard, only qualified personnel with appropriate service documentation should service the pulse eximeter.

### 0.3 Battery warnings

### MARNING A

- Improper operation may cause the internal lithium ion battery to be hot, ignited or exploded, and it may lead to the decrease of the battery capacity. It is necessary to read the operation manual carefully and pay more attention to warning message.
- Do not open the battery compartment. Only the qualified service personnel authorized by the manufacturer can open the battery compartment and replace the battery, and batteries of same model and specification should be replaced.
- · Be careful when connecting the battery with polarity.
- Do not use the battery near fire or environmental temperature exceeds 60 °C.
   Do not heat or splash the battery or throw it into fire or water.
- Do not destroy the battery. Do not pierce battery with a sharp object such as a needle. Do not hit with a hammer, step on or throw or drop the battery. Do

not disassemble or modify the battery. The battery can heat, smoke, deformation or burning.

- When leakage or foul smell is found, stop using the battery immediately. If your skin or cloth comes into contact with leaked liquid, cleanse it with clean water at once. If the leaked liquid splashes into your eyes, do not wipe them. Irrigate them with clean water first and go to see a doctor immediately.
- Properly dispose of or recycle the depleted battery according to local regulations

### 0.4 General precautions on environment

Do not keep or operate the pulse oximeter under the environment listed below.

	Avoid placing in an area exposed to moisture. Do not touch the pulse oximeter with wet hand.		Avoid exposure to direct sunlight
	Avoid placing in an area where high variation of temperature exists. Operating temperature ranges from 5 °C ~ 40 °C. Operating humidity ranges from 30 % ~ 85 %.	<b>2</b> 7	Avoid in the vicinity of electric heater.
S. M.	Avoid placing in an area where there is an excessive humidity rise or ventilation problem.		Avoid placing in an area where there is an excessive shock or vibration.
	Avoid placing in an area where chemicals are stored or where there is in danger of gas leakage.		Avoid dust and especially metal material enters into the pulse oximeter.
(0,0,0)	Do not disjoint or disassemble the pulse oximeter. Bistos Co., Ltd. does not have liability of it.		Power off when the pulse oximeter is not fully ready to operate. Otherwise, the pulse oximeter could be damaged.

### 1 System basics

### 1.1 Intended use

The BT-710 pulse oximeters acquire the blood oxygen saturation (SpO<sub>2</sub>) and pulse rate (PR). The signals are converted into digital data and processed, examines the data for alarm

conditions and display the data. The pulse oximeter also provides operating control for the user. The pulse oximeter intend to use in hospital clinical area such as general ward, to provide additional information to the medical and nursing staff about the blood oxygen saturation of the patient. The BT-710 pulse oximeters are intended to be used only under regular supervision of clinical personnel. It is suitable for adult and pediatric, neonate. The intended locations of use are hospitals and clinics.

- 1) Intended patient population
  - Adult (>18 years adults) and Pediatrics (30 days < and <18 years) and Neonate (0 days < and < 30days)</li>
- 2) Intended user profile
  - Doctor, physicians or nursing staff who is qualified personnel
  - Basic experiences or knowledge on medical field, especially on patient monitoring
  - Trained or requested to read IEU before use
- 3) Environment of use
  - Hospital and clinic
    - Requirements: Stable power source

### 4) Scope of application

This pulse oximeter is suitable for bedside monitoring of patient. This pulse oximeter enables the monitoring of blood oxygen saturation (SpO<sub>2</sub>) and pulse rate (PR). It is equipped with a replaceable built-in battery to provide convenience for the patient movement in hospital.

### 5) Indications and contraindications

### Blood oxygen saturation (SpO<sub>2</sub>)

### Indication:

- Monitoring effectives of oxygen therapy
- A reading is needed to facilitate the completion of an early warning score to inform clinical assessment
- Sedation or anesthesia
- Transport of patients who are unwell and require oxygenation assessment
- Haemodynamic instability (e.g. cardiac failure or Myocardial Infarction)
- Respiratory illness e.g. asthma, chronic obstructive pulmonary disease
- Monitoring during administration of respiratory depressant drugs, e.g. opiate epidural or patient-controlled analgesia.
- Assessing oxygen saturation during physical activity e.g. in pulmonary rehabilitation

### Contraindications

 Pulse oximetry does not give an indication of haemoglobin so if the patient is profoundly anaemic then their oxygen saturation may by normal but they may still be hypoxic

Source: NHS. "Clinical Procedure\_ Procedure for Pulse Oximetry/SPO2". Wirral Community NHS Trust. Sep. 2013

### 1.2 Operating principle

Refer to the chapters for every physiological parameter from chapter 5.

### 1.3 System configurations

Basic configuration of BT-710

- Main body with 4.3" touch screen and built-in lithium-ion battery
- Adult SpO<sub>2</sub> sensor probe
- AC/DC adapter

### 1.4 Product outlook







Figure 1-1: Front view

Figure 1-2: Rear view

Figure 1-4: Bottom view

### 1.5 Description of pulse oximeter



Figure 1-5: Front view

Name		Description	
1	DC power indicator	Turned on when the pulse oximeter is being powered by the adapter.	
2	Display area	Display the waveform and measured value	
3	<b>し</b>	- Power On: Press down the key more than 2 seconds.	

		the system will display the alarm message "The system	
4	[Alarm reset]	To reset the alarm condition.	
5	Setting]	Enter to the setting mode. Press again to close the setting mode.	



Figure 1-6: Top view

	Name	Description	
6	SpO <sub>2</sub>	SpO <sub>2</sub> sensor probe interface	



Figure 1-7: Bottom view

	Name	Description	
7	SD card interface	For software upgrade	
8	Power adapter	5V, 2A adapter	
9	Lanyard eyelet	For convenient hand held	

### 1.6 Understanding the display



Figure 1-8: Standard display

	Description		
1	Current alarm message. When an alarm occurs, this area will displayed yellow red depending on the alarm type.		
2	SpO <sub>2</sub> value. Display the measured SpO <sub>2</sub> value.		
3	SpO <sub>2</sub> upper alarm limit. Display the user set upper alarm limit		
4	SpO2 lower alarm limit. Display the user set lower alarm limit		
5	SpO <sub>2</sub> waveform. Display the measured SpO <sub>2</sub> waveform. The waveform is not normalized.		
6	Pulse rate value. Display the measured pulse rate per minute.		
7	Sleep mode. Touch this area makes the pulse oximeter to enter the sleep mode. To exit "sleep mode", press [Power] or [Alarm reset] or [Setting]		
8	Perfusion Index. Display the measured perfusion index.		
9	Pulse rhythm strength.		
10	Battery Status		
11	Patient type.		
12	The unit of SpO <sub>2</sub> .		

### 1.7 Essential performance

This device Pulse oximeter provides patient vital signs such as pulse rate, blood oxygen saturation and perfusion index by placing the sensors to the appropriate site of patient. The device is composed with display, control circuit and panel, and input part for SpO<sub>2</sub> sensors. It detects SpO<sub>2</sub> and PR using specific sensors. The detected analog signal amplified and converted to digital. This converted data feed to the CPU and converted to the display format as number and waveform. This device is incorporated with alarm system. The alarm generated when the detected signal range is beyond the user set alarm limits.

### 2 Preparing for operations

### 2.1 Installation

To ensure normal working of the pulse oximeter, read this chapter before use, and install as required.

### **⚠** WARNING

All analog and digital devices connected to the pulse oximeter must be certified by IEC standards (e.g. IEC 60950 Data processing equipment standard and IEC 60601-1 Medical equipment standard). Furthermore, all configurations shall comply with valid version of IEC 60601-1 standard. The personnel connecting additional devices to the input / output signal ports are responsible for the compliance with IEC 60601-1 standard. If there is any question, please contact Bistos

- The copyright of pulse oximeter software belongs to Bistos. Without permission, any organization or individual shall not interpolate, copy or exchange by any means or form.
- When the pulse oximeter is combined with other devices, it must comply with IEC 60601-1:2005 + A1:2012, and should not be connected with multi-socket outlet or extension cord.

Prior to installation, the operator must ensure that the following space, power, environmental requirements are met.

### 2.1.1 Unpack and check

BT-710 pulse oximeter was inspected rigorously at the factory before delivery, in order to avoid being hit when transported, carried out careful packaging. Before unpacking, carefully inspect the package. If any damage, please immediately contact the Bistos. Unpack in the correct way, carefully remove the pulse oximeter and accessories from the box and check with the packing list. Check if there is any mechanical damage, and the all listed components are completely packed. If you have questions, please contact the marketing department of Bistos or agency.

Please keep the packing box and materials for use in future transporting or storage.

- 2.1.2 Power requirements
- DC power supply adapter (Model: UE10WCP1-050200SPA)

Input: A.C. 100 V ~ 240 V. 50/60 Hz

Output: D.C. 5 V. 2.0 A

- Built-in rechargeable lithium-ion battery: D.C. 3.7 V, 3000mAh (Model: JHY605085)
- 2.1.3 Environmental requirements

The storage, transport and use of the pulse oximeter must meet the following environmental requirements

Operating	Ambient temperature	5℃~40℃			
environment	Relative humidity	30 % ~ 85 % (Non-condensing)			
environment	Atmospheric pressure	700 ~ 1060 mbar (hPa)			
Transportation	Prevent severe shock, vibration, rain and snow splashing during				
Transportation	transport.				
Storage	The packaged pulse oximeter should be stored in well-ventilated room with ambient temperature -20 °C ~ 60 °C, relative humidity 0 ~ 95 % (Non-condensing), atmospheric pressure 700 ~ 1060 mbar(hPa), and without corrosive gases.				

The operating environment of the pulse oximeter should avoid noise, vibration, dust, corrosive or flammable and explosive materials. In order to allow air flowing smoothly and achieve good heat dissipation, at least 2 inches (5cm) clearance should be kept around the device. When the device is moved from one environment to another, the device may have condensation due to the differences in temperature or humidity. In this case, wait until the condensation disappears before using the device.



### MARNING

Ensure that the pulse eximeter is used under specified environment. Fail to do this, the technical specifications declared in this manual may not be met and it may result in damage to equipment and other unforeseen consequences.

### 2.2 Connecting to power



· Do not try to open the pulse oximeter when the power is connecting.

Connect to power adapter in the following steps:

- Make sure that the AC power supply meets the following specifications: A.C. 100V-240V. 50/60Hz
- Use the power adapter provided with the pulse oximeter. Plug the power adapter into the power connector of the pulse eximeter, and plug the other end of the power adapter into the mains (low voltage power supply network facilities) power outlet.

### NOTE

 The operation of pulse eximeter after the supply mains has been interrupted and is restored after a period of time that is longer than 30s.

### 3 Basic operations

### 3.1 Turn on

- 3.1.1 Check the pulse oximeter
- Before turn on the pulse oximeter, check whether there is mechanical damage to the pulse oximeter, and whether the external cables and accessories are connected correctly.
- Plug the power adapter into the AC power outlet. If using battery power, make sure the battery is fully charged.
- Check all the functions required to make sure that the pulse oximeter operates properly.



If the pulse oximeter is damaged, or fails to work normally, do not use it.

Please contact the maintenance personnel or Bistos immediately.

### 3.1.2 Start the pulse oximeter

If finish to check the pulse oximeter, it is ready to start the pulse oximeter.

Press the [Power] key and the system enters the main interface within seconds.

- Once the power is supplied, the system performs a power-on self test to check the functions before start-up. If any fatal error occurs during boot up, the system will alarm. If this case persists, please stop to using the pulse oximeter and contact the maintenance personnel or Bistos.
- Check all available pulse oximeter functions to ensure that the pulse oximeter operate properly.

- If the pulse oximeter equipped with a battery, charge the battery after each use to ensure sufficient power.
- After unpacking and use the pulse oximeter at first time, the pulse oximeter should be powered with adapter.

### 3.1.3 Connect the sensors

Connect the SpO<sub>2</sub> sensor probe to the pulse oximeter and the monitoring site of patient.

### 3.1.4 Start monitoring

Start monitoring in the following steps:

- Check if the patient cable is connected properly.
- Check if the settings of the pulse oximeter are corrects, such as patient type.
- For the details of parameter measurement or monitoring, see the appropriate section.
- The operator can operate according to their own habits, standing in front. left or right of the monitor, easy to observe and operate the monitor.

### 3.2 Turn off

Turn off the pulse oximeter in the following steps

- Disconnect the sensor probe connected to the patient.
- Press and hold the [III [Power] key for 2 seconds to pop up the 3 seconds countdown window, and the pulse oximeter turns off in 3 seconds.

### **ACAUTION**

• If the pulse oximeter is not turned off properly, you can simply disconnect the power to shutdown forcibly. But the forced shutdown may cause data loss, and it is not recommended.

### 4 Setup the pulse oximeter


In the main screen, Press the key enter the setup menu.



- - Select "Patient Type", you can choose "Adult" or "Pediatric" or "Neonate". According to patient type, the settable PR alarm limit and default are changed.
- ✓ Select "Alarm Delay", you can choose "off", "1s", "2s", "3s", "4s", "5s", "6s", "7s" or "8s". It sets the time from the occurrence of a triggering event either in the patient, for physiological alarm conditions to when the alarm system determines that an alarm condition exists. Technical alarm conditions are not delayed.
- ✓ Select "Backlight", and enter value ( Range: 0-5 ). It can set brightness of the background according to needs of use. "5" is the bright-
- ✓ Select "Alarm Volume", which is alarm tone volume and settable among (Range: 0-9), "0" is alarm sound off.
- ✓ Select "Alarm Reminder Signal", you can choose "On", or "Off". It is periodic signal that reminds the operator that the alarm system is in an alarm signal inactivation for both physiological alarm conditions

- and technical alarm conditions
- ✓ Select "Alarm Reminder Interval", when the "Alarm Reminder Signal" sets "On", you can choose alarm interval to "1min", "2min" or "3min".
- ✓ Select ▶ key to display the following interface.
- ✓ Select "Year", "Month", "Day", "Hour", "Minute" and "Second", and set the current date.
- ✓ Select "Date Format", and set among the below date. formats in accordance with custom
  - "YYYY-MM-DD": Year- Month-Day "MM-DD-YYYY": Month -Day-Year.
  - "DD-MM-YYYY": Day-Month-Year.
- ✓ Select key to display the following interface.

✓ Select "Displays", set "SpO2".



- "English": The interface language of the pulse oximeter is English. "Türkce": The interface language of the pulse oximeter is Turkish. "Español": The interface language of the pulse oximeter is Spanish.
- "Français": The interface language of the pulse eximeter is French. "Polski": The interface language of the pulse oximeter is Polish. "Italiano": The interface language of the monitor is Italian. "Deutsch": The interface language of the monitor is German.
- ✓ Select "Default", clears all the setting except the date and time and go back to the initial factory setup.
- ✓ Select key to display the following interface.



### 5.1 Overview

Circulary

UH 95 SMN 360

MP:55

Blood oxygen saturation (SpO<sub>2</sub>) is the percentage of oxyhemoglobin (HbO2) capacity bound by oxygen in the blood in the total hemoglobin (Hb) capacity that can be combined, that is, the concentration of oxygen in the blood.

The principle for monitoring the pulse SpO2 is to fix the probe fingerstall on the patient's finger, use the finger as a transparent container for hemoglobin, use 660nm wavelength red light and 905nm near-infrared light as the incident light, maximum output power is 300 mW, measure the light transmission intensity through the tissue bed, and calculate the concentration of hemoglobin and SpO2.

The passing lights depend on a variety of factors, most of which are constant. However, one of these factors, the arterial blood flow, changes with time, as it is pulsating. By measuring the light absorbed during pulsating, it is possible to obtain the arterial blood SpO2. Detection pulsation can give a "plethysmography" wave and pulse rate signal.



The main screen displays "SpO<sub>2</sub>" value and "plethysmography" wave. While the monitor gets and displays the physiological signal from a patient, the pulse tone which is the auditory information signal. The volume of this pitch tone is user adjustable.

This pulse oximeter applies to measure  $SpO_2$  of adults (>18 years), pediatric (<18 years,>30 days), neonate (<30 days). Contact  $SpO_2$  probe to Patient's finger to get " $SpO_2$ " value and "olethysmography" wave.

SpO<sub>2</sub> function of this pulse eximeter has been calibrated in factory.

### NOTE

· Information about wavelength range can be especially useful to clinicians.

### 5.2 Safety information



- Please use SpO<sub>2</sub> sensor supplied from Bistos, operate in accordance with the Manual, and observe all warnings and precautions.
- Before monitoring, check whether the sensor probe is normal. When SpO<sub>2</sub> sensor probe is unplugged from the socket, the screen will display "SpO2 No Sensor connected" error message, and trigger an audible and visual alarm simultaneously.
- If the sensor or sensor packaging has signs of damage, do not use this SpO<sub>2</sub> sensor: return it to the manufacturer.
- If there is carboxyhemoglobin, methemoglobin or dye diluted chemical, the SpO<sub>2</sub> value will have deviation.
- When the patient has a tendency to hypoxia, use the oximeter to analyze blood samples in order to fully grasp the patient's condition.
- Do not put the sensor on limbs with arterial duct or intravenous tube.
- Do not intertwine electrosurgical equipment cable with the patient cable.
- Avoid using the pulse oximeter and sensors while using the NMR equipment, in order to avoid severe burns to the patient as a result of induced currents.
- During long time continuous monitoring of a patient, check the position of SpO<sub>2</sub> sensor once every 2 hours, and move properly when the skin changes or every four hours. Some patients may require more frequent inspection, such as patients with perfusion disorders or sensitive skin, because persistent and prolonged monitoring may increase unpredictable skin changes, such as allergies, redness, blistering or pressure necrosis.
- Before using, verify compatibility between the monitor and probe, otherwise it
  may cause injury to the patient.
- Functional testers cannot be used to evaluate the accuracy of pulse oximetry.

### NOTE

- Do not put the oxygen probe and blood pressure cuff on the same limb, because blood flow occlusion during blood pressure measurement will affect the SpO, readings.
- The pulse oximeter cannot be used to verify the accuracy of SpO<sub>2</sub> probe and SpO<sub>2</sub> equipment.

### 5.3 Monitoring steps

- (1) Act the appropriate SpO<sub>2</sub> sensor according to the patient.
- (2) Turn on the pulse eximeter, and connect the SpO<sub>2</sub> sensor probe to the pulse eximeter.
- (3) Clean the measurement site, such as finger with nail polish.
- (4) Put the SpO<sub>2</sub> sensor probe on the patient's finger.
- (5) Select the appropriate alarm settings.
- (6) Start monitoring.

### NOTE

- When you turn on the pulse oximeter, plug in SpO<sub>2</sub> probe and connect patient's finger (or toe), monitor displays SpO<sub>2</sub> wave, "SpO<sub>2</sub> Pulse Search" displayed in the technical alarm area until the monitor measured SpO<sub>2</sub> value and pulse rate. "SpO<sub>2</sub> Search Timeout" displayed in the technical alarm area until the pulse oximeter measured pulse rate. Check the sensor mounting position, whether the sensor is damaged or sensor type. Reconnect the sensor or use new sensor.
- The alarm message for signal inadequacy such as "SpO2 No Sensor connected", "SpO2 Sensor Off", "SpO2 Search Timeout", "SpO2 Signal Unstable", "SpO2 Sensor Error" and low priority alarm and indication for signal instability like "---" are meaning that the value of SpO2 or PR might be inaccurate.

### 5.4 Setting SpO<sub>2</sub>

Select the main screen to enter the SpO2 set interface





Select "Wave Speed" and set wave speed to "12.5mm/s" or "25mm/s". The faster speed, the smoother wave.

- Select "Wave Mode", and set the wave drawing mode to "Scan" or "Fill". This displays a wave as line graph or filled line graph.
- Select "Average Time", and set the average time to "2-4s", "4-8s", "8-16s". According to the selected option, the module displays the average of data collected within that time. The shorter the averaging time is, the quicker the monitor responds to changes in the patient's oxygen saturation level.
- Select "Pulse Volume", and set the volume for auditory information signal for pulse rate ( Range: 0-9 ).
- Select "Sensitivity", and you can choose "1", "2", "3". It is responding speed for changes of parameter setup. It is faster when setting value is higher.
- Select key to display the following interface.



- Select "Patient Type", you can choose "Adult" or "Pediatric" or "Neonate". According to patient type, the settable PR alarm limit and default are changed.
- Select "SpO2 Low Limit", and set the SpO2 alarm low limit value (Range: 0-99), Adult/Pediatric/Neonate default: 90. The input cannot enter higher than the SpO2 alarm upper limit.
- Select "SpO2 High Limit", and set the SpO2 alarm upper limit value (Range: 1-100), Adult/Pediatric/Neonate default: 100/100/95. The input cannot enter less than the SpO2 alarm low limit.
- Select "PR Low Limit", and set the RP alarm low limit value (Adult Range: 15-299, Pediatric/Neonate Range: 15-349), Adult/Pediatric/Neonate default: 50/75/100. The input cannot enter higher than the PR alarm upper limit.
- Select "PR High Limit", and set the PR alarm upper limit value (Adult Range: 16-300, Pediatric/Neonate Range: 16-350), Adult/Pediatric/Neonate default: 120/160/200. The input cannot enter less than the PR alarm low limit.
- Select "PI Low Limit", and set the PI alarm low limit value (Range: 0.00-19.90). Default: 0.00.
- Select "PI High Limit", and set the PI alarm upper limit value (Range: 0.10-20.00). Default: 20.00.
- Select "Alarm Level", you can choose "Mid" or "High" of priority level for the physiological alarm.
- Alarm, Select "Default", the alarm parameter is set to the default value.

### NOTE

 Setting values including alarm settings are stored and does not changed even when the power supply interrupted.

### 5.5 Measuring influencing factors

During operation, the following factors can affect the accuracy of SpO<sub>2</sub> measurement:

✓ High-frequency radio wave interference, such as interference generated by the host

21

system or interference from electrosurgery instrument connected to the system.

- ✓ Intravenous dye.
- ✓ Too frequent movement of the patient.
- ✓ External light radiation.
- ✓ Sensor is improperly installed or improperly in contact with the patient.
- ✓ Sensor temperature.
- ✓ The sensor is placed on limbs with blood pressure cuff, arterial duct or lumen tube.
- ✓ Concentration of non-functional hemoglobin such as carboxyhemoglobin (COHb) and methemoglobin (MetHb).
- ✓ Shock, anemia, hypothermia, and the application of vasoconstrictor drugs may reduce the
  arterial blood flow to a level that cannot be measured.
- ✓ The measurement also depends on the absorption of specific wavelengths of light by oxyhemoglobin and reduced hemoglobin. If there is any other substance that absorbs the same wavelength, the measurement may have false or low SpO₂ values, such as: carbon hemoglobin, methemoglobin, methylene blue, and indigo carmine.
- ✓ SpO₂ probe described in Annex is recommended.
- ✓ Operating environment limit: Operating temperature range: 5 ~ 40 °C, Humidity range: 30%~85% (non-condensing) Atmospheric pressure: 700hPa ~ 1060hPa.

#### 5.6 Technical description

- The SpO2 sensor probe (specified in Chapter 11) material which contacts patients or other staff has passed the biocompatibility test and meet the requirements of ISO 10993-1
- Fluke's index 2XL Oxygen Analyzer can be used to check the function of the pulse oximeter and can be used to assess the accuracy of the pulse rate but cannot be used to assess the accuracy of blood oxygen.
- Functional testers cannot be used to evaluate the accuracy of pulse oximetry probe and pulse oximetry.
- Measure the maximum temperature between the oxygen probe and the tissue contact surface: Measured as described in Annex BB of ISO 80601-2-61,the temperature is less than 41°C
- The PR accuracy shall be stated as the root-mean-square (rms) difference between paired pulse rate data recorded pulse oximeter equipment and with a reference method (E.G, an electronic pulse simulator, ECG heart rate, palpated pulse, thoracic auscultator or a second pulse equipment which has been qualified by comparison to one of these references).

#### 6 Review

Use the trend screen to recall all the historical patient data in a list, including monitoring time (in 1-minute intervals), and SpO<sub>2</sub> and PR values. The most recent measurements display at the top of the list.



Figure 6-1: SpO<sub>2</sub> Trend (1)



Figure 6-2: SpO<sub>2</sub> Trend (2)

To view historical data:

- ✓ Press the [Setting] key twice to display the Trend screen.
- ✓ Select or key to page up and down to view patient data.
- ✓ Press key to exit the Trend screen and display the Main screen or press the 

  [Setting] key to see the trend graph.

- ✓ Select 'SpO2', 'PR' or 'PI' to see.
- ✓ Press key to exit the Trend screen and display the Main screen.

#### 7 Alarm

Alarm means that the pulse oximeter prompts the medical staff through sound and light when the abnormal changes in vital signs are monitored or the pulse oximeter has a failure or is unable to monitor the patient successfully. For better viewing those alarms, it is recommended to look straight at the display when using the monitor.

#### ⚠ WARNING

In any single area (e.g. intensive care unit), it can be a potential danger if there
are the same or similar devices using different alarm preset.

After setting, the alarm and other parameters of the monitor won't be lost when the system is power off, unless modified manually. Connect the power again and turn on the monitor, it will resume normal working, and the alarm and other parameters remain unchanged.

BT-710 Operation manual

23

#### 7.1 Alarm types

According to the nature of the alarm, the alarms of the monitor can be divided into physiological alarms, technical alarms, alarm reminder signal and prompt messages.

- ✓ Physiological alarms
  - A physiological alarm is usually triggered when a physiological parameter of the patient exceeds the alarm limit or the patient has physiological abnormalities. The information of physiological alarm is displayed in the physiological alarm area on too of the screen.
- √ Technical alarms

Technical alarm is also known as a system error message, which is caused by improper operation or system failure resulting in system malfunction or monitoring result distorted. The information of technical alarm is displayed in the technical alarm area on top of the screen.

- ✓ Alarm Reminder Signal
  - Alarm Reminder Signal is periodic signal that reminds the operator that the alarm system is in an alarm signal inactivation for both physiological alarm conditions and technical alarm conditions.
- ✓ Prompt messages

Strictly speaking, the prompt messages are not alarms. The monitor also will display some information associated with system status in addition to the physiological alarms and technical alarms, and generally such information do not involve the patient's vital signs. The prompt messages generally appear in the technical alarm area and parameters area.

#### 7.2 Alarm condition priorities

According to the severity of the alarm conditions, the physiological alarms of the monitor can be divided into high priority, medium priority.

- ✓ High priority alarms
  - The patient is in critical condition that is life-threatening, and should be immediately rescued, or the monitor has a serious mechanical failure or malfunction, causing it unable to detect the patient's critical state and endangering the patient's life.
- ✓ Medium priority alarms

The patient's physiological signs are abnormal and appropriate measures or treatment should be taken immediately, or although it won't endanger the patient's life, the mechanical failure or deactivation of the monitor will affect the normal monitoring of key physiological parameters.

- ✓ Low priority alarms
  - The certain monitoring function is invalid due to mechanical failure or deactivation, but it won't endanger the patient's life.

The levels of some physiological alarms can be modified. But, the priority of all technical alarms and some physiological alarms have been set in the monitor at the factory and cannot be modified by the user.

2020.03

#### 7.3 Alarm mode

When an alarm occurs, the monitor uses the following audible or visual alarm to prompt the

- ✓ Visual alarm
- ✓ Audible alarm
- ✓ Alarm info
- ✓ Parameter flashing

Of which, the visual alarm, audible alarm, and alarm information distinguish the alarm levels in a different manner respectively.

#### 7 3 1 Visual alarm

When an alarm occurs, the alarm indicator will flash in different colors and frequencies to prompt the alarm priority.

- ✓ High priority alarm: Red
- ✓ Medium and Low priority alarm: Yellow

#### 7.3.2 Audible alarm

An audible alarm is that the monitor prompts the alarm priorities with alarm tone characteristics when an alarm occurs.

Audible alarm pattern:

- ✓ Medium priority alarm: Beep-beep-beep
- ✓ Low priority alarm: Beep-

#### 7.3.3 Alarm information

Alarm information displayed on the physiological or technical alarm area of the monitor indicates the corresponding alarm information when an alarm occurs. The system will distinguish the alarm priority with different background colors:

- ✓ High priority alarm: Red
- ✓ Medium priority alarm: Yellow
- ✓ Low priority alarm: Yellow

The following flags in front of physiological alarms are used to distinguish the alarm priorities.

- ✓ High priority alarm: \*\*\*
- ✓ Medium priority alarm: \*\*
- ✓ Low priority alarm: \*

#### 7.3.4 Parameter flashing

When the physiological parameter values in the parameter area will flash once per second, and the indicator for the upper limit and lower limit of the parameter area will also flash at the same frequency, it indicating that the parameter exceeds the upper limit or lower limit.

#### 7.4 Alarm state

#### 7 4 1 Alarm reset

Select button, and you can temporarily turn off the alarm sound of currently occurring physiological alarms of the pulse oximeter, but the alarm information is still retained. For technical alarms, clear the alarm state, display alarm prompt information, the alarm state icon area displays the icon. When a new physiological alarm or technical alarm occurs, the alarm reset is automatically canceled.

#### 7.4.2 Alarm sound off

When "Alarm Volume" sets "0", the alarm state area on the screen shows the con. If "Alarm Volume" sets bigger than "0", the system will cancel alarm sound off state.

#### ▲ WARNING

 When the alarm is off, and the alarm reminder signal is on, the system will have alarm reminder tone.

#### 7.5 Alarm information

This chapter lists some important physiological and technical alarm information, and some alarms are not necessarily listed.

Note that in this chapter: P column indicates the default alarm priority: H indicates high priority, M indicates middle priority, L indicates low priority, and "\*" indicates configurable priority by the user.

Corresponding countermeasures are listed for each alarm message. If you operate in accordance with the countermeasures but the problem persists, contact your service personnel.

✓ Physiological alarm information.

Source	SpO <sub>2</sub>				
		Default		Causes and countermeasures	
Parameter	Alarm message	Value (Adult/ Pediatric/ Neonate)	Р		
SpO₂ High Limit	SpO₂ Too High	100/100/95		SpO <sub>2</sub> value is higher than the settled upper alarm limit or lower than the settled lower alarm limit. Check the patient's physiological	
SpO <sub>2</sub> Low Limit	SpO <sub>2</sub> Too Low	90	М	condition, and check if the patient categoriand alarm limit settings are appropriate for the patient.	
PR High Limit	PR Too High	120/160/200	*	PR value is higher than the settled upper alarm limit or lower than the settled lower alarm limit. Check the patient's physiological	
PR Low Limit	PR Too Low	50/75/100		condition, and check if the patient category and alarm limit settings are appropriate for the patient.	

PI High Limit	PI Too High	0.00		PI value is higher than the settled upper alarm limit or lower than the settled lower alarm limit. Check the patient's physiological condition, and when the PI alerts the clinician to consider another monitoring site	
PI Low limit	PI Too Low	20.00			
-	No Pulse	-	Н	The physiological signal shows a potentially life-threatening drop in oxygen saturation. Check the patient's physiological condition.	

#### ✓ Technical alarm information

Source	Alarm	P	Causes and countermeasures
	SpO <sub>2</sub> No Sensor connected	L	Indicates the condition that SpO <sub>2</sub> sensor cable is unplugged from the socket.
SpO₂	SpO <sub>2</sub> Search Timeout	L	When you turn on the monitor, plug in SpO <sub>2</sub> probe and connect patient's finger (or toe), monitor displays SpO <sub>2</sub> wave, "SpO <sub>2</sub> valse Search' displayed in the technical alarm area until the monitor measured SpO <sub>2</sub> value and pulse rate. "SpO <sub>2</sub> Search Timeout' displayed in the technical alarm area until the monitor measured pulse rate. Check the sensor mounting position, whether the sensor is damaged or sensor type. Reconnect the sensor or use new sensor.
	SpO <sub>2</sub> Sensor Off	L	SpO <sub>2</sub> sensor falls off from the patient or monitor, malfunctions, or sensor other than specified in this Manual is used. Check the sensor mounting position, whether the sensor is damaged or sensor type. Reconnect the sensor or use new sensor.
	Battery Low	М	Connect to DC power supply, and charge the battery, and power with the battery as needed after fully charged.
Battery	System will shutdown	N/A	Press down the keys more than 2 seconds and the system will display the alarm message "The system will shut down 3 seconds"

#### 8 Battery

#### 8.1 Overview

The pulse oximeter has a built-in rechargeable battery to ensure that the pulse oximeter can also be used normally in case of patient transfer or power failure. When the oximeter is connected to a DC power source, it will charge the battery no matter whether the pulse oximeter is turned on or not. In the case of power failure, the system will automatically use the battery to power the pulse oximeter to avoid interrupting the pulse oximeter working. The battery icon on the screen indicates the battery status:



Battery is working properly and is fully charged.



Battery is working properly and the green part indicates the battery power.



Battery power is low, and requires charging immediately, or else the pulse oximeter will turn off automatically.



Battery is not installed.



Battery is properly installed and being charged.

The battery power can only maintain for some time. Low battery voltage will trigger a high level technical alarm "Battery Low": in this case, connect the pulse oximeter to DC power and charge the battery.

#### 8.2 Battery usage guide

Battery life depends on the frequency and time of use. If the battery maintenance and storage are proper, the lithium battery life is three years. If you do not use the battery properly, its life may be shortened. It is recommended to replace the lithium battery once every three years.

In order to ensure the maximum capacity of the battery, please note the following usage guide:

- ✓ Before using the battery, please read this manual and labels on the battery surface. carefully.
- ✓ Do not drop the battery.
- ✓ If it won't be used for a long time (over three months), please store the battery properly. Charge the battery to 50%, and wrap the battery with non-conductive material in order to avoid direct contact with metal, resulting in damage. Keep the battery in a cool dry place.
- ✓ Check the battery performance once every two years. Before servicing the pulse oximeter or you suspect that the battery is the fault, also check the battery performance.

#### ♠ WARNING

- Keep the battery out of the reach of children.
- Use only the designated battery.
- · If the battery is damaged or leaks, replace it immediately. Do not use a defective battery for the monitor.
- Do not open the battery compartment. Only the qualified service personnel authorized by the manufacturer can open the battery compartment and replace or check the battery, and when it needs replacement, the replacement should be same model with established battery.
- Do not disassemble the battery, throw it in fire, or short-circuit it. Battery fire, explosion and leakage may lead to personal injury. Do not touch the leaking battery with bare hand directly.

#### 8.3 Checking battery performance

Please refer to the following steps to check the battery performance:

- ✓ Disconnect the pulse eximeter from the patient and stop the monitoring or measurement.
- Connect DC power to the monitor, and charge the battery for more than 4 hours uninterruptedly.
- ✓ Disconnect the DC power and power the pulse oximeter with battery until the pulse oximeter is turned off.
- ✓ Battery duration reflects the battery performance.

If the battery operating time is significantly shorter than the time stated in specifications, please contact our service personnel for replacing the battery.



 Do not open the battery compartment. Only the qualified service personnel authorized by the manufacturer can open the battery compartment and replace or check the battery, and when it needs replacement, the replacement should be same model with established battery.

#### 8.4 Battery recycling

If the battery has visible damage or cannot store power, it should be replaced and recycled properly. Follow the appropriate regulations to dispose of used batteries.

## ⚠ WARNING

 Do not disassemble the battery, throw it in fire, or short-circuit it. Battery fire, explosion and leakage may lead to personal injury. Do not touch the leaking battery with naked hand directly.

#### 9 Caring and cleaning

#### 9.1 Overview

In the using process, please make sure that there is no dust on or near your device. To prevent damage, please use the diluted detergents and disinfectants specified in this manual, and use the lowest possible concentration. For the damage or accident caused by using other materials or methods, our company does not assume any responsibility.

#### 9.2 Cleaning

The device should be cleaned regularly. In the heavily polluted environment, increase the frequency of cleaning. Before cleaning, please consult the hospital about device cleaning requirements.

Below are available cleaning agents:

- ✓ Diluted ammonia
- ✓ Diluted sodium hypochlorite (washing bleach)
- ✓ Diluted formaldehyde

- ✓ Hydrogen peroxide (3%)
- √ Ethanol (70%)
- √ Isopropanol (70%)

Before cleaning:

- ✓ Turn off the pulse oximeter, disconnect the power cord.
- ✓ Use a soft cotton ball to adsorb appropriate amount of cleaning agent and wipe the display screen.
- ✓ Use a soft lint-free cloth to adsorb appropriate amount of cleaning agent and wipe the surface of the device.
- ✓ If necessary, use a clean, dry, lint-free cloth to remove any excess detergent.
- ✓ Dry the device naturally in a ventilated cool environment.

#### **M** WARNING

- Before cleaning the pulse oximeter or sensor, turn off the power and disconnect the DC power.
- The pulse oximeter should be kept clean. It is recommended to regularly clean
  the enclosure surface and the display screen. Cleaning the enclosure with nonetching cleaner such as soap and water.

#### **ACAUTION**

- · To avoid damaging the pulse oximeter:
- Do not use strong solvents such as acetone.
- Most cleaners must be diluted before use. Diluting should be according to the manufacturer's instructions
- Do not use abrasive materials (such as steel wool).
- Do not allow any liquid entering into the enclosure, and never immerse any part of the device into liquid.
- Do not leave any cleaning solution on the surface of any part of the device.

#### NOTE

- Wipe the pulse oximeter and sensor surface with medical alcohol, dry it naturally or with clean, dry, lint-free cloth.
- Bistos is not liable for effectiveness of using these chemicals for infectious disease control. Please consult the infectious disease control officers or experts of the hospital for advice.

#### 9.3 Disinfection

In order to avoid damage to the product, we recommend that the product is disinfected only when it is deemed necessary by the hospital maintenance procedures. We also recommend

30

that the instrument to be disinfected must first be cleaned.

#### **ACAUTION**

 To prevent damage to the pulse oximeter, do not disinfect the pulse oximeter with gas (EtO) or formaldehyde.

#### 10 Maintenance

#### MARNING

If the hospitals or institutions using this instrument can't implement a satisfactory maintenance schedule, it will result in device failure and may endanger human health.

#### 10.1 Checking

Check the following basic items before using the pulse oximeter:

- ✓ Check for any mechanical damage.
- ✓ Check all wires and accessories
- Check all instrument functions that may be used for pulse oximeter and ensure that the instrument is in good working condition.

If the instrument function has any sign of damage, do not use this pulse oximeter. Please contact the hospital's professional maintenance personnel or our customer service personnel

Every 6-12 months or after each repair, a comprehensive examination must be performed by trained and qualified technical service personnel, including functional safety checks. The specific inspection items are as follows:

- ✓ Environment and power meet the requirements.
- ✓ Device and accessories have no mechanical damage.
- ✓ The power supply has no wear, and the insulation is good.
- ✓ Specified accessories are used.
- ✓ Alarm system is functioning correctly.
- ✓ Battery performance meets the requirements.
- ✓ Monitoring functions are in good working condition.
- ✓ Ground impedance and leakage current meet the requirements.

All checks that require disassembling the instrument must be performed by qualified service personnel. Safety and maintenance checks may also be carried out by the Bistos personnel.

#### 10.2 Trouble shootings

Problem	Solution	
Device not power on	Check the battery. If the battery is low, it should be promptly charged.	
SpO <sub>2</sub> waveforms or	Is the red light on the finger sensor flashing? If not, there	

values do not displayed on the screen	might be poor contact. Check the patient cable and the connector.	
	Is the patient's arm under pressure? Never take blood pressure and SpO <sub>2</sub> measurements on the same arm	
	Is the environmental temperature too low? Never expose the patient's arm to cold air since this can affect readings.	
	Has all patient nail polish, especially blue or purple, been removed?	
SpO <sub>2</sub> values turn on and off during SpO <sub>2</sub> monitoring	During long term monitoring, patient movement might result in $SpO_2$ interruptions. Keep the patient stabilized. $SpO_2$ interruptions due to patient hand motion are normal.	

#### 10.3 Maintenance plan

The following tasks can only be done by qualified service personnel of Bistos. When the following maintenance is needed, please contact your service representative. Before testing or maintenance clean and disinfect the device

Inspection / Maintenance Item	Frequency			
Check the safety according to IEC 60601-1	At least once every two years, after replacing the power supply or the pulse oximeter falls down.			
Check all monitoring or measur- ing functions not listed	At least once every two years, or when you suspect that the measured value is not accurate.			

#### 11 Accessories

#### **M** WARNING

- Use the accessories specified in this manual. Using other accessories may damage the monitor, or cannot reach the safety and performance claimed in this manual.
- The operating and storage environment of the monitor should meet the requirements of the accessories. Please refer to the manual of the accessories for these requirements.
- If the packaging or accessories have any sign of damage, do not use such accessories.
- For SpO<sub>2</sub> Sensor, the normal life time is two years. Please replace in time.

#### Standard accessories are as follows:

No.	Description	QTY	Type-number
1	Adult Finger Clip SpO₂ Sensor Probe	1	Manufacturer: Unimed Medical Supplies,Inc
			U403-01

## 12 Specifications

## 12.1 Safety specifications

#### 12.1.1 Product category

In accordance with classification specified in the European Medical Device Directive 93/42/EEC, this pulse oximeter is Class IIb device. The pulse oximeter is classified as follows in accordance with IEC 60601-1:

Category Name	Specification	
Type of electric shock protection	Class II and internally powered equipment	
Electric shock protection grade	Type BF applied part	
Explosion protection grade	Common equipment, no explosion protection	
Liquid inlet protection grade	IPX2	
Operating mode	Continuous mode	
Movement	Hand-held equipment	

#### 12 1 2 Power

Power		
Adapter   Input: AC 100 ~ 240V (50/60 Hz) Output: DC 5V / 2.0A		
	3.7V Li-ion battery 3000 mA	
Rechargeable Battery	Operating Time(When it fully charged): 5 hours	
	Charging Time(Fully): 4 hours	

#### 12.2 Hardware specifications

Physical Characteris	tics
Dimensions	Main Unit: 84(W) X 158.5(H) X 34.5(D)
Weight	< 1.5 Kg for standard configuration
Display	
Туре	Color TFT touch screen LCD
Size	4.3"
Audio	
	Alarm tone (45 ~ 85 dB)
Speaker	Pulse tone
	Alarm sound meet the IEC 60601-1-8 standard require-
	ments

#### 12.3 Functional specifications

SpO <sub>2</sub>	
Standards compliant	ISO 80601-2-61:2011
Display range	0% ~ 100%
SpO <sub>2</sub> display resolution	1%

SpO <sub>2</sub> accuracy	±2% (at the range 70%~100%)(adult/pediatric mode) ±3% (at the range 70%~100%) (neonate mode)		
	not define when lower than 70%;		
	Measurement accuracy verification		
	The SpO <sub>2</sub> accuracy has been ve	erified in human experi-	
	ments by comparing with arte	rial blood sample refer-	
	ence measured with a CO-oximeter. Pulse oximeter		
	measurements are statistically		
	two-thirds of the measuremer	•	
	within the specified accuracy r	ange compared to CO-	
	oximeter measurements.		
	The accuracy of the oximeter has been validated by a clinical trial involving 12 healthy adult subjects - 4		
	women and 8 men. Among	, ,	
	subjects, light skin are 5 subjects, dark skin are 3 subjects, the age from 21 to 28.		
	Over the range of 70% to 100%	6. overall accuracy was	
	determined by calculating the		
	across all samples and is 1.44%.		
SpO <sub>2</sub> alarm limit range	Upper alarm limit	1%~100%	
	Lower alarm limit	0%~99%	
SpO <sub>2</sub> alarm signal	No delay when the "Alarm Delay" sets "off".		
generation delay	If the "Alarm Delay" sets more than "1s", the alarm		
	generation is delayed by the set time.		
SpO <sub>2</sub> value refresh	The data update period is 1s/time		
period			
	Low sensitivity (when	Within 6.8s	
	sensitivity sets 1)	Maril: 0.4	
Data averaging	Intermediate sensitivity	Within 3.4s	
	(when sensitivity sets 2) Advanced sensitivity (when	Within 3.4s	
	sensitivity sets 3)	WILIIII 3.45	
	Select "Alarm Delay", you can choose "off", "1s", "2s",		
Alarm condition delay	"3s", "4s", "5s", "6s", "7s" or "8s".		
1	33 , 43 , 33 , 03 , 73 0	. 55 .	

PR	
Display range	25~250bpm
Resolution	±1 bpm
Accuracy	±2% or ±2bpm,whichever is greater

BT-710 Operation manual

34

#### 13 Manufacturer's declaration on EMC

BT-710 needs special precautions regarding EMC (Electromagnetic compatibility) and needs to be used according to the EMC information provided in this user manual. Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect the BT-710 and should be kept at least 1 m away from the equipment.

#### NOTE

- Using unqualified accessories, sensors and cables will increase the electromagnetic emission and reduce the electromagnetic immunity of the device.
- Do not put the device close to other devices or stack together. When necessary, observe the device closely to ensure that it runs normally in the environment.
- The device requires special EMC protection, and it is necessary to install and maintain it in the environment that meets the following EMC information.
- Even if other devices comply with CISPR emission requirements, they may also cause interference to this device.
- When the input signal amplitude is smaller than the minimum amplitude specified in the technical specifications, it may result in inaccurate measurements.
- Mobile communication devices or wireless network devices may have an impact on the device

#### 13.1 Electromagnetic emissions

The BT-710 is intended for use in the electromagnetic environment specified below.			
The customer or the user of the BT-710 should assure that it is used in such an environment.			
Emissions test	Compliance	Electromagnetic environment - guidance	
RF emissions CISPR 11	Group 1	The BT-710 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF emissions CISPR 11	Class A	The BT-710 is suitable for use in all establishments other than domestic, and may be used in domestic establish-	
Harmonic emissions IEC 61000-3-2	Class A	ments and those directly connected to the public low- voltage power supply network that supplies buildings use	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	for domestic purposes, provided the following warning is heeded:  Warning: This BT-710 is intended for use by healthcare professionals only. This equipment/ system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the BT-710 or shielding the location.	

# 13.2 Recommended separation distances between portable and mobile RF communications equipment and BT-710

The BT-710 is Intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the BT-710 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BT-710 as recommended below, according to the maximum output power of the communications equipment.

Rated maxi- mum output	Separation distance according to frequency of transmitter [m]		
power of transmitter [W]	150 kHz to 80 MHz $d = 3.5\sqrt{p}$	80 MHz to 800 MHz $d = 3.5\sqrt{p}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3}\right] \sqrt{p}$
0.01	0.35	0.35	0.23
0.1	1.11	1.11	0.74
1	3.5	3.5	2.34
10	11.07	11.07	7.38
100	35	35	23.24

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1) At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

#### 13.3 Electromagnetic immunity

The BT-710 is intended for use in the electromagnetic environment specified below.

The customer or the user of the BT-710 should assure that it is used in such an environment.

THE CUSCOMET OF C	The customer of the user of the B1-710 should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance	
Electrostatic discharge (ESD)	±8 kV Contact	±8 kV Contact	Floors should be wood, concrete or ceramic tile. If floors are covered	
, ,	±15 kV air	±15 kV air	with synthetic material, the	
IEC 61000-4- 2:2009			relative humidity should be at least 30 %.	
Electrical fast	±2 kV for power	±2 kV for power	Mains power quality should be	
transient/burst	supply lines	supply lines	that of a typical commercial or	
IEC 61000-4-	±1 kV for in- put/output lines	±1 kV for in- put/output lines	hospital environment.	
4:2004	(>3m)	(>3m)		
Surge	±1 kV differential	±1 kV differential	Mains power quality should be	
	mode	mode	that of a typical commercial or	
IEC 61000-4-	±2 kV common	±2 kV common	hospital environment.	
5:2006	mode	mode		
Voltage dips,	< 5 % UT (> 95 %	< 5 % UT (> 95 %	Mains power quality should be	
short interrup-	dip in Uτ) for 0.5	dip in Uτ) for 0.5	that of a typical commercial or	
tions and	cycles	cycle	hospital environment. If the user	

voltage variations on power supply input lines IEC 61000-4- 11:2004	40 % <i>U</i> τ (60 % dip in <i>U</i> τ ) for 5 cycles  70 % <i>U</i> τ (30 % dip in <i>U</i> τ) for 25 cycles  <5 % <i>U</i> τ (> 95 % dip in <i>U</i> τ) for 5 s	40 % <i>U</i> τ (60 % dip in <i>U</i> τ ) for 5 cycles  70 % <i>U</i> τ (30 % dip in <i>U</i> τ) for 25 cycles  <5 % <i>U</i> τ (> 95 % dip in <i>U</i> τ) for 5 s	of the BT-550 image intensifier requires continued operation during power mains interruptions, it is recommended that the BT-710 be powered from an uninterruptible power supply.	
Power frequen- cy (50 Hz and 60 Hz) magnet- ic field IEC 61000-4- 8:2010	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environ- ment.	
NOTE UT is the a.c	NOTE UT is the a.c. mains voltage prior to application of the test level.			

The BT-710 is intended for use in the electromagnetic environment specified below.

The customer or the user of the BT-710 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level
Conducted RF	3 Vrms	3 Vrms
IEC 61000-4-6:2009	150 kHz to 80 MHz	
Radiated RF	3 V/m	3 V/m
IFC 61000-4-3	80 MHz to 2.5 GHz	

#### Electromagnetic environment - guidance

Portable mobile RF communications equipment should be used no closer to any part of the BT-550,

including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

#### Recommended separation distance

 $d = 1.2, \sqrt{p} (d = 3.5, \sqrt{p})$ 

 $d = 1.2\sqrt{p}$  (Resp:  $d = 3.5\sqrt{p}$ ) 80 to 800MHz

 $d = 1.2 \sqrt{p}$  800M to 2.5GHz

where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and *d* is the recommended separation distance in meters (m).

Field strengths from fixed RF transmitters, as deter-mined by an electromagnetic site survey <sup>a</sup>, should be less than the compliance level in each frequency range. <sup>b</sup>

Interference may occur in the vicinity of equipment marked with the following symbol:

0.0

NOTE 1) At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

\* Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BT-710 is used exceeds the applicable RF compliance level above, the BT-550 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the BT-710.

Over the frequency range 150 kHz to 80MHz, field strengths should be less than 3 V/m.

#### **Product Warranty**

Product Name	Pulse oximeter
Model Name	BT-710
Serial No.	
Warranty Period	2 Years
Date of Purchase	
Customer	Hospital: Address: Name: Telephone:
Sales Agency	
Manufacture	Bistos Co., Ltd.

Thank you for purchasing BT-710.

#### Service Telephone and Fax. Numbers

Telephone: +82 31 750 0340 Fax: +82 31 750 0344

Bistos Co., Ltd.

7<sup>th</sup> FL., A Bldg., Woolim Lions Valley 5-cha, 302, Galmachi-ro, Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea

> www.bistos.co.kr bistos@bistos.co.kr

EC REP Obelis s.a

Bd. Général Wahis 53 1030 Brussels / BELGIUM Telephone: +32 2. 732.59.54

Fax: +49 2. 732.60.03





<sup>\*</sup> This product is manufactured and passed through strict quality control and inspection.

<sup>\*\*</sup> Compensation standard concerning repair, replacement, refund of the product complies with "Framework Act on Consumers" noticed by Fair Trade Commission of Republic of Korea.



# BT-720 Patient monitor Operation Manual



Keep this manual for future reference

P/N: 720-ENG-OPM-EUR-R03

# **Proprietary Material**

Information and descriptions contained in this manual are the property of Bistos Co., Ltd. and may not be copied, reproduced, disseminated, or distributed without express written permission from Bistos Co., Ltd.

Information furnished by Bistos Co., Ltd is believed to be accurate and reliable. However, no responsibility is assumed by Bistos for its use, or any infringements of patents or other rights of third parties that may result from its use. No license is granted by implication or otherwise under any patent or patent rights of Bistos Co., Ltd.

The information contained herein is subjects to change without notice.

## Prepared by:

Bistos Co., Ltd. 7<sup>th</sup> FL., A Bldg., Woolim Lions Valley 5-cha, 302, Galmachi-ro, Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea

Telephone: +82 31 750 0340 Fax: +82 31 750 0344

Revision R03 October, 2020

Printed in Korea Copyright © Bistos Corporation 2018. All rights reserved.

# **Contents**

0.	Safety information	4
	0.1 General precautions, warnings and cautions	
	0.2 Shock hazards	7
	0.3 Battery warnings	
	0.4 General precautions on environment	9
1.	System basics	9
	1.1 Intended use	
	1.2 Operating principle	
	1.3 System configurations	
	1.4 Product outlook	
	1.5 Description of monitor	
	1.6 Understanding the display	
2.	·	
۷.	2.1 Installation	
	2.2 Connecting to power	
3.	Basic operations	14
	3.1 Turn on	14
	3.2 Turn off	15
4.	Setup the monitor	15
5.		
	5.1 Overview	
	5.2 Safety information	
	5.3 Monitoring steps	
	5.4 Setting SpO <sub>2</sub>	18
	5.6 Technical description	
6.	NIBP	19
	6.1 Overview	
	6.2 Safety information	20
	6.3 Measurement limits	
	6.4 Measurement procedure	
	6.5 Setting NIBP	22
	6.6 Clean and disinfection method of NIBP cuff	23
7.	Review	24
8.	Alarm	25
9.	Battery	26
	9.1 Overview	
	9.2 Battery usage guide	26
	9.3 Checking battery performance	27
	9.4 Battery recycling	27

10.	Cai	ring and cleaning	27
		Overview	
	10.2	Cleaning	27
	10.3	Disinfection	28
11.		aintenance	
		Checking	
	11.2	Maintenance plan	29
12.	Ac	cessories	· <b>2</b> 9
13.		ecifications	
		Safety specification	
	13.2	Hardware specifications	30
	13.3	Functional specification	31
14.	Co	mmon faults and maintenance	32
15.	Ma	anufacturer's declaration on EMC	33
	15.1	Electromagnetic emissions	33
		Recommended separation distances between portable and mobile RF con	
	e q u	ipment and BT-720	3 4
	15.3	Electromagnetic immunity	35
Pro	duct	Warranty	37

# **O Safety information**

Before using BT-720 Patient monitor, read this entire manual and be fully understood the following safety information to prevent injury of patient and user.

## **Symbols Used**

The following symbols identify all instructions that are important to safety. Failure to follow these instructions can lead to injury or damage to the patient monitor. When used in conjunction with the following words, the symbols indicate:

wing words, the symbols indicate:		
<b>MARNING</b>	Can lead to serious injury or death.	
<b>A</b> CAUTION	Can lead to minor injury or product/property damage	
following symbols	are placed on product, label, packaging and this manual in order to stand for	

The following symbo the information abo	ols are placed on product, label, packaging and this manual in order to stand fo ut:
	Used to identify safety information. Be well-known this information thoroughly before using BT-720.
	Used to identify safety information.  Be well-known this information thoroughly before using BT-720
IPX1	Indicates the protection level against the ingress of liquid.  IPX1 is protection against some water drops falling vertically.  It correspond the device, patient monitor.
IPX2	Indicates the protection level against the ingress of liquid.  IPX2 is protection from some water drops when the device is tilted up to and including 15°.  It correspond the accessories for SpO2.
<b>(3)</b>	Refer to operation manual. Read manual before placing the device.
	Indicates DC power supply.  Indicates the device is in the battery operation mode.
Э́Д	Indicates nurse call interface.
器	Indicates network interface.
<b>←</b>	Indicates USB interface.
<b>♦</b> ••	Indicates power adapter polarity.
سا	Indicates the production date.
***	Indicates the manufacturer.
SN	Indicates the serial number of the device.
EC REP	Indicates the authorized representative in the European Community of manufacturer.
<b>∤</b>	Indicates a defibrillation-proof type BF applied part.
- <b>W</b>	Indicates a defibrillation-proof type CF applied part.
	Indicates CLASS II equipment.(Adapter)
53	Indicates the date after which the medical device is not to be used.

<del>**</del>	Indicates to keep the device dry.
I	Indicates the medical device that can be broken or damaged if not handled carefully.
<u> </u>	Indicates to keep upright
XI6	Indicates the maximum stacking limit.
1	Indicates the temperature limitation for operation, transport and storage.
<b>%</b>	Indicates the humidity limitation for operation, transport and storage.
<b>\$•</b> \$	Indicates the range of atmospheric pressure to which the medical device can be safely exposed.
LATEX	Indicates the device contains natural rubber latex.(Accessories)
	Indicates the packing material is recyclable.
X	Indicates to not dispose the device together with unsorted municipal waste(for EU only). The solid bar symbol indicates that mains adapter is put on the market after 13 August 2005.

## 0.1 General precautions, warnings and cautions

- Examine the patient monitor and any accessories periodically to ensure that the cables, adapter cords and instruments do not have visible evidence of damage that may affect patient safety or performance. The recommended inspection interval is once per week or less. Do not use the patient monitor if there is any visible sign of damage.
- Only the DC power adapter supplied with the BT-720 is approved for use with the device.
- Do not attempt to service the BT-720 patient monitor. Only qualified service personnel by Bistos Co. Ltd. should attempt any needed internal servicing.
- Perform periodic safety testing to insure proper patient safety. This should include leakage current measurement and insulation testing. The recommended testing interval is once per year.
- If the hospital or healthcare institutions using this device fail to implement a satisfactory
  maintenance schedule, it will result in device failure and may endanger the patient's
  safety.
- Use the patient monitor under the conditions specified in this operation manual. Beyond
  the conditions, the patient monitor may not function properly and the measurement
  results may not accurate and may result in device failure or endangering the patient's
  safety.
- Do not operate the BT-720 patient monitor if it fails to pass the power on self-test procedure.
- During the operation, do not disconnect any cable.
- The BT-720 patient monitor is intended to be used by clinical professionals or trained doctors, nurses or laboratory assistant.
- Do not service and maintain or clean the device including accessories while in use with a
  patient.
- Using the device to one patient at a time.



## WARNING

- Thoroughly read and understand the manual prior to use of the BT-720. Failure to
  do so could result in personal injury or equipment damage.
- The device is intended for clinical patient monitoring, and only trained and qualified doctors and nurses should use the device.
- The alarm volume, upper and lower alarm limits should be set according to the
  actual situation of the using environment. Do not just rely on audio alarm system
  while monitoring the patient, because too low alarm volume or muted alarm may
  result in notice failure of alarm situation and endanger the patient's safety. Please
  pay close attention to the actual clinical status of the patient.
- Use only the power adapter supplied with monitor.
- Position the monitor where it is easy to de-energize the monitor when needed.
- Do not open the enclosure to avoid an electric shock. Any repair and upgrade of monitor should be done by service personnel trained and authorized by Bistos. Co., Ltd.
- When handling packaging materials, abide by local laws and regulations or hospital waste disposal regulations. Keep the packaging materials away from children.
- Do not use in the presence of flammable anesthetics to prevent explosion or fire.
- Install the power lines and cables of accessories carefully to avoid patient entanglement or suffocation, cables tangled or electrical interference.
- When the monitor is used together with electrosurgical devices, the user (a doctor or a nurse) should ensure the safety of the patient and instrument.
- The physiological wave, physiological parameters and alarm information displayed on the monitor are only for the doctor's reference and should not be directly used as the basis for clinical treatment.
- · This is not a therapeutic device.
- For patients with pacemakers, the cardio tachometer may count the pacemaker pulse in case of a cardiac arrest or arrhythmias. Never rely solely on the cardio tachometer alarm. Closely monitor the patients with pacemaker. For the inhibition of the device on pacemaker, refers to this manual.
- Use of accessories other than those listed and approved for use with this product may result in increased emissions or decreased immunity.
- Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this manual. In addition, portable and mobile RF communications equipment can affect medical electrical equipment.
- The equipment shall not be used adjacent to other devices unless verification of normal operation in the configuration in which it is to be used can be achieved.
- Keep matches, and all other sources of ignition, out of the room in which the
  patient monitor is located. Textiles, oils, and other combustibles are easily ignited
  and burn with great intensity in air enriched with oxygen. Personal injury or
  equipment damage could occur.
- A fire and explosion hazard exists when performing cleaning or maintenance procedures in an oxygen-enriched environment.
- The patient monitor has been validated with the accessories and options listed in this manual and found to comply with all relevant safety and performance requirements applicable to the device. It is therefore the responsibility of the person or organization who makes an unauthorized modification, or incorporates an unapproved attachment to the device.
- An operator may only perform maintenance procedures specifically described in

this manual.

Do not remove the covers of a BT-720 yourself to avoid damage to the equipment and unexpected electrical shock. Only qualified Bistos service engineer must repair or replace components.



- Please install or carry the instrument properly to prevent damage due to falling, collision, strong vibration or other mechanical force.
- Avoid instrument splashed by water.
- Avoid high temperatures, the instrument should be used within a temperature range of 5 °C ~ 40 °C ₀
- Avoid using instrument in the environment such as pressure is too high, poor ventilation, dusty, or contain salt, sulfur gas and chemical.
- Before using the monitor, check the monitor and accessories if there is damage that may affect patient safety. If there is obvious damage or aging, replace the parts before use. The replacement should be made with same parts of original parts.
- Before powering on the device, make sure that the power used by the device complies with the supply voltage and frequency requirements on the equipment label or in the Operator's Manual.
- Equipment should be tested at least once a year, the test should be done and recorded by trained, have security testing knowledge and experienced personnel. If there are any problems in the tests, they must be repaired.
- When the instrument and accessories are about to exceed the useful life (expected service life: 5 years), it must be treated in accordance with relevant local laws and regulations or the hospital's rules and regulations.
- Do not connect to other equipment or network which not specified in the instruction for use, in risk of external high voltage.
- Do not connect any equipment or accessories that are not approved by the manufacturer or according to IEC 60601-1 to the monitor. The operation or use of non-approved equipment or accessories with the monitor is not tested or supported, and monitor operation and safety are not guaranteed in such a case.
- Any non-medical equipment (such as the external printer) is not allowed to be used within the patient vicinity (1.5m/6ft.).
- Parts and accessories used must meet the requirements of the applicable safety standards, and/or the system configuration must meet the requirements of the medical electrical systems standard.
- Ensure that the conductive parts of electrodes and associated connectors, including neutral electrodes, do not come in contact with earth or any other conducting objects.
- Protection of ME EQUIPMENT against effects of discharge of a cardiac defibrillator depends on use of proper cables.

## 0.2 Shock hazards



## WARNING

- Unplug the monitor from its power source prior to cleaning or maintenance to prevent personal injury or equipment damage.
- Some chemical cleaning agents may be conductive and leave a residue that may permit a build-up of conductive dust or dirt. Do not allow cleaning agents to contact electrical

- components and do not spray cleaning solutions onto any of these surfaces. Personal injury or equipment damage could occur.
- Do not expose the unit to excessive moisture that would allow for liquid pooling. Personal injury or equipment damage could occur.
- Do not touch the patient and signal input/output parts simultaneously
- Due to the risk of electrical shock hazard, only qualified personnel with appropriate service documentation should service the monitor.

## 0.3 Battery warnings



# **WARNING**

- Improper operation may cause the internal lithium ion battery to be hot, ignited or exploded, and it may lead to the decrease of the battery capacity. It is necessary to read the operation manual carefully and pay more attention to warning message.
- Do not open the battery compartment. Only the qualified service personnel authorized by the manufacturer can open the battery compartment and replace the battery, and batteries of same model and specification should be replaced.
- Be careful when connecting the battery with polarity.
- Do not use the battery near fire or environmental temperature exceeds 60  $^{\circ}$ C. Do not heat or splash the battery or throw it into fire or water.
- Do not destroy the battery. Do not pierce battery with a sharp object such as a needle. Do not hit with a hammer, step on or throw or drop the battery. Do not disassemble or modify the battery. The battery can heat, smoke, deformation or burning.
- When leakage or foul smell is found, stop using the battery immediately. If your skin or cloth comes into contact with leaked liquid, cleanse it with clean water at once. If the leaked liquid splashes into your eyes, do not wipe them. Irrigate them with clean water first and go to see a doctor immediately.
- Properly dispose of or recycle the depleted battery according to local regulations.

## 0.4 General precautions on environment

Do not keep or operate the equipment under the environment listed below.

	Avoid placing in an area exposed to moisture. Do not touch the equipment with wet hand.		Avoid exposure to direct sunlight
	Avoid placing in an area where high variation of temperature exists. Operating temperature ranges from $5^{\circ}$ ~ $40^{\circ}$ C. Operating humidity ranges from $30^{\circ}$ ~ $85^{\circ}$ M.		Avoid in the vicinity of electric heater.
Suggest and the suggest and th	Avoid placing in an area where there is an excessive humidity rise or ventilation problem.		Avoid placing in an area where there is an excessive shock or vibration.
	Avoid placing in an area where chemicals are stored or where there is in danger of gas leakage.	100	Avoid dust and especially metal material enter into the equipment
<del>OOZ</del> is	Do not disjoint or disassemble the device. Bistos Co., Ltd. does not have liability of it.		Power off when the equipment is not fully ready to operate. Otherwise, the equipment could be damaged

## 1 System basics

## 1.1 Intended use

The BT-720 Patient Monitors acquire the physiological signals for non-invasive blood pressure (NIBP), pulse rate (PR) and blood oxygen saturation (SpO<sub>2</sub>). The signals are converted into digital data and processed, examines the data for alarm conditions and display the data. The monitor also provides operating control for the user. The patient monitor intend to use in hospital clinical area such as intensive care units, operating room, emergency department, to provide additional information to the medical and nursing staff about the physiological condition of the patient. The BT-720 patient monitors are intended to be used only under regular supervision of clinical personnel. It is suitable for adult and pediatric, neonate. The intended locations of use are hospitals and clinics.

#### 1) Intended patient population

Adult (>18 years adults) and Pediatrics (30 days < and <18 years) and Neonate (0 days < and <30days)</li>

#### 2) Intended user profile

- Doctor, physicians or nursing staff who is qualified personnel
- Basic experiences or knowledge on medical field, especially on patient monitoring
- Trained or requested to read IFU before use

#### 3) Environment of use

- Hospital and clinic
- Requirements: Stable power source

#### 4) Scope of application

This monitor is suitable for bedside monitoring of patient. This monitor enables blood oxygen saturation (SpO<sub>2</sub>), pulse rate (PR) and monitoring. It is equipped with a replaceable built-in battery to provide convenience for the patient movement in hospital.

### 5) Indications and contraindications

## Blood oxygen saturation (SpO<sub>2</sub>)

## Indication:

- Monitoring effectives of oxygen therapy
- A reading is needed to facilitate the completion of an early warning score to inform clinical assessment
- Sedation or anesthesia
- Transport of patients who are unwell and require oxygenation assessment
- Haemodynamic instability (e.g. cardiac failure or Myocardial Infarction)
- Respiratory illness e.g. asthma, chronic obstructive pulmonary disease
- Monitoring during administration of respiratory depressant drugs, e.g. opiate epidural or patient-controlled analgesia.
- Assessing oxygen saturation during physical activity e.g. in pulmonary rehabilitation

#### Contraindications

 Pulse oximetry does not give an indication of haemoglobin so if the patient is profoundly anaemic then their oxygen saturation may by normal but they may still be hypoxic

Source: NHS. "Clinical Procedure\_ Procedure for Pulse Oximetry/SPO2". Wirral Community NHS Trust. Sep, 2013

#### Non-invasive blood pressure (NIBP)

## Indication:

- To determine a patient's blood pressure
- Screen for hypertension
- Following the effect of anti-hypertensive treatments in a patient to optimize their management
- Assessing a person's suitability for a spot or certain occupations

- Estimation of cardiovascular risk
- Determining for the risk of various medical procedure
- Figuring out whether a patient is clinically deteriorating or is at risk.

#### Contraindications

- Oscillometric blood pressure devices may not be accurate in patients with weak or thready pulse
- In patients with heart beats below 50 beats/minutes, even if the rhythm is regular, some of the semi-automatic devices are unable to reduce their deflation rate sufficiently so that too rapid a falling in cuff pressure results in underestimation of systolic blood pressure and overestimation of diastolic blood pressure.
- Do not apply to limb with AV fistula, significant injury or burn, or lymph node removal post mastectomy.

Source: [1] NHS. "Clinical Procedure\_ Procedure for Blood Pressure Monitoring". Wirral Community NHS Trust. Dec, 2013

[2] Clinical Quality& Patient Safety Unit, QAS. *Clinical Practice Procedures:*Assessment/Non-invasive blood pressure. Queensland Government, 2016. https://www.ambulance.qld.gov.au/clinical.html

#### 1.2 Operating principle

Refer to the chapters for every physiological parameter from chapter 5 to chapter 6.

## 1.3 System configurations

Basic configuration of BT-720

- Main body with 4.3" touch screen and built-in lithium-ion battery
- Adult SpO2 probe and extension cable
- AC/DC adapter

## Options of BT-720

Non-invasive blood pressure cuff

Picture	Name	Description	Qty
	Adult SpO2 sensor (standard)	SpO2 sensor for adult	1ea
	SpO₂ extension cord (Standard)	Cord to connect the SpO2 sensor and main body	1ea
	Adult NIBP cuff (Option)	Measures NIBP for adult	1ea
	NIBP extension tube (Option)	Tube to connect the NIBP cuff and main body	1ea
3	Adapter (Standard)	For power supply	1ea

## 1.4 Product outlook



Figure 1-1: Front view




Figure1-2 : Rear view

## 1.5 Description of monitor




Figure 1-3: Front view

rigures 3. Front view			
	Name	Description	
1	NIBP tube interface. This is option. When you purchase the mon only with SpO2, this interface will be blocked.		
2	SpO2	SpO2 cable interface	
3	Display area Display the waveform and measured value		
4	Alarm indicator	Indicates the priority of physiological alarm and technical alarms in different colors and flashing frequencies.  - High priority: Red, fast flashing (1.4 ~ 2.8 Hz)  - Medium priority: Yellow, slow flashing (0.4 ~ 0.8 Hz)  - Low priority: Yellow, constant on	
5	[Setting]	Enter to the setting mode. Press again to close the setting mode.	
6	[NIBP]	Start and stop the non-invasive blood pressure measurement manuall This is option. When you purchase only with SpO2 this button will be deleted.	
7	[Alarm reset]	To reset the alarm condition.	

8	DC power indicator	Turned on when the monitor is being powered by the adapter.
9	Battery indicator	<ul> <li>On: The battery is being charged or has been fully charged.</li> <li>Off: The battery has not been installed.</li> <li>Flashing: The monitor is being powered by the battery.</li> </ul>
10	[Power]	<ul> <li>Power On: Press down the key more than 2 seconds.</li> <li>Power Off: Press down the keys more than 2 seconds and the system will display the alarm message "The system will shut down 3 seconds".</li> </ul>




Figure1-4: Rear view

	Name	Description
1	Network port	For CMS
2	Type A USB	For software upgrade
3	Type B USB	Reserved for data transfer
4	Power adapter	15V, 2.4A adapter

# 1.6 Understanding the display




Figure 1-5: Standard display

	Description
	Current alarm message. When an alarm occurs, this area will displayed yellow or white
	depending on the alarm type.
2	SpO <sub>2</sub> value. Display the measured SpO <sub>2</sub> value.
3	SpO₂ upper alarm limit. Display the user set upper alarm limit
4	SpO <sub>2</sub> waveform. Display the measured SpO <sub>2</sub> waveform.
5	SpO₂ lower alarm limit. Display the user set lower alarm limit

6	Perfusion Index. Display the measured perfusion index.
7	Sweep speed. Display the user set SpO₂ waveform speed.
8	RTIP. Display the current measured cuff pressure value.
9	MAP. Display the average pressure value.
10	Diastolic blood pressure value.
11	Systolic blood pressure value.
12	The unit of NIBP.
13	Battery status.
14	Patient type.
15	Display current time.
16	The unit of SpO <sub>2</sub> .

#### 1.7 Essential performance

This device Patient Monitor provides various patient vital signs such as pulse rate, blood oxygen saturation and non-invasive blood pressure by placing the sensors to the appropriate site of patient. The device is composed with display, control circuit and panel, and input part for various sensors. It detects SpO2, PR and NIBP using specific sensors and cuff. The detected analog signal amplifies and converted to digital. This concerted data feed to the CPU and converted to the display format as number and waveform. This device is incorporated with alarm system. The alarm generated when the detected signal range is beyond the user set alarm limits.

## 2 Preparing for operations

## 2.1 Installation

To ensure normal working of the monitor, read this chapter before use, and install as required.



## WARNING

- All analog and digital devices connected to the monitor must be certified by IEC standards (e.g. IEC 60950 Data processing equipment standard and IEC 60601-1 Medical equipment standard). Furthermore, all configurations shall comply with valid version of IEC 60601-1 standard. The personnel connecting additional devices to the input / output signal ports are responsible for the compliance with IEC 60601-1 standard. If there is any question, please contact Bistos.
- If the patient cable interface and network interface are connected with multiple devices, the total electric leakage current cannot exceed the allowable value.
- The copyright of monitor software belongs to Bistos. Without permission, any
  organization or individual shall not interpolate, copy or exchange by any means or form.
- When the monitor is combined with other devices, it must comply with IEC 60601-1:2005 + A1:2012, and should not be connected with multi-socket outlet or extension cord.
- Do not connect the device on other equipment or network, to which a signal input/output part may be connected.

Prior to installation, the operator must ensure that the following space, power, environmental requirements are met.

#### 2.1.1 Unpack and check

BT-720 patient monitor was inspected rigorously at the factory before delivery, in order to avoid being hit when transported, carried out careful packaging. Before unpacking, carefully inspect the package. If any damage, please immediately contact the Bistos. Unpack in the correct way, carefully remove the monitor and accessories from the box and check with the packing list. Check if there is any mechanical damage, the all listed are completely packed. If you have questions, please contact the marketing department of Bistos or agency.

Please keep the packing box and materials for use in future transporting or storage.

## 2.1.2 Placement requirements

Equipment installation must meet:

 Ensure that the operating floor and the monitor have enough space for connecting the accessory wires.

## 2.1.3 Power requirements

DC power supply adapter

Input: A.C. 100 V ~ 240 V, 50/60 Hz

Output: D.C. 15 V, 2.5 A

Built-in rechargeable lithium-ion battery: D.C. 11.1 V, 4400 mAh

## 2.1.4 Environmental requirements

The storage, transport and use of the monitor must meet the following environmental requirements.

Operating	Ambient temperature	5°C ~ 40°C
Operating environment	Relative humidity	30 % ~ 85 % (Non-condensing)
	Atmospheric pressure	700 ~ 1060 mbar (hPa)
Transportation	Prevent severe shock, vibration, rain and snow splashing during transport.	
Storage	The packaged monitor should be stored in well-ventilated room with ambient temperature -20 $^{\circ}$ C $^{\circ}$ 60 $^{\circ}$ C, relative humidity 0 $^{\circ}$ 95 % (Noncondensing), atmospheric pressure 700 $^{\circ}$ 1060 mbar(hPa), and without corrosive gases.	

The operating environment of the monitor should avoid noise, vibration, dust, corrosive or flammable and explosive materials. In order to allow air flowing smoothly and achieve good heat dissipation, at least 2 inches (5cm) clearance should be kept around the device.

When the device is moved from one environment to another, the device may have condensation due to the differences in temperature or humidity. In this case, wait until the condensation disappears before using the device.



## WARNING

Ensure that the monitor is used under specified environment. Fail to do this, the
technical specifications declared in this manual may not be met and it may result in
damage to equipment and other unforeseen consequences.

## 2.2 Connecting to power



## WARNING

- Do not try to open the monitor when the power is connecting.
- During the operation, do not disconnect any cable.

#### Connect to power adapter in the following steps:

- Make sure that the AC power supply meets the following specifications: a.c.100V-240V, 50/60Hz.
- Use the power adapter provided with the monitor. Plug the power adapter into the power connector of the monitor, and plug the other end of the power adapter into the mains (low voltage power supply network facilities) power outlet.

## 3 Basic operations

## 3.1 Turn on

## 3.1.1 Check the monitor

- Before turn on the monitor, check whether there is mechanical damage to the monitor, and whether the external cables and accessories are connected correctly.
- Plug the power adapter into the AC power outlet. If using battery power, make sure the battery is fully charged.
- Check all the functions required for patient monitoring to make sure that the monitor operates properly.



#### WARNING

If the monitor is damaged, or fails to work normally, do not use it for patient monitoring.
 Please contact the maintenance personnel or Bistos immediately.

## 3.1.2 Start the monitor

If finish to check the monitor, it is ready to start the monitor.

Press the [Power] key and the system enters the main interface within seconds.

- If any fatal error occurs during self-test, the system will alarm. If this case persists, please stop to using the monitor and contact the maintenance personnel or Bistos.
- Check all available monitor functions to ensure that the monitor operate properly.
- If the monitor equipped with a battery, charge the battery after each use to ensure sufficient power.
- After unpacking, when use the monitor first time, the monitor should be powered with adapter.

#### 3.1.3 Connect the sensors

Connect the required sensor to the monitor and the monitoring site of patient.

## 3.1.4 Start monitoring

Start monitoring in the following steps:

- Check if the patient cable and the sensor are connected properly.
- Check if the settings of the monitor are corrects, such as patient type.
- For the details of parameter measurement or monitoring, see the appropriate section.
- The operator can operate according to their own habits, standing in front, left or right of the monitor, easy to observe and operate the monitor.

## 3.2 Turn off

Turn off the monitor in the following steps

- Disconnect the cables and sensors connected to the patient.
- Press and hold the [Power] key for 2 seconds to pop up the 3 seconds countdown window, and the monitor turns off in 3 seconds.



 If the monitor is not turned off properly, you can simply disconnect the power to shutdown forcibly. But the forced shutdown may cause data loss, and it is not recommended.

# 4 Setup the monitor

In the main screen, Press the [Setting] key enter the setup menu.



- Select "Patient Type", you can choose "Adult" or "Pediatric" or "Neonate".
- Select "Alarm Delay", you can choose off, 1s, 2s, 3s, 4s, 45s, 5s, 6s, 7s or 8s.
- Select "Alarm Volume", and enter value (Range: 0-9).

- Select "Backlight", and enter value (Range: 0-5), 5 is the brightest.
- Select "Alarm Reminder Signal", you can choose "On", or "Off".
- Select "Alarm Reminder Interval", you can choose "1min", "2min" or "3min".



- Select "Date Format", and set the date format in accordance with custom
  - "YYYY-MM-DD": Year- Month-Day.
  - "MM-DD-YYYY": Month -Day-Year.
  - "DD-MM-YYYY": Day-Month-Year.
- Date (YMD)": Set the year, month, and day.
- Fig. "Time (24H)": Set the hour, minute and second.



- Select "Screen Setup", set "SpO2+NIBP".
- Select "Language", and select the option as needed:
  - "English": The interface language of the monitor is English.
  - "Türkçe": The interface language of the pulse oximeter is Turkish.
  - "Español": The interface language of the pulse oximeter is Spanish.
  - "Français": The interface language of the pulse oximeter is French.  $\,$
  - "Polski": The interface language of the pulse oximeter is Polish.
  - "Italiano": The interface language of the monitor is Italian.
  - "Deutsch": The interface language of the monitor is German.
  - "Magyar": The interface language of the monitor is Hungary.
  - "Portuguese": The interface language of the monitor is Portugal.
- Select "Network", you can choose "On", or "Off".
- Select "Default", back to the initial state.
- Select "Network Setup", set up the network.



- Select "Passwords", set Passwords.
- Select "SpO2", set "S-089".

# 5 SpO<sub>2</sub>

## 5.1 Overview

Blood oxygen saturation ( $SpO_2$ ) is the percentage of oxyhemoglobin (HbO2) capacity bound by oxygen in the blood in the total hemoglobin (Hb) capacity that can be combined, that is, the concentration of oxygen in the blood.

The principle for monitoring the pulse  $SpO_2$  is to fix the probe fingerstall on the patient's finger or toe, use the finger (or toe) as a transparent container for hemoglobin, use 660nm wavelength red light and 950nm near-infrared light as the incident light, maximum output power is 300 mW, measure the light transmission intensity through the tissue bed, and calculate the concentration of hemoglobin and  $SpO_2$ . The passing lights depend on a variety of factors, most of which are constant. However, one of these factors, the arterial blood flow, changes with time, as it is pulsating. By measuring the light absorbed during pulsating, it is possible to obtain the arterial blood  $SpO_2$ . Detection pulsation can give a "plethysmography" wave and pulse rate signal.

The main screen displays "SpO<sub>2</sub>" value and "plethysmography" wave.

This monitor applies to measure  $SpO_2$  of adults (>18 years) and pediatric (<18 years,>30 days), neonate (<30 days). Contact  $SpO_2$  probe to Patient's finger (or toe) to get " $SpO_2$ " value and "plethysmography" wave.

SpO<sub>2</sub> function of this monitor has been calibrated in factory.

The monitor is defibrillation proof, so the monitor operates normally after defibrillation.

## 5.2 Safety information



#### WARNING

- Please use SpO2 sensor supplied from Bistos, operate in accordance with the Manual, and observe all warnings and precautions.
- Before monitoring, check whether the sensor cable is normal. When SpO₂ sensor cable is unplugged from the socket, the screen will display "SpO2 Sensor Off" error message, and trigger an audible and visual alarm simultaneously.
- If the sensor or sensor packaging has signs of damage, do not use this SpO<sub>2</sub> sensor; return it to the manufacturer.
- If there is carboxyhemoglobin, methemoglobin or dye diluted chemical, the SpO<sub>2</sub> value will have deviation.
- When the patient has a tendency to hypoxia, use the oximeter to analyze blood samples in order to fully grasp the patient's condition.
- Do not put the sensor on limbs with arterial duct or intravenous tube.
- Do not intertwine electrosurgical equipment cable with the sensor cable.
- Avoid using the monitor and sensors while using the NMR equipment, in order to avoid severe burns to the patient as a result of induced currents.
- During long time continuous monitoring of a patient, check the position of SpO<sub>2</sub> sensor once every 2 hours, and move properly when the skin changes or every four hours.
   Some patients may require more frequent inspection, such as patients with perfusion disorders or sensitive skin, because persistent and prolonged monitoring may increase unpredictable skin changes, such as allergies, redness, blistering or pressure necrosis.
- Before using, verify compatibility between the monitor, probe and cable, otherwise it
  may cause injury to the patient.
- Functional testers cannot be used to evaluate the accuracy of pulse oximetry and pulse oximetry.
- SpO2 low alarm limit cannot be less than 85.

## NOTE

- Do not put the oxygen probe and blood pressure cuff on the same limb, because blood flow occlusion during blood pressure measurement will affect the SpO<sub>2</sub> readings.
- The monitor cannot be used to verify the accuracy of SpO<sub>2</sub> probe and SpO<sub>2</sub> equipment.

## 5.3 Monitoring steps


- . Select the appropriate SpO<sub>2</sub> sensor according to the patient.
- 2. Turn on the monitor, and connect the SpO<sub>2</sub> lead wire to the monitor.
- 3. Clean the measurement site, such as finger with nail polish.
- 4. Put the SpO<sub>2</sub> sensor probe on the patient's finger or toe.
- 5. Select the appropriate alarm settings.
- 6. Start monitoring.

#### NOTE

Turn on the monitor, plug in SpO<sub>2</sub> probe and connect patient's finger (or toe), monitor displays SpO<sub>2</sub> wave, "SpO2 Pulse Search" displayed in the technical alarm area until the monitor measured SpO<sub>2</sub> value and pulse rate. "SpO2 Search Timeout" displayed in the technical alarm area until the monitor measured pulse rate. Check the sensor mounting position, whether the sensor is damaged or sensor type. Reconnect the sensor or use new sensor.

## 5.4 Setting SpO<sub>2</sub>

Select SpO<sub>2</sub> parameter area to enter the SpO<sub>2</sub> set interface



- Select "Wave Speed" and set wave speed to "12.5mm/s" or "25mm/s". The faster speed, the smoother wave.
- Select "Wave Mode", and set the wave drawing mode to "Scan" or "Fill".
- Select "Average Time", and set the average time to "2-4s", "4-8s", "8-16s".
- The user can set the pulse volume. The pulse volume can be set to 0, 1, 2, 3, 4, 5, 6, 7, 8, or 9.
- Select key to display the following interface.



- Select "SpO2 Low Limit", and enter value (Range: 0-99), Adult/Pediatric/Neonate Default: 90.
- Select "SpO2 High Limit", and enter value (Range: 1-100) , Adult/Pediatric/Neonate Default: 100/100/95.
- Select "Alarm Level", you can choose "Mid" or "High".
- Select "PR Low Limit", and enter value (Adult Range: 15-299, Pediatric/Neonate Range: 15-

- 349) , Adult/ Pediatric / Neonate default: 50/75/100
- Select "PR High Limit", and enter value (Adult Range: 16-300, Pediatric/Neonate Range: 16-350), Adult/Pediatric/Neonate default: 120/160/200.
- Alarm ,Select "Default", the alarm parameter is set to the default value
- Select "PI Low Limit", and enter value (Range: 0.00-19.90). Default: 0.00.
- Select "PI High Limit", and enter value (Range: 0.10-20.00). Default: 20.00.

#### 5.5 Measuring influencing factors

During operation, the following factors can affect the accuracy of SpO<sub>2</sub> measurement:

- High-frequency radio wave interference, such as interference generated by the host system or interference from electrosurgery instrument connected to the system.
- Intravenous dye.
- Too frequent movement of the patient.
- External light radiation.
- Sensor is improperly installed or improperly in contact with the patient.
- Sensor temperature.
- The sensor is placed on limbs with blood pressure cuff, arterial duct or lumen tube.
- Concentration of non-functional hemoglobin such as carboxyhemoglobin (COHb) and methemoglobin (MetHb).
- Shock, anemia, hypothermia, and the application of vasoconstrictor drugs may reduce the arterial blood flow to a level that cannot be measured.
- The measurement also depends on the absorption of specific wavelengths of light by oxyhemoglobin and reduced hemoglobin. If there is any other substance that absorbs the same wavelength, the measurement may have false or low SpO<sub>2</sub> values, such as: carbon hemoglobin, methemoglobin, methylene blue, and indigo carmine.
- SpO<sub>2</sub> probe described in Annex is recommended.
- Operating environment limit: Operating temperature range: 5 ~ 40 °C, Humidity range: 30%~85% (non-condensing) Atmospheric pressure: 700hPa ~ 1060hPa.

#### 5.6 Technical description

- Accessories have passed the biocompatibility test and meet the requirements of ISO 10993-1.
- Fluke's index 2XL Oxygen Analyzer can be used to check the function of the monitor and can be used to assess the accuracy of the pulse rate but cannot be used to assess the accuracy of blood oxygen.
- Functional testers cannot be used to evaluate the accuracy of pulse oximetry probe and pulse oximetry.
- Measure the maximum temperature between the oxygen probe and the tissue contact surface: Measured as described in Annex BB of ISO 80601-2-61, the temperature is less than 41° C.

#### 6 NIBP (Option)

#### 6.1 Overview

The monitor uses oscillometric method to measure noninvasive blood pressure (NIBP).

The oscillometric method for measuring blood pressure is to inflate a cuff with a certain amount of pressure until the arterial blood flow has been completely blocked. As applied pressure decreases, the arterial blood flow which was completely occluded gradually opened, and completely opened. Then, the pulsation of the arterial vascular wall will generate a shock wave in the cuff. SBP, MAP, and DBP are obtained by measuring and analyzing cuff pressure oscillations when deflating.

- Produce first most clear signal reflect SBP
- Oscillation amplitude reaches the peak reflect MAP
- When the cuff pressure is suddenly lowered reflect DBP

Measuring mode: manual mode and automatic mode. Each mode shows systolic, mean and diastolic blood pressure.

- Manual mode
  - Using Manual mode start to measures by hand
- Automatic mode measures

Use manual mode to open automatic mode, then the measure will automatically turn to automatic mode after a certain time. During measurement, any error will stop the current automatic measurement, but not affect next automatic measurement unless the time interval less than 30s. If the time interval less than 30s, should delay the next automatic measurement, keep the interval more than 30s.

The time interval can be choose In Automatic mode as 1, 2, 3, 4, 5, 10, 15, 30, 60, 90, 120, 180, 240, 480 minutes.

The monitor is defibrillation proof, so the monitor operates normally after defibrillation.

#### 6.2 Safety information



#### WARNING

- Do not carry out non-invasive blood pressure measurement on patients with sickle cell disease and skin damage or any expected damage. Do not measure NIBP on traumatic body part. This may cause further injury.
- When pediatric and neonate patients are measured, in order to ensure the cuff pressure
  does not exceed its maximum measurement range of patient types (Adult mode:
  300mmHg and Pediatric mode: 240mmHg, Neonate mode: 150mmHg), you must ensure
  that you have selected the correct patient type (see patient information menu settings).
  Using the wrong type of pattern is likely to endanger the patient to patient safety, as
  higher blood pressure levels for adults does not apply to pediatric and neonate.
- For patients with severe coagulation disorder, determine if the automatic blood pressure measurement is carried out according to the clinical evaluation, since the friction of body and cuff may produce hematoma.
- Do not install a cuff on the limbs with intravenous infusion or duct, because it may lead
  to tissue damage around the duct when the cuff is inflated and makes the infusion slow
  down or be blocked.
- The inflatable tube connecting the blood pressure cuff and the monitor should be smooth without entanglement. The pressure generated by being kinked connection tubing may cause blood flow interference.
- For patients with severe thrombotic disorders, determine whether to carry out automatic blood pressure measurement according to the clinical situations, since the limb bundled with a cuff may produce hematoma.
- Measure blood pressure frequently will affect the distribution of blood flow, May endanger the safety of patients.
- Check the patient's physiological condition before measure blood pressure, in order to
  ensure that long time measure will not damage the circulation of patients
- For mastectomy patients, applying the NIBP cuff on the surgery side arm can cause lymphedema. Measure blood pressure on opposite side arm.
- Pressurization of the cuff can temporarily cause loss of function of simultaneously used monitoring device on the same limb.
- Measurement results may be affected by posture and mental state of the patient.
- If there are doubts on the measurement results, please use other blood pressure measurements and compare, if necessary, contact the Equipment Division.

#### 6.3 Measurement limits

According to the patient's condition, the oscillometric method has some limitations. This measurement is to look for the regular pulse waves generated by arterial pressure. If the patient's condition makes this detection method difficult, the measured value becomes unreliable, and pressure measurement time increases. The user should be aware that the following conditions may interfere with measurement method, making the pressure measurement unreliable or extend the time. In this case, the patient's condition does not allow measurement.

#### Patient movement

If the patient is talking, moving, shaking or cramping, the measurement will be unreliable or even impossible, as these may interfere with the detection of arterial pressure pulse, and extend the pressure measurement time.

#### Arrhythmia

If the patient shows arrhythmia which results in irregular heartbeat, the measurement will be unreliable and even cannot be done, and the pressure measurement time will be extended.

> Use of an artificial heart-lung machine

If a patient is connected to an artificial heart-lung machine, the measurement will be impossible.

Pressure changes

If the arterial pressure pulse is being analyzed to obtain a measured value at a certain time and the blood pressure of the patient changes rapidly, the measurement will be unreliable or impossible.

Severe shock

If the patient is in severe shock or hypothermia, the pressure measurement will not be reliable, because the decrease of blood flow to the periphery would cause decrease in arterial pulsation.

Limit heart rate

If the heart rate is below 40bpm (beats / min) or above 240bpm (beats / min), the blood pressure measurement is impossible.

Obese patients

A thick layer of fat around a limb blocks the arterial oscillation so that it cannot reach the cuff. The accuracy is lower than normal.

Environmental Requirements

Measuring blood pressure should meet the environment range as follow:

ambient humidity 30% ~ 85%, no condensing,

ambient temperature  $5 \sim 40 \,^{\circ}\text{C}$ ,

Atmospheric pressure: 700hPa ~ 1060hPa.

NIBP performance and measurement accuracy will be affected beyond the range.

#### 6.4 Measurement procedure

#### 6.4.1 Prepare the measurement

- 1. Turn on the monitor, and check if it works properly.
- 2. Verify the patient category, and make changes if improper. Depending on the current patient type, the patient type is selected in the patient information interface.
- 3. Connect the blood pressure cuff extension tube to the monitor.
- 4. Select the cuff in accordance with the following method, make sure that the cuff is completely deflated, and then tie it to the upper arm or thigh of the patient.
  - Determine the limb circumference of the patient.
  - Select the appropriate cuff (marked with appropriate limb circumference). Cuff width should be 40% of the limb circumference (50% for neonate) or 2/3 of the upper arm length. The length of the inflated part of the cuff should be sufficient for 50%~80% around the limb.
  - Place the cuff on the upper arm or thigh of the patient, and ensure that the marking "ARTERIA" is located just above the appropriate artery. Make sure that the cuff does not

wrap too tight around the limb, or it may cause distal discoloration or even ischemia.

#### 6.4.2 Patient posture requirements during measurement

- 1. Sit comfortable or lie down relaxedly.
- 2. No crossing legs.
- 3. Back and elbow should be supported.
- 4. The center of NIBP cuff and the right atrium are at in the same level.
- 5. Remind patients, no talking during measurement and try to relax.

#### NOTE

- When have doubt about blood pressure measuring result, re-measure after the
  patient sit-in about 5 minutes. If still have doubt, replace the blood pressure
  measuring equipment and measure again.
- The operator should be in the position where he/she can readily operate the sphygmomanometer.

#### 6.4.3 Start/stop measurement


Use the [NIBP start/stop] key on the monitor panel to start / stop the blood pressure measurement.

#### 6.4.4 Correcting measurement results

The position of limb blood pressure measurement should be in the same horizontal position of the patient's heart. Otherwise, correct the measurement results with the following correction method.

- If the cuff is above the heart level position, increase 0.75mmHg (0.10kPa) per centimeter of gap to the measured results.
- If the cuff is below the heart level position, subtract 0.75mmHg (0.10kPa) per centimeter of gap from the measured results.
- If the patient is obese or clothes are too thick, subtract 5mmHg ~ 10mmHg (0.65kPa ~ 1.3kPa) from the measured results.

#### 6.5 Setting NIBP



- Select "Patient Type", you can choose "Adult" or "Pediatric" or "neonate".
- Select "Interval", you can choose Manual, 1min, 2min, 3min, 4min, 5min, 10min, 15min, 30min, 60min, 90min, 120min, 180min, 240min or 480min.
- Select "Unit", and select the unit "mmHg" or "kPa".
- Select "Initial Pressure", and set the appropriate cuff pressure value. When the patient is adult, the pressure can be select from "140", "160", "180". The default cuff pressure value is "160".
- Select "Initial pressure", and set the appropriate cuff pressure value. When the patient is pediatric, the pressure can be select from "140", "160". The default cuff pressure value is "140".
- Select "Initial pressure", and set the appropriate cuff pressure value. When the patient is neonate, the pressure can be select from "100", "120". The default cuff pressure value is "100".

- Select "Start Measure", starts the blood pressure measurement; select "Stop Measure", stop the blood pressure measurement.
- Select "Continue Measure", continued measure pressure.
- Select "Reset", and restore the inflation pressure of the blood pressure pump to currently configured initial settings. When the blood pressure pump is not working properly, but no warning is given, you can reset the blood pressure pump, and automatically restores the blood pressure pump.
- > Select key to display the following interface.



- Select "SYS Low Limit", and enter value (Adult Range: 30-279, Pediatric Range: 30-229, Neonate Rang: 30-144). Adult/ Pediatric /Neonate default: 90/70/40.
- Select "SYS High Limit", and enter value (Adult Range: 31-280, Pediatric Range: 31-230, Neonate Rang: 30-145). Adult/ Pediatric /Neonate default: 160/120/90.
- Select "DIA Low Limit", and enter value (Adult Range: 10-219, Pediatric Range: 10-164, Neonate Rang: 10-104). Adult/ Pediatric /Neonate default: 50/40/20.
- Select "DIA High Limit", and enter value (Adult Range: 11-220, Pediatric Range: 11-165, Neonate Rang: 11-105). Adult/ Pediatric /Neonate default: 90/70/60..
- Select "Alarm Level", you can choose "Mid" or "High".
- Select "Default", the alarm parameter is set to the default value.



- Select "Calibration", Users can not calibrate NIBP. If calibration is required, please contact your service representative. Cuff pressure sensor should be checked and calibrated at least once every two years by qualified professional service personnel.
- Select "Leak Test", the purpose of leakage test is to detect if the sealing of the air passage is in good condition. If the leakage test passes, the alarm area displays "Leakage test Stopped". If not passed, the alarm area displays "Cuff leak" message. NIBP leakage test shall be at least once every two years or when you think that the reading is not accurate.

#### 6.6 Clean and disinfection method of NIBP cuff

If necessary, NIBP cuff and NIBP extension tube can be cleaned and disinfected together without separated

#### 6.6.1 Cleaning method

- Prepare enzyme cleaning agent, distilled water and 10% solvent, respectively in different spray bottle.
- Sprinkle cleaning agent on NIBP cuff, connector and extension tube, keep 1 minute for the dry stains.
- Use a soft cloth to wipe smooth face. Use soft hair brush to brush visible stain and irregular surface
- 4. Rinsed with copious amounts of distilled water.

#### NOTE

- Please be especially careful to clean the air ball and control valve of whole air system.
   Do not allow any liquid entering into reversing valve and saturated valve.
- Don't use a soft cotton ball and fiber to clean this accessory because they will stick on the cuff and extension tube.

#### 6.6.2 Disinfection method

- Sprinkle bleach solution (Formula: the proportion of water and bleaching powder to 1:10) then keep 5 minutes
- 2. Wipe off excess bleach solution and elute with distilled water again
- 3. Natural dry cuff

#### 7 Review

Use the Trend screen to recall all the historical patient data in a list, including monitoring time (in 1-minute intervals), and SpO<sub>2</sub>, PR, SYS, DIA and MAP values. The most recent measurements display at the top of the list.



Figure 7-1: SpO2 trend (1)

#### To view historical data:

- Press the [Setting] key twice to display the Trend screen.
- Select or key to page up and down to view patient data.
- Press key to exit the Trend screen and display the Main screen or press the (Setting) key to see the trend graph.





Figure 7-2: SpO2 trend (2)

- Select 'SpO2', 'PR' or 'PI' to see.
- Press key to exit the Trend screen and display the Main screen or press the [Setting] key to see the NIBP Trend.



Select or key to page up and down to view patient data.

#### 8 Alarm

Select button, and you can temporarily turn off the alarm sound of currently occurring physiological alarms of the monitor, but the alarm information is still retained. For technical alarms, clear the alarm state, display alarm prompt information, the alarm state icon area displays the

clear the alarm state, display alarm prompt information, the alarm state icon area displays the icon. When a new physiological alarm or technical alarm occurs, the alarm reset is automatically canceled.

Physiological alarm information.

Source	Alarm message	
SpO <sub>2</sub>	SpO <sub>2</sub> Too High	
	SpO <sub>2</sub> Too Low	
	PR Too High	
	PR Too Low	
	Systolic Too High	
	Systolic Too Low	
NIBP	Mean pressure Too High	
	Mean pressure Too Low	
	Diastolic Too High	
	Diastolic Too Low	

Technical alarm information.

Source	Alarm message		
6.0	SpO <sub>2</sub> Communication Stop		
	SpO <sub>2</sub> Communication Error		
	SpO <sub>2</sub> No Sensor		
SpO <sub>2</sub>	SpO <sub>2</sub> Sensor Off		
	SpO <sub>2</sub> Search Timeout		
	SpO <sub>2</sub> Search Pulse		
NIBP	NIBP Communication Stop		
	NIBP self-check error		
	Cuff type error		
	Cuff loose or no cuff		
	Cuff leak		
	Air pressure error		

NIBP signal weak
NIBP over range
NIBP signal unstable
NIBP over pressure
NIBP signal saturated
NIBP system error
Measurement timeout

#### 9 Battery

#### 9.1 Overview

The monitor has a built-in rechargeable battery to ensure that the monitor can also be used normally in case of patient transfer or power failure. When the monitor is connected to an DC power source, it will charge the battery no matter whether the monitor is turned on or not. In the case of power failure, the system will automatically use the battery to power the monitor to avoid interrupting the monitor

The battery icon on the screen indicates the battery status:

	Battery is working properly and is fully charged.
	Battery is working properly and the green part indicates the battery power.
	Battery power is low, and requires charging immediately, or else the pulse oximeter will turn off automatically.
	Battery is not installed.
4	Battery is properly installed and being charged.

The battery power can only maintain for some time. Low battery voltage will trigger a high level technical alarm "Battery Low"; in this case, connect the monitor to DC power and charge the battery.

#### 9.2 Battery usage guide

Battery life depends on the frequency and time of use. If the battery maintenance and storage are proper, the lithium battery life is three years. If you do not use the battery properly, its life may be shortened. It is recommended to replace the lithium battery once every three years.

In order to ensure the maximum capacity of the battery, please note the following usage guide:

- Do not drop the battery.
- Check the battery performance once every two years. Before servicing the monitor or you suspect that the battery is the fault source, also check the battery performance.



#### WARNING

- Keep the battery out of the reach of children. Use only the designated battery.
- If the battery is damaged or leaks, replace it immediately. Do not use a defective battery for the monitor.
- Do not open the battery compartment. Only the qualified service personnel authorized by the manufacturer can open the battery compartment and replace or check the battery, and when it needs replacement, the replacement should be same model with established battery.
- Do not disassemble the battery, throw it in fire, or short-circuit it. Battery fire, explosion and leakage may lead to personal injury. Do not touch the leaking battery with bare hand directly.

#### 9.3 Checking battery performance

Please refer to the following steps to check the battery performance:

- Disconnect the monitor from the patient and stop all monitoring or measurement.
- Connect DC power to the monitor, and charge battery for more than 4 hours uninterruptedly.
- $\triangleright$ Disconnect the DC power and power the monitor with battery until the monitor is turned off.
- Battery duration reflects the battery performance.

If the battery operating time is significantly shorter than the time stated in specifications, please contact our service personnel for replacing the battery.



#### WARNING

Do not open the battery compartment. Only the qualified service personnel authorized by the manufacturer can open the battery compartment and replace or check the battery, and when it needs replacement, the replacement should be same model with established battery.

#### 9.4 Battery recycling

If the battery has visible damage or cannot store power, it should be replaced and recycled properly. Follow the appropriate regulations to dispose of used batteries.



#### 🔼 WARNING

Do not disassemble the battery, throw it in fire, or short-circuit it. Battery fire, explosion and leakage may lead to personal injury. Do not touch the leaking battery with naked hand directly.

#### 10 Caring and cleaning

#### 10.1 Overview

In the using process, please make sure that there is no dust on or near your device. To prevent damage, please use the diluted detergents and disinfectants specified in this Manual, and use the lowest possible concentration. For the damage or accident caused by using other materials or methods, our company does not assume any responsibility.

#### 10.2 Cleaning

The device should be cleaned regularly. In the heavily polluted environment, increase the frequency of cleaning. Before cleaning, please consult the hospital about device cleaning requirements. Below are available cleaning agents:

- $\triangleright$ Diluted ammonia
- $\triangleright$ Diluted sodium hypochlorite (washing bleach)
- Diluted formaldehyde
- $\triangleright$ Hydrogen peroxide (3%)
- Ethanol (70%)
- Isopropanol (70%)

#### Before cleaning:

- Turn off the monitor and disconnect the power.
- Use a soft cotton ball to adsorb appropriate amount of cleaning agent and wipe the display screen.
- $\triangleright$ Use a soft lint-free cloth to adsorb appropriate amount of cleaning agent and wipe the surface of the device.
- If necessary, use a clean, dry, lint-free cloth to remove any excess detergent.
- Dry the device naturally in a ventilated cool environment.



#### WARNING

- Before cleaning the monitor or sensor, turn off the power and disconnect the DC power.
- The monitor should be kept clean. It is recommended to regularly clean the enclosure surface and the display screen. Cleaning the enclosure with non-etching cleaner such as soap and water.



#### CAUTION

- To avoid damaging the monitor:
  - Do not use strong solvents such as acetone.
  - Most cleaners must be diluted before use. Diluting should be according to the manufacturer's instructions.
  - > Do not use abrasive materials (such as steel wool).
  - Do not allow any liquid entering into the enclosure, and never immerse any part of the device into liquid.
  - Do not leave any cleaning solution on the surface of any part of the device.

#### NOTE

- Wipe the monitor and sensor surface with medical alcohol, dry it naturally or with clean, dry. lint-free cloth.
- Bistos is not liable for effectiveness of using these chemicals for infectious disease control. Please consult the infectious disease control officers or experts of the hospital for advice.

#### 10.3 Disinfection

In order to avoid damage to the product, we recommend that the product is disinfected only when it is deemed necessary by the hospital maintenance procedures. We also recommend that the instrument to be disinfected must first be cleaned.



#### CAUTION

 To prevent damage to the monitor, do not disinfect the monitor with gas (EtO) or formaldehyde.

#### 11 Maintenance



#### WARNING

 If the hospitals or institutions using this instrument can't implement a satisfactory maintenance schedule, it will result in device failure and may endanger human health.

#### 11.1 Checking

Check the following basic items before using the monitor:

- Check for any mechanical damage.
- Check all exposed wires, insertions and accessories.
- Check all instrument functions that may be used for patient monitoring and ensure that the instrument is in good working condition.

If the instrument function has any sign of damage, do not use this monitor for any patient monitoring. Please contact the hospital's professional maintenance personnel or our customer service personnel. Every 6-12 months or after each repair, a comprehensive examination must be performed by trained and qualified technical service personnel, including functional safety checks; the specific inspection items are as follows:

- Environment and power meet the requirements.
- Device and accessories have no mechanical damage.
- The power supply has no wear, and the insulation is good.
- Specified accessories are used.
- Alarm system is functioning correctly.
- Battery performance meets the requirements.
- Monitoring functions are in good working condition.
- Ground impedance and leakage current meet the requirements.

If the instrument function has any sign of damage, do not use this monitor for any patient monitoring. Please contact the hospital's professional maintenance personnel or our customer service personnel. All checks that require disassembling the instrument must be performed by qualified service personnel. Safety and maintenance checks may also be carried out by the Company's personnel.

#### 11.2 Maintenance plan

The following tasks can only be done by qualified service personnel of Bistos. When the following maintenance is needed, please contact your service representative. Before testing or maintenance, clean and disinfect the device.

Inspection / Maintenance Item	Frequency
Check the safety according to IEC	At least once every two years, after replacing the
60601-1	power supply or the monitor falls down.
Check all monitoring or measuring	At least once every two years, or when you suspect
functions not listed	that the measured value is not accurate.
NIBP leakage test	At least once every two years, or follow hospital
MBP leakage test	regulations
NIDD colibration	At least once every two years, or follow hospital
NIBP calibration	regulations

#### 12 Accessories



#### WARNING

- Use the accessories specified in this manual. Using other accessories may damage the monitor, or cannot reach the safety and performance claimed in this manual.
- The operating and storage environment of the monitor should meet the requirements of the accessories. Please refer to the manual of the accessories for these requirements.
- Disposable accessories can only be used once, because repeated use can cause performance degradation.
- If the packaging or accessories have any sign of damage, do not use such accessories.
- For  $SpO_2$  Sensor and Blood Pressure Cuff, the normal life time is two years. Please replace in time.

#### Standard accessories are as follows:

No.	Description	QTY	Type-number
4 4 4 5 2 2 2 6 2	Adult Finger Clin CnO Conser	1	Manufacturer:
1	1 Adult Finger Clip SpO <sub>2</sub> Sensor		Unimed Medical Supplies,Inc
2	SpO2 extension cable		U403-01
3 Power Adapter			Manufacturer:
	1	DONGGUAN SHILONG GUHUA ELECTRONIC	
		CO., LTD	
			UE36LCP1-150240SPA

Optional accessories list is as follows:

No.	Description	QTY	Type-number
	Adult Non-Invasive blood pressure cuff	1	Manufacturer:
1			Shenzhen Med-link Electronics Tech Co.,Ltd
blood pressure curr		Y000A1	
2 NIBP extension tube			Manufacturer:
	NIDD and and an Archae		XIAMEN CONJOIN ELECTRONICS TECHNOLOGY CO.,
	1	LTD	
			CJP37-C12B1

#### 13 Specifications

#### 13.1 Safety specifications

#### 13.1.1 Product category

In accordance with classification specified in the European Medical Device Directive 93/42/EEC, this monitor is Class IIb device. The monitor is classified as follows in accordance with IEC 60601-1:

Category Name	Specification
Type of electric shock protection	Class II and internally powered equipment When you question the integrity of the external protective earthing or protective ground conductor parameter of the equipment, the device must be powered by the internal power supply (battery).
Electric shock protection grade	Type CF applied part (defibrillation proof)
Explosion protection grade	Common equipment, no explosion protection
Liquid inlet protection grade	IPX1
Operating mode	Continuous mode
Movement	Portable equipment

#### 13.1.2 Power

Power	
Adapter	Input: AC 100 ~ 240V (50/60 Hz)
	Output: DC 15V / 2.4A
Rechargeable Battery	11.1V Li-ion battery 4400 mA
	Operating Time(When it fully charged): 5 hours
	Charging Time(Fully): 4 hours

13.2 Hardware specifications

Physical Characteristics	
Dimensions	Main Unit: 254(W) X 90(H) X 185(D)
Weight	< 1.5 Kg

Display	
Type	Color TFT touch screen LCD
Size	4.3", 480 x 272 pixels

LED	
Alarm Indicator	Yellow & Red
Adapter power indicator	1 green

Battery status indicator	1 green
--------------------------	---------

Audio	
	Alarm sound (45 ~ 85 dB), key pressing sound
Speaker	PR sound
	Alarm sound meet the IEC 60601-1-8 standard requirements

Alarm signal	
Alarm delay	Off, 1s, 2s, 3s, 4s, 5s, 6s, 7s, 8s, depending on the setup

Data storage	
Trend	168 hours. Resolution: 1 min

Environment		
	Operating	Transport and storage
Temperature	5 ℃~40 ℃ (41°F~104°F)	-20℃~60℃(-4°F~140°F)
Relative Humidity	30~85% RH, Non-condensing	0~95 % RH, Non-condensing
Atmospheric	70kPa~106kPa	70kPa~106kPa
pressure		

#### **13.3 Functional specifications**

#### 13.3.1 SpO<sub>2</sub>

SpO <sub>2</sub>			
Standards compliant	ISO 80601-2-61:2011		
Display range	0% ~ 100%		
SpO <sub>2</sub> display resolution	1%		
SaO₂ accuracy	±2% (70%~100%) (adult/pediatric mode); ±3% (70%~100%) (neonate mode); not define when lower than 70%;		
CaO2 alama limit ranga	Upper alarm limit	1%~100%	
SpO2 alarm limit range	Lower alarm limit	0%~99%	
SpO <sub>2</sub> alerting signal generates a delay	No delay		
SpO <sub>2</sub> value refresh period	1s/time		
	Low sensitivity	6~8s	
Average period	Intermediate sensitivity	4 ~ 6s	
	Advanced sensitivity	2~4s	
	Low sensitivity	<8s	
Alarm condition delay period	Intermediate sensitivity	<6s	
	Advanced sensitivity	<2s	
Alarm sign generates delay period	0s		

PR	
Measuring range	25~250bpm
Resolution	1 bpm
Accuracy	±2% or ±2bpm,whichever is greater

#### 13.3.2 NIBP

NIBP	
MIDI	

	-30:2009			
Automatic oscillometric method				
Manual, automatic				
100, 000 times				
1/2/3/4/5/10/15/30/60/90/120/180/240/480min				
20~40s				
		Adult	Pediatric	Neonate
Systolic bloo	d pressure	40-270	40-200	40-130
Mean blood	pressure	20-230	20-175	20-100
Diastolic blo	od pressure	10-210	10-162	10-90
Maximum average error: ±5mmHg				
Maximum standard deviation: 8mmHg				
1mmHg				
Default Pressure setting range				
Adult	160mmHg	140mmHg,	160mmHg, 18	0mmHg
Pediatric	140mmHg	140mmHg, 160mmHg,		
Neonate	100mmHg	100mmHg,	120mmHg,	
Adult: 300mmHg				
Pediatric: 240mmHg				
Neonate: 150mmHg				
Adult: 320~330mmHg Pediatric: 265~275mmHg				
			Neonate: 160~165mmHg	
±3mmHg				
	Manual, aut 100, 000 tim 1/2/3/4/5/1 20~40s Systolic bloo Mean blood Diastolic blo Maximum at 1mmHg Adult Pediatric Neonate Adult: 300m Pediatric: 24 Neonate: 15 Adult: 320~3 Pediatric: 26 Neonate: 16	Manual, automatic  100, 000 times  1/2/3/4/5/10/15/30/60/90  20~40s  Systolic blood pressure Mean blood pressure Diastolic blood pressure Maximum average error: ±5 Maximum standard deviation 1mmHg  Default Adult 160mmHg Pediatric 140mmHg Neonate 100mmHg Adult: 300mmHg Pediatric: 240mmHg Neonate: 150mmHg Adult: 320~330mmHg Pediatric: 265~275mmHg Neonate: 160~165mmHg	Manual, automatic  100, 000 times  1/2/3/4/5/10/15/30/60/90/120/180/240,  20~40s  Adult  Systolic blood pressure 40-270  Mean blood pressure 10-210  Maximum average error: ±5mmHg  Maximum standard deviation: 8mmHg  1mmHg  Default Pressure se  Adult 160mmHg 140mmHg,  Pediatric 140mmHg 140mmHg,  Neonate 100mmHg 100mmHg,  Adult: 300mmHg  Pediatric: 240mmHg  Neonate: 150mmHg  Adult: 320~330mmHg  Pediatric: 265~275mmHg  Neonate: 160~165mmHg	Manual, automatic  100, 000 times  1/2/3/4/5/10/15/30/60/90/120/180/240/480min  20~40s  Adult Pediatric  Systolic blood pressure 40-270 40-200  Mean blood pressure 10-210 10-162  Maximum average error: ±5mmHg  Maximum standard deviation: 8mmHg  1mmHg  Default Pressure setting range  Adult 160mmHg 140mmHg, 160mmHg, 180  Pediatric 140mmHg 140mmHg, 160mmHg, 180  Reonate 100mmHg 100mmHg, 120mmHg,  Adult: 300mmHg  Pediatric: 240mmHg  Neonate: 150mmHg  Adult: 320~330mmHg  Pediatric: 265~275mmHg  Neonate: 160~165mmHg

Electrical characteristics	
Supply voltage	10V~14V DC
Maximum power consumption	3.6w
Quiescent current	50mA
Maximum current during	180mA
measurement	
Maximum current during inflation	300mA

#### 14 Common faults and maintenance

The following table shows the common faults on the operation, and the solution.

Faults	Solution	
Not power on	Check the battery. If the battery is low, please contact to Bistos.	
Blank Screen	Check the screen and screen line.	
The system time is not correct	Set up error, can be reset through the system User     Maintenance menu.	
	The button battery on main control board is run out, please contact to Bistos.	
No SpO₂ waveform or value	Is the red light on the finger sensor flashing? If not there might be poor contact. Check the extension cable and the connector.	
	2. Is the patient's arm under pressure? Never take blood pressure and SpO2 measurements on the same arm.	
	3. Is the environmental temperature too low? Never expose the	

	patient's arm to cold air since this can affect the readings.		
	4. Has all patient nail polish, especially blue or purple been		
	removed?		
SpO2 value turn on	During long term monitoring, patient movement might result in		
and off during	SpO2 interruption. Keep the patient stabilized. SpO2 interruptions		
monitoring	due to patient hand motion are normal.		
Blood pressure	1. Check whether the pump is broken.		
measurement does	2. Check whether the trachea is broken.		
not start	3. Check whether the blood pressure plate is normal.		
Blood pressure	Check whether the blood pressure cuff is leak.		
started, but couldn't	2. Check whether the NIBP extension tube and machine connect		
measure the value	is well.		
	3. Check whether the deflating valve on blood pressure plate is		
	normal.		
	4. Check whether the pressure sensor is normal.		

If the above doesn't solve the problem, please contact Bistos after-sales department or dealers.

#### 15 Manufacturer's declaration on EMC

BT-720 needs special precautions regarding EMC (Electromagnetic compatibility) and needs to be used according to the EMC information provided in this user manual. Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect the BT-720 and should be kept at least 1 m away from the equipment.

#### NOTE

- Using unqualified accessories, sensors and cables will increase the electromagnetic emission and reduce the electromagnetic immunity of the device.
- Do not put the device close to other devices or stack together. When necessary, observe the device closely to ensure that it runs normally in the environment.
- The device requires special EMC protection, and it is necessary to install and maintain it in the environment that meets the following EMC information.
- Even if other devices comply with CISPR emission requirements, they may also cause interference to this device.
- When the input signal amplitude is smaller than the minimum amplitude specified in the technical specifications, it may result in inaccurate measurements.
- Mobile communication devices or wireless network devices may have an impact on the device.

#### 15.1 Electromagnetic emissions

The BT-720 is intended for use in the electromagnetic environment specified below.		
The customer or the user of the BT-720 should assure that it is used in such a		720 should assure that it is used in such an environment.
<b>Emissions test</b>	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The BT-720 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The BT-720 is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following

		warning is heeded:
Harmonic emissions IEC 61000-3-2	Class A	Warning: This BT-720 is intended for use by healthcare professionals only. This equipment/ system may cause radio interference or may disrupt the operation of nearby
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the BT-720 or shielding the location.

## 15.2 Recommended separation distances between portable and mobile RF communications equipment and BT-720

The BT-720 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the BT-720 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BT-720 as recommended below, according to the maximum output power of the communications equipment.

Rated	Separation distance according to frequency of transmitter [m]		
maximum output power of transmitter [W]	150 kHz to 80 MHz $d = 3.5\sqrt{p}$	80 MHz to 800 MHz $d=3.5\sqrt{p}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3}\right] \sqrt{p}$
0.01	0.35	0.35	0.23
0.1	1.11	1.11	0.74
1	3.5	3.5	2.34
10	11.07	11.07	7.38
100	35	35	23.24

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1) At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

#### 15.3 Electromagnetic immunity

The BT-720 is intended for use in the electromagnetic environment specified below. The customer or the user of the BT-720 should assure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment -guidance
Electrostatic discharge (ESD)	±8 kV Contact	±8 kV Contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with
IEC 61000-4- 2:2009	±12 KV all	±12 KA 9IL	synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst	±2 kV for power supply lines ±1 kV for input/output	±2 kV for power supply lines ±1 kV for input/output lines (>3m)	Mains power quality should be that of a typical commercial or hospital environment.
4:2004	lines (>3m)		·
Surge IEC 61000-4-	±1 kV differential mode ±2 kV common	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or
5:2006	mode		hospital environment.
Voltage dips, short interruptions and voltage variations	< 5 % <i>U</i> τ (> 95 % dip in <i>U</i> τ) for 0.5 cycles	< 5 % <i>U</i> τ (> 95 % dip in <i>U</i> τ) for 0.5 cycle	Mains power quality should be that of a typical commercial or
on power supply input lines	40 % <i>U</i> τ (60 % dip in <i>U</i> τ ) for 5 cycles	40 % Uτ (60 % dip in $U$ τ ) for 5 cycles	hospital environment. If the user of the BT-720 image intensifier
IEC 61000-4- 11:2004	70 % <i>U</i> τ (30 % dip in <i>U</i> τ) for 25	70 % <i>U</i> τ (30 % dip in <i>U</i> τ) for 25 cycles	requires continued operation during power mains interruptions, it is
	cycles <5 % <i>U</i> τ (> 95 %	<5 % <i>U</i> τ (> 95 % dip in <i>U</i> τ ) for 5 s	recommended that the BT-720 be powered from an uninterruptible
	dip in $U_T$ ) for 5 s		power supply.
Power frequency (50 Hz and 60 Hz) magnetic field	3 A/m	3 A/m	Power frequency magnetic fields should be at levels
IEC 61000-4- 8:2010			characteristic of a typical location in a typical commercial or hospital environment.

NOTE  $U_T$  is the a.c. mains voltage prior to application of the test level.

The BT-720 is intended for use in the electromagnetic environment specified below.

The customer or the user of the BT-720 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level
Conducted RF IEC 61000-4-6:2009	3 Vrms 150 kHz to 80 MHz	3 Vrms
Radiated RF	3 V/m	3 V/m
IEC 61000-4-3	80 MHz to 2.5 GHz	

Electromagnetic environment - guidance

Portable mobile RF communications equipment should be used no closer to any part of the BT-720, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

Recommended separation distance

$$d = 1.2\sqrt{p} (d = 3.5\sqrt{p})$$

$$d - 1.2\sqrt{p}$$
 (Resp:  $d - 3.5\sqrt{p}$ ) 80 to 800MHz

$$d - 1.2\sqrt{p}$$
 800M to 2.5GHz

where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).

Field strengths from fixed RF transmitters, as deter-mined by an electromagnetic site survey <sup>a</sup>, should be less than the compliance level in each frequency range. <sup>b</sup>

Interference may occur in the vicinity of equipment marked with the following symbol:



NOTE 1) At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>&</sup>lt;sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BT-720 is used exceeds the applicable RF compliance level above, the BT-720 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the BT-720.

<sup>&</sup>lt;sup>b</sup> Over the frequency range 150 kHz to 80MHz, field strengths should be less than 3 V/m.

#### **Product Warranty**

Product Name	Patient Monitor
Model Name	BT-720
Serial No.	
Warranty Period	2 Years
Date of Purchase	
Customer	Hospital: Address: Name: Telephone:
Sales Agency	
Manufacture	Bistos Co., Ltd.

<sup>\*</sup> Thank you for purchasing BT-720.

#### Service Telephone and Fax. Numbers

Telephone: +82 31 750 0340 Fax: +82 31 750 0344

Bistos Co., Ltd. 7<sup>th</sup> FL., A Bldg., Woolim Lions Valley 5-cha, 302, Galmachi-ro, Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea

> www.bistos.co.kr bistos@bistos.co.kr

EC REP Obelis s.a

Bd. Général Wahis 53 1030 Brussels / BELGIUM Telephone: +32 2. 732.59.54 Fax: +49 2. 732.60.03





<sup>\*</sup> This product is manufactured and passed through strict quality control and inspection.

<sup>\*</sup> Compensation standard concerning repair, replacement, refund of the product complies with "Framework Act on Consumers" noticed by Fair Trade Commission of Republic of Korea.



# BT-720 VITAL SIGN MONITOR



## 4.3" Vital Sign Monitor

SpO<sub>2</sub>, Pulse, NIBP, Masimo SpO<sub>2</sub> Touch screen / Portable design / PC view software







SpO2





SpO2+NIBP

# \*\* Technical Specification

Model			BT-720
Category			Vital Sign Monitor
	Dispaly		4.3" Touch LCD
	Measu remert Range		0-100 %
	Accuracy (70-100%)	Adult / Pedlatrlc	±2 %
C <sub>2</sub> O		Neonate	±3 %
SpO <sub>2</sub>	Accuracy (0-69%)		Unspecfied
	Perfusio	n Index	0.05-20 %
	Pulse Rate F	Range (bpm)	25-250
	Met	hod	Automatic Osdlometric
	Operation Mode		Manual / Auto / STAT
	Parameter		Systoic, Diastolic, Mean
	0 1 11 10	Adult	40-270
	Systolic Range (mmHg)	Pediatric	40-200
NIBP**		Neonate	40-130
NIDE	Dlastolic Range (mmHg)	Adult	10-210
		Pediatric	10-160
		Neonate	10-90
		Adult	20-230
	Mean Range (mmHg)	Pediatric	20-175
	(IIIIIIII)	Neonate	20-100
SpO <sub>2</sub> -Masimo*		Masimo SpO <sub>2</sub>	
Type (capacity)		Li-ino (4400 mAh)	
Battery	Run Time		8 hour
	Charging Time		8 hour
	PC Software Interface		RJ45, SD card slot
			2 year








Trend Graphic

Trend Tabular

PC Viewer Software



# BT-770 Patient monitor Operation Manual



Keep this manual for future reference

P/N: 770-ENG-OPM-EUR-R00

## **Proprietary Material**

Information and descriptions contained in this manual are the property of Bistos Co., Ltd. and may not be copied, reproduced, disseminated, or distributed without express written permission from Bistos Co., Ltd.

Information furnished by Bistos Co., Ltd is believed to be accurate and reliable. However, no responsibility is assumed by Bistos for its use, or any infringements of patents or other rights of third parties that may result from its use. No license is granted by implication or otherwise under any patent or patent rights of Bistos Co., Ltd.

The information contained herein is subjects to change without notice.

#### Prepared by:

Bistos Co., Ltd. 7<sup>th</sup> FL., A Bldg., Woolim Lions Valley 5-cha, 302, Galmachi-ro, Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea

Telephone: +82 31 750 0340 Fax: +82 31 750 0344

Revision R00 September, 2018

Printed in Korea Copyright © Bistos Corporation 2018. All rights reserved.

## **Contents**

	Safety information ······	
	0.1 General precautions, warnings and cautions	9
	0.2 Shock hazards ······	···· 13
	0.3 Battery warnings	
	0.4 General precautions on environment	··· 15
	System basics ······	
	1.1 Intended use ·····	
	1.2 Operating principle ·····	
	1.3 System configurations ·····	
	1.4 Product outlook ·····	
	1.5 Description of monitor ·····	
	1.6 Understanding the display ·····	
	1.7 Smart Hotkeys ·····	
	1.8 Essential performance ······	··· 28
	Preparing for operation ······	
	2.1 Installation ·····	
	2.2 Connecting to power ·····	···· 31
	Basic operations ······	
	3.1 Turn on	
	3.2 Turn off	
	3.3 Basic operations ·····	
	3.4 Operation mode ·····	
	3.5 Measurement setup ······	
	3.6 Freezing waves ·····	
	3.7 Other common setup ······	37
4.	Patient information management ······	41
	4.1 Patient setup menu ·····	
	4.2 Admitting a patient	
	4.3 Patient information	
	4.4 Discharging a patient ·····	
	4.5 Clear alarms ·····	
	4.6 Clear trend ·····	
	4.7 Clear NIBP trend ·····	45
	Display format ·····	
	5.1 Selecting user interface ······	
	5.2 Display description	··· 47
	Alarm	
	6.1 Alarm types ·····	
	6.2 Alarm condition priorities	50
	6.3 Alarm mode	

	6.4 Alarm states ·····	
	6.5 Alarm setup	56
	6.6 Latch alarm ·····	58
	6.7 Manual event ·····	59
	6.8 Alarm record ·····	59
7.	ECG	
	7.1 Overview	
	7.2 Safety information	
	7.3 Monitoring steps ·····	
	7.4 ECG display ·····	
	7.5 ECG setup ·····	
	7.6 Alarm setup ·····	68
8.	RESP	
	8.1 Overview ·····	
	8.2 Safety information	
	8.3 Placing electrodes for respiration monitoring	70
	8.4 Respiration display ·····	
	8.5 Respiration setup ·····	
	8.6 Alarm setup ·····	······ 73
9.	PR	
	9.1 Overview	
	9.2 Display	
	9.3 Setting PR sound ·····	
	9.4 Alarm setup ·····	74
10	D. SpO <sub>2</sub> ·····	
	10.1 Overview	75
	10.2 Safety information	
	10.3 Monitoring steps ·····	77
	10.4 Display ·····	78
	10.5 Setting SpO <sub>2</sub> ·····	
	10.6 Measuring influencing factors ······	80
	10.7 Alarm setup ·····	80
	10.8 Technical description	81
11	1. NIBP ·····	82
	11.1 Overview ·····	82
	11.2 Safety information	83
	11.3 Measurement limits ·····	
	11.4 Measurement procedure ·····	85
	11.5 NIBP display ·····	87
	11.6 Setting inflation pressure ·····	
	11.7 NIBP reset ·····	87
	11.8 Clean and disinfection method of NIBP cuff	88
	11.9 Alarm setup ·····	

12	. TEMP	89
	12.1 Overview ·····	
	12.2 Safety information	
	12.3 Measurement steps ······	
	12.4 Measurement requirements ·····	
	12.5 Temperature display ·····	
	12.6 Setting temperature unit	
	12.7 Alarm setup ·····	· 91
13	. Review	
	13.1 Reviewing trend chart ······	
	13.2 Reviewing trend table ······	
	13.3 NIBP measurement review	· 94
14	. Battery ·····	
	14.1 Overview	
	14.2 Battery usage guide ·····	
	14.3 Checking battery performance	
	14.4 Battery recycling ·····	97
15	. Caring and cleaning ·····	
	15.1 Overview	. 98
	15.2 Cleaning ·····	· 98
	15.3 Disinfection ·····	. 99
16	. Maintenance ·····	
	16.1 Checking ·····	
	16.2 Viewing software version information	
	16.3 Maintenance plan ·····	· 101
	16.4 ECG calibration	· 102
17	. Accessories ······	103
18		104
	18.1 Safety specification ·····	
	18.2 Hardware specifications ······	
	18.3 Functional specification	· 106
19	. Alarm information ······	
	19.1 Physiological alarm ·····	· 111
	19.2 Technical alarm ·····	· 113
20	. Default parameter configuration	115
21	. Common faults and maintenance ·······	117
22	. Manufacturer's declaration on EMC	
	22.2 Recommended separation distances between portable and mobile RF communic equipment and BT-770	· 120
	22.3 Electromagnetic immunity ·····	· 121

Product Warranty ······ 123
Figure 1-1: Front view21
Figure 1-2: Side view
Figure 1-3: Rear view
Figure 1-4: Front view23
Figure 1-5: Side view
Figure 1-6: Rear view25
Figure 1-7: Standard display 27
Figure 3-1: "Settings" menu
Figure 4-1: "Patient" menu
Figure 4-2: "Quick Admit" menu
Figure 4-3: "Patient Info" menu
Figure 5-1: Standard display47
Figure 5-2: Big ECG format
Figure 5-3: Big font format
Figure 5-4: ECG 7-Lead full screen format ······· 49
Figure 7-1: 3-Lead placement method ······ 63
Figure 7-2: 5-Lead placement method ······ 64
Figure 7-3: ECG wave in standard display format
Figure 7-4: ECG parameter in standard display format
Figure 7-5: "ECG Setup" menu
Figure 8-1: 5-lead respiration electrode placement 70
Figure 8-2: Respiration wave
Figure 8-3: Respiration parameter display71
Figure 8-4: "RESP setup" menu
Figure 9-1: PR parameter display 74
Figure 10-1: SpO <sub>2</sub> parameter display ······ 78
Figure 10-2: SpO <sub>2</sub> wave
Figure 10-3: "SpO <sub>2</sub> Setup" menu
Figure 11-1: NIBP parameter display87
Figure 12-1: TEMP parameter display90
Figure 13-1: Trend chart92
Figure 13-2: "Trend" table93
Figure 13-3: NIBP measurement review 94

# **O** Safety information

Before using BT-770 Patient monitor, read this entire manual and be fully understood the following safety information to prevent injury of patient and user.

#### **Symbols Used**

The following symbols identify all instructions that are important to safety. Failure to follow these instructions can lead to injury or damage to the patient monitor. When used in conjunction with the following words, the symbols indicate:

<b>MARNING</b>	Can lead to serious injury or death.
$\triangle$ CAUTION	Can lead to minor injury or product/property damage

The following symbols are placed on product, label, packaging and this manual in order to stand for the information about:



Used to identify safety information.

Be well-known this information thoroughly before using BT-770.



Used to identify safety information.

Be well-known this information thoroughly before using BT-770

Indicates the protection level against the ingress of liquid.

IPX1 is protection against some water drops falling vertically.

It correspond the device, patient monitor and accessory, temperature

sensor

Indicates the protection level against the ingress of liquid.

IPX2 is protection from some water drops when the device is tilted up to and including 15°.

It correspond the accessories for SpO2 and ECG.



Refer to operation manual. Read manual before placing the device.


\_\_\_\_ Indicates DC power supply.

Indicates the device is in the battery operation mode.

→ ☐ Indicates nurse call interface.

吕 Indicates network interface.

2018.09



Bistos Co., Ltd.

device can be safely exposed.

P/N: 770-ENG-OPM-EUR-R00

Indicates the range of atmospheric pressure to which the medical



Indicates the device contains natural rubber latex. (Accessories)



Indicates the packing material is recyclable.



Indicates to not dispose the device together with unsorted municipal waste(for EU only). The solid bar symbol indicates that mains adapter is put on the market after 13 August 2005.

## 0.1 General precautions, warnings and cautions

- Examine the patient monitor and any accessories periodically to ensure that the cables, adapter cords and instruments do not have visible evidence of damage that may affect patient safety or performance. The recommended inspection interval is once per week or less. Do not use the patient monitor if there is any visible sign of damage.
- Only the DC power adapter supplied with the BT-770 is approved for use with the device.
- Do not attempt to service the BT-770 patient monitor. Only qualified service personnel by Bistos Co. Ltd. should attempt any needed internal servicing.
- Perform periodic safety testing to insure proper patient safety. This should include leakage current measurement and insulation testing. The recommended testing interval is once per year.
- If the hospital or healthcare institutions using this device fail to implement a satisfactory maintenance schedule, it will result in device failure and may endanger the patient's safety.
- Use the patient monitor under the conditions specified in this operation manual. Beyond the conditions, the patient monitor may not function properly and the measurement results may not accurate and may result in device failure or endangering the patient's safety.
- Do not operate the BT-770 patient monitor if it fails to pass the power on self-test procedure.
- During the operation, do not disconnect any cable.
- The BT-770 patient monitor is intended to be used by clinical professionals or trained doctors, nurses or laboratory assistant.
- Do not service and maintain or clean the device including accessories while in use with a patient.
- Using the device to one patient at a time.

#### WARNING

- Thoroughly read and understand the manual prior to use of the BT-770. Failure to do so could result in personal injury or equipment damage.
- The device is intended for clinical patient monitoring, and only trained and qualified doctors and nurses should use the device.
- The alarm volume, upper and lower alarm limits should be set according to the actual situation of the using environment. Do not just rely on audio alarm

system while monitoring the patient, because too low alarm volume or muted alarm may result in notice failure of alarm situation and endanger the patient's safety. Please pay close attention to the actual clinical status of the patient.

- Use only the power adapter supplied with monitor.
- Position the monitor where it is easy to de-energize the monitor when needed.
- Do not open the enclosure to avoid an electric shock. Any repair and upgrade
  of monitor should be done by service personnel trained and authorized by
  Bistos. Co., Ltd.
- When handling packaging materials, abide by local laws and regulations or hospital waste disposal regulations. Keep the packaging materials away from children.
- Do not use in the presence of flammable anesthetics to prevent explosion or fire.
- Install the power lines and cables of accessories carefully to avoid patient entanglement or suffocation, cables tangled or electrical interference.
- When the monitor is used together with electrosurgical devices, the user (a doctor or a nurse) should ensure the safety of the patient and instrument.
- The physiological wave, physiological parameters and alarm information displayed on the monitor are only for the doctor's reference and should not be directly used as the basis for clinical treatment.
- This is not a therapeutic device.
- For patients with pacemakers, the cardio tachometer may count the
  pacemaker pulse in case of a cardiac arrest or arrhythmias. Never rely solely on
  the cardio tachometer alarm. Closely monitor the patients with pacemaker. For
  the inhibition of the device on pacemaker, refers to this manual.
- Use of accessories other than those listed and approved for use with this
  product may result in increased emissions or decreased immunity.
- Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this manual. In addition, portable and mobile RF communications equipment can affect medical electrical equipment.
- The equipment shall not be used adjacent to other devices unless verification
  of normal operation in the configuration in which it is to be used can be
  achieved.
- Keep matches, and all other sources of ignition, out of the room in which the
  patient monitor is located. Textiles, oils, and other combustibles are easily
  ignited and burn with great intensity in air enriched with oxygen. Personal

injury or equipment damage could occur.

- A fire and explosion hazard exists when performing cleaning or maintenance procedures in an oxygen-enriched environment.
- The patient monitor has been validated with the accessories and options listed in this manual and found to comply with all relevant safety and performance requirements applicable to the device. It is therefore the responsibility of the person or organization who makes an unauthorized modification, or incorporates an unapproved attachment to the device.
- An operator may only perform maintenance procedures specifically described in this manual.
- Do not remove the covers of a BT-770 yourself to avoid damage to the equipment and unexpected electrical shock. Only qualified Bistos service engineer must repair or replace components.

### $\triangle$ CAUTION

- Please install or carry the instrument properly to prevent damage due to falling, collision, strong vibration or other mechanical force.
- Avoid instrument splashed by water.
- Avoid high temperatures, the instrument should be used within a temperature range of 5  $^{\circ}$  C  $^{\sim}$  40  $^{\circ}$  C  $_{\circ}$
- Avoid using instrument in the environment such as pressure is too high, poor ventilation, dusty, or contain salt, sulfur gas and chemical.
- Before using the monitor, check the monitor and accessories if there is damage that may affect patient safety. If there is obvious damage or aging, replace the parts before use. The replacement should be made with same parts of original parts.
- Before powering on the device, make sure that the power used by the device complies with the supply voltage and frequency requirements on the equipment label or in the Operator's Manual.
- Equipment should be tested at least once a year, the test should be done and recorded by trained, have security testing knowledge and experienced personnel. If there are any problems in the tests, they must be repaired.
- When the instrument and accessories are about to exceed the useful life (expected service life: 5 years), it must be treated in accordance with relevant local laws and regulations or the hospital's rules and regulations.
- Do not connect to other equipment or network which not specified in the

instruction for use, in risk of external high voltage.

- Do not connect any equipment or accessories that are not approved by the
  manufacturer or according to IEC 60601-1 to the monitor. The operation or use
  of non-approved equipment or accessories with the monitor is not tested or
  supported, and monitor operation and safety are not guaranteed in such a
  case.
- Any non-medical equipment (such as the external printer) is not allowed to be used within the patient vicinity (1.5m/6ft.).
- Parts and accessories used must meet the requirements of the applicable safety standards, and/or the system configuration must meet the requirements of the medical electrical systems standard.
- Ensure that the conductive parts of electrodes and associated connectors, including neutral electrodes, do not come in contact with earth or any other conducting objects.
- Protection of ME EQUIPMENT against effects of discharge of a cardiac defibrillator depends on use of proper cables.

## **0.2 Shock hazards**

#### WARNING

- Unplug the monitor from its power source prior to cleaning or maintenance to prevent personal injury or equipment damage.
- Some chemical cleaning agents may be conductive and leave a residue that may permit a build-up of conductive dust or dirt. Do not allow cleaning agents to contact electrical components and do not spray cleaning solutions onto any of these surfaces. Personal injury or equipment damage could occur.
- Do not expose the unit to excessive moisture that would allow for liquid pooling. Personal injury or equipment damage could occur.
- Do not touch the patient and signal input/output parts simultaneously
- Due to the risk of electrical shock hazard, only qualified personnel with appropriate service documentation should service the monitor.

## 0.3 Battery warnings

#### WARNING

- Improper operation may cause the internal lithium ion battery to be hot, ignited or exploded, and it may lead to the decrease of the battery capacity. It is necessary to read the operation manual carefully and pay more attention to warning message.
- Do not open the battery compartment. Only the qualified service personnel authorized by the manufacturer can open the battery compartment and replace the battery, and batteries of same model and specification should be replaced.
- Be careful when connecting the battery with polarity.
- Do not use the battery near fire or environmental temperature exceeds 60  $^{\circ}$ C. Do not heat or splash the battery or throw it into fire or water.
- Do not destroy the battery. Do not pierce battery with a sharp object such as a needle. Do not hit with a hammer, step on or throw or drop the battery. Do not disassemble or modify the battery. The battery can heat, smoke, deformation or burning.
- When leakage or foul smell is found, stop using the battery immediately. If your skin or cloth comes into contact with leaked liquid, cleanse it with clean water at once. If the leaked liquid splashes into your eyes, do not wipe them. Irrigate them with clean water first and go to see a doctor immediately.
- Properly dispose of or recycle the depleted battery according to local regulations.

# **0.4 General precautions on environment**

Do not keep or operate the equipment under the environment listed below.

	Avoid placing in an area exposed to moisture. Do not touch the equipment with wet hand.	Avoid exposure to direct sunlight
	Avoid placing in an area where high variation of temperature exists.  Operating temperature ranges from 5°C ~ 40°C.  Operating humidity ranges from 30% ~ 85 %.	Avoid in the vicinity of electric heater.
	Avoid placing in an area where there is an excessive humidity rise or ventilation problem.	Avoid placing in an area where there is an excessive shock or vibration.
	Avoid placing in an area where chemicals are stored or where there is in danger of gas leakage.	Avoid dust and especially metal material enter into the equipment
00%	Do not disjoint or disassemble the device. Bistos Co., Ltd. does not have liability of it.	Power off when the equipment is not fully ready to operate. Otherwise, the equipment could be damaged.

# 1 System basics

# 1.1 Intended use

The BT-770 Patient Monitors acquire the physiological signals such as ECG, respiratory rate, non-invasive blood pressure (NIBP), blood oxygen saturation ( $SpO_2$ ) and temperature. The signals are converted into digital data and processed, examines the data for alarm conditions and display them. The monitor also provides an operation control panel for users. The patient monitor intend to use in hospital clinical area such as intensive care units, cardiac care units, operating room, emergency department, to provide additional information to the medical and nursing staff about the physiological condition of the adult. The BT-770 patient monitors are intended to be used only under regular supervision of clinical personnel. It is suitable for adult and pediatric. The intended locations of use are hospitals and clinics.

#### 1) Intended patient population

Pediatrics (30 days < and <18 years) and adult (>18 years adults)

#### 2) Intended user profile

- Doctor, physicians or nursing staff who is qualified personnel
- Basic experiences or knowledge on medical field, especially on patient monitoring
- Trained or requested to read IFU before use

#### 3) Environment of use

- Hospital and clinic
- Requirements: Stable power source

#### 4) Scope of application

This monitor is suitable for bedside monitoring of adults and pediatric. This monitor enables ECG, respiration (RESP), pulse rate (PR), blood oxygen saturation (SpO $_2$ ), noninvasive blood pressure (NIBP) and temperature (TEMP) monitoring. It is equipped with a replaceable built-in battery to provide convenience for the patient movement in hospital.

#### 5) Indications and contraindications

#### Blood oxygen saturation (SpO<sub>2</sub>)

#### Indication:

- Monitoring effectives of oxygen therapy
- A reading is needed to facilitate the completion of an early warning score to inform clinical assessment
- Sedation or anesthesia
- Transport of patients who are unwell and require oxygenation assessment
- Haemodynamic instability (e.g. cardiac failure or Myocardial Infarction)
- Respiratory illness e.g. asthma, chronic obstructive pulmonary disease
- Monitoring during administration of respiratory depressant drugs, e.g. opiate epidural or patient-controlled analgesia.
- Assessing oxygen saturation during physical activity e.g. in pulmonary rehabilitation

#### Contraindications

 Pulse oximetry does not give an indication of haemoglobin so if the patient is profoundly anaemic then their oxygen saturation may by normal but they may still be hypoxic

Source: NHS. "Clinical Procedure\_ Procedure for Pulse Oximetry/SPO2". Wirral Community NHS Trust. Sep, 2013

#### Non-invasive blood pressure (NIBP)

#### Indication:

- To determine a patient's blood pressure
- Screen for hypertension
- Following the effect of anti-hypertensive treatments in a patient to optimize their management
- Assessing a person's suitability for a spot or certain occupations
- Estimation of cardiovascular risk
- Determining for the risk of various medical procedure
- Figuring out whether a patient is clinically deteriorating or is at risk.

#### Contraindications

- Oscillometric blood pressure devices may not be accurate in patients with weak or thready pulse
- In patients with heart beats below 50 beats/minutes, even if the rhythm is regular, some of the semi-automatic devices are unable to reduce their deflation rate sufficiently so that too rapid a falling in cuff pressure results in underestimation of systolic blood pressure and overestimation of diastolic blood pressure.
- Do not apply to limb with AV fistula, significant injury or burn, or lymph node removal post mastectomy.

Source: [1] NHS. "Clinical Procedure\_ Procedure for Blood Pressure Monitoring". Wirral Community NHS Trust. Dec, 2013

[2] Clinical Quality& Patient Safety Unit, QAS. Clinical Practice Procedures: Assessment/Non-invasive blood pressure. Queensland Government, 2016. https://www.ambulance.qld.gov.au/clinical.html

#### Electrocardiography (ECG)

#### Indication:

- The electrocardiogram (ECG) has proven to be among the most useful diagnostic test in clinical medicine. It is routinely used in the evaluation of patients to detect myocardial injury, ischemia and the presence of prior infarction, in the assessment of patients with electrolyte abnormalities, drug toxicities and implanted defibrillators and pacemakers.
- In addition to its usefulness in the evaluation of ischemic coronary disease, the ECG, in conjunction with ambulatory ECG monitoring, is of particular use in the diagnosis of disorders of the cardiac rhythm and in the evaluation of syncope. Other common uses of the ECG include the assessment of metabolic disorders and side effects of pharmacotherapy, as the evaluation of primary and secondary cardiomyopathic processes, among others.

#### Contraindications

No absolute contraindications to performing an ECG exist, other than patient refusal.
 Some patients may have allergies or, more commonly, sensitivities to the adhesive used to affix the leads; in these cases, hypoallergenic alternatives are available from various manufacturers.

Source: Tarek, A. "Electrocardiography", < Medscape >, Apr 17, 2017

#### Temperature (TEMP)

#### Indication:

- To obtain the baseline temperature to enable comparisons to be made with future recordings
- To enable close observation in resolving hypothermia/hyperthermia
- To observe and monitor patients for changes indicating an infection
- To monitor the effect of treatment for antimicrobial therapy for infection
- Using before and during a blood transfusion to monitor for signs of a reaction

#### Contraindications

- No known contraindications

# 1.2 Operating principle

Refer to the chapters for every physiological parameter from chapter 7 to chapter 12.

# 1.3 System configurations

Basic configuration of BT-770

- Main body with 12 inch touch screen and built-in lithium-ion battery
- ECG cable and electrode
- Adult SpO2 probe and extension cable
- Non-invasive blood pressure cuff
- Temperature probe
- AC/DC adapter

## Options of BT-770

• External plug-in printer

Picture	Name	Description	Qty
Taxan	ECG cable and lead wire (standard)	Measures ECG	1ea
ECG Electrodes	ECG electrode (standard)	Electrode for ECG measurement	1ea
100 C	Adult SpO2 sensor (standard)	SpO2 sensor for adult	1ea
	SpO₂ extension cord (Standard)	Cord to connect the SpO2 sensor and main body	1ea
and the second s	Adult NIBP cuff (standard)	Measures NIBP for adult	1ea
	NIBP extension tube (standard)	Tube to connect the NIBP cuff and main body	1ea
And the second s	Temperature sensor (Standard)	Measures the body temperature	1ea
	Adapter (Standard)	For power supply	1ea

# 1.4 Product outlook




Figure 1-1: Front view



Figure1-2: Side view



Figure1-3: Rear view

# 1.5 Description of monitor




Figure 1-4: Front view

	Name	Description
1	Alarm indicator	Indicates the priority of physiological alarm and technical alarms in different colors and flashing frequencies.  - High priority: Red, fast flashing (1.4 - 2.8 Hz)  - Medium priority: Yellow, slow flashing (0.4 - 0.8 Hz)  - Low priority: Yellow, constant on
2	Display area	Display the waveform and measured value
3	(Power)	<ul> <li>Power On: Press down the key more than 2 seconds.</li> <li>Power Off: Press down the keys more than 2 seconds and the system will display the alarm message "The system will shut down 3 seconds".</li> </ul>
4	Battery indicator	<ul> <li>On: The battery is being charged or has been fully charged.</li> <li>Off: The battery has not been installed.</li> <li>Flashing: The monitor is being powered by the battery.</li> </ul>
5	DC power indicator	Turned on when the monitor is being powered by the adapter.
6	[Alarm reset]	To reset the alarm condition.

7	[Alarm pause]	To pause the alarm sound. Alarm pause time can be set as 1, 2, 3, 4, 5, 10, 15 minutes, and permanent. Default setting is 2 minutes.
8	[NIBP start/stop]	Start and stop the non-invasive blood pressure measurement.
9	Control knob	Rotate: move the cursor.  Press: select the menu or execute a command.
10	[Setting]	Enter to the setting mode. Press again to close the setting mode.
11	₩ [Freeze]	Freeze/unfreeze the waveform.




Figure1-5: Side view

	Name	Description
1	SpO2	SpO2 cable interface
2	T1	Temperature probe interface
3	T2	Temperature probe interface
4	ECG	ECG cable interface
5	NIBP	NIBP cuff interface




Figure1-6: Rear view

	Name	Description
1	Handle	Handle for main body transport
2	Speaker holes	For alarm and synchronizing sound
3	Air outlet	Heat dissipation
4	Battery cover	Battery compartment cover

5	Bracket	To wall mount the monitor
6	ID label	Identify the monitor information
7	Air intake	For ventilation
8	Auxiliary output interface	Nurse call
9	Network port	For CMS
10	USB port	For trend or software upgrade
11	Power adapter	15V, 2.4A adapter

# 1.6 Understanding the display




Figure 1-7: Standard display

	Name	Description
1	Information area	Include patient information, alarm status icon, physiological and technical alarms. In DEOM mode, it displays "DEMO".
2	Waveform area	Mainly display the waves of physiological parameters with name of the parameter on the left side.
3	Parameter area	Show the corresponding parameter measured value and current upper and lower alarm limits of each parameter module. The parameters are shown in fixed position, that is, from top to bottom and from left to right:  - ECG - NIBP - SpO <sub>2</sub> and PR - TEMP - RESP
4	Information Tip Area	Display the network status, battery status, automatic identification screen brightness icon.
5	Hot key icons	Shows the hotkeys, which are frequently used for some common operations.
6	Date and Time area	Display the current date and time.

# 1.7 Smart Hotkeys

Smart hotkeys are graphic hotkeys displayed at the bottom of the main screen of the monitor, and enable the user to use specific features conveniently.

Кеу	Name	Description
溪	[Pause]	Alarm pause
	[Pat. Set]	Patient information setting
	[NIBP]	NIBP measurement start and stop
	[Event]	Manual event mark
	[Displays]	Change the display format
	[Freeze]	Freeze the waveform
	[Trend]	Trend display
	[Print]	Print key
	[Settings]	Setup menu
	[Volume]	Volume setup key
<b>6</b>	[Unlock]	Touch screen lock key

# 1.8 Essential performance

This device Multi-parameter Patient Monitor provides various patient vital signs such as pulse rate, ECG, respiration, blood oxygen saturation, blood pressure and temperature by placing or inserting the various sensors to the appropriate site of patient. The device is composed with display, control circuit and panel, and input part for various sensors. It detects ECG, SpO2, NIBP, etc. using ECG cable and specific probes and sensors. The detected analog signal amplifies and converted to digital. This concerted data feed to the CPU and converted to the display format as number and waveform. This device is incorporated with alarm system. The alarm generated when the detected signal range is beyond the user set alarm limits.

# 2 Preparing for operations

# 2.1 Installation

To ensure normal working of the monitor, read this chapter before use, and install as required.

#### WARNING

- All analog and digital devices connected to the monitor must be certified by IEC standards (e.g. IEC 60950 Data processing equipment standard and IEC 60601-1 Medical equipment standard). Furthermore, all configurations shall comply with valid version of IEC 60601-1 standard. The personnel connecting additional devices to the input / output signal ports are responsible for the compliance with IEC 60601-1 standard. If there is any question, please contact Bistos.
- If the patient cable interface and network interface are connected with multiple devices, the total electric leakage current cannot exceed the allowable value.
- The copyright of monitor software belongs to our company. Without permission, any organization or individual shall not interpolate, copy or exchange by any means or form.
- When the monitor is combined with other devices, it must comply with IEC 60601-1:2005 + A1:2012, and should not be connected with multiple socket outlet or extension cord.
- Do not connect the device on other equipment or network, to which a signal input/output part may be connected.

Prior to installation, the operator must ensure that the following space, power, environmental requirements are met.

#### 2.1.1 Unpack and check

BT-770 patient monitor was inspected rigorously at the factory before delivery, in order to avoid being hit when transported, carried out careful packaging. Before unpacking, carefully inspect the package. If any damage, please immediately contact the Bistos. Unpack in the correct way, carefully remove the monitor and accessories from the box and check with the packing list. Check if there is any mechanical damage, the all listed are completely packed. If you have questions, please contact the marketing department of Bistos or agency.

Please keep the packing box and materials for use in future transporting or storage.

## 2.1.2 Placement requirements

Equipment installation must meet:

- The left and right side of the monitor should have space more than 100 cm from the wall
- Back on the monitor should have space more than 50 cm.
- Ensure that the operating floor and the monitor have enough space for connecting the accessory wires.

### 2.1.3 Power requirements

DC power supply adapter

Input: A.C. 100 V - 240 V, 50/60 Hz

Output: D.C. 15 V, 2.5 A

Built-in rechargeable lithium-ion battery: D.C. 11.1 V, 4400 mAh

## 2.1.4 Environmental requirements

The storage, transport and use of the monitor must meet the following environmental requirements.

Operating	Ambient temperature	5°C ~ 40 °C	
environment	Relative humidity	30 % ~ 85 % (Non-condensing)	
	Atmospheric pressure	700 ~ 1060 mbar (hPa)	
Transportation	Prevent severe shock, vibration, rain and snow splashing during transport.		
	The packaged monitor should be stored in well-ventilated room with		
Storage	ambient temperature -20 °C $\sim$ 60°C, relative humidity 0 $\sim$ 95 % (Noncondensing), atmospheric pressure 700 $\sim$ 1060 mbar(hPa), and without		
	corrosive gases.		

The operating environment of the monitor should avoid noise, vibration, dust, corrosive or flammable and explosive materials. In order to allow air flowing smoothly and achieve good heat dissipation, at least 2 inches (5cm) clearance should be kept around the device.

When the device is moved from one environment to another, the device may have condensation due to the differences in temperature or humidity. In this case, wait until the condensation disappears before using the device.

#### WARNING

Ensure that the monitor is used under specified environment. Fail to do this, the technical specifications declared in this manual may not be met and it may result in damage to equipment and other unforeseen consequences.

# 2.2 Connecting to power

### **⚠** WARNING

- Do not try to open the monitor when the power is connecting.
- During the operation, do not disconnect any cable.

Connect to power adapter in the following steps:

- Make sure that the AC power supply meets the following specifications: a.c.100V-240V, 50/60Hz.
- Use the power adapter provided with the monitor. Plug the power adapter into the power connector of the monitor, and plug the other end of the power adapter into the mains (low voltage power supply network facilities) power outlet with protective earth.

# 3 Basic operations

## 3.1 Turn on

#### 3.1.1 Check the monitor

- Before turn on the monitor, check whether there is mechanical damage to the monitor, and whether the external cables and accessories are connected correctly.
- Plug the power adapter into the AC power outlet. If using battery power, make sure the battery is fully charged.
- Check all the functions required for patient monitoring to make sure that the monitor operates properly.



#### WARNING

If the monitor is damaged, or fails to work normally, do not use it for patient monitoring. Please contact the maintenance personnel or Bistos immediately.

#### 3.1.2 Start the monitor

If finish to check the monitor, it is ready to start the monitor.

Press the [Power] key, the yellow warning lights flash once and the system enter the program reading interface; finally the system makes a "tick" sound, the boot screen disappears, and the system enters the main interface.

- If any fatal error occurs during self-test, the system will alarm. If this case persists, please stop to using the monitor and contact the maintenance personnel or Bistos.
- Check all available monitor functions to ensure that the monitor operate properly.
- If the monitor equipped with a battery, charge the battery after each use to ensure sufficient power.
- After unpacking, when use the monitor first time, the monitor should be powered with adapter.

#### 3.1.3 Connect the sensors

Connect the required sensor to the monitor and the monitoring site of patient.

## 3.1.4 Start monitoring

Start monitoring in the following steps:

- Check if the patient cable and the sensor are connected properly.
- Check if the settings of the monitor are corrects, such as patient type.
- For the details of parameter measurement or monitoring, see the appropriate section.
- The operator can operate according to their own habits, standing in front, left or right of the monitor, easy to observe and operate the monitor.

# 3.2 Turn off

Turn off the monitor in the following steps

- Disconnect the cables and sensors connected to the patient.
- Press and hold the [Power] key for 2 seconds to pop up the 3 seconds countdown window, and the monitor turns off in 3 seconds.

# **A**CAUTION

• If the monitor is not turned off properly, you can simply disconnect the power to shutdown forcibly. But the forced shutdown may cause data loss, and it is not recommended.

# 3.3 Basic operations

## 3.3.1 Using the control knob

Control knob can be used to perform the following operations:

- Rotate: Rotate control knob clockwise or counter clockwise to move the cursor.
- Press: Press control knob to perform an action, such as access to a menu or execute a command.

Control knob is the main control means. On the interface or the menu, the green highlighted box that moves with the knob turning is called the cursor. By turning the control knob, you can position the cursor in order to perform the desired operation.

## 3.3.2 Using keys

The monitor has three types of keys:

- Soft keys: Within the display these keys allow quick access to certain menus or performing certain actions, including:
  - Parameter hotkeys: Select a parameter area and enter the appropriate parameter setup menu, including drug calculation and time setup.
  - Wave hotkeys: Select a wave area and enter the appropriate parameter setup.
  - Smart hotkeys: The shortcut keys that the user can operate quickly are displayed at bottom of the screen. Refer to '1.7 Smart Hotkeys'.
- Hard keys: The physical keys on the monitor, such as the [Alarm pause] key on the front panel.
- Popup keys: Menu keys relevant to the tasks that automatically appear on the monitor screen when need, such as, the confirmation key popped up when you need to confirm the change.

## 3.3.3 Using the touch screen

Click on the touch screen to quickly and easily perform specific operation.

## 3.3.4 Using soft keyboard

If you choose a menu which needs to enter characters, the system will display the soft keyboard on the screen. If you finish entering, press [Enter] key to confirm that you have finished entering and close the soft keyboard.

## 3.3.5 Using menu

Select the (Settings) smart key on the monitor or press the (Setting) key on the monitor panel to open the "Settings" mode as shown below. You can set-up the monitor.



Figure 3-1: "Setting" menu

The style of other menus is basically similar to the "Settings", and generally consists of the following components:

- Menu title: A title of the current menu.
- Close menu: Close the current menu. Exit the current menu or close the current menu and return to the previous menu.
- Main display area: Display options, buttons or prompt messages. The symbol ">>"
   indicates that selecting this option can enter the corresponding submenu.
- Confirmation key area: Some menus contain a confirmation key area to confirm the menu operations, including confirmation and cancel key.

# 3.4 Operation mode

The monitor has 2 operating modes, of which the demo mode is protected by a password.

1. Monitoring mode (operating mode)

This is the daily operating mode of patient monitoring; you can change some settings in accordance with the patients, such as alarm limits. However, when the patient is discharged, the monitor will restore these settings to default according to pre-set configuration.

#### 2. Demo mode

This mode is protected by a password for demonstration purpose only.

- Enter the demo mode:

  - Select "Demo Mode>>" → enter the password and confirm, and the monitor enters the demo mode.
- Exit demo mode:

  - > Select "Exit Demo >>" and the monitor exits the demo mode.

## **MARNING**

 The demo mode is mainly used to show the monitor's performance and for user training. In actual clinical use, the demo function is prohibited in order to avoid mistaking the displayed waves and parameters as those of the patient, thus affecting patient monitoring, and delaying diagnosis and treatment.

# 3.5 Measurement setup

This section only describes the general settings of measuring wave in monitor mode; for other specific settings of each parameter, please refer to the appropriate section.

Select the wave area of a parameter to enter the appropriate setup menu. The setup menu defines the specific wave setup of the parameter, such as wave gain and wave speed. You may set the waves of different parameters as needed.

# 3.6 Freezing waves

In the patient monitoring process, you can freeze the wave on the screen, review and carefully observe the patient's condition during this time. Freeze / unfreeze the wave as follows:

Select Freeze hotkey or press the [Freeze] key on the monitor panel to freeze the displayed wave of the monitor.

Select Freeze hotkey or press the [Freeze] key on the monitor panel again to release the freezing state.

# 3.7 Other common setup

The common setup of the monitor is the general setup that defines how the monitor works, for example: alarm volume setting. They may affect the setup of multiple measurements or display interfaces.

# 3.7.1 Defining the monitor

When install the monitor or change the usage occasion, the monitor should be defined as follows:

- Select "User Maintenance >>" →enter the password and confirm → "User Maintenance" menu.
  - Select "Device Name": Enter device name through the soft keyboard on the screen.
  - Select "Department": Enter the sector and department using the device through the soft keyboard on the screen.
  - Select "Bed Number": Enter the bed number through the soft keyboard on the screen.

### 3.7.2 Language setup

Set the monitor language in the following steps:

- Select [Settings]Smart Hotkey or press [Setting] key on the monitor panel → "Settings".
- Select "User Maintenance >>" →enter the password and confirm →"User Maintenance" menu.
- Select "Language", and select the option as needed:
  - "English": The interface language of the monitor is English.
  - "Türkçe": The interface language of the monitor is Turkish.
  - "Español": The interface language of the monitor is Spanish.
  - "Français": The interface language of the monitor is French.

### 3.7.3 Date and time

Set the monitor time in the following steps:

- Select [Settings]Smart Hotkey or press [Setting] key on the monitor panel → "Settings";
- ➤ Select "User Maintenance >>" →enter the password and confirm →"User Maintenance" menu.
- Select "Time Setup >>" → enter "Time Setup>>" menu.
- Or you can enter the "Time Setup" directly by touching the time display area on the display.
- "Date (YYYY-MM-DD)": Set the year, month, and day.
- "Time (24H)": Set the hour, minute and second.
- Select "Date Format", and set the date format in accordance with custom
  - "YYYY-MM-DD": Year- Month-Day.
  - "MM-DD-YYYY": Month -Day-Year.
  - "DD-MM-YYYY": Day-Month-Year.
- "Time Format", set the time format is 24H.

#### 3.7.4 Volume control

- 1. Alarm Volume
  - Select (Volume) smart hotkey → "Volume Setup" menu.
  - Select "Alarm Volume": Set alarm volume from 1 to 9.
- 2. QRS Volume
  - Select (Volume) smart hotkey → "Volume Setup" menu.
  - > Select "QRS Volume": Set QRS volume from 0 to 9. 0 means off.
- 3. Pulse Volume

  - > Select "Pulse Volume": Set pulse volume from 0 to 9. 0 means off.
- 4. Touch Volume
  - Select (Volume) smart hotkey → "Volume Setup" menu.
  - > Select "Touch Volume": Set touch volume from 0 to 9. 0 means off.
- 5. Key Volume

  - Select "Key Volume": Set key volume from 0 to 9. 0 means off.

#### 3.7.5 Setting parameter unit

You can select a preferred unit through the following operations

- Select  $\bigcirc$  [Settings]Smart Hotkey or press  $\bigcirc$  [Setting] key on the monitor panel  $\rightarrow$  "Settings".
- Select "User Maintenance >>" → enter the password and confirm →"User Maintenance" menu.
- ➤ Select "Unit Setup >>" →"Unit Setup" menu.
  - Select "Height Unit", and select the unit "cm" / "inch" as needed.
  - Select "Weight Unit", and select the unit "kg" / "lb" as needed.

- "ST Unit" fixed as "mV", is not optional.
- Select "Pressure Unit", and select the unit "mmHg" / "kPa" as needed.
- Select "TEMP Unit", and select the unit " $^{\circ}$ C" / " $^{\circ}$ F"as needed.

# 4 Patient information management

Connect the patient to the monitor, and the monitor will display and store the physiological data of the patient, so the patient can be monitored without admitting the patient. However, admitting the patient correctly is very important.

If the monitor has admitted the patient, it is recommended to operate the monitor to discharge the current patient before connecting to (not admitted) the next patient. Otherwise, the data of the previous patient will be stored in the data of the current patient.

#### WARNING

- Whether the patient is admitted or not, the system will give a default value to "Patient Type" and "Pace Maker", "Patient Type" default "Adult", "Pace Maker" default "No", and the user must confirm that the default value is appropriate for the patient being monitored.
- For patients with pacemakers, "Pace Maker" must be set to "Yes". Otherwise, the pacing pulse will be treated as normal QRS wave group, and the system is unable to detect the alarm status of "ECG Signal weak".
- For patients without a pacemaker, "Pace Maker" must be set to "No". Otherwise, the system is unable to detect the arrhythmias (including PVCs count) related to ventricular premature beats, and fails to perform ST segment analysis.

# 4.1 Patient setup menu

You can manage the patient through the "Patient" menu. To enter "Patient" menu, operate as follows:

Select (Settings)Smart Hotkey or press (Setting) key on the m	onitor panel → "Settings"
→ "Patient >>" → "Patient" menu;	

Or

Select [Pat. Set] Smart Hotkey to enter "Patient" menu, as shown in Fig. 4-1.



Figure 4-1 "Patient" menu

# 4.2 Admitting a patient

Admit a patient as follows:

In "Patient" menu, select "Quick Admit"  $\rightarrow$  "Warning" message  $\rightarrow$  "OK"  $\rightarrow$  "Quick Admit" menu, as shown in Figure 4-2.



Figure 4-2 "Quick Admit" menu

- > Select "Patient Type", and set the patient category as needed: "Adult" and "Pediatric".
- > Select "Pace Maker", and set whether the patient wears a pacemaker according to the patient condition: "Yes" or "No".
- After setting, select "OK" to save the current setup or select "Cancel" and do not save the current setup.

# 4.3 Patient information

To edit patient information, operate as follows:

In the "Patient" menu, select "Patient Info". The "Patient Info" menu as shown in Figure 4-3 will be displayed.




Figure 4-3. "Patient Info" menu

- 1. Select "Last Name", and enter patient's surname through the soft keyboard (Letters: not more than 20 characters).
- 2. Select "First Name", and enter patient name through the soft keyboard ( Letters: not more than 20 characters ) .
- 3. Select "Patient ID", and enter the patient ID through the soft keyboard (Letters: not more than 20 characters).
- 4. Select "Case Number", and enter the case number through the soft keyboard ( Letters: not more than 20 characters ) .
- 5. Select "Gender", and set the patient's gender.
- 6. Select "Patient Type", and set the patient category as needed: Adult and Pediatric.
- 7. Select "Pace Maker", and set whether the patient wears a pacemaker.

- 8. Select "Height(cm)", and set the patient's height via the pop-up keyboard on the screen(Range: 0 ~ 250).
- 9. Select "Weight (kg)", and set the patient's weight via the pop-up keyboard on the screen(Range: 0 ~ 350).
- 10. Select "Blood Type", and set the patient's blood type: A, B, AB or O.
- 11. Select "Admission Date (MM-DD-YYYY)", and set the date of admitting the patient.
- 12. Select "Birthday (MM-DD-YYYY)", and set the birth date of the patient.

After setting, select "OK" to save the current setting or select "Cancel" and do not save the current setting.

# 4.4 Discharging a patient

To discharge a patient, operate as follows:

In the "Patient" menu, select "Discharge Patient"  $\rightarrow$  "Warning" message  $\rightarrow$  "OK" to finish the operation of discharging a patient.

After the patient is discharged, all the information of the patient stored in the monitor will be cleared. Therefore, discharge the patient only when needed.

# 4.5 Clear alarms

To clear alarms, operate as follows:

In the "Patient" menu, select "Clear Alarms"  $\rightarrow$  "Warning" message  $\rightarrow$  "OK" to finish the operation of clear alarms.

After the alarm is cleared, all the information of alarms stored in the monitor will be cleared. Therefore, clear alarm only when needed.

# 4.6 Clear trend

To clear trend, operate as follows:

In the "Patient" menu, select "Clear Tabular Trend"  $\rightarrow$  "Warning" message  $\rightarrow$  "OK" to finish the operation of clear tabular trend.

After the tabular trend was cleared, all the information of tabular trend stored in the monitor will be cleared. Therefore, clear tabular trend only when needed.

P/N : 770-ENG-OPM-EUR-R00 Bistos Co., Ltd. 2018.09

# 4.7 Clear NIBP trend

To clear NIBP trend, operate as follows:

In the "Patient" menu, select "Clear NIBP Trend"  $\rightarrow$  "Warning" message  $\rightarrow$  "OK" to finish the operation of clear NIBP trend.

After the NIBP trend was cleared, all the information of NIBP trend stored in the monitor will be cleared. Therefore, clear NIBP trend only when needed.

# **5** Display format

The monitor has four display format, which are "Normal Screen", "Big ECG Screen", "Big font Screen", and "ECG 7-Lead Full-Screen". The user can select the display format according to needs, and get different screen information.

# 5.1 Selecting user interface

Select the user interface as follows:

- Select □ [Displays] smart hotkey → Screen Select;
- Select the display format according to needs:
  - "Normal Screen": Standard interface.
  - "Big ECG Screen": Big ECG interface.
  - "Big font Screen": Big font interface.
  - "ECG 7-Lead Full-Screen": ECG 7-Lead Full interface.

# 5.2 Display description

# 5.2.1 Normal display format




Figure 5-1: Standard Display

The normal display provides the parameter wave being monitored and the parameters displayed in the parameter area. This is the basic display of the monitor. In this display mode all parameters, two ECG waves, one blood oxygen saturation percentage wave, one respiratory wave are displayed.

## 5.2.2 Big ECG format

The big ECG format is as shown in Figure 5-2.




Figure 5-2: Big ECG format

# 5.2.3 Big font format

The big font format is as shown in Figure 5-3.



Figure 5-3: Big font format

## 5.2.4 ECG 7-Lead full screen format

The ECG 7-Lead full screen format is as shown in Figure 5-4.




Figure 5-4: ECG 7-Lead full screen format

# 6 Alarm

Alarm means that the monitor prompts the medical staff through sound and light when the abnormal changes in vital signs are monitored or the monitor has a failure or is unable to monitor the patient successfully.

#### WARNING

- In any single region (e.g. ICU), it has potential danger if the same or similar devices use different alarm setup.
- After setting, the alarm and other parameters of the monitor won't be lost when the system is power off, unless modified manually. Connect the power again and turn on the monitor, it will resume normal working, and the alarm and other parameters remain unchanged.

# 6.1 Alarm types

According to the nature of the alarm, the alarms of the monitor can be divided into physiological alarms, technical alarms and prompt messages.

#### Physiological alarms

A physiological alarm is usually triggered when a physiological parameter of the patient exceeds the alarm limit or the patient has physiological abnormalities. The information of physiological alarm is displayed in the physiological alarm area on top of the screen.

#### Technical alarms

Technical alarm is also known as a system error message, which is caused by improper operation or system failure resulting in system malfunction or monitoring result distorted. The information of technical alarm is displayed in the technical alarm area on top of the screen.

#### Prompt messages

Strictly speaking, the prompt messages are not alarms. The monitor also will display some information associated with system status in addition to the physiological alarms and technical alarms, and generally such information do not involve the patient's vital signs. The prompt messages generally appear in the technical alarm area and parameters area.

# 6.2 Alarm condition priorities

According to the severity of the alarm conditions, the physiological alarms of the monitor can be divided into high priority, medium priority and low priority.

#### High priority alarms

The patient is in critical condition that is life-threatening, and should be immediately rescued, or the monitor has a serious mechanical failure or malfunction, causing it unable to detect the patient's critical state and endangering the patient's life.

#### Medium priority alarms

The patient's physiological signs are abnormal and appropriate measures or treatment should be taken immediately, or although it won't endanger the patient's life, the mechanical failure or disoperation of the monitor will affect the normal monitoring of key physiological parameters.

#### Low priority alarms

The patient's physiological signs are abnormal and appropriate measures or treatment may need to be taken, or certain monitoring function is invalid due to mechanical failure or disoperation, but it won't endanger the patient's life.

The priority of all technical alarms and some physiological alarms have been set in the monitor at the factory and cannot be modified by the user. The levels of some physiological alarms can be modified.

## 6.3 Alarm mode

When an alarm occurs, the monitor uses the following audible or visual alarm to prompt the user:

- Visual alarm
- Audible alarm
- Alarm info
- Parameter flashing

Of which, the visual alarm, audible alarm, and alarm information distinguish the alarm levels in a different manner respectively.

#### 6.3.1 Visual alarm

When an alarm occurs, the alarm indicator will flash in different colors and frequencies to prompt the alarm priority.

- ➤ High priority alarm: Red, fast flashes.
- Medium priority alarm: Yellow, slow flashes.
- Low priority alarm: Yellow, lit without flashing.

### 6.3.2 Audible alarm

An audible alarm is that the monitor prompts the alarm priorities with different sound characteristics when an alarm occurs.

Medium priority alarm: Beep-beep-beep

Low priority alarm: Beep

#### 6.3.3 Alarm information

Alarm information displayed on the physiological or technical alarm area of the monitor indicates the corresponding alarm information when an alarm occurs. The system will distinguish the alarm priority with different background colors:

➤ High priority alarm: Red

Medium priority alarm: Yellow

Low priority alarm: Yellow

The following flags in front of physiological alarms are used to distinguish the alarm priorities.

High priority alarm: \*\*\*

Medium priority alarm: \*\*

Low priority alarm: \*

### 6.3.4 Parameter flashing

When the physiological parameter values in the parameter area will flash once per second, and the upper limit and lower limit of the parameter will also flash at the same frequency, it indicating that the parameter exceeds the upper limit or lower limit.

# 6.4 Alarm states

In addition to the above alarm modes, you can also set the monitor to the following three alarm states as needed, and display different alarm icons on the screen:

- > Alarm Reset
- > Alarm sound off
- > Alarm pause
- > Alarm off

#### 6.4.1 Alarm reset

Select button, and you can temporarily turn off the alarm sound of currently occurring physiological alarms of the monitor, but the alarm information is still retained. For technical alarms, clear the alarm state, display alarm prompt information, the alarm state icon area displays the icon. When a new physiological alarm or technical alarm occurs, the alarm reset is automatically canceled.

#### 6.4.2 Alarm sound off

The alarm sound can be turned off through the following operations:

- Select "User Maintenance >>"→enter the password and confirm →"User Maintenance" menu.
- Select "Alarm Param >>"→" Alarm Param" menu.
- > Set "Minimum Alarm Volume" to "0". "Minimum Alarm Volume" range from 0 to 9, the default value is 1.
- Select (Volume) smart hotkey → "Volume Setup" menu.
- > Set "Alarm Volume" to "0".

When the alarm sound is turned off, the alarm state area on the screen shows the icon.

If "Minimum Alarm Volume" is larger than 0, the system will cancel alarm sound off state.

#### WARNING

• When the alarm is off, and the alarm reminder signal is on, the system will have alarm reminder tone.

#### 6.4.3 Alarm pause

Press the [Pause] smart hotkey or [Alarm pause] key on the monitor panel to temporarily stop the alarm of the monitor in the following steps:

- Pause]smart hotkey will appear magnified and reverse colored icon.
- > The light alarm and audible alarm of the physiological alarms are suspended, and the alarm information is not displayed.
- > The remaining time of alarm pause is displayed in the physiological alarm area.
- Alarm parameters and upper / lower limit stop flashing.
- > The audible alarm and light alarm of technical alarms are suspended, but the alarm message is still displayed.

After the alarm pause is finished, the monitor will automatically cancel the alarm pause state. During the alarm pause, you can also press [Pause] smart hotkey or [Alarm pause] key on the monitor panel to cancel the alarm pause manually.

You can set the alarm pause time as follows:

- Select (Settings]Smart Hotkey or press (Setting) key on the monitor panel → "Settings".
- Select "User Maintenance >>" → enter the password and confirm →"User Maintenance" menu.
- Select "Alarm Param >>"→" Alarm Param" menu.
- Select "Alarm Pause Time", and set the alarm pause time.
  - "1min" /"2min" /"3min" /"4min" /"5min" /"10min" /"15min" "Permanent". By default, the alarm pause time is 2 minutes.
  - "Permanent" means alarm off.
  - It is recommended that the SpO2 alarm pause time shall not more than two minutes.

#### 6.4.4 Alarm off

As shown in 6.4.3, if the "Alarm Pause Time" is set to "Permanent", press the [Pause] smart hotkey or [Alarm pause] key on the monitor panel, and the monitor will turn off the alarm. In this case, except the alarm prompt characteristics maintained in alarm pause state:

- 🌣 [Pause] smart hotkey will appear magnified 🔀 icon.
- The physiological alarm area displays "Alarm Pause".

You can press the [Pause] smart hotkey again to manually cancel the alarm off.

If the monitor is in the alarm state of suspension or high priority technical alarm is triggered, the alarm and the alarm off pause are automatically canceled.

#### WARNING

When the alarm volume is set to '0' or the alarm pause time is set to permanent, the monitor does not sound an alarm when an alarm occurs. Therefore, the operator should use this feature carefully.

# 6.5 Alarm setup

# 6.5.1 Setting the alarm delay time

To limit alarm of continuous measurement parameter, you can set the alarm delay time. If the alarm condition disappears during the delay period, the monitor will not generate an alarm. In "Alarm Param" menu, select "Alarm Delay" time and "ST Alarm Delay" time.

The specific operation is as follows:

- Select (Settings) smart hotkey or press (Setting) key on the monitor panel → "Settings".
- Select "User Maintenance >>" → enter the password and confirm → "User Maintenance" menu.
- Select "Alarm Param >>" → "Alarm Param" menu.
- > Select "Alarm Delay", and set the alarm delay time as needed:
  - "Off": Turn off the alarm delay.
  - "1s" / "2s" / "3s" / "4s" / "5s" / "6s" / "7s" / "8s": Alarm delay time is 1 sec, 2 sec, 3 sec, 4 sec, 5 sec, 6 sec, 7 sec or 8 sec. By default, the alarm delay time is 4 seconds.
- Select "ST Alarm Delay", and set the ST alarm delay time as needed.
  - "Off": ST alarm delay is off
  - "10s" / "20s" / "30s" / "45s" / "1min" / "2min" / "3min": ST alarm delay time is 10 sec, 20 sec, 30 sec, 45 sec, 1 min, 2 min or 3 min. By default, the ST alarm delay time is 20 seconds.

### 6.5.2 Setting the alarm reminder signal and alarm reminder interval

The alarm reminder signal can be turned on or off. When the alarm is off and the alarm sound is off, and then the alarm reminder signal is on. You can set the alarm reminder interval as needed: "1min" / "2 min" / "3 min".

The specific operation is as follows:

- Select "User Maintenance >>" →enter the password and confirm → "User Maintenance" menu.
- ➤ Select "Alarm Param >>" → "Alarm Param" menu.
- Select "Alarm Reminder Signal", and set the alarm reminder signal as needed:

- "On": The alarm Reminder Signal is on.
- "Off": The alarm Reminder Signal is off.
- Select "Alarm Reminder Interval", and set the alarm reminder interval as needed:
  - "1min" / "2 min" / "3 min": Alarm reminder interval is 1 min, 2 min or 3 min. By default, the alarm reminder interval is 3 min.

#### 6.5.3 Setting a parameter alarm

You can set the parameter alarm for every alarm separately. For SpO<sub>2</sub>, as an example, select "Alarm Setup >>" in the "Settings" menu and select "SpO<sub>2</sub>" and enter the SpO<sub>2</sub> alarm setup menu.

- 1. Turn on / off alarm
- > Select "Alarm Switch" and set the alarm switch as follows:
  - "On": Turn on SpO<sub>2</sub> alarm; when the parameter alarm occurs, the monitor will prompt according to the set alarm level.
  - "Off": Turn off  $SpO_2$  alarm;  $\bowtie$  icon is displayed in the parameter area, and the monitor won't prompt the parameter alarm.
- 2. Set the alarm priority
- Select "Alarm Level", and set the alarm priority as follows:
  - "Mid": Set the alarm priority to medium.
  - "High": Set the alarm priority to high.

Note: Regulatory requirements, the parameter (ECG, blood oxygen saturation, blood pressure) can set the alarm priority high and mid.

#### 3. Set the alarm limit

In any cases, the alarm system only allows setting the values within the effective range of the system, and the upper alarm limit must be higher than the lower alarm limit.

- ➤ Select "SpO₂ Low Limit" and set the lower limit of SpO₂ alarm.
- ➤ Select "SpO₂ High Limit" and set the upper limit of SpO₂ alarm.
- > Select "PR Low Limit" and set the lower limit of PR alarm.
- Select "PR High Limit" and set the upper limit of PR alarm.

Туре	Adults		Pediatric	
	Range	Default	Range	Default
SpO2 Low Limit	0-99	90	0-99	90
SpO2 High Limit	1-100	100	1-100	100
PR Low Limit	15-299	50	15-349	75
PR High Limit	16-300	120	16-350	160

- 4. Restore default alarm setup
- > Select "Default", and restore the alarm setup to the factory setup.

#### NOTE

- When setting the upper and lower alarm limits, confirm the patient category and set its range according to the clinical need. If the setting exceeds the alarm limits, the alarm system will fail easily.
- When the alarm limit is turned on, and the upper and lower alarm limits are manually set, the monitor will display the upper and lower alarm limits continuously, and the initial alarm preset value will not be provided additionally.

# 6.6 Latch alarm

Physiological alarms can be set to "Latching" or "No latching".

- Latching": Even if the cause of physiological alarm is cleared, the system will still be "latched", that is, continue to display the alarm information corresponding to physiological alarm, the alarm sound also continues, but the alarm mode has the following changes:
  - Parameters and upper or lower alarm limit are no longer flashing.
  - Display the time that the latest alarm was triggered after the alarm message in the physiological alarm area.
- "No latching": After the causes of physiological alarm are cleared, the system will no longer prompt the physiological alarm.

The default alarm of the system is no-latching alarm. You can set the alarm as latching or no-latching in the following steps.

- Select (Settings) smart hotkey or press (Setting) key on the monitor panel → "Settings".
- Select "User Maintenance >>" →enter the password and confirm → "User Maintenance" menu.
- ➤ Select "Alarm Param >>" → "Alarm Param" menu.
- Select "Latching Alarm", and set the alarm as needed:
  - "Latching": Latching alarm.
  - "No latching": Non-latching alarm.

# 6.7 Manual event

In the patient monitoring process, some events may have an impact on the patient, resulting in changes of some monitoring waves or parameters. In order to assist in the analysis of these effects, you can manually record these events through the [Event] smart hotkey, and then view it in the event review, refer to 15.4 Event Review for detailed operation.

# 6.8 Alarm record

When the monitor's machine alarm system is powered down, all alarm records are not saved.

Physiological alarm can store 200 alarm records, if full of 200, the latest alarm records will replace the beginning of the record;

Technical alarm can store 100 alarm records, if full of 100, the latest alarm records will replace the beginning of the record.

# **7 ECG**

### 7.1 Overview

Electrocardiogram (ECG) is produced by the continuous electrical activity of the patient's heart, and displayed with wave and numeric on the monitor in order to accurately assess the physiological state of the patient at the time. The ECG cable should be connected properly, so as to obtain a correct measurement value and normal display. This monitor can simultaneously display 7 ECG waves.

Patient cable consists of two parts.

- Wires connected to the monitor
- ECG electrodes connected to the patient

Connect to the monitor with five lead ECG cable, and ECG can display two different waves by adjusting the two leads. You can use the control knob to change the lead name on the left of the ECG wave on the screen and select the lead to be monitored.

The parameters displayed in the parameter area of the monitor include heart rate (HR), ST segment measurements and arrhythmia counts per minute. All these parameters can be used as alarm parameters.

The monitor is designed for defibrillation proof, so the monitor operates normally after defibrillation.

#### NOTE

• In the factory setup, ECG wave display in the first two waves from top in the wave area in the normal display format.

# 7.2 Safety information

#### WARNING

- To monitor ECG signal, ECG cable and ECG electrodes specified in this manual must be used.
- When connecting the electrodes or patient cable, make sure that the patient is absolutely not connected with any other conductive parts or in contact with the ground. In particular, make sure that all the ECG electrodes, including the neutral electrodes, are attached to the patient and prevent them from contact with the conductive parts or ground.
- When using electrosurgical (ES) equipment, users should put ECG electrodes at middle of the ES earthing plate and ES knives to prevent from burns. Cables of ES equipment cannot be wrapped with ECG cables together.
- During use of ES equipment, don't put electrodes near the earthing plate of such equipment, otherwise ECG signals will be much disturbed.
- For patients who wear a pace maker, pacing pulse analysis must be turned on. Otherwise, the pacing pulse may be counted as a normal QRS wave, make the ECG signal too weak to detect the alarm.
- Periodically check the skin that the electrode is placed at. If there is any sign of allergy or irritation, replace the electrode or change the placement position.
- Electrosurgical (ESU) device interference, defibrillator discharge
  - When the patient needs defibrillation, do not use non-defibrillator type ECG cables. For defibrillation protection, please use the accessories specified by manufacturer. (Refer to Chapter 17. Accessories)
  - During defibrillation, the operating personnel shall not touch the patient, tables and instrument.
  - During defibrillation, the ECG cable connected with the patient's body may be damaged. Check if the function is normal again before using these cables.
  - The monitor will recover within 10 seconds after defibrillation and will not lose any stored data. During electrosurgery or defibrillation, the measurement accuracy may be temporarily reduced. This does not affect the safety of the patient or the instrument.
- Do not expose the monitor to X-ray or strong magnetic fields (e.g. MRI).

# 7.3 Monitoring steps

### 7.3.1 Preparation

Before placing the electrode, prepare the patient's skin in the following steps.

- Skin preparation: Since the skin is a poor conductor, it is very important to treat the patient's skin for electrode placement appropriately to make good contact between the electrode and the skin. Select the flat position with less muscles for the electrode placement, and refer to the method below for treatment of the skin:
  - Remove the body hair at the position for electrode placement.
  - Gently rub the skin at the position for electrode placement to remove dead skin cells.
  - Wash the skin thoroughly with soap and water (do not use ether and pure alcohol, as this will increase the skin's impedance).
  - Dry the skin completely before placing the electrode.
- Install the spring clip or stud prior to the placement of the electrodes.
- > Place the electrode on the patient.
- Connect the ECG cable and ECG interface.

#### **⚠** WARNING

 Check if the lead is adequately attached and do not have any damage before monitoring. When the ECG cable is unplugged, the screen will display "ECG Lead Off" prompt, and trigger an audible and visual alarm.

#### 7.3.2 Selecting lead

- ightharpoonup Select the ECG parameter area or wave area ightharpoonup "ECG Setup" menu.
- Select "Other Setup >>" → "ECG Other Setup" menu.
- Select "Lead Type", and select the ECG lead as needed.
  - "3-Lead": 3-lead; ECG wave options: I, II, III.
  - "5-Lead": 5-lead; ECG wave options: I, II, III, AVR, AVL, AVF, V.

## 7.3.3 Lead name and corresponding color

The lead names in European standard and U.S. standard (represented with R, L, N, F, C in European Standard, and represented with RA, LA, RL, LL, V in the U.S. standard) are shown in Table 9 -1.

Table 7-1: Lead Name in European Standard and American Standard

European Standard (EN)		American Standard (AHA)		
Lead Name	Color	Lead Name	Color	
R	Red	RA	White	
L	Yellow	LA	Black	
F	Green	LL	Red	
N	Black	RL	Green	
С	White	V	Brown	

# 7.3.4 Installing the electrodes

#### > 3-lead

The electrode placement position of 3-lead is shown in Fig. 7-1.

- R/RA electrode: placed below the clavicle, near the right shoulder.
- L/LA electrode: placed below the clavicle, near the left shoulder.
- F/LL electrode: placed on the left abdomen.




Figure 7-1: 3-Lead placement method

#### 5-lead

The electrode placement position of 5-lead is shown in Fig. 9-2:

- R/RA electrode: placed below the clavicle, near the right shoulder.
- L/LA electrode: placed below the clavicle, near the left shoulder.
- N/RL electrode: placed on the right abdomen.
- F/LL electrode: placed on the left abdomen.
- C/V electrode: placed on the chest wall.




Figure 7-2: 5-Lead placement method

#### NOTE


- To ensure the patient safety, all leads must be connected to the patient.
- If the electrodes are attached correctly, but the ECG wave is not accurate, then replace the lead.
- Interference from ungrounded instrument near the patient and ESU may cause waveform problem.

#### 7.3.5 Checking the pacemaker

Before ECG monitoring, it is very important to set the pace maker state of the patient properly. If the patient has a pacemaker, set "Pace Maker" to "Yes", and the icon displays in the patient information area. When the system detects a pacing signal, the "| "symbol will be marked in the top of the ECG wave.

You can change the pacing state in the following method:

Select the patient information area to pop up the "Patient Info" menu, or select the [Pat. Sat] smart hotkey and select "Patient Info" menu.



Diagnostic, Monitor, Surgery will not affect rejection of pacemaker pulses.

## **⚠** WARNING

• For patients with pacemakers, the cardio tachometer may count the pacemaker pulse in case of a cardiac arrest or arrhythmias. Never rely solely on the cardio tachometer alarm. Closely monitor the patients with pacemaker. For the inhibition of the device on pacemaker, refers to this manual.

# 7.4 ECG display

#### ECG wave display

The monitor displays two ECG waves on the normal screen. Fig. 7-3 below is the monitoring interface of 5-lead, and is for reference purposes only. The graphics displayed on your monitor may be slightly different.




Figure 7-3: ECG wave in normal display format

In addition, when "Pace Maker" is set to "Yes", and the patient wears a pacemaker, the "  $\mid$  "symbol will be marked in the top of the ECG wave.

### > ECG parameter display

The ECG parameter area of the monitor in the normal screen is shown in Fig. 9-4:




Figure 7-4: ECG parameter in standard display format

# 7.5 ECG setup

Select the ECG parameter area or wave ECG area or select the [Settings] smart hotkey and "Parameter Setup >>" and "ECG Setup >>" to pop up the "ECG Setup" menu, which is as shown below. You can set the ECG through the "ECG Setup" menu.




Figure 7-5: "ECG Setup" menu

- Select "HR Source", and set the heart rate source.
  - "Auto": Automatically select HR source.
  - "ECG": Select ECG monitoring as the HR source.
  - "SpO2": Select SpO2 monitoring as the HR source.
- Select "Calculated Channel" and select heart rate calculated channel.
  - "Auto": Automatically select Heart Rate calculated Channel.
  - "I": Select the first ECG waveform as the Heart Rate calculated Channel.
  - "II": Select the second ECG waveform as the Heart Rate calculated Channel.
  - "V": Select the third ECG waveform as the Heart Rate calculated Channel.
- > Select "ECG1" / "ECG2" to set the display wave channel. Select "ECG1" / "ECG2", and set the names of upper ECG wave and lower ECG wave on the screen.
  - ECG1/ECG2 should not be the source of the same wave, source waveform can be set I/II/III/AVR/AVL/AVF/V.
- Select "Wave Gain", and set the ECG wave gain. When the wave is shorter, increase the

wave gain factor appropriately; when the wave is high or the peak cannot be displayed, reduce the wave gain appropriately, gain can be set  $Auto/\times0.25/\times0.5/\times1/\times2$ . When Big ECG interface is selected, gain can be set  $Auto/\times0.25/\times0.5/\times1/\times2/\times4$ .

- Select "Wave Speed", and set the wave speed. The wave speed is "12.5mm/s"/"25mm/s"/"50mm/s". The default is 25mm/s.
- > Select "Filter", and set the filter mode:
  - "Diagnostic": Diagnostic mode
  - "Monitor": Monitor mode
  - "Surgery": Surgery mode
- > Select "Arrhythmia Analysis>>", and set the alarm switch, alarm level, alarm record.
- Select "ST Analysis>>", and set the ST analysis, ST Channel, ST Alarm Setup.
  - "ST analysis": You can set "Off" or "On".
  - "ST Channel": You can set "1", "2", "3".
  - "ST Alarm Setup": and set the alarm switch, alarm level, alarm record.
- > Select "Other Setup>>", and set the QRS Volume, Lead Type, Notch Filter, Pace Maker.

# 7.6 Alarm setup

Select "Alarm Setup >>"  $\rightarrow$  "Alarm Setup" interface to set ECG related alarms; see 6.5 Alarm Setup for the setting method.

# 8 RESP

# 8.1 Overview

Thoracic electrical bio impedance is a method used for measuring the respiration. When the patient is breathing, the thoracic impedance between two ECG electrodes changes due to thoracic activity. The monitor generates a respiratory wave on the screen by measuring the changing impedance value. The monitor calculates the respiration rate (RR) according to the wave cycle.

# 8.2 Safety information

#### NOTE

Respiration monitoring does not apply to patient with large range of activities, as this may lead to false alarm.

#### WARNING

- Do not use anti-electric knife ECG cable for respiration monitoring.
- Respiration measurement cannot identify the apnea because it will alarm if the next respiration is not detected in the predetermined period after last respiration, and therefore it cannot be used for diagnostic purpose.

# 8.3 Placing electrodes for respiration monitoring

Since the skin is a poor conductor, it is very important to treat the patient's skin for electrode placement appropriately to get better respiration signals. Refer to 7.3.1.

Respiration measurement uses standard ECG cable and electrode placement method. You can use different ECG cables (3-lead or 5-lead). Respiratory signal is measured between two ECG electrodes. If standard ECG electrode position is used, the two electrodes are R (right arm) and L (left arm) electrodes of I lead or R (right arm) and F (left leg) electrode of II lead.

#### NOTE

• For optimal respiration wave, R and L electrodes should be placed horizontally if I lead is selected for respiration measurement. R and F electrodes should be places diagonally if II lead is selected for respiration measurement.

Figure 8-1 below shows the placement of 5-lead electrodes.




Figure 8-1: 5-lead respiration electrode placement

# 8.3.1 Adjusting position of respiration electrode

If you want to measure ECG and respiration simultaneously, you may need to adjust the position of the two electrodes for respiration measurement.

#### NOTE

 Adjusting the standard position of ECG electrodes will lead to changes in the ECG wave, and may affect the ST and arrhythmia analysis.

## 8.3.2 Cardiomotility superimpose

The effect of cardiomotility on the respiratory wave is called cardiomotility superimposing. When the respiration electrodes collect impedance changes caused by rhythmic blood flow, this will happen. Placing the respiration electrodes correctly will reduce this effect. The liver and ventricle should avoid the connection of respiration electrode, so that the heart or pulsating flow won't generate artifact.

# 8.4 Respiration display

Respiration wave is displayed as shown in Figure 8-2.




Figure 8-2: Respiration wave

Respiration parameters are displayed as shown in Figure 10-3.




Figure 8-3: Respiration parameter display

# 8.5 Respiration setup

Select the RESP parameter zone or respiration wave area  $\rightarrow$  "RESP Setup" menu, which is shown below. You can set respiration through "RESP Setup" menu.



Figure 8-4: "RESP Setup" menu

# 8.5.1 Setting apnea time

Apnea alarm is a high level alarm for monitoring the apnea. In "RESP Setup" menu, set "Apnea Delay" to an appropriate value and set the apnea alarm time. When the apnea time of the patient is longer than the set time, the monitor will trigger an alarm. Set time can be set "20s" / "25s" / "30s" / "35s" / "40s" / "45s" / "50s" / "55s" / "60s", default apnea alarm time is 20s.

#### 8.5.2 Adjusting wave gain

In "RESP Setup" menu, select "Wave Gain", and set the wave gain: the greater gain, the higher wave amplitude. Gain can be set "x0.5" / "x1" / "x2", the default is "x1".

# 8.5.3 Setting sweep speed

In "RESP Setup" menu, select "Wave speed", and set the sweep speed: the faster sweep speed, the smoother wave. The wave speed is "6.25mm/s" / "12.5mm/s" / "25mm/s". The default is "12.5mm/s".

### 8.5.4 Setting calculated channel

In "RESP Setup" menu, select "Calculated Channel", and set the calculated channel. The calculated channels are "RA-LA", "RA-LL", "LA-RL" and "LL-RL". The default is "RA-LA".

# 8.5.5 Setting sensitivity

In "RESP Setup" menu, select "sensitivity", and set the sensitivity. The sensitivity is "1" / "2" / "3" / "4" / "5". The default is "2".

# 8.6 Alarm setup

Select "Alarm Setup >>"→ "Alarm Setup" interface to set respiration related alarms; see 6.5 Alarm Setup for the setting method.

# 9 PR

### 9.1 Overview

The mechanical activity of the heart causes arterial pulsation, and PR (pulse rate) value can be obtained by measuring this pulsation. PR value can be obtained through SpO<sub>2</sub> measurement.

The average calculation of the heart rate is the direct average method. The refresh rate is every 1 second.

# 9.2 Display

The color of PR parameter area is same as SpO<sub>2</sub> parameter color of PR source, as shown in Fig. 9-1:




Figure 9-1: PR parameter display

# 9.3 Setting PR sound

Select SpO<sub>2</sub> parameter area or SpO<sub>2</sub> wave area  $\rightarrow$  "SpO<sub>2</sub> Setup" menu;

Select "Pulse Volume" to set "Pulse Volume" to 0~9. Select 0 to turn off the pulse volume, and select 9 to set the maximum volume.

#### NOTE

HR sound has higher priority than PR sound. When HR makes a sound, PR won't.
 When HR sound set to 0, PR can make a sound.

# 9.4 Alarm setup

Select PR parameter area  $\rightarrow$  "SpO<sub>2</sub> Setup" menu  $\rightarrow$  "Alarm Setup >>" to enter the "Alarm Setup" interface, and set PR alarm switch, alarm level and upper/lower alarm limit. See 6.5 Alarm Setup for detailed setting method.

# **10** SpO<sub>2</sub>

### 10.1 Overview

Blood oxygen saturation (SpO<sub>2</sub>) is the percentage of oxyhemoglobin (HbO2) capacity bound by oxygen in the blood in the total hemoglobin (Hb) capacity that can be combined, that is, the concentration of oxygen in the blood.

The principle for monitoring the pulse SpO<sub>2</sub> is to fix the probe fingerstall on the patient's finger or toe, use the finger (or toe) as a transparent container for hemoglobin, use 660nm wavelength red light and 950nm near-infrared light as the incident light, maximum output power is 300 mW, measure the light transmission intensity through the tissue bed, and calculate the concentration of hemoglobin and SpO<sub>2</sub>.

The passing lights depend on a variety of factors, most of which are constant. However, one of these factors, the arterial blood flow, changes with time, as it is pulsating. By measuring the light absorbed during pulsating, it is possible to obtain the arterial blood SpO<sub>2</sub>. Detection pulsation can give a "plethysmography" wave and pulse rate signal.

The main screen displays "SpO<sub>2</sub>" value and "plethysmography" wave.

This monitor applies to measure SpO<sub>2</sub> of adults (>18 years) and pediatric (<18 years,>30 days). Contact SpO<sub>2</sub> probe to Patient's finger (or toe) to get "SpO<sub>2</sub>" value and "plethysmography" wave.

SpO<sub>2</sub> function of this monitor has been calibrated in factory.

The monitor is defibrillation proof, so the monitor operates normally after defibrillation.

# 10.2 Safety information

## WARNING

- Please use SpO2 sensor specified in this Manual, operate in accordance with the Manual, and observe all warnings and precautions.
- Before monitoring, check whether the sensor cable is normal. When SpO2 sensor cable is unplugged from the socket, the screen will display "SpO2 Sensor Off" error message, and trigger an audible and visual alarm simultaneously.
- If the sensor or sensor packaging has signs of damage, do not use this SpO2 sensor; return it to the manufacturer.
- If there is carboxyhemoglobin, methemoglobin or dye diluted chemical, the SpO2 value will have deviation.
- When the patient has a tendency to hypoxia, use the oximeter to analyze blood samples in order to fully grasp the patient's condition.

- Do not put the sensor on limbs with arterial duct or intravenous tube.
- Do not intertwine electrosurgical equipment cable with the sensor cable.
- Avoid using the monitor and sensors while using the NMR equipment, in order to avoid severe burns to the patient as a result of induced currents.
- During long time continuous monitoring of a patient, check the position of SpO2 sensor once every 2 hours, and move properly when the skin changes or every four hours. Some patients may require more frequent inspection, such as patients with perfusion disorders or sensitive skin, because persistent and prolonged monitoring may increase unpredictable skin changes, such as allergies, redness, blistering or pressure necrosis.
- When the measured pulse rate is not complete, then the "---".
- Before using, verify compatibility between the monitor, probe and cable, otherwise it may cause injury to the patient.
- Functional testers cannot be used to evaluate the accuracy of pulse oximetry and pulse oximetry.
- SpO2 low alarm limit cannot be less than 85.

#### NOTE

- Do not put the oxygen probe and blood pressure cuff on the same limb, because blood flow occlusion during blood pressure measurement will affect the SpO<sub>2</sub> readings.
- The monitor cannot be used to verify the accuracy of SpO<sub>2</sub> probe and SpO<sub>2</sub> equipment.

# 10.3 Monitoring steps

- 1. Select the appropriate SpO<sub>2</sub> sensor according to the patient.
- 2. Turn on the monitor, and connect the  $SpO_2$  lead wire to the monitor.
- 3. Clean the measurement site, such as finger with nail polish.
- 4. Put the SpO<sub>2</sub> sensor probe on the patient's body.
- 5. Select the appropriate alarm settings.
- 6. Start monitoring.

#### NOTE

• Turn on the monitor, plug in SpO<sub>2</sub> probe and connect patient's finger (or toe), monitor displays SpO<sub>2</sub> wave, "SpO<sub>2</sub> Pulse Search" displayed in the technical alarm area until the monitor measured SpO<sub>2</sub> value and pulse rate. "SpO<sub>2</sub> Search Timeout" displayed in the technical alarm area until the monitor measured pulse rate. Check the sensor mounting position, whether the sensor is damaged or sensor type. Reconnect the sensor or use new sensor.

# 10.4 Display

SpO2 parameter area is as shown in figure 10-1.




Figure 10-1: SpO<sub>2</sub> parameter display

SpO2 wave is as shown in figure 10-2.




Figure 10-2: SpO<sub>2</sub> wave

# 10.5 Setting SpO<sub>2</sub>

Select SpO<sub>2</sub> parameter area or SpO<sub>2</sub> wave area  $\rightarrow$  "SpO<sub>2</sub> Setup" menu, which is shown below. You can set SpO<sub>2</sub> through "SpO<sub>2</sub> Setup" menu.




Figure 10-3: "SpO2 Setup" menu

#### 10.5.1 Setting wave speed

➤ Select "Wave Speed" and set wave speed to "12.5mm/s" or "25mm/s"; the faster speed, the smoother wave.

#### 10.5.2 Setting wave mode

- > Select "Wave Mode", and set the wave drawing mode
  - "Scan": Scan mode.
  - "Fill": Fill mode.

#### 10.5.3 Setting average time

Select "Average Time", and set the average time to "2-4s", "4-8s", "8-16s".

#### 10.5.4 Pulse volume

The user can set the pulse volume. The pulse volume can be set to 0, 1, 2, 3, 4, 5, 6, 7, 8, or 9. By default, the pulse volume is set to 3.

# 10.6 Measuring influencing factors

During operation, the following factors can affect the accuracy of SpO<sub>2</sub> measurement:

- ➤ High-frequency radio wave interference, such as interference generated by the host system or interference from electrosurgery instrument connected to the system.
- Intravenous dye.
- > Too frequent movement of the patient.
- > External light radiation.
- > Sensor is improperly installed or improperly in contact with the patient.
- Sensor temperature.
- > The sensor is placed on limbs with blood pressure cuff, arterial duct or lumen tube.
- Concentration of non-functional hemoglobin such as carboxyhemoglobin (COHb) and methemoglobin (MetHb).
- Shock, anemia, hypothermia, and the application of vasoconstrictor drugs may reduce the arterial blood flow to a level that cannot be measured.
- ➤ The measurement also depends on the absorption of specific wavelengths of light by oxyhemoglobin and reduced hemoglobin. If there is any other substance that absorbs the same wavelength, the measurement may have false or low SpO<sub>2</sub> values, such as: carbon hemoglobin, methemoglobin, methylene blue, and indigo carmine.
- SpO<sub>2</sub> probe described in Annex is recommended.
- ➤ Operating environment limit: Operating temperature range: 5 ~ 40 °C, Humidity range:30%~85% (non-condensing) Atmospheric pressure: 700hPa ~ 1060hPa.

# 10.7 Alarm setup

In "SpO<sub>2</sub> Setup" menu, select "Alarm Setup >>" to enter "Alarm Setup" interface, and set SpO<sub>2</sub> alarm switch, alarm level, upper and lower alarm limit. See 6.5 Alarm Setup for detailed setting method.

# 10.8 Technical description

- Accessories have passed the biocompatibility test and meet the requirements of ISO 10993-1.
- Fluke's index 2XL Oxygen Analyzer can be used to check the function of the monitor and can be used to assess the accuracy of the pulse rate but cannot be used to assess the accuracy of blood oxygen.
- Functional testers cannot be used to evaluate the accuracy of pulse oximetry probe and pulse oximetry.
- Measure the maximum temperature between the oxygen probe and the tissue contact surface: Measured as described in Annex BB of ISO 80601-2-61, the temperature is less than 41 °C.

# **11** NIBP

### 11.1 Overview

The monitor uses oscillometric method to measure noninvasive blood pressure (NIBP).

The oscillometric method for measuring blood pressure is to inflate a cuff with a certain amount of pressure until the arterial blood flow has been completely blocked. As applied pressure decreases, the arterial blood flow which was completely occluded gradually opened, and completely opened. Then, the pulsation of the arterial vascular wall will generate a shock wave in the cuff. SBP, MAP, and DBP are obtained by measuring and analyzing cuff pressure oscillations when deflating.

- Produce first most clear signal reflect SBP
- Oscillation amplitude reaches the peak reflect MAP
- When the cuff pressure is suddenly lowered reflect DBP

Measuring mode: manual, cycle, and continuous. Each mode shows systolic, mean and diastolic blood pressure.

#### Manual mode

Using Manual mode start to measures by hand

#### Automatic mode measures

Use manual mode to open automatic mode, then the measure will automatically turn to automatic mode after a certain time. During measurement, any error will stop the current automatic measurement, but not affect next automatic measurement unless the time interval less than 30s. If the time interval less than 30s, should delay the next automatic measurement, keep the interval more than 30s.

The time interval can be choose In Automatic mode as 1, 2, 3, 4, 5, 10, 15, 30, 60, 90, 120, 180, 240, 480 minutes.

#### Continuous mode

Choose continuous mode, 5 seconds after complete a measurement start the next measurement, continue 5 minutes then stop. During measurement, any error will stop the continuous measurement. If the first measurement time is over 4 minutes and 40 seconds but less than 5 minutes, the continuous mode will stop before 5 minutes, if the first measurement time is over 5 minutes, the continuous mode will stop after 5 minutes.

The monitor is defibrillation proof, so the monitor operates normally after defibrillation.

# 11.2 Safety information

#### WARNING

- Do not carry out non-invasive blood pressure measurement on patients with sickle cell disease and skin damage or any expected damage. Do not measure NIBP on traumatic body part. This may cause further injury.
- When pediatric patients are measured, in order to ensure the cuff pressure does not exceed its maximum measurement range of patient types (Adult mode: 300mmHg and Pediatric mode: 240mmHg), you must ensure that you have selected the correct patient type (see patient information menu settings). Using the wrong type of pattern is likely to endanger the patient to patient safety, as higher blood pressure levels for adults does not apply to pediatric.
- For patients with severe coagulation disorder, determine if the automatic blood pressure measurement is carried out according to the clinical evaluation, since the friction of body and cuff may produce hematoma.
- Do not install a cuff on the limbs with intravenous infusion or duct, because it may lead to tissue damage around the duct when the cuff is inflated and makes the infusion slow down or be blocked.
- The inflatable tube connecting the blood pressure cuff and the monitor should be smooth without entanglement. The pressure generated by being kinked connection tubing may cause blood flow interference.
- For patients with severe thrombotic disorders, determine whether to carry out automatic blood pressure measurement according to the clinical situations, since the limb bundled with a cuff may produce hematoma.
- Measure blood pressure frequently will affect the distribution of blood flow, May endanger the safety of patients.
- Check the patient's physiological condition before measure blood pressure, in order to ensure that long time measure will not damage the circulation of patients
- For mastectomy patients, applying the NIBP cuff on the surgery side arm can cause lymphedema. Measure blood pressure on opposite side arm.
- Pressurization of the cuff can temporarily cause loss of function of simultaneously used monitoring device on the same limb.
- Measurement results may be affected by posture and mental state of the patient.
- If there are doubts on the measurement results, please use other blood pressure measurements and compare, if necessary, contact the Equipment Division.

# 11.3 Measurement limits

According to the patient's condition, the oscillometric method has some limitations. This measurement is to look for the regular pulse waves generated by arterial pressure. If the patient's condition makes this detection method difficult, the measured value becomes unreliable, and pressure measurement time increases. The user should be aware that the following conditions may interfere with measurement method, making the pressure measurement unreliable or extend the time. In this case, the patient's condition does not allow measurement.

#### Patient movement

If the patient is talking, moving, shaking or cramping, the measurement will be unreliable or even impossible, as these may interfere with the detection of arterial pressure pulse, and extend the pressure measurement time.

#### Arrhythmia

If the patient shows arrhythmia which results in irregular heartbeat, the measurement will be unreliable and even cannot be done, and the pressure measurement time will be extended.

#### Use of an artificial heart-lung machine

If a patient is connected to an artificial heart-lung machine, the measurement will be impossible.

#### Pressure changes

If the arterial pressure pulse is being analyzed to obtain a measured value at a certain time and the blood pressure of the patient changes rapidly, the measurement will be unreliable or impossible.

#### Severe shock

If the patient is in severe shock or hypothermia, the pressure measurement will not be reliable, because the decrease of blood flow to the periphery would cause decrease in arterial pulsation.

#### Limit heart rate

If the heart rate is below 40bpm (beats / min) or above 240bpm (beats / min), the blood pressure measurement is impossible.

#### Obese patients

A thick layer of fat around a limb blocks the arterial oscillation so that it cannot reach the cuff. The accuracy is lower than normal.

#### Environmental Requirements

Measuring blood pressure should meet the environment range as follow:

ambient humidity 30% ~ 85%, no condensing,

ambient temperature 5  $\sim$  40  $^{\circ}$ C,

Atmospheric pressure: 700hPa ~ 1060hPa.

NIBP performance and measurement accuracy will be affected beyond the range.

# 11.4 Measurement procedure

### 11.4.1 Prepare the measurement

- 1. Turn on the monitor, and check if it works properly.
- 2. Verify the patient category, and make changes if improper. Depending on the current patient type, the patient type is selected in the patient information interface.
- 3. Connect the blood pressure cuff extension tube to the monitor.
- 4. Select the cuff in accordance with the following method, make sure that the cuff is completely deflated, and then tie it to the upper arm or thigh of the patient.
  - Determine the limb circumference of the patient.
  - ➤ Select the appropriate cuff (marked with appropriate limb circumference). Cuff width should be 40% of the limb circumference or 2/3 of the upper arm length. The length of the inflated part of the cuff should be sufficient for 50%~80% around the limb.
  - Place the cuff on the upper arm or thigh of the patient, and ensure that the marking φ is located just above the appropriate artery. Make sure that the cuff does not wrap too tight around the limb, or it may cause distal discoloration or even ischemia.

#### 11.4.2 Patient posture requirements during measurement

- Sit comfortable or lie down relaxedly.
- 2. No crossing legs.
- 3. Back and elbow should be supported.
- 4. The center of NIBP cuff and the right atrium are at in the same level.
- 5. Remind patients, no talking during measurement and try to relax.

#### NOTE

- When have doubt about blood pressure measuring result, re-measure after the patient sit-in about 5 minutes. If still have doubt, replace the blood pressure measuring equipment and measure again.
- The operator should be in the position where he/she can readily operate the sphygmomanometer.

### 11.4.3 Start/stop measurement

Use the [NIBP start/stop] key on the monitor panel or [NIBP] smart hotkey on the display to start / stop the blood pressure measurement.

### 11.4.4 Correcting measurement results

The position of limb blood pressure measurement should be in the same horizontal position of the patient's heart. Otherwise, correct the measurement results with the following correction method.

- ➤ If the cuff is above the heart level position, increase 0.75mmHg (0.10kPa) per centimeter of gap to the measured results.
- ➤ If the cuff is below the heart level position, subtract 0.75mmHg (0.10kPa) per centimeter of gap from the measured results.
- ➢ If the patient is obese or clothes are too thick, subtract 5mmHg ~ 10mmHg (0.65kPa ~ 1.3kPa) from the measured results.

#### 11.5 NIBP display

NIBP measurement has no waveform display, and only displays NIBP measurement results in the parameter area, as shown in Fig. 11-1. The figure below is for reference only. The graphics displayed on the monitor may be slightly different.




Figure 11-1: NIBP parameter display

## 11.6 Setting inflation pressure

If necessary, you can manually set the initial cuff inflation pressure as follows.

- ➤ Select the NIBP parameter area → "NIBP Setup" menu;
- ➤ Select "Initial Pressure", and set the appropriate cuff pressure value. When the patient is adult, the pressure can be select from "140", "160", "180". The default cuff pressure value is "160".
- Select "Initial pressure", and set the appropriate cuff pressure value. When the patient is pediatric, the pressure can be select from "140", "160". The default cuff pressure value is "140".

## 11.7 NIBP reset

Select NIBP parameter area  $\rightarrow$  "NIBP Setup" menu  $\rightarrow$  Select "Reset", and restore the inflation pressure of the blood pressure pump to currently configured initial settings. When the blood pressure pump is not working properly, but no warning is given, you can reset the blood pressure pump, and automatically restores the blood pressure pump.

### 11.8 Clean and disinfection method of NIBP cuff

If necessary, NIBP cuff and NIBP extension tube can be cleaned and disinfected together without separated

#### 11.8.1 Cleaning method

- Prepare enzyme cleaning agent, distilled water and 10% solvent, respectively in different spray bottle.
- 2. Sprinkle cleaning agent on NIBP cuff, connector and extension tube, keep 1 minute for the dry stains.
- 3. Use a soft cloth to wipe smooth face. Use soft hair brush to brush visible stain and irregular surface
- 4. Rinsed with copious amounts of distilled water.

#### NOTE

- Please be especially careful to clean the air ball and control valve of whole air system. Do not allow any liquid entering into reversing valve and saturated valve.
- Don't use a soft cotton ball and fiber to clean this accessory because they will stick on the cuff and extension tube.

#### 11.8.2 Disinfection method

- 1. Sprinkle bleach solution (Formula: the proportion of water and bleaching powder to 1:10) then keep 5 minutes
- 2. Wipe off excess bleach solution and elute with distilled water again
- 3. Natural dry cuff

## 11.9 Alarm setup

In "NIBP Setup" menu, select "Alarm Setup >>" to enter "Alarm Setup" interface, and set NIBP alarm switch, alarm level, upper and lower alarm limit. See 6.5 Alarm Setup for detailed setting method.

# **12 TEMP**

#### 12.1 Overview

The monitor has two temperature measurement channels; the temperature sensor will measure the body temperature, and calculate the difference between the body temperature data.

The monitor is designed for defibrillation proof, so the monitor operates normally after defibrillation.

### 12.2 Safety information



#### WARNING

- Before monitoring, check if the probe cable is normal. Unplug the temperature probe cable from the jack, the screen will display "TEMP1 Sensor Off" and "TEMP2 Sensor Off" prompt and make an alarm sound.
- Calibrate the temperature measuring instrument at least once every two years (or according to hospital procedures). When calibration is required, please contact Bistos.

# 12.3 Measurement steps

Please refer to the following steps:

- 1. Turn on the monitor and check if it works normally.
- Select the appropriate temperature probe according to the patient category and measurement needs.
- Insert the probe lead wire into the temperature probe interface.
- Attach the probe to the patient properly.
- 5. Make sure that the alarm settings apply to the patient.

When measuring body temperature, temperature probe can be attached to body surface such as the neck, armpits, ears and other locations.

#### 12.4 Measurement requirements

The normal measuring range of monitor is 0°C~50°C, and the accuracy is consistent in this range.

The environmental temperature range for body temperature measuring is  $5^{\circ}C^{\sim}40^{\circ}C$ . Get the right temperatures for the shortest measurement time is 40s, and the measuring interval is 1s.



#### WARNING

Please measure the body temperature in the specified environment temperature range, or else it may be dangerous.

#### 12.5 Temperature display

The monitor can display the body temperature of two channels (T1 and T2) and the alarm limits, difference between the two temperature (TD) and temperature units. Select Temp parameter area and open the [Temp Setup] menu.

Temperature display area is as shown below:




Figure 12-1: TEMP parameter display

#### WARNING

The operator, prior to use, need to check the compatibility of the probe and thermometer. If the temperature value displayed by the monitor has significant difference from the body temperature under normal condition, please check if the probe resistance of the monitor matches the resistance set in the monitor system; if not, please replace a probe with appropriate resistance or adjust the monitor and select the appropriate resistance. Incompatible probe will affect the critical properties.

## 12.6 Setting temperature unit

-You can define your favorite temperature unit as follows:

Select TEMP parameter area → "TEMP Setup" menu.

In the "TEMP Setup" menu, set "Unit" to "  $^{\circ}C$  " or "  $^{\circ}F$  "

# 12.7 Alarm setup

In "TEMP Setup" menu, select "Alarm Setup >>" to enter "Alarm Setup" interface, and set TEMP alarm switch, alarm level, upper and lower alarm limit. See 6.5 Alarm Setup for detailed setting method.

## 12.8 Technical description

Accessories have passed the biocompatibility test and meet the requirements of ISO 10993-1.

# 13 Review

The monitor provides up to 168 hours trend data review of all monitoring parameters, 1000 groups of NIBP measurement data and 200 physiological alarm events, 100 technical alarm events. The user can select trend chart or trend table to view trend change; or view the latest wave.

## 13.1 Reviewing trend chart

Select [Trend] smart hotkey to enter [Trend] menu, and select [Graphic] to enter the following window.




Figure 13-1: Trend chart

- In the trend chart, use the following method to select the parameter to be reviewed:
  - Select the parameter box, rotate the shuttle to select the parameters to be reviewed, press the shuttle, and set the parameter box as the parameter to be reviewed.
- > Browse the trend chart in the following method:
  - Select and to move the trend cursor.
  - Select and b to turn pages to left or right and move the trend chart.
  - The cursor top displays the current time corresponding to the current cursor position, and the left of the trend chart window displays the parameter values of the time, which will change automatically with the move of trend cursor.

# 13.2 Reviewing trend table

Select [Trend] smart hotkey to enter "Trend" menu, select "Tabular" and enter the following window.



Figure 13-2: "Trend" table

- > Browse the trend table in the following method:
  - Select and to turn pages to left or right and move the trend table to observe the target parameters.
  - Select and to turn pages up or down and move the trend table to observe more data.

## 13.3 NIBP measurement review

Select Trend smart hotkey to enter "Trend" menu, and select "NIBP" to enter the following window




Figure 13-3: NIBP measurement review

This window shows the measurement time of noninvasive blood pressure, systolic blood pressure "SYS", diastolic blood pressure "DIA", mean blood pressure "MAP" and pulse rate "PR". The monitor can store 1000 sets of NIBP measurements in total.

- > NIBP viewing method is as follows:
  - Select and to turn pages up or down and move the trend table to observe more data.

# 14 Battery

#### 14.1 Overview

The monitor has a built-in rechargeable battery to ensure that the monitor can also be used normally in case of patient transfer or power failure. When the monitor is connected to an DC power source, it will charge the battery no matter whether the monitor is turned on or not. In the case of power failure, the system will automatically use the battery to power the monitor to avoid interrupting the monitor working.

The battery icon on the screen indicates the battery status:



Battery is working properly and is fully charged.



Battery is working properly and the green part indicates the battery power.



Battery power is low, and requires charging immediately, or else the pulse oximeter will turn off automatically.



Battery is not installed.



Battery is properly installed and being charged.

The battery power can only maintain for some time. Low battery voltage will trigger a high level technical alarm "Battery Low"; in this case, connect the monitor to DC power and charge the battery.

#### 14.2 Battery usage guide

Battery life depends on the frequency and time of use. If the battery maintenance and storage are proper, the lithium battery life is three years. If you do not use the battery properly, its life may be shortened. It is recommended to replace the lithium battery once every three years.

In order to ensure the maximum capacity of the battery, please note the following usage guide:

- Do not drop the battery.
- Check the battery performance once every two years. Before servicing the monitor or you suspect that the battery is the fault source, also check the battery performance.

#### WARNING

- Keep the battery out of the reach of children.
- Use only the designated battery.
- If the battery is damaged or leaks, replace it immediately. Do not use a defective battery for the monitor.
- Do not open the battery compartment. Only the qualified service personnel authorized by the manufacturer can open the battery compartment and replace or check the battery, and when it needs replacement, the replacement should be same model with established battery.
- Do not disassemble the battery, throw it in fire, or short-circuit it. Battery fire, explosion and leakage may lead to personal injury. Do not touch the leaking battery with bare hand directly.

## 14.3 Checking battery performance

Please refer to the following steps to check the battery performance:

- Disconnect the monitor from the patient and stop all monitoring or measurement.
- Connect DC power to the monitor, and charge the battery for more than 4 hours uninterruptedly.
- Disconnect the DC power and power the monitor with battery until the monitor is turned off.
- Battery duration reflects the battery performance.

If the battery operating time is significantly shorter than the time stated in specifications, please contact our service personnel for replacing the battery.



#### WARNING

Do not open the battery compartment. Only the qualified service personnel authorized by the manufacturer can open the battery compartment and replace or check the battery, and when it needs replacement, the replacement should be same model with established battery.

# 14.4 Battery recycling

If the battery has visible damage or cannot store power, it should be replaced and recycled properly. Follow the appropriate regulations to dispose of used batteries.



#### WARNING

Do not disassemble the battery, throw it in fire, or short-circuit it. Battery fire, explosion and leakage may lead to personal injury. Do not touch the leaking battery with bare hand directly.

# 15 Caring and cleaning

#### 15.1 Overview

In the using process, please make sure that there is no dust on or near your device. To prevent damage, please use the diluted detergents and disinfectants specified in this Manual, and use the lowest possible concentration. For the damage or accident caused by using other materials or methods, our company does not assume any responsibility.

#### 15.2 Cleaning

The device should be cleaned regularly. In the heavily polluted environment, increase the frequency of cleaning. Before cleaning, please consult the hospital about device cleaning requirements.

Below are available cleaning agents:

- Diluted ammonia
- Diluted sodium hypochlorite (washing bleach)
- Diluted formaldehyde
- Hydrogen peroxide (3%)
- Ethanol (70%)
- Isopropanol (70%)

#### Before cleaning:

- Turn off the monitor and disconnect the power.
- Use a soft cotton ball to adsorb appropriate amount of cleaning agent and wipe the display screen.
- Use a soft lint-free cloth to adsorb appropriate amount of cleaning agent and wipe the surface of the device.
- If necessary, use a clean, dry, lint-free cloth to remove any excess detergent.
- Dry the device naturally in a ventilated cool environment.

#### WARNING

Before cleaning the monitor or sensor, turn off the power and disconnect the DC power.

The monitor should be kept clean. It is recommended to regularly clean the
enclosure surface and the display screen. Cleaning the enclosure with nonetching cleaner such as soap and water.

#### **A**CAUTION

- To avoid damaging the monitor:
  - Do not use strong solvents such as acetone.
  - Most cleaners must be diluted before use. Diluting should be according to the manufacturer's instructions.
  - Do not use abrasive materials (such as steel wool).
  - > Do not allow any liquid entering into the enclosure, and never immerse any part of the device into liquid.
  - Do not leave any cleaning solution on the surface of any part of the device.

#### NOTE

- Wipe the monitor and sensor surface with medical alcohol, dry it naturally or with clean, dry, lint-free cloth.
- Bistos is not liable for effectiveness of using these chemicals for infectious disease control. Please consult the infectious disease control officers or experts of the hospital for advice.

## 15.3 Disinfection

In order to avoid damage to the product, we recommend that the product is disinfected only when it is deemed necessary by the hospital maintenance procedures. We also recommend that the instrument to be disinfected must first be cleaned.

## $\triangle$ CAUTION

• To prevent damage to the monitor, do not disinfect the monitor with gas (EtO) or formaldehyde.

# 16 Maintenance

#### $\wedge$

#### WARNING

 If the hospitals or institutions using this instrument can't implement a satisfactory maintenance schedule, it will result in device failure and may endanger human health.

## 16.1 Checking

Check the following basic items before each using the monitor:

- Check for any mechanical damage.
- Check all exposed wires, insertions and accessories.
- Check all instrument functions that may be used for patient monitoring and ensure that the instrument is in good working condition.

If the instrument function has any sign of damage, do not use this monitor for any patient monitoring. Please contact the hospital's professional maintenance personnel or our customer service personnel.

Every 6-12 months or after each repair, a comprehensive examination must be performed by trained and qualified technical service personnel, including functional safety checks; the specific inspection items are as follows:

- Environment and power meet the requirements.
- Device and accessories have no mechanical damage.
- The power supply has no wear, and the insulation is good.
- Specified accessories are used.
- Alarm system is functioning correctly.
- > Battery performance meets the requirements.
- Monitoring functions are in good working condition.
- Ground impedance and leakage current meet the requirements.

If the instrument function has any sign of damage, do not use this monitor for any patient monitoring. Please contact the hospital's professional maintenance personnel or our customer service personnel.

All checks that require disassembling the instrument must be performed by qualified service personnel. Safety and maintenance checks may also be carried out by the Company's personnel.

# 16.2 Viewing software version information

You can view the software version through the following steps:

- Select Settings Smart Hotkey → "Settings" Menu;
- ➤ Select "Monitor Info>>" → "Monitor Info" menu;
- > "Monitor Info" menu displays the software version information of the monitor.

## 16.3 Maintenance plan

The following tasks can only be done by qualified service personnel of Bistos. When the following maintenance is needed, please contact your service representative. Before testing or maintenance, clean and disinfect the device.

Inspection / Maintenance Item	Frequency
Check the safety according to IEC	At least once every two years, after replacing the
60601-1	power supply or the monitor falls down.
Check all monitoring or measuring	At least once every two years, or when you suspect
functions not listed	that the measured value is not accurate.
NIDD lookogo tost	At least once every two years, or follow hospital
NIBP leakage test	regulations
NIBP calibration	At least once every two years, or follow hospital
NIBP CAIIDIALION	regulations

#### 16.4 ECG calibration

In the process of using the monitor, the displayed ECG signals may be inaccurate due to hardware or software problems, mainly shown as waveform amplitude becoming larger or smaller. At this moment, you need to calibrate ECG.

Prepare the following instruments before testing:

- ECG simulator
- ECG cable
- Vernier caliper

The calibration method is as follows:

- > Connect the ECG cable to the monitor.
- Connect the ECG electrodes to the ECG simulator.
- Select Settings Smart Hotkey → "Settings" Menu;
- Select "User Maintenance >>" →enter the password and confirm → "User Maintenance" menu.
- Select "Module Maintenance >>" → "Module Maintenance" menu.
- Select "ECG"  $\rightarrow$  "ECG Maintenance" menu, and select "Calibration" to calibrate the ECG.
- Measure the wave amplitude with a caliper; in different filtering modes,  $\times 0.25$  is  $2.5 \pm 5\%$  (mm),  $\times 0.5$  is  $5.0 \pm \%$  5 (mm),  $\times 1$  is  $10.0 \pm \%$  5 (mm), and  $\times 2$  is  $20.0 \pm \%$  5 (mm). Comparing the amplitude of the square wave with the ruler, the error range should be within 5%.
- When calibration is complete, select "Stop Calibration" to exit.

# 17 Accessories

#### WARNING

- Use the accessories specified in this manual. Using other accessories may damage the monitor, or cannot reach the safety and performance claimed in this manual.
- The operating and storage environment of the monitor should meet the requirements of the accessories. Please refer to the manual of the accessories for these requirements.
- Disposable accessories can only be used once, because repeated use can cause performance degradation.
- If the packaging or accessories have any sign of damage, do not use such accessories.
- For ECG Cables, SpO<sub>2</sub> Sensor, Blood Pressure Cuff and Temperature Sensor, the normal life time is two years. Please replace in time.

#### Standard accessories are as follows:

No.	Description	QTY	Type-number
1	ECG Cables and lead-wires ECG electrodes(5)	1	Manufacturer: Shenzhen Launch Electronics Tech CO., Ltd  98ME01AC009(AHA standard) or 98ME01EC009(IEC standard)
2	Adult Finger Clip SpO <sub>2</sub> Sensor	1	Manufacturer: Unimed Medical Supplies,Inc
3	SpO2 extension cable		U403-01
4	Adult Non-Invasive blood pressure cuff	1	Manufacturer: Shenzhen Med-link Electronics Tech Co.,Ltd Y000A1
5	NIBP extension tube	1	Manufacturer: XIAMEN CONJOIN ELECTRONICS TECHNOLOGY CO., LTD CJP37-C12B1
6	Temperature Sensor	1	Manufacturer: Shenzhen taijia electronic Co., Ltd SPT4520010N
7	Power Adapter	1	Manufacturer: DONGGUAN SHILONG GUHUA ELECTRONIC CO., LTD UE36LCP1-150240SPA

# 18 **Specifications**

# 18.1 Safety specifications

## 18.1.1 Product category

In accordance with classification specified in the European Medical Device Directive 93/42/EEC, this monitor is Class IIb device. The monitor is classified as follows in accordance with IEC 60601-1:

Category Name	Specification
Type of electric shock protection	Class II and internally powered equipment When you question the integrity of the external protective earthing or protective ground conductor parameter of the equipment, the device must be powered by the internal power supply (battery).
Electric shock protection grade	Type CF applied part (defibrillation proof)
Explosion protection grade	Common equipment, no explosion protection
Liquid inlet protection grade	IPX1
Operating mode	Continuous mode
Movement	Portable equipment

#### 18.1.2 Power

Power	
Adapter	Input: AC 100 - 240V (50/60 Hz)
	Output: DC 15V / 2.4A
Rechargeable Battery	11.1V Li-ion battery 4400 mA
	Operating Time(When it fully charged): 5 hours
	Charging Time(Fully): 4 hours

# 18.2 Hardware specifications

<b>Physical Characteristics</b>	
Dimensions	Main Unit: 320(W) X 250(H) X 65(D)
Weight	<= 2.8 Kg for standard configuration

Display	
Туре	Color TFT touch screen LCD
Size and resolution	12", 800 X 600 pixels

Audio	
Speaker	Alarm sound (45 ~ 85 dB), key pressing sound
	QRS sound, PR sound
	Alarm sound meet the IEC 60601-1-8 standard requirements

Alarm signal	
Alarm delay	Off, 1s, 2s, 3s, 4s, 5s, 6s, 7s, 8s, depending on the setup Default is 4s
Pause duration	1min, 2min, 3min, 4min, 5min, 10min, 15min or permanent, depending on the setup Default is 2 minutes.

Data storage	
Trend	168 hours. Resolution: 1 min
Alarm event	200 physiological alarm events, 100 technical alarm events
NIBP measurement result	1000 groups

Environment		
	Operation	Transport and storage
Temperature	5~ 40°C (41°F~104°F)	-20 ~ 60°C (-4°F~140°F)
Humidity	30~ 85% non-condensing	0 ~ 95 % non-condensing
Atmospheric pressure	70~106 kPa	70~106 kPa

# 18.3 Functional specifications

# 18.3.1 ECG/TEMP/RESP

ECG			
Lead Type	3 lead	1, 11, 111	
	5 lead	I, II, III, aVR, aVL, aVF, V	
Disales es está de e	Auto, 2.5mm/mV(x0.2	Auto, 2.5mm/mV(x0.25), 5 mm/mV(x0.5),	
Display sensitivity	10mm/mV(x1.0), 20mi	m/mV(x2.0), 40mm/mV(x4.0)	
Wave sweep speed	12.5 mm/s, 25 mm/s, 5	50 mm/s	
	Diagnostic mode	0.05 - 100 Hz	
Band width	Monitor mode	0.5 - 40 Hz	
	Surgery mode	1 - 25 Hz	
CMRR	>100 dB		
Notch	50/60 Hz notch filter can be set to on or off		
Differential input	> 5 MΩ		
impedance			
Electrode polarization	± 400 mV		
voltage range			
Baseline recovery time	<5s after defibrillation (in monitor and surgery mode)		
Calibration signal	1 mV (peak – peak), accuracy ± 3%		
Lead-off detection	Measuring electrode: < 0.1 uA		
current	Drive electrode: < 1uA		

Pacing pulse	
	For PACE MAKER pulses that meet the criteria below, pacing pulse will be marked on the screen.
Pulse identification	Detection range(Amplitude): ± 2 mV ~ ± 700 mV
	Pulse width: 0.2ms ~ 2.0 ms
Average HR	Calculate from 15s data
Interval of HR refreshing	Calculate once every second
HR change response time	Time from 80 bpm to 120 bpm: ≤ 10 sec
nk change response time	Time from 80 bpm to 40 bpm: ≤ 10 sec
	For T-wave with 100ms QRS wave, 350ms QT period, 180ms
Tall T-wave suppression	duration and 1.2mV amplitude, the HR calculation will not be
	affected
Without overshoot	Amplitudes (ap) from ±2 mV to ±700 mV and pulse widths
rejection of pacemaker	from 0.1 ms to 2.0 ms.
pulses	
Tall T-wave rejection	2mV
capability	

HR		
Measuring range	Adult: 15 ~ 300 bpm Pediatric: 15 ~ 350 bpm	
Resolution	1 bpm	
Heart rate measurement error	± 1 bpm or ± 1%, whichever is greater	
Hoost water manageming	Ventricular bigeminy	80 ± 1 bpm
Heart rate measuring	Slow alternating ventricular bigeminy	60 ± 1 bpm
accuracy and response to	Rapid alternating ventricular bigeminy	120 ± 1 bpm
irregular rhythm	Bidirectional systoles	90 ± 2 bpm
	1 mV, 206 bpm Ventricular tachycardia	<10 s
	0.5 mV, 206 bpm Ventricular tachycardia	<10 s
Time to alarm for tachycardia	2 mV, 206 bpm Ventricular tachycardia	<10 s
	2 mV, 195 bpm Ventricular tachycardia	<5 s
	1 mV, 195 bpm Ventricular tachycardia	<5 s
	4 mV, 195 bpm Ventricular tachycardia	<5 s

HR Alarm	
HR upper limit	Adult: 16 ~ 300, 1 bpm step
nk upper iiiiiit	Pediatric: 16 ~ 350, 1 bpm step
HR lower limit	Adult: 15 ~ 299, 1 bpm step
HK lower little	Pediatric: 15 ~ 349, 1 bpm step

TEMP	
Standard compliance	ISO 80601-2-56:2009
Measurement method	Thermistor
Operating mode	Direct mode
Measuring range	0 °C ~ 50.0 °C (32 °F ~ 122.0 °F)
Resolution	0.1 ℃
Measurement accuracy	± 0.3 °C
Number of channel	2

TEMP Alarm	
T1/T2 upper limit	0.1 °C ~ 50.0 °C , 0.1 °C /°F step
T1/T2 lower limit	0 °C ~ 49.9 °C, 0.1 °C /°F step
TD upper limit	0 °C ~ 50.0 °C, 0.1 °C/°F step

RESP	
Measurement method	Thoracic electrical bio impedance method
Measuring range	Lead RA-LA, RA-LL,LA-RL,LL-RL
Wave gain	X0.5, x1, x2
Respiratory impedance	0.2 ~ 3 Ω
range	
Base line impedance	500 ~ 2 000 Ω
Scan speed	6.25 mm/s, 12.5 mm/s, 25 mm/s
Measurement accuracy	± 2 rpm
Measurement range	0 ~ 120 rpm

RR Alarm	
DD upper limit	Adult: 7 ~ 120
RR upper limit	Pediatric: 7 ~ 150
DD lower limit	Adult: 6 ~ 119
RR lower limit	Pediatric: 6 ~ 149

### 18.3.2 NIBP

NIBP					
Standards compliant	IEC 80601-2-30:2009/A1:2013				
Measurement method	Automatic oscillometric method				
Operating mode	Manual, au	tomatic, conti	nuo	us	
Useful life	100, 000 tir	nes			
Measurement interval in automatic mode	1/2/3/4/5/	10/15/30/60/9	90/1	20/180/240/480	min
Typical measurement time	20~40s				
				Adult	Pediatric
Normal mode measuring	Systolic blo	od pressure		40-270	40-200
range (mmHg)	Mean blood pressure			20-230	20-175
	Diastolic blo	ood pressure		10-210	10-162
Measurement accuracy	Maximum average error: ±5mmHg				
ivicasurement accuracy	Maximum standard deviation: 8mmHg				
Resolution	1mmHg				
		Default Pressure setting range			nge
Initial inflation pressure	Adult	160mmHg	14	10mmHg, 160mm	Hg, 180mmHg
	Pediatric 140mmHg 140mmHg, 160mmHg,			Hg,	
Overpressure protection	Adult: 300mmHg				
Overpressure protection point (software)	Pediatric: 240mmHg				
Overpressure protection	Adult: 320~330mmHg				
point (hardware)	Pediatric: 265~275mmHg				
Static Pressure accuracy	±3mmHg				

NIBP Alarm			
		Adult	Pediatric
NUDD upper limit (mm l/g)	SYS	31 ~ 280	31 ~ 230
NIBP upper limit (mmHg) 1 mmHg step	MAP	11 ~ 240	11 ~ 175
	DIA	11~220	11~165
NUDD love r limit (mmllg)	SYS	30 ~ 279	30 ~ 229
NIBP lower limit (mmHg) 1 mmHg step	MAP	10 ~ 239	10 ~ 174
	DIA	10~219	10 ~ 164

NIBP Electrical characteristics	
Supply voltage	10V~14V DC
Maximum power	3.6w
consumption	
Quiescent current	50mA
Maximum current during	180mA
measurement	
Maximum current during	300mA
inflation	

#### 18.3.3 SpO<sub>2</sub>

SpO <sub>2</sub>	
Standards compliant	ISO 80601-2-61:2011
	•••

#### Measurement accuracy verification

The  $SpO_2$  accuracy has been verified in human experiments by Comparing with arterial blood sample reference measured with a CO-oximeter. Pulse oximeter measurements are statistically distributed and about two-thirds of the measurements are expected to come within the specified accuracy range compared to CO- oximeter measurements.

The accuracy of the oximeter has been validated by a clinical trial involving 12 healthy adult subjects - 4 women and 8 men. Among them medium skin are 4 subjects, light skin are 5 subjects, dark skin are 3 subjects, the age from 21 to 28.

Overall accuracy was determined by calculating the root mean square error across all samples and is 1.56%".

Display range	0% ~ 100%		
SpO <sub>2</sub> display resolution	1%		
SaO <sub>2</sub> checking accuracy	±2% (70%~100%);		
	not define when lower thar	า 70% ;	
SnO2 alarm limit range	Upper alarm limit	1%~100%	
SpO2 alarm limit range	Lower alarm limit	0%~99%	
SpO₂ alerting signal	No delay		
generates a delay			
SpO <sub>2</sub> value refresh period	1s/time		
	Low sensitivity	6~8s	
Average period	Intermediate sensitivity	4 ~ 6s	
	Advanced sensitivity	2~4s	
Alama aanditian dalam	Low sensitivity	<8s	
Alarm condition delay period	Intermediate sensitivity	<6s	
period	Advanced sensitivity	<4s	
Alarm sign generates	Os		
delay period			

PR	
Measuring range	25~250bpm
Resolution	1% bpm
Accuracy	±2% or ±2bpm, whichever is greater

PR alarm		
Upper limit	Adult: 16 ~ 300	
Upper limit	Pediatric: 16 ~ 350	
Lougarlimit	Adult: 15 ~ 299	
Lower limit	Pediatric: 15 ~ 349	

# 19 Alarm information

This chapter lists some important physiological and technical alarm information, and some alarms are not necessarily listed.

Note that in this chapter: P column indicates the default alarm priority: H indicates high priority, M indicates middle priority, L indicates low priority, and "\*" indicates priority set by the user.

Corresponding countermeasures are listed for each alarm message. If you operate in accordance with the countermeasures but the problem persists, contact your service personnel.

### 19.1 Physiological alarms

Source	Alarm message	Р	Causes and countermeasures
	HR Too High		HR value is higher than the upper alarm limit or
			lower than the lower alarm limit. Check the
	IID Ta a Lavo	M*	patient's physiological condition, and check if the
	HR Too Low		patient category and alarm limit settings are
			appropriate for the patient.
			PVCs value is higher than the upper alarm limit or
			lower than the lower alarm limit. Check the
	PVCS Too High	M*	patient's physiological condition, and check if the
			patient category and alarm limit settings are
			appropriate for the patient.
	Asystole	Н	The patient has arrhythmia. Check the patient's
	VF/VTA	Н	condition, electrodes, cables and lead wires.
	R on T	M*	condition, electrodes, cables and lead wires.
ECG	Frequent PVC	M*	
100	Couplet PVC		
	Single PVC	M*	
	PVC Bigeminy	M*	
	PVC Trigeminy	M*	
	Tachycardia	M*	
	Bradycardia	M*	
	Miss Beat	M*	
	Pacemaker Not Capture	Н	Pacemaker works abnormally; check the
	Pacemaker Not work	Н	pacemaker.
		Н	The patient ECG signal is too weak, and the
	ECG Signal weak		system can't analyze. Check the patient's
			condition, electrodes, cables and leads.
	ST-I Too High	M*	ST value is higher than the upper alarm limit or
	ST-I Too Low	IVI "	

Source	Alarm message	Р	Causes and countermeasures
	ST-II Too High		lower than the lower alarm limit. Check the
	ST-II Too Low		patient's physiological condition, and check if the
	ST-III Too High		patient category and alarm limit settings are
	ST-III Too Low		appropriate for the patient.
	RR Too High		Patient PR value is higher than the upper alarm
			limit or lower than the lower alarm limit. Check
		M*	the patient's physiological condition, and check if
	RR Too Low		the patient category and alarm limit settings are
Resp			appropriate for the patient.
			The patient's respiratory signal is too weak, and
	Apnea(RESP)	Н	the system can't analyze. Check the patient's
	,		condition, electrodes, cables and leads.
	RESP ARTIFACT	H*	Respiration heartbeat interference
	T1 Too High		T1/T2 value is higher than the upper alarm limit
	T1 Too Low		or lower than the lower alarm limit. Check the
	T2 Too High		patient's physiological condition, and check if the
	T2 Too Low		patient category and alarm limit settings are
			appropriate for the patient.
Temp		− M*	TD value is higher than the upper alarm limit or
			lower than the lower alarm limit. Check the
	TD Too High		patient's physiological condition, and check if the
	TD 100 mgm		patient category and alarm limit settings are
			appropriate for the patient.
	SpO₂ Too High		SpO <sub>2</sub> value is higher than the upper alarm limit or
			lower than the lower alarm limit. Check the
			patient's physiological condition, and check if the
	SpO₂ Too Low		patient category and alarm limit settings are
			appropriate for the patient.
SpO <sub>2</sub>	PR Too High	M*	PR value is higher than the upper alarm limit or
			lower than the lower alarm limit. Check the
	DD To a La		patient's physiological condition, and check if the
	PR Too Low		patient category and alarm limit settings are
			appropriate for the patient.
	NIBP signal weak		NIBP value is higher than the upper alarm limit or
	NIBP-Sys Too High		lower than the lower alarm limit. Check the
NIBP	NIBP-Sys Too Low		patient's physiological condition, and check if the
INIDE	NIBP-Mean Too High	IVI	patient category and alarm limit settings are
	NIBP-Mean Too Low		appropriate for the patient.
	NIBP-Dia Too High		

# 19.2 Technical alarms

Source	Alarm message	Р	Causes and countermeasures	
			Connect to AC power supply, and charge the	
System	Battery Low	Н	battery, and power with the battery as needed	
			after fully charged.	
ECG	ECG Comm. Stop	Н	ECG module failure, or communication failure	
	ECG Comm. Error	Н	between the module and the host; please restart	
	ECG Config Error	Н	the device.	
	ECG Selfcheck Error	Н		
	ECG Lead Off	M*	The electrodes are not connected to the patier	
	ECG YY OFF (YY is a lead		firmly or fall off, or lead wires and the main cable	
	name)	M*	fall off. Check the connection of electrodes and	
			lead wires.	
Tomn	TEMP1 Sensor Off	L	The temperature sensor falls off from the patient.	
Temp	TEMP2 Sensor Off	L	Check the sensor connection.	
	SpO <sub>2</sub> Comm. Stop	Н	SpO₂ module failure, or communication failure	
	SpO <sub>2</sub> Comm. Error		between the module and the host; please restart	
		Н	the device.	
	SpO <sub>2</sub> No Sensor	L	SpO <sub>2</sub> sensor falls off from the patient or monitor,	
	SpO <sub>2</sub> Sensor Off		malfunctions, or sensor other than specified in	
		L	this Manual is used. Check the sensor mounting	
SpO <sub>2</sub>	SpO₂ Search Timeout		position, whether the sensor is damaged or	
- p - <u>z</u>	, -	L	sensor type. Reconnect the sensor or use new	
			sensor.	
	SpO₂ Search Pulse		Sensor signal is poor or too weak. Check the	
	• -		patient's condition, and place the sensor in a	
		L	suitable position. If the failure persists, replace	
			the sensor.	
	NIBP Comm. Stop	Н	NIBP module failure, or communication failure	
	NIBP Comm. Error	Н	between the module and the host; please restart	
	NIBP Selfcheck error	Н	the device.	
	NIBP CFG Error	Н	the device.	
	NIBP system error	Н	If failure occurs during measurement, the system	
NIBP			can't analyze and calculate. Check the patient's	
INIDI	Measurement timeout	L	condition, check the connections or replace the	
			cuff, and then re-test.	
			The used cuff does not match the set patient	
	Cuff type error	L	category. Verify the patient category and replace	
	Can type ciroi	-	the cuff.	
			NIBP cuff isn't placed or connected properly, or	
	Cuff loose or no cuff L		there is gas leak.	
AUDO	C (() - 1			
NIBP	Cuff leak	L	Check cuff and inflation tube.	
	Air pressure error	L	Ambient atmospheric pressure is abnormal.	
	,		Confirm that the environment complies with the	

Source	Alarm message	Р	Causes and countermeasures
			monitor's specifications, and check whether
			there are special reasons affecting ambient
			pressure.
	NIDD aver your	1	The measured blood pressure of the patient
	NIBP over range	L	exceeds the measuring range.
			Patient's pulse may be weak or cuff is too loose.
	NIDD signal week	١.	Check the condition of the patient, and place the
	NIBP signal weak	L	cuff in a suitable position. If the failure persists,
			replace the cuff.
	NIBP signal unstable		Excessive movement may result in too much
		L	motion artifact or interference in the signal
			during measurement.
	NIBP signal saturated	1.	Motion signal amplitude is too large due to
		L	movement and other reasons.
	NIBP over pressure		Cuff overpressure, and gas blockage may occur;
			check the gas path, and then re-measure.
	Modulo reset failed		NIBP module reset error; check the gas path is
	Module reset failed	L	blocked, and then restart the measurement.

# 20 Default parameter configuration

This chapter lists the important factory default settings of different departments in monitor configuration mode. Users can not change the default configuration, but can modify the settings as required and save as user-defined configuration.

	required and so	Oution			lts
Module	Option			Adult	Pediatric
	Alarm level			Mid	Mid
	Alarm record			Off	Off
	Lead type			5-lead	5-lead
	Calculation of	channel		Auto	Auto
	Power frequ	ency supp	ression	On	On
	Alarm limits			50~120 on	75~160 on
		ST segme	ent analysis	Off	Off
	ST	Alarm le	vel	Mid	Mid
	segment	Alarm re	cord	Off	Off
	analysis	Alarm lin	nits	-0.2~0.2 on	-0.2~0.2 on
ECG		Alarm le	vel	Mid	Mid
		Alarm re	cord	Off	Off
	Arrhythmia	Alarm lin	nits	0~10 on	0~10 on
	analysis	ARR	Alarm switch	On	On
		alarm	Alarm level	Mid	Mid
		settings	Alarm record	Off	Off
	Gain			x1	x1
	Wave velocit	ty		25.0mm/s	25.0mm/s
	Filter mode			Monitor	Monitor
	Wave color			Green	Green
	Wave style			Color scale	Color scale
	Alarm level Alarm record		Mid	Mid	
				Off	Off
	Pressure unit			mmHg	mmHg
	Measurement mode			Adult	Pediatric
NIBP	Interval			Manual	Manual
NIDP	Display color	r		White	White
	Pre-inflation value			150	100
	Systolic bloo	d pressure	limit	90~160 on	70~120 on
	Mean blood	pressure l	imit	60~110 on	50~90 on
	Diastolic blo	od pressui	re limit	50~90 on	40~70 on
	Alarm level			Mid	Mid
	Alarm record	t		Off	Off
520	Alarm limits			90~100 on	90~100 on
SpO <sub>2</sub>	Wave velocit	ty		25.0	25.0
	Wave color			Cyan	Cyan
	Wave style			Line	Line
	Alarm level			Mid	Mid
RESP	Alarm record	t		Off	Off
	Apnea alarm	)		20 sec	20 sec

Module	Oution	Module defaults		
	Option	Adult	Pediatric	
	Alarm limits	8~30 on	8~30 on	
	Gain	x1	x1	
	Wave velocity	12.5	12.5	
	Wave color	Yellow	Yellow	
	Wave style	Line	Line	
	Alarm source	SpO <sub>2</sub>	SpO <sub>2</sub>	
PR	Alarm level	Mid	Mid	
PK	Alarm record	Off	Off	
	Alarm limits	50~120 on	75~160 on	
	Alarm level	Mid	Mid	
	Alarm record	Off	Off	
	Display color	White	White	
TEMP	Temperature unit	$\mathbb{C}$	$\mathbb{C}$	
	T1 alarm limits	36.0~39.0 on	36.0~39.0 on	
	T2 alarm limits	36.0~39.0 on	36.0~39.0 on	
	TD alarm limits	0.0~2.0 on	0.0~2.0 on	

# **21** Common faults and maintenance

The following table shows the common faults on the operation, and the solution.

Faults	Solution	
Blank Screen	Connects the monitor to check the screen and screen line	
	whether normal.	
The system time is not correct	1. Set up error, can be reset through the system User	
	Maintenance menu.	
	2. The button battery on main control board is run out,	
	please change the button battery.	
No ECG waveform	See the ECG cable and lead-wires whether in good	
	condition, disconnected or electrode rusting result in connection fail.	
	2. Look at whether the ECG cable and lead type are	
	consistent.	
Unable to do ST analysis	1. Check the ECG Setup→ ST Analysis→ ST Analysis is	
	set to "On".	
	2. Check the ECG Setup→ Other Setup→ Paced whether	
	be set to "On". If Paced is set to "On", means the	
	patient have a pacemaker; in this case the machine is	
	not doing the ST analysis.	
No SpO₂ waveform or value	Check whether the SPO2 Sensor is connected and in good condition.	
Blood pressure does not start	1. Check whether the pump is broken.	
	2. Check whether the trachea is broken.	
	3. Check whether the blood pressure plate is normal.	
Blood pressure started, but	1. Check whether the blood pressure cuff is leakage.	
couldn't measure the value	2. Check whether the NIBP extension tube and machine	
	connect is well.	
	3. Check whether the deflating valve on blood pressure	
	plate is normal.	
	4. Check whether the pressure sensor is normal.	

If the above doesn't solve the problem, please contact Bistos after-sales department or dealers.

# 22 Manufacturer's declaration on EMC

BT-770 needs special precautions regarding EMC (Electromagnetic compatibility) and needs to be used according to the EMC information provided in this user manual. Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect the BT-770 and should be kept at least 1 m away from the equipment.

#### NOTE

- Using unqualified accessories, sensors and cables will increase the electromagnetic emission and reduce the electromagnetic immunity of the device.
- Do not put the device close to other devices or stack together. When necessary, observe the device closely to ensure that it runs normally in the environment.
- The device requires special EMC protection, and it is necessary to install and maintain it in the environment that meets the following EMC information.
- Even if other devices comply with CISPR emission requirements, they may also cause interference to this device.
- When the input signal amplitude is smaller than the minimum amplitude specified in the technical specifications, it may result in inaccurate measurements.
- Mobile communication devices or wireless network devices may have an impact on the device.

# 22.1 Electromagnetic emissions

The BT-770 is intended for use in the electromagnetic environment specified below. The customer or the user of the BT-770 should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance	
RF emissions CISPR 11	Group 1	The BT-770 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF emissions CISPR 11	Class A	The BT-770 is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings	
Harmonic emissions IEC 61000-3-2	Class A	used for domestic purposes, provided the following warning is heeded:  Warning: This BT-770 is intended for use by healthcare professionals only. This equipment/	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the BT-770 or shielding the location.	

# 22.2 Recommended separation distances between portable and mobile RF communications equipment and BT-770

The BT-770 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the BT-770 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BT-770 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum	Separation dis	stance according to frequenc [m]	y of transmitter
output power of transmitter [W]	150 kHz to 80 MHz $d = 3.5\sqrt{p}$	80 MHz to 800 MHz $d=3.5\sqrt{p}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3}\right] \sqrt{p}$
0.01	0.35	0.35	0.23
0.1	1.11	1.11	0.74
1	3.5	3.5	2.34
10	11.07	11.07	7.38
100	35	35	23.24

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1) At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

P/N : 770-ENG-OPM-EUR-R00 Bistos Co., Ltd. 2018.09

# 22.3 Electromagnetic immunity

The BT-770 is intended for use in the electromagnetic environment specified below. The customer or the user of the BT-770 should assure that it is used in such an environment.

Immunity test	IEC 60601	Compliance level	Electromagnetic
,	Test level		environment -guidance
Electrostatic	±8 kV Contact	±8 kV Contact	Floors should be wood,
discharge (ESD)			concrete or ceramic tile. If
	±15 kV air	±15 kV air	floors are covered with
IEC 61000-4-2:2009			synthetic material, the
			relative humidity should
			be at least 30 %.
Electrical fast	±2 kV for power	±2 kV for power	Mains power quality
transient/burst	supply lines	supply lines	should be that of a typical
	±1 kV for	±1 kV for	commercial or hospital
IEC 61000-4-4:2004	input/output lines	input/output lines	environment.
	(>3m)	(>3m)	
Surge	±1 kV differential	±1 kV differential	Mains power quality
	mode	mode	should be that of a typical
IEC 61000-4-5:2006	±2 kV common	±2 kV common	commercial or hospital
	mode	mode	environment.
Voltage dips, short	< 5 % <i>U</i> т (> 95 %	< 5 % <i>U</i> т (> 95 %	Mains power quality
interruptions and	dip in $U_T$ ) for 0.5	dip in $U_T$ ) for 0.5	should be that of a typical
voltage variations on	cycles	cycle	commercial or hospital
power supply input			environment. If the user
lines	40 % <i>U</i> т (60 % dip	40 % <i>U</i> т (60 % dip	of the BT-770 requires
	in <i>U</i> τ ) for 5 cycles	in <i>U</i> τ ) for 5 cycles	continued operation
IEC 61000-4-11:2004			during power mains
	70 % <i>U</i> т (30 % dip	70 % <i>U</i> т (30 % dip	interruptions, it is
	in <i>U</i> τ) for 25	in <i>U</i> τ) for 25	recommended that the
	cycles	cycles	BT-770 be powered from
			an uninterruptible power
	<5 % <i>U</i> т (> 95 %	<5 % <i>U</i> т (> 95 %	supply.
	dip in <i>U</i> τ ) for 5 s	dip in <i>U</i> τ ) for 5 s	
Power frequency (50	3 A/m	3 A/m	Power frequency
Hz and 60 Hz)			magnetic fields should be
magnetic field			at levels characteristic of
			a typical location in a
IEC 61000-4-8:2010			typical commercial or
			hospital environment.

NOTE  $U_T$  is the a.c. mains voltage prior to application of the test level.

The BT-770 is intended for use in the electromagnetic environment specified below. The customer or the user of the BT-770 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level
Conducted RF	3 Vrms	3 Vrms
IEC 61000-4-6:2009	150 kHz to 80 MHz	
Radiated RF	3 V/m	3 V/m
IEC 61000-4-3	80 MHz to 2.5 GHz	

#### Electromagnetic environment - guidance

Portable mobile RF communications equipment should be used no closer to any part of the BT-770, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

#### **Recommended separation distance**

$$d - 1.2\sqrt{p} \ (d - 3.5\sqrt{p})$$

$$d - 1.2\sqrt{p}$$
 (Resp:  $d - 3.5\sqrt{p}$ ) 80 to 800MHz

$$d - 1.2\sqrt{p}$$
 800M to 2.5GHz

where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and *d* is the recommended separation distance in meters (m).

Field strengths from fixed RF transmitters, as deter-mined by an electromagnetic site survey <sup>a</sup>, should be less than the compliance level in each frequency range. <sup>b</sup>

Interference may occur in the vicinity of equipment marked with the following symbol:



NOTE 1) At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- <sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BT-770 is used exceeds the applicable RF compliance level above, the BT-770 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the BT-770.
- <sup>b</sup> Over the frequency range 150 kHz to 80MHz, field strengths should be less than 3 V/m.

## **Product Warranty**

Product Name	Patient Monitor
Model Name	BT-770
Serial No.	
Warranty Period	2 Years
Date of Purchase	
Customer	Hospital:
	Address:
	Name:
	Telephone:
Sales Agency	
Manufacture	Bistos Co., Ltd.

- \* Thank you for purchasing BT-770.
- \* This product is manufactured and passed through strict quality control and inspection.
- Compensation standard concerning repair, replacement, refund of the product complies
   with "Framework Act on Consumers" noticed by Fair Trade Commission of Republic of
   Korea.

#### Service Telephone and Fax. Numbers

Telephone: +82 31 750 0340 Fax: +82 31 750 0344

Bistos Co., Ltd.

7<sup>th</sup> FL., A Bldg., Woolim Lions Valley 5-cha, 302,
Galmachi-ro, Jungwon-gu, Seongnam-si,
Gyeonggi-do, Korea

www.bistos.co.kr bistos@bistos.co.kr

EC REP Obelis s.a

Bd. Général Wahis 53 1030 Brussels / BELGIUM

Telephone: +32 2. 732.59.54 Fax: +49 2. 732.60.03







# EC CERTIFICATE Full Quality Assurance System

Certificate No.: 243269-2017-CE-KOR-NA-PS Rev. 5.0

Project No.: PRJC-533956-2015-MSL-KOR

Valid Until: 01 September 2023

This is to certify that the quality system of:

# Bistos Co., Ltd.

7th Fl., A Bldg., Woolim Lions Valley 5-cha, 302, Galmachi-ro, Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea

For design, production and final product inspection/testing of:

Monitoring devices of vital physiological parameters and Utilising non-ionizing radiation

Has been assessed with respect to:

The conformity assessment procedure described in Annex II excluding section 4 of Council Directive 93/42/EEC on Medical Devices, as amended

and found to comply

Place and date: Høvik, 30th April 2021

**Check Validity** 

For the issuing office: Notified Body 2460 DNV Product Assurance AS



Hazem Tinawi Technical Reviewer



Certificate No.: 243269-2017-CE-KOR-NA-PS Rev. 5.0 Place and date: Høvik, 30th April 2021

Further details of the product(s) and conditions for certification are given overleaf.

#### **Jurisdiction**

Application of Council Directive 93/42/EEC of 14 June 1993, adopted as "Forskrift om Medisinsk Utstyr" by the Norwegian Ministry of Health and Care Services.

#### Certificate history:

Revision	Description	Issue Date
0.0	Replaces certificate EU1308401, Rev2.0 (NB 0470) following transfer of Notified Body functions to DNV GL Nemko Presafe AS (NB 2460)	01 September 2017
1.0	EU Rep change	13 April 2018
2.0	Re-certification for Fetal monitor and Neonatal Phototherapy unit (BT-300, BT-350, FM-20, Biocare FM-1, BT-400) Scope extension for pulse oximeter and patient monitor (BT-710, BT-720, BT-740, BT-770) The accessories (Feotal Doppler system probe and Cardiotocograph transducers) are removed (AY-DOP-300, AY-DOP-350, AY-UC-300, AY-UC-350)	01 September 2018
3.0	Editorial change	13 February 2020
4.0	Scope extension to new model (BT-780)	26 April 2021
5.0	Editorial change in model name (typo error)	30 <sup>th</sup> April 2021

#### Products covered by this Certificate:

Product Description	Product Name	Class
Fetal monitor	<ul> <li>BT-300</li> <li>BT-350</li> <li>FM-20</li> <li>Biocare FM-1</li> </ul>	lla
Neonatal Phototherapy unit	■ BT-400	lla
Pulse Oximeter	■ BT-710	IIb
Patient Monitor	<ul> <li>BT-720</li> <li>BT-740</li> <li>BT-770</li> <li>BT-780</li> </ul>	IIb



Certificate No.: 243269-2017-CE-KOR-NA-PS Rev. 5.0 Place and date: Høvik, 30th April 2021

The complete list of devices is filed with the Notified Body

#### Sites covered by this certificate

Site Name	Address
Bistos Co., Ltd.	7th Fl., A Bldg., Woolim Lions Valley 5-cha, 302, Galmachi-ro, Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea

#### **EU Representative**

OBELIS S.A, Bd. General Wahis, 53, 1030 Brussels, Belgium





Certificate No.: 243269-2017-CE-KOR-NA-PS Rev. 5.0

Place and date: Høvik, 30th April 2021

#### Terms and conditions

The certificate is subject to the following terms and conditions:

- Any producer (see 2001/95/EC for a precise definition) is liable for damage caused by a
  defect in his product(s), in accordance with directive 85/374/EEC, as amended, concerning
  liability of defective products.
- The certificate is only valid for the products and/or manufacturing premises listed above.
- The Manufacturer shall fulfil the obligations arising out of the quality system as approved and uphold it so that it remains adequate and efficient.
- The Manufacturer shall inform the Notified Body of any intended updating of the quality system and the Notified Body will assess the changes and decide if the certificate remains valid.
- Periodical audits will be held, in order to verify that the Manufacturer maintains and applies
  the quality system. the Notified Body reserves the right, on a spot basis or based on
  suspicion, to pay unannounced visits.

The following may render this Certificate invalid:

- Changes in the quality system affecting production.
- Periodical audits not held within the allowed time window.

## Conformity declaration and marking of product

When meeting with the terms and conditions above, the producer may draw up an EC declaration of conformity and legally affix the CE mark followed by the Notified Body identification number.

End of Certificate



# Management System Certificate

Certificate No.: 243275-2017-AQ-KOR-NA-PS Rev 4.0

Initial Certification Date: 12 August 2004

Valid Until: 09 September 2024

This is to certify that the quality system of:

# Bistos Co., Ltd.

7th Fl., A Bldg., Woolim Lions Valley 5-cha, 302, Galmachi-ro, Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea

has been found to conform to the Quality Management System standard:

ISO 13485:2016/NS-EN ISO 13485:2016

This certificate is valid for the following scope:

Design and Development, Manufacturing, Sales, Distribution, and Servicing of Ultrasound Doppler system, Fetal monitor, Phototherapy, Patient Monitor, Pulse Oximeter, Incubator, Head-worn light, Infant Warmer and Electric Breast Pump.

Place and date: Høvik, 23 June 2021

**Check Validity** 



NORWEGIAN ACCREDITATION For the issuing office: DNV Product Assurance AS

Tone Elise Kolpus Lead Auditor

MSYS 018

ficate is subject to terms and conditions as set out in the Certification Agreement. Failure to comply may render this Certificate invalid.



Certificate No.: 243275-2017-AQ-KOR-NA-PS Rev. 4.0 Place and date: Høvik, 22 June 2021

Site Name	Address	Site Specific Scope
Head Office	7th Fl., A Bldg., Woolim Lions Valley 5-cha, 302, Galmachi-ro, Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea	Design and Development, Sales, Distribution, and Servicing of Ultrasound Doppler system, Fetal monitor, Phototherapy, Patient Monitor, Pulse Oximeter, Incubator, Head-worn light, Infant Warmer and Electric Breast Pump.
Factory	116~122ho, Jungang Induspia 3, 27, Sagimakgol-ro 105beon-gil, Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea	Manufacturing of Ultrasound Doppler system, Fetal monitor, Phototherapy, Patient Monitor, Pulse Oximeter, Incubator, Head-worn light, Infant Warmer and Electric Breast Pump.