

EMDP	100000
Ethyl Glucuronide (ETG)	
Ethyl Glucuronide	500
Fentanyl(FTY)	
Norfentanyl	20
Fentanyl	100
Synthetic Cannabis (K2)	
JWH-018 Pentanoic Acid	25
JWH-073 Butanoic Acid	25
JWH-018 N-4-hydroxypentyl	2000
JWH-018 (Spice Cannabinoid)	1000
JWH-018 4-Hydroxypentyl metabolite-D5 (indole-D5)	1000
JWH-073 (Spice Cannabinoid)	2000
JWH-073 3-Hydroxybutyl metabolite	1000
JWH-073 3-Hydroxybutyl metabolite-D5 (indole-D5)	1000
JWH-019 6-hydroxypentyl	1000
JWH-122 N-4-hydroxypentyl	2000
JWH-210 5-Hydroxypentyl metabolite	5000
AM2201 4-Hydroxypentyl metabolite	1000
Ketamine (KET)	
Ketamine	1000
Methadone	50000
Pethidine	12500
Methylamphetamine	12500
Methoxyphenamine	12500
Promethazine	25000
Phencyclidine	25000
Methcathinone (MC)	
Methcathinone	1000
Methylene dioxymethamphetamine (MDMA)	
3,4-Methylenedioxyamphetamine HCl (MDMA)	500
3,4-Methylenedioxyamphetamine HCl (MDA)	3000
3,4-Methylenedioxyethylamphetamine (MDE)	300
Methamphetamine (MET)	
D(+)-Methamphetamine	1000

D-Amphetamine	50000
Chloroquine	50000
(+/-)-Ephedrine	50000
(-)-Methamphetamine	25000
(+/-)3,4-methylenedioxyamphetamine(MDMA)	2000
b-Phenylethylamine	50000
Trimethobenzamide	10000
Morphine (MOP)	
Morphine	300
Codeine	300
Ethyl Morphine	300
Hydrocodone	5000
Hydromorphone	5000
Morphine-3-b-d-glucuronide	1000
Thebaine	30000
Methaqualone (MQL)	
Methaqualone	300
Methadone (MTD)	
Methadone	300
Doxylamine	50000
Opiate (OPI)	
Morphine	2000
Codeine	2000
Ethylmorphine	5000
Hydrocodone	12500
Hydromorphone	5000
Levorphanol	75000
s-Monoacetylmorphine	5000
Morphine 3-b-D-glucuronide	2000
Norcodeine	12500
Normorphone	50000
Oxycodone	25000
Oxymorphone	25000
Thebaine	100000
Oxycodone (OXY)	
Oxycodone	100
Dihydrocodeine	20000

Codeine	100000
Hydromorphone	100000
Morphine	>100,000
Acetylmorphine	>100,000
Buprenorphine	>100,000
Ethylmorphine	>100,000
Phencyclidine (PCP)	
Phencyclidine	25
4-Hydroxyphencyclidine	12500
Propoxyphene (PPX)	
d-Propoxyphene	300
d-Norpropoxyphene	300
Tricyclic Antidepressants (TCA)	
Notriptyline	1000
Nordoxepine	1000
Trimipramine	3000
Amitriptyline	1500
Promazine	1500
Desipramine	200
Imipramine	400
Clomipramine	12500
Doxepine	2000
Maprotiline	2000
Promethazine	25000
Marijuana (THC)	
11-nor-D9-THC-9-COOH	50
11-nor-D8-THC-9-COOH	30
11-hydroxy-D9-Tetrahydrocannabinol	2500
D8- Tetrahydrocannabinol	7500
D9- Tetrahydrocannabinol	10000
Cannabinol	10000
Cannabidiol	100000
Tramadol (TRA)	
Tramadol	100
(+/-) Chlorpheniramine	500000
Dipehnhydramine	250000
Pheniramine	>500,000
PCM	>250,000

Cross-Reactivity

Considering the complexity of clinical urine specimens and the possibility that various urine specimens contain potentially interfering substances, we simulated above situations by adding the potentially interfering substances to a certain concentration as specimen. The following components show no cross-reactivity when tested with One Step Multi-Drug Test Panel at a concentration of 100 µg/ml.

Non Crossing-Reacting Compounds

Acetaminophen	Epinephrine HCl	Metoprolol tartrate
Aciclovir	Esomeprazole	Mifepristone
Afrin	Estroven	Montelukast
Aleve	Fenofibrate	Mosapride Citrate
Amiodarone HCl	Fluvoxamine	Narcotine
Amlodipine Mesylate	Fuel	Nifedipine
Amoxicillin	Gabapentin	Nikethamide
Ampicillin	Glibenclamide	Nimodipine
Aripiprazole	Gliclazide	Omeprazole
Aspirin	Glipizide	Penfluridol
Atorvastatin	Glucosamine Chondroitin	Penicillin V Potassium
Atropine	Glucose	Pioglitazone HCl
Captopril	Haloperidol	Piracetam
Carbamazepine	Heartburn Relief	Pravastatin sodium
Cefaclor	Hydrochlorothiazide	Propylthiouracil
Cefradine	I Caps	Rifampicin
Cephalexin	Isosorbide dinitrate	Sildenafil citrate
Ciprofloxacin	Ketoconazole	Simvastatin
Clarithromycin	Levofloxacin	Spirolactone
Clopidogrel bisulfate	Levonorgestrel	Tetracycline
Clozapine	Levothyroxine sodium	Trazodone HCl
Cortisone	Lidocaine HCl	Triamterene
CVS	Lisinopril	Vitamin B1
Dextromethorphan HBr	Lithium carbonate	Vitamin B2
Diclofenac sodium	Loratadine	Vitamin C
Domperidone	Magnesium	Zencore Plus2
Enalapril maleate	Mega-T Plus	













From the results above, it is clear that Multi-Drug Test panel resists well against interference from these substances.

BIBLIOGRAPHY OF SUGGESTED READING

- Baselt, R.C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications, Davis, CA, 1982.
- Ellenhorn, M.J. and Barceloux, D. G Medical Toxicology. Elsevier Science Publishing Company, Inc., New York, 1988
- Gilman, A. G., and Goodman, L. S. The Pharmacological Fluids, in Martin WR(ed): Drug Addiction I, New York, Spring – Verlag, 1977.
- Harvey, R.A., Champe, P.C. Lippincotts Illustrated Reviews. Pharmacology. 91-95, 1992.

- Hawks RL, CN Chiang. Testing for drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monography 73, 1986
- Hofmann F.E., A Handbook on Drug and Alcohol Abuse: The Biomedical Aspects, New York, Oxford University Press, 1983.
- McBay, A. J. Clin. Chem. 33,33B-40B, 1987.

MEANING OF SYMBOLS ON PACKAGE

	Consult instructions for use		Tests per kit		Authorized Representative
	For in vitro diagnostic use only		Use by		Do not reuse
	Store between 2-30°C		Lot Number		Do not use if package is damaged
	Manufacturer		Keep dry		Keep away from sunlight

	Hangzhou Aichek Medical Technology CO.,Ltd. Jinxing Cun, Yuhang Community, Yuhang District (Future Sci-Tech City), Hangzhou, Zhejiang, P.R. China
	SUNGO Europe B.V. Olympisch Stadion 24, 1076DE Amsterdam, Netherlands

Version: 02
Effective Date: 2021-11-20