Prospektüs

Prospectus

FASTLAK

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FASTLAK

Instruction for Use

Surgical, Poly(glycolide(90%)-co-lactide(10%))
Rapid Absorbable PGLA Suture
Synthetic, Multifilament, Braided, Undyed

DESCRIPTION

FASTLAK suture is synthetic absorbable sterile surgical suture which is flexible strand composed of a copolymer synthesized from 90% glycolide and 10 % L-lactide (derived from glycolic and lactic acids). The empirical formula of the copolymer is $(C_4H_6O_2)_m(C_3H_4O_2)_n$. Multifilament suture is consisting of elementary fibres which are assembled by braiding.

FASTLAK is gamma degraded or having lower molecular weight form of PGL to supply rapid hydrolysis/absorption of PGL.

FASTLAK suture is prepared by coating suture material with a mixture composed of equal parts calcium stearate ($C_{18}H_{35}O_2$) and 30% glycolide and 70% lactide ($C_3H_4O_2$, $C_4H_6O_2$) where $m+n=3.7$. When FASTLAK suture is introduced into a living organism it is absorbed by that organism cause no undue tissue irritation. Copolymer and the coating with calcium stearate have been found to be nonantigenic, nonpyrogenic and elicit only a mild tissue reaction during absorption.

FASTLAK suture has been presented in USP 8/0 to 2 (metric sizes 0.4 and 5) gauge sizes in a variety of lengths, as non-needed or attached to stainless steel needles of varying types and sizes. The suture is available in the undyed (natural, beige) form.

FASTLAK sutures meet United States Pharmacopoeia (U.S.P.) requirements as described in the U.S.P. "Monograph for Absorbable Surgical Sutures" and European Pharmacopoeia (E.P.) requirements as described in the E.P. "Monograph for Sutures, Sterile Synthetic Absorbable Braided, 01/2008:0667" with the exception of an occasional slight oversize in some gauges and the exception of knot tensile strength. The knot tensile strength of FASTLAK meets USP requirements for collagen sutures as well as the requirements of the European Pharmacopoeia for "Chorda Resorbilis Sterilis".

INDICATIONS

FASTLAK suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures but not for use in cardiovascular and neurological tissues.

ACTIONS

FASTLAK suture elicits a minimal acute inflammatory reaction in tissue and ingrowth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of suture occur by means of hydrolysis gradually and decrease the strength in the body. After hydrolysis it is executed from the body as carbon dioxide and water. Absorption begins as a loss of tensile strength followed by a loss of mass.

Animal and in-vitro hydrolysis studies indicate that FASTLAK suture provides 50% tissue support during at least one week post implantation. All of the original tensile strength (Break Strength Retention-BSR) is lost between 10-14 days post implantation and suture absorption is essentially complete 42 days.

Contraindications

This suture, being absorbable, should not be used where extended approximation of tissue is required.

WARNINGS

Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing for wound closure, as risk of wound dehiscence may vary with the site of application and the tissue material used. Physicians should consider the in vivo performance (under ACTIONS section) when selecting a suture for use in patients.

The use of this suture may be inappropriate in elderly malnourished or debilitated patients, or in patients suffering from conditions which may delay wound healing.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation.

As FASTLAK is an absorbable suture material, the use of supplemental non-absorbable sutures should be considered by the surgeon in the closure of the sites which may undergo expansion, stretching, or distension, or which may require additional support.

Do not re-sterilize. Discard open packages and unused sutures.
Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

PRECAUTIONS

Skin sutures which must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed. Under some circumstances, notably orthopedic procedures, immobilization of joints by external support may be employed at the discretion of the surgeon.

Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply, as suture extrusion and delayed absorption may occur.

In handling FASTLAK or any other suture materials, care should be taken to avoid damage to needle and suture. Avoid crushing or ripping damage due to application of surgical instruments such as forceps or needle holders.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks.

FASTLAK sutures, which are treated to enhance handling characteristics require the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstances and experiences of the surgeon.

Avoid prolonged exposure to elevated temperatures.

Dispose of contaminated and unused products are in accordance with local and facility requirements. Discard used needles in "sharps" containers.

ADVERSE REACTIONS

Adverse reactions associates with the use of FASTLAK sutures include transitory local irritation at the wound site, acute inflammatory tissue reaction, erythema and induration during the absorption process of subcuticular sutures. Like all foreign bodies FASTLAK may irritate an existing infection.

HOW SUPPLIED

FASTLAK suture is available sterile, as braided undyed (natural) strands in sizes 8/0 to 2 (metric sizes 0.4-5), in a variety of lengths, as non-needed or attached to stainless steel needles of varying types and sizes. FASTLAK suture is for single use only.

STORAGE

Store below 25°C and keep away from sunlight. Protect from humidity. Do not use after expiry date.

SYMBOLS USED ON LABELLING

Do not request

Do not sterilize

Do not use if pack age is damaged

Manufacturer

Date of Manufacture, Year

Exp. Date, Year-Month

Sterile R: Irradiation

Batch Number

Undyed, Absorbable, Braided, Coated

Catalogue Number

Store below 25°C

Keep away from sunlight

Protect from humidity

Recyclable pack

Attention, See instruction for use

2292

IU-FL rev-06, 15-01-2018

Issue date: 11.09.2012

"EASSI (The European Association of Surgical Suture Industry) has developed a system of symbols which is designed to describe various suture product characteristics in an intuitive, pictorial manner. The use of symbols is permitted by the Medical Device Directive (MDD 93/42/EEC) and enables companies to provide information to the customer without having to provide multilingual translations."

BOZ TIBBİ MALZEME SANAYİ VE TİCARET A.Ş.
Sağlık Mah. Sağlık Sokakı No:33/55/HİTye-Cankaya, ANKARA/TÜRKİYE
Tel: +90 (312) 254 03 40 (5 hat) Faks:+90 (312) 254 03 50
web: www.boztibbi.com e-mail: boz@boztibbi.com

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GB

FASTSORB

Surgical, Polyglycolic Acid
Rapid Absorbable Suture
Synthetic, Multifilament, Braided, Undyed

Instruction for Use**DESCRIPTION**

FASTSORB Surgical Degraded Polyglycolic acid (PGA) suture; synthetic, sterile, absorbable, is homopolymer of glycolic acid. The empirical formula of polyglycolic acid is $(-O-CH_2-CO-O-CH_2-CO-)_n$. The polymer chain is degraded by gamma radiation. FASTSORB is coated with the mixture of calcium stearate $(C_{18}H_{35}O_2Ca)$ and polycaprolactone $(C_{12}H_{22}O_5)_n$. The suture and coating material are inert, noncollagenous, nonantigenic, and non-pyrogenic. FASTSORB suture is undyed in natural beige color.

INDICATIONS

FASTSORB suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures but not for use in cardiovascular and neurological tissues. FASTSORB suture is for single use only.

ACTIONS

PGA suture elicits a minimal acute inflammatory reaction in tissue and ingrowth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of PGA suture occurs by means of hydrolysis gradually and decreases the strength in the body. After hydrolysis it's executed from the body as carbon dioxide and water. Absorption begins as a loss of tensile strength followed by a loss of mass. PGA suture retains approximately 75% of the original tensile strength (Break Strength Retention-BSR) in 7 days. It loses its strength in 42 days.

CONTRAINDICATIONS

The sutures being absorbable should not be used where extended approximation of tissue is required. The sutures should not be used for cardiovascular and neurosurgical tissue.

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Under some circumstances, notably orthopedic procedures, immobilization of joint by external support may be employed at the discretion of the surgeon.

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PRECAUTIONS

In handling FASTSORB or any other suture materials, care should be taken to avoid damage to needle and suture. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point.

Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks.

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ADVERSE REACTIONS

Adverse effects associated with the use of the device include wound dehiscence, failure to provide adequate wound support in closure of the sites where expansion, stretching, or distension occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients, or in patients suffering from conditions which may delay wound healing, infection, minimal acute inflammatory tissue reaction localized irritation when skin sutures are left in place for greater than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calculus formation in urinary and biliary tracts when prolonged contact with salt solution such as urine and bile occurs, and transitory local irritation at the wound site.

HOW SUPPLIED

FASTSORB sutures are available as sterile, braided, undyed (beige) strands in sizes 8/0 thru 6 (metric sizes 0.4 – 8.0) in a variety of lengths, with a variety of needles. FASTSORB sutures are available in one, two and three dozen boxes. FASTSORB suture is supplied sterile.

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
AR

Fastsorb

Surgical Suture, Polyglycolicacid (PGA)
Fast Absorbable, Synthetic, Braided, Undyed

Surgical Sutures
Absorbable

Boz®

Definition

Synthetic absorbable, coated, short term wound support and short term mass absorption.

Material

Polyglycolic Acid (PGA)

Coating

Calcium Stearate and polycaprolactone

Structure

Braided

Colour

Beige (undyed)

Size

EP : 0.4 to 8

USP : 8/0 to 6

Tensile Strength Retention

7 Days 75 %

Mass Absorption

42 days

Sterilization Method

Gamma Irradiation

Main Indication

Gynaecology (episiotomies)

Urology

Ophtalmic surgery (conjunctive suturing)

Odontology (oral mucosa)









Pediatric surgery

Skin closer - intracutaneous, subcutaneous.

Ligatures




Fastsorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
12 mm							
3/8 Reverse Cutting 330 Micron			FS1545	6/0 0,7	75	12	Undyed
13 mm							
1/2 Round Body, Taper Point 330 Micron			FS1563	6/0 0,7	75	12	Undyed
1/2 Round Body, Taper Point 380 Micron			FS1886	5/0 1,0	75	12	Undyed
3/8 Reverse Cutting 330 Micron			FS5912	6/0 0,7	75	12	Undyed
3/8 Reverse Cutting 380 Micron			FS1929	5/0 1,0	75	12	Undyed
3/8 Reverse Cutting Premium 380 Micron			FS1941	5/0 1,0	45	12	Undyed
16 mm							
1/2 Round Body, Taper Point			FS1989	5/0 1,0	75	12	Undyed
			FS2503	4/0 1,5	75	12	Undyed
3/8 Reverse Cutting			FS1673	6/0 0,7	75	12	Undyed
			FS2019	5/0 1,0	75	12	Undyed
			FS2539	4/0 1,5	75	12	Undyed
3/8 Inside Cutting			FS2557	4/0 1,5	75	12	Undyed
			FS3038	3/0 2,0	75	12	Undyed
19 mm							
3/8 Reverse Cutting			FS2639	4/0 1,5	75	12	Undyed
			FS3108	3/0 2,0	75	12	Undyed
20 mm							
1/2 Round Body, Taper Point			FS2228	5/0 1,0	75	12	Undyed

Fastsorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm							
1/2 Round Body, Taper Point			FS2683	4/0 1,5	75	12	Undyed
			FS3141	3/0 2,0	75	12	Undyed
3/8 Reverse Cutting			FS2250	5/0 1,0	75	12	Undyed
22 mm							
1/2 Round Body, Taper Point			FS2745	4/0 1,5	75	12	Undyed
25 mm							
1/2 Round Body, Taper Point			FS2787	4/0 1,5	75	12	Undyed
			FS3251	3/0 2,0	75	12	Undyed
3/8 Reverse Cutting			FS3273	3/0 2,0	75	12	Undyed
			FS3756	2/0 3,0	75	12	Undyed
26 mm							
1/2 Round Body, Taper Point			FS3785	2/0 3,0	75	12	Undyed
			FS4214	0 3,5	75	12	Undyed
3/8 Reverse Cutting			FS3821	2/0 3,0	75	12	Undyed
			FS4236	0 3,5	75	12	Undyed
			FS3345	3/0 2,0	75	12	Undyed
30 mm							
1/2 Round Body, Taper Point			FS3368	3/0 2,0	75	12	Undyed
			FS3857	2/0 3,0	75	12	Undyed
			FS4243	0 3,5	75	12	Undyed
			FS4607	1 4,0	75	12	Undyed





Fastorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm						
3/8 Reverse Cutting		FS2885	4/0 1,5	75	12	Undyed
		FS3396	3/0 2,0	75	12	Undyed
		FS4628	1 4,0	75	12	Undyed
35 mm						
1/2 Round Body, Taper Point		FS3915	2/0 3,0	75	12	Undyed
		FS4299	0 3,5	75	12	Undyed
		FS4305	0 3,5	90	12	Undyed
		FS4643	1 4,0	75	12	Undyed
		FS4648	1 4,0	90	12	Undyed
3/8 Reverse Cutting		FS3422	3/0 2,0	75	12	Undyed
		FS3932	2/0 3,0	75	12	Undyed
37 mm						
1/2 Round Body, Taper Point		FS3954	2/0 3,0	90	12	Undyed
40 mm						
1/2 Round Body, Taper Point		FS3981	2/0 3,0	75	12	Undyed
		FS3988	2/0 3,0	90	12	Undyed
		FS4367	0 3,5	75	12	Undyed
		FS4374	0 3,5	90	12	Undyed

Fastorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
1/2 Round Body, Taper Point 	FS4712	1 4,0	75	12	Undyed
	FS4718	1 4,0	90	12	Undyed
	FS5001	2 5,0	75	12	Undyed
3/8 Reverse Cutting 	FS4005	2/0 3,0	75	12	Undyed
45 mm					
1/2 Round Body, Taper Point 	FS4421	0 3,5	75	12	Undyed
	FS4771	1 4,0	75	12	Undyed
	FS5034	2 5,0	75	12	Undyed
Straight Reverse Cutting 	FS2911	4/0 1,5	75	12	Undyed
50 mm					
1/2 Round Body, Taper Point 	FS4037	2/0 3,0	75	12	Undyed
	FS4453	0 3,5	75	12	Undyed
	FS4811	1 4,0	75	12	Undyed
	FS5065	2 5,0	75	12	Undyed
65 mm					
1/2 Round Body, Taper Point 	FS4864	1 4,0	100	12	Undyed

Fastsorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
150 cm					
FASTSORB Non-Needle Suture	FS1373	6/0 0,7	150	12	Undyed
	FS2325	4/0 1,5	150	12	Undyed
	FS2934	3/0 2,0	150	12	Undyed
	FS3510	2/0 3,0	150	12	Undyed
	FS4113	0 3,5	150	12	Undyed
	FS4531	1 4,0	150	12	Undyed
	FS4907	2 5,0	150	12	Undyed

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GB

GLIKOSORB

Surgical, Polyglycolic Acid
Absorbable PGA Suture
Synthetic, Multifilament, Braided, Undyed/Violet

Instruction for Use**DESCRIPTION**

GLIKOSORB Surgical Polyglycolic acid (PGA) suture; synthetic, sterile, absorbable, is homopolymer of glycolic acid. The empirical formula of polyglycolic acid is $(-O-CH_2-CO-O-CH_2-CO-)_n$. GLIKOSORB is coated with the mixture of calcium stearate $(C_{18}H_{35}O_4)_2$ and polycaprolactone $(C_{12}H_{22}O_5)_n$. The suture and coating material are inert, noncollagenous, nonantigenic, and non-pyrogenic. Coloured violet (D&C Violet No. 2 C.I. # 60725 / Solvent Violet 13) to enhance visibility in tissue.

INDICATIONS

GLIKOSORB suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures but not for use in cardiovascular and neurological tissues. GLIKOSORB suture is for single use only.

ACTIONS

PGA suture elicits a minimal acute inflammatory reaction in tissue and ingrowth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of PGA suture occurs by means of hydrolysis gradually and decreases the strength in the body. After hydrolysis it's executed from the body as carbon dioxide and water. Absorption begins as a loss of tensile strength followed by a loss of mass. PGA suture retains approximately 75% of the original tensile strength (Break Strength Retention-BSR) at two weeks post implantation. At three weeks, approximately 50% of the original tensile strength (Break Strength Retention-BSR) is retained for sizes 6-0 and larger and approximately 40% of its original tensile strength (Break Strength Retention-BSR) is retained for sizes 7-0 and smaller. All of the original tensile strength (Break Strength Retention-BSR) is lost between four and five weeks post implantation. In 60-90 days it is totally absorbed.

CONTRAINDICATIONS

The sutures being absorbable should not be used where extended approximation of tissue is required. The sutures should not be used for cardiovascular and neurological tissue.

WARNINGS

Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Physicians should consider the in vivo performance (under ACTIONS section) when selecting a suture for use in patients.

Skin sutures which must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated.

Under some circumstances, notably orthopedic procedures, immobilization of joint by external support may be employed at the discretion of the surgeon. Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply, as suture extrusion and delayed absorption may occur.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation. Do not re-sterilize. Discard open packages and unused sutures.

Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

PRECAUTIONS

In handling PGA or any other suture materials, care should be taken to avoid damage to needle and suture. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

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PGA sutures, which are treated to enhance handling characteristics requires the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstances and experiences of the surgeon.

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ADVERSE REACTIONS

Adverse effects associated with the use of the device include wound dehiscence, failure to provide adequate wound support in closure of the sites where expansion, stretching, or distension occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients, or in patients suffering from conditions which may delay wound healing, infection, minimal acute inflammatory tissue reaction localized irritation when skin sutures are left in place for greater than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calculus formation in urinary and biliary tracts when prolonged contact with salt solution such as urine and bile occurs, and transitory local irritation at wound site.

HOW SUPPLIED

GLIKOSORB sutures are available braided, dyed (violet) and undyed strands in sizes 8/0 thru 6 (metric sizes 0.4 – 8) in a variety of lengths, with a variety of needles. GLIKOSORB sutures are available in one, two and three dozen boxes. GLIKOSORB suture is supplied sterile.

STORAGE

Store below 25 °C and keep away from sunlight. Protect from humidity. Do not use after expiry date.

SYMBOLS USED ON LABELLING

 Do not reuse	 Undyed, Absorbable, Braided, Coated
 Do not re-sterilize	 Catalogue Number
 Do not use if pack age is damaged	 Store below 25°C
 Manufacturer	 Keep away from sunlight
YYYY  Date of Manufacture, Year	 Protect from humidity
 YYYY-MM Expiry Date, Year-Month	 Recyclable pack
 STERILE EO Sterile EO: Ethylene oxide	 Attention, See instruction for use
 LOT Batch Number	
 Dyed, Absorbable, Braided, Coated	2292

IFU-GS-rev-06-15-01-2018

Issue date: 11.09.2012

"EASSI(The European Association of Surgical Suture Industry) has developed a system of symbols which is designed to describe various suture product characteristics in an intuitive, pictorial manner. The use of symbols is permitted by the Medical Device Directive (MDD 93/42/EEC) and enables companies to provide information to the customer without having to provide multilingual translations."

İBOZ TİBBİ MALZEME SANAYİ VE TİCARET A.Ş.
Sağlık Mah. Sakarlık Sokak No:33/5/İhlye Çankaya/ANKARA/TÜRKİYE
Tel: +90 (312) 254 03 40 (5 hat) Faks:+90 (312) 254 03 30
web: www.boztibbi.com e-mail: boz@boztibbi.com

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Boz®

Definition

Synthetic absorbable, coated, mid term wound support and mid term mass absorption.

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













16 mm

1/2 Round Body, Taper Point			GS1990	5/0 1,0	75	12	Violet
			GS1992	5/0 1,0	75	12	Undyed
			GS2503	4/0 1,5	75	12	Undyed
			GS2504	4/0 1,5	75	12	Violet




Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
16 mm							
1/2 Round Body, Taper Point			GS3005	3/0 2,0	75	12	Violet
3/8 Reverse Cutting			GS1674	6/0 0,7	75	12	Violet
			GS2014	5/0 1,0	45	12	Undyed
			GS2540	4/0 1,5	75	12	Violet
			GS3025	3/0 2,0	75	12	Violet
3/8 Inside Cutting			GS2140	5/0 1,0	45	12	Violet
			GS2144	5/0 1,0	75	12	Violet
			GS2557	4/0 1,5	75	12	Undyed
			GS3039	3/0 2,0	75	12	Violet
17 mm							
1/2 Round Body, Taper Point			GS2584	4/0 1,5	75	12	Violet
18 mm							
1/2 Round Body, Taper Point			GS5525	5/0 1,0	75	12	Violet
			GS2613	4/0 1,5	75	12	Violet
			GS3624	2/0 3,0	75	12	Violet
3/8 Round Body, Taper Point			GS3085	3/0 2,0	75	12	Violet
19 mm							
3/8 Reverse Cutting			GS1706	6/0 0,7	75	12	Violet
			GS2640	4/0 1,5	75	12	Violet
3/8 Reverse Cutting Premium			GS2211	5/0 1,0	45	12	Violet


Glikosorb



Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color	
20 mm							
1/2 Round Body, Taper Point			GS2229	5/0 1,0	75	12	Violet
			GS2685	4/0 1,5	75	12	Violet
			GS3142	3/0 2,0	75	12	Violet
1/2 Tapercutting			GS2673	4/0 1,5	75	12	Violet
3/8 Reverse Cutting			GS2250	5/0 1,0	75	12	Undyed
			GS2716	4/0 1,5	45	12	Undyed
			GS2719	4/0 1,5	75	12	Violet
			GS3164	3/0 2,0	45	12	Violet
			GS3168	3/0 2,0	75	12	Violet
			GS3663	2/0 3,0	75	12	Violet
22 mm							
1/2 Round Body, Taper Point			GS2747	4/0 1,5	75	12	Violet
			GS3196	3/0 2,0	75	12	Violet
			GS3671	2/0 3,0	75	12	Violet
			GS4575	1 4,0	75	12	Violet
3/8 Reverse Cutting			GS2763	4/0 1,5	75	12	Violet
			GS3214	3/0 2,0	90	12	Violet

Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
24 mm					
1/2 Reverse Cutting  	GS5524	4/0 1,5	75	12	Violet
	GS5523	3/0 2,0	75	12	Violet
	GS3692	2/0 3,0	75	12	Violet

25 mm					
1/2 Round Body, Taper Point  	GS2788	4/0 1,5	75	12	Violet
	GS3252	3/0 2,0	75	12	Violet
	GS3746	2/0 3,0	100	12	Violet
	GS3724	2/0 3,0	75	12	Violet
1/2 Tapercutting  	GS3718	2/0 3,0	75	12	Violet
3/8 Reverse Cutting  	GS2804	4/0 1,5	75	12	Violet
	GS3757	2/0 3,0	75	12	Violet
3/8 Inside Cutting  	GS2812	4/0 1,5	45	12	Violet
	GS3284	3/0 2,0	45	12	Violet
	GS3763	2/0 3,0	75	12	Violet

26 mm					
1/2 Round Body, Taper Point  	GS2292	5/0 1,0	75	12	Violet
	GS2828	4/0 1,5	75	12	Violet
	GS3313	3/0 2,0	75	12	Violet
	GS3318	3/0 2,0	90	12	Violet
	GS3786	2/0 3,0	75	12	Violet

Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
26 mm					
1/2 Round Body, Taper Point 	GS4215	0 3,5	75	12	Violet
	GS4592	1 4,0	90	12	Violet
1/2 Reverse Cutting 	GS4230	0 3,5	75	12	Violet
3/8 Round Body, Taper Point 	GS2841	4/0 1,5	75	12	Violet
3/8 Reverse Cutting 	GS3345	3/0 2,0	75	12	Undyed
	GS3346	3/0 2,0	75	12	Violet
	GS3821	2/0 3,0	75	12	Undyed
	GS3822	2/0 3,0	75	12	Violet
30 mm					
1/2 Round Body, Taper Point 	GS3369	3/0 2,0	75	12	Violet
	GS3374	3/0 2,0	90	12	Violet
	GS3858	2/0 3,0	75	12	Violet
	GS4244	0 3,5	75	12	Violet
	GS4608	1 4,0	75	12	Violet
1/2 Reverse Cutting 	GS3882	2/0 3,0	75	12	Violet
	GS4266	0 3,5	75	12	Violet
1/2 Tapercutting 	GS3843	2/0 3,0	75	12	Violet


Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm						
3/8 Reverse Cutting 		GS2887	4/0 1,5	75	12	Violet
		GS3397	3/0 2,0	75	12	Violet
		GS3891	2/0 3,0	75	12	Violet
		GS4273	0 3,5	75	12	Violet
32 mm						
1/2 Round Body, Taper Point 		GS5508	3/0 2,0	75	12	Violet
35 mm						
1/2 Round Body, Taper Point 		GS3411	3/0 2,0	75	12	Violet
		GS3916	2/0 3,0	75	12	Violet
		GS3922	2/0 3,0	90	12	Violet
		GS4300	0 3,5	75	12	Violet
		GS4644	1 4,0	75	12	Violet
3/8 Reverse Cutting 		GS3423	3/0 2,0	75	12	Violet
		GS3933	2/0 3,0	75	12	Violet
36 mm						
1/2 Round Body, Taper Point 		GS5509	3/0 2,0	75	12	Violet
1/2 Reverse Cutting 		GS4675	1 4,0	75	12	Violet





Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
37 mm					
1/2 Round Body, Taper Point 	GS3950	2/0 3,0	75	12	Violet
	GS4337	0 3,5	75	12	Violet
	GS4681	1 4,0	75	12	Violet
1/2 Reverse Cutting 	GS5507	2/0 3,0	75	12	Violet
	GS4354	0 3,5	90	12	Violet
	GS4994	2 5,0	75	12	Violet
40 mm					
1/2 Round Body, Taper Point 	GS3982	2/0 3,0	75	12	Violet
	GS4368	0 3,5	75	12	Violet
	GS4375	0 3,5	90	12	Violet
	GS4713	1 4,0	75	12	Violet
	GS4719	1 4,0	90	12	Violet
	GS5003	2 5,0	75	12	Violet
	GS5006	2 5,0	90	12	Violet
1/2 Reverse Cutting 	GS4397	0 3,5	75	12	Violet
	GS4745	1 4,0	75	12	Violet
	GS4749	1 4,0	90	12	Violet
	GS5019	2 5,0	75	12	Violet





Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm						
1/2 Reverse Cutting		GS5023	2 5,0	90	12	Violet
1/2 Tapercutting		GS4363	0 3,5	75	12	Violet
3/8 Reverse Cutting		GS4756	1 4,0	75	12	Violet
		GS4408	0 3,5	75	12	Violet
45 mm						
1/2 Round Body,Taper Point		GS4422	0 3,5	75	12	Violet
		GS4772	1 4,0	75	12	Violet
		GS5035	2 5,0	75	12	Violet
		GS5040	2 5,0	90	12	Violet
1/2 Tapercutting		GS4766	1 4,0	75	12	Violet
		GS5738	2 5,0	75	12	Violet
3/8 Reverse Cutting		GS4432	0 3,5	75	12	Violet
		GS4787	1 4,0	75	12	Violet
Straight Reverse Cutting		GS2909	4/0 1,5	45	12	Violet
48 mm						
1/2 Round Body,Taper Point		GS4796	1 4,0	75	12	Violet

Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
50 mm					
1/2 Round Body, Taper Point 	GS4454	0 3,5	75	12	Violet
	GS4459	0 3,5	90	12	Violet
	GS4812	1 4,0	75	12	Violet
	GS4818	1 4,0	90	12	Violet
	GS5066	2 5,0	75	12	Violet
	GS5070	2 5,0	90	12	Violet
1/2 Reverse Cutting 	GS4471	0 3,5	75	12	Violet
60 mm					
Straight Reverse Cutting 	GS3492	3/0 2,0	75	12	Violet
	GS4074	2/0 3,0	75	12	Violet
63 mm					
3/8 Blunt 	GS4858	1 4,0	75	12	Violet
64 mm					
3/8 Blunt 	GS4500	0 3,5	90	12	Violet
	GS4861	1 4,0	90	12	Violet

Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
6x45 cm					
GLIKOSORB Non-Needle Suture	GS2950	3/0 2,0	6x45	12	Violet
	GS3543	2/0 3,0	6x45	12	Violet
	GS4137	0 3,5	6x45	12	Violet
	GS4563	1 4,0	6x45	12	Violet
10x50 cm					
GLIKOSORB Non-Needle Suture	GS2954	3/0 2,0	10x50	12	Violet
	GS3548	2/0 3,0	10x75	12	Violet
12x45 cm					
GLIKOSORB Non-Needle Suture	GS2344	4/0 1,5	12x45	12	Violet
	GS2958	3/0 2,0	12x45	12	Violet
	GS3552	2/0 3,0	12x45	12	Violet
	GS4145	0 3,5	12x45	12	Violet
	GS4554	1 4,0	12x45	12	Violet
17x45 cm					
GLIKOSORB Non-Needle Suture	GS2966	3/0 2,0	17x45	12	Violet
	GS3558	2/0 3,0	17x45	12	Violet
150 cm					
GLIKOSORB Non-Needle Suture	GS1737	5/0 1,0	150	12	Violet
	GS2326	4/0 1,5	150	12	Violet
	GS2935	3/0 2,0	150	12	Violet

Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
150 cm					
GLIKOSORB Non-Needle Suture	GS3511	2/0 3,0	150	12	Violet
	GS4114	0 3,5	150	12	Violet
	GS4532	1 4,0	150	12	Violet

250 cm					
GLIKOSORB Non-Needle Suture	GS2332	4/0 1,5	250	12	Violet
	GS2942	3/0 2,0	250	12	Violet
	GS3518	2/0 3,0	250	12	Violet
	GS4121	0 3,5	250	12	Violet
	GS4541	1 4,0	250	12	Violet