U11 Plus GL
LabStrip
URINALYSIS

IVD

Urinalysis

Test strip for semi-quantitative urine analysis from fresh urine. Use only with **LabUMat 2** automated urine chemistry analyzer or for visual reading. The **LabStrip U11 Plus GL** is an *in vitro* diagnostic medical device for professional laboratory use only in conformity with the Directive 98/79/EC.

150 pcs urine test strips for rapid determination of Bilirubin, Urobilinogen, Ketones (Acetoacetic Acid), Ascorbic acid, Glucose, Protein (Albumin), Blood, pH value, Nitrite, Leucocytes and Specific Gravity in urine. Refer to the carton and label for specific parameter combination on the product you are using.

Summary and Explanation:

Screening test for recognition of liver disease, biliary and hepatic obstructions, diabetes, and haemolytic diseases, urological, and nephrological diseases associated with haematuria or haemoglobinuria, diseases of the kidneys and urinary tract, pathological shifts in the pH value, as well as for investigation of the sediment.

Clinical Utility:

Bilirubin: Intended to measure the levels of bilirubin conjugates in urine. Measurements of urinary bilirubin and conjugates are used in the diagnosis and treatment of certain liver and bile diseases.

Urobilinogen: Intended to detect and estimate urobilinogen (a bile pigment degradation product of red blood cell haemoglobin) in urine. Estimations obtained by this device are used in the diagnosis and treatment of liver diseases and haemolytic (red blood cells) disorders.

Ketones: Intended to detect ketones in urine. Nitrite identification is used for in the diagnosis and treatment of acidosis (a condition characterized by abnormally high acidity of body fluids) or ketosis (a condition characterized by increased production of ketone bodies) and for monitoring patients with diabetes.

Ascorbic acid: Intended to measure the level of ascorbic acid (vitamin C) in urine.

Glucose: Intended to measure glucosuria (glucose in urine). Urinary glucose measurements are used in the diagnosis and treatment of carbohydrate metabolism disorders including diabetes mellitus and hyperglycemia.

Protein: Intended to identify protein in urine. Identification of urinary protein is used in diagnosis and treatment of renal diseases.

Blood: Intended to detect occult blood in urine. Occult blood indicates serious urological or kidney diseases.

pH: Intended to estimate the pH of urine. Estimations of pH are used to evaluate the acidity or alkalinity of urine as it relates to numerous renal and metabolic disorders and in the monitoring of patients with certain diets. Persisting high pH values indicate urinary tract infections.

Nitrite: Intended to identify nitrite in urine. Nitrite identification is used in the diagnosis and treatment of urinary tract infections of bacterial origin.

Leucocytes: Intended to detect leucocytes in urine. Leucocytes indicate inflammatory diseases of the kidney and the urinary tract, and suggest need of further investigation.

Specific Gravity: Intended to provide an estimation of renal ability of urine concentration or urine dilution. The specific gravity of urine varies in accordance with the drinking quantity as well as different disorders.

Principle of the Procedure:

Bilirubin: A red azo compound is obtained in the presence of acid by coupling of bilirubin with diazonium salt.

Urobilinogen: The test is based on the coupling of urobilinogen with a stabilised diazonium salt to a red azo compound.

Ketones: Acetone and acetoacetic acid react with sodium nitroprusside in alkaline solution to give violet coloured complex (Legal's test).

Ascorbic acid: The detection is based on the decolouration of Tillman's reagent.

Glucose: The detection is based on the glucoseoxidase-peroxidase-chromogen reaction. Apart from glucose no other compound in urine is known to give a positive reaction.

Protein: The test is based on the „protein error“ principle of the indicator. The test is especially sensitive in the presence of albumin. Other proteins are indicated with less sensitivity.

Blood: The detection is based on the pseudoperoxidative activity of haemoglobin and myoglobin, which catalyze the oxidation of an indicator by organic hydroperoxide and chromogene producing a green colour.

pH: The test paper contains indicators which clearly change colour between pH 5 and pH 9 (from orange to green to turquoise).

Nitrite: The colour test is based on the principle of Griess reaction. Any degree of pink-orange colouration should be interpreted as a positive nitrite test suggestive of $\geq 10^3$ organisms/ml urine.

Leucocytes: The test is based on the esterase activity of granulocytes. This enzyme splits heterocyclic carboxylates. The component released reacts with a diazonium salt producing a violet colour.

Specific Gravity: The test is based on a colour change of the reagent from blue-green to greenish yellow depending on the concentration of ions in the urine. The test permits the determination of urine density between 1.000 and 1.030.

Kit Components:

Each kit contains everything needed to perform 150 tests:

- 150 pcs **LabStrip U11 Plus GL** test strips,
- Label with colour scale for visual reading,
- 1 pc Calibration card for checking and setting of **LabUMat 2** automated urine chemistry analyzer,
- 1 pc Instructions for use.

Other required appliances for urine analysis:

- **LabUMat 2** automated urine chemistry analyzer with instructions for use.
- Clean, detergent free and dry container for urine collection.

Reagents Composition:

The reagents in the individual test fields are formulated to contain:

| | | |
|--------------------------|--------------------------------------|--------|
| Bilirubin: | Diazonium salt | 3.1 % |
| Urobilinogen: | Diazonium salt | 3.6 % |
| Ketones: | Sodium nitroprusside | 2.0 % |
| Ascorbic acid: | 2,6-dichloro-phenol-indophenol | 0.7 % |
| Glucose: | Glucose oxidase | 2.1 % |
| | Peroxidase | 0.9 % |
| | O-Tolidine hydrochloride | 5.0 % |
| Protein: | Tetra-bromophenol blue | 0.2 % |
| Blood: | Isopropylbenzol-hydroperoxide | 21.0 % |
| | Tetramethylbenzidine-dihydrochloride | 2.0 % |
| pH: | Bromthymol blue | 10.0 % |
| | Methyl red | 2.0 % |
| Nitrite: | Sulfanilic acid | 1.9 % |
| | Tetrahydrobenzol[h]quinolon-3-ol | 1.5 % |
| Leucocytes: | Carboxylic acid ester | 0.4 % |
| | Diazonium salt | 0.2 % |
| Specific Gravity: | Bromthymol blue | 2.8 % |

Concentrations given are based on reagent composition (w/w) at time of manufacture and may vary within manufacturing tolerances.

Warning!

- Every item in the package can be handled as household waste. As reactive materials are present at very low quantity, the product does not come under the scope of the relevant EU regulations for dangerous materials.
- Keep away from swallowing, touching with skin or mucous membranes of chemicals!
- For in vitro diagnostic use only!
- The vial cap containing a non toxic, molecule-sieve based drying agent blocks the test strips from air-humidity. In case of swallowing the drying agent chemicals, drink substantial amount of liquid.
- If you have any questions, please turn to your local distributor!

Arrangements for urine analysis:

Warning!

Use only **LabUMat 2** automated urine chemistry analyzer for **LabStrip U11 Plus GL** test strip urine analysis.

If you open a new **LabStrip U11 Plus GL** test strip package you can find a calibration card for checking and setting of **LabUMat 2** automated urine chemistry analyzer. You can find the instructions for settings in the user manual of the meter. Adjustment of the meter is not required until the start of next strip vial.

Follow the user manual of **LabUMat 2** automated urine chemistry analyzer.

Specimen Collection and Preparation:

- Collect urine in a clean, dry container which allows complete immersion of all the fields on the test strip.
- Do not add preservatives.
- Test the specimen as soon as possible, with the sample well mixed but not centrifuged.
- The use of fresh morning urine is recommended for optimal nitrite tests, as well as for the valid determination of bilirubin and urobilinogen, since these compounds are unstable when exposed to light and room temperature (+15 to +25 °C).
- If immediate testing is not possible, the sample should be stored in the refrigerator (+2 to +8 °C) and then brought to room temperature (+15 to +25 °C) before used in the test.
- Non-preserved urine at room temperature may undergo pH changes due to microbial proliferation, which may interfere with protein determination.
- If cleanly voided specimens are not collected from females, positive results for leucocytes may be found due to contamination from outside the urinary tract.
- Skin cleansers containing chlorhexidine may affect protein test results if specimen positive contamination occurs.

Procedure and Notes:

- Use only fresh, mixed, not centrifuged urine. First morning urine is recommended. Perform the urine analysis in 4 hours after sample collection! Keep urine away from light.
- Collect specimen in clean, rinsed containers, free of detergents. Do not use preservatives!
- After removing the required number of strips, close immediately the container securely with the vial cap containing drying agent.
- Do not touch test areas of the reagent strip.

In case of instrumental reading

- Read carefully the instructions for use of **LabUMat 2** automated urine chemistry analyzer.
- In case of reading with **LabUMat 2** automated urine chemistry analyzer, the meter feeds the strips and the samples are pipetted automatically. Then the result is displayed after 1 minute incubation period.

In case of visual reading

- The reagents have two handles. Take the end that is closer to the tube cap in the original package. You hold it right if the dark (green) SG pad is closest to your hand.
- Immerse the test strip into the urine for approx. 2 sec, so that all reagent areas are covered.
- Remove excess urine from the strip by wiping the edge of the strip on the edge of the urine container or on absorbent paper.
- To prevent interference from adjacent test areas, incubate the strip in horizontal position.
- Compare the reagent areas on the strip with the corresponding color charts on the container about 60 sec (between 60-120 sec for the Leucocyte test) after immersion. Do not read later than 2 minutes.
- Colours on the colour chart are representing the nominal values of the test fields. Actual results are located around nominal values.
- The white field between the specific weight and leucocyte test field is for instrument measurements and is used to compensate for the intrinsic colour of urine.

Warning!

- Do not take the strip out of the meter during reading procedure.
- Prior to measuring always make sure that the process is performed according to the instructions for use of the **LabUMat 2** automated urine chemistry analyzer.
- Do not perform urine analysis at temperatures below +15 °C or above +35 °C.
- Strips must be kept away from heat and direct sunlight.
- Do NOT ever reuse test strips. After measurement completed eliminate the strip carefully.
- Until usage, store the test strips in original packages. Strips in each vial should not be mixed.
- Diagnoses and therapies can not be derived from one single test result only, instead should be based on all available medical diagnoses.
- Never use the strip if more than 5 minutes spent from the moment of its removing from the vial.

Biological risk!

Handle all specimens and strips as if they contained infectious agents. When the assay procedure is completed, dispose of specimens and strips carefully follow the relevant, local instructions.

- In rare occasions, the varying test conditions, due to the heterogeneity of the different urine (for reasons of different levels of activators, inhibitors, or different ion concentrations) may cause variation in the intensity and contrast of the colors.
- Not all cases of interference with every component of any medicines are known. The color reaction of the pads might change. We, therefore, recommend another test at the end of any medication with drugs.
- Always follow the general working instructions for laboratories as well.
- The test strips do NOT contain toxic materials!

Results:

Results are determined visually by direct comparison of reacted test fields with the colour chart on the container label. Visual colour charts represent nominal test values for each test field-actual values may vary around the nominal values.

The leucocyte and blood (erythrocyte) tests are not quantitative determinations, but serve as screening methods for the presence of Leucocytes and blood (erythrocytes) in urine. Microscopic examination of specimens with a positive Leucocyte or blood test result should be performed if quantitative results are required.

Ascorbic acid may interfere with the glucose, nitrite, bilirubin and blood test results (see Limitations below). If a positive Ascorbic acid result is found, either repeat the test at least 10 hours one day after discontinuation of Vitamin C administration or use a photometric test unaffected by Ascorbic acid.

When using **LabUMat 2** automated urine chemistry analyzer please refer the user manual of those instrument.

System operation:

Each test strip has 11 measuring zones. These zones contain sensitive reagents. Colour of the test field is changing as a result of chemical reaction of urine contact. The **LabUMat 2** automated urine chemistry analyzer detects the colouration and displays the result.

Limitations of the Procedure:

Note: Diagnostic or therapeutic decisions should not be based on any single result or method.

Bilirubin: The reaction is unaffected by pH of urine. False low or negative results may be simulated by large amounts of vitamine C or nitrite or by longer exposure of the sample to direct light. Increased concentration of urobilinogen can reinforce the sensitivity of the field. Different urine constituents (e.g. urine indicane) can lead to atypical coloration. For metabolites of drugs see urobilinogen.

Urobilinogen: The reaction is unaffected by pH of urine. Higher concentration of formaldehyde or exposure of the urine to light for a longer period of time may lead to lowered or falsely negative results. Beetroot (excreted pigments) or metabolites of drugs which give a colour at low pH (phenazopyridine, azo dyes, p-aminobenzoic acid or other medicaments which have a red intrinsic coloration in acidic medium) may produce false positive results. Prolonged exposure to light is to be avoided.

Ketones: Phthalein compounds and derivatives of anthrachinone interfere by producing a red coloration in the alkaline range which may mask the coloration of ketones.

Ascorbic acid: As ascorbic acid already in low concentrations can disturb various test fields, especially the glucose assay in low concentrations, the test must be repeated if the ascorbic acid reaction is positive, however, at the earliest 10 hours after the last vitamin C intake (medication, fruit and vegetables).

Glucose: High concentrations of ascorbic acid in urines with a low glucose concentration (up to 100 mg/dl (5.5 mmol/l)) may inhibit the reaction and lead to lower or false negative results. Repeat the test 10 hours one day after stopping the intake of vitamin C. Pay attention to the ascorbic acid field. In addition an inhibitory effect is produced by gentisic acid, a pH value of <5 and high specific gravity. False positive reactions can also be produced by a residue of peroxide containing cleansing agents or others.

Protein (albumin): Falsely positive results are possible in high alkaline urine samples (pH >9) and in the presence of high specific gravity, after infusions with polyvinylpyrrolidone (blood substitute) after intake of medicaments containing quinine and also by disinfectant residues containing quaternary ammonium groups in the urine sampling vessel.

Blood: Microhaematuria does not affect the colour of urine and is only detectable by microscopic or chemical tests. From a level approx. 25 Ery/ μ l and above, even at high concentrations of ascorbic acid normally no negative results are observed. Falsely positive reactions can also be produced by a residue of peroxide containing cleansing agents, activities of microbial oxidase due to infections of the urogenital tract or by formaline. For establishing an individual diagnosis, it is therefore indispensable to take into consideration also the clinical manifestations.

The number of erythrocytes which are detected by sediment analysis may be lower than the result of the test strip, because lysed cells are not detected by sediment analysis.

pH: No interferences are known.

Nitrite: Before testing the patient should ingest vegetable-rich meals, reduce fluid intake and discontinue antibiotic and vitamin C therapy 3 days prior to the test. False positive results may occur in stale urine samples, in which nitrite has been formed by contamination of the specimen and in urines containing dyes (derivatives of pyridinium, beetroot).

A negative result even in the presence of bacteriuria can have the following reasons: bacteria not containing nitrate reductase, antibiotic treatment diet with low nitrate content, high diuresis, high content of ascorbic acid or insufficient incubation of the urine in the bladder.

Leucocytes: Strongly coloured compounds (e.g. nitrofurantoin) may disturb the colour of the reaction. High concentrations of glucose, oxalic acid, drugs containing cephalixin, cephalothine or tetracycline can lead to weakened reaction. False-positive reactions may be caused by contamination of vaginal secretion. The number of leucocytes which are detected by sediment analysis may be lower than the result of the strip, because lysed cells are not detected by sediment analysis. Partial cytolysis intensifies the colour response, particularly in the region of the maximum analytical sensitivity. Leucocyte esterase results may be positive in the absence of observable cells if the leucocytes have lysed. False-positive reactions may be caused by formaldehyde (preservative). Protein concentrations above 5 g/l or a high specific gravity may diminish the colour response. Bacteria, trichomonas and erythrocytes do not, however, react with the test field.

Specific Gravity: Highly acidic urines (pH <6) yield slightly elevated results whereas highly alkaline urines (pH >8) yield diminished results. Glucose and urea do not interfere.

Expected values:

Bilirubin: Normally, no bilirubin is detectable in urine. Concentration of 0.5 mg/dl and more lead to a colour of red-orange peach and indicate the early stage of a liver disease. The colour fields correspond to the following bilirubin concentrations: neg. (negative), 1 (+), 3 (++), 6 (+++) mg/dl or neg. (negative), 17 (+), 50 (++) , 100 (+++) μ mol/l. Values of 0.5-1 mg/dl bilirubin are indicated.

Urobilinogen: The normal concentration of urobilinogen in urine goes from 0.1-1.8 mg/dl (1.7-30 μ mol/l). Concentrations of >2 mg/dl (35 μ mol/l) are considered to be pathological. Absolute absence of urobilinogen in the urine, which is also pathological, cannot be demonstrated using test strips. The colour fields correspond to the following urobilinogen concentrations: norm. (normal), 2 (+), 4 (++) , 8 (+++) , 12 (++++) mg/dl or norm. (normal), 35 (+), 70 (++) , 140 (+++) , 200 (++++) μ mol/l.

Ketones: Normally, the urine is free of ketones. Detectable concentrations of ketones can originate from physiological stress (fasting, pregnancy, excessive sport). Phenylketones in higher concentrations will produce variable colours. β -Hydroxybutyric acid is not detected. The colour fields correspond to the following acetoacetic acid values: neg. (negative), 15 (+), 50 (++) , 150 (+++) mg/dl or neg. (negative), 1.5 (+), 5 (++) , 15 (++) mmol/l. Values of 5 mg/dl acetoacetic acid or 50 mg/dl acetone are indicated.

Ascorbic acid: In the presence of ascorbic acid a colour change takes place from grey blue to orange. The colour fields correspond to the following ascorbic acid concentrations: neg. (negative), 20 (+), 40 (++) mg/dl or neg. (negative), 1.14 (+), 2.28 (++) mmol/l. Values of 5-10 mg/dl or 0.6-1.1 mmol/l are indicated.

Glucose: Normally, glucose cannot be detected in the urine although small amounts are secreted also by the healthy kidney. Changes in the colorations less than 50 mg/dl (2.8 mmol/l) are to be considered normal. The colour fields correspond to the following ranges of glucose concentrations: norm. (normal), 50 (+), 150 (++) , 500 (+++) , 1000 (++++) mg/dl or norm. (normal), 2.8 (+), 8 (++) , 28 (+++) , 56 (++++) mmol/l. Values of 40 mg/dl glucose are indicated.

Protein (albumin): Normally, no protein is detectable in the urine of healthy subjects. Colour intensities for the protein (albumin) field which match the intensity of the 0.3 g/l colour field are to be classed as pathological. The colour fields correspond to the following ranges of albumin concentrations: neg. (negative), 30 (+), 100 (++) , 500 (+++) mg/dl or neg. (negative), 0.3 (+), 1.0 (++) , 5.0 (++) g/l. Values of approx. 15 mg/dl albumine are indicated.

Blood: Intact erythrocytes are reported by punctual colorations on the test pad, haemoglobin and myoglobin are reported by a homogeneous green coloration. The color fields correspond to the following values: neg. (negative), approx. 5-10 (+), approx. 50 (++) , approx. 300 (++) Ery/ μ l. Values of approx. 5 erythrocytes/ μ l are indicated.

pH: Normal fresh urine of healthy people varies between 5 to 6 pH values. The colour fields correspond to the following pH values: 5, 6, 7, 8, 9.

Nitrite: Negative results do not exclude significant bacteriuria (insufficient incubation, urinary tract infections due to bacteria not containing nitrate reductase). Red or blue borders or edges which may be present must not be interpreted as a positive result. Values of 0.05-0.1 mg/dl nitrite are indicated.

Leucocytes: Urines of healthy subjects do not contain any leucocytes. Positive results, even when constantly varying from negative to 25, are to be considered as clinically relevant. Discolorations which no longer match the negative colour field are to be classed as positive. The colour fields correspond to the following values: neg. (negative), approx.25 (+), approx. 75 (++) , approx. 500 (++) leucocytes/ μ l. Values of 10-20 leucocytes/ μ l are indicated.

Specific Gravity: Normal fresh urine of healthy people varies between 1.015 to 1.025 values. The colour scale has optimized at a pH of the urine of 6. The colour fields correspond to the values 1.000, 1.005, 1.010, 1.015, 1.020, 1.025, 1.030.

Measuring Ranges:

| Parameters | neg. | trace* | + | ++ | +++ | ++++ | | | |
|---------------------------|-------|--------|----------|--------|---------|------|------|-----|---|
| Bilirubin (mg/dL) | neg. | 0,5 | 1 | 3 | 6 | | | | |
| Urobilinogen (mg/dL) | norm. | - | 2 | 4 | 8 | 12 | | | |
| Ketones (mg/dL) | neg. | 5 | 15 | 50 | 150 | | | | |
| Ascorbid acid* (mg/dL) | neg. | - | 20 | 40 | 100 | | | | |
| Glucose (mg/dL) | norm. | 25 | 50 | 150 | 500 | 1000 | | | |
| Protein (mg/dL) | neg. | 15 | 30 | 100 | 500 | | | | |
| Blood (Ery/ μ L) | neg. | - | ca. 5-10 | ca. 50 | ca. 300 | | | | |
| pH* | 5 | 5,5 | 6 | 6,5 | 7 | 7,5 | 8 | 8,5 | 9 |
| Nitrite | neg. | - | pos. | | | | | | |
| Leukocytes (Leu/ μ L) | neg. | - | ca. 25 | ca. 75 | ca. 500 | | | | |
| Compensation filed | | | | | | | | | |
| Specific Gravity | 1000 | 1005 | 1010 | 1015 | 1020 | 1025 | 1030 | | |

* **LabUMat2** is able to detect trace categories in case of Bilirubin, Ketones, Glucose and Protein, and half unit in case of pH and 100 mg/dL in case of Asc. acid.

Storage and Stability:

Keep diagnostic test strips in tightly closed original tubes in a dry, dark and cool place (between +2 to +30 °C). After removing the required number of strips, close immediately the container securely. Do not remove the desiccant from the original cap

The strips must be kept away from moisture, direct sunlight, elevated temperature and chemical fumes. Under proper conditions test strips are stable up to the stated expiry date even after opening. Do not touch the test pads.

System quality control:

Each laboratory should evaluate according to its own working standards.

During quality control, the legal requirements and guidelines shall be observed regarding frequency of tests, target values and target ranges, as well as the documentation of the results.

Corrective actions shall be taken if the results are out of the target range.

For periodically checking your strips and system we offer the QC the Dipper (Quantimetrix Corporation, REF: 1440-01), the Dropper (Quantimetrix Corporation, REF:1440-02) or DipAndSpin (Quantimetrix Corporation, REF:1470-01) urine dipstick control solutions.

Dip the test strip into control solution instead of urine sample!

Read the chapter „Checking of the system“ in the instructions for use of the urine analyzer.

Specific Performance Characteristics:

The performance characteristics of **LabStrip U11 Plus GL** test strips are based on both clinical and analytical studies. Sensitivity is dependent upon the colour perception of the reader, the presence or absence of interfering specimens, and the lighting conditions for visual reading.

Each colour block on the chart corresponds to a range of analyte concentrations.

Bilirubin: In 90 % of urines tested, bilirubin concentrations of 0.5 mg/dl produced a positive result. After a longer reaction time an unspecific yellow colour may develop, which can create positive interference.

Urobilinogen: Based on the work of Kutter [10], a concentration of 1 mg/dl urobilinogen will yield a positive result. The test is sufficiently sensitive so that normal specimens will produce a slight pink colour.

Ketones: In 90 % of urines tested, acetoacetic acid at 8 mg/dl produced a positive result. The test field reacts less sensitively with acetone. Hydroxybutyric acid is not detected.

Ascorbic acid: In 90 % of urines tested, ascorbic acid at 20 mg/dl yielded a positive result.

Glucose: The maximum sensitivity is 20 mg/dl. The test field response is adjusted so that pathological glucose concentrations of 30 mg/dl (Fine) [11] are recognized. Sugars other than glucose and other reducing substances do not react in this test. Possible interference by ascorbic acid can be detected by the adjoining test field which reacts with the ascorbic acid.

Protein: In 90 % of urines tested, albumin concentrations of 12 mg/dl produced a positive result. The test is more sensitive to albumin than to globulin, Bence-Jones proteins and mucoproteins. A negative result does not exclude the presence of these other proteins.

Blood: The test permits the differentiation of intact erythrocytes from hemoglobin or myoglobin. Erythrocytes react as diffuse patches on the field. The practical sensitivity of the test lies between 5 and 10 Ery/ μ l.

A study of 625 fresh urine specimens, comparing results with those using another test strip for blood demonstrated a clinical specificity of 90.2 % and sensitivity of 81 %.

pH: pH values are determined to within 1 unit over the range from 5 to 9. Readings are not affected by variations in the urine buffer concentration.

Nitrite: The maximum sensitivity is 0.05 mg/dl, which is equivalent to about 100,000 bacteria/ml. In early morning urine, 90 % of all infections are detected by the nitrite test. Although most uropathogenic bacteria are able to reduce nitrate to nitrite (e. g. Klebsiella, E. coli, Proteus, Aerobacter, Citrobacter, etc.), results depend on the number of bacteria, the nitrate content, and the retention time of the urine.

Leucocytes: In 90 % of urines tested, concentrations of 20 leucocytes/ μ l produced positive results. Any pink coloration of the test field should be considered clinically significant. When result from this method were compared with those using another test strip for leucocytes for 822 fresh urines, clinical specificity of 80 % and sensitivity of 89.2 % were determined.

Specific Gravity: In 86 % of 102 urines tested, the finds based on the colour chart were found to be within the acceptable range of + or – one colour increment compared to the findings of the reference refractometer.

Intra-assay

Within run precision was determined by using 10 replicas of two levels (normal, abnormal) of control urine. The negative and positive values were correctly identified 100 % of time for all the parameters.

Inter-assay

Between run precision was determined by using two levels (normal, abnormal) of control urine in 10 independent assays and with three different lots of reagent over a 6 months period. The negative and positive values were correctly identified 100 % of time for all the parameters.

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Literature/Literatur

- [1] **Legal, E. A.:** New Acetone Reaction and its Applicability for the Examination of Urine. Chem. Centr. 15: 652 (1983)
- [2] **Chertack, M. und Sherrick, J.:** Evaluation of Nitroprusside Dip Test for Ketone Bodies. J. A. M. A. 167: 1621 (1958)
- [3] **Roe, J. H.:** Chemical Determination of ascorbic, dehydroascorbic and diketogulonic Acids. Methods of Biochemical Analysis, Vol 1: 115 (1954) ed. by d. Glick, Interscience Publisher, New York
- [4] **Comer, J.:** Semiquantitative Specific Test Paper for Glucose in Urine. Anal. Chem. 28: 1748 (1956)
- [5] **Appel, W., Nurck, C. und Merkle, U.:** A Rapid Test for Urinary glucose with an Ascorbic Acid Zone. Medical laboratory 6: 29–39 (1979)
- [6] **Sorenson, S.:** The Measurement of the Hydrogen Ion Concentration and Its Importance for Enzymatic process. Biochem. Z. 21: 131 (1909)
- [7] **Vonderschmitt, D. und Scholer, A.:** Teststreifen für Screening-Untersuchungen zum semiquantitativen Nachweis von Proteinurinen. J. Clin. Chem. Biochem. 19: 997 (1981)
- [8] **Leonards, J.:** Simple Test for Hematuria compared with Established Tests. J. A. M. A. 179: 807 (1962)
- [9] **Weltmann, O.:** Method for the Simple Detection of Urinary Tract Infections. Wien. Med. Wschr. 72: 618 (1922)
- [10] **Kutter, D. und Humbel, R.:** Quantitative Assay of Urinary Urobilinogen with p-Methoxybenzene Diazoniumfluoroborate. Clin. chim. Acta 45: 61–66 (1922)
- [11] **Fine, J.:** Glucose Content of Normal Urine. Brit. Med. J. 1: 1209–1214 (1965)

ANA-9901GL-1

Manufacturer/Hergestellt von:

77 ELEKTRONIKA Kft.
1116 Budapest, Fehérvári út 98.
HUNGARY/UNGARN
Tel: + 36 (1) 2061480
Fax: + 36 (1) 2061481
E-mail: sales@e77.hu
www.e77.hu

Markings/Symbole

In vitro diagnostic medical device
In vitro Diagnostikum

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