Generic Name : Ferritin Kit (Particle-enhanced Immunoturbidimetric Assay Method)

Abbreviated name: FER

Order Information

Cat. No.	Package size
FER0202	R1 2×18 mL+ R2 2×10 mL
FER0203	R1 1×12 mL + R2 1×7 mL
FER0204	R1 1×20 mL + R2 1×12 mL
FER0302	R1 2×18 mL+ R2 2×10 mL+ Cal 4×2 mL
FER0303	R1 1×12 mL + R2 1×7 mL + Cal 4×2 mL
FER0304	R1 1×20 mL + R2 1×12 mL + Cal 4×2 mL
FER1302	R1 1×12 mL + R2 1×7 mL + Cal 4×2 mL

Intended use

In vitro test for the quantitative determination of FER concentration in serum or plasma on photometric systems.

Summary

Ferritin is the iron storage protein with a molecular weight of approximately 450KD. Ferritin is mainly detectable in liver, spleen and bone marrow and is also found in human serum with small amounts. Serum ferritin concentration usually reflects body iron stores and is considered as one of the most reliable indicators of iron status of patients.

Whereas elevated concentration of ferritin may occur due to cell destruction, liver disease, tumor cell production and chronic inflammation, low level ferritin is always indicative of an iron deficiency.

Method

Particle-enhanced Immunoturbidimetric Assay Method

Reaction Principle

When Latex bounded with anti-ferritin antibodies reacts with the antigen in the sample, agglutination occurs. This agglutination is detected as an absorbance change which is proportional to ferritin concentration in the sample. The actual concentration is then determined by a calibration curve prepared from calibrators of known concentration.

Reagents

Components and concentration

R1	Tris buffer	20 mmol/L
	Preservative	0.5 g/L
R2	Latex coated with anti-human FER antibody	0.15w/v%
	Preservative	1 g/L
Cal	Human FER	/

Warnings and precautions

1. For in vitro diagnostic use only.
2. All products derived from human source have been tested HBsAg, HIV antibody, and HCV antibody negative. However, as no testing method can rule out the potential risk of infection with absolutely certainty, the material should be considered potential infectious.
3. Take appropriate precautions for the use of laboratory reagents. Do not swallow and Avoid contact with skin and mucous membranes.
4. When the reagent accidentally enters the eyes and mouth, or contact with the skin, immediately wash with plenty of water. If necessary, visit the doctor for further medical treatment.
5. Disposal of all waste material should be in accordance with local guidelines.
6. Do not allow to mix reagents with different lots.
7. Avoid direct exposure to sunlight and freezing. The results can not be assured when stored at inappropriate condition.
8. Do not use the product after its expiration date, as the reliability of the measured values cannot be assured after this date.
9. Material safety data sheet is available for professional user on request.

Reagent Preparation

R1: Ready to use

R2: Ready to use

Calibrator Preparation

Calibrators are ready to use.

Storage and stability

Up to expiration date indicated on the label, when stored unopened at 2-8°C and protected from light.

Once opened, the R1 and R2 reagent is stable for 28 days when refrigerated on the analyzer or refrigerator.

Once opened, the calibrator is stable for 14 days at 2-8°C and protected from light.

Contamination of the reagents must be avoided.

Do not freeze the reagents.

Reagent blank absorbency

The absorbance of FER reagent blank at 570 nm should be >0.5 A.

Materials required but not provided

1. Calibrators: FER Calibrator, Cat.No.OSF181 or within FER kit.
2. Controls: FER Control, Cat.No.OQM161.
3. NaCl solution 9 g/L.
4. General laboratory equipments.

Specimen collection and preparation

a) Human Serum or heparinized plasma is suitable for samples.

b)Samples Stability:

7 days at 20-25°C

7 days at 4-8°C

1 year at -20°C (only freeze once)

Assay procedure

	Blank	Sample
Reagent 1	100 µL	100 µL
Dist. water	5 µL	—
Sample	—	5 µL
Mix and incubate at 37°C 4 min, then add:		
Reagent 2	50 µL	50 µL
Mix thoroughly at 37°C, and read the absorbance (A Blank), then read the absorbance (A Sample) again 3 min. later.		
$\Delta A = [A \text{ Sample}] - [A \text{ Blank}]$		

Calibration

1. It is recommended to use Mindray calibrators and 9 g/L NaCl for multipoint calibration.
2. Traceability of the calibrator: This method has been standardized against WHO Standard 94/572.
3. Calibration frequency: 28 days
It is recommended to calibrate when new working reagent is made.
After reagent lot changed.
As required following quality control procedures.

Quality control

Two level control materials should be analyzed with each batch of samples. In addition, the controls should be run with each new calibration, with each new reagent cartridge, and after specific maintenance or troubleshooting procedures as detailed in the appropriate system manual. We recommend using the supporting controls to verify the performance of the measurement procedure; other suitable control material can be used in addition.

Each laboratory should establish its own internal quality control scheme and procedures for corrective action if control can not recover within the acceptable tolerances.

Calculation

The analyzer calculates the FER concentration of each sample automatically after calibration.

Conversion factor: $\mu\text{g} / \text{mL} \times 10^{-3} = \text{ng} / \text{mL}$

Reference Intervals

Each laboratory should establish its own reference intervals based upon its patient population. The reference intervals measured at 37°C listed below

were taken from literature:

Sample Type	S.I. Units
Serum/Plasma	Male: 30-400 ng/mL Female: 15-150 ng/mL

Performance Characteristics

Representative performance data obtained from Mindray system (Mindray BS series analyzers / Mindray FER Reagent) is given below. Results may vary if a different instrument, an individual laboratory or a manual procedure is used.

Limitations-interference

The following substances were tested for interference with this methodology. Criterion: Recovery within ± 10 % of initial value.

Substance	Level Tested	Observed Effect
Intralipid	150 mg/dL	NSI*
Bilirubin	40 mg/dL	NSI*
Hemoglobin	500 mg/dL	NSI

* NSI: No Significant Interference (within $\pm 10\%$)

Linearity Range

The Mindray System provides the following analytical ranges:

Sample Type	S.I. Units
Serum/Plasma	10~1000 ng/mL

If the value of sample exceeds 1000 ng/mL, the sample should be diluted with 9 g/L NaCl solution (e.g. 1+9) and the result should be multiplied by 10.

Analytic Sensitivity/Limit of Detection

The lowest measurable FER concentration that can be distinguished from zero is 10 ng/mL with 99.7% confidence.

Precision

Precision performance using the CLSI Approved Guideline EP5-A2 to assay serum control appears in the table below. U: ng/mL

Type of Imprecision	Level 1			Level 2		
	Mean	SD	CV (%)	Mean	SD	CV (%)
Within-run		1.48	1.38		2.02	0.70
Between-run	107.31	0.88	0.82	289.19	1.76	0.61
Between-day		2.31	2.16		1.43	0.50
Within-device		2.89	2.69		3.04	1.05

Method Comparison

A comparison between Mindray System (Mindray BS series analyzers /Mindray FER Reagent) (y) and Hitachi System (Hitachi analyzer/FER Reagent) (x) using 40 samples gave following correlation (ng/mL): $y=0.9941x+0.7063$, $R^2=0.9883$

Details of the comparison experiments are available on request.

References

1. Williams WJ, Beutler E, Ersler AJ, Lichtman MA, eds. Hematology, 4th ed. New York: McGraw-Hill, 1990.
2. San Diego Declaration, Erythropoietin use and response in end-stage renal disease. The American Society of Nephrology, Annual meeting, San Diego. J Am Soc Nephrol 1995; 3:35.
3. Wick M, Pinggera W, Lehmann P. Iron Metabolism. Anemias. Clinical Aspects and Laboratory, 5th ed. Vienna/New York: Springer-Verlag, 2003.
4. Cook, J.D., Lipschitz, D.A., Laughton, M.B.B., Miles, E.M. & Finch, C.A.: Serum ferritin as a measure of iron stores in normal subjects. Am. J. Clin. Nutr. 27:680, 1974.
5. Marcus, D.M. & Zinberg, N.: Isolation of ferritin from human mammary and pancreatic carcinomas by means of antibody immunoadsorbents. Arch. Biochem. Biophys. 162:493, 1974.

Graphical symbols

In Vitro Diagnostic
medical device

Batch Code

European
Conformity

Authorized representative in
the European Community

Use By

Consult
Instructions for use

Biological
Risks

Temperature
Limit

Manufacturer

Catalogue
number

© 2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China

Website: www.mindray.com

E-mail Address: service@mindray.com

Tel: +86-755-81888998

Fax: +86-755-26582680

EC-Representative: Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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1. Williams WJ, Beutler E, Ersler AJ, Lichtman MA, eds. Hematology, 4th ed.

New york: McGraw-Hill,1990.

- 2.San Diego Declaration, Erythropoietin use and response in end-stage renal disease. The American Society of Nephrology, Annual meeting, San Diego. J Am Soc Nephrol 1995; 3:35.
- 3.Wick M, Pinggera W, Lehmann P. Iron Metabolism. Anemias.Clinical Aspects and Laboratory, 5th ed. Vienna/New York: Springer-Verlag, 2003.
- 4.Cook, J.D., Lipschitz, D.A., Laughton, M.B.B., Miles, E.M. & Finch, C.A.: Serum ferritin as a measure of iron stores in normal subjects. Am. J. Clin. Nutr. 27:680, 1974.
- 5.Marcus, D.M. & Zinberg, N.: Isolation of ferritin from human mammary and pancreatic carcinomas by means of antibody immonoadsorbents. Arch. Biochem. Biophys. 162:493, 1974.

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