# **GI-MA**

**For use on IMMULITE® 2000 systems**

**SIEMENS**

# IMMULITE® 2000 GI-MA

## English

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of CA19-9 in serum, as an aid in patient management.

Catalog Number: **L2KGI2** (200 tests)  
Test Code: **GIM** Color: **Dark Pink**

CA19-9 antigen levels in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. Results reported by the laboratory to the physician must include the identity of the assay used to measure CA19-9 antigen levels. Values obtained with different assays cannot be used interchangeably. Before changing assays, the laboratory must confirm baseline values for patients being serially monitored.

## Summary and Explanation

The GI-MA monoclonal antibody reacts specifically with the sialylated Lewis<sup>a</sup> blood group antigen (CA19-9). The antigen is expressed as a monosialoganglioside on secreted mucus glycoproteins in colorectal and pancreatic tumors. It can also be detected on mucin in patients' sera.<sup>1</sup> Sources include normal pancreas; bile duct; and gastric, colic, endometric and salivary epithelia. In healthy individuals, the CA19-9 antigen circulates at low levels, normally about 37 U/mL.<sup>16</sup> Elevated levels have been seen in benign inflammatory diseases of the hepatobiliary tract.<sup>2</sup> Studies have also shown that mucin bearing this antigen was more frequently detected in the sera of patients with pancreatic cancer than any other gastrointestinal carcinoma, including colorectal.<sup>3</sup> It is also found in bile duct, ovarian mucinous cystadenocarcinomas and uterine adenocarcinomas.

Pancreatic cancer in the nonjaundiced patient can be difficult to diagnose, even with the aid of modern scanning

techniques. The detection and quantification of mucin can be used in the serological diagnosis of the disease.

The use of assays for detecting the CA19-9 antigen in pancreatic cancer is well documented, with 68–94% sensitivity and 76–100% specificity being reported using approximately 37 U/mL as a cutoff.<sup>3-13</sup> It has been found in 18% of colon cancer patients and 67% of hepatobiliary cancer patients.<sup>3</sup> The CA19-9 antigen has also been found in the sera of cystic fibrosis patients,<sup>15</sup> and has been used in the serological diagnosis of the disease. Although useful as an adjunct in cancer diagnosis, it should not be considered a screening or diagnostic test when used alone.

## Principle of the Procedure

IMMULITE 2000 GI-MA is a solid-phase, two-site sequential chemiluminescent immunometric assay.

**Incubation Cycles:** 2 × 30 minutes

## Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 GI-MA has not been tested with all possible variations of tube types.

**Volume Required:** 50 µL serum

**Storage:** 1 day at 2–8°C, or 2 months at –20°C.<sup>16</sup>

## Warnings and Precautions

For *in vitro* use only.

**CAUTION:** This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

---

<b>H412</b> <b>P273, P501</b>	Harmful to aquatic life with long lasting effects. Avoid release to the environment. Dispose of contents and container in accordance with all local, regional, and national regulations. <b>Contains:</b> sodium azide; GI-MA Adjustors
----------------------------------	--

---

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

### GI-MA Bead Pack (L2GI12)

With barcode. 200 beads, coated with murine monoclonal anti-CA19-9 antibody. Stable at 2–8°C until expiration date.

**L2KGI2:** 1 pack

### GI-MA Reagent Wedge (L2GIA2)

With barcode. 11.5 mL protein-based buffer. 11.5 mL alkaline phosphatase

(bovine calf intestine) conjugated to murine monoclonal anti-CA19-9 antibody in buffer, with preservative. Stable at 2–8°C until expiration date.

**L2KGI2:** 1 wedge

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

### GI-MA Adjustors (LGIL, LGIH)

Two vials (Low and High) of CA19-9 in a lyophilized nonhuman serum matrix, with preservative. Reconstitute each vial with **4.0 mL** distilled or deionized water. Let stand for 30 minutes, then mix by *gentle* swirling or inversion. Stable at 2–8°C for 30 days, after reconstitution, or for 6 months (aliquotted) at –20°C.

**L2KGI2:** 1 set

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

## Kit Components Supplied Separately

### Multi-Diluent 2 (L2M2Z, L2M2Z4)

For the on-board dilution of patient samples. One vial, concentrated (ready-to-use), nonhuman protein/buffer matrix, with preservative. Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C.

**L2M2Z:** 25 mL    **L2M2Z4:** 55 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

**L2M2Z:** 3 labels    **L2M2Z4:** 5 labels

**L2SUBM:** Chemiluminescent Substrate

**L2PWSM:** Probe Wash

**L2KPM:** Probe Cleaning Kit

**LRXT:** Reaction Tubes (disposable)

**L2ZT:** 250 Sample Diluent Test Tubes (16 × 100 mm)

**L2ZC:** 250 Sample Diluent Tube Caps

Also required

Distilled or deionized water; test tubes; controls

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

**Recommended Adjustment Interval:**  
2 weeks

**Quality Control Samples:** Follow government regulations or accreditation requirements for quality control frequency.

Use controls or serum pools with at least two levels (low and high) of CA19-9.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined by an appropriate internal laboratory quality control scheme.

## Expected Values

Based on its relationship to IMMULITE GI-MA (see Method Comparison), the assay can be expected to have essentially the same reference ranges.

In a preliminary study, serum samples from 70 adults, both males and females, in apparent good health were processed by the IMMULITE GI-MA assay. The results ranged from nondetectable to 33 U/mL (absolute upper value reported), with a median of 5.1 U/mL.

Additionally, a reference range study for IMMULITE GI-MA was performed using serum samples from adult volunteers, including both men and nonpregnant women ranging from approximately 20 to 70 years of age (central 95%: 22–64 years, median: 40 years). The subjects were in apparent good health, based on a questionnaire.

Blood samples were collected in France, Germany, the Netherlands and Portugal. Results were generated by an independent laboratory in the Netherlands

using IMMULITE kits. (The samples were collected in plain glass tubes without anticoagulants, gel barriers, or clot-promoting additives, and assayed in singlicate.)

Medians and 95th percentiles for relevant subgroups are tabulated below. Percentiles were determined nonparametrically.

Group	Median	95%ile	Units	<i>n</i>
Males	3.0	18.3	U/mL	428
Females	4.2	18.7	U/mL	370
Combined	3.5	18.4	U/mL	798

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

## Limitations

CA19-9 is not expressed in about 3–7% of the population, both in healthy and ill subjects, who belong to the rare blood group of Lewis.<sup>a-b-</sup>

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

## Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in U/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

**Calibration Range:** Up to 1000 U/mL

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

**Analytical Sensitivity:** 1.0 U/mL

**High-Dose Hook Effect:** None up to 50,000 U/mL

**Precision:** Samples were processed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Samples spiked 1 to 19 with three CA19-9 solutions (400, 1000 and 2000 U/mL) were assayed. (See "Recovery" table for representative data.)

**Specificity:** The antibody is highly specific for CA19-9. (See "Specificity" table.)

**Anticoagulants:** No significant effect.

**Bilirubin:** Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

**Biotin:** Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

**Hemolysis:** Presence of hemoglobin in concentrations up to 512 mg/dL has no effect on results, within the precision of the assay.

**Lipemia:** Presence of triglycerides in concentrations up to 3000 mg/dL has no effect on results, within the precision of the assay.

**Method Comparison:** The assay was compared to IMMULITE GI-MA on 92 samples. (Concentration range: nondetectable to approximately 1000 U/mL. See graph.) By linear regression:

$(IML\ 2000) = 1.09 (IML) - 0.43\ U/mL$   
 $r = 0.998$

Means:

102 U/mL (IMMULITE 2000)  
93.2 U/mL (IMMULITE)

## References

1) Magnani JL, Stepkowski Z, Koprowski H,

Ginsburg V. Identification of the gastrointestinal and pancreatic cancer-associated antigen detected by monoclonal antibody 19-9 in the sera of patients as a mucin. *Cancer Res* 1983;43:5489-92. 2) Frebourg T, Bercoff E, Manchon N, et al. The evaluation of CA19-9 antigen level in the early detection of pancreatic cancer. *Cancer* 1988;62:2287. 3) Herlyn M, Sears HF, Stepkowski Z, Koprowski H. Monoclonal antibody detection of a circulating tumor-associated antigen. 1. Presence of antigen in the sera of patients with colorectal, gastric and pancreatic carcinoma. *J Clin Immunol* 1982;2:135-40. 4) Cooper EH, Forbes MA, Taylor M. An evaluation of DUPAN-2 in pancreatic cancer and gastrointestinal disease. *Br J Cancer* 1990;62:1004-5. 5) Rhodes JM, Ching C-K. Serum diagnostic tests for pancreatic cancer. In: Neoptolemos J, editor. *Bailliere's clinical gastroenterology*. London: Bailliere Tindall, 1990: 833-52. 6) Steinberg W. The clinical utility of the CA19-9 tumor-associated antigen. *Am J Gastroenterol* 1990;85:350-5. 7) Molina LM, Diez M, Cava MT, Maestro ML, Ortega MD, Mendiz JG, et al. Tumor markers in pancreatic cancer: a comparative clinical study between CEA, CA19-9 and CA50. *Int J Biol Markers* 1990;5:127-32. 8) Ferrara C, Basso D, Fabris C, Malesci A, Fogar P, Megiato T, et al. Comparison of two newly identified tumor markers (CAR-3 and DUPAN-2) with CA19-9 in patients with pancreatic cancer. *Tumori* 1991;77:56-60. 9) Kawa S, Kato M, Oguchi H, Kobayashi T, Furuta S, Kanai M. Elevated serum levels of Dupan-2 in pancreatic cancer patients negative for Lewis blood group phenotype. *Br J Cancer* 1991;64:899-902. 10) Kobayashi T, Kawa S, Tokoo M, Oguchi H, Kiyosawa K, Furuta S, et al. Comparative study of CA-50 (time-resolved fluorimmunoassay), Span-1, and CA19-9 in the diagnosis of pancreatic cancer. *Scand J Gastroenterol* 1991;26:787-97. 11) Lucarotti ME, Habib NA, Kelly SB, Rothnie ND, Nelson O, Lindholm L, et al. Clinical evaluation of combined use of CEA, CA19-9 and CA50 in the serum of patients with pancreatic carcinoma. *Eur J Surg Oncol* 1991;17:51-3. 12) Masson P, Palsson B, Andren-Sandberg A. Evaluation of CEA, CA19-9, CA-50, CA-195 and TATI with special reference to pancreatic disorders. *Int J Pancreatol* 1991;8:333-44. 13) Motoo Y, Kawakami H, Watanabe H, Satomura Y, Ohta H, Okai T, et al. Serum sialyl-Tn antigen levels in patients with digestive cancers. *Oncology* 1991;48:321-6. 14) Niederau C, Grendell JH. Diagnosis of pancreatic carcinoma. Imaging techniques and tumor markers. *Pancreas* 1992;7:66-86. 15) Roberts DD, Rose MC, Wang WT, et al. Isolation and characterization of mucin from the serum of cystic fibrosis patients. *Am J Respir Cell Mol Biol* 1990;2:373. 16) Burtis CA, Ashwood ER, editors. *Tietz textbook of clinical chemistry*. 2nd ed. Philadelphia: W.B. Saunders, 1994: 920.

## Technical Assistance

Available outside the United States only.  
For technical assistance, contact your  
National Distributor.

[www.siemens.com/diagnostics](http://www.siemens.com/diagnostics)

The Quality System of Siemens Healthcare  
Diagnostics Products Ltd. is certified to ISO  
13485.

---

## Tables and Graphs

---

### Intraassay Precision (U/mL)

	Mean <sup>1</sup>	SD <sup>2</sup>	CV <sup>3</sup>
1	7	0.6	9.2%
2	30	2.5	8.2%
3	384	25.1	6.5%
4	616	46.5	7.6%

### Interassay Precision (U/mL)

	Mean <sup>1</sup>	SD <sup>2</sup>	CV <sup>3</sup>
1	7	0.7	10.0%
2	30	2.9	9.4%
3	384	26.3	6.8%
4	616	50.8	8.3%

### Linearity (U/mL)

	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	8 in 8 <sup>5</sup>	21	—	—
	4 in 8	10	11	91%
	2 in 8	4.7	5.3	89%
	1 in 8	1.8	2.6	69%
2	8 in 8	20	—	—
	4 in 8	9.1	10	91%
	2 in 8	4.3	5.0	86%
	1 in 8	1.8	2.5	72%
3	8 in 8	52	—	—
	4 in 8	27	26	104%
	2 in 8	12	13	92%
	1 in 8	5.6	6.5	86%
4	8 in 8	64	—	—
	4 in 8	29	32	91%
	2 in 8	14	16	88%
	1 in 8	6.1	8.0	76%
5	8 in 8	129	—	—
	4 in 8	59	65	91%
	2 in 8	28	32	88%
	1 in 8	13	16	81%
6	8 in 8	240	—	—
	4 in 8	110	120	92%
	2 in 8	51	60	85%
	1 in 8	28	30	93%

### Specificity


Compound <sup>1</sup>	Amount Added <sup>2</sup>	% Cross reactivity <sup>3</sup>
CA125	1000 U/mL	ND <sup>4</sup>
CA15-3	720 U/mL	ND
CEA	1000 ng/mL	ND

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2015 Siemens Healthcare Diagnostics. All rights reserved.