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SNOL 58/350 LSP11

Economical low temperature oven SNOL 58/350 is designed for thermal treatment such as drying, heating, thermal testing, aging and similar purposes in forced airflow environment up to 350°C. Forced air convection allows a homogeneous temperature distribution during all process, which ensures optimal results; good technical parameters ensure high quality results. This oven can be used in scientific laboratories, educational institutions, medicine and industry.

Basic model

- ✓ Control panel is placed on the left side;
- ✓ Ventilation motor is placed on the back side;
- ✓ Forced horizontal air circulation;
- ✓ Mild steel chamber;
- ✓ 3 shelves;
- ✓ Outside casing – metal sheet, powder painted grey;
- ✓ Insulation – rock wool (complete lack of asbestos);
- ✓ Door opening to the right, inner door – stainless steel, silicon joint;
- ✓ Microprocessor – temperature controller;
- ✓ Low electric power usage for increased energy efficiency;
- ✓ Short heating up / cooling down period.

Options

- ✓ Additional shelves;
- ✓ Reinforced shelves;
- ✓ Metal tray;
- ✓ Digital timer;
- ✓ Buzzer;
- ✓ Over-temperature protection;
- ✓ Data recorder;
- ✓ Data communication/USB;
- ✓ Calibration of temperature measurement system;
- ✓ Oven exterior made from stainless steel;
- ✓ Table for supporting the oven;
- ✓ Additional 1 year warranty.

Specifications

Technical data	SNOL 58/350 LSP11
Volume, L	58
Maximum temperature, °C	350
Continuous operating temperature, °C	T+10-350
Power	2
Rated supply voltage, V	230
Number of phases	1
Rated frequency, Hz	50
Chamber material	Mild steel
Maximum heating-up time (without charge), min.	40
Temperature stability	1
Temperature uniformity	3°C
Airflow	horizontal
Shelves (number)/max.	3/7
Chamber width, mm	390
Chamber depth, mm	375
Chamber height, mm	360
Overall width, mm	670
Overall depth, mm	615
Overall height, mm	580
Mass (Netto), kg	40

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*LOW TEMPERATURE LABORATORY RESISTIVE
ELECTRIC FURNACE*

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TABLE OF CONTENTS

1. PURPOSE	26
2. SYMBOLS AND MARKINGS	26
3. PACKAGE CONTAINS	27
4. SAFETY REGULATIONS	27
5. PARTS AND TECHNICAL SPECIFICATIONS	28
5.1. PARTS OF THE ELECTRIC FURNACE	28
5.2. TECHNICAL SPECIFICATIONS	30
6. REQUIREMENTS FOR SAFE USAGE	30
7. UNPACKING, CHECKING EQUIPMENT, DELIVERY	31
8. ASSEMBLY	31
9. PREPARATION FOR USE	32
10. OPERATION	33
11. STORAGE	33
12. TRANSPORTATION	33
13. WARRANTY	34
14. TROUBLESHOOTING	34
15. MAINTENANCE	35
16. ELECTRICAL PRINCIPLE SCHEME	36
17. THE PACKAGING PARTS REMOVAL	42
18. FURNACE CONFORMITY CERTIFICATION	43

1. PURPOSE



Low temperature laboratory electric Furnace **SNOL XX/350** (further called as "Furnace") is designed for heating or drying materials in the ambient air temperatures from T+10 °C to 350 °C in static conditions (T – temperature of the room where Furnace is being operated).

2. SYMBOLS AND MARKINGS

Markings

	alternating current
	earth/grounding contact
	danger of electrical shock
	caution, hot surface

EXAMPLE OF PRODUCT LABEL:

Identification on the label	Meaning
Model: SNOL 58/350	Type of Furnace, 58 – total chamber volume (58 dm ³); 350 – nominal temperature (350°C).
IP 20	Degrees of protection provided by enclosures EN60529
No 0000	Product serial number
230V	Voltage of power supply, V
~50/60 Hz	Frequency of electrical current, Hz
2,0 kW	Nominal power, kW
40 kg	Weight of the product, kg
2022 01	Date of manufacture (year, month)
	Mark of conformity with EU requirements
	Mark of conformity with Russia requirements

Warning symbols

	This symbol precedes any important safety instructions. Personal and property safety could be put at risk if these instructions are not followed. Before switching on the low temperature Furnace, it is very important to read and follow all the instructions of this manual.
	This symbol precedes manual requirements which must be followed with great care in order to ensure the Furnace's correct operation, to prevent damage and to protect the safety operators.
	This symbol precedes important information.

3. PACKAGE CONTAINS

The user is supplied with:

Name	Quantity pcs.
Electrical Furnace SNOL XX/350 , pc.	1
Feet, pc.	4
Instruction manual, pc	1
Temperature controller user manual, pc.	1
Shelves, pc.	3

4. SAFETY REGULATIONS

The Furnace must be used in accordance with the requirements of the user's technical safety operation rules for electric equipment and local legislation.

The operator working with the electric Furnace must be informed about the operational guidelines for electric equipment up to 1000V, acquainted with safety regulations for the operation of an electric Furnace, and knowledgeable about the product's operating principles and proper use.


The Furnace must be earthed and should only be connected to a mains electrical supply that has an earth/grounding contact. The Furnace must be connected to a mains electrical supply of the correct rating, check the label on the rear of the Furnace for voltage, frequency and power usage.

	It is prohibited to operate the electric Furnace if the earth/grounding is insecure.
	It is prohibited to operate the electric Furnace if any of external protective shields are removed.
	If the electric Furnace is operated for a longer period of time, the external surfaces and the door may get hot. Wear heat-resistant gloves when opening the door of the hot electric Furnace, loading or unloading it, or touching the hot external surfaces.
	Combustible or explosive materials or materials that can become combustible or cause explosion during thermal treatment must not be loaded into the low temperature Furnace.
	Do not load unknown materials into the Furnace.
	The gas emitted during thermal treatment has to be directed outside through a ventiduct (required by existing regulations).
	The room where the Furnace is operated must have appropriate ventilation.
	It is forbidden to operate the Furnace in a closed, unventilated room.
	A room in which several Furnaces are in operation must have a dedicated means of ventilation.
	In the event of obvious malfunctions, the Furnace must be immediately turned off. Maintenance can be performed by SnolTherm service operators, an authorized specialist or a user who strictly follows SnolTherm written maintenance instructions.
	In the event the Furnace ceases normal operations, disconnect it from the power supply (the mains) and take measures to resolve the problem. The manufacturer is not responsible for damages arising from the use of the Furnace for purposes for which it was not intended; the user bears all responsibility. Any operations that put at risk the Furnace's safety are forbidden. General safety and accident prevention rules should be applied. Only Furnaces that are technically sound should be used.
	Failure to comply with instructions could result in negative consequences for the user.

The manufacturer 's advice cannot cover all issues. Therefore, the user assumes responsibility for the reaction and the impact that heat has on materials and the risks related to conducting tests. Our advice will help to reduce the dangers faced by users as well as protect the Furnace and the tested materials from damages.

5. PARTS AND TECHNICAL SPECIFICATIONS

5.1. Parts of the electric Furnace

The Furnace is composed of housing, chamber, door, control panel.

The Furnace's housing made of sheet steel. The working chamber (stainless or mild steel) is installed on the housing's right side. Shelves can be placed in the chamber. Evenly distribute charge on the shelves. A tape heating element is used for heating.

Thermal insulation material is located between the housing and the working chamber and in the door.

The control panel is located on left side.


*There are holes for ventilation in the housing. In Furnace model 58/350, a circulating ventilator is located in the chamber.

The Furnace has a 230 V power supply, 50/60 Hz frequency. A digital temperature controller and a thermocouple temperature sensor mounted in the chamber work together to control and adjust the temperature. The controller shows the set and present temperature.

Electric scheme is attached (see page 36).

ADDITIONAL ACCESSORIES MAY BE INSTALLED IN THE FURNACE:

1. Process observation window;
2. Protection from overheating, which prevents the Furnace from exceeding the maximum temperature. If the maximum temperature is exceeded, the safety device turns the heating off automatically. If this happens, allow the Furnace to cool down by 50-100°C below overheating safety device activation temperature. Turn the Furnace on. To reset safety device to the inactive mode, find the red button at the front of the Furnace (Pic.2) and press it (In each case, the mounting position of the button may vary). If the safety device is going to run again contact the manufacturer or to the company which sold the product (company addresses at the end of the passport);

Picture 2

3. Timer is designed to delay execution of the program. If timer GLX72Q41, is mounted in the Furnace then, when the main switch is on, turn on the fixed control of the timer. For that reason, press the first, counting from left to right, button of the timer. The symbol of closed contacts  1 on the display will appear and the Furnace will start working normally.
If you want to use the timer to turn the Furnace on / off at the certain time, by following the User Manual of the timer, set up the Furnace on/off program. After putting in the charge turn the Furnace on, so that it starts to heat up, and set the timer from fixed to programming control. For this reason, press first button on the timer. The note **“FIX”** should fade from timer’s display;
4. Buzzer is to alarm about the end of the program;
5. RS-485 / USB interface and software for data recording, viewing and configuring the temperature controller running your thermal treatment process. The software is designed for Windows operating system. Computer software allows simply run, review and display charts on thermal process temperatures and other settings;
6. Data recorder Eurotherm 6100 is for basic visualization and recording requirements. The 6100 has full color display and utilises touch screen technology for clear and intuitive configuration and operation. It further supports a USB port as standard to enable the use of a mouse, keyboard or bar code scanner. Data can be moved manually or automatically archived to multiple locations: removable media, network servers or the Eurotherm Review database on a PC. The recorder can easily be integrated into a larger system and data files can be transferred across the network;

5.2. Technical specifications

		Measure	SNOL 58/350	SNOL 67/350
Nominal power		kW	2,0±9%	2,0±9%
Voltage		V	230	
Frequency		Hz	50/60	
Ingress Protection class (EN60529)		IP20	IP20	
Nominal chamber temperature		°C	350	
Number of phases			1	
Atmosphere in the chamber			Air	
Temperature fluctuations**		°C	±1	±2
Limits of automatic temperature control		°C	T+10 ÷ 350	T+10 ÷ 350
Temperature distribution ***		°C	±2,5	±14
Heating time ****		min	40	40
Dimensions of the chamber:	width	mm	390	390
	depth	mm	380	445
	height	mm	360	390
Number of shelves (standard/max)		pcs.	3/7	3/7
Non-concentrated load of one shelf		kg	15	15
Total limiting load		kg	30	30
Overall dimensions of the Furnace	width	mm	675	675
	depth	mm	630	625
	height	mm	585	600
Weight		kg	40	37

T – temperature of the room where Furnace is being operated

**Temperature fluctuations at steady mode without load. The value of the parameter is given after the temperature controller has been auto-tuned with a typical (see temperature controller instruction manual).

***Temperature distribution in the working space without load and shelves when the low temperature laboratory electric Furnace is heated to the nominal temperature

****Heating time for reaching the nominal temperature without load.

6. REQUIREMENTS FOR SAFE USAGE

The electric Furnace should be placed in a closed, ventilated room that meets the following requirements:

	- the electric Furnace should be placed on a base that is horizontal (an allowed unevenness for a 1m base is ± 1mm.), hard, made of incombustible material; - height above sea level – up to 2000 m; - ambient temperature from +5 to +40°C; - ambient relative air humidity should not exceed 80% at the temperature of +31°C, and should not exceed 50% at the temperature of +40°C; - the voltage fluctuation in the mains (power supply) should not exceed ± 10% of the nominal value; - the environment should not pose an explosion risk; it also should not contain a large amount of electrically conductive dust, water vapor, aggressive gas; - it is recommended that the Furnace be placed under an exhaust ventilation hood; - the electric Furnace should not be subjected to vibration and shocks; - it is prohibited to exceed the maximum rated temperature, otherwise the service life of heating elements will be shortened and the thermocouple could be damaged.
---	---


7. UNPACKING, CHECKING EQUIPMENT, DELIVERY

Please check the product in its original packaging, if there are any damages of transportation. Inform the carrier immediately if transportation damage has occurred.

Carefully flip the packaging with the Furnace in it. Unscrew wood screws that hold pallet to the cover. Take off the cover.

Remove transportation screws from the bottom of the pallet (Pic. 2).

Take the protective materials out of the chamber.


Picture 2

	 Attention ! Sliding or tilting of the unit. Damage to the unit. Risk of injury by lifting heavy loads. Do not lift or transport the unit using the door handle, the door or the lower housing. Forbidden to raise or to carry one person, if the stove weighs more than 30 kg. Products, which weight from 30 to 60 kg, have lift and transport two people. Products, which weight in excess of 60 kg, lifting and transporting to be used auxiliaries (lifts, etc.).
--	--

If you need to return, Use the original packaging for safe transportation. Dispose of the transport packaging (See page 42).

8. ASSEMBLY

Unpack the electric Furnace and place it in the location prepared for its operation. Gently flip the Furnace and at the bottom part (Pic.2) screw in the feet (Pic.3). Stabilize Furnace by twisting feet. Put in the shelves to the chamber.

Picture 2

Picture 3

Voltage specified in technical specification chart (see page 30) must correspond to the rated voltage of the supply network.


Before connecting to the mains, check the serviceability of the cord, plugs and sockets with a grounding contact.

Plug the cord to the mains.

Do not position the furnace in a manner that would make it difficult to operate your main power disconnect switch.

Prepare the place of use for the Furnace.


-The supporting surface, where the Furnace will be placed on, must be smooth and non-combustible.

-To prevent damage to the floor covering or the risk of fire if a hot object falls during removal from the chamber, the floor covering in the working area of the Furnace must be non-combustible.

-The Furnace must be installed at a safe distance from the surrounding surfaces:

- The minimum distance from the rear and side surfaces of the Furnace to flammable surrounding surfaces is $L = 0.5 \text{ m}$;
- The minimum distance from the top surface of the Furnace to flammable surrounding surfaces is $L1 = 1 \text{ m}$;
- The minimum distance from the rear and side surfaces of the Furnace to non-combustible surrounding surfaces can be reduced to $L = 0.2 \text{ m}$;
- The minimum distance from the rear surface of the Furnace with a ventilation hole or chimney in the rear wall to the surrounding surfaces $L = 1 \text{ m}$;

In individual cases, safety distances must be adapted to the environment and larger distances must be selected.

The room where the Furnace is installed must have a good ventilation system. Working with the processes that emit vapors, fumes, etc. in addition, we recommend installing hoods or connecting the furnace chimneys to the ventilation system as shown in the figure.

9. PREPARATION FOR USE

The Furnace should be dried before using it for the first time or if it has been out of operation for a long period of time and stored in humid conditions. It should be dried as follows:

- without load, heat the Furnace until it reaches a temperature ranging from $100-150^{\circ}\text{C}$ and maintain this temperature range for 2-3 hours;
- heat up to the nominal temperature. Maintain for 1-2 hours.

The dried Furnace can now be operated.

While drying, some 'smoke' may appear but it has no effect on further operation of the electric Furnace.

10. OPERATION

Open the door of the electric Furnace.


Put the load on the shelves at least 1/10 of the space away from the side walls to ensure a better circulation of air.

When putting or taking off the charge do not damage thermocouple, mounted in the chamber.

The data indicated in technical specifications gives to the working space.

A, B, C – dimensions of the chamber;
 a, b, c – spaces between the walls and the load;
 $a = A \times 0,1$
 $b = B \times 0,1$
 $c = C \times 0,1$
 Working volume (space)
 $V = (A-2a) \times (B-2b) \times (C-2c)$

Close the electric Furnace's door.

Turn on the switch; a pilot lamp will light up, indicating that the power is on.

Not leaving the Furnace at high temperatures for longer than necessary will prolong its life.

After work switch off the Furnace. The switching-off will be indicated by the pilot lamp that fades. When program ends the bell sounds.

11. STORAGE

Store the electric Furnace in its packing in a location that meets the following requirements:

- heated and naturally ventilated room with an ambient temperature ranging from +5 to +35°C;
- atmosphere with sulphurous gas concentration not exceeding 0,13 mg/m³ and chlorine salts concentration not exceeding 0,3 mg/m² per 24 hours;
- the ambient relative air humidity should not exceed 80% at the temperature of 25°C.

During loading and unloading of the Furnace, protect it from any shocks.

12. TRANSPORTATION

The electric Furnace should be transported only in closed means of transport in which the temperature ranges from - 50 to + 50°C; at the temperature of +25°C, the ambient relative air humidity must not exceed 80%.

You can transport the electric Furnace by any means of closed transport, but protect it from any rough, jarring movements. The packaging does not protect the electric Furnace from the effects of improper handling.

13. WARRANTY

The manufacturer guarantees that the electric Furnace meets company standards.

The Furnace is guaranteed for 12 months from the date of purchase provided that the user follows the instructions regarding the Furnace's storage, transportation and usage. However, the guarantee does not extend more than 24 months beyond the Furnace's manufacture date.

Manufacture defects that appear during the warranty period shall be repaired at the manufacturer's expense. The manufacturer's address:

Umega Group, AB, SnolTherm ,
Narkûnai LT-28104 Utena, Lithuania;
Tel +370 389 54586; Fax +370 389 61223;
E-mail: info@snoltherm.com ;
www.snol.com.

14. TROUBLESHOOTING

Problem	Possible Cause	Solution
The Furnace doesn't turn on.	There is no rated power supply voltage; The fuse is burnt out.	To check the power supply; Change the fuse.
The Furnace doesn't heat up.	The temperature controller is not switched on; Damaged heating element; Damaged solid-state relay; Tripped a protective device.	Switch on temperature controller; Change heating element; Change solid-state relay; The back wall to activate the device. *
Heating time is longer than expected.	Low voltage of the power supply; Door does not maintain a hermetic seal.	Check the power supply's voltage; Tighten the door's closure.
Temperature fluctuation or large over-shoot of the set temperature.	Temperature Controller requires control terms tuning.	Perform an Auto-tune (AT) on the temperature controller (look at temperature controller instructions).

	*If after switching on the protective device it triggers again, the Furnace must not be operated. In that case the manufacturer or the manufacturer 's representative should be addressed for the repairs of the Furnace.
---	---

	In order to continually improve our products and their usage, supplements to this publication may be introduced without notice.
---	---

ASSESSMENT TECHNIQUE OF THE LABORATORY ELECTRIC FURNACE'S PARAMETERS

TESTED PARAMETERS:

1. Verification of nominal temperature and temperature stability without charge,
2. Measurement of temperature distribution

TEST METHODS

All tests are performed when the electric Furnace is heated to the nominal temperature set point and operates in the established temperature mode.

Before measuring the temperature stability, it is necessary to select the PID parameters of the temperature controller (see User Manual of temperature controller).

1. Verification of nominal temperature and temperature stability is performed without charge.

In the cold Furnace, the control thermocouple is fastened as close to the working thermocouple as possible. The nominal temperature is set at the temperature controller. When the Furnace is heated and the thermodynamic equilibrium is settled, the control points of the thermocouple are recorded in the interval of not less than 1 hour, every 60 seconds (but not at longer intervals). The arithmetic mean of these values is the measured value of the nominal temperature.

The maximum (positive and negative) temperature deviations from the average temperature are temperature stability.

2. The temperature distribution is measured by thermocouples with the recorder. The thermocouples are positioned at the corners of the chamber, from surfaces at a distance of 0.1 edge lengths, and in the center. At least five thermocouples are used. The nominal thermal temperature is set. When the Furnace is heated and the thermodynamic equilibrium stabilizes, the temperature is captured. The result is calculated according to the formula:

$$T = \pm \frac{T_{\max} - T_{\min}}{2}$$

Where:

T - temperature distribution, °C;

$T_{\max}$ - is the maximum temperature obtained during measurement at one of the points, °C;

$T_{\min}$ - is the minimum temperature obtained during measurement at one of the points, °C;

15. MAINTENANCE

 Attention !	The electric Furnace should be disconnected from the mains (power supply) and cooled.
---	--


Change of fuse	Press the holder of the fuse, turn anticlockwise and take out from the socket. Change a burnt-out fuse. Screw in the holder of the fuse in the reverse order. Ensure the fuse holder is firmly tightened.
Once every six months	Once every six months visually examine the power lead and inside of the chamber for damage.
Cleaning	When an operation is completed, clean the external surfaces of the electric Furnace, except for the markings, with a water-dampened piece of cloth. Clean out the scale from the heating chamber with a water-dampened piece of cloth.
Cleaning of the fan openings	Clean the vents with a vacuum cleaner.

You are responsible for the correct boot sterilization.

Before cleaning or sterilization (in addition to what is recommended by the manufacturer), the user must make sure that his chosen method of sterilization and cleaning products will not damage the equipment.

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
If recycling is not possible, all packing parts can also be disposed of with normal waste.

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18. FURNACE CONFORMITY CERTIFICATION

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Low temperature laboratory Furnace model was tested and meets the company standards.

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Control department seal and signature

**ATITIKTIES DEKLARACIJA
ДЕКЛАРАЦИЯ СООТВЕТСТВИЯ
DECLARATION OF CONFORMITY
KONFORMITÄTSEKTLÄRUNG**

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Low – temperature laboratory electric Furnaces **SNOL 0,1÷5000/30÷600** are designed for laboratory analytical and thermal treatment of different materials at temperatures at 30°C to 600°C.

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We herewith declare under our sole responsibility that the following product is in conformity with the requirements of the European Union Low Voltage Directive 2014/35 / EU - Standards LST EN 61010-1: 2011, LST EN 61010-2-010: 2015, Safety of Machinery

for the European Union Requirements for Directive 2006/42 / EC - Standard LST EN ISO 12100: 2011, European Union Electromagnetic Compatibility Directive 2014/30 / EU - Standard LST EN 61000-3-2:2014, LST EN 61000-6-1:2007, LST EN 61000-6-1:2007/P:2008 and complies with the requirements of LST EN 55011: 2010 GOST 12.2.091-2012, GOST IEC 61010-2 -010-2013 - TP TC 004/2011, to standard GOST P 51522.1-2011 - CU TR 020/2011.

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Bureau Veritas Certification Holding SAS – UK Branch certifies that the Management System of the above organisation has been audited and found to be in accordance with the requirements of the management system standards detailed below

ISO 9001:2015

Scope of certification

Production and design of laboratory and industrial electric and gas stoves, prefabricated shelves, mobile racks, euro pallet stands, metal workspaces, metal cabinets; production of high temperature thermal insulation materials

Original Cycle Start Date: **2003-07-01**

Expiry date of previous cycle: **NA**

Certification / Recertification Audit date:

Certification / Recertification cycle start date: **2021-07-01**

Subject to the continued satisfactory operation of the organization's Management System, this certificate expires on: **2024-06-30**

Certificate No.: **LT006620**

Rev: **1**

Issue date: **2023-03-02**

0008

Certification Body Address: 5th Floor, 66 Prescott Street, London, E1 8HG, United Kingdom

Local Office: Ukmergės g. 369A, LT-12142 Vilnius, Lithuania

Further clarifications regarding the scope and validity of this certificate, and the applicability of the management system requirements, please call: +370 5 233 79 75

**ATITIKTIES DEKLARACIJA
ДЕКЛАРАЦИЯ СООТВЕТСТВИЯ
DECLARATION OF CONFORMITY
KONFORMITÄTSEKTLÄRUNG**

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Low – temperature laboratory electric furnaces **SNOL 0,1÷5000/30÷600** are designed for laboratory analytical and thermal treatment of different materials at temperatures at 30°C to 600°C.

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We herewith declare under our sole responsibility that the following product is in conformity with the requirements of the European Union Low Voltage Directive 2014/35 / EU - Standards LST EN 61010-1: 2011, LST EN 61010-2-010: 2015, Safety of Machinery for the European Union Requirements for Directive 2006/42 / EC - Standard LST EN ISO 12100: 2011, European Union Electromagnetic Compatibility Directive 2014/30 / EU - Standard LST EN 61000-3-2:2014, LST EN 61000-6-1:2007, LST EN 61000-6-1:2007/P:2008 and complies with the requirements of LST EN 55011: 2010 GOST 12.2.091-2012, GOST IEC 61010-2 - 010-2013 - TP TC 004/2011, to standard GOST P 51522.1-2011 - CU TR 020/2011.

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www.snol.com

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SnolTherm GmbH

Registered office: Asbacher Straße 27 a, 53577

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Permanent establishment: Winchesterstr. 2,

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www.snol.ua

Letter of Authorization

Snoltherm, UAB as Lithuanian producer of SNOL laboratory and industrial furnaces and ovens hereby confirms that the company:

KIRANTONI“ SRL
Miron Costin street, nr. 18-132, Chisinau
MD-2068, Republic of Moldova

Has the right to sell SNOL production furnaces, ovens and spare parts as well as to provide after sales service for these machines in the territory of **Moldova**. This authorization letter is valid until 2024.12.31. Any other letter of authorization or any other powers of attorney issued previously to the abovementioned person are null and void from [2024.01.08].

Andrius Ignatavičius
Sales Director

Snoltherm, UAB
Plento str. 3, Narkūnai, LT-28104, Utena dis.,
Lithuania