**Biosigma SpA**

**a Dominique Dutscher Company**

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) |

Fax ++39 0426 302228 | SMS ++39 348 4077376 |

E-mail info@biosigmaeu.com | http://www.biosigma.com |

**M521D**

26/05/2022

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
### FOR THE STANDARDIZATION OF MICROSCOPIC URINALYSIS


FAST - READ 102® system improves the standardization of microscopic urinalysis and provides precision and reproducibility such as to guarantee constant readings which will not be influenced by variations of techniques among different operators. FAST - READ 102® is made of a slide, protected by an optically transparent film, with 10 independent chambers containing a standard volume of 7 µL. After dispensing the sample on the slide application area by means of a capillary mechanism, the sediment is homogeneously distributed in the reading chamber. Each chamber is fitted with its own system for the collection of excess urine to prevent any possible contamination. Furthermore FAST - READ 102® S slide allows easy identification of the sample by using the numbers printed on the sides.

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
*After completing the preparation of the urinary sediment by normal centrifugation,*

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*Use the pipette to mix the urinary sediment by repeatedly pressing the bulb; then*

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*HOW TO USE THE NUMBERING SYSTEM  
Example: For the identification of samples from number 131 to 140 mark number 1 on the left*

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

insert a pipette into the test tube and dispense one drop of STAIN.

dispense a drop of sample onto the appropriate area of the slide.

side and number 3 on the right side of the slide.

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### GB - FAST READ 102° CELL COUNT METHOD

Fast Read 102° is a plastic device composed of 10 counting chambers. Each chamber contains a GRID with 10 SQUARES, subdivided into 16 small squares (called SECTORS) .

After filling the counting chamber with the sample, proceed to the counting of cells distributed in N squares. Considering that the grid contains 10 squares and each square has a dimension of 1 x 1 mm, a depth of 0.1 mm and a volume of 0.1 ul, the formula for determining the concentration of cells (cells / ml) is:

$$[\text{Cells/ml}] = \frac{(\sum \text{cells counted in N squares})}{N} \times \text{dilution factor} \times 10^4$$

$10^4$  = conversion from 0.1ul to 1 ml

Attention to the cells at the edges, you should only count those on two sides to avoid the risk of overestimates or underestimates.

In the example shown above we perform the cell counting on a sample concentrated 10 times:

N = 5 (number of squares considered for counting)

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*The device is manufactured under license Hycor Biomedical Inc. (EU Patent Number 326349)*

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**MATERIALS**: The materials utilized to produce the in vitro labware devices are selected on the base of the properties indicated for the individual destination use. The device BVS100 is manufactured in METHACRYLATE (PMMA).

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### **IN VITRO DIAGNOSTIC MEDICAL DEVICE CE**

**For professional use only**

#### **RECOMMENDATIONS FOR USE**

(as for information)

- *Use only for laboratory analysis.*
- *Do not hold the device to the flame or heat sources that may cause damage.*
- *Do not use after expiry or if box is already open.*
- *Do not vary the intended purpose of the product.*
- *Keep out of reach of children.*
- *Store in a dry place, away from sources of heat and direct sunlight.*
- *Waste disposal: utilize the provided D.P.I and discharge under the laws in force.*