A large, white, circular MRI machine with a central opening for the patient. The machine has a complex, multi-segmented design. The text "GE HealthCare" is visible on the left side of the gantry, and "SIGNA Changegun" is on the right. There are also some smaller labels like "ECO" and "4.7T" near the top and bottom of the central opening. The machine is shown from a slightly elevated, side-on perspective.


**ALL INFORMATION IN  
THIS DRAWING SET IS  
PRELIMINARY, SUBJECT  
TO CHANGE!**

[illegible]

CONFIDENTIAL  
SECRET DE AFACERI

# EQUIPMENT LAYOUT

| ITEM                                 | DESCRIPTION                              | DIMENSIONS<br>LxWxH (mm) | WEIGHT<br>(kg) |
|--------------------------------------|------------------------------------------|--------------------------|----------------|
| 1                                    | MAGNET                                   | 1938x2450x2400           | 4553           |
| 2                                    | FIXED PATIENT TABLE (PT)                 | 2158x572x965             | 136            |
| 3                                    | MAGNET RUNDOWN UNIT (MRU)                | 287x206x172              | 3.2            |
| 4                                    | FURNITURE FOR COILS (NOT SUPPLIED BY GE) | -                        | -              |
| 5                                    | REAR PEDESTAL                            | 853x510x956              | 66             |
| 6                                    | INTEGRATED COOLING CABINET (ICC)         | 902x801x2048             | 545            |
| 7                                    | INTEGRATED SYSTEM CABINET (ISC)          | 1275x974.5x2000          | 1116           |
| 8                                    | MAGNET MONITOR (MON)                     | 381x260x127              | 4.5            |
| 9                                    | MAIN DISCONNECT PANEL (MDP)              | 216x508x914              | 43.5           |
| 10                                   | OPERATOR CONSOLE (GOC)                   | -                        | 59.3           |
| 11                                   | PENETRATION PANEL (PP)                   | 255x243x723              | -              |
| WALL - ACCORDING TO RECEIVED DRAWING |                                          |                          |                |
| WALL TO DEMOLISH                     |                                          |                          |                |
| RF SHIELD - 100 dB ATTENUATION       |                                          |                          |                |
| EXAM ROOM HEIGHT                     |                                          |                          |                |
| FINISHED FLOOR TO SLAB HEIGHT        |                                          |                          | 3.50 m         |
| FALSE CEILING HEIGHT                 |                                          |                          | 2.70 m         |

## NOTES

- A) Magnet access 2300x2500 mm
- B) Removable panels

## MAGNETIC PROXIMITY LIMITS

| Gauss (mT) Limit   | Equipment                                                                                                                                                                                                                                                                                                                                           |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 0.5 gauss (0.05mT) | Nuclear camera                                                                                                                                                                                                                                                                                                                                      |
| 1 gauss (0.1mT)    | Positron Emission Tomography scanner, Linear Accelerator, Cyclotrons, Accurate measuring scale, Analog image intensifiers, Bone Densitometers, Video display (tube), CT scanner, Ultrasound, Ultrasound, Electron microscope                                                                                                                        |
| 3 gauss (0.3mT)    | Power transformers, Main electrical distribution transformers                                                                                                                                                                                                                                                                                       |
| 5 gauss (0.5mT)    | Cardiac pacemakers, Neurostimulators, Biostimulation devices                                                                                                                                                                                                                                                                                        |
| 10 gauss (1mT)     | Magnetic computer media, Line printers, VCRs, Film processor, X-ray tubes, Emergency generators, Commercial laundry equipment, Food preparation area, Water cooling equipment, HVAC equipment, Major mechanical equipment room, Credit cards, watches, and clocks, Air conditioning equipment, Fuel storage tanks, Motors greater than 5 horsepower |
| 50 gauss (5mT)     | Metal detector for screening, LCD panels, Telephones                                                                                                                                                                                                                                                                                                |
| No Limit           | Digital Detectors                                                                                                                                                                                                                                                                                                                                   |


The customer must provide detail defining ferrous material below the magnet to the Project Manager so the GE Healthcare MRI Siting and Shielding team can review for compliance.

### STEEL MASS LIMITS TO MAGNET ISOCENTER (3x3 m [10x10 ft] AREA UNDER MAGNET)


| kg/m <sup>2</sup> | lb/ft <sup>2</sup> | Distance Below Top Surface Of Floor |       |
|-------------------|--------------------|-------------------------------------|-------|
|                   |                    | mm                                  | in    |
| 0                 | 0                  | 0 - 73                              | 0-3   |
| 9.8               | 2                  | 73 - 124                            | 3-5   |
| 14.7              | 3                  | 124 - 251                           | 5-10  |
| 39.2              | 8                  | 251 - 327                           | 10-13 |
| 98.0              | 20                 | 327+                                | 13+   |

The actual field strength can be affected by Magnetic shielding, Earth's magnetic field, other magnetic fields and stationary or moving metal. This information must be used to evaluate potential site interaction of GE Healthcare equipment with other non-GE Healthcare equipment. Magnetic shielding can be installed to prevent interaction between the magnet and nearby sensitive devices. The GE Healthcare Project Manager of Installation (PMI) can work with the customer to coordinate the magnetic shielding site evaluation. The customer is responsible for installation of all magnetic shielding.


## MAGNETIC FRINGE FIELD TOP VIEW

## MAGNETIC FRINGE FIELD FRONT VIEW

## MAGNETIC FRINGE FIELD SIDE VIEW

CONFIDENTIAL

SECRET DE AFACERI

EQUIPMENT LAYOUT TOP VIEW  
WITH MAGNETIC FIELD

| ITEM | DESCRIPTION                              | DIMENSIONS<br>LxWxH (mm) | WEIGHT<br>(kg) |
|------|------------------------------------------|--------------------------|----------------|
| 1    | MAGNET                                   | 1938x2450x2400           | 4553           |
| 2    | FIXED PATIENT TABLE (PT)                 | 2158x572x865             | 136            |
| 3    | MAGNET RUNDOWN UNIT (MRU)                | 287x206x172              | 3.2            |
| 4    | FURNITURE FOR COILS (NOT SUPPLIED BY GE) |                          |                |
| 5    | REAR PEDESTAL                            | 853x510x956              | 66             |
| 6    | INTEGRATED COOLING CABINET (ICC)         | 902x801x2048             | 545            |
| 7    | INTEGRATED SYSTEM CABINET (ISC)          | 1275x974.5x2000          | 1116           |
| 8    | MAGNET MONITOR (MON)                     | 381x260x127              | 4.5            |
| 9    | MAIN DISCONNECT PANEL (MDP)              | 216x508x914              | 43.5           |
| 10   | OPERATOR CONSOLE (GOC)                   |                          | 59.3           |
| 11   | PENETRATION PANEL (PP)                   | 255x243x723              | -              |

WALL - ACCORDING TO RECEIVED DRAWING

WALL TO DEMOLISH


RF SHIELD - 100 dB ATTENUATION


EXAM ROOM HEIGHT

FINISHED FLOOR TO SLAB HEIGHT

FALSE CEILING HEIGHT

CONFIDENTIAL  
SECRET DE AFACERI


CONFIDENTIAL  
SECRET DE AFACERI

EQUIPMENT LAYOUT FRONT VIEW  
WITH MAGNETIC FIELD

| ITEM | DESCRIPTION                              | DIMENSIONS<br>LxWxH (mm) | WEIGHT<br>(kg) |
|------|------------------------------------------|--------------------------|----------------|
| 1    | MAGNET                                   | 1938x2450x2400           | 4553           |
| 2    | FIXED PATIENT TABLE (PT)                 | 2158x572x965             | 136            |
| 3    | MAGNET RUNDOWN UNIT (MRU)                | 287x206x172              | 3.2            |
| 4    | FURNITURE FOR COILS (NOT SUPPLIED BY GE) | -                        | -              |
| 5    | REAR PEDESTAL                            | 853x510x956              | 66             |
| 6    | INTEGRATED COOLING CABINET (ICC)         | 902x801x2048             | 545            |
| 7    | INTEGRATED SYSTEM CABINET (ISC)          | 1275x5974.5x2000         | 1116           |
| 8    | MAGNET MONITOR (MON)                     | 381x260x127              | 4.5            |
| 9    | MAIN DISCONNECT PANEL (MDP)              | 216x508x914              | 43.5           |
| 10   | OPERATOR CONSOLE (GOC)                   | -                        | 59.3           |
| 11   | PENETRATION PANEL (PP)                   | 255x243x723              | -              |

- STRUCTURE - ACCORDING TO RECEIVED DRAWING
- WALL - ACCORDING TO RECEIVED DRAWING
- WALL TO DEMOLISH
- RF SHIELD - 100 dB ATTENUATION


- NOTES
- A) Note: Define RF shield's inset according to provisions made by the RF Shield vendor.
  - B) Warning! 5 Gauss line outside the Magnet room limits.

SECRET DE AFACERI


CONFIDENTIAL

FLOOR - ELECTRICAL LAYOUT

| ITEM         | QTY | DESCRIPTION                                                                                     |
|--------------|-----|-------------------------------------------------------------------------------------------------|
| 1            |     | VibroAcoustic Dampening Kit (see Floor Structural Detail)                                       |
| 2            |     | Patient table (PT) rebar free area                                                              |
| 3            |     | Vertical wall duct 200x60 from false ceiling to horizontal wall duct                            |
| 4            |     | 200x60 horizontal wall duct                                                                     |
| 5            | 1   | Magnet Rundown Unit (MRU)                                                                       |
| 6            | 1   | Main Disconnect Panel (MDP)                                                                     |
| Basic system |     |                                                                                                 |
| 2            |     | System emergency off (EO1, 2), (recommended height 1.50m-1.85m above floor)                     |
| 10           |     | Electrical outlets 230V 10/16A +G                                                               |
| 1            |     | Electrical outlet for Service 230V 10/16A +G and service lamp above false ceiling (recommended) |
| 1            |     | Service lamp in PEN closet (recommended)                                                        |
| 1            |     | Door interlock switch                                                                           |
| 3            |     | R4S network socket                                                                              |
| 2            |     | Emergency exhaust fan switch                                                                    |
| 1            |     | Electrical outlet 230V 10/16A +G, on uninterruptible power supply (if available)                |
| Wall duct    |     |                                                                                                 |


MAGNET ON VIBROACOUSTIC DAMPENING KIT

SCALE 1:25

MAGNET ROOM FLOOR SPECIFICATIONS


Magnet, Enclosure, and Patient Table areas must be flat and level within 3 mm (0.125 in) within the shaded area shown.

The finished floor must support the weight of all components (e.g., patient table, gradient coil replacement cart) throughout operation and service life.

SCALE 1:50

TABLE FRAME ANCHOR MOUNTING REQUIREMENTS

NOT TO SCALE


1. The RF Shield vendor must design and install the table frame anchor bolt.
2. The table frame anchor hole must be drilled after the Magnet is delivered.
3. The table frame anchor must not contact floor rebar or other structural steel.
4. The table frame anchor must electrically contact the RF shield at point of entry.
5. The table frame anchors must have the following properties:
  - a. Anchors must be two-part assembly (male/female).
  - b. Female side must be expansion- or epoxy-type.
  - c. Male side must be a bolt or threaded rod with appropriate-sized nut (bolt or rod must be removable - not epoxied or cemented in place).
  - d. Anchors must be electrically conductive.
  - e. Anchors must be non-magnetic.
  - f. Anchors must not induce galvanic corrosion with the RF shield.
  - g. Anchors must be commercially procured.
  - h. The anchor rod hole clearance in the table anchor base is 11 mm (0.43 in). The anchor rod diameter must be sized appropriately.
  - i. Anchors must meet the following clamping force: 2669 N (600 lbf).
  - j. The anchor rod must extend 40 mm ± 1.5 mm (1.57 in ± 0.05 in) above the finished floor.
  - k. The anchor rod must be at least 1.5 mm (0.06 in) in total length (length above the floor plus embedded length).
  - l. The RF Shield vendor must perform a pull test on the anchor (equal to the clamping force). Results must be provided to the GE Healthcare Project Manager of Installation (PMI).

CONFIDENTIAL  
SECRET DE AFACERI

**CONFIDENTIAL**  
**SECRET DE AFACERI**


# CEILING LAYOUT

| ITEM | DESCRIPTION                                                                                                  |
|------|--------------------------------------------------------------------------------------------------------------|
| 1    | He exhaust pipe (See Cryogenics Layout for details)                                                          |
| 2    | Magnet curtain kit                                                                                           |
| 3    | Horizontal cable ladder 450                                                                                  |
| 4    | Vertical cable ladder 200                                                                                    |
| 5    | Cable tray 300x60, above false ceiling                                                                       |
| 6    | Opening in the false ceiling 200x60 and vertical wall duct 200x60 from false ceiling to horizontal wall duct |
| 7    | Main Disconnect Panel (MDP)                                                                                  |

## CABLE TRAYS IN EQUIPMENT ROOM

CABLE TRAY DETAIL


| Minimum Cable Tray Width |     |             | For                             |
|--------------------------|-----|-------------|---------------------------------|
| From                     | To  | Width       |                                 |
| ISC                      | ICC | 229 [9 in]  | Electrical, Water               |
| ICC                      | PP  | 300 [12 in] | Electrical, Water, Helium lines |


NOT TO SCALE

## MINIMUM MAGNET CEILING HEIGHT (TOP VIEW)

Shaded area within solid lines indicates floor to ceiling minimum height of 2500 mm (8'-2.5"). If the ceiling height is between 2500 mm (8'-2.5") and 2667 mm (8'-9"), the flexible main lead extension for low ceiling height (2.5M Low Ceiling Kit-Passive, M7000GGM) is required for ramping the magnet.

## CABLE TRAYS REQUIREMENTS IN MAGNET ROOM


Cable Tray Requirements (90° Magnet Interface)

- 1 - Ceiling
- 2 - Finished Floor
- 3 - Magnet isocenter
- 4 - Minimum cable tray height for Groups A,B,C,D required at back of Magnet: 2449 mm [96.4 in].
- 5 - Maximum height from floor to top of tray (anywhere in Magnet room): 3251 mm [128 in].
- 6 - Minimum distance from top of cable tray to ceiling: 254 mm [10 in].
- 7 - Minimum distance on either side of an obstruction: 254 mm [10 in], minimum distance from top of cable tray to obstruction: 178 mm [7 in].
- 8 - Gradient cable support to isocenter: 1067 mm [42 in].
- 9 - Other cable support to isocenter: 864 ±12 mm [34 ±0.5 in].
- 10 - Minimum distance between trays: 12 mm [0.5 in].
- 11 - Non-ferrous cable support
- 12 - The center of the gradient cable group is 89 mm from the inside edge of the tray, in line with the magnet center.
- 13 - The center of the gradient cable group is 89 mm [3.5 in] from magnet center.
- 14 - Outer gradient cable is 450 mm [18 in] from edge of non-ferrous cable support.
- 15 - Distance between non-ferrous cable support: at least 305 mm [12 in].

CONFIDENTIAL  
SECRET DE AFACERI

CABLE CONCEALMENT

SCALE 1:25

CONFIDENTIAL  
SECRET DE AFACERI

## POWER REQUIREMENTS

### SPECIFICATIONS OF MAIN POWER INPUT

|                                   |                                                |
|-----------------------------------|------------------------------------------------|
| POWER SUPPLY                      | 380/400/415/480V ±7.5%/-10%, THREE-PHASE + GND |
| FREQUENCIES                       | 50/60Hz ± 3Hz                                  |
| POWER FACTOR                      | 0.9                                            |
| TOTAL SYSTEM 50 MILLISECOND POWER | 80KVA                                          |
| TOTAL SYSTEM CONTINUOUS POWER     | 64KVA                                          |

- Governing electrical codes may require a neutral wire. If present, neutral must be terminated in MDP.
- Power input must be separated from any others which may generate transients (elevators, air conditioning, radiology rooms equipped with high speed film changers, ...).
- Total harmonic distortion less than 2.5%. Phase imbalance must not exceed 2%.
- Lock-out/tag-out: The Main Disconnect Panel (MDP) shall provide an external single point lock-out/tag-out feature for the entire system and a means to externally lock-out/tag-out each output breaker independently. Each lock-out/tag-out feature shall accommodate a standard sized lock hasp.

### SPECIFICATIONS OF OPTIONAL BACK-UP POWER SUPPLY

MAGNET MONITOR REQUIRES A 110/220 VAC, 50/60 HZ, 3.0 A FACILITY SUPPLIED OUTLET. POWER AT THE OUTLET MUST BE CONTINUOUSLY AVAILABLE.

#### FOR CRYOCOOLER COMPRESSOR

|                   |                                                                                  |
|-------------------|----------------------------------------------------------------------------------|
| POWER INPUT       | 380/400/415/480V, THREE-PHASE + G                                                |
| POWER REQUIREMENT | MIN 9KVA                                                                         |
| POWER CONSUMPTION | MAX 7.2KW / STEADY STATE 6.5KW at 50Hz<br>MAX 8.3KW / STEADY STATE 7.5KW at 60Hz |
| FREQUENCY         | 50/60Hz ± 3Hz                                                                    |

Power to Cryocooler Compressor must be removed when emergency off circuit is actuated.

### CABLES

- Power and cable installation must comply with the distribution diagram.
- Size of the Main power input cable is determined by the customer, taking its length and admissible voltage drops into consideration.
- All cables must be isolated and flexible, cable color codes must comply with standards for electrical installation.
- The cables from signaling and remote control (V.Emergency Off Buttons,...) will go to Main Panel with a pigtail length of 1.5m (60in), and will be connected during installation.
- Each conductor will be identified and isolated (screw connector).

### GROUND SYSTEM


- The equipotential link will be by means of an equipotential bar.
- The grounding point of MDP is directly connected to the building's ground by an isolated copper cable.
- The impedance of the earth bar should be less than or equal to 2 ohms.

CONFIDENTIAL  
SECRET DE AFACERI


## TYPICAL MAGNET ROOM GROUNDING

### GROUNDING REQUIREMENTS

- All power lines into the RF shielded room require an RF filter.
- All electrical devices (for example, outlets, light fixtures, and so on) must have a ground wire from device power source and be grounded to the RF Shield at the RF Common Ground Stud.
- Resistance between any two grounded devices must not exceed 0.1 ohm to ensure equal potential ground system within the Magnet Room.
- Do not ground non-MR equipment to the MR ground system.
- If needed, electrical devices can be grounded at the System Cabinet rear panel.
- The common ground stud must be installed near the penetration point(s) of the GE equipment, into the RF shield between the Equipment Room and Magnet Room.
- For additional information refer to RF Shielded Room manual 5850260-1EN

## FACILITY SUPPLIED WIRING

This diagram displays minimum power requirements for GE equipment and should be used as a guide to determine appropriate wire sizes per local regulatory requirements.

CC Integrated Cooling Cabinet  
SC Integrated System Cabinet  
MDP Main Disconnect Panel  
PDU Power Distribution Unit (Inside ISC)  
PP Penetration Panel

**Notes:**

- (1) Refer to Power Distribution detail for more information
- (2) Size incoming wires from GE equipment according to CB sizes listed on
- (3) Power Distribution detail
- (4) Refer to Lighting Requirements detail
- (5) A cable is supplied by GE but may be extended if needed.

Cable SUPPLIED by CUSTOMER  
 Cable SUPPLIED by GE - Customer Installed  
 Equipment SUPPLIED by CUSTOMER  
 Equipment SUPPLIED by GE

Accepts following range of standard stranded conductors.  
All wire types, color and sizing to be selected in  
accordance with governing electrical code(s).

| GE MDP M70022MB (380V-415V with RCD) |            |           |            |          |  |
|--------------------------------------|------------|-----------|------------|----------|--|
| Item                                 | Phase      |           | Ground     |          |  |
|                                      | AWG/ktcmil | sq mm     | AWG/ktcmil | sq mm    |  |
| A                                    | 2 - 3/0    | 35 - 95   | 4 - 250    | 16 - 125 |  |
| B                                    | 6 - 3/0    | 10 - 95   | 4 - 250    | 16 - 125 |  |
| C                                    | 10 - 1     | 6 - 35    | 10 - 1     | 6 - 35   |  |
| D                                    | 12 - 6     | 2.5 - 10  | 6 - 4      | 13 - 21  |  |
| E                                    | 14 - 10    | 2.5 - 6   | 14 - 10    | 2.5 - 6  |  |
| F                                    | 22 - 17    | 0.5 - 2.5 |            |          |  |

**NOTES:**

- 1) Cryocooler Compressor (CIV) must operate 24 hours per day, 7 days per week.
- 2) Runners E0059, E3500 and N3527 are GE supplied cables. All other wiring is customer supplied and installed.
- 3) Two remote Emergency Off Maintained Buttons are supplied with the MDP. Emergency Off removes power from all outputs when activated.
- 4) All MDP output circuits drop out on loss of power. The Cryocooler (CIV) circuit will automatically restart upon restoration of power.
- 5) All circuit breakers have a short circuit current interrupting rating of 25,000 amps.
- 6) MDP is NRTL labeled.
- 7) All feeder circuits require dedicated ground wires.

## Power Distribution

Rev 0 | Date 29/NOV/2023

**CHAMPION.DWG**


## SIGNA CHAMPION

1000

11

---

INTERCONNECTIONS

| CABLES ROUTING FOR OPTIONS            |      |                          |               |    |
|---------------------------------------|------|--------------------------|---------------|----|
| Option                                | From | To                       | Cable Length  |    |
|                                       |      |                          | m             | ft |
| Magnetic Resonance Elastography (MRE) | MRE  | Magnet Isocenter         | Nominal 7.31  | 24 |
|                                       | MRE  | ISC cabinet              | Maximum 10.06 | 33 |
|                                       | MRE  | Ethernet Hub In ISC      | 15.24         | 50 |
|                                       | MRE  | Customer Supplied Outlet | 15.24         | 50 |
|                                       |      |                          | 60 Hz         | 20 |
|                                       |      |                          | 50 Hz         | 25 |

| GRADIENT CABLE LENGTH OPTIONS |                          |          |
|-------------------------------|--------------------------|----------|
| Length Identifier             | Available lengths m (ft) | Proposed |
| L9 (Magnet Room)              | 4.5 (14.8)               | -        |
|                               | 6.5 (21.3)               | -        |
|                               | 8.5 (27.9)               | -        |
|                               | 10.5 (34.4)              | -        |
|                               | 12.5 (41.0)              | -        |

GENERAL NOTE: PMI must validate proposed selectables and send confirmation to OTR. Refer to MyProjects if Proposed field is blank.

GOLDSEAL/SILVER PREFERRED NOTE: Cable lengths listed may differ from what is shipped with the system. Contact the Goldseal group for actual lengths to be delivered.

ROOM MOVE NOTE: Cable lengths listed may differ from what is included with reinstalled system. Contact the local field engineer for actual lengths to be delivered.

## LIGHTING REQUIREMENTS

- All lighting fixtures and associated components must meet all RF shielded room and RF grounding requirements (e.g., track lighting is not recommended due to possible RF noise).
- All removable lighting fixtures and associated components must be non-magnetic.
- All lighting must use direct current (the DC must have less than 5% ripple).
- 300 lux must be provided at the front of the magnet for patient access and above the magnet for servicing.
- Fluorescent lighting must not be used in the magnet room.
- Lighting must be adjusted using a discrete switch or a variable DC lighting controller.
- SCR dimmers or rheostats must not be used.
- DC LED lighting may be used if the DC power converter and RF sources are all located outside the magnet room RF shield.
- NOTE: LED lighting could cause image quality issues due to RF interference. Make sure a MR-compatible LED lighting solution is chosen.
- Battery chargers (e.g., used for emergency lighting) must be located outside the magnet room.
- LED Lighting or short filament length incandescent bulbs are recommended.
- Linear lamps are not recommended due to the high burnout rate.

## CONNECTIVITY REQUIREMENTS

Your new GE Healthcare imaging modality will require local and remote connectivity to enable our full range of digital support:

- Local connectivity - This allows your system to connect to local devices such as PACS and modality worklist. We will require network information to configure the system(s), and a live ethernet port(s) prior to the delivery of the system(s).
- Remote connectivity - Your GE Healthcare service warranty includes InSite™ (applicable to InSite capable products), a powerful broadband-based service which enables digital tools that can help guard your hospital against equipment downtime and revenue loss by quickly connecting you to a GE Healthcare expert.

Depending on product family and software version, imaging systems can be connected in one of the following methods:

1. TLS over TCP Port 443 (Preferred method for new products) via:
  - a. DNS resolution
  - b. Customer-provided Proxy or
  - c. GE Proxy (Available in some regions)
2. Site-to-Site IPsec VPN tunnel

Please provide the GE project manager with the contact information for the resource that can provide information required to set up these connections. GEHC will send out communication to these contacts, which will include the project's Connectivity requirements, and a Connectivity form. This form will need to be completed and returned to GEHC prior to delivery of the system to ensure the system is tested and connectivity is enabled prior to the completion of the installation.

CONFIDENTIAL  
SECRET DE AFACERI

## TEMPERATURE AND HUMIDITY SPECIFICATIONS

### IN-USE CONDITIONS

|                       | MAGNET ROOM | CONTROL ROOM | EQUIPMENT ROOM |
|-----------------------|-------------|--------------|----------------|
| Temperature           | 15-21°C     | 15-32°C      | 15-32°C        |
|                       | 59-69.8°F   | 59-89.6°F    | 59-89.6°F      |
| Temperature gradient  | ≤ 3°C/h     | ≤ 3°C/h      | ≤ 3°C/h        |
| [1]                   | ≤ 5°F/h     | ≤ 5°F/h      | ≤ 5°F/h        |
| Relative humidity     | 30% to 60%  | 30% to 75%   | 30% to 75%     |
| Humidity gradient [2] | ≤ 5%/h      | ≤ 5%/h       | ≤ 5%/h         |

#### NOTE

Maximum ambient temperature is derated by 1°C [33.8°F] per 300 m [984 ft] above 950 m [3112 ft] (not to exceed 2600 m [8530 ft]).

(1) Operating temperature gradient limits shall be between -3°C/hr (-5°F/hr) and 3°C/hr (5°F/hr), when averaged over 1 hour.

(2) Operating humidity gradient limits shall be between -5% RH/hour and 5% RH/hour, when averaged over 1 hour.

### AIR RENEWAL

According to local standards.

### STORAGE CONDITIONS

|                       |                |
|-----------------------|----------------|
| Temperature           | -30°C to 55°C  |
|                       | -22°F to 131°F |
| Temperature gradient  | ≤ 20°C/h       |
|                       | ≤ 68°F/h       |
| Relative humidity [1] | 10% to 90%     |
| Humidity gradient     | ≤ 30%/h        |

(1) non-condensing

Material should not be stored for more than 90 days.

NOTE  
In case of using air conditioning systems that have a risk of water leakage it is recommended not to install it above electric equipment or to take measures to protect the equipment from dripping water.

## HEAT DISSIPATION DETAILS

| DESCRIPTION                         | ROOM      | IDLE |        | AVERAGE |        | MAX  |        |
|-------------------------------------|-----------|------|--------|---------|--------|------|--------|
|                                     |           | W    | BTU/hr | W       | BTU/hr | W    | BTU/hr |
| Magnet (MAG) and Patient Table (PT) | Magnet    | 561  | 1915   | 1200    | 4095   | 2400 | 8189   |
| Penetration Panel (PP)              | Equipment | 0    | 0      | 0       | 0      | 0    | 0      |
| Main Disconnect Panel (MDP)         | Equipment | 132  | 450    | 132     | 450    | 264  | 901    |
| Integrated System Cabinet (ISC)     | Equipment | 280  | 955    | 2130    | 7270   | 7100 | 24232  |
| Integrated Cooling Cabinet (ICC)    | Equipment | 0    | 0      | 500     | 1706   | 1000 | 3412   |
| Cryocooler Compressor (CRY)         | Equipment | 500  | 1706   | 500     | 1706   | 500  | 1706   |
| Magnet Monitor (MON)                | Equipment | 60   | 205    | 60      | 205    | 60   | 205    |
| Operator Workspace equipment (OW)   | Control   | 1450 | 4947   | 1450    | 4947   | 1450 | 4947   |

**CONFIDENTIAL**  
**SECRET DEAFACER**

## MAGNET ROOM VENTING REQUIREMENTS


### HVAC VENT REQUIREMENTS

- HVAC vendor must comply with Magnet room temperature and humidity specifications and RF shielding specifications.
- RF Shield vendor must install open pipe or honeycomb HVAC waveguides.
- All serviceable parts in the Magnet room (e.g.: diffusers) must be non-magnetic.
- Waveguides must be nonmagnetic and electrically isolated.
- Incoming air must contain at least 5% air from outside the Magnet room (inside or outside the facility) to displace residual helium.

### EMERGENCY VENT REQUIREMENT


- Exhaust vent system is supplied by the customer.
- All items within the RF enclosure must be non-magnetic.
- The exhaust vent system must be tested and operational before the magnet is installed.
- The exhaust intake vent must be located near the magnet cryogenic vent at the highest point on the finished or drop ceiling.
- The Magnet room exhaust fan and exhaust intake vent must have a capacity of at least **1200 CFM (34 m³/min)** with a minimum of **12 room air exchanges per hour**.
- The exhaust fan must be placed above RF shielding located outside 10 gauss (1mT) and with appropriate waveguide.
- The system must have a manual exhaust fan switch near the Operator Workspace and in the Magnet room near the door (the switches must be connected in parallel).
- All system components must be accessible for customer inspection, cleaning and maintenance

## EMERGENCY EXHAUST VENT


# HVAC - CHILLED WATER LAYOUT

| ITEM | DESCRIPTION                                                         |
|------|---------------------------------------------------------------------|
| 1    | Air diffusers (number and position to be defined)                   |
| 2    | A/C in/out waveguide (number and position to be defined)            |
| 3    | Emergency air exhaust waveguide (number and position to be defined) |
| 4    | Integrated Cooling Cabinet (ICC)                                    |
| 5    | Water backup                                                        |
| 6    | Facility chiller                                                    |
|      | Opening in the false ceiling / wall                                 |
|      | Opening in the RF cage                                              |

**CONFIDENTIAL**  
**SECRET DE AFACERI**

CHILLED WATER BLOCK DIAGRAM


LIQUID COOLANT SPECIFICATIONS

| PARAMETER                                                     | REQUIREMENTS                                                                                                                                                                                                                                                                                                           |
|---------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Availability                                                  | Continuous                                                                                                                                                                                                                                                                                                             |
| Antifreeze                                                    | 0-40% propylene glycol water (PGW) or ethylene glycol water (EGW)                                                                                                                                                                                                                                                      |
| Minimum Flow                                                  | 50 L/min [13.2 gpm]                                                                                                                                                                                                                                                                                                    |
| Maximum Flow                                                  | 80 L/min [21.1 gpm]                                                                                                                                                                                                                                                                                                    |
| Typical Flow                                                  | 60 L/min [15.8 gpm]                                                                                                                                                                                                                                                                                                    |
| Maximum Pressure Drop in ICC at Minimum Flow                  | 1.8 bar [26.1 psi] with 40% PGW, 1021 kg/m <sup>3</sup> [63.7 lbs/ft <sup>3</sup> ] density<br>1.3 bar [19.3 psi] with pure water, 994 kg/m <sup>3</sup> [62.1 lbs/ft <sup>3</sup> ] density                                                                                                                           |
| Maximum Pressure Drop in ICC at Maximum Flow                  | 3.4 bar [49.3 psi] with 40% PGW, 1021 kg/m <sup>3</sup> [63.7 lbs/ft <sup>3</sup> ] density<br>2.5 bar [36.5 psi] with pure water, 994 kg/m <sup>3</sup> [62.1 lbs/ft <sup>3</sup> ] density                                                                                                                           |
| Maximum inlet pressure to ICC                                 | 5.52 bar [80 psi]                                                                                                                                                                                                                                                                                                      |
| Chiller size                                                  | Minimum 36 kW                                                                                                                                                                                                                                                                                                          |
| Condensation protection                                       | Facility plumbing to the ICC must be properly routed and insulated to prevent equipment damage or safety hazards                                                                                                                                                                                                       |
| Minimum continuous heat load                                  | 7.5 kW                                                                                                                                                                                                                                                                                                                 |
| Inlet temperature to ICC                                      | 5 to 15°C [41 to 59° F]                                                                                                                                                                                                                                                                                                |
| Customer supplied feeder hose (from main water supply to ICC) | 25.4 mm [1 in] minimum hose inside diameter.<br>If the hose is longer than 10 m [32.8 ft] and shorter than 30 m [98.4 ft], a hose with a 31.75 mm [1.25 in] minimum ID is recommended. If using 31.75 mm [1.25 in] hose, an adaptor is required to reduce the inner diameter to 25.4 mm [1 in] for the ICC connection. |
| Water quality                                                 | Refer to pre-installation manual for detailed specifications                                                                                                                                                                                                                                                           |


CITY WATER BACKUP SPECIFICATIONS FOR COMPRESSOR

INLET WATER FLOW/TEMPERATURE FOR CRYOCOOLER COMPRESSOR  
\*RECOMMENDED TEMPERATURE AND WATER FLOW TO CRY ENTRY

|                       | MIN                    | MAX                   | IDEAL |
|-----------------------|------------------------|-----------------------|-------|
| INLET TEMP (°C)       | 4                      | 28                    | 12-16 |
| INLET FLOW (l/min)    | 4                      | 10                    | 7-8   |
| INLET PRESSURE (kPa)  | 200                    | 690                   |       |
| TEMP RISE             | 26°C at 4 l/min flow   | 10°C at 10 l/min flow |       |
| HEAT DISSIPATION (kW) | 7.2 kW                 |                       |       |
| PRESSURE DROP         | 60 kPa at 8 l/min flow |                       |       |

MINIMUM RECOMMENDED CHILLER CONFIGURATION


- FLOW METER  
VALVE  
FILTER
- Equipment supplied by GEHC  
Customer/contractor supplied and installed


CONFIDENTIAL  
SECRET DE AFACERI

CONFIDENTIAL  
SECRET DE AFACERI


CRYOGENICS TOP VIEW

CRYOGENICS SIDE VIEW

CRYOGENICS FRONT VIEW

**SECRET DE AFACERI**

1. The 203 mm [8 in] OD vent material must be one of the following materials with the wall thickness indicated:

- NOTE**  
All welds on the pipe must be ground down to a smooth 203 mm [8 in] diameter so that it can be clamped to the Ventglas with enough force.


4. If required, bellows pipe less than 300 mm [12 in] in length may be used as a thermal expansion joint.

5. The vent pipe must withstand the maximum pressure listed in the Pre-Installation Manual.
6. Waveguide vent material must match the outside diameter of the magnet flanged vent adapter.


### KEY COMPONENTS:

- RF waveguide extended from wall to magnet adapter
- Must be all same material and all welded or brazed
- Support system must withstand 8229 N [1850 lbs]
- GE venglas must be installed in vertical section directly over magnet
- Vent cap must prevent ingress of weather elements due to strong wind (for example, rain, snow, hail, sand, and so on) and foreign material debris (for example, leaves, bird nests, and so on).

**NOTE:** for the Dielectric Break outside the Magnet room please see the Typical Cryogen Vent Pipe Detail

## CRYOGENIC VENTING ( EXTERIOR )

## NOTES

- (1) Restricted area: minimum distance between vent pipe and ground is 3660mm [144 in]. Barriers are required.
- (2) Public area: barriers are not required if height is > 5000mm [197 in].
- (3) The bottom of the 90° elbow must be at least 914mm [36 in] above the roof deck (or higher if at risk of being blocked by drifting snow, sand, or other potential obstructions.)

**NOT TO SCALE**

# MAGNET CRYOGENIC VENT SYSTEM PRESSURE DROP MATRIX

| Outer dia. of pipe (D)           | Distance of vent system component magnet | Pressure drop for straight pipe |       | Std sweep 45° elbow |       | Long sweep 45° elbow |       | Std sweep 90° elbow |       | Long sweep 90° elbow |       | 90° miter |       |
|----------------------------------|------------------------------------------|---------------------------------|-------|---------------------|-------|----------------------|-------|---------------------|-------|----------------------|-------|-----------|-------|
|                                  |                                          | psi/ft                          |       | kPa/m               |       | psi                  |       | kPa                 |       | psi                  |       | kPa       |       |
|                                  |                                          | ft                              | m     | ft                  | m     | ft                   | m     | ft                  | m     | ft                   | m     | ft        | m     |
| 6 in. (inside dia. $\pm 14mm$ )  | 10                                       | 3.048                           | 0.146 | 3.311               | 0.156 | 5.544                | 0.256 | 3.696               | 0.169 | 10.394               | 0.476 | 6.93      | 0.314 |
|                                  | 20                                       | 6.096                           | 0.292 | 6.622               | 0.312 | 11.088               | 0.512 | 7.392               | 0.336 | 20.788               | 0.936 | 13.86     | 0.628 |
|                                  | 30                                       | 9.144                           | 0.437 | 9.933               | 0.468 | 16.626               | 0.768 | 11.088              | 0.512 | 31.576               | 1.408 | 20.788    | 0.936 |
|                                  | 40                                       | 12.192                          | 0.583 | 13.244              | 0.624 | 22.176               | 1.024 | 14.784              | 0.672 | 42.768               | 1.920 | 27.712    | 1.248 |
|                                  | 50                                       | 15.240                          | 0.729 | 16.556              | 0.780 | 28.224               | 1.272 | 17.976              | 0.816 | 53.760               | 2.376 | 35.616    | 1.584 |
|                                  | 60                                       | 18.288                          | 0.875 | 19.864              | 0.940 | 34.208               | 1.520 | 22.960              | 0.976 | 64.800               | 2.880 | 43.200    | 1.920 |
|                                  | 70                                       | 21.336                          | 1.021 | 23.008              | 1.096 | 40.192               | 1.768 | 28.000              | 1.216 | 74.880               | 3.376 | 51.840    | 2.304 |
|                                  | 80                                       | 24.384                          | 1.167 | 26.152              | 1.260 | 46.176               | 2.016 | 34.000              | 1.464 | 80.864               | 3.664 | 58.880    | 2.640 |
|                                  | 90                                       | 27.432                          | 1.313 | 29.304              | 1.416 | 52.160               | 2.264 | 40.000              | 1.712 | 86.848               | 3.952 | 65.920    | 2.976 |
|                                  | 100                                      | 30.480                          | 1.459 | 32.448              | 1.576 | 58.144               | 2.512 | 46.000              | 1.960 | 92.832               | 4.240 | 72.960    | 3.312 |
| 8 in. (inside dia. $\pm 19mm$ )  | 10                                       | 3.048                           | 0.146 | 3.311               | 0.156 | 5.544                | 0.256 | 3.696               | 0.169 | 10.394               | 0.476 | 6.93      | 0.314 |
|                                  | 20                                       | 6.096                           | 0.292 | 6.622               | 0.312 | 11.088               | 0.512 | 7.392               | 0.336 | 20.788               | 0.936 | 13.86     | 0.628 |
|                                  | 30                                       | 9.144                           | 0.437 | 9.933               | 0.468 | 16.626               | 0.768 | 11.088              | 0.512 | 31.576               | 1.408 | 20.788    | 0.936 |
|                                  | 40                                       | 12.192                          | 0.583 | 13.244              | 0.624 | 22.176               | 1.024 | 14.784              | 0.672 | 42.768               | 1.920 | 27.712    | 1.248 |
|                                  | 50                                       | 15.240                          | 0.729 | 16.556              | 0.780 | 28.224               | 1.272 | 17.976              | 0.816 | 53.760               | 2.376 | 35.616    | 1.584 |
|                                  | 60                                       | 18.288                          | 0.875 | 19.864              | 0.940 | 34.208               | 1.520 | 22.960              | 0.976 | 64.800               | 2.880 | 43.200    | 1.920 |
|                                  | 70                                       | 21.336                          | 1.021 | 23.008              | 1.096 | 40.192               | 1.768 | 28.000              | 1.216 | 74.880               | 3.376 | 51.840    | 2.304 |
|                                  | 80                                       | 24.384                          | 1.167 | 26.152              | 1.260 | 46.176               | 2.016 | 34.000              | 1.464 | 80.864               | 3.664 | 58.880    | 2.640 |
|                                  | 90                                       | 27.432                          | 1.313 | 29.304              | 1.416 | 52.160               | 2.264 | 40.000              | 1.712 | 86.848               | 3.952 | 65.920    | 2.976 |
|                                  | 100                                      | 30.480                          | 1.459 | 32.448              | 1.576 | 58.144               | 2.512 | 46.000              | 1.960 | 92.832               | 4.240 | 72.960    | 3.312 |
| 10 in. (inside dia. $\pm 24mm$ ) | 10                                       | 3.048                           | 0.146 | 3.311               | 0.156 | 5.544                | 0.256 | 3.696               | 0.169 | 10.394               | 0.476 | 6.93      | 0.314 |
|                                  | 20                                       | 6.096                           | 0.292 | 6.622               | 0.312 | 11.088               | 0.512 | 7.392               | 0.336 | 20.788               | 0.936 | 13.86     | 0.628 |
|                                  | 30                                       | 9.144                           | 0.437 | 9.933               | 0.468 | 16.626               | 0.768 | 11.088              | 0.512 | 31.576               | 1.408 | 20.788    | 0.936 |
|                                  | 40                                       | 12.192                          | 0.583 | 13.244              | 0.624 | 22.176               | 1.024 | 14.784              | 0.672 | 42.768               | 1.920 | 27.712    | 1.248 |
|                                  | 50                                       | 15.240                          | 0.729 | 16.556              | 0.780 | 28.224               | 1.272 | 17.976              | 0.816 | 53.760               | 2.376 | 35.616    | 1.584 |
|                                  | 60                                       | 18.288                          | 0.875 | 19.864              | 0.940 | 34.208               | 1.520 | 22.960              | 0.976 | 64.800               | 2.880 | 43.200    | 1.920 |
|                                  | 70                                       | 21.336                          | 1.021 | 23.008              | 1.096 | 40.192               | 1.768 | 28.000              | 1.216 | 74.880               | 3.376 | 51.840    | 2.304 |
|                                  | 80                                       | 24.384                          | 1.167 | 26.152              | 1.260 | 46.176               | 2.016 | 34.000              | 1.464 | 80.864               | 3.664 | 58.880    | 2.640 |
|                                  | 90                                       | 27.432                          | 1.313 | 29.304              | 1.416 | 52.160               | 2.264 | 40.000              | 1.712 | 86.848               | 3.952 | 65.920    | 2.976 |
|                                  | 100                                      | 30.480                          | 1.459 | 32.448              | 1.576 | 58.144               | 2.512 | 46.000              | 1.960 | 92.832               | 4.240 | 72.960    | 3.312 |

## Notes

- Elbows with angles greater than 90 deg must not be used
- Data in Table is based on the following facts and assumptions:
  - Initial flow conditions at magnet interface
  - EM energy (13MJ) is dumped to He during quench and rises He temperature to 10 Kelvin
  - Gas temperature starting at 10 Kelvin and increase with length determined by thermal energy balance
  - 90% He is assumed to be evacuated within 30 sec. None left after quench.
  - Absolute roughness is assumed to be 0.25 mm.
  - $R/D = 1.0$  for standard sweep elbows,  $R/D = 1.5$  for long sweep elbows where D = outer diameter of pipe; R = radius of bend
- The total pressure drop of the entire cryogenic vent system must be less than 20 psi (137.9 kPa). The calculation starts at the magnet vent interface and ends at the termination point outside the building.

SECRET DE AFACERI


CONFIDENTIAL

## PENETRATION PANEL WALL OPENINGS


### View from equipment room side.

Refer to Equipment Layout page for site specific orientation of Penetration Panel.  
Integrated System Cabinet must be covered on both sides for safety.


#### OPTION 2: PENETRATION PANEL ABOVE ISC

#### OPTION 1: PENETRATION PANEL, EITHER SIDE OF ISC

## PENETRATION PANEL CLOSET

CONFIDENTIAL


SECRET DE AFACERI

## ACOUSTICS SPECIFICATIONS

Acoustic and vibroacoustic information is provided for site planning and architectural design activities. It is the customer's responsibility to hire a qualified acoustic engineer for solutions to further attenuate this transmitted noise and vibration, if required. The actual room noise level may vary based on room design, optional equipment, and usage:

Control Room: 62dBA  
Equipment Room: 80dBA  
Magnet Room: 122dBA\*  
(maximum sound pressure level at magnet bore isocenter)


\* Frequency: 20 Hz to 20kHz

| FREQUENCY (Hz) | AMPLITUDE (m/s²) |
|----------------|------------------|
| 2              | 0.05             |
| 10             | 0.05             |
| 20             | 0.15             |
| 30             | 0.2              |
| 35             | 0.2              |
| 50             | 0.2              |

## SOUND PRESSURE SPECTRAL DISTRIBUTION

1/3 Band Relative SPL


## VIBRATION SPECIFICATIONS


Excessive vibration can affect MR image quality. Vibration testing must be performed early in the site planning process to ensure vibration is minimized. Both steady state vibration (exhaust fans, air conditioners, pumps, etc.) and transient vibrations (traffic, pedestrians, door slamming, etc.) must be assessed. The magnet cannot be directly isolated from vibration. Any vibration issue must be resolved at the source.

Transient vibration levels above the specified limits in the MR Site Vibration Test Guidelines must be analyzed. Any transient vibration that causes vibration to exceed the steady-state level must be mitigated.

MAGNET STEADY-STATE VIBRATION SPECIFICATIONS

## VIBRATION TRANSMITTED THROUGH VIBROACOUSTIC MAT

CONFIDENTIAL

SECRET DE AFACERI


## DELIVERY

### ROUTING


- The customer is solely liable for routing of components from dock to final site.
- GE must be able to move system components in or out with no need to uncrate or disassemble any of the components. The entire passageway must be cleared, adequately lighted and free from dust.
- The floor and its surfacing must be able to withstand the live load of components and handling equipment.
- Floor surfacing must be continuous.
- The customer must protect any fragile flooring surfaces.

### SPECIFICATIONS FOR MAGNET ROUTING

- Floor must be able to withstand a moving load of 4553kg [10,038 lb]
- Recommended opening height: 2.5m [98.5in], width: 2.5m [98.5in]. If recommended dimensions cannot be met refer to pre-installation manual for detailed specifications.
- Maximum slope: 30°

otes: Width can be reduced by 56 [2.2 in] if right center enclosure is removed.  
Short lifting beams are available for small rooms. L=2000 [78.7 in], available as GE part no. 5796140.  
PATH WITH 90 DEGREE TURN  
STRAIGHT PATH

### STORAGE CONDITIONS


- System components except the magnet should be stored in a cleaned room.
- Temperature = -30 to 55°C [-22 to 131°F], relative humidity 10-90% non condensing.
- Phantoms and coolant must be kept from freezing.
- The magnet will be delivered after GE validation of the site.

### INSTALLATION AND DELIVERY ACCEPTANCE


- A survey of the site established by the customer and GE will make the decision for the delivery time.
- This survey of the site (a form is made available by GE) is only to check if the apparent conditions of the site allow the equipment to be delivered.
- If the site is not ready, GE can delay the delivery time.

## MINIMUM MAGNET CEILING HEIGHT (TOP VIEW)

Shaded area within solid lines indicates floor to ceiling minimum height of 2500 mm (8'-2.5"). If the ceiling height is between 2500 mm (8'-2.5") and 2667 mm (8'-9"), the flexible main lead extension for low ceiling height (2.5M Low Ceiling Kit-Passive, M7000GM) is required for ramping the magnet.

## GRADIENT COIL REPLACEMENT



| EQUIPMENT                                                         | DIMENSIONS LxWxH |              | WEIGHT |      | NOTE                                                                                                                                             |
|-------------------------------------------------------------------|------------------|--------------|--------|------|--------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                   | mm               | in           | kg     | lbs  |                                                                                                                                                  |
| Replacement VRMw gradient coil assembly on a shipping cradle/cart | 910x2444x1499    | 35.8x96.2x59 | 1449   | 3194 | Initial gradient coil assembly is shipped installed in the magnet. Shipping/installation cart is used to install replacement coil assembly only. |
|                                                                   |                  |              |        |      |                                                                                                                                                  |

The weight bearing structure of the site should support any additional weight of the main replacement parts occurring during maintenance of the magnet, throughout the whole lifecycle of the MR.

SECRET DE AFACERI


CONFIDENTIAL

MAGNET ENCLOSURE




NOTE:  
Enclosure dimensions are for reference only, NOT FOR SITE PLANNING USE.  
● Center of gravity  
DETAIL NOT TO SCALE

LOW HEIGHT FIXED TABLE


INTEGRATED COOLING CABINET (ICC)



| PICC WATER SUPPLY AND RETURN CONNECTIONS |                                                    |
|------------------------------------------|----------------------------------------------------|
| Item                                     | Description                                        |
| 1                                        | Facility return                                    |
| 2                                        | Facility supply                                    |
| 3                                        | City Water Return for Cryocooler compressor backup |
| 4                                        | City Water Supply for Cryocooler compressor backup |

INTEGRATED SYSTEM CABINET


GLOBAL OPERATOR CABINET (GOC)

MAIN DISCONNECT PANEL


MAGNET MONITOR AND MAGNET RUNDOWN UNIT (MRU)

SCALE 1:20

PENETRATION PANEL

Equipment Dimensions (2)

SIGNA CHAMPION

EN-MRI-TYP-SIGNA\_CHAMPION.DWG

1:20 | Rev 0 | Date 29/NOV/2023

26/27

## DISCLAIMER

### GENERAL SPECIFICATIONS

- GE is not responsible for the installation of developers and associated equipment, lighting, cassette trays and protective screens or derivatives not mentioned in the order.
- The final study contains recommendations for the location of GE equipment and associated devices, electrical wiring and room arrangements. When preparing the study, every effort has been made to consider every aspect of the actual equipment expected to be installed.
- The layout of the equipment offered by GE, the dimensions given for the premises, the details provided for the pre-installation work and electrical power supply are given according to the information noted during on-site study and the wishes expressed by the customer.
- The room dimensions used to create the equipment layout may originate from a previous layout and may not be accurate as they may not have been verified on site. GE cannot take any responsibility for errors due to lack of information.
- Dimensions apply to finished surfaces of the room.
- Actual configuration may differ from options presented in some typical views or tables.
- If this set of final drawings has been approved by the customer, any subsequent modification of the site must be subject to further investigation by GE about the feasibility of installing the equipment. Any reservations must be noted.
- The equipment layout indicates the placement and interconnection of the indicated equipment components. There may be local requirements that could impact the placement of these components. It remains the customer's responsibility to ensure that the site and final equipment placement complies with all applicable local requirements.
- All work required to install GE equipment must be carried out in compliance with the building regulations and the safety standards of legal force in the country concerned.
- These drawings are not to be used for actual construction purposes. The company cannot take responsibility for any damage resulting therefrom.

### CUSTOMER RESPONSIBILITIES

- It is the responsibility of the customer to prepare the site in accordance with the specifications stated in the final study. A detailed site readiness checklist is provided by GE. It is the responsibility of the customer to ensure all requirements are fulfilled and that the site conforms to all specifications defined in the checklist and final study. The GE Project Manager of Installation (PMI) will work in cooperation with the customer to follow up and ensure that actions in the checklist are complete, and if necessary, will aid in the rescheduling of the delivery and installation date.
- Prior to installation, a structural engineer of record must ensure that the floor and ceiling is designed in such a way that the loads of the installed system can be securely borne and transferred. The layout of additional structural elements, dimensioning and the selection of appropriate installation methods are the sole responsibility of the structural engineer. Execution of load bearing structures supporting equipment on the ceiling, floor or walls are the customer's responsibility.

THE UNDERSIGNED, HEREBY CERTIFIES THAT I HAVE READ AND APPROVED THE PLANS IN THIS DOCUMENT.

| DATE | NAME | SIGNATURE |
|------|------|-----------|
|      |      |           |

SIGNA CHAMPION | EN-MRI-TYP-SIGNA. CHAN


## CUSTOMER SITE READINESS REQUIREMENTS

| REQUIRED MANUALS FOR SYSTEM PRE-INSTALLATION                                                                                                                                              |                     |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|
| Description                                                                                                                                                                               | Document Number*    |
| Product specific Pre-Installation Manual                                                                                                                                                  | Refer to cover page |
| Magnet Room Venting                                                                                                                                                                       | 5850263             |
| RF Shielded Room Pre-Installation Requirements for MR systems                                                                                                                             | 5850260             |
| IEC Electromagnetic Compatibility                                                                                                                                                         | 5850261             |
| Acoustic Room Details                                                                                                                                                                     | 5850262             |
| Magnet Venting Conformance Assessment Form                                                                                                                                                | 2705036             |
| *documents can be accessed in multiple languages at <a href="https://customer-doc.cloud.gehealthcare.com/#/cdp/dashboard">https://customer-doc.cloud.gehealthcare.com/#/cdp/dashboard</a> |                     |

- A mandatory component of this drawing set is the GE Healthcare Pre-installation manual. Failure to reference the Pre-installation manual will result in incomplete documentation required for site design and preparation.
- The items on the GE Healthcare Site Readiness Checklist **DOC1809666** are **REQUIRED** to facilitate equipment delivery to the site. Equipment will not be delivered if these requirements are not satisfied.
- Any deviation from these drawings must be communicated in writing to and reviewed by your local GE Healthcare installation project manager prior to making changes.
- Make arrangements for any rigging, special handling, or facility modifications that must be made to deliver the equipment to the installation site. If desired, your local GE Healthcare installation project manager can supply a reference list of rigging contractors.
- New construction requires the following:
  1. Secure area for equipment,
  2. Power for drills and other test equipment,
  3. Restrooms.
- Provide for refuse removal and disposal (e.g. crates, cartons, packing)
- It is required to minimize vibrations within the scan room. It is the customer's responsibility to contract a vibration consultant/engineer to implement site design modifications to meet the GE vibration specification. Refer to the system Pre-installation manual for vibration specifications.

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

[illegible]

ing and shedding of possibly exsiccating rainforest

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