
The challenge of clots in POC blood gas testing

In critical care, accurate blood gas analyzer measurements are important for making informed treatment decisions fast. But even when the most rigorous preanalytical precautions are followed,

blood clots can still be introduced into an analyzer via a sample, and may potentially affect measurement, calibration or QC. Macro Clot Check and Micro Clot Check work in tandem to reduce this risk.

MACRO CLOT CHECK uses advanced flow monitoring to detect macro clots that disrupt flow and can affect sensor performance.

MICRO CLOT CHECK uses continuous sensor response analysis to detect and remove micro clots that do not affect flow but which can attach on or near sensors and potentially affect sensor performance.

Advanced Clot Checks and Removal

Macro Clot Check and Micro Clot Check combine to deliver Advanced Clot Checks and Removal on the ABL90 FLEX and FLEX PLUS. Micro Clot Check integrates three novel clot detection methods.

1 BACKLOG CLOT DETECTION

Detects if a micro clot has attached to a sensor after a measurement.

2 IDLE CLOT DETECTION

Detects clots that affect sensors over time, typically small clots on or near sensors.

3 AIR/LIQUID CLOT DETECTION

Detects clots by comparing actual and expected sensor response when measuring on air or liquid.

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