CompoDock

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Table of Contents

1 Important information

1.1	Organization of the Operating Instructions	3
1.2	How to use the Operating Instructions	3
1.3	Brief description	4
1.4	Intended use	5
1.4.1	Fields of application	5
1.5	Standard accessories	5
1.6	Classification	5
1.7	Maintenance and cleaning	6
1.8	Repair	6
1.9	Safety precautions	6
1.10	Manufacturer's address	10

2 Design

2.1	Front view	11
2.2	Rear view / ports	12

3 Sealing principle

4 Operation

4.1	CompoDock Counter	17
4.2	Establishing a sterile connection	18
4.3	Turning the device off	23

5 Accessories (options)

5.1	CompoDock Counter	24
5.2	Side tables	24
5.3	Bag trays (included in the scope of delivery)	24
5.4	Accessories required for use with DockMaster Net	24

6 Cleaning / disinfection

6.1	Surface cleaning	26
6.2	Surface disinfection	26
6.3	Cleaning the electrodes	27
6.4	Cleaning the tube guides	27

7 Errors and alarms

7.1	General notes	28
-----	---------------------	----

8 Initial start-up

9 Description of the device

9.1	Specifications	34
9.1.1	Dimensions, weight and housing material	34
9.1.2	Electrical safety	34
9.1.3	Guidance and manufacturer's declaration on EMC (IEC 60601-1-2:2001)	35
9.1.4	Electrical supply	38
9.1.5	Type label	38
9.1.6	Validated tubing systems	38
9.1.7	Symbols	39
9.2	Storage / transport	39
9.3	Replacement of the main fuses	40
9.4	Environmental compatibility and recycling	42
9.5	Warranty	42

10 Certificates

10.1	EC declaration of conformity	43
10.2	UL classification	44

1 Important information

1.1 Organization of the Operating Instructions

Editorial information	The current version of these Operating Instructions 5/06.09 is valid for machine software V 2.4X and higher.
Changes	Changes to the Operating Instructions will take effect in new editions. In general, this manual is subject to change without notice.

1.2 How to use the Operating Instructions

These Operating Instructions are part of the accompanying documents and thus an integral part of the CompoDock. They contain any information necessary for the use of the desktop sealing device.

The Operating Instructions must be carefully studied before attempting to operate the CompoDock.

Explanation of the Note and Caution symbols used:

Note

Informs the operator that if the steps are not followed as described, a specific function will be executed incorrectly or will not be executed at all, or will not produce the desired effect.

Caution

Advises the operator against certain procedures or actions that could cause damage to the equipment or may have adverse effects on individuals.

Chapter "Accessories (Options)" describes how to operate the options and contains ordering information.

1.3 Brief description

Sealing device	The CompoDock is a high-frequency connection device that permits the sterile connection of two sterile, standard medical PVC tubing sets, such as blood bags, leukocyte depletion filters or apheresis sets.
Tube specifications	Outer diameter: 3.9 mm – 4.6 mm Internal diameter: 2.8 mm – 3.1 mm Wall thickness: 0.5 mm – 0.9 mm All medical PVC tubes complying with the DIN / ISO standard 3826 and having the above dimensions can be used. Sterile sealing of other tube brands or types has been validated by Fresenius Kabi on sample tubes from Baxter, Terumo, Macopharma, and Pall. Since Fresenius Kabi has not concluded any quality agreement with the companies mentioned, we cannot make any general statement on the application of the tubes we use with regard to the systems produced by these companies. For that reason, every operator is himself responsible for validating the sealing quality with the bag systems used in his blood bank, on startup of the CompoDock. It is permitted to test the docks by means of the usual pressure and/or tensile strength tests used in the blood banks. Alternatively, it is also possible to send the docks to Fresenius for a check of their mechanical integrity. As the validations carried out by NPBI¹ and Fresenius² have shown, a mechanically integer dock justifies the assumption that the docking process has been completed sterile as specified above and that, thus, the closed system has been maintained.
Tube combinations	The following tube combinations can be docked: – wet – wet – wet – dry – dry – dry Tubes may only be docked when they have reached room temperature (18–22 °C).
Monitoring	The display provides the following information: – Processing state – Instructions – Error messages – Test programs for service technicians
Function	Two tubes are connected to one another in two steps: – High-frequency sealing generates defined closed tube ends. – The tube ends are heated and connected. 1.Product development 2.Product redevelopment

Method	This sterile connection of blood lines is called "sterile docking", thus permitting sealing under "total containment". Since it radiates heat, thus not having any direct contact with the PVC material of the tubes, the heating element in the CompoDock can be used for an indefinite number of sterile connections.
Operation	Two keys permit easy operation. Operation of the sealing device also requires that all instructions, safety notes and regulations specified in these Operating Instructions are complied with. Any national standards and guidelines must be observed.
Bag trays	The CompoDock provides the possibility to mount bag trays to the tube clamps or folding side tables to the housing. The advantage of the bag trays is that the tube ends must only be about 8 cm long for docking because the bag trays are moved with the tube clamps.

1.4 Intended use

1.4.1 Fields of application

The CompoDock is used for the production of special cellular blood components, e.g. in blood banks where tube connections are needed that provide protection against microbial contamination.

The CompoDock is a high-frequency connection device that permits sterile connection of medical PVC tubes to one another.

Here, sterile connection is achieved when germs from the environment are prevented from being carried into the closed systems. The closed system is maintained after the docking process.

1.5 Standard accessories

CompoDock Counter 100

1.6 Classification

The CompoDock reflects the latest state of technology and complies with the requirements of IEC 61010-1:90 + A1:92 + A2:95 and with the EMC Guideline 89/336/EEC and EN 60601-1-2: 2001.

The CompoDock is protected against imitation by US patents 4.737.214 and EP 0.208.004.

1.7 Maintenance and cleaning

We generally recommend to conclude a maintenance agreement to prevent expensive repair work and to ensure an optimum service life of the device.

Initiate the maintenance measures once a year or when prompted to do so by a service message on the display.

This service message is displayed after 10,000 sealing processes.

Please contact the nearest authorized Fresenius Kabi Technical Service. Service should be performed within a month's time or within a number of 2000 sealing processes.

Cleaning, above all of the tube clamps, is essentially required for maintaining full operability of the device and must, therefore, be carried out regularly and carefully.

1.8 Repair

Assembly, extensions, adjustments, modifications or repairs may only be carried out by the manufacturer or persons authorized by him.

1.9 Safety precautions

Caution

The CompoDock is laboratory equipment. It is not permitted to dock tubes that are in connection with the patient.

Caution

It is always possible that the processed blood is infectious. It must therefore always be treated as being potentially infected.

Caution

If a contamination is detected during the sealing process, the complete system must be discarded even if the dock has not yet been opened.

It cannot be ensured that the inner connection site of two tubes is completely tight.

Caution

Before using the CompoDock, the operator has to make sure of the operability and the proper state of the device and has to observe any other enclosed information and safety-relevant servicing instructions.

Caution

The CompoDock may only be operated and used by persons disposing of the appropriate training, knowledge and experience required.

Caution

When used as intended, the CompoDock may only be used together with CE-certified tubing systems which have been checked for proper functioning and approved by the operator.

Caution

Check each dock for mechanical integrity. A dock of low physical quality may not be considered to be ready-to-use. This is not a typical CompoDock phenomenon, but a general phenomenon of sterile sealing, i.e. sterility of each individual blood component cannot be ensured even if produced properly.

Caution

Improper dock!

- While being docked, the tubes must not be exposed to external forces.
(E.g. do not pull on the tubes.)

This might result in incorrect positioning of the tubes and, thus, in an improper and leaky dock.

Caution

Improper dock!

- Contaminated sealing electrodes may affect the sealing process adversely. This may result in high-frequency arc formation and improper docks.
- Contaminated tube guides result in incorrect positioning of the tubings and, thus, in improper docks.

Contaminated sealing electrodes and tube guides must be cleaned immediately.

Caution

Risk of injury due to high frequencies.

Risk of burning!

- While sealing, keep your fingers away from the sealing electrodes.
 - Do not reach into the unit while carrying out the sealing process.
-

Caution

Risk of injury!

There is the danger of contusing fingers between the heating element cover and the bag trays.

- Do not reach into the unit while carrying out the sealing process.
-

Caution

Device failure caused by fluid entering into the device.

Damage to electronic components.

Contamination of the mechanics.

- Never allow any fluid to enter into the device.
-

Caution

Electric shock!

- If fluid or blood has entered in the interior of the device, immediately disconnect the power plug. Have the device checked by an authorized service technician.

Never use damaged or defective devices.

Use only devices which permit safe operation.

Use devices only after it has been ensured that they are used for their intended application.

**Note**

Repeatedly switching the CompoDock on and off may shorten the life of the heating element substantially. It is recommended to switch on the CompoDock only once a day.

**Note**

Place the bags or other heavy parts of the device that are connected to the tubes onto the bag trays or the side tables (accessories).

**Note**

Connect or disconnect the CompoDock Counter only when the device is turned off.

**Note**

Any wet or cold tubes that are carrying condensate must not be docked immediately. The wet or cold surface causes formation of high-frequency arcs and improper docks. Tubings should only be docked when they have reached room temperature (18–22 °C).

**Note**

To achieve a perfect sealing quality, it is essential for both tubes to be inserted properly in the tube guides and the central block.

Carefully check the tubes for correct positioning before continuing the sealing process.

There may be improper docks.

**Note**

The maximum load on the bag trays must not exceed 2 kg.

**Note**

When placing bag systems onto the bag trays, ensure that the bags do not project beyond the stop. Otherwise, the tubing would impermissibly become subject to tensile forces.

1.10 Manufacturer's address

Please address any inquiries to:

Germany

Fresenius Kabi AG

D-61346 Bad Homburg


Phone: +49 (0) 6172 / 608-0

Local service:


2 Design

2.1 Front view

- 1 Display
- 2 **Start** key
- 3 **Stop** key
- 4 Heater cover
- 5 Stop of bag tray
- 6 Bag tray
- 7 Right clamping lever
- 8 Right tube clamp
- 9 Positioning cover
- 10 Central block
- 11 Left tube clamp
- 12 Left clamping lever
- 13 Power LED
- 14 CompoDock Counter compartment

2.2 Rear view / ports


- 1 Power input module
- 2 Main fuses (hidden behind the flap)
- 3 ON/OFF main power switch
- 4 Connection to power supply (for device plug)
- 5 Ventilation slots
- 6 Scanner connection
- 7 Network connection

3 Sealing principle

The steps below describe the CompoDock sealing process.


Positioning and securing

Carefully insert the tubes to be connected in the tube guides of the tube clamps and in the central block. Then firmly press down the positioning cover to support positioning. When pressed down, the positioning cover actuates a switch, so that the sealing process can be started.


Clamping

When the clamping levers are closed, two clamping jaw pairs squeeze off the tubes. The tube ends remain closed until completion of the sealing process.


Sealing off

Depending on the service setting, the sealing electrodes are closed either automatically or after the **Start** key has been pressed, thus producing a defined dock end at each tube. In the same step, the remaining parts of the tubes are provided with a larger dock. This ensures the operator's safety while the sterile sealing process is in progress.


Disconnecting

While the sealing process is in progress, the tube clamps move outwards, thus disconnecting the docked tube ends from the remaining parts in the central block.


Positioning

The two tube clamps move to the heating position.


**Heating**

A heating element is positioned between the two tube ends. From a very short distance, this heating element radiates heat of a temperature of 460 °C, thus heating both tube ends.

Establishing the tube connection


Once the tube ends are heated sufficiently, the heating element returns to its initial position. The tube clamps move towards one another and connect the tubes. The connection is completed after a preset cooling time.

Removing the tube

After having opened the clamping levers, remove the tube from the tube clamps and open the internal connection site of the tubes.

The sealed tubes are separated by a thin skin in the interior of the sealing site. This thin skin must be opened immediately after completion of the sealing process.

4 Operation

4.1 CompoDock Counter

Operation of the CompoDock requires the CompoDock Counter. The number of sealing processes available (e.g. 2500 possible docks) is programmed in the CompoDock Counter.

The number of docks is decremented by one for each completed sealing process. The number of remaining sealing processes is displayed in the standby mode.

Once the available number of docks is used up, a new Counter is necessary.

Note

Connect or disconnect the CompoDock Counter only when the device is turned off.

Caution

Loss of available docks in the CompoDock Counter.

The Counter may become unusable.

- Never attempt to prevent the Counter from decrementing.
-

Note

In the CompoDock, only Counter modules may be used that have been procured from Fresenius Kabi. Any transfer of CompoDock Counters to third parties is not permitted. Any intervention in the CompoDock Counters or in the electronics of the CompoDock for the purpose of manipulating the Counter reading is not permitted and will be taken to a criminal and civil court and may result in heavy sentences and claims for damages.

4.2 Establishing a sterile connection

Caution

The minimum length of the tube ends to be docked must be 8 cm. We recommend a tube length of 10-12 cm, as this will make insertion of the tube easier and more user-friendly.

It is not possible to reliably connect shorter tubes to one another, because the tube would then be subject to tensile forces during the docking process. Such tensile forces may result in an improper dock.

Caution

Never use tube stripping tongs to strip across a dock.



Considerable mechanical damage to the dock may occur, resulting in leakage and thus loss of sterility.

Note

The **Stop** key is provided as emergency stop button. If a problem is detected during the sealing process, press the **Stop** key. The process will then be stopped. Once the problem has been corrected, continue the process by pressing the **Start** key.

- Insert the CompoDock Counter.
- Connect the device to the power supply and turn it on with the main power switch.


- Open the clamping jaws with the clamping levers.

```
* COMPODOCK V2.422 *
Fresenius Kabi
System check|/|
```

The green LED on the information display lights, showing that the device is ready for operation.
The system check is performed.


Searching Dockmaster

The device attempts to establish a connection to DockMaster. If no DockMaster software is found, the device will switch to its normal operating mode.


- Place the bag system onto the bag trays and ensure it does not project at the stops.
With projecting bags the length of tubing available for sealing is reduced which will result in improper docks.
- Insert tubes and
- slightly pull on the tube when inserting it and press it into the guide.

* xxxx docks left
Insert tubes and
press tube cover

- Check the tube for proper positioning by carefully moving it to and fro in the guide.
Correct insertion of the tube in the guide is essential for proper positioning of the tubes, thus having a direct influence on the dock quality. If the tubes fail to be inserted correctly in the guide, this may more easily lead to improper docks.

- Press the positioning cover fully and firmly down and release it again. The positioning cover remains in the closed position. An audible signal will be emitted.
- Close the clamping levers within 10 sec. after the beep.
The tubes are clamped and fixed with two clamping jaws.
- The sealing process starts automatically after the tube clamps have been closed. The positioning cover opens.
If the Autostart function is deactivated in the Service menu, the sealing process only starts after the **Start** key has been pressed.

* xxxx docks left
close clamps
to start processing

Process start check

Note

Do not open the heater cover.

This would disturb the sealing process and result in improper docks.

Processing ...
Sealing tube ends
Please wait...

The sealing process is in progress.

The display will show the following messages:

1. **Sealing tube ends**
2. **Tubes to heat pos.**
3. **Heating tube ends**
4. **Docking**
5. **Cooling down dock**

Dock completed
Open clamps and
Take out the tubes

After completion of the sealing process, an audible signal is emitted.

- Open both clamping levers.

Dock completed
Take out the tubes
and close cover


- Take the tube in your hands at its right and left ends and pull it up and out of the blocks.

● Establishing a passage

Note

Open the internal connection site of the tubes immediately. Otherwise the connection site cools down and fails to open properly.

- Take the tube very closely at the right and left ends of the dock between thumb and forefinger of both of your hands and press it together.

The passage is established if fluid is transported.

Note

The sealing seams at the remaining tubes are not as strong.

- Prevent exerting too much force onto the remaining tubes.
- If necessary, produce a strong sealing seam with the help of a commercial tube sealing device.

Dock completed
Press **START**

- Remove both remaining tubes from the central block;
- Press the **Start** key.

Dock completed

Returning motors

Both tube clamps return to the standby position. A new sealing process can be started.

● Checking the sealing seam

Caution

Check each dock for mechanical integrity.


Sterility of each individual blood component cannot be ensured even if produced properly.

Companies and institutions must dispose of a properly functioning pharmaceutical quality assurance system.


- in Germany e.g. PharmBetrV § 1a
 - in other countries, observe the local regulations!
 - A sterility test and, thus, the validation of the quality of blood preparation is subject to the operator's scope of duties. Random checks for sterility and hemolysis must be carried out on the basis of valid regulations. Fresenius principally recommend a check of all docks.
-

Quality controls of docks that are currently produced must be carried out at regular intervals.

Visual check

Proper connection:
Symmetrical dock

Improper connection:
Asymmetrical dock, possibly caused by incorrect insertion of the tubes

4.3 Turning the device off

Note

The CompoDock may only be turned off if the clamping jaws are in the standby position (initial position).

If turned off earlier, the number of docks left on the Compodock Counter will be decremented by one each. If turned off earlier several times, the complete number of available docks on the Compodock Counter will be lost.

Turning power off

The tube clamps are in the standby position:

- Use the main power switch to turn the device off.

The tube clamps are not in the standby position:

- Press the Start key.

The tube clamps return to the standby position

- Use the main power switch to turn the device off.

5 Accessories (options)

Caution

Use only accessories / options recommended by the manufacturer.

Assembly, extensions, adjustments, modifications or repairs may only be carried out by the manufacturer or persons authorized by him.

5.1 CompoDock Counter

CompoDock Counter 500	Part number: 902 870 1 Available docks: 500 sealing processes
CompoDock Counter 1000	Part number: 902 869 1 Available docks: 1000 sealing processes
CompoDock Counter 2500	Part number: 902 868 1 Available docks: 2500 sealing processes

5.2 Side tables

Side table, right	Part number: M61 110 1
Side table, left	Part number: M61 117 1

5.3 Bag trays (included in the scope of delivery)

Bag tray, right	Part number: M64 774 1
Bag tray, left	Part number: M64 773 1

5.4 Accessories required for use with DockMaster Net

Software DockMasterNET, CD + User Manual	Part number: 902 892 1
Scanner for CompoDock	Part number: 902 893 1

Scanner holder for M67 893 1	Part number: M60 814 1
D-SUB null-modem cable, 9-way, socket/socket, 3m	Part number: M67 235 1
D-SUB extension, 1:1, 9- way, plug/socket, 5m	Part number: M67 236 1
Interface converter, USB to 4xRS-232	Part number: M67 237 1
Interface converter, USB to 2xRS-232	Part number: M67 238 1
USB cable 2.0, FullSpeed A plug/B plug, 5m	Part number: M67 239 1
USB 2.0 active extension, A plug/A socket, 5m	Part number: M67 240 1

6 Cleaning / disinfection

Caution

For cleaning and disinfecting the surfaces, turn off the device and disconnect it from power supply.

6.1 Surface cleaning

Caution

When cleaning the CompoDock, ensure that fluid is prevented from entering inside.

Prevent spraying inside the device or onto its surfaces.

Clean the surfaces of the device with a moist cloth and a neutral cleaning agent.

If fluid (blood, blood components or cleaning agent) has entered into the device, turn off the device and notify the service technician.

Caution

Risk of burning

Never clean the heating element.

Caution

It is always possible that the processed blood is infectious. It must therefore always be treated as being potentially infected.

Caution

Only use mild and neutral cleaning agents.

Do not use any abrasive and corroding agents nor any products that solve plastic or lubricants.

6.2 Surface disinfection

Note

Use only alcohol-based disinfectants, e.g. Freka®-NOL from Fresenius.

Wipe disinfection

For disinfection and cleaning, wipe off the CompoDock with a soft cloth dampened with disinfectant.

6.3 Cleaning the electrodes

Daily checks

The sealing electrodes must be checked every day. If contaminated, the electrodes must be cleaned, because contaminations may adversely affect the sealing process.

To facilitate cleaning of the sealing electrodes, simultaneously press the **Start** and **Stop** keys. The tube clamps move to the 'cleaning position'. Press the **Start** key again to open and close the electrodes.

Use a cotton applicator dampened with disinfectant (Freka[®]-NOL) for cleaning. Allow the electrodes to dry after cleaning.

Press the **Stop** key to return the blocks to their original standby position.

6.4 Cleaning the tube guides

Daily checks

Check the tube holders on either side every day.

If contaminated, the tube holders must be cleaned, because contaminations support incorrect positioning of the tubes. This may result in reduced dock quality and/or to repeated improper docks.

7 Errors and alarms

7.1 General notes

If the CompoDock detects any irregularity in the sealing process, an alarm message will be displayed. Usually, this alarm message consists of two parts:

- **Cause of the error**

Shows where the error has occurred

- **Result of the error**

There may be three possible results of an error:

- The error does not have any effect on the sealing result.
The sealing process is continued automatically or can be restarted by pressing the **Start** key.
- The error might have an effect on the final result of the sealing process.
The operator must press the **Stop** key to stop the sealing process.
- This error involves a technical problem that cannot be eliminated by the operator.
Turn off the device with the main power switch.

If further problems occur, please contact the nearest authorized Fresenius Kabi Technical Service.

Display message	Error cause	Possible solution
Clamp opened	The tube clamps are open.	Close the tube clamps.
Left clamp not open	The left clamping lever is closed.	Open the left clamping lever.
Left clamp not closed	The left clamping lever is open.	Close the left clamping lever.
Right clamp not open	The right clamping lever is closed.	Open the right clamping lever.
Right clamp not closed	The right clamping lever is open.	Close the right clamping lever.
Both clamps not open	Both clamping levers are closed.	Open both clamping levers.
Both clamps not closed	Both clamping levers are open.	Close both clamping levers.
Heater cover not open	The heater cover is closed.	Open the heater cover.
Heater cover not closed	The heater cover is open.	Close the heater cover.

Display message	Error cause	Possible solution
Counter is damaged or not placed	The CompoDock Counter is defective or inserted improperly.	Turn the device off. Insert the Counter properly. If necessary, replace the Counter by a new one.
Counter data error	The CompoDock Counter is defective.	Turn the device off. Replace the Counter by a new one.
No docks left in counter	The available number of docks is used up.	Insert a new Counter.
Tube(s) not placed	Tubes are missing or are inserted improperly.	Insert the tubes properly.
Dock cooling time interrupted	The cooling time is interrupted. The tubes were removed too early.	Remove the tube connection only after completion of the sealing process.
Seal electrodes contact	Sealing electrode contact problems	Clean the sealing electrodes.
Release tube cover	The tube cover is pressed down.	Open the tube cover.
Open clamps Contamination risk	The tube clamps were opened while the sealing process was in progress.	Do not interfere with the docking process currently in progress.

Display message	Error cause	Possible solution
Maximum time opening seal electrodes	Time error The sealing electrodes open too slowly.	Correct any possible mechanical clamping. Notify the service technician.
Minimum time opening seal electrodes	Time error The sealing electrodes open too fast.	
Maximum time closing seal electrodes	Time error The sealing electrodes close too slowly.	
Minimum time closing seal electrodes	Time error The sealing electrodes close too fast.	
Maximum time moving to outer position	Movement in X1 direction too slow	
Minimum time moving to outer position	Movement to X1 position too fast	
Maximum time moving to heating position	Movement to X2 heating position Maximum running time exceeded	
Minimum time moving to heating position	Movement to X2 heating position Minimum running time fallen below	
Maximum time moving to docking position	Movement to X3 docking position Maximum running time exceeded	
Maximum time moving forward	Movement to Y0 position Maximum running time exceeded	
Minimum time moving forward	Movement to Y0 position Minimum running time fallen below	
Maximum time moving backward	Movement to Y1 position Maximum running time exceeded	
Minimum time moving backward	Movement to Y1 position Minimum running time fallen below	
Cannot reach start position	No movement to initial position.	
Maximum seal time	Maximum sealing time exceeded	
Xmotor seal position wrong	Error in X-motor sealing position	
Ymotor seal position wrong	Error in Y-motor sealing position	
Heater position wrong	Error in heating element position	
Cannot lower heater	The temperature of the heating element cannot be reduced.	
Cannot raise heater	The temperature of the heating element cannot be increased.	

Display message	Error cause	Possible solution
Counter failure 127	Counter monitoring	Notify the service technician.
24 volt supply too low	The 24-V voltage is too low.	
24 volt supply too high	The 24-V voltage is too high.	
5 volt supply too low	The 5-V voltage is too low.	
5 volt supply too high	The 5-V voltage is too high.	
Blower fail	The fan is defective / occluded.	
Heater failure	The heating element is defective.	
Seal electrodes position wrong	Error in the position of the sealing electrodes	
Tech. Maintenance required	Maintenance date reached	
XMOTOR movement detected	Position error	
Dock position fail	Error in dock position	
System watchdog fail	Error in watchdog system	
Controller memory fail (RAM)	RAM memory error	
Controller memory fail (EPROM)	ROM/EPROM memory error	
System settings fail	EEPROM system parameter error	

8 Initial start-up

The initial start-up and the installation may only be performed by an authorized Fresenius Kabi service technician.

The technician can still change the sealing parameters to match the tube combinations used.

Caution

When the CompoDock is installed, the air opening on the rear must always stay uncovered, so that the ventilated air can flow out undisturbed.

Note

Repeatedly switching the CompoDock on and off may shorten the life of the heating element substantially. It is recommended to switch on the CompoDock only once a day.

Unpacking

- Take the CompoDock out of its packaging box.
- Keep the packaging box for a possible future transport.

Place of installation

Place the CompoDock on solid base.

Variations in temperature

An appropriate acclimatization time must be ensured in case of great variations in temperature.

Scope of delivery


Prior to the initial start-up, check that all parts have been accounted for.

Power cable 1x
CompoDock Counter 1x
Bag trays 2x
Operating Instructions

General external checks

Check the CompoDock on the outside for possible damage in transit or other damage.

Connections

The rear panel of the CompoDock is provided with a connection for a grounded power cable.

Connect the power cable to the supply connection on the rear of the device and to a receptacle outlet.

CompoDock Counter

Insert the CompoDock Counter into the counter compartment.

Turning the device power on

- Switch the device on using the main power switch on its rear.

The green Power LED on the information display lights, showing that the device is ready for operation.

The system check is performed.

```
* COMPODOCK V2.422 *  
Fresenius Kabi  
System check    | / |
```

The display shows the following messages:

This is followed by an audible alarm to check the audible alarm generator.

```
Searching Dockmaster
```

The device attempts to establish a connection to DockMaster.

If no DockMaster software is found, the device will switch to its normal operating mode.

```
*   xxxx docks left  
  
Insert tubes and  
Close clamps
```

The CompoDock is ready for operation.

9 Description of the device

9.1 Specifications

9.1.1 Dimensions, weight and housing material

Dimensions	Height: 22 cm
	Width: 34 cm
	With bag trays: 70 cm
	With side tables: 75 cm
	Depth: 45 cm
Weight	Approx. 21.5 kg incl. bag tray
Housing material	1.4305 ground stainless steel sheet and USt37 sheet steel, 1.5 mm thick; completely powder-coated

9.1.2 Electrical safety

Complies with the following requirements	UL 61010-1 Second Edition: 2008 IEC 61010-1 Second Edition: 2008 CAN/CSA-C22.2 No. 610101-1 Second Edition EMC Guideline 89/336/EEC, EN 61000 and EN 55011
Type of protection against electric shock	Safety class I
Degree of protection against ingress of liquids	IPX0

9.1.3 Guidance and manufacturer's declaration on EMC (IEC 60601-1-2:2001)

● Electromagnetic emissions

Guidance and manufacturer's declaration – electromagnetic emissions		
The CompoDock is intended for use in the electromagnetic environment specified below. The customer or the user of the CompoDock should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 2	The CompoDock must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
RF emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	The CompoDock is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

● Electromagnetic immunity

Guidance and manufacturer's declaration – electromagnetic immunity			
The CompoDock is intended for use in the electromagnetic environment specified below. The customer or the user of the CompoDock should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC-61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines Not applicable	Mains power quality should be that of a typical commercial and/or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial and/or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial and/or hospital environment. If the operator of the CompoDock demands continued function even if power failures occur, it is recommended to ensure supply of the Compodock from an uninterruptible power supply or from a battery.

Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: U_T is the a.c. mains voltage prior to application of the test level.			
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 V _{rms} 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 V 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the CompoDock, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ 150 kHz to < 80 MHz $d = 1.2 \sqrt{P}$ 80 MHz to < 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range ^b.  Interference may occur in the vicinity of equipment marked with the following symbol.
Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CompoDock is used exceeds the applicable RF compliance level above, the CompoDock should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the CompoDock. b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

● **Recommended separation distances between portable and mobile RF telecommunication devices and the CompoDock**

Recommended separation distances between portable and mobile RF telecommunication devices and the CompoDock			
The CompoDock is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the operator of the CompoDock can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CompoDock as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to < 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to < 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

9.1.4 Electrical supply


● CompoDock

Input voltage	100 - 120 / 200 - 240 Volt AC autoranging 50 - 60 Hz
Fuses	T6.3 A 250 V (2x)
Power consumption	300 watts max.
RF frequency	40.68 MHz. $\pm$ 10 kHz
Nominal RF output	90 watts

9.1.5 Type label

The type label shown is only an example.

For the current data of the device, consult the device type label.

● Operating conditions

Operating temperature range	15 °C to 35 °C
Relative humidity	35% to 85%
Atmospheric pressure	700 hPa to 1060 hPa
Position of use	horizontal

9.1.6 Validated tubing systems

The following medical PVC tubes have been validated by Fresenius Kabi.

Part number	Specification
M60024	3,1 mm x 4.2 mm
M60025	3,0 mm x 4.55 mm

Sterile sealing of other tube brands or types has been validated by Fresenius Kabi on sample tubes from Baxter, Terumo, Macopharma, and Pall. Since Fresenius Kabi has not concluded any quality agreement with the companies mentioned, we cannot make any general statement on the application of the tubes we use with regard to the systems produced by these companies.

9.1.7 Symbols

Caution, consult accompanying documents; General danger

IPX0

Protection against ingress of liquids: IPX0

E-code xxx

Equipment code of the device

The CE mark documents the manufacturer's declaration that the requirements of the EC directives applied by him have been complied with.

Conformity of this device with legal protection requirements (89/336/EEC) is documented in the specifications and by the CE mark.

Device emits non-ionizing radiation

Identification of electric and electronic devices

40,68 MHz

Radio frequency oscillator

Keep hands out of machinery

9.2 Storage / transport

Environmental conditions

Temperature -20 to +45 °C
Humidity 20 – 85 %
Atmospheric pressure 500 – 1060 hPa

Storage

The CompoDock should be stored in a well ventilated room.

Transport outside buildings

The CompoDock must be transported in its original packaging or in the transport packaging which is available as an accessory part.

Transport packaging:

Description	Part number
Packaging box	M61 125 1
Foam insert for bottom	M61 126 1
Foam insert for cover	M61 127 1
Handle for packaging	M60 625 1 (required 2x)
Plastic bag	M62 416 1

Transport inside buildings

The CompoDock may only be held on the circumferential lower edge for carrying (see figure)

Caution

Incorrect transport may damage accessory components (clamps, tubing covers, heater door).

Never lift the CompoDock by the accessory components.

9.3 Replacement of the main fuses

Caution

Danger to life caused by electrical voltage

Touching live parts will cause an electric shock.

- Switch off the device using the main switch.
 - Disconnect the power plug.
-

Caution


The device may be damaged.

When replacing a fuse observe the specifications regarding the fuses. (see **Electrical supply**, page en--38).

Step 1

The two main fuses are accessible from outside and positioned in the power input module behind a protective cover.


- Open the protective cover by means of a screwdriver.

Step 2

- Use a screwdriver to remove the fuse module from the plugged connections.

Step 3

- Pull out the fuse module.

Step 4

- Use a screwdriver to prise the defective fuse from the fuse holder.

**Note**

Always replace both fuses in case of a fuse defect.

**Installation of the new
fuses**

Install the new fuses in reverse order.

9.4 Environmental compatibility and recycling

Only environmentally compatible and recyclable materials are used for the manufacture of the CompoDock. Prior to the disposal of the CompoDock, any potential risk of infection must be eliminated by appropriate disinfection. The respective regulations for the disposal of electronic scrap should be followed for the disposal of electronic boards. Disposal of batteries should always conform to applicable regulations. Further information regarding disposal is available on request.

9.5 Warranty

Warranty claims are, among others, excluded in the following cases: improper use of the device, non-respect of the notes in the Operating Instructions, opening of the device or damage to the test seal within the warranty period.

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We / Wir Fresenius Kabi AG D - 61346 Bad Homburg, Germany
Fresenius Kabi AG Transfusion Technology Division