

# SuctionPro72™ Closed Ventilation Suction System



**Reduce Infection. Reduce Patient Stay. Reduce Costs.**

# SuctionPro72™



## SuctionPro 72™

The Portex® SuctionPro72™ Closed Ventilation Suction System is a single patient use suctioning device for the removal of secretions from the tracheobronchial tree of ventilator dependent adult patients. Intended for 72-hour use.

## Key Features

- 3-day recommended duration of use
- Clear pathway evacuation port
- Lockable thumb valve end cap
- Sterile, single patient use
- Clear T piece for visualisation of the pathway
- Soft but strong catheter sleeve
- MDI Adaptor for integrated inhaler capability
- Patient labels now coloured by day for easy identification
- Trac-Wedge™ device to aid in disconnection of the catheter from the patient's endotracheal or tracheostomy tube
- Swivel connector to reduce torque to patient in some packs



Fig. 1



Fig. 3



Fig. 4



Fig. 2

## Instructions for Use

### To Lavage

**Fig. 1** Hold T-piece in one hand and advance catheter approximately 10cm into the airway. Instill saline solution through the irrigation inlet.

### To Suction:

**Note:** Patients may benefit from pre-oxygenation with 100% oxygen.

**Fig. 2** Make sure the suction control valve lock is in the 'OPEN' position. Advance the catheter to the desired depth whilst holding the patient end steady. If resistance is met, withdraw the catheter 2-3cm before applying suction.

**Fig. 3** Grasp the control valve and apply backwards-sliding pressure on the blue thumb actuator to suction.

**Note:** Maximum suction is achieved by sliding the actuator fully back.

Withdraw the catheter slowly with suction activated in a straight motion to avoid kinking until blue mark is just fully visible in catheter sleeve.

**Fig. 4** Ensure that the catheter tip is out of the breathing path and in-line with the saline port. Begin to clean catheter tip with saline. The saline should be administered through the irrigation inlet whilst vacuum is applied making sure that the tip and area surrounding it is fully flushed with saline. Release control valve actuator and turn the valve lock to the 'CLOSE' position when finished.

## Preparation

- Before attaching the system to the patient turn on the suction, make sure the lock is turned to the 'OPEN' position and check the operation of the control valve by sliding back the actuator. Once in the fully back position release and make sure that the device shuts correctly
- Attach male connector of the SuctionPro72™ Closed Ventilation Suction System device to suction tubing
- Attach the suction tubing to the male connector of the SuctionPro72™ closed ventilation suction system
- Attach ventilator circuit to dual-swivel or Tpiece adaptor
- Attach the dual-swivel or Tpiece adaptor to the tracheal or tracheostomy tube connector

Parameter	Value	Parameter	Value	Parameter	Value
Temp (oral)	37.5	HR (bpm)	78	RR (bpm)	18
Temp (axillary)	37.2	SpO2 (%)	98	PEEP (cm H2O)	5
Temp (rectal)	37.8	MAP (mmHg)	85	FiO2 (%)	21
Temp (nasal)	37.1	SBP (mmHg)	120	Flow (L/min)	30
Temp (esophageal)	37.4	DBP (mmHg)	70	Volume (ml)	500
Temp (skin)	36.8	PPV (%)	3	Compliance (ml/cm H2O)	20
Temp (core)	37.3	AP (mmHg)	60	Resistance (cm H2O/L/min)	1.5
Temp (periphery)	36.5	CI (L/min/m2)	2.2	Power (W)	50
Temp (limb)	36.2	SV (ml)	50	Work of Breathing (J/min)	100
Temp (ear)	37.0	Stroke Volume (ml)	70	Work of Breathing Index (J/min/L)	0.5
Temp (tympanic)	37.4	Cardiac Output (L/min)	5.5	Work of Breathing Index (J/min/L)	0.5

Day Label

# SuctionPro72™

Comprehensive product range available in single and dual lumen configurations, with coloured day labels. MDI adaptor in non dual swivel options. Each Portex® SuctionPro72™ suction system offers a wide range of options to enhance patient care outcomes and accommodate clinical practices. Available in a case of 20 units.

## IDENTIFICATION MATRIX

	Part No	10 FR	12 FR	14 FR	16 FR	300mm Length	570mm Length	Flex Tube	Coudé Tip	Dual Swivel	MDI adaptor
SINGLE LUMEN	Z110-10	•					•				•
	Z110-12		•				•				•
	Z110-14			•			•				•
	Z110-16				•		•				•
	Z115-10	•				•					•
	Z115-12		•			•					•
	Z115-14			•		•					•
	Z115-16				•	•					•
	Z116-14			•		•		•			•
	Z118-14			•			•	•			•
	Z120-10	•					•		•		•
	Z120-12		•				•		•		•
	Z120-14			•			•		•		•
	Z120-16				•		•		•		•
	Z130-14*			•			•				•
	Z130-16*				•		•				•
	Z135-14*			•		•					•
	Z135-16*				•	•					•
	Z150-10	•					•				•
	Z150-12		•				•				•
Z150-14			•			•				•	
Z150-16				•		•				•	
Z155-10	•					•				•	
Z155-12		•				•				•	
Z155-14			•			•				•	
Z155-16				•		•				•	
Z156-14			•			•		•		•	
Z160-14			•				•		•	•	
Z160-16				•			•		•	•	
DOUBLE LUMEN	Z210-12		•				•				•
	Z210-14			•			•				•
	Z210-16				•		•				•
	Z215-12		•			•					•
	Z215-14			•		•					•
	Z216-14			•		•		•			•
	Z250-12		•				•			•	
	Z250-14			•			•			•	
	Z250-16				•		•			•	
	Z255-12		•			•				•	
	Z255-14			•		•				•	
	Z256-14			•		•		•		•	

\*not CE marked

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

**Smiths Medical International Ltd.**  
 Ashford, UK TN25 4BF  
 Phone: +44 (0) 845 8500445  
[www.smiths-medical.com](http://www.smiths-medical.com)

Find your local contact information at: [www.smiths-medical.com/customer-support](http://www.smiths-medical.com/customer-support)

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Portex, SuctionPro 72 and the Smiths Medical and Portex design marks are registered trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2013 Smiths Medical. All rights reserved. RE194257GB-082013

**CE Rx**  
 0473 ONLY

**smiths medical**

MHYTCA-1030

# Respiratory Care Solutions Improving Quality of Life



## Improving quality of life...



**The need to rehabilitate patients effectively after respiratory disorders, is of extreme importance, not only to decrease patient recovery time for improved hospital efficiency, but also for the well-being and quality of life of the patient.**

Introducing a Pulmonary Rehabilitation Programme (PRP) can help patients with lung disease achieve the highest possible level of functioning. PRPs have been shown to improve quality of life, mitigate symptoms, improve exercise tolerance, and lower the number of hospital admissions.

Two key elements of a PRP are Bronchial Hygiene Therapy (BHT) and Lung Expansion Therapies.

Smiths Medical has developed a comprehensive range of respiratory care products designed for both hospital and home use.

These products help to rehabilitate patients with the aim of improving their physical and social performance. By focusing on the rehabilitation of patients and continued lung training, hospitals can potentially see cost savings due to reduced hospital stays and home rehabilitation.

#### References:

1. Wiersgalla Susan, RRT, RCP, North Memorial Medical Center, Robbinsdale, MN. Abstract presented at the 48th International Respiratory Congress for the AARC Annual Convention and Exhibition on October 5th, 2002 in Tampa, Florida.
2. Steen HJ, Redmond AOB, O'Neill D, Beattie F. Acta Paediatr Scand. Evaluation of the PEP mask in cystic fibrosis. 1991; 80:51-56.
3. Tyrell JC, Hiller EJ, Martin J. Face mask physiotherapy in cystic fibrosis. Archives of Dis in Child 1986; 61: 598-611.
4. Mahlmeister MJ, Fink JB, Hoffman GL, Fifer LF. "Positive-expiratory-pressure mask therapy: Theoretical and Practical Considerations and a Review of the Literature", Respiratory Care, 1991;36:1218-1230.
5. Guell R. Breath, Home- Based Rehabilitation 2008,5 pg-37





Figure 1




Figure 2




Figure 3  
Post-Surgical Atelectasis Improvement in CABG Patients.  $p < .001$

## Rehabilitation

For hospitalised patients, the correct rehabilitation can save the hospital both time and money, as well as making the experience more tolerable for the patient.

Many hospitals today, use Intermittent Positive Pressure Breathing (IPPB) to treat and reverse atelectasis as well as being used in re-expanding lung parenchyma. IPPB is a form of assisted ventilation that triggers a positive pressure breath to the patient when the machine senses either effort by the patient, or a negative pressure of an inspiratory breath. Once a certain pressure is reached it then returns to zero. However, IPPB starts with a negative pressure and does not deliver positive expiratory pressure during exhalation. (See figure 1)

EzPAP® is a simple and effective "In-Hospital" method for delivering positive airway pressure throughout the

breathing cycle. Pressure does decrease during inhalation, but it always remains positive, helping to open airways and re-inflate collapsed alveoli. (See figure 2). EzPAP® also provides the additional benefit of Positive Expiratory Pressure therapy during exhalation offering further rehabilitation to the patient. (See figures 4,5 &6)

The combination of therapies offers a safe and effective alternative to IPPB for hospitalised patients. EzPAP® is easy to use and inexpensive in comparison to IPPB. Wiersgella has also demonstrated that EzPAP®, when used post operatively on coronary artery bypass graft patients, shows measurable improvements in atelectasis levels. <sup>1</sup> (figure 3)

Positive Expiratory Pressure therapy can also be achieved using the acapella® duet from Smiths Medical. The acapella® duet is a vibratory PEP device that

## from Hospital...

has a built-in port for aerosolised medication via a small volume nebuliser. Medication such as bronchodilator or mucolytic treatments can be delivered simultaneously with PEP therapy, cutting treatment time to one session, saving time and effort for both clinicians and patients.

Products that encourage PEP therapy are usually well tolerated by most patients and due to the option of self-administering, are an ideal choice for patients who are out of ICU. In contrast to traditional chest physiotherapy (CPT) options, PEP therapy offers effective secretion removal at a low cost. By being independent of daily assistance from clinicians, the compliance level is enhanced and this in turn can also reduce related health-care costs.<sup>2,3</sup>





Figure 4  
Collateral Ventilation




Figure 5 Pursed-lips breathing (or use of a fixed orifice resistor such as a PEP device) creates back pressure that splints the airway open during exhalation.<sup>17</sup>




Figure 6 Elevated intrathoracic pressure can compress unstable airways during exhalation.<sup>17</sup>

acapella® duet and acapella® choice, are inexpensive, easy to use standalone options, providing vibratory PEP therapy to remove secretions. They can be used in any position and offer high and low flow rates for the patient. These devices, with the option of a mouthpiece or mask, are a fully versatile product for both patient and hospital.

TheraPEP® is an alternative product offering PEP therapy. This therapy can also be self-administered in half the time of CPT.<sup>4</sup> TheraPEP® can accommodate virtually any lung capacity and allows inhalation and exhalation without removal from the mouth. With a 22mm ID connector to allow small volume nebulisers or MDI spacers and the option of mouthpiece or mask, this is another great choice for secretion clearance and atelectasis reversal in the hospital setting or at home.

## to Home

For post-surgical patients, medical devices that help improve their lungs are an important function. Incentive spirometers encourage patients to take slow and deep breaths to expand the lungs. Smiths Medical offers both the Coach® 2 incentive spirometer and the CliniFLO® low-flow incentive spirometer. Coach® 2 combines a one-way valve, highly visible piston and easy to understand graphics indicating correct inspiratory flow rate to help patients perform and monitor their own post-surgical breathing exercises without the need for direct supervision. CliniFLO® is ideal for generic, paediatric or weakened patients due to flow settings as low as 100ml/sec.

Once a patient has been trained to self-administer their chosen therapy in a clinical setting, they are able to continue this therapy at home. acapella® duet, acapella® choice, TheraPEP®, Coach® 2 and CliniFLO® all offer the

versatility to be used in a clinical and/or home setting to provide continued therapy. Each is lightweight, easy to use and transportable. Studies have indicated that home-based rehabilitation programmes not only provide similar benefits to hospital rehabilitation programmes but also can reduce the use of medication and the number of hospitalisations.<sup>5</sup>

Continued lung exercises for patients, whether it is post-surgical in the hospital or at home are key to facilitating patient recovery and therefore an improved quality of life. The respiratory care range from Smiths Medical tailors to each patient's need whilst enabling the clinician to save time and money.




## EzPAP®



### Features and benefits:

- 22mm OD patient end to accommodate mouthpiece or 3 mask sizes.
- Scalloped ambient air inlet with gaps to resist occlusion.
- Pressure port with detachable cap, allows connection to a gauge for easy monitoring.
- Complete procedural kit.
- Disposable manometer.
- Easy to use, no extensive training.
- Can be used in conjunction with aerosol medication (e.g. nebuliser) via 22mm connection.

## acapella® Vibratory PEP Therapy System



### Features and benefits:

- Convenient built-in nebuliser port, standard sized to fit most medication nebulisers (acapella® duet only).
- $\geq 10\text{L}/\text{min}$  expiratory flow requirement.
- Streamlined body design offers easy grip (acapella® duet only).
- Adjustable frequency and flow resistance settings.
- Clear colouring aids in visual recognition of cleanliness (acapella® duet only).
- Tethered cap to reduce risk of contamination when in resting position (acapella® duet only).
- Easily disassembled for heat disinfection by boiling, autoclaving and dishwasher (top shelf only)
- Functional in any position – Trendelenburg, standing or sitting.
- Distal 22mm OD fitting allows nebuliser connection via tee adaptor (acapella® choice only).
- One-way inspiratory valve allows inhalation without removal from the mouth.
- Proximal 22mm OD connection allows use with mouthpiece or mask.

## TheraPEP®



### Features and benefits:

- Six Fixed Orifice Options.
- Built-in durable pressure indicator.
- 22mm OD patient end.
- Inspiratory valve.
- Can accommodate virtually any patient's lung capacity.
- Resists breakage, unlike fragile, costly manometers.
- Provides immediate, visual 360° feedback of prescribed pressure.
- May be used with a mask or mouthpiece, or Nebulizer.
- Allows inhalation and exhalation without removing from mouth

## Coach® 2 and CliniFLO®



### Features and benefits:

#### Coach® 2

- Ensures patients inhale, rather than exhale into the unit.
- Easily adjustable for each patient's use.
- Can be seen by patients emerging from the effects of anaesthesia.
- Easy to train.
- Stays with the patient for maximum compliance.
- Accommodates patients requiring supplemental oxygen.
- Saves space.
- Colourful deep-sea characters, games, puzzles and stickers.

#### CliniFLO®

- Can accommodate virtually any patient.
- Reduces the possibility of contamination.
- Provides immediate patient feedback.
- Easy to train.
- Reduces the chance that the setting will be changed inadvertently.

## EzPAP®

### ORDERING INFORMATION

Product Code	PRODUCT DESCRIPTION	UNITS / CASE
23-0747	EzPAP® System with Mouthpiece	10
23-0757	EzPAP® System with Disposable Manometer and Mouthpiece	10
23-1747	EzPAP® System with Paediatric Mask	1
23-2747	EzPAP® System with Medium Mask	1
23-3747	EzPAP® System with Large Mask	1
23-6000*	EzPAP® Kit: One Pressure Gauge (With Gauge Protector), Three EzPAP® Units (with Mouthpiece), Three Gauge Guards and Ten 22 mm ID Adaptors in a Durable Plastic Box.	1

Each system includes EzPAP® one 7 ft oxygen tube, one pressure port cap, and one of the above

\* not CE marked

## acapella® Vibratory PEP Therapy System

### ORDERING INFORMATION

PART No.	PRODUCT DESCRIPTION	UNITS/CASE
27-9000	acapella® duet kit, includes acapella® duet, mouthpiece, Portex® SVN, oxygen tubing, collapsible flex tubing	10
27-9001	acapella® duet and mouthpiece only	10
007760	Portex® updraft medication nebulizer	50
27-7000	acapella® choice with mouthpiece	10
21-1530	acapella® DH with mouthpiece	10
21-3530	acapella® DH with Paediatric mask	1
21-5530	acapella® DH with medium mask	1
21-7530	acapella® DH with large mask	1
21-1015	acapella® DM with mouthpiece	10
21-3015	acapella® DM with Paediatric mask	1
21-5015	acapella® DM with medium mask	1
21-7015	acapella® DM with large mask	1

### ACCESSORIES

Product Code	PRODUCT DESCRIPTION	UNITS/CASE
27-0050	Replacement mouthpiece; fits all acapella® family designs	50

Note: Product is for single patient use only.

For further information please visit:  
[www.smiths-medical.com/respiratorycare](http://www.smiths-medical.com/respiratorycare)

## TheraPEP®

### ORDERING INFORMATION

Part No.	PRODUCT DESCRIPTION	UNITS/CASE
20-1112	TheraPEP® System with Mouthpiece	10
20-3112	TheraPEP® System with Paediatric Mask	1
20-5112	TheraPEP® System with Small Mask	1
20-7112	TheraPEP® System with Large Mask	1

### ACCESSORIES

Part No.	PRODUCT DESCRIPTION	UNITS/CASE
20-3115	Paediatric Mask	1
20-5115	Small Adult Mask	1
20-7115	Large Adult Mask	1
20-0005	Connector, Straight, 22 mm I. D.	10
20-0010	TheraPEP® Pressure Port	10
20-0022	TheraPEP® Pressure Port, Tubing, Indicator	10
20-0050	TheraPEP® Mouthpiece 22 mm I. D.	50
20-0120	TheraPEP® Pressure Port and Resistor	10
20-1110	TheraPEP® Pressure Port, Resistor and Mouthpiece	10

All TheraPEP® Systems Include: Pressure Port, Resistor, 22 mm ID Straight Connector, Tubing and Pressure Indicator.

## Coach® 2 Incentive Spirometers

### ORDERING INFORMATION

PART No.	PRODUCT DESCRIPTION	Volume
22-4000	Coach®2 One way valve	4000ml
22-4001	Coach®2	4000ml
22-2500	Coach®2 One way valve	2500ml
22-2501	Coach®2	2500ml
22-2000	Coach®2 Kids One way valve	2000ml

## CliniFLO® Low-Flow Incentive Spirometers

### ORDERING INFORMATION

PART No.	PRODUCT DESCRIPTION
22-1200	CliniFlo®

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

**Smiths Medical International Ltd.**  
 Ashford, UK TN25 4BF  
 Phone: +44 (0) 845 8500445  
[www.smiths-medical.com](http://www.smiths-medical.com)

Find your local contact information at: [www.smiths-medical.com/customer-support](http://www.smiths-medical.com/customer-support)

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Portex, EzPAP, Coach, CliniFLO, acapella, TheraPEP, and the Smiths Medical and Portex design marks are registered trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2013 Smiths Medical. All rights reserved. RE194286GB-082013

**CE**  
**Rx**  
**0473 ONLY**

**smiths medical**

MROCCA-0012

# Certificate of Registration



This is to certify that the quality management system of

## SMITHS MEDICAL CZECH REPUBLIC a.s.

Olomoucká 306, 753 01 Hranice, Czech Republic

has been assessed and registered by AMTAC Certification Services Limited as conforming to the requirements of:

### EN ISO 13485:2012

The quality management system is applicable to:

The design, assembly, manufacture, packaging and supply of:

Obstetrics and Gynaecology Devices and Accessories,  
Interventional Imaging Devices and Accessories,  
Oxygen & Humidity Management Devices,  
Pain Management Devices and Accessories,  
Invasive Patient Pressure Monitoring Devices and Accessories,  
Tracheostomy Devices,  
Disposable Infusion Devices,  
Suction Catheters Devices,  
Intubation Systems Devices.

<b>Certificate Number:</b>	<b>1201-04 B</b>
<b>Initial Certification Date:</b>	<b>10 January 2014</b>
<b>Certificate Effective Date:</b>	<b>22 May 2017</b>
<b>Certificate Expiry Date:</b>	<b>28 February 2019</b>



*Brian Johnson*

*AMTAC Certification Services Limited, Milton Keynes, UK*

*This certificate is the property of AMTAC Certification Services Ltd*



061

In the issuance of this certificate, Intertek assumes no liability to any party other than to the Client, and then only in accordance with the agreed upon Certification Agreement. This certificate's validity is subject to the organization maintaining their system in accordance with Intertek's requirements for systems certification. Validity may be confirmed via email at [certificate.validation@intertek.com](mailto:certificate.validation@intertek.com) or by scanning the code to the right with a smartphone. AMTAC Certification Services Limited is owned by AMTAC Certification Services Holdings Limited, which is a wholly owned subsidiary of Intertek UK Holdings Limited. AMTAC Certification Services Limited is an accredited body registered under UKAS with the identification number of 061. In the issuance of this certificate, AMTAC assumes no liability to any party other than to the Client and then only in accordance with the agreed Terms and Conditions.



The certificate remains the property of Intertek, to whom it must be returned upon request.

# Certificat de Înregistrare

## Intertek

Se certifică prin prezenta că sistemul de management al calității al

### SMITHS MEDICAL CZECH REPUBLIC a.s.

Olomoucká 306, 753 01 Hranice, Republica Cehă

a fost evaluat și înregistrat de AMTAC Certification Services Limited ca fiind conform cerințelor:

### EN ISO 13485:2012

Sistemul de management al calității este aplicabil pentru:

Proiectarea, asamblarea, fabricarea, ambalarea și furnizarea de:

- Dispozitive și Accesorii pentru Obstetrică și Ginecologie,
- Dispozitive și Accesorii Intervenție Imagistică,
- Dispozitive de Management al Oxigenului și Umidității,
- Dispozitive și Accesorii de Management al Durerii,
- Dispozitive și Accesorii Invasive de Monitorizare a Tensiunii Pacientului,
- Dispozitive Traheotomie,
- Dispozitive de Unică Folosință pentru Injecții,
- Dispozitive Catetere Aspirare,
- Sisteme de Dispozitive de Intubare.

**Certificat Număr:** 1201-04 B  
**Data Inițială a Certificării:** 10 Ianuarie, 2014  
**Data Efectivă a Certificatului:** 22 Mai 2017  
**Data Expirării Certificatului:** 28 Februarie 2019



Brian Johnson

AMTAC Certification Services Limited, Milton Keynes, UK

Acest Certificat este proprietatea AMTAC Certification Services Ltd



061

În emiterea acestui certificat, Intertek nu-și asumă nicio răspundere față de altă Parte în afară de Client, și aceasta, numai în conformitate cu Acordul de Certificare agreed. Validitatea acestui certificat depinde de păstrarea de către organizația a sistemului în conformitate cu cerințele Intertek pentru sistemele de certificare. Validitatea poate fi confirmată prin email la [certification@intertek.com](mailto:certification@intertek.com) sau prin scanarea codului din dreapta cu un. AMTAC Certification Services Limited este deținută de AMTAC Certification Services Holdings Limited, care este o succursală deținută integral de Intertek UK Holdings Limited. AMTAC Certification Services Limited este un organism acreditat înregistrat la UKAS cu numărul de identificare 051. În emiterea acestui certificat, AMTAC nu-și asumă nicio responsabilitate față de nicio parte, altă decât Clientul, și aceasta numai în conformitate cu Termenii și Condițiile agreeate. Certificatul rămâne proprietatea Intertek, căreia îi trebuie returnat la cerere. CT-ISO 9001:2008-UKAS-EN-LT-P-04.jan.12



\*\*\*\*\*

Subsemnata **MUSUROIA MIRELA**, traducator autorizat de Ministerul Justitiei, certific exactitatea acestei traduceri cu textul înscrisului original în limba engleza, ce a fost vizat de mine.

Traducător autorizat  
Nr. 2769/2015



# Certificate of Registration



This is to certify that the quality management system of

## SMITHS MEDICAL INTERNATIONAL LIMITED

Boundary Road, Hythe, Kent, CT21 6JL, UK

has been assessed and registered by AMTAC Certification Services Limited as conforming to the requirements of:

### EN ISO 13485:2012

The quality management system is applicable to:

Design of:

Breathing Systems, Drainage Devices, Feeding Devices, Filtration Devices, Infusion Disposables, Intubation Systems, Obstetrics and Gynaecology Devices, Oxygen and Humidity Management Devices, Pressure Monitoring Devices, Respiratory Mechanics Devices, Resuscitation Devices, Suction Catheters, Tracheostomy Tubes, Vascular Access Devices, Cardio Thoracic Catheters and Drapes

Additional Site:

Human Resources and Training, Shipping, Demand Planning, Post Market Surveillance, Market Intelligence, E-Business, International Customer Services, Business Development, Registrations, Finance, Wallace Women's Healthcare

<b>Certificate Number:</b>	<b>053-01 B</b>
<b>Initial Certification Date:</b>	<b>20 October 2005</b>
<b>Certificate Effective Date:</b>	<b>23 July 2016</b>
<b>Certificate Expiry Date:</b>	<b>22 July 2019</b>



061



*Brian Johnson*

*AMTAC Certification Services Limited, Milton Keynes, UK*

*This certificate is the property of AMTAC Certification Services Ltd*

In the issuance of this certificate, Intertek assumes no liability to any party other than to the Client, and then only in accordance with the agreed upon Certification Agreement. This certificate's validity is subject to the organization maintaining their system in accordance with Intertek's requirements for systems certification. Validity may be confirmed via email at [certificate.validation@intertek.com](mailto:certificate.validation@intertek.com) or by scanning the code to the right with a smartphone. AMTAC Certification Services Limited is owned by AMTAC Certification Services Holdings Limited, which is a wholly owned subsidiary of Intertek UK Holdings Limited. AMTAC Certification Services Limited is an accredited body registered under UKAS with the identification number of 061. In the issuance of this certificate, AMTAC assumes no liability to any party other than to the Client and then only in accordance with the agreed Terms and Conditions.



The certificate remains the property of Intertek, to whom it must be returned upon request.

# Certificat de Înregistrare

## Intertek

Prin prezenta se certifică faptul că sistemul de management al calității al

### SMITHS MEDICAL INTERNATIONAL LIMITED

Boundary Road, Hythe, Kent, CT21 6JL, UK

a fost evaluat și înregistrat de AMTAC Certification Services Limited ca fiind conform cerințelor:

### EN ISO 13485:2012

Sistemul de management al calității este aplicabil pentru:

Proiectarea:

Sistemelor de Respirat, Dispozitivelor de Drenaj, Dispozitivelor de Alimentare, Dispozitivelor de Filtrare, Infuzoare de Unică Folosință, Sisteme de Intubare, Dispozitive de Obstetrică și Ginecologie, Dispozitive de Management al Oxigenului și Umidității, Dispozitive de Monitorizare a Presiunii, Dispozitive Mecanice de Respirație, Dispozitive de Resuscitare, Catetere Sucțiune, Tuburi Traheotomie, Dispozitive Acces Vascular, Catetere Cardio Toracice și Comprese Chirurgicale

Poziții Suplimentare:

Resurse Umane și Pregătire Profesională, Transport, Cerere Planificare, Supraveghere Post Piață, Cunoașterea Pieței, E-Business, Servicii Client Internațional, Dezvoltare Afacere, Înregistrări, Finanțe, Wallace Women's Healthcare

Certificat Număr: Inițial Data	053-01 B
Certificării: Dată Efectivă	20 Octombrie 2005
Certificat:	23 Iulie 2016
Data Expirării Certificatului:	22 Iulie 2019

Brian Johnson

AMTAC Certification Services Limited, Milton Keynes, UK  
Acest certificat este proprietatea AMTAC Certification Services Ltd



061

În emiterea acestui certificat, Intertek nu-și asumă nicio răspundere față de altă Parte în afară de Client, și aceasta, numai în conformitate cu Acordul de Certificare agreeat. Validitatea acestui certificat depinde de păstrarea de către organizație a sistemului în conformitate cu cerințele Intertek pentru sistemele de certificare. Validitatea poate fi confirmată prin email la [certificate.validation@intertek.com](mailto:certificate.validation@intertek.com) sau prin scanarea codului din dreapta cu un. AMTAC Certification Services Limited este deținută de AMTAC Certification Services Holdings Limited, care este o sucursală deținută integral de Intertek UK Holdings Limited. AMTAC Certification Services Limited este un organism acreditat înregistrat la UKAS cu numărul de identificare 061. În emiterea acestui certificat, AMTAC nu-și asumă nicio responsabilitate față de nicio parte, alta decât Clientul, și aceasta numai în conformitate cu Termenii și Condițiile agreeate. Certificatul rămâne proprietatea Intertek, careia îi trebuie returnat la cerere.



Intertek Intertek Intertek Intertek Intertek