INSTRUCTIONS FOR USE

PipelineTM Flex
Embolization Device
Shield TechnologyTM

TABLE OF CONTENTS

Pipeline™ Flex Embolization Device with Shield Technology™

English	4
Français	5
Deutsch	7
Italiano	8
Español	10
Svenska	12
Nederlands	13
Português	15
Suomi	16
Dansk	18
Ελληνικά	20
Čeština	21
Magyar	23
Русский	24
Polski	26
Türkçe	28
Norsk	29
Slovenčina	30
Română	32
Български	33
한국어	35
عربي	37
Symbol Glossary	38

Instructions for Use

Pipeline™ Flex Embolization Device with Shield Technology™

CAUTION

- This device should be used only by physicians trained in percutaneous, intravascular techniques and procedures at medical facilities with the appropriate fluoroscopic equipment.
- The Pipeline™ Flex embolization device with Shield Technology™ should be used by physicians who have received appropriate training for this device.
- Carefully inspect the sterile package and the Pipeline™ Flex embolization device with Shield Technology™ prior to use to verify that neither has been damaged during shipment. Do not use kinked or damaged components.
- The Pipeline™ Flex embolization device with Shield Technology™ is not to be used after the expiration date imprinted on the product label.
- Do not reprocess or resterilize. Reprocessing and resterilization increase the risks of patient infection and may compromise device performance.

DESCRIPTION

The Pipeline™ Flex embolization device with Shield Technology™ consists of a permanent implant combined with a guide-wire-based delivery system. The Pipeline™ Flex embolization device with Shield Technology™ implant is a braided, multi-alloy, mesh cylinder woven from platinum/tungsten and cobalt-chromium-nickel alloy wires. A photograph of Pipeline™ Flex embolization device with Shield Technology™ is shown in Figure 1a and the design of the distal delivery system is shown in Figure 1b. The woven wires of the device provide approximately 30% metal coverage of the arterial wall surface area. The implant is designed for placement in a parent vessel across the neck of an intracranial aneurysm (IA). The expanded or un-constrained diameter is 0.25 mm larger than the labeled diameter.

Shield Technology™ is a synthetic polymer surface treatment and is not derived from any animal or human sources.


The tip coil is made of platinum-tungsten alloy, the proximal bumper is platinum-iridium alloy, and the tip, distal and proximal solder joints are composed of tin-silver. The protective sleeves are designed to protect the distal portion of the braid while the Pipeline™ Flex embolization device with Shield Technology™ is advanced through the micro catheter. The proximal bumper and resheathing pad allow the user to push the Pipeline™ Flex embolization device with Shield Technology™ out of the micro catheter when the delivery system is advanced. The resheathing pad also allows the user to resheath the Pipeline™ Flex embolization device with Shield Technology™ back into the micro catheter. The resheathing marker provides the user fluoroscopic visualization for the limit of resheathing the Pipeline™ Flex embolization device with Shield Technology™.

The Pipeline™ Flex embolization device with Shield Technology™ implant is mounted on a 304 stainless steel micro guide-wire approximately 200 cm long and compressed inside an introducer sheath. The Pipeline™ Flex embolization device with Shield Technology™ is designed to be delivered only through a compatible micro catheter of 0.027 inch (0.69 mm) inside diameter at least 135 cm in length.

Figure 1a: The Pipeline™ Flex embolization device with Shield Technology™

Figure 1b: The Pipeline™ Flex embolization device with Shield Technology™

- | | | |
|----------------------|----------------------|------------------------|
| 1. Tip Coil | 5. Fluorosafe Marker | 9. Resheathing Pad |
| 2. Proximal Bumper | 6. Delivery Wire | 10. Resheathing Marker |
| 3. Introducer Sheath | 7. Distal Marker | |
| 4. Braid | 8. PTFE Sleeves | |

INDICATIONS FOR USE

The Pipeline™ Flex embolization device with Shield Technology™ is intended for endovascular embolization of cerebral aneurysms.

CONTRAINDICATIONS

- Patients who have not received antiplatelet agents prior to the procedure.
- Patients with active bacterial infection.
- Patients in whom antiplatelet and/or anticoagulation therapy is contraindicated.
- The Pipeline™ Flex embolization device with Shield Technology™ should not be used alone as sole therapy for acutely ruptured aneurysms.

POTENTIAL COMPLICATIONS

Potential complications include, but are not limited to, the following:

- Adverse reaction to antiplatelet/anticoagulation agents or contrast media
- Ischemia
- Blindness
- Mass effect
- Coma
- Neurological deficits
- Death
- Parent Artery Stenosis
- Device fracture
- Perforation
- Device migration or misplacement
- Perforator occlusion
- Dissection of the parent artery
- Post-procedure bleeding
- Distal embolization including to a previously uninvolved territory
- Ruptured or perforated aneurysm
- Embolism
- Seizure
- Groin injury
- Stroke
- Headache
- Thromboembolism
- Hemorrhage
- Transient Ischemic Attack (TIA)
- Hematoma or hemorrhage at the puncture site
- Vasospasm
- Hydrocephalus
- Vessel occlusion
- Infection
- Vessel perforation
- Intracerebral bleeding
- Vision impairment
- Chronic pain
- AV fistula
- Pseudoaneurysm

PRECAUTIONS

- Do not use in patients in whom the angiography demonstrates the anatomy is not appropriate for endovascular treatment, due to conditions such as severe intracranial vessel tortuosity or stenosis.
- Do not attempt to re-position the Pipeline™ Flex embolization device with Shield Technology™ after deployment past the resheathing marker.
- Placement of multiple Pipeline™ Flex embolization devices with Shield Technology™ may increase the risk of ischemic complications.
- The appropriate anti-platelet and anti-coagulation therapy should be administered in accordance with standard medical practice.
- A thrombosed aneurysm may aggravate pre-existing, or cause new, symptoms of mass effect and may require medical therapy.
- Do not place the Pipeline™ Flex embolization device with Shield Technology™ in patients whom a pre-existing stent is in place in the parent artery at the target aneurysm location.
- Use of implants with labeled diameter larger than the parent vessel diameter may result in decreased effectiveness and additional safety risks due to incomplete foreshortening resulting in an implant longer than anticipated.
- Person with known allergy to cobalt/chromium alloy (including major elements of cobalt, chromium, nickel, molybdenum or tungsten) may suffer an allergic reaction to the Pipeline™ Flex embolization device with Shield Technology™ implant.
- Person with known allergy to tin, silver, stainless steel or silicone elastomer may suffer an allergic reaction to the Pipeline™ Flex embolization device with Shield Technology™ delivery system.

COMPATIBILITY

The Pipeline™ Flex embolization device with Shield Technology™ is designed to be delivered through a compatible 0.027 (0.69 mm) inside diameter micro catheter. Unconstrained diameter of the Pipeline™ Flex embolization device with Shield Technology™ is 0.25 mm greater than the labeled diameter (on the packaging). Do not use the Pipeline™ Flex embolization device with Shield Technology™ in vessel diameters that are larger than the labeled diameter.

MAGNETIC RESONANCE IMAGING

Non-clinical testing has demonstrated that the Pipeline™ Flex embolization device with Shield Technology™ is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 3 Tesla or less.
- Spatial gradient field of 720 Gauss/cm or less.
- Maximum whole-body-averaged specific absorption rate (SAR) of 4.0 W/kg for 15 minutes of scanning

In non-clinical testing, the Pipeline™ Flex embolization device with Shield Technology™ produced a temperature rise of less than 0.6°C at a maximum whole body averaged specific absorption rate (SAR) of 4.0 W/kg assessed by calorimetry for 15 minutes of MR scanning in a 3 Tesla MR 750 GE Signa 20.0 System MR Scanner.

The Pipeline™ Flex embolization device with Shield Technology™ may create local field inhomogeneity and susceptibility artifacts which may degrade the diagnostic quality of the MRI images. Based on the non-clinical testing of the 5.0 mm device using standard views, the worst case maximum artifact was < 4 mm when subjected to 3.0 Tesla. Local field artifact from the Pipeline™ Flex embolization device with Shield Technology™ may decrease the accuracy of MR angiogram in assessing vessel luminal patency.

MR image quality may be compromised if the area is in the exact same area or relatively close to the position of the Pipeline™ Flex embolization device with Shield Technology™. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant.

DIRECTIONS FOR USE

- Using standard interventional radiographic technique, place the micro catheter tip at least 20 mm past the distal edge of the aneurysm. Gently retract the micro catheter to reduce slack in the micro catheter prior to inserting Pipeline™ Flex embolization device with Shield Technology™.
 - Note:** It is recommended to use a heparinized saline drip to continuously flush micro catheter during Pipeline™ Flex embolization device with Shield Technology™ use.
- Choose a Pipeline™ Flex embolization device with Shield Technology™ with labeled diameter that approximates the target vessel diameter.
 - Select an appropriately sized Pipeline™ Flex embolization device with Shield Technology™ such that its fully expanded diameter is equivalent to that of the largest target vessel. An incorrectly sized Pipeline™ Flex embolization device with Shield Technology™ may result in inadequate device placement, incomplete opening, or migration.
 - The Pipeline™ Flex embolization device with Shield Technology™ foreshortens substantially (50-60%) during deployment. Take device foreshortening into account when deploying the Pipeline™ Flex embolization device with Shield Technology™.
- Choose a Pipeline™ Flex embolization device with Shield Technology™ with labeled length that is at least 6 mm longer than the aneurysm neck.
- Remove packaging hoop from the pouch and pull the distal end of the introducer sheath from the blue clip on the packaging hoop.
- Carefully remove system from packaging hoop until the delivery wire is exposed.
- Partially insert introducer sheath into the rotating hemostatic valve (RHV) at the catheter hub and close the RHV. Confirm back flush of the heparinized saline at the proximal end of the introducer sheath prior to advancing the Pipeline™ Flex embolization device with Shield Technology™ into the micro catheter.
- Advance introducer sheath into the RHV; visually confirm the tip of the sheath is seated deeply in the hub of the micro catheter.
- Secure introducer sheath to the hub by locking down the RHV tightly.
- Advance the proximal end of delivery wire until it aligns with proximal end of introducer sheath.
- Remove the introducer sheath.
 - Note:** The delivery wire has a fluorosafe marker no further than 125 cm from the distal end.
 - Caution:** The fluorosafe marker is only compatible with micro catheters with a minimum length of 135 cm.
- Advance the Pipeline™ Flex embolization device with Shield Technology™ into the micro catheter by pushing the delivery wire until the tip of the delivery wire aligns with the tip of the micro catheter.
 - Caution:** If high forces or excessive friction is encountered during delivery, discontinue delivery of the device and identify the cause of the resistance, remove device and micro catheter simultaneously. Advancement of the Pipeline™ Flex embolization device with Shield Technology™ against resistance may result in damage or patient injury.
 - Caution:** The presence of other indwelling endovascular stents may interfere with proper deployment and function of the Pipeline™ Flex embolization device with Shield Technology™.
- Once the tip of delivery system and micro catheter are aligned, verify that the Pipeline™ Flex embolization device with Shield Technology™ is in the desired location. Distal end of Pipeline™ Flex embolization device with Shield Technology™ should be placed at least 3 mm past the distal edge of the aneurysm.
- Begin to deliver the Pipeline™ Flex embolization device with Shield Technology™ using a combination of unsheathing the Pipeline™ Flex embolization device with Shield Technology™ and pushing the delivery wire simultaneously.

WARNING


Pushing delivery wire without retracting the micro catheter at the same time will cause the open end braid to move distally in the vessel. This may cause damage to the braid or vessel

- After the distal end of the Pipeline™ Flex embolization device with Shield Technology™ has successfully expanded, deploy the remainder of Pipeline™ Flex embolization device with Shield Technology™ by pushing the delivery wire and/or unsheathing the Pipeline™ Flex embolization device with Shield Technology™. Resheathing and/or manipulation of the micro catheter, by locking down the delivery wire and moving both as a system, may facilitate expansion of the Pipeline™ Flex embolization device with Shield Technology™.
 - Caution:** Under fluoroscopy, carefully monitor the tip coil during Pipeline™ Flex embolization device with Shield Technology™ deployment.
 - Caution:** If the delivery wire cannot be retracted into the micro catheter, carefully remove the delivery wire and micro catheter simultaneously.
- Resheathing Instructions: During deployment of the Pipeline™ Flex embolization device with Shield Technology™ resheathing can be performed by advancing the micro catheter while pulling delivery wire.
 - The Pipeline™ Flex embolization device with Shield Technology™ can be resheathed until the resheathing marker has reached the distal marker of the micro catheter (see Figure 2 below).
 - The Pipeline™ Flex embolization device with Shield Technology™ is fully resheathed when the distal marker is retracted completely inside the micro catheter. The system is designed to allow for 2 full cycles of resheathing of the Pipeline™ Flex embolization device with Shield Technology™.

WARNING

Resheathing the Pipeline™ Flex embolization device with Shield Technology™ more than 2 full cycles may cause damage to the distal or proximal ends of the braid.

Figure 2. Pipeline™ Flex embolization device with Shield Technology™ (Resheathing schematic as seen under fluoroscopy, image not to scale).

1. Proximal End of the device
2. Micro Catheter
3. Micro Catheter Distal Marker
4. Resheathing Marker
5. Proximal Bumper
6. Delivery Wire
7. Resheathing Limit
8. Device Detached

16. After the entire Pipeline™ Flex embolization device with Shield Technology™ is deployed, advance the micro catheter through the device making sure not to dislodge the braid. When the micro catheter tip is distal to the Pipeline™ Flex embolization device with Shield Technology™, retract the delivery wire into the micro catheter tip.

Caution: If the catheter cannot be advanced through the Pipeline™ Flex embolization device with Shield Technology™, carefully remove the delivery wire through the Pipeline™ Flex embolization device with Shield Technology™ construct.

17. Carefully inspect the deployed Pipeline™ Flex embolization device with Shield Technology™ under fluoroscopy to confirm that it is completely apposed to the vessel wall and not kinked. If the device is not fully apposed or is kinked, consider using a balloon catheter, micro catheter, or guide-wire to fully open it.

DISPOSAL: The implant and or delivery system should be disposed of or returned to the manufacturer per institutional guidelines.

QUESTIONS AND ANSWERS

Q If excessive friction is experienced during the insertion of delivery system at any time during the delivery of Pipeline™ Flex embolization device with Shield Technology™, what should I do?

A Carefully remove the entire system simultaneously (micro catheter and delivery system).

Q Can I retrieve the Pipeline™ Flex embolization device with Shield Technology™ if the distal end of the Pipeline™ Flex embolization device with Shield Technology™ has expanded at an undesirable location?

A Yes. A partially deployed Pipeline™ Flex embolization device with Shield Technology™ can be resheathed per resheathing instructions, step 15 in the Directions for Use.

Q Can I retrieve a fully deployed Pipeline™ Flex embolization device with Shield Technology™?

A Once fully deployed, the Pipeline™ Flex embolization device with Shield Technology™ cannot be removed. A second Pipeline™ Flex embolization device with Shield Technology™ can be deployed if needed.

Q Can I place a second Pipeline™ Flex embolization device with Shield Technology™ inside another Pipeline™ Flex embolization device with Shield Technology™?

A Yes. A second Pipeline™ Flex embolization device with Shield Technology™ can be placed inside another Pipeline™ Flex embolization device with Shield Technology™. After placing the first Pipeline™ Flex embolization device with Shield Technology™, advance the micro catheter over the delivery wire while keeping the delivery core wire across the Pipeline™ Flex embolization device with Shield Technology™. Position the micro catheter at the desired location and retrieve the delivery wire. Select a new appropriate Pipeline™ Flex embolization device with Shield Technology™ and deploy it as normal.

Caution: Placement of multiple Pipeline™ Flex embolization devices with Shield Technology™ may increase the risk of ischemic complications.

Q If there is a difference between the proximal and distal diameter, which Pipeline™ Flex embolization device with Shield Technology™ diameter do I choose?

A Choose a Pipeline™ Flex embolization device with Shield Technology™ that matches larger (typically proximal) vessel diameter to ensure proper anchoring.

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Micro Therapeutics, Inc.