TrailBlazer™

Support Catheter

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Instructions for use

Device description

The TrailBlazer™ support catheter is an over-the-wire (OTW), single lumen catheter (1) with an atraumatic tapered tip (3). The catheter system is offered in 9 models that are compatible with 0.36 mm (0.014 in), 0.46 mm (0.018 in), and 0.89 mm (0.035 in) guidewires. Depending on the model, the support catheter has a working length of 65 cm, 90 cm, 135 cm, or 150 cm. For specific sizing information, see Table 1 for model specifications.


Figure 1. TrailBlazer support catheter

As shown in Figure 1, the distal catheter shaft has 3 radiopaque markers (2) that aid in positioning the catheter. The distal 40 cm portion of the catheter (4) has a hydrophilic coating. The proximal lumen (6) begins within the manifold (5), transitions through the catheter shaft, and terminates at the distal tip. This lumen is used to pass the catheter over a guidewire. The maximum diameters of compatible guidewires are listed in the Table 1 specifications. The catheter guidewire size and length are printed on the strain relief (7).

Table 1. Size specifications

MODEL	SC-014-135	SC-014-150	SC-018-090	SC-018-135	SC-018-150	SC-035-065	SC-035-090	SC-035-135	SC-035-150
Maximum guidewire diameter mm (in)	0.36 (0.014)	0.36 (0.014)	0.46 (0.018)	0.46 (0.018)	0.46 (0.018)	0.89 (0.035)	0.89 (0.035)	0.89 (0.035)	0.89 (0.035)
Working length cm	135	150	90	135	150	65	90	135	150
Minimum guidewire length cm	180	180	150	180	180	150	150	180	180
Marker band spacing mm	15	15	15	15	15	50	50	50	50
Proximal shaft diameter, mm (in)	0.99 (0.039)	0.99 (0.039)	1.12 (0.044)	1.12 (0.044)	1.12 (0.044)	1.60 (0.063)	1.60 (0.063)	1.60 (0.063)	1.60 (0.063)
Distal shaft diameter mm (in)	0.69 (0.027)	0.69 (0.027)	0.79 (0.031)	0.79 (0.031)	0.79 (0.031)	1.30 (0.051)	1.30 (0.051)	1.30 (0.051)	1.30 (0.051)
Tip outside diameter mm (in)	0.51 (0.020)	0.51 (0.020)	0.58 (0.023)	0.58 (0.023)	0.58 (0.023)	1.04 (0.041)	1.04 (0.041)	1.04 (0.041)	1.04 (0.041)
Minimum guide catheter, mm (Fr)	1.7 (5)	1.7 (5)	1.7 (5)	1.7 (5)	1.7 (5)	2.3 (6)	2.3 (6)	2.3 (6)	2.3 (6)
Guide catheter/sheath min I.D. mm (in)	1.14 (0.045)	1.14 (0.045)	1.14 (0.045)	1.14 (0.045)	1.14 (0.045)	1.68 (0.066)	1.68 (0.066)	1.68 (0.066)	1.68 (0.066)
Minimum introducer sheath, mm (Fr)	1.4 (4)	1.4 (4)	1.4 (4)	1.4 (4)	1.4 (4)	1.7 (5)	1.7 (5)	1.7 (5)	1.7 (5)

Indications for use

TrailBlazer support catheters are percutaneous, single lumen catheters designed for use in the peripheral vascular system. TrailBlazer support catheters are intended to guide and support a guidewire during access of the vasculature, to allow for wire exchanges, and to provide a conduit for the delivery of saline solutions or diagnostic contrast agents.

Contraindications

There are no known contraindications.

Warnings

- This device is supplied **sterile** for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risks of patient infection and compromised device performance.
- If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the device or the vessel lumen. Carefully withdraw the support catheter.
- The TrailBlazer support catheter should be used by physicians with adequate training in the use of the device.
- Always use fluoroscopy or other visualization when manipulating the catheter.

Precautions

- Federal law restricts this device to sale by or on the order of a physician.
- Carefully inspect the sterile package and the device before use. Do not use the device if the packaging or the device is damaged.
- Do not wipe the catheter surface with dry gauze.
- Do not attempt to pass the TrailBlazer support catheter through a smaller-sized introducer sheath than indicated on the label. Refer to Table 1 for size specifications.
- Ensure that the guidewire diameter and length are compatible with the TrailBlazer support catheter. See Table 1 for size specifications.
- To avoid kinking, advance the catheter slowly, in small increments, until the proximal end of the guidewire emerges from the catheter.
- Do not exceed the maximum recommended infusion pressure of 2068 kPa (300 psi).
- Do not advance the catheter into a vessel with a diameter smaller than the outer diameter of the catheter.
- After use, this product may be a potential biohazard. Handle and dispose of the device in accordance with acceptable medical practice and applicable local, state, and federal laws and regulations.

Potential adverse effects

The following complications may result, but may not be limited to the following conditions:

- Access site complications
- Arterial dissection
- Arterial spasm
- Arterial thrombosis
- Catheter fracture with tip separation and distal embolization
- Death
- Distal embolization (air, blood clots, or plaque)
- Intimal disruption
- Local or systemic infection including sepsis
- Perforation and vessel rupture
- Surgical intervention

How supplied

Warning: This device is supplied **sterile** for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risks of patient infection and compromised device performance.

Storage

Store the support catheter away from sunlight at room temperature in a dry place. Do not store catheters where they are directly exposed to organic solvents, ionizing radiation, or ultraviolet light.

Recommended items

Prepare the following items using sterile technique:

- 10 cc syringe filled with sterile heparinized saline
- Three-way stopcock
- Appropriately sized exchange guidewire (refer to Table 1 for size specifications)
- Appropriately sized hemostatic introducer sheath (refer to Table 1 for size specifications)

Preparation

Caution: Carefully inspect the sterile package and the device before use. Do not use the catheter if the packaging or the device is damaged.

- Carefully remove the catheter from the inner package.
- The TrailBlazer support catheter is coated with a hydrophilic coating. Before inserting the catheter, activate the coating by immersing the catheter in normal saline, for approximately 30 to 60 seconds, or by wiping the catheter shaft with a saturated gauze sponge.

- Caution:** Do not wipe the catheter surface with dry gauze.
- To displace air, connect a 10 cc syringe filled with heparinized saline to the luer lock on the manifold and flush out the lumen.

Directions for use

- Introduce the TrailBlazer support catheter percutaneously using an appropriately sized introducer sheath.

Caution: Do not attempt to pass the TrailBlazer support catheter through a smaller-sized introducer sheath than indicated on the label. Refer to Table 1 for size specifications.
- Place the prepared catheter over a previously inserted guidewire.

Caution: Ensure that the guidewire diameter and length are compatible with the TrailBlazer support catheter. See Table 1 for size specifications.

- Caution:** To avoid kinking, advance the catheter slowly, in small increments, until the proximal end of the guidewire emerges from the catheter.
- Carefully advance the TrailBlazer support catheter to the desired vasculature, using fluoroscopic imaging to visualize the marker bands.
Warning: If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the device or the vessel lumen. Carefully withdraw the support catheter.
Warning: Always use fluoroscopy or other visualization when manipulating the catheter.
Caution: Do not advance the catheter into a vessel with a diameter smaller than the outer diameter of the catheter.
 - Carefully remove the catheter, making sure to maintain guidewire position.
Note: To use the catheter for infusion, withdraw the guidewire and attach an infusion device (syringe, power injector, and so forth). Do not exceed 2068 kPa (300 psi) inlet infusion pressure.
Caution: Do not exceed the maximum recommended infusion pressure of 2068 kPa (300 psi).

Table 2. TrailBlazer support catheter flow rates (ml/second)

Model	Guidewire diameter mm (in)	Length cm	Saline		Contrast*	
			1034 kPa (150 psi)	2068 kPa (300 psi)	1034 kPa (150 psi)	2068 kPa (300 psi)
SC-014-135	0.36 (0.014)	135	1.2	1.7	0.5	0.8
SC-014-150	0.36 (0.014)	150	1.2	1.7	0.5	0.8
SC-018-090	0.46 (0.018)	90	2.2	3.1	1.1	1.9
SC-018-135	0.46 (0.018)	135	2.0	3.0	0.9	1.6
SC-018-150	0.46 (0.018)	150	1.9	2.8	0.8	1.4
SC-035-065	0.89 (0.035)	65	7.8	11.8	5.4	8.9
SC-035-090	0.89 (0.035)	90	7.3	10.6	4.3	7.2
SC-035-135	0.89 (0.035)	135	6.0	8.7	3.2	5.5
SC-035-150	0.89 (0.035)	150	6.0	8.7	3.2	5.2

*A mix of 75% Isovue 370 contrast media and 25% saline equivalent solution

Disposal

Caution: After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state, and national laws and regulations.


Warranty Disclaimer

This product is subject to standard Medtronic warranty terms.

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Figur 1. TrailBlazer supportkateter

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© 2017 Medtronic. All rights reserved. Medtronic and Medtronic logo are trademarks of Medtronic. ™* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
ev3 Inc.
Medtronic Australasia Pty Ltd