Mac-Lab User Manual

AltiX BT22 Edition

5222001-1EN
Revision 2.0
US English
© 2022 General Electric Company

SECRET DE
AFACERI

1.1.6 Clinical benefits

Mac-Lab™ displays and records clinical data (ECG waveforms, heart rate, pulse oximetry, respiration rate, CO2, temperature, and hemodynamic measurements) and provides the clinical tools needed to help physicians deliver patient diagnoses.

CardioLab™ displays and records clinical data (ECG waveforms, intracardiac signals, stimulus data, ablation data, pulse oximetry, respiration rate, CO2, temperature, invasive and noninvasive blood pressures) and provides the clinical tools needed to help physicians deliver patient diagnoses.

ComboLab displays and records clinical data (ECG waveforms, heart rate, intracardiac signals, stimulus data, ablation data, pulse oximetry, respiration rate, CO2, temperature, hemodynamic measurements, invasive and noninvasive blood pressures) and provides the clinical tools needed to help physicians deliver patient diagnoses.

The MLCL Client Software displays clinical data (ECG waveforms, heart rate, intracardiac signals, stimulus data, ablation data, pulse oximetry, respiration rate, CO2, temperature, hemodynamic measurements, invasive and noninvasive blood pressures) and provides the clinical tools needed to help physicians deliver patient diagnoses.

1.1.7 Patient benefit

The Mac-Lab, CardioLab, and ComboLab recording systems maintain a comprehensive archive for performed procedures that may include waveform data, reports, and a time-based procedure log. The hemodynamic or electrophysiology study information is archived and available to provide continuity in patient care should additional treatments be necessary if indicated by a healthcare professional. The study information may also be transmitted to an Electronic Medical Records system. All study data stored or transmitted by the Mac-Lab, CardioLab, or ComboLab is encrypted to maintain the security of patient medical records and to ensure private health information is protected.

1.1.8 Intended Use/Intended Purpose

1.1.8.1 Types of systems

The following system names are used:

- **Mac-Lab:** A hemodynamic recording system typically used in catheterization laboratories.
- **CardioLab:** An electrophysiology recording system built on a common platform with Mac-Lab.
- **ComboLab:** A single system that contains the features of both the Mac-Lab and the CardioLab systems.
- **MLCL Client Software:** The Mac-Lab and CardioLab acquisition software that is made available for installation on stand-alone systems to operate as a review workstation.

1.1.8.2 Mac-Lab system

The Mac-Lab system is intended for recording hemodynamic clinical data, which can then be displayed, filtered, digitized, amplified, measured, calculated, and/or transmitted for storage, analysis, and viewing at distributed locations.

3.2.1 Control room components

Acquisition computer

The acquisition computer provides the ability to record the patient's real time waveform data. It uses the Microsoft Windows operating system and has the following drives:

- DVD drive: used for installation and service of the computer.
- SD drive: used to store backup copies of patient procedure data.

The required media for the SD card drive is an SDHC or SDXC card with capacity of at least 16GB, and speed of at least Class 6 (Class 10 and higher is also acceptable).

Desk (optional)

Provides a workspace for control room components.

HD Hub

A small enclosure that houses the required components for HD video (three HD video distribution amplifiers and one HD switcher). The HD Hub distributes video and communication signals to the control and procedure rooms.

Isolation transformer

Provides isolated power for components.

Monitors

For Mac-Lab systems, two monitors are standard to display *Real-Time* and *Review* windows. For CardioLab and ComboLab systems, an optional third monitor displays images acquired from the X-ray and ultrasound systems or can be used as a workspace integrator display.

Printer (optional)

Used to print snapshots, images, and reports after a procedure.

NOTE

Connect the printer only to the hospital network, not to the acquisition computer.

Uninterruptible power supply (UPS)

The UPS provides emergency power to the acquisition computer in the event of a power outage. It will sustain power to the computer for a minimum of 20 seconds. When the UPS battery is low, it will beep. If this happens, do not attempt to replace the battery. Contact a GE representative.

NOTE

The UPS provides uninterrupted power to the acquisition computer. Only the first monitor will be powered during a power outage.

NOTE

When loss of power occurs, close the study, exit the Mac-Lab/CardioLab software, and shut down the acquisition computer. When power is restored, turn on the acquisition computer, start the Mac-Lab/CardioLab software, and continue the study

(see 5.5.1.4 Continue existing study on page 71). See 17.6.12 Force restart the app on page 229.

Video splitter

Used to split the video signal from the *Real-Time*, *Review*, or *Image* monitors for display on additional monitors.

Video switch

Used to provide one switched video output to a remote monitor.

Wired barcode scanner (optional)

Used to scan in supplies and medication used in a procedure.

Wireless barcode scanner (optional)

Used to scan in supplies and medication used in a procedure.

The wireless barcode scanner must be in proximity of the acquisition computer and/or the barcode scanner base to scan items into the procedure. It is recommended to place the barcode scanner back in the cradle when not in use to ensure that the wireless barcode scanner is charged. The wireless barcode scanner used with the Mac-Lab/CardioLab system includes a label to ensure that the proper scanner is used in the event that there are multiple wireless barcode scanning devices in the lab environment.

3.2.2 Procedure room components

Analog IO box (optional)

Distributes ECG and blood pressure signals from the cath amplifier (PDM) to other equipment and provides 4 analog input ports.

CO₂ module (optional)

The CO₂ module acquires the patient's respiration data for the system. The CO₂ module is powered by the PDM Base Station Plus or PDM Slim Connect.

Large display monitor (LDM) (optional)

The LDM is a large widescreen flat-panel LCD monitor used for simultaneous display of several image sources, including the acquisition system.

Monitors

One remote monitor is standard. Switched video is optional between the *Real-Time* and *Review* windows. A second monitor is also optional.

PDM

The PDM acquires patient data for the Mac-Lab system. The PDM is powered by the PDM Base Station Plus or PDM Slim Connect unit.

To view configuration information about the cath amplifier used to acquire patient data, select **Administration > Cath Amplifier Configuration**.

6.4.8 Help menu

To view the online version of the Mac-Lab user manual, select **Help > Mac-Lab Help**.

Select **Help > About Mac-Lab** to view information regarding the Mac-Lab system.

6.5 Real-Time and Review window toolbar functions

The following are additional toolbar functions for the *Real-Time* and *Review* windows.

Sweep/Scroll Speed

To adjust the waveform's sweep/scroll speed (in mm/sec) in the *Real-Time* or *Review* window, highlight the respective window and perform one of the following:

- Press the **Sweep (Up)** or **Sweep (Down)** keys on the keyboard.
- Select the - and + icons in the toolbar.
- Select a sweep speed from the drop-down menu.
- Select **Real-Time > Change Speed** or (**Review > Change Speed**) and the appropriate sweep speed from the listed options (5, 10, 25, 50, 100, 200 or 400). 800 or 1600 are available in the *Review* window.

To add a custom sweep speed to the drop-down list:

1. Highlight the sweep speed and type the new sweep speed number in the white area.
2. Deselect the area or press **TAB**.

NOTE

The new sweep speed will not be added to the **Change Speed** menu option and only the last entered custom speed will be saved for the next study.

6.5.1 Study Configuration

See 6.6 *Study Configuration – Real-Time* window on page 105 or 6.7 *Study Configuration – Review* window on page 114.

6.5.2 Time and amplitude calipers

Time calipers are used to measure the time interval of a signal in milliseconds (ms) or beats per minute (BPM). Amplitude calipers measure the amplitude of ECG signals in millivolts (mV).

NOTE

A maximum of 10 calipers can be displayed in the *Real-Time* or *Review* window at one time.

Amplitude and time calipers can be displayed at the same time.

To perform an amplitude or time measurement in either the *Real-Time* or *Review* window:

1. Click the **Time Calipers** or **Amplitude Calipers** icon in the respective windows.
2. Click a signal in the respective window to drop the first caliper leg.

Icon	See section
--	5.5.1.10 Resolve backups on page 74
--	5.5.1.11 Purge oldest study data on page 74
	5.5.1.12 Restore study on page 74

5.5.1.1 New study

Real-time data acquisition begins when you create a new study.

1. In the *Navigator* window, select the **New Study** icon or select **Study > New**.
2. If User Tracking is set to non-authenticated or authenticated with user selection, enter or select a name in the *Choose Name for User Tracking* window. See 7.1.11.2 *Choose the user on page 123* for more information.

A new study opens and the *Patient Information* window appears.

3. Enter the last name, first name, and MRN of the patient in the appropriate fields. All other fields, including the **Contact** or **Worklist** tabs, are not required to begin a new study and can be completed at any time during the study.

NOTE

To duplicate patient information if the patient returns to the catheterization laboratory for an additional study, highlight the patient name in the catalog and select **New Study**. This patient's demographic information is transferred to the new study and is available for editing.

4. Select **OK**.

NOTE

For emergency cases, you can bypass the *Patient Information* window by pressing **Shift + Enter** on the keyboard. The *Patient Information* window disappears and patient information can be edited at a later time. The following information is entered in the *Patient Information* window when it is bypassed:

- **Last name:** = current time (for example, 02-04-56).
- **First name:** = No name.
- **MRN** = a generated number.
- **Study date:** = current date and time.

5.5.1.2 Patient search

This action is only available if using a CVIS (Cardiovascular Information System) that supports the patient search feature, such as the Centricity Cardio Workflow (CCW).

See 5.7 *Patient information from CVIS on page 79* for information about using a CVIS that supports active patient search.

1. From the *Navigator* window, select the **Patient Search** icon or select **Study > Patient Search**. The *Search for Patient* window appears.
2. Search by **Patient** or **Date of Admission**.

To search by **Patient**:

5.6.2 Export custom data

After you change or import any of the items in the table in [5.6 Import and export custom data on page 77](#), create a backup.

1. Select **Utilities > Export [Lists, Report Formats, Registries, or Custom Form Layouts]**.
2. Click **Next** at the *Data Export Wizard*.
3. Navigate to the export location (for example, a removable disk or network location).

NOTE

If necessary, create a new folder that indicates the item exported and the date (for example, `Report_Formats_15MAR2005`).

The export can take up to several minutes. A confirmation displays when the export completes.

4. Click **OK** to finish.

5.7 Patient information from CVIS

The system allows retrieval of patient demographics, account and order information from the Cardiovascular Information System (CVIS) database.

NOTE

The system DICOM configurations must be specifically configured to connect to the CVIS database. See [18.1.7 Connectivity on page 234](#).

See [5.5.1.2 Patient search on page 69](#) and [5.5.1.3 Schedule search on page 70](#) for details about getting patient information from a CVIS that supports those features.

5.7.1 Active patient search

A CVIS that supports active patient search allows entry of patient demographics and exam history prior to beginning a case study. The Mac-Lab case can also be started from the CVIS. At that time, the study and patient data is transferred to the Mac-Lab. At completion of the study on the Mac-Lab, the study data can be transferred back to the CVIS.

To transfer study data from a CVIS that supports active patient search to the Mac-Lab system, refer to the CVIS operator manual.

The Centricity Cardio Workflow (CCW) supports active patient search.

To use active patient search when the CVIS client is not installed on the system, do the following steps:

1. At the *Navigator* window, select **Study > Get Active Patient**.
The *Get Active Patient* window appears and is populated with the study data from the CVIS.
2. Select **Start New Study**.
The *Patient Information* window appears with the demographics information from the CVIS.
3. Select **OK** and begin the study.

5.7.2 Switch between systems

To switch between the CVIS system and the Mac-Lab system, push the **Action + TAB** keys on the keyboard. Push the **Action** key and tab to the appropriate system.

medications to the master lists. If the CVIS option is not purchased, see [13.1.2.2 Customizing master lists - systems without CVIS on page 173](#).

To select a medication item:

1. Press the **Action + B** keys or select **Study > Medications**.
2. Select the desired medication category in the **Root** directory.
3. Select the medication item by double-clicking the item or highlight the item and press the **Enter** key.

The selected medication is displayed in the *Medications used list* window.

4. Type the dosage in the **Amount** column.
5. If required, select appropriate staff member in the **Ordered By** drop-down list.
6. If required, select appropriate staff member in the **Given By** drop-down list.
7. If appropriate, type the time the medication was stopped in the *Time Stop* window.
8. Press the **Enter** key to confirm or click the **Check** icon next to the time column.
9. To delete a medication item, highlight the item and click the **Delete** icon.

NOTE

To quickly search for a medication, click the **Name** heading in the *Medication pick list* window to alphabetize the list. Begin typing the item to be searched.

The *Medication pick list* category window can be displayed or hidden by clicking the **Show Tree List** icon.

NOTE

An optional barcode scanner can be used to quickly scan in medication items. Scan the barcode on the medication package. The barcode must be pre-entered into the master list before the scanner will recognize the item.

NOTE

The **Amount**, **Subtotal**, and **Selection Total** fields of Medication items in the *Medications used list* will be formatted according to the [F.2 Do Not Use List on page 301](#).

NOTE

The **Administered Unit** for a medication indicates whether the medication is dispensed per unit of time, by patient weight, or both. The **Subtotal** for a medication will be blank if a medication is dispensed per unit of time and the **Time Stop** has not been set. Similarly, the **Subtotal** will be blank for medications dispensed by weight if the patient's weight has not been entered. When the **Subtotal** is blank, the **Selection Total** will indicate that the total is unavailable due to missing data. When all required medication data has been entered, the **Subtotal** and **Selection Total** will indicate these values for the selected medication.

13.1.1.6 Radiology

The **Radiology** tab provides a location to enter information regarding patient radiation exposure. This information can be printed as a report section.

To access the **Radiology** list:

1. Press the **Action + N** keys or select **Study > Radiology**.

2. The **Exam Record DAP**, **Exam Fluoro DAP**, **Exam Total DAP**, **Total Fluoro Time** and **Total runs** will automatically be populated with information from select GE X-ray systems if the X-ray interface option is present. **Exam Total Dose** will also be automatically populated if the x-ray system provides this information.

It is possible to manually populate the **Interventional Fluoro time**, **Diagnostic Fluoro time**, **Cine ID**, **Archive ID** and **CD ID**. **Exam Total Dose** may be manually populated. The **Exam Total DAP** will automatically be calculated if the manual calculation is entered in the **Exam Record DAP** and **Exam Fluoro DAP** information.

To populate those fields, place the cursor in the appropriate field and manually enter the information. Make sure to manually calculate the **Total Fluoro Time**. The **Exam Total DAP** will automatically be calculated if the manual calculation is entered in the **Exam Record DAP** and **Exam Fluoro DAP** information.

3. The information of each cine run will appear on the bottom half of the *Radiology* window if the X-ray interface option is present.
4. Click **Refresh** to refresh the current values displayed by the X-ray system.
5. Click **Close** when finished.

13.1.1.7 Complications

The **Complications** tab contains information on complications that may have occurred during a procedure.

To select a complication item:

1. Press the **Action + M** keys or select **Study > Complications**.
2. Select the desired **Complication** category in the **Root** directory.
3. Select the complication item by double-clicking the item or highlight the item and press the **Enter** key.
4. The selected complication is displayed in the *Complications used list* window.
5. To delete a complication, highlight that item and click the **Delete** icon.

NOTE

To quickly search for a complication, click the **Name** heading in the list item window to alphabetize the list. Begin typing the item to be searched.

The **Complications pick list** category window can be displayed or hidden by clicking the **Show Tree List** icon.

13.1.1.8 Contrast

The **Contrast** tab contains information on contrast agents that were used during a procedure.

To select a contrast item:

1. Select the **Contrast** tab in the *Master Lists* window or select **Study > Contrast**.
2. Select the desired contrast category in the **Root** directory.
3. Select the contrast item by double-clicking the item or highlight the item and press the **Enter** key.
4. The selected contrast is displayed in the *Contrast used list* window.
5. Type the dosage in the **Amount** column.

5. Select the **Display hospital name** check box, if applicable.
6. Select **Printer Setup** to choose the print setup options.

18.1.6 Hospital Information

Select **Hospital Info** in the *System Settings* window.

These settings store the name of the site, address of the site, and lab location. This information is used on window prints, snapshots, and is exported to cardiovascular information systems. It is not automatically imported into reports.

Hospital information is stored locally on standalone systems. In networked environments, hospital information is stored on the INW server and on each acquisition system to allow unique lab locations to be specified. Review workstations do not have their own **Hospital Info** settings.

For networked environments, the system state determines where hospital information is taken from, when it is used.

System state	Hospital Info settings taken from
At <i>Navigator</i> window with System , Host , or Backup catalog selected	Acquisition system
At <i>Navigator</i> window with Network catalog selected	INW server
In study that has Study Location of acquisition system	Acquisition system
In study that has Study Location of INW server	INW server

Configure Hospital Info on an acquisition system

1. From the *Navigator* window, select the **System** catalog.
2. Select **Administration > System Settings** to open the *System Settings* window.
3. Select **Hospital Info** in the *System Settings* window.
4. Enter hospital information in the fields provided. **Lab location** refers to the location of the acquisition system (e.g., **Lab 2**).

Configure Hospital Info on an INW server

1. From the *Navigator* window, select the **Network** catalog.
2. Select **Administration > System Settings** to open the *System Settings* window.
3. Select **Hospital Info** in the *System Settings* window.
4. Enter hospital information in the fields provided. **Lab location** refers to the general location for all of the labs in the network (e.g., **3rd Floor of North Tower**).

18.1.7 Connectivity

Select **Connectivity** in the *System Settings* window.

18.1.7.1 Local system

DICOM settings

Use the information in the following table to populate the fields on the **Local System** tab.

Field	Description
AE Title	Enter the name of the system. You can change the name based on the number of systems (e.g., two or more systems can be renamed MLab A , MLab B , etc.). NOTE AE Title is not necessarily the same as the machine or host name. AE Title refers to the DICOM name, not the machine name.
Listen Port	The default is 1225 for normal transport and 2762 for secure transport. This value is typically only changed if it causes a conflict.
DICOM server security The DICOM server can run in one of three security modes: None , Legacy , or Enhanced . Choose the option that accommodates all remote DICOM systems that use the Mac-Lab/CardioLab system's DICOM services. For example, if the X-ray system supports Legacy mode but the ultrasound system does not support secure transport, select None . The remote DICOM system must be configured to use the same security mode. If the configurations do not match, the systems cannot communicate with each other. Certificate validation may require configuring both the Mac-Lab/CardioLab system and the remote DICOM system with the other system's certificate data. Refer to the Mac-Lab/CardioLab Privacy and Security User Manual (PN 5222004-1EN) for details.	
Server Security > None	Transport security is disabled. Select this option if any remote DICOM host that uses the Mac-Lab/CardioLab system's DICOM services does not support secure DICOM transport.
Server Security > Legacy	This is the secure transport compatibility mode. Select this option if any remote DICOM host that uses the Mac-Lab/CardioLab system's DICOM services supports only the retired DICOM transport methods.
Server Security > Enhanced	This transport mode utilizes certificate validation. Select this option if every remote DICOM host that uses the Mac-Lab/CardioLab system's DICOM services supports the current DICOM transport methods and certificate validation.
Maximum search responses to accept	Enter the appropriate number. The default is 150. The maximum search responses returned is 999. This setting prevents the system from being overrun by broad queries. For example, if this value is set to 150, a search for all patients on the CVIS system displays no more than 150 results.

When starting a case study using demographics provided by the CVIS, the following fields in the *Patient Information* window are disabled by default in order to ensure consistency of data among the CVIS, Mac-Lab, and X-ray systems: **Last Name**, **Middle Name**, **First Name**, **MRN**, and **Date of Birth**. When a case is started from a scheduled study, the **Requested Procedure Identifier** and **Scheduled Procedure Identifier** fields on the **Worklist** tab are also disabled.

If the demographic fields listed previously need to be changed, uncheck the **Lock required patient demographics fields from Modality Worklist Host** check box at any time before or during a case.

After demographic information in the *Patient Information* window is changed, the system displays a message that states that these fields are required to ensure consistency with associated systems, such as the x-ray or CVIS. When prompted, you may choose to apply these changes or not as follows:

- Select **Yes** to apply changes made to any field in the *Patient Information* window.
- Select **No** to apply changes to any field except **Last Name**, **Middle Name**, **First Name**, **MRN** and **Date of Birth**.
- Select **Cancel** to return to the *Patient Information* window without changing anything.

It is recommended to return to the *System Settings* window and select the **Lock required patient demographics fields from Modality Worklist Host** check box after the changes have been made.

Chapter 1 Introduction

1.1 Manual information

1.1.1 Purpose

This manual is intended solely as a reference for the operation of the Mac-Lab System, including the Mac-Lab application of the MLCL Client Software. It is not intended as a replacement for thorough and complete training on the system capabilities and operation by a qualified instructor; it is intended as a supplement to such training.

It is important that all of the instructions in this manual be followed. However, these instructions in no way supersede current medical practices regarding patient care and safety.

1.1.2 Product references

The name of the product described in this manual is Mac-Lab. It will also be referred to as "the system" throughout this manual.

1.1.3 Printed copies

Contact your local GE representative to request a paper copy of this manual.

1.1.4 Intended audience

This manual is written for clinical professionals who use, maintain and/or troubleshoot the Mac-Lab. Clinical professionals are expected to have a working knowledge of appropriate medical procedures, practices and terminology used in the treatment of patients.

1.1.4.1 Operators

- Operators can record and calculate clinical data from pediatric or adult patients as they undergo various procedures. Clinical data may be manually entered or acquired via interfaced GE amplifier modules or other interfaced systems.
- It is recommended that Mac-Lab/CardioLab operators be trained to become familiar with the capabilities and operations of the Mac-Lab/CardioLab system. The following training options are available:
 - Classroom training
 - On-site customer training
 - Remote training

1.1.4.2 System administrator

- Coordinating the overall MacLab/CardioLab system is the responsibility of the system administrator. This person has access to additional system supervisory functions, including operational setup parameters and system backup functions.

6.4.1.1 Case event

The Mac-Lab allows you to link a case to a primary physician and to track the supplies and medications used by that physician during the case. This information is primarily used for statistical purposes.

The following is the workflow recommended when using the **Case Event** function.

1. Add the primary physician to the staff master list (**Study > Staff**) prior to the start of the case event. If desired, add the supplies and medications. See [13.1 Master lists on page 167](#) for instructions on how to select staff, supplies and medication.
2. After the lists have been chosen, begin logging time for the case study by selecting **Study > Case Event**.

The *Case Event* window opens.

3. Select the **Case Type** and **Primary Physician**. Select **OK**.

An entry appears in the *Log* window that indicates that the case has started. The system links any medication or supply used by the primary physician to that case.

4. To change the primary physician or to track supplies or medications for an additional diagnostic or interventional case:
 - a. Select **Study > Case Event**.

The logging of the current case stops and a new case begins. The log is updated with a new case start entry.

- b. The *Case Event* window appears. Select a different **Primary Physician** or **Case Type**. Every supply or medication used after you select **OK** is linked to that physician.

NOTE

Add the physician to the staff list if that physician has not previously been used in the case and does not appear in the **Primary Physician** field.

To change the physician or case type, double-click in the **Summary** column of the *Log* window at **Case Event**. The *Case Event* window appears. Choose a different physician or case type.

To launch a new case event using a macro instead of selecting **Study > Case Event**, see [Chapter 14 Macros on page 179](#) for instructions.

6.4.1.2 Intervention

The *Intervention* window provides the user the ability to specify between coronary, graft and peripheral interventional events. Each *Intervention* window provides specific attributes and treatment options and gives the user the ability to tie inventory supplies to that intervention event. Each of these events can be opened by a macro command, a menu command or through the **Intervention** icon in the *Log* window. The stopwatch is linked to an event each time it is started or stopped. The user also has the ability to have multiple windows open simultaneously to provide better documentation.

Table 6-3 Intervention types and commands

Intervention type	Menu command
Coronary intervention	• Select Study > Intervention > Coronary. • Click the Create new intervention event icon in the <i>Log</i> window and click Coronary.

Table 6-3 Intervention types and commands (Table continued)

Intervention type	Menu command
Graft intervention	- Select Study > Intervention > Graft. - Click the Create new intervention event icon in the <i>Log</i> window and click Graft.
Peripheral intervention	- Select Study > Intervention > Peripheral. - Click the Create new intervention event icon in the <i>Log</i> window and click Peripheral.

6.4.1.2.1 Intervention attempts

To begin an intervention attempt, select **Study > Intervention** and the appropriate event (**Coronary**, **Graft** or **Peripheral**), or select the **Intervention** icon in the *Log* window.

To create a macro that opens an intervention event, see [Chapter 14 Macros on page 179](#).

NOTE

Multiple intervention events can be performed at the same time utilizing separate windows.

The *Intervention Event* window appears.

1. Select a lesion number and corresponding vessel segment using the drop-down menu.
2. Select the **Lesion Attributes**, **Treatment Attributes**, and **Supplies**.
3. Select a treatment type from the **Treatment Modalities** section of the *Intervention* window by selecting the **Insert Item** icon and selecting the treatment type from the drop-down menu.

Supplies can be added at any time during or after the event.

4. Do one of the following to begin the attempt:
 - Select the **Start New Attempt** button. A **Timer** appears.
 - Select the **Add New Attempt** button or the **Deploy** button. A new attempt is created but no timer appears.

NOTE

Each segment vessel combination is unique to a lesion. There cannot be identical segment vessel combinations on another lesion.

5. Add, edit, or delete treatment attributes by selecting the **Insert Item**, **Edit Item/Folder**, or **Delete** icons in the **Treatment Modalities** field. Deleting a treatment deletes all of the supplies and attempts associated with that treatment. Individual attempts may be deleted from the *Log* window.
6. Add, edit, or delete supplies that have already been added to the study by selecting the **Insert Item**, **Edit Item/Folder**, or **Delete** icons in the **Supply** field. To add a supply, select the **Add Supply** button.
7. Select the check box next to **Save vitals** to save the vitals to the *Log* window.

NOTE

If another attempt is made on the same lesion for the same treatment, the **Attempt Number** automatically increments (for example, **Attempt 2 of 2**).

6.4.1.2.2 Treatment selections

Treatment selections are the same in all three *Intervention* windows:

7.3.2 ST Segment and the Log window

An ST Segment Event is placed in the *Log* window with each ST Segment reference obtained.

The ST Segment event contains the time of the event and an automatically generated segment number (example ST Segment: 1) description. Any saved segment may be brought into the *ST Segment* window by selecting the event time. This will also update the signal data in the *Review* window. Double-clicking on the event summary will bring up the *ST Segment Event Description* dialog box, in which the segment description can be changed.

An ST Segment section may be added to any report format. This section will display the time of the ST Segment, as well as the ST Segment values for all surface leads. The report will separate the segment measurements by phase. See [Chapter 16 Reports on page 199](#).

7.4 Image Viewer

Use Image Viewer to view and annotate snapshots and images in a study.

During a study, you can annotate a captured image to highlight unique characteristics of the captured data. Annotations can be saved with the image as part of the study. The image and annotations can also be saved to a file or copied to the clipboard.

7.4.1 Annotation

7.4.1.1 Add annotation objects

Add lines, rectangles, ellipses, polylines, polygons, pointers, freehand drawings, highlights, text, and notes to the image.

1. Select an annotation type in the annotation toolbar, then left-click and drag.
2. When adding a text or note annotation, press **Ctrl + Enter** to create a new line of text below the current line.
3. To complete adding a text or note annotation, press **Enter** or left-click outside of the object.

7.4.1.2 Select annotation objects

Copy, paste, move, resize, delete, and change properties of annotations.

- To select an individual object, move the cursor over the object until the cursor turns into a circle, then left-click.
- To select multiple objects:
 1. Hold down the shift key.
 2. Move the cursor over an object until the cursor turns into a circle, then left-click.
 3. Repeat the previous step for each object you want to select.
 4. Release the shift key.
 - Alternatively, draw a focus box around the objects by holding down the left-mouse button while dragging.
- To select all objects, right-click in the image and choose **Select All**.

6.3.4 FFR and DFR equalization

When using FFR or DFR in a study, the **P1** and **P2** channels must be equalized before you do an FFR or DFR measurement. When the catheter and pressure wire are in place, right-click **FFR/DFR Value Display** to do the equalize function.

- The FFR or DFR display must be turned on. Select either **Real-Time > FFR Display On** or **Real-Time > DFR Display On** to turn on the display or toggle between them..
- To turn on the FFR or DFR display with a macro, see [Chapter 14 Macros on page 179](#) for instructions.

After you select **Equalize**, the **FFR/DFR Value Display** shows that the system is **Equalizing**. This process may fail if the catheter and pressure wire are too far apart. If this happens, the FFR/DFR Value Display again shows **Equalize**. The largest allowable offset between the two channels is 10 mmHg, and the largest allowable time difference is 100 ms. Adjust the position of the catheter and pressure wire and select **Equalize** again.

When equalization is complete, the system calibrates the **P1** and **P2** channels so that mean, amplitude, and time are aligned. The system inserts a note into the study with the mean calibration value, which is called the offset. The system also displays the Pd/Pa ratio in the pressure status area. You can now do FFR or DFR measurements. See [12.3.9 Do an FFR measurement on page 159](#) or [12.3.10 Do a DFR measurement on page 160](#) for more information.

The calibration resets when the FFR or DFR channel is turned off, the study is closed, or if either **P1** or **P2** is zeroed. Calibration is maintained when you toggle between FFR and DFR without turning off the channel. This allows the system to change between FFR and DFR measurements without recalibrating.

- Follow FFR or DFR accessory preparation and checks recommended by the accessory manufacturer. Ensure accessories are not damaged prior to use.
- To equalize FFR or DFR (equalize **P1** and **P2** channels) using a macro instead of right-clicking the **FFR/DFR Value Display**, see [Chapter 14 Macros on page 179](#) for instructions.

6.3.5 Shortcut to turn channels on and off

6.3.5.1 Turn off a channel

Channels can be turned on or off quickly using the following shortcut:

If the setting remains unchecked, any operator will be able to change these required demographics in every subsequent study done on this system.

INNOVA CENTRAL/Touch Screen settings

For INNOVA CENTRAL/Touch Screen configuration, refer to [Appendix D INNOVA CENTRAL/Touch Screen](#) on page 295.

18.1.7.2 Remote hosts

X-Ray System Settings

Select the *Incoming Remote Hosts* tab to access the **X-Ray System Settings**.

Use the information in the following table to populate the **X-Ray System Settings** fields.

Field	Description
AE Title	Enter the DICOM name of the X-ray system in the laboratory. NOTE AE Title is not necessarily the same as the machine or host name. AE Title refers to the DICOM name, not the machine name.
Hostname or IP Address	Enter the host name or IP address of the x-ray system in the laboratory.
Port Number	Enter the appropriate connection port.
Modality	The default is XA . Contact your service representative to change this field.
Text Encoding	Determines what character set to use to package data sent to the X-ray system in the laboratory. To change the default, pick one that the X-ray system supports from the drop-down menu. Refer to the DICOM Conformance Statement for the X-ray system for a list of its supported character sets.
Client Security Outgoing connections to the remote host are made in one of three security modes: None , Legacy , and Enhanced . The remote DICOM system must be configured to use the same security mode. If the configurations do not match, the systems cannot communicate with each other. Certificate validation may require configuring both the Mac-Lab/CardioLab system and the remote DICOM system with the other system's certificate data. Refer to the Mac-Lab/CardioLab Privacy and Security User Manual (PN 5222004-1EN) for details.	
Client Security > None	Transport security is disabled. Select this option if the remote DICOM host does not support secure DICOM transport.
Client Security > Legacy	This is the secure transport compatibility mode. Select this option if the remote DICOM host supports only the retired DICOM transport methods.
Client Security > Enhanced	This transport mode utilizes certificate validation. Select this option if the remote DICOM host supports current DICOM transport methods and certificate validation.
Automatic Exam Start	Select this option to enable the Mac-Lab/CardioLab system to initiate a study on the Innova system. NOTE Only select this option if the Innova system supports the automatic exam start feature.

1.1.8.3 CardioLab system

The CardioLab system is intended for recording electrophysiology clinical data, which can then be displayed, filtered, digitized, amplified, measured, calculated, and/or transmitted for storage, analysis, and viewing at distributed locations.

1.1.8.4 ComboLab system

The ComboLab system is the combination of both the Mac-Lab and CardioLab systems intended for recording hemodynamic and electrophysiology clinical data, respectively, which can then be displayed, filtered, digitized, amplified, measured, calculated, and/or transmitted for storage, analysis, and viewing at distributed locations.

1.1.8.5 MLCL Client Software

The MLCL Client Software is intended for recording, documenting, and/or reviewing clinical data for hemodynamic and electrophysiology procedures, which can then be displayed, filtered, digitized, amplified, measured, calculated, and/or transmitted for storage, analysis, and viewing at distributed locations.

1.1.9 Compliance

The system complies with the following standards:

- IEC/EN 60601-1 3rd ed (Amendment 1) Medical Electrical Equipment, Part 1 General Requirements for Basic Safety and Essential Performance
- IEC/EN 60601-1-2 4th ed Electromagnetic disturbances - Requirements and tests
- IEC/EN 60601-1-4 Programmable electrical medical systems

1.1.10 Electromagnetic Disturbances (EMD)

The Mac-Lab, CardioLab, and ComboLab systems are intended for use as specified within the Indications for Use.

WARNING

The Mac-Lab, CardioLab, and ComboLab systems are not intended to be used together with high-frequency surgical equipment. The combined use may result in decreased immunity of the Mac-Lab, CardioLab, and ComboLab system.

Electromagnetic disturbances do not impact the basic safety of the Mac-Lab system. The Mac-Lab system does not have essential performance, as defined in IEC 60601-1.

Electromagnetic disturbances do not impact the basic safety of the CardioLab system.

Electromagnetic disturbances could result in lost or degraded display and accuracy of intracardiac waveforms, which is the essential performance of the CardioLab system. If this occurs, intracardiac signal noise is possible.

- Refer to the troubleshooting section of the Mac-Lab/CardioLab/Centricity Cardiology INW Service Manual (PN 5222010-1EN) for instructions for resolving the different types of intracardiac signal noise.

2. Hover the mouse pointer over a measurement event in the *Cath Measurement* window to show a tool tip that contains the time of the measurement event and the complete event summary.
3. Select a measurement event or phase event in the *Cath Measurement* window to update the *Review* window to the signal data for that event.

7.8 Workspace Integrator display

Your acquisition system may be equipped with Workspace Integrator. This feature provides an additional display on which you can simultaneously view a third-party application. When you launch the Mac-Lab/CardioLab application, the application fills all displays except for the Workspace Integrator display. The display that doesn't show the Mac-Lab/CardioLab application is the Workspace Integrator display. Workspace Integrator also allows you to use the same keyboard and mouse to interact with both applications.

This section refers to the second application as a third-party application. However, the second application may be a GE Healthcare application, such as the CCW Client.

7.8.1 Set up the Workspace Integrator display

Some third-party applications require additional configuration to launch them. To simplify launch of these applications, you may want to pin the application to the **Start** menu, pin the application to the taskbar, or put a shortcut to the application on the desktop.

For example: If you use the Remote Desktop Connection tool to access your third-party application, do the following to pin Remote Desktop Connection to the **Start** menu or taskbar:

1. Select the **Start** button on the taskbar.
2. Select **Windows Accessories**.
3. Right-click **Remote Desktop Connection**.
 - If you want to pin this application to the **Start** menu, select **Pin to Start**.
 - If you want to pin this application to the task bar, select **More > Pin to taskbar**.

You may need to work with your site network administrator to prepare your system to launch the third-party application.

If you are using Citrix™ with Workspace Integrator, you may need to contact your site network administrator to install and configure Citrix to work with third-party applications.

Make sure to launch, close, and re-launch the application running on the Windows Integrator display to confirm that it opens on the correct display. Also confirm that the application's dialogs only open on the Windows Integrator display. It is the customer's responsibility to set up the Workspace Integrator's application, and any dialogs, so that they appear only on the Workspace Integrator display.

7.8.2 Use the Workspace Integrator display

1. Launch your third-party application.

You may need to provide application credentials or other information to access the third-party application.

2. Launch the Mac-Lab/CardioLab application.

The Mac-Lab/CardioLab application is in all of the displays except one. The third-party application is in the remaining display.

6.4.1.4 Export data

To manually export HL7 data to a third party device, select **Study > Export Data**.

NOTE

To export study reports selectively in PDF format, you must include the # at the beginning of the report format name when you generate the report. Reports that are generated in this format are included in the exported data when active translations include the **Selective Reports in PDF** option.

6.4.1.5 Close study

To close the study, select **Study > Close Study**.

NOTE

If the **Data export on study close** option is set, wait at least 5 minutes to ensure completion of the export. During the wait period, do not log off Windows or turn off the machine.

6.4.2 Measurements menu

6.4.2.1 Measure single active pressure – *Real-Time* window

1. Activate the pressure to measure by pressing the **P1**, **P2**, **P3**, or **P4** key on the keyboard or click the **P1**, **P2**, **P3**, or **P4** channel on the *Real-Time* window.
2. Select **Measurements > Active Pressure**, click the **Active Pressure Measurement** icon on the *Real-Time* window toolbar, or press the **Meas** key on the keyboard.

The pressure measurement is saved to the *Review* window and the corresponding numeric result is displayed in the *Log* window.

6.4.2.2 Measure single active pressure – *Review* window

1. In the left side of the *Review* window, activate the pressure to measure by clicking the desired pressure label.
2. Select **Measurements > Active Pressure**, click the **Active Pressure Measurement** icon on the *Review* window toolbar, or press the **Meas** key on the keyboard.

The pressure measurement is displayed in the *Review* window and the corresponding numeric result is displayed in the *Log* window.

6.4.2.3 Measure all pressures – *Real-Time* window

1. Activate a pressure by pressing the **P1**, **P2**, **P3**, or **P4** key on the keyboard, or click the **P1**, **P2**, **P3**, or **P4** channel on the *Real-Time* window.
2. Select **Measurements > All Pressures**, click the **All Pressure Measurements** icon in the *Real-Time* window or press the **Shift + Meas** keys on the keyboard.

The pressure measurements are saved to the *Review* window and the corresponding numeric result(s) is displayed in the *Log* window.

6.4.2.4 Measure all pressures – *Review* window

1. In the left side of the *Review* window, activate the pressure to measure by clicking a pressure label.

2. Select **Measurements > All Pressures**, click the **All Pressure Measurements** icon in the *Review* window, or press the **Shift + Meas** keys on the keyboard.

The pressure measurements are displayed in the *Review* window and the coLogrresponding numeric result(s) is displayed in the window.

6.4.2.5 Relabel and remeasure a pressure

To relabel and remeasure a pressure if the label was incorrect in the *Real-Time* window when a measurement was made, perform the following steps:

1. Make the inaccurately labeled pressure signal active in the *Review* window by clicking the pressure label. The label will turn orange and a rectangle will be placed around it.
2. Relabel the pressure correctly using the methods described in [12.2.1 Change the pressure site label on page 154](#).
3. Delete the erroneous pressure event in the *Log* window by right-clicking on the time and select **Delete > Yes**.
4. With the appropriate pressure still active in the *Review* window, remeasure the pressure either using menu commands, icons, or the keyboard.

6.4.2.6 FFR or DFR

See [12.3.9 Do an FFR measurement on page 159](#) or [12.3.10 Do a DFR measurement on page 160](#).

6.4.2.7 Pullback

See [12.4.1 Initiate a pullback recording on page 161](#).

6.4.2.8 Valve

See [6.5.18 Valve window – Review window on page 104](#).

6.4.2.9 Gradient (P-P)

See [12.6 Peak-to-Peak window on page 165](#).

6.4.2.10 TDCO

See [11.1.2.1 TDCO setup on page 149](#).

6.4.2.11 Manual CO

See [11.1.1 Manual cardiac output on page 149](#).

6.4.2.12 Blood analysis

See [11.1.3 Fick cardiac output on page 151](#).

6.4.2.13 Resistance

To display resistance results, select **Measurements > Resistances**.

6.4.2.14 Stroke work

To display stroke work results, select **Measurements > Stroke Work Results**.

6.4.2.15 A wave analysis

To display A wave analysis, select **Measurements > A wave Analysis On** or execute a macro command. See [Chapter 14 Macros on page 179](#) for instructions on building a macro.

6.4.2.16 Stat vitals

WARNING

AUDIBLE INDICATORS MALFUNCTION

The audible indicators are for reference only. The audible indicators are not designed for use as a patient alarm.

WARNING

VISUAL INDICATORS MALFUNCTION

The visual indicators are for reference only. The visual indicators are not designed for use as a patient alarm.

To perform a stat vitals measurement, select **Measurements > Stat Vitals** or press the **Action + L** keys on the keyboard.

NOTE

The current heart rate, SpO₂, respiration rate, temperature, blood pressure (either invasive or NBP as indicated in the measurement configuration (see [6.6.7.2 Measurements on page 114](#)) and optional inspired and/or expired CO₂ are recorded in the Log window.

6.4.2.17 Start/Stop NBP

NOTE

The NBP key toggles the action of the blood pressure cuff (inflate/deflate) from any acquisition or review system. To avoid confusion, initiate NBP from one designated system only. Ensure all users understand that an NBP request may be initiated from any acquisition or review system and once initiated, will be cancelled if the NBP key is pressed again.

NOTE

NBP inflation pressures are displayed on both the acquisition and review systems, but cuff inflation error message windows are only displayed on the acquisition system.

1. To do a stat NBP measurement, select **Measurements > Start/Stop NBP** or press the **NBP** key on the keyboard. The cuff inflates and the measurement, heart rate, SpO₂, respiration rate and temperature are recorded into the Log window.
2. To stop a cuff in the middle of an inflation, select **Measurements > Start/Stop NBP** or press the **NBP** key on the keyboard.

6.4.2.18 Auto-log vitals

To automatically log vitals at time intervals defined in the *Measurement Configuration* window (see [6.6.7.2 Measurements on page 114](#)), select **Measurements > Auto-Log Vitals**.

NOTE

The current heart rate, SpO₂, respiration rate, temperature, blood pressure (either invasive or NBP as indicated in the Measurement Configuration (see 6.6.7.2 Measurements on page 114) and optional inspired and/or expired CO₂ are recorded in the Log window.

6.4.2.19 Auto-log NBP (optional)

To automatically log non-invasive blood pressure measurements with vitals, select **Measurements > Auto-Log NBP**. Set the vitals configuration to **Auto-Log NBP** when using NBP as the vital logging pressure source.

Auto-Log Vitals is most useful when you use an invasive blood pressure as the blood pressure channel.

Double-click the vitals event in the Log window to open a vitals event window. You can enter a level of consciousness and a level of pain in the vitals event window. Because the auto-log NBP records a set of vitals with each reading, this information may be redundant.

6.4.2.20 Vital configuration

Select **Measurements > Vital Configuration** to open the *Vital Configuration* tab in the *System Settings* window. See Vital Configuration.

6.4.2.21 Audible indicators

NOTE

The following instructions only apply when audible indicators are enabled from the *Vital Configuration* tab (see Vital Configuration).

WARNING**AUDIBLE INDICATORS MALFUNCTION**

The audible indicators are for reference only. The audible indicators are not designed for use as a patient alarm.

1. Perform one of the following to silence the audible indicator:
 - Click the **Turn On Or Off Audible Indicator** icon in the toolbar.
 - Select **Measurements > Silence Audible Indicator** from the *Real-Time* or *Review* screen.
 - Press **Ctrl + D** on the keyboard.
2. Click **Yes** at the confirmation window.

NOTE

A confirmation window is not displayed if the icon is selected from the toolbar.

6.4.2.22 Data extraction

All channels of data can be exported from the system for a user-selectable time. The data is extracted in either binary or text.

On large data sets, the data extraction operation may take several minutes to complete. During this time, the application does not respond to other user inputs. If you attempt to interact with either the application or the *Data Extraction* window, the system displays an error message and the data extraction bar stops advancing. If this happens, ignore the message. The operation will complete successfully even if the data extraction bar stops advancing.

Manual de utilizare Mac-Lab

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5222001-1RO
Revizia 1.0
Romanian | Limba română
© 2022 General Electric Company

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