PMMA Intraocular Lens

DESCRIPTION / INTENDED USE

S_{IDA}-LENS PMMA Intraocular Lens (IOL) manufactured by Sidapharm, is an optical device for the replacement of natural crystalline lens in human eyes.

The IOL comprises of two parts; the central clear optic (which acts as visual zone) and the peripheral haptic (which helps in lens anchorage). It is made from clinical grade Polymethyl Methacrylate (PMMA) with an added UV-absorbing compound.

To allow the surgeon's flexibility in meeting individual patient requirements, Sidapharm offers several distinct types of lenses, which range from 6.00 to +35.00 diopter* (please refer to IOLs box for each one's specific characteristics).

SUPPLY FORMAT

S_{IDA}-LENS PMMA lenses are supplied in a lens carrying case made of Polystyrene contained in a heat-sealed medical-grade pouch, sterilized using Ethylene oxide (ETO). The contents of the pouch are sterile unless the package is opened or damaged.

INDICATIONS

These IOLs can be used for visual correction of aphakia in adult patients whose cataract lens has been removed by extra-capsular cataract extraction. The lens is intended to be placed in the capsular bag.

PRECAUTIONS

- Do not keep the lens in direct sunlight or at temperature greater than 45°C. Keep from freezing.
- Do not use if the sterile packaging has been opened or damaged.
- Only skilled surgeons with experience in either viewing and/or assisting numerous surgical implantation and successfully completed at least one course on IOL implantation, should attempt implantation of these lenses.
- Pouch should be opened only under sterile conditions.
- Do not soak or rinse lens in solutions other than sterile BSS or equivalent of such.
- Do not attempt to autoclave and resterilize this lens.
- Do not autoclave this lens.
- Handle the lens carefully. Locking forceps or needle holders should never be used to pick up lenses.
- Do not reshape the supporting structures (haptics).
- The patient must be advised that the doctor or the medical center should be informed of any side effects not referred in this information.
- Do not use lens if it is accidentally dropped. The implications of reuse other than adverse reactions are not clinically well known.

CONTRA-INDICATIONS

Patients with the following conditions are not suitable candidates for IOL operation. Otherwise, there is risk for their eyesight.

- Patient with eyes having fore shortened anterior segments, e.g. microphthalmic or certain forms of chronic angle - closure glaucoma
- For anterior chamber implantation, previous history of, or predisposition to retinal detachment
- Patients in whom the IOL may interfere with the ability to observe, diagnose or treat posterior segment disease
- Recurrent anterior or posterior segment inflammation of unknown aetiology
- Abnormality of the iris which would preclude adequate fixation of the lens, such as aniridia, hemi-iridectomy or severe atrophy.

WARNINGS

- The effectiveness of UV-absorbing lenses in reducing the incidence of retinal disorders has not been established.
- The safety of the use of the Neodymium-YAG laser and IOLs with UV-absorbing materials has not been established. The physician is urged to use extreme caution in such cases where a patient with UV-absorbing IOLs is treated with a Neodymium-YAG laser.
- The compression force exerted on the eye tissue by the lens is not established. The physician should have knowledge in the selection of the type of lens depending on the eye dimensions.
- "In-the-bag lenses" should be used only when capsulorhexis is performed and the posterior capsule is in good condition.
- Care should be taken to avoid breakage of haptic while inserting lens through the scleral tunnel of small incision.

ADVERSE REACTIONS

The reuse of IOL may cause:

- Cross contamination which will lead to microbial infection.
- Improper mechanical properties.
- Repeated use may affect aesthetic appearance of lens which will lead to visual defects,

INSTRUCTIONS FOR USE

Do not attempt to expand the haptics, flex them out of the plane of the lens, twist or torque the lens, in order to avoid breakage.

- In a sterile environment, peel apart to open the pouch and remove the lens case.
- Unscrew the lens case with care and grasp the lens optic using smooth edged forceps and lift.
- Soak or rinse lens in sterile balanced salt solution or sterile normal saline solution.
- The lens, due to some static charge produced while opening the lens case, may stick to the cap in rare instances.
- Ensure that the lens is in good condition for optic and haptic surfaces for adherence of any particles.
- Prior to implantation examine the lens for power, type, proper configuration and optical surfaces.
- During insertion handle lenses by haptic portion only.
- Carton box contains extra labels. These are for convenience in maintaining and reporting records of implantable lenses during clinical investigation. One of the labels could be affixed in patient's case sheet, for future reference.
- Rinse the lenses in sterile B.S. Solution before implantation, to minimize the static discharge which enhances the microbial contamination.

SPECIFIC INSTRUCTIONS

- The need for a secondary iridectomy for pupillary block may be prevented by one or more iridectomies at the time of IOL implantation, especially in cases of anterior chamber IOL implantation
- Physicians considering lens implantation in patients with cataracts associated with macular degeneration or retinal degeneration, glaucoma, or corneal diseases, should be aware that the patient may not achieve improvement in central visual acuity caused by these pathologies
- Secondary glaucoma has been reported occasionally in patients with pre-existing glaucoma who received posterior chamber lens implants. The intraocular pressure of implant patients with glaucoma should be carefully monitored post-operatively
- The safety and effectiveness of a posterior chamber lens if placed in the anterior chamber has not been established. Such implantation has been shown to be unsafe in some cases
- The "A" constant value is only an estimate value. It is recommended that surgeon derives his/her own value based on his/her clinical experience.

Applicable Models: SDANT (anterior chamber)
SDHS, SDALP (posterior chamber)

* Out-of-range dioptries may be produced upon request

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