Instruction for use

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«Chlamydia IgM-ИФА»

**A SOLID-PHASE ENZYME IMMUNOASSAY
FOR THE QUALITATIVE DETERMINATION
OF IgM ANTIBODIES TO CHLAMYDIA SPP.
IN HUMAN SERUM OR PLASMA**

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
CONTENT

1. INTENDED USE	10
2. SUMMARY AND EXPLANATION	10
3. PRINCIPLE OF THE TEST	11
4. WARNINGS AND PRECAUTIONS	11
5. KIT COMPONENTS	12
6. SPECIMEN COLLECTION AND STORAGE	13
7. TEST PROCEDURE	13
8. QUALITY CONTROL	15
9. CALCULATION OF RESULTS	15
10. EXPECTED VALUES	15
11. PERFORMANCE CHARACTERISTICS	15

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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUALITATIVE DETERMINATION OF IgM ANTIBODIES TO *CHLAMYDIA TRACHOMATIS* IN HUMAN SERUM OR PLASMA

1. INTENDED USE

A solid-phase enzyme immunoassay for the qualitative determination of IgM antibodies to *Chlamydia trachomatis* in serum or plasma.

For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 46 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Chlamydia spp. (*trachomatis*, *psittaci*, *pneumoniae*, *pecorum* and others) are obligate intracellular parasites and one of the most common human and animal pathogens. E.g., *C. trachomatis* has been recognized as a frequent cause of sexually transmitted diseases, which may lead to serious deterioration of reproductive function both in men and women (chronic pelvic inflammation and infertility). This species may also cause conjunctivitis and keratitis in adults and *ophthalmia neonatorum* in children born to an infected mother. *C. psittaci* is transmitted from birds to humans and causes an atypical pneumonia (ornithosis); *C. pneumoniae* is the causative agents of both acute pneumonia and chronic bronchitis. Due to intracellular 'depot' of the microorganism, *Chlamydia* infections may be asymptomatic and recurrent.

Detection of specific antibodies to *Chlamydia spp.* is a valuable tool to detect infection and to monitor its course and treatment.

IgM antibodies to *Chlamydia* are the first to appear – they become detectable within 1-2 weeks after infection and indicate a primary infection.

IgA antibodies to *Chlamydia* appear later on (usually – on day 10-14) and may be found both in serum and on mucosal surfaces (secretory IgA) and in high levels are indicative of chronic infection. IgA antibodies to *Chlamydia* react well to treatment and are used to monitor its effectiveness.

IgG antibodies to *Chlamydia* appear usually on week 3 after infection and indicate the level of immune response to the infection; their low levels may also be anamnestic and indicative of an earlier infection.

When primary *Chlamydia* infection is suspected, the detection of IgM is highly diagnostic. In reinfection, IgM level may be rarely detected while IgG and IgA levels rise quickly, often in one to two weeks.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated with murine monoclonal antibodies to human IgM. Total IgM from the specimen binds to anti-IgM fixed on the microwell surface. Unbound material is removed by washing procedure. *Chlamydia trachomatis*, labelled with peroxidase enzyme, is then added into the microwells and binds to anti-*Chlamydia trachomatis* IgM if present in total IgM fixed. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Positivity index (PI, %) is calculated by the formula (see Calculations). Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

4.1. For professional use only.

4.2. This kit is intended for *in vitro* diagnostic use only.

4.3. INFECTION HAZARD: There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

4.4. Avoid contact with stop solution containing 5.0% H₂SO₄. It may cause skin irritation and burns.

4.5. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

4.6. Do not use the kit beyond the expiration date.

4.7. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

4.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

4.9. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

4.10. Do not mix reagents from different lots.

4.11. Replace caps on reagents immediately. Do not swap caps.

4.12. Do not pipette reagents by mouth.

4.13. Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

4.14. Material Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

4.15. Material Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/ diluted components
1	<i>Chlamydia</i> IgM EIA strips, 8x12 wells	1	pcs		until exp.date
2	CONTROL - CONTROL + Control sera (0.5 ml and 0.2 ml, resp.)	2	pcs	colourless; red	2 months
3	CONJUGATE concentrate, 1.4 ml	1	pcs	blue	Concentrate - until exp.date Diluted - 1 day at 2-8 °C
4	DILUTION buffer, 14 ml	1	pcs	blue	until exp.date
5	DILUTION buffer 14 ml	1	pcs	purple	until exp.date
6	SUBSTRATE solution, 14 ml	1	pcs	colourless	until exp.date
7	WASHING solution concentrate 26X, 22 ml	1	pcs	colourless	Concentrate - until exp.date Diluted was- hing solution - 45 days at 2-8 °C or 15 days at RT
8	STOP solution, 14 ml	1	pcs	colourless	until exp.date
9	PLATE SEALING tape	2	pcs		N/A
10	INSTRUCTION Chlamydia IgM EIA	1	pcs		N/A
11	QC data sheet Chlamydia IgM EIA	1	pcs		N/A

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 90–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 10–250 µl;
- Dry thermostat for 37 °C ± 2 °C
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0-3.0.

5.3. Storage and stability of the Kit

Store the whole kit at 2 to 8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED) to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at 2-8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

7. TEST PROCEDURE**7.1. Reagent preparation**

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18 to +25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 26X by 26 dilutions in distilled water.

7.2. Procedural Note:

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 4 wells for control samples CONTROL - CONTROL + (3 and 1 wells resp.) and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Pipet 90 µl of EIA Sample buffer into each well.
3	Pipet 10 µl of control samples CONTROL - CONTROL + and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
4	Incubate 30 minutes at 37 °C.
5	Prepare washing solution by 26X dilution of washing solution concentrate BUF WASH 26X with distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
6	Prepare working conjugate solution by dilution of conjugate concentrate 11 fold by conjugate dilution buffer. ATTENTION: working conjugate solution is unstable and should not be stored! Prepare the volume required for actual assay run.
7	Dispense 100 µl of working conjugate solution into the wells. Cover the wells by plate adhesive tape.
8	Incubate 30 minutes at +18...+25°C.
9	Wash the strips 5 times.
10	Dispense 100 µl of SUBS TMB into the wells.
11	Incubate 10-20 minutes at +18...+25°C.
12	Dispense 100 µl of STOP into the wells.
13	Measure OD (optical density) at 450 nm.
14	Set photometer blank on air.

8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications:

For the assay to be valid, the following requirements should be met:

1. OD450 for CONTROL+ should be ≥ 1.5 AU.
2. OD450 for CONTROL- should not be more than 0.15 AU for all replicates.

9. CALCULATION OF RESULTS

1. Calculate the mean absorbance values (OD450) for CONTROL- in triplicates and each pair of samples.

2. Calculate the cut-off value: (mean OD450 for CONTROL-) + 0.2.

3. Calculate Positivity Index (PI) for each sample:

$PI = \text{mean OD450}(\text{sample}) / \text{Cut-off}$.

10. EXPECTED VALUES

If PI value is greater than 1.1, the result is POSITIVE.

If PI value is less than 0.9, the result is NEGATIVE.

If PI value is between 0.9 and 1.1, the result is EQUIVOCAL. Such samples should be retested. If the result is equivocal again, a new sample should be obtained 2-4 weeks later and tested again. If the result remains equivocal, the sample should be considered negative. Therapeutic consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutic measures. NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

11. PERFORMANCE CHARACTERISTICS

11.1. Specificity

Specificity of the test was evaluated on 94 serum specimens found negative. All tested specimens were found negative. Based on these data, specificity of the test is 100%.

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Russian Diagnostic
Manufacturers Association

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Russian Association
of Medical Laboratory
Diagnosticians

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