## Evolut™ PRO+ System

Evolut™ PRO+ Transcatheter Aortic Valve • Транскатетърна аортна клапа Evolut™ PRO+ • Katetrizačně implantovaná aortální chlopeč Evolut™ PRO+ • Evolut™ PRO+ transkateter-aortaklap • Evolut™ PRO+ Transkatheter-Aortenklappe • Διακαθετηριακή αορτική βαλβίδα Evolut™ PRO+ • Válvula aórtica transcáteter Evolut™ PRO+ • Transkateteraalne aordiklapp Evolut™ PRO+ • Katetrilla implantoitava Evolut™ PRO+ -aorttaläppä • Valve aortique transcathéter Evolut™ PRO+ • Transkateterski aortni zalistak Evolut™ PRO+ • Evolut™ PRO+ transzkátéteres aortabillentyű • Katup Aorta Transkateter Evolut™ PRO+ • Valvola aortica transcátetera Evolut™ PRO+ • Evolut™ PRO+ транскатетерлік қолқа клапаны • „Evolut™ PRO+“ transkateterinis aortos vožtuvas • Evolut™ PRO+ transkatetra aortas vārstulis • Транскатетерски аортен залисток Evolut™ PRO+ • Evolut™ PRO+-transkatheter-aortaklep • Evolut™ PRO+ transkateter aortaklaff • Zastawka aortalna do implantacji przezcewnikowej Evolut™ PRO+ • Válvula aórtica transcáteter Evolut™ PRO+ • Válvula aórtica transcáteter Evolut™ PRO+ • Valvā aorticā transcateter Evolut™ PRO+ • Аортальный клапан для транскатетерной установки Evolut™ PRO+ • Transkatétrová náhrada aortálnej chlopne Evolut™ PRO+ • Transkatetrská aortna zaklopka Evolut™ PRO+ • Evolut™ PRO+ transkateterski aortni zalistak • Evolut™ PRO+ kateterburen aortaklaff • Evolut™ PRO+ Transkateter Aortik Kapak • Аортальный клапан для транскатетерной установки Evolut™ PRO+ • Van động mạch chủ qua đường ống thông Evolut™ PRO+

Evolut™ PRO+ Delivery Catheter System • Система на доставящ катетър Evolut™ PRO+ • Systém zaváděcího katétru Evolut™ PRO+ • Evolut™ PRO+ fremføringskatetersystem • Evolut™ PRO+ Kathetereinführsystem • Σύστημα καθετήρα παράδοσης Evolut™ PRO+ • Sistema de catéter de liberación Evolut™ PRO+ • Sisestuskateetri süsteem Evolut™ PRO+ • Evolut™ PRO+ -sisäänvientikatetrijärjestelmä • Système de pose Evolut™ PRO+ • Kateterski sustav za uvodenje Evolut™ PRO+ • Evolut™ PRO+ bevezetőkatéter-rendszer • Sistem Kateter Penyaluran Evolut™ PRO+ • Sistema del catetere di rilascio Evolut™ PRO+ • Evolut™ PRO+ жеткізу катетері жүйесі • „Evolut™ PRO+“ įvedimo kateterio sistema • Evolut™ PRO+ ievadkatetra sistēma • Катетерски систем за испорачување Evolut™ PRO+ • Evolut™ PRO+-plaatsingskathetersysteem • Evolut™ PRO+-innføringskatetersystem • System ceownika wprowadzającego Evolut™ PRO+ • Sistema de cateter introdutor Evolut™ PRO+ • Sistema de cateter introdutor Evolut™ PRO+ • Sistem de poziționare cu cateter Evolut™ PRO+ • Система катетерной доставки Evolut™ PRO+ • Systém aplikačného katétra Evolut™ PRO+ • Katetrski sistem za dovajanje Evolut™ PRO+ • Evolut™ PRO+ sistem katetera za plasiranje • Evolut™ PRO+ införingskatetersystem • Evolut™ PRO+ Taşima Kateteri Sistemi • Катетерна система доставки Evolut™ PRO+ • Hệ thống ống thông chuyển van Evolut™ PRO+

Evolut™ PRO+ Loading System • Система за зареждане Evolut™ PRO+ • Zasouvací systém Evolut™ PRO+ • Evolut™ PRO+ indføringsssystem • Evolut™ PRO+ Ladesystem • Σύστημα φόρτωσης Evolut™ PRO+ • Sistema de montaje Evolut™ PRO+ • Laadimissüsteem Evolut™ PRO+ • Evolut™ PRO+ -latausjärjestelmä • Système de chargement Evolut™ PRO+ • Sustav za potiskivanje Evolut™ PRO+ • Evolut™ PRO+ betöltőrendszer • Sistem Pemuatan Evolut™ PRO+ • Sistema di caricamento Evolut™ PRO+ • Evolut™ PRO+ жүктеу жүйесі • „Evolut™ PRO+“ įstatymo sistema • Evolut™ PRO+ ielādes sistēma • Систем за веведување Evolut™ PRO+ • Evolut™ PRO+-laadsysteem • Evolut™ PRO+-ladesystem • System ładowania Evolut™ PRO+ • Sistema de carregamento Evolut™ PRO+ • Sistema de carregamento Evolut™ PRO+ • Sistem de încărcare Evolut™ PRO+ • Приспособление для загрузки клапана Evolut™ PRO+ • Zavádzací systém Evolut™ PRO+ • Sistem za vstavljanje Evolut™ PRO+ • Evolut™ PRO+ sistem za ubacivanje • Evolut™ PRO+ inmatningssystem • Evolut™ PRO+ Yükleme Sistemi • Система завантаження Evolut™ PRO+ • Hệ thống nạp Evolut™ PRO+

Instructions for Use • Инструкции за употреба • Návod k použití • Brugsanvisning • Gebrauchsanweisung • Οδηγίες χρήσης • Instrucciones de uso • Kasutusjuhend • Käyttöohjeet • Mode d'emploi • Upute za upotrebu • Használati útmutató • Petunjuk Penggunaan • Istruzioni per l'uso • Қолдану жөнiндегi нұсқаулар • Naudojimo instrukcijos • Lietošanas pamācība • Упатство за употреба • Gebruiksaanwijzing • Bruksanvisning • Instrukcja użytkowania • Instruções de utilização • Instruções de utilização • Instrucțiuni de utilizare • Инструкция по применению • Pokyny na používanie • Navodila za uporabo • Uputstva za upotrebu • Bruksanvisning • Kullanım Talimatları • Інструкція з використання • Hướng dẫn sử dụng

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
- EN 1 Catheter tip guide tube  
2 Inflow cone  
3 Backplate

- 4 Outflow cone  
5 Capsule guide tube  
6 Locking collar

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## 1 Device description

The Medtronic Evolut PRO+ system is a recapturable transcatheter aortic valve replacement system, which includes the Evolut PRO+ transcatheter aortic valve, the delivery catheter system, and the loading system.

### 1.1 Evolut PRO+ transcatheter aortic valve (bioprosthesis)

The bioprosthesis is designed to replace the native or surgical bioprosthetic aortic heart valve without open heart surgery and without concomitant surgical removal of the failed valve. The support frame is manufactured from Nitinol, which has multilevel, self-expanding properties and is radiopaque. The bioprosthesis is manufactured by suturing valve leaflets and an inner skirt from porcine pericardium into a tri-leaflet configuration. The bioprosthesis has a porcine pericardial tissue outer skirt (wrap), which is 1.5 cells in height and is sutured to the inflow section of the bioprosthesis.

The bioprosthesis is processed with alpha-amino oleic acid (AOA), which is a compound derived from oleic acid, a naturally occurring long-chain fatty acid. AOA is an antimineralization treatment shown to reduce both early and late valvular calcification. The bioprosthesis is available for a range of aortic annulus diameters (refer to *Table 1*).

**Table 1.** Patient anatomical criteria

Bioprosthesis model	Size	Aortic annulus diameter	Aortic annulus perimeter ( $\pi \times$ aortic annulus diameter)
EVPROPLUS-23	23 mm	18 mm to 20 mm	56.5 mm to 62.8 mm
EVPROPLUS-26	26 mm	20 mm to 23 mm	62.8 mm to 72.3 mm
EVPROPLUS-29	29 mm	23 mm to 26 mm	72.3 mm to 81.7 mm
EVPROPLUS-34	34 mm	26 mm to 30 mm	81.7 mm to 94.2 mm

**Table 2.** Patient access criteria

Bioprosthesis model	Delivery catheter system model	Capsule OD	Arterial access	Subclavian access with an IMA graft	Direct aortic access
EVPROPLUS-23, EVPROPLUS-26, EVPROPLUS-29	D-EVPROP23-29	18 Fr (6 mm)	$\geq 5$ mm	$\geq 5.5$ mm	$\geq 60$ mm from basal plane
EVPROPLUS-34	D-EVPROP34	22 Fr (7.33 mm)	$\geq 6$ mm	$\geq 6.5$ mm	

### 1.2 Delivery catheter system (catheter)

The catheter comes in different models. Refer to *Table 3* for system compatibility. Refer to *Figure 2* and *Figure 3* for catheter components.

The catheter facilitates the placement of the bioprosthesis within the annulus of the aortic valve. The catheter assembly is flexible and compatible with a 0.035 in (0.889 mm) guidewire. The distal (deployment) end of the system features an atraumatic, radiopaque catheter tip and a capsule that covers and maintains the bioprosthesis in a crimped position. The capsule includes a distal flare to enable the bioprosthesis to be partially or fully recaptured after partial deployment. A stability layer is fixed at the handle and extends down the outside of the catheter shaft. It provides a barrier between the retractable catheter and the introducer sheath and vessel walls, thus enabling the catheter to retract freely. The Evolut PRO+ inline sheath is assembled over the stability layer, which functions as a hemostatic introducer sheath and minimizes the access site size to the capsule diameter. The 23-29 mm catheter model is compatible with sheaths that can accommodate an 18 Fr (6.00 mm) device. The 34 mm catheter model is compatible with sheaths that can accommodate a 22 Fr (7.33 mm) device.

The delivery catheter system consists of a catheter with an integrated handle to provide the user with accurate and controlled deployment. The handle is on the proximal end of the catheter and is used to load, deploy, recapture, and reposition the bioprosthesis. The handle features a gray front grip used to stabilize the system. The deployment knob turns to deploy the bioprosthesis precisely. Arrows on the deployment knob indicate the direction of rotation required to deploy the bioprosthesis. If desired, the deployment knob can be turned in the opposite direction to partially or fully recapture the bioprosthesis if the radiopaque capsule marker band has not yet reached the distal end of the radiopaque paddle attachment. Once the radiopaque capsule marker band reaches the distal end of the radiopaque paddle attachment, it is at the point of no recapture. The deployment knob also features a trigger, which can be engaged to make macro adjustments to the capsule position. A blue hand rest connects to the deployment knob. The end of the handle features a tip-retrieval mechanism, which can be used to withdraw the catheter tip to meet the capsule after the device has been fully deployed.

The catheter packaging contains an integrated loading bath and a removable tray with 3 rinsing bowls for loading and rinsing the bioprosthesis. The integrated loading bath features a mirror, which aids in accurate placement of the bioprosthesis frame paddles during loading. In addition to these features, the device packaging is swiveled and secured to facilitate the bioprosthesis loading procedure (*Figure 4* and *Figure 5*).

### 1.3 Loading system

The loading system compresses the bioprosthesis into the catheter. The loading system comes in different models.

Refer to *Table 3* for system compatibility. Refer to *Figure 6* for components.

**Table 3.** System compatibility

Bioprosthesis model	Compatible loading system models	Compatible catheter models
EVPROPLUS-23	L-EVPROP23-29	D-EVPROP23-29
EVPROPLUS-26		
EVPROPLUS-29		
EVPROPLUS-34	L-EVPROP34	D-EVPROP34

**Note:** Verify the correct delivery system and loading system are used for the selected valve. Valves are designed for compatibility with specific delivery and loading systems.

## 2 Warnings and precautions

Carefully read all warnings, precautions, and instructions for use for all components of the system before use. Failure to read and follow all instructions or failure to observe all stated warnings could cause serious injury or death to the patient.

### 2.1 Warnings

- This procedure should only be performed where emergency aortic valve surgery can be performed promptly.
- Carefully verify that the relevant patient anatomical parameters are within the specifications set forth in *Section 1.1*.
- The heart team should consider the need for revascularization when selecting a treatment option for aortic valve replacement (AVR).
- Do not freeze the bioprosthesis. Freezing causes damage to the bioprosthesis. Check the freeze indicator located inside the labeled carton. If the product has been exposed to a freeze-thaw condition, the indicator vial will break, which causes dye to escape and stain the paper backing. Do not use the bioprosthesis if this or any other sign of freezing is observed.
- Package contents are supplied sterile. Do not use if package is damaged.
- The outside of the bioprosthesis container and the outer surfaces of the catheter pouch and the loading system outer pouch are nonsterile and must not be placed in the sterile field.
- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- If a bioprosthesis and catheter have been removed from a patient, dispose of both the bioprosthesis and catheter; do not attempt to reuse either component. Both the bioprosthesis and catheter must be replaced with new sterile components.
- Package contents should not be opened until implantation and sizing are certain.
- Do not handle or manipulate the bioprosthesis with sharp or pointed objects.
- Do not use products if the use-by date printed on the product packaging label has elapsed.
- Store the bioprosthesis at room temperature.
- Carefully remove the catheter from the packaging to prevent it from kinking.
- Perform the **Bioprosthesis rinsing procedure** (*Section 7.3*) before loading the bioprosthesis into the catheter.
- After the bioprosthesis has been loaded into the catheter, maintain the distal end of the catheter with the bioprosthesis in cold (0°C to 8°C [32°F to 46°F]) sterile saline until implantation.
- Use the deployment knob to deploy and recapture the bioprosthesis. Do not use the trigger for deploying or recapturing because it could cause inaccurate placement of the bioprosthesis.
- System failure could occur if an incorrect combination of devices is used. Refer to the system compatibility table in *Section 1.3*.
- Mechanical failure of the delivery catheter system and/or accessories may result in patient complications.
- Calcific degeneration could cause accelerated structural deterioration of the valve in the following individuals:
  - Children, adolescents, or young adults
  - Patients with altered calcium metabolism (for example, chronic renal failure, or hyperparathyroidism)
- The safety and effectiveness of the Evolut PRO+ bioprosthesis implanted within a failed pre-existing transcatheter bioprosthesis have not been demonstrated.

### 2.2 Patient selection considerations

The safety and effectiveness of the Medtronic Evolut PRO+ system have not been evaluated in the pediatric population.

The safety and effectiveness of the bioprosthesis for aortic valve replacement have not been evaluated in the following patient populations:

- Patients with a pre-existing prosthetic heart valve with a rigid support structure in either the mitral or pulmonic position if either the pre-existing prosthetic heart valve could affect the implantation or function of the bioprosthesis or the implantation of the bioprosthesis could affect the function of the pre-existing prosthetic heart valve
- Patients with liver failure (Child-Pugh Class C)
- Patients with cardiogenic shock manifested by low cardiac output, vasopressor dependence, or mechanical hemodynamic support
- Patients who are pregnant or breastfeeding
- Patients with a congenital unicuspid valve
- Patients with prohibitive left ventricular outflow tract calcification
- Patients with sinus of Valsalva anatomy that would prevent adequate coronary perfusion
- Patients with significant aortopathy requiring ascending aortic replacement
- Patients who do not meet the criteria for symptomatic severe native aortic stenosis as defined by the following:
  - Aortic valve area of  $\leq 1.0 \text{ cm}^2$  (or aortic valve area index of  $\leq 0.6 \text{ cm}^2/\text{m}^2$ ) by the continuity equation, AND mean gradient  $\geq 40 \text{ mmHg}$  or maximal aortic valve velocity  $\geq 4.0 \text{ m/s}$  by resting echocardiogram, low-dose dobutamine, or exercise stress echocardiography
- Patients who do not meet the criteria for asymptomatic severe native aortic stenosis as defined by the following:
  - Very severe aortic stenosis with an aortic valve area of  $\leq 1.0 \text{ cm}^2$  (or aortic valve area index of  $\leq 0.6 \text{ cm}^2/\text{m}^2$ ), AND maximal aortic velocity  $\geq 5.0 \text{ m/s}$ , or mean gradient  $\geq 60 \text{ mmHg}$  by transthoracic echocardiography at rest, OR
  - Aortic valve area of  $\leq 1.0 \text{ cm}^2$  (or aortic valve area index of  $\leq 0.6 \text{ cm}^2/\text{m}^2$ ), AND a mean gradient  $\geq 40 \text{ mmHg}$ , or maximal aortic valve velocity  $\geq 4.0 \text{ m/s}$  by transthoracic echocardiography at rest, AND an exercise tolerance test that demonstrates a limited exercise capacity, abnormal BP response, or arrhythmia, OR
  - Aortic valve area of  $\leq 1.0 \text{ cm}^2$  (or aortic valve area index of  $\leq 0.6 \text{ cm}^2/\text{m}^2$ ), AND mean gradient  $\geq 40 \text{ mmHg}$ , or maximal aortic valve velocity  $\geq 4.0 \text{ m/s}$  by transthoracic echocardiography at rest, AND a left ventricular ejection fraction  $< 50\%$

The safety and effectiveness of the bioprosthesis for aortic valve replacement have not been evaluated in patient populations presenting with the following conditions:

- Blood dyscrasias as defined: leukopenia ( $\text{WBC} < 1000 \text{ cells}/\text{mm}^3$ ), thrombocytopenia (platelet count  $< 50\,000 \text{ cells}/\text{mm}^3$ ), history of bleeding diathesis or coagulopathy, or hypercoagulable states
- Mixed native aortic valve disease (aortic stenosis and aortic regurgitation with predominant aortic regurgitation [3–4+])
- Moderate to severe (3–4+) or severe (4+) mitral or severe (4+) tricuspid regurgitation
- Hypertrophic obstructive cardiomyopathy
- New or untreated echocardiographic evidence of intracardiac mass, thrombus, or vegetation
- Moderate to severe mitral stenosis
- Severe ventricular dysfunction with left ventricular ejection fraction (LVEF)  $< 20\%$
- Symptomatic carotid or vertebral artery disease

Current heart valve replacement guidelines recommend bioprosthetic valves for patients at 65 years of age or above. Consider the applicable clinical guidelines for valve replacement when selecting the appropriate valve type for each patient.

### 2.3 Implantation precautions

- Before catheter insertion, remove the loading stylet.
- Before catheter insertion, the heart team determines whether to perform balloon aortic valvuloplasty (BAV) of the native aortic valve to facilitate placement of the bioprosthesis.
- While the catheter is in the patient, ensure the guidewire is extending from the proximal end of the catheter. Do not remove the guidewire from the catheter while the catheter is inserted in the patient.
- There will be some resistance when the catheter is advanced through the vasculature. If there is a significant increase in resistance, stop advancement and investigate the cause of the resistance (for example, magnify the area of resistance) before proceeding. Do not force passage. Forcing passage could increase the risk of vascular complications (for example, vessel dissection or rupture).
- For direct aortic access procedures, use caution when advancing and maintaining the introducer sheath position to avoid tissue and vascular injury and unacceptable access site bleeding.

## 2.4 Repositioning precautions

- If the radiopaque capsule marker band has not yet reached the distal end of the radiopaque paddle attachment, the bioprosthesis can be recaptured or repositioned. During deployment, the deployment knob provides a tactile indication as a notification before the point of no recapture.
- Once the radiopaque capsule marker band reaches the distal end of the radiopaque paddle attachment (point of no recapture), retrieval of the bioprosthesis from the patient (for example, use of the catheter) is not recommended. Retrieval after the point of no recapture may cause mechanical failure of the delivery catheter system, aortic root damage, coronary artery damage, myocardial damage, vascular complications, prosthetic valve dysfunction (including device malposition), embolization, stroke, and/or emergent surgery.
- During deployment, the bioprosthesis can be advanced or withdrawn as long as annular contact has not been made. Once annular contact is made, the bioprosthesis cannot be advanced in the retrograde direction; recapture until the bioprosthesis is free from annular contact, and then reposition in the retrograde direction. If necessary, and the radiopaque capsule marker band has not yet reached the distal end of the radiopaque paddle attachment, the bioprosthesis can be withdrawn (repositioned) in the antegrade direction. However, use caution when moving the bioprosthesis in the antegrade direction.

**Caution:** Use the handle of the delivery system to reposition the bioprosthesis. **Do not** use the outer catheter shaft.

- Physicians should use judgment when considering repositioning a fully deployed bioprosthesis (for example, using a snare, balloon, and/or forceps). Repositioning the bioprosthesis is not recommended except in cases where imminent serious harm or death is possible (for example, coronary occlusion). Repositioning of a deployed valve may cause aortic root damage, coronary artery damage, myocardial damage, vascular complications, prosthetic valve dysfunction (including device malposition), embolization, stroke, and/or emergent surgery.
- Do not attempt to retrieve or to recapture a bioprosthesis if any one of the outflow struts is protruding from the capsule. If any one of the outflow struts has deployed from the capsule, the bioprosthesis must be released from the catheter before the catheter can be withdrawn.

## 2.5 Postimplantation precautions

- Ensure the capsule is closed before catheter removal.
- When using a separate introducer sheath, if increased resistance is encountered when removing the catheter through the introducer sheath, do not force passage. Increased resistance may indicate a problem and forced passage may result in damage to the device and/or harm to the patient. If the cause of resistance cannot be determined or corrected, remove the catheter and introducer sheath as a single unit over the guidewire, and inspect the catheter and confirm that it is complete.
- Exercise care when crossing an implanted bioprosthesis with the catheter or other devices.
- The safety and performance of the Evolut PRO+ bioprosthesis implanted within a transcatheter bioprosthesis have not been demonstrated. However, in the event that the Evolut PRO+ bioprosthesis must be implanted within a transcatheter bioprosthesis to improve valve function, valve size and patient anatomy must be considered before implantation of the Evolut PRO+ bioprosthesis to ensure patient safety (for example, to avoid coronary obstruction).
- Clinical long-term durability has not been established for the bioprosthesis. Evaluate bioprosthesis performance as needed during patient follow-up.

## 2.6 Post-implant balloon dilatation considerations

If valve function or sealing is impaired due to excessive calcification or incomplete expansion, a post-implant balloon dilatation (PID) of the bioprosthesis may improve valve function and sealing. If the heart team determines that balloon dilatation is appropriate, consider all of the following factors when selecting the dilatation parameters to ensure patient safety:

- Balloon model
- Balloon size
- Balloon position
- Inflation pressure
- Patient anatomy

Two primary factors must be considered when selecting a maximum balloon diameter for post-implant balloon dilatation:

- To mitigate trauma to the annulus
  - A compliant or semi-compliant balloon (for example, B. Braun Z-Med I™\* / Z-Med II™\*, InterValve V8™\*) should not exceed the diameter of the native aortic annulus. For TAV in SAV, the balloon should not exceed the inner diameter of the surgical bioprosthetic valve.
  - A non-compliant balloon (for example, Bard TRUE™\* Dilatation) should be at least 1 mm smaller than the diameter of the native aortic annulus. For TAV in SAV, the balloon should be at least 1 mm smaller than the inner diameter of the surgical bioprosthetic valve.

- To mitigate trauma to the Evolut TAV bioprosthetic leaflets
  - The maximum balloon size chosen for dilatation using a compliant or semi-compliant balloon should not exceed the TAV waist diameter beyond the level set forth in *Table 4* with an applied inflation pressure of no greater than 2 atm.
  - The maximum balloon size chosen for dilatation using a non-compliant balloon should not exceed 1 mm more than the TAV waist diameter with an applied inflation pressure of no greater than 2 atm (see *Table 4*).

**Table 4.** Post-implant balloon dilatation sizing

Evolut PRO+ size	23 mm			26 mm				29 mm				34 mm				
Native annulus (SAV inner) diameter (in mm)	18	19	20	20	21	22	23	23	24	25	26	26	27	28	29	30
TAV waist diameter (in mm)	20	20	20	22	22	22	22	23	23	23	23	24	24	24	24	24
Maximum balloon diameter (in mm) for compliant and semi-compliant balloons at 2 atm	18	19	20	20	21	22	23	23	24	25	26	26	27	28	28	28
Maximum balloon diameter (in mm) for non-compliant balloons at 2 atm	17	18	19	19	20	21	22	22	23	24	24	25	25	25	25	25

**Caution:** Overexpansion of the narrowest portion (waist) of the Evolut PRO+ TAV beyond the levels set forth in *Table 4* has been demonstrated through bench data to cause damage to the bioprosthetic leaflets. Complaints of damage to the bioprosthetic leaflets during post-implant balloon dilatation have been reported in some clinical cases, resulting in moderate to severe aortic insufficiency, which may be detected acutely or during follow-up.

It is important to note that the mechanical compliance properties of the selected balloon influence the dilatation dynamics.

Balloons should not be inflated beyond 2 atm of applied pressure.

Compliant and semi-compliant (softer) balloons will more readily conform to the hourglass profile of the TAV bioprosthesis at lower pressures, but must be inflated at pressures that preserve the hourglass profile of the TAV.

Conversely, non-compliant (stiffer) balloons will achieve the nominal diameter during inflation irrespective of the underlying annulus or TAV resistance and should be downsized (see *Table 4*).

For additional instructions on the use of balloon catheter devices refer to the specific balloon catheter manufacturer’s labeling.

In the event that larger balloon diameters than those listed in *Table 4* are required to expand the Evolut PRO+ TAV due to clinically important residual aortic regurgitation or stenosis, using “bailout” intraventricular balloon positioning when performing PID avoids expansion of the narrowest portion (waist) of the Evolut PRO+ TAV. This can mitigate the risk of leaflet damage. Dilatation with intraventricular balloon positioning should be performed with caution in the setting of a smaller ventricle cavity, presence of LVOT calcification, or wire positioning that interferes with mitral valve function, in order to avoid any unintended balloon interaction with anatomy. The balloon’s length and diameter, along with the individual patient anatomy, must be considered. Care should also be taken not to exceed the annular diameters when performing PID with intraventricular balloon positioning (see *Table 4*).

In the event that a bailout PID with intraventricular balloon positioning is performed, the nominal diameter of the balloon should not exceed the annular diameter when using compliant or semi-compliant balloons; the nominal diameter of the balloon should be at least 1 mm smaller than the annular diameter when using non-compliant balloons.

### 3 Potential adverse events

Adverse events that may be associated with the use of the Evolut PRO+ system include, but may not be limited to, the following:

- Death
- Myocardial infarction, cardiac arrest, cardiogenic shock, cardiac tamponade
- Coronary occlusion, obstruction, or vessel spasm (including acute coronary closure)
- Cardiovascular injury (including rupture, perforation, tissue erosion, or dissection of vessels, ascending aorta trauma, ventricle, myocardium, or valvular structures that may require intervention)
- Emergent surgical or transcatheter intervention (for example, coronary artery bypass, heart valve replacement, valve explant, percutaneous coronary intervention [PCI], balloon valvuloplasty)
- Prosthetic valve dysfunction (regurgitation or stenosis) due to fracture; bending (out-of-round configuration) of the valve frame; underexpansion of the valve frame; calcification; pannus; leaflet wear, tear, prolapse, or retraction; poor valve coaptation; suture breaks or disruption; leaks; mal-sizing (prosthesis-patient mismatch); malposition (either too high or too low)/malplacement
- Prosthetic valve migration/embolization
- Prosthetic valve endocarditis
- Prosthetic valve thrombosis
- Delivery catheter system malfunction resulting in the need for additional re-crossing of the aortic valve and prolonged procedural time

- Delivery catheter system component migration/embolization
- Stroke (ischemic or hemorrhagic), transient ischemic attack (TIA), or other neurological deficits
- Individual organ (for example, cardiac, respiratory, renal [including acute kidney failure]) or multi-organ insufficiency or failure
- Major or minor bleeding that may require transfusion or intervention (including life-threatening or disabling bleeding)
- Vascular access-related complications (for example, dissection, perforation, pain, bleeding, hematoma, pseudoaneurysm, irreversible nerve injury, compartment syndrome, arteriovenous fistula, stenosis)
- Mitral valve regurgitation or injury
- Conduction system disturbances (for example, atrioventricular node block, left-bundle branch block, asystole), which may require a permanent pacemaker
- Infection (including septicemia)
- Hypotension or hypertension
- Hemolysis
- Peripheral ischemia
- Bowel ischemia

General surgical risks applicable to transcatheter aortic valve implantation:

- Abnormal lab values (including electrolyte imbalance)
- Allergic reaction to antiplatelet agents, contrast medium, or anesthesia
- Exposure to radiation through fluoroscopy and angiography
- Permanent disability

## 4 Patient selection and treatment

### 4.1 Indications

The Evolut PRO+ system is indicated for patients presenting with severe native aortic valve stenosis. For patients presenting with severe native bicuspid aortic valve stenosis, the Evolut PRO+ system is indicated for patients who are at intermediate or greater risk for surgical aortic valve replacement (AVR) where intermediate risk is defined as Society of Thoracic Surgeons operative risk score  $\geq 4\%$  or documented heart team agreement of risk for AVR due to frailty or comorbidities. For patients presenting at low risk for AVR ( $< 4\%$ ), the system is indicated for patients  $\geq 70$  years of age with an LVEF  $> 30\%$ .

The Evolut PRO+ system is also indicated for patients with a stenosed, insufficient, or combined surgical bioprosthetic valve failure necessitating valve replacement who are at high or greater risk for surgical aortic valve replacement (AVR) where high risk is defined as Society of Thoracic Surgeons operative risk score  $\geq 8\%$  or documented heart team agreement of risk for AVR due to frailty or comorbidities.

Patients must present with anatomical dimensions as described in *Section 1.1*.

### 4.2 Contraindications

The Evolut PRO+ system is contraindicated for patients presenting with any of the following conditions:

- Known hypersensitivity or contraindication to aspirin, heparin (HIT/HITTS) and bivalirudin, ticlopidine, clopidogrel, Nitinol (Titanium or Nickel), or sensitivity to contrast media, which cannot be adequately premedicated
- Preexisting mechanical heart valve in aortic position
- Ongoing sepsis, including active endocarditis

### 4.3 Access route considerations

Patients must present with transarterial access vessels with diameters that are  $\geq 5.0$  mm when using Model D-EVPROP23-29 or  $\geq 6.0$  mm when using Model D-EVPROP34, or patients must present with an ascending aortic (direct aortic) access site  $\geq 60$  mm from the basal plane. Patients must also present with anatomical dimensions as described in *Section 1.1*.

The risks of long-term anticoagulation and/or antiplatelet therapy should be considered. Implantation of the bioprosthesis is not recommended in patients with aortic root angulation (the angle between the plane of the aortic valve annulus and the horizontal plane/vertebrae) of  $> 30^\circ$  for right subclavian/axillary access or  $> 70^\circ$  for femoral and left subclavian/axillary access.

For subclavian access, patients with a patent Left Internal Mammary Artery (LIMA) graft must present with access vessel diameters that are either  $\geq 5.5$  mm when using Model D-EVPROP23-29 or  $\geq 6.5$  mm when using Model D-EVPROP34. Use caution in patients with a preexisting patent Left Internal Mammary Artery/Right Internal Mammary Artery (LIMA/RIMA) graft. For direct aortic access, ensure the access site and trajectory are free of patent RIMA or a pre-existing patent RIMA graft.

For direct aortic and subclavian access procedures, care must be exercised when using the tip-retrieval mechanism to ensure adequate clearance to avoid advancement of the catheter tip through the bioprosthesis leaflets during device closure.

For direct aortic access procedures, use a separate introducer sheath; do not use the Evolut PRO+ inline sheath. Maintain the Evolut PRO+ inline sheath at the proximal end of the catheter throughout the procedure.

For transfemoral access, use caution in patients who present with multiplanar curvature of the aorta, acute angulation of the aortic arch, an ascending aortic aneurysm, or severe calcification in the aorta and/or vasculature. If  $\geq 2$  of these factors are present, consider an alternative access route to prevent vascular complications.

#### 4.4 Valve-in-valve considerations

Implanting the Evolut PRO+ bioprosthesis in a degenerated surgical aortic bioprosthesis (valve-in-valve) should generally be avoided in patients with any of the following conditions:

- The degenerated bioprosthesis presents with a significant concomitant perivalvular leak (between prosthesis and native annulus), is not securely fixed in the native annulus, or is not structurally intact (for example, wireform frame fracture).
- The degenerated bioprosthesis presents with a partially detached leaflet that in the aortic position may obstruct a coronary ostium.

#### 4.5 Bicuspid patient considerations

If the patient presents with a bicuspid aortic valve, the heart team should consider the age of the patient and the need for ascending aorta intervention when determining the appropriate treatment option for the patient.

### 5 How supplied

#### 5.1 Bioprosthesis

The bioprosthesis is sterilized with and stored in glutaraldehyde solution. The bioprosthesis is supplied sterile and nonpyrogenic in a jar with a screw cap. A freeze indicator is placed inside the labeled carton. If the freeze indicator has been activated, do not use the bioprosthesis.

#### 5.2 Catheter

The catheter is packaged in a single-pouch configuration. The catheter is sterile if the package is undamaged and unopened. The catheter is sterilized by ethylene oxide.

#### 5.3 Loading system

The loading system is packaged in a double-pouch configuration. The loading system is sterile if the packages are undamaged and unopened. The loading system is sterilized by ethylene oxide.

### 6 Additional materials (not included with the system)

- Sterile saline
- Sterile syringe, 10 mL (10 cm<sup>3</sup>)
- Extra-support guidewire, 0.035 in (0.889 mm)  $\times$  a minimum length of 240 cm
- Valvuloplasty balloon
- Standard procedure supplies

### 7 Instructions for use

**Caution:** Do not allow the bioprosthesis to dry. Maintain tissue moisture with irrigation or immersion.

#### 7.1 Inspection prior to use and swivel tray setup

1. Carefully inspect the package before opening.

**Caution:** Do not use after the Use By date or if the integrity of the sterile package has been compromised (for example, damaged package).

**Caution:** Do not use the bioprosthesis if the freeze indicator has been activated.

2. Remove the product from the protective package.
3. Visually check that the product is free of defects. Do not use if any defects are noted.
4. Remove the locking clip attached to the rinsing bowls.
5. Remove the rinsing bowls from the integrated loading bath.
6. Remove the locking clips that connect the distal and proximal trays.
7. Lift the tray connector from the distal tray, and swivel the distal tray 180° counterclockwise.

- Clip the tray tab on the distal tray to the tray tab holder on the proximal tray.
- Fill the integrated loading bath with cold, sterile saline (0°C to 8°C [32°F to 46°F]).

## 7.2 Preparation of the catheter and loading system

- Attach a 10 mL syringe filled with sterile saline to the capsule flush port on the proximal end of the handle. Leave the syringe in place until loading is completed.
- Carefully lift the distal end of the catheter to a near vertical orientation. To prevent kinking, do not bend the catheter severely.
- Open the capsule and expose the paddle attachment.  
**Note:** Use the deployment knob to open the capsule completely until the paddle attachment is fully exposed.
- With the capsule held vertically, flush the capsule flush port. Verify that no catheter leakage is observed during any of the flushing steps. If leakage is observed, use a new system.
- Submerge the capsule completely in the cold saline bath while flushing the capsule flush port. Continue flushing the capsule until it is completely submerged in the bath to prevent air from entering the catheter.  
**Note:** After the bioprosthesis has been loaded into the capsule, the capsule flush port can no longer be flushed.
- Secure a locking clip in the clip holder to angle the catheter tip into the integrated loading bath.
- Place the loading system components in the integrated loading bath.

## 7.3 Bioprosthesis rinsing procedure

- Fill each of the 3 rinsing bowls (provided within the packaging) with approximately 500 mL of fresh, sterile saline at ambient temperature (15°C to 25°C [59°F to 77°F]).  
**Caution:** Do not handle or manipulate the bioprosthesis with sharp or pointed objects. Use atraumatic forceps only.
- Confirm the integrity of the primary bioprosthesis container. Open the container and remove the bioprosthesis. Let any remaining solution drain from the bioprosthesis completely.  
**Note:** Retain the container with the original solution. It may be needed to store and return a rejected bioprosthesis.
- Compare the serial number on the container with the serial number on the tag attached to the bioprosthesis.  
**Caution:** If the serial numbers do not match, do not use the bioprosthesis.
- Carefully remove the serial number tag from the bioprosthesis and retain the tag.
- Immerse the entire bioprosthesis in a sterile rinsing bowl.
- Gently agitate the bioprosthesis by hand for 15 seconds to remove the glutaraldehyde from the bioprosthesis.
- Repeat *Step 5* and *Step 6* in one of the remaining rinsing bowls.
- Leave the bioprosthesis submerged in sterile saline in the third rinsing bowl until it is ready to be loaded.

## 7.4 Bioprosthesis loading procedure

Perform the bioprosthesis loading procedure while the distal end of the catheter is immersed in the integrated loading bath filled with cold, sterile saline (0°C to 8°C [32°F to 46°F]). The bioprosthesis should remain immersed in saline during the loading process to minimize the introduction of air into the loaded system.

**Note:** Confirm the loading system and catheter sizes are compatible with the bioprosthesis size (*Table 3*).

**Note:** Refer to *Figure 6* for loading system components.

**Caution:** Rapid capsule advancement can contribute to difficulties with loading the valve. Slowly advancing the capsule helps facilitate successful loading.

- Submerge and cool the bioprosthesis in the integrated loading bath filled with cold, sterile saline.
- Ensure that the capsule guide tube is fully open (unlocked) with the locking collar at the proximal end of the capsule guide tube.
- Advance the capsule guide tube over the catheter shaft toward the handle and across the catheter tip.
- Once the catheter tip has been crossed, fully advance the locking collar to the distal end of the capsule guide tube until it is closed (locked).
- Continue to advance the capsule guide tube over the catheter shaft towards the handle until it contacts the distal end of the capsule.  
**Caution:** Do not attempt to advance the capsule guide tube over the capsule; this will prevent the capsule flare from expanding fully and prevent proper loading.
- Ensure that the backplate has been inserted into the inflow cone and the exposed part of the backplate is facing up.
- Insert the inflow portion of the bioprosthesis frame into the inflow cone. Ensure that the bioprosthesis frame paddle marked with a "C" is facing up and that the paddles are aligned with the paddle attachment pockets.
- Secure the outflow cone onto the inflow cone until it locks.

9. Insert the catheter tip guide tube completely into the distal end of the inflow cone. Inspect the outflow struts of the bioprosthesis and, if needed, manually manipulate so that they are evenly spaced and the bioprosthesis frame paddles are approximately 180° apart.
10. Insert the distal catheter tip into the catheter tip guide tube.  
**Note:** Allow the loading tool to rest on the loading bath floor to ensure coaxial alignment with the catheter to assist in seating the bioprosthesis frame paddles within the paddle attachment pockets.
11. Retract the catheter tip guide tube to set the bioprosthesis frame paddles into the paddle attachment pockets.  
**Note:** If the bioprosthesis frame paddles do not seat properly within the paddle attachment pockets upon retracting the catheter tip guide tube, slightly manipulate the position of the loading tool until paddle seating is achieved.  
**Note:** If necessary, it is acceptable to manually compress the bioprosthesis frame paddles with fingertips to help seat the paddles within the paddle attachment pockets.  
**Note:** Ensure both bioprosthesis frame paddles are completely seated within the paddle attachment pockets before continuing to the next step.
12. Hold the loading tool stationary with one hand, and with the other hand manually advance the capsule guide tube so that the distal section covers the paddle attachment pockets and the top portion of the outflow struts. Use the mirror to ensure that both bioprosthesis frame paddles are positioned correctly in the paddle attachment pockets and the outflow struts are within the distal tip of the capsule guide tube.
13. Advance the capsule to cover the bioprosthesis frame paddles, pausing when the capsule covers the proximal half of the paddles to confirm the paddles are both still properly seated before advancing further. Use the mirror to ensure that both paddles are captured in the capsule.  
**Caution:** Do not advance the capsule over the bioprosthesis frame paddles unless they are fully seated in the center of the paddle attachment pockets. Advancing the capsule before the paddles are fully seated could damage the capsule and result in emboli.
14. Advance the capsule to capture the bioprosthesis outflow struts. Use the mirror to ensure that all bioprosthesis outflow struts are symmetrical and captured in the capsule.
15. Continue to advance the capsule until the distal end of the capsule guide tube covers the distal end of the commissure pad of the bioprosthesis. The capsule guide tube should completely cover the commissure pad.
16. Remove the backplate and the catheter tip guide tube from the outflow cone.
17. While holding the capsule guide tube stationary, advance the inflow cone to crimp the inflow portion of the bioprosthesis frame until the outflow cone contacts the capsule guide tube. During this step, the outflow cone contacts the locking collar component and moves the locking collar to the proximal end of the capsule guide tube.  
**Note:** The capsule guide tube will be in the unlocked configuration after this step.  
**Note:** Ensure the bioprosthesis frame axis is visually aligned (coaxial) with the inflow cone axis during the insertion of the bioprosthesis into the inflow cone. Complete the insertion of the bioprosthesis into the inflow cone in one uninterrupted movement.
18. Advance the capsule over the bioprosthesis until the capsule comes within 5 mm of the catheter tip.
19. Remove the capsule guide tube together with the outflow cone and inflow cone from the catheter.
20. Advance the capsule to close the gap between the capsule and catheter tip completely.  
**Caution:** Stop advancing the capsule once the gap to the catheter tip is closed. Advancing the capsule farther could damage the capsule.
21. Slightly rotate the deployment knob in the direction of the arrows to relieve stress. Ensure that the capsule does not separate from the catheter tip.  
**Note:** After the bioprosthesis has been loaded into the capsule, the capsule flush port can no longer be flushed.
22. Visually and tactilely inspect the capsule for a misloaded bioprosthesis. The capsule should be straight, smooth, and free of any bends, protrusions, or discolorations. If any of these conditions are felt or observed, the bioprosthesis is likely to be misloaded.  
**Note:** If a misload is detected, do not attempt to reload the bioprosthesis. Discard the entire system. The valve, catheter, loading system, loading tray, and saline must all be replaced with new sterile components.
23. Attach a 10 mL syringe filled with sterile saline to the stability layer flush port on the distal end of the handle and flush.
24. Remove the loading stylet from the guidewire lumen at the capsule.
25. Attach a 10 mL syringe filled with sterile saline to the wire lumen flush port on the proximal end of the handle and flush.
26. Attach a 10 mL syringe filled with sterile saline to the Evolut PRO+ inline sheath flush port and flush.
27. Before inserting into a patient, visually inspect the loaded bioprosthesis under fluoroscopy.  
**Note:** If a misload is detected, do not attempt to reload the bioprosthesis. Discard the entire system. The valve, catheter, loading system, loading tray, and saline must all be replaced with new sterile components.
28. Leave the bioprosthesis submerged in sterile saline until implantation.

## 7.5 Bioprosthesis implantation

**Note:** Use systemic anticoagulation during the implantation procedure based on hospital or physician protocol. If heparin is contraindicated, consider an alternative anticoagulant.

### 7.5.1 Bioprosthesis deployment

1. Prepare the vascular access site according to standard practice.

**Note:** The need for predilatation of the native aortic valve is determined by the heart team. If needed, predilate the native aortic valve with an appropriate diameter valvuloplasty balloon.

**Note:** Balloon predilatation of a stenotic surgical aortic bioprosthesis has not been evaluated. In cases where there is severe stenosis, predilatation of the surgical aortic bioprosthesis may be performed, and the steps used are identical to native valve predilatation.

2. Backload the catheter onto the guidewire. Insert the catheter tip and capsule through the access site, while maintaining the Evolut PRO+ inline sheath tip against the proximal end of the capsule. Then, insert the Evolut PRO+ inline sheath through the access site, maintaining contact with the capsule. Maintain guidewire position across the aortic valve.

**Note:** For direct aortic access procedures, use a separate introducer sheath; do not use the Evolut PRO+ inline sheath. Maintain the Evolut PRO+ inline sheath at the proximal end of the catheter throughout the procedure.

3. Under fluoroscopic guidance, advance the catheter over the guidewire to the aortic annulus. **Do not** rotate the catheter as it is advanced; rotating the handle does not rotate the capsule.

**Caution:** There will be some resistance when the catheter is advanced through the vasculature. If there is a significant increase in resistance, stop advancement and investigate the cause of the resistance (for example, magnify the area of resistance) before proceeding. Do not force passage. Forcing passage could increase the risk of vascular complications (for example, vessel dissection or rupture).

**Caution:** Persistent force on the catheter can cause the catheter to kink, which could increase the risk of vascular complications (for example, vessel dissection or rupture).

4. Position the catheter so that the bioprosthesis is at the recommended target depth of 3 mm relative to the valve annulus. If the implant depth is <1 mm or >5 mm, consider recapture (*Section 7.5.2*).

**Caution:** Bioprosthesis implant depth <1 mm may contribute to an increased risk of prosthetic valve migration. Bioprosthesis implant depth >5 mm may contribute to an increased risk of conduction disturbances, which may require a permanent pacemaker.

**Note:** For surgical bioprosthetic valves, consider the features of the valve when determining placement of the bioprosthesis.

**Note:** Physicians should consider patient anatomy when determining implant depth.

5. To deploy the bioprosthesis, rotate the deployment knob in the direction of the arrows. The capsule retracts and exposes the bioprosthesis. Continue deploying the bioprosthesis in a controlled manner, adjusting valve position as necessary and noting the position of the radiopaque capsule marker band and paddle attachment.

**Warning:** Use the deployment knob to deploy and recapture the bioprosthesis. Do not use the trigger for deploying or recapturing because it could cause inaccurate placement of the bioprosthesis.

**Note:** Consider using controlled pacing (90 to 120 bpm) because it may increase valve stability during this stage of deployment, especially in patients with larger anatomies.

**Note:** Slight antegrade repositioning of a partially deployed bioprosthesis (before the radiopaque capsule marker band reaches the distal end of the radiopaque paddle attachment) can be achieved by carefully withdrawing the catheter.

**Caution:** Use the catheter handle to reposition the bioprosthesis. **Do not** use the outer catheter shaft.

6. Before the radiopaque capsule marker band reaches the distal end of the radiopaque paddle attachment, evaluate the bioprosthesis position.

**Note:** When the bioprosthesis is approximately 2/3 deployed, the deployment knob provides a tactile indication as a notification before the point of no recapture. Once the radiopaque capsule marker band reaches the distal end of the radiopaque paddle attachment, it is at the point of no recapture.

7. Either complete bioprosthesis deployment or initiate bioprosthesis recapture.

### 7.5.2 Bioprosthesis recapture (optional)

The bioprosthesis is recapturable during deployment before the radiopaque capsule marker band reaches the distal end of the radiopaque paddle attachment. Deployment of the bioprosthesis can be attempted 3 times. If the bioprosthesis is recaptured a third time, it must be removed from the patient.

1. Rotate the deployment knob in the opposite direction of the arrows to recapture the bioprosthesis. A partially recaptured bioprosthesis can be repositioned or fully recaptured.

**Warning:** Use the deployment knob to deploy and recapture the bioprosthesis. Do not use the trigger for deploying or recapturing because it could cause inaccurate placement of the bioprosthesis.

2. To fully recapture the bioprosthesis, continue rotating the deployment knob until the gap between the capsule and catheter tip is closed.  
**Caution:** Stop advancing the capsule once the gap to the catheter tip is closed. Advancing the capsule farther could damage the capsule.
3. Reposition the recaptured bioprosthesis at the recommended target depth of 3 mm relative to the valve annulus. If the implant depth is <1 mm or >5 mm, consider recapture.  
**Caution:** Bioprosthesis implant depth <1 mm may contribute to an increased risk of prosthetic valve migration. Bioprosthesis implant depth >5 mm may contribute to an increased risk of conduction disturbances, which may require a permanent pacemaker.  
**Note:** For surgical bioprosthetic valves, consider the features of the valve when determining placement of the bioprosthesis.  
**Note:** Physicians should consider patient anatomy when determining implant depth.
4. Redeploy the bioprosthesis (*Section 7.5.1, Step 5 and Step 6*).
5. Either complete bioprosthesis redeployment or initiate bioprosthesis recapture. If the bioprosthesis has been recaptured 3 times, withdraw the recaptured bioprosthesis.

### 7.5.3 Catheter removal

1. Under fluoroscopic guidance, confirm that the catheter tip is coaxial with the inflow portion of the bioprosthesis.
2. For transfemoral access, withdraw the catheter until the catheter tip is positioned in the descending aorta. For direct aortic access and subclavian access, withdraw the catheter until the catheter tip is close to the distal tip of the introducer sheath.
3. Under fluoroscopic guidance, close the catheter capsule before withdrawal.  
**Caution:** Close the capsule until it is aligned with the catheter tip. Do not overcapture the catheter tip, because it could interfere with catheter withdrawal through the introducer sheath or cause vessel trauma upon removal.  
**Caution:** Ensure the capsule is closed before catheter removal.  
**Caution:** When using a separate introducer sheath, if increased resistance is encountered when removing the catheter through the introducer sheath, do not force passage. Increased resistance may indicate a problem and forced passage may result in damage to the device and/or harm to the patient. If the cause of resistance cannot be determined or corrected, remove the catheter and introducer sheath as a single unit over the guidewire, and inspect the catheter and confirm that it is complete.
4. Withdraw the catheter until the capsule meets the distal end of the Evolut PRO+ inline sheath.  
**Note:** For direct aortic access procedures, maintain the Evolut PRO+ inline sheath at the proximal end of the catheter.
5. Withdraw the catheter and Evolut PRO+ inline sheath together, and dispose of the device in accordance with local regulations and hospital procedures.
6. Perform routine aortogram to assess the bioprosthesis for proper expansion and function.

**Caution:** Overexpansion of the narrowest portion (waist) of the Evolut PRO+ TAV beyond the levels set forth in *Table 4* has been demonstrated through bench data to cause damage to the bioprosthetic leaflets. Complaints of damage to the bioprosthetic leaflets during post-implant balloon dilatation have been reported in some clinical cases, resulting in moderate to severe aortic insufficiency, which may be detected acutely or during follow-up.

## 8 Training

Medtronic provides detailed operator training. For training material and in-service support, contact a Medtronic representative.

**Caution:** Implantation using the Evolut PRO+ system should be performed only by physicians who have received Evolut PRO+ training.

## 9 MRI safety information

Nonclinical testing and modeling have demonstrated that the Evolut PRO+ bioprosthesis is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3.0 T
- Maximum spatial field gradient of 2500 gauss/cm (25 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg (Normal Operating Mode)

Based on nonclinical testing and modeling, under the scan conditions defined above, the Evolut PRO+ bioprosthesis is expected to produce a maximum in vivo temperature rise of less than 4.0°C after 15 minutes of continuous scanning. Based on nonclinical data, the image artifact caused by the device will extend no greater than 7 mm from the Evolut PRO+ bioprosthesis when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

Scanning under the conditions defined above may be performed immediately after implantation.

For deployment of the Evolut PRO+ bioprosthesis inside of a failed surgical bioprosthetic aortic valve, consult the MRI labeling pertaining to the failed valve for additional image artifact information. The presence of other implants or medical circumstances of the patient may require lower limits on some or all of the above parameters.

## **10 Disclaimer of warranty**

**The following disclaimer of warranty applies to customers outside the United States:**

**ALTHOUGH THE MEDTRONIC EVOLUT PRO+ SYSTEM, EVOLUT PRO+ TRANSCATHETER AORTIC VALVE, EVOLUT PRO+ DELIVERY CATHETER SYSTEM, AND EVOLUT PRO+ LOADING SYSTEM, HEREAFTER REFERRED TO AS “PRODUCT”, HAVE BEEN CAREFULLY DESIGNED, MANUFACTURED, AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.**

The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this DISCLAIMER OF WARRANTY is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this DISCLAIMER OF WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid.

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Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks ("TM") belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries.
Medtronic CoreValve LLC
Medtronic B.V.