24 March 2009

Mr. Jeff Wang
LumiQuick Diagnostics, Inc.
2946 Scott Blvd.
Santa Clara, CA 95054

Dear Mr. Jeff Wang:

I am writing to inform you that today, we have notified by registered mail the Competent Authority in the following countries:

Austria	Bulgaria	Cyprus	Czech Republic	Denmark	Estonia
Finland	France	Germany	Greece	Hungary	Iceland
Ireland	Italy	Latvia	Liechtenstein	Lithuania	Luxembourg
Malta	The Netherlands	Norway	Poland	Portugal	Switzerland
Romania	Slovakia	Slovenia	Spain	Sweden	
United Kingdom					

With this notification, LumiQuick Diagnostics, Inc. has met the requirements of the In-vitro Diagnostics Directive, 98/79/EC for the following devices:

- Adeno/Rota Virus
- Cardiac Marker
- Dengue IgG/IgM Combo (registered only in Italy and The Netherlands)
- Drugs of Abuse
- Fecal Occult Blood (registered only in Italy and The Netherlands)
- H. Pylori Ab/Ag
- HCG
- Legionella (registered only in Italy and The Netherlands)
- LH (registered only in Italy and The Netherlands)
- Strep A (registered only in Italy and The Netherlands)

As of today and without any further notice from the respective Competent Authorities, LumiQuick Diagnostics, Inc. can consider the respective devices and Authorized Representative as officially registered.

If you have any questions, please do not hesitate to contact me.

Yours sincerely,

Rene van de Zande
President & CEO
Emergo Europe

Declaration of Conformity

PRODUCT IDENTIFICATION		
Product name	Model/number	
H. Pylori Ab/Ag Test Devices		
QuickProfile H. Pylori Antigen Test Card	71020	
QuickProfile H. Pylori Antibody Test Card Whole Blood	71024	
QuickProfile H. Pylori Antibody Test Card-Serum	71046	
QuickProfile H. Pylori Antigen Test Strip	71061	
QuickProfile H. Pylori Antibody Serum Test Strip	71064	
QuickProfile H. Pylori Antibody WB Test Strip	71086	
MANUFACTURER		
Name of company	Address	Representative
LumiQuick Diagnostics, Inc.	2946 Scott Blvd. Santa Clara, CA 95054 USA	Jeff Wang
AUTHORIZED REPRESENTATIVE		
Name of company	Address	Telephone/email
Emergo Europe	Prinsessegracht 20 2514 AP The Hague, Netherlands	+31.70.345.8570 - phone +31.70.346.7299 - fax europe@emergogroup.com
CONFORMITY ASSESSMENT		
Device classification	Route to compliance	Standards applied
Class: Self-Certify	Annex III of IVDD 98/79/EC Council Directive	ISO 13485:2003

LumiQuick Diagnostics, Inc. declares that the above mentioned products meet the provision of the Council Directive 98/79/EC for In Vitro Diagnostic Medical Devices and Directive 98/79/EC as transposed in the national laws of the Member States.

COMPANY REPRESENTATIVE: Jeff Wang

TITLE: Quality Systems Manager

SIGNATURE:

DATE: 28/04/2017

Declaration of Conformity

PRODUCT IDENTIFICATION		
Product name	Model/number	
Infectious diseases		
HIV I&II Test Strip	71001	
HIV I&II Test Card	71002	
HCV Antibody Test Card	71030	
HBsAg Test Card	71004	
MANUFACTURER		
Name of company	Address	Representative
LumiQuick Diagnostics, Inc.	2946 Scott Blvd. Santa Clara, CA 95054 USA	Jeff C. Wang
AUTHORIZED REPRESENTATIVE		
Name of company	Address	Telephone/email
Emergo Europe	Molenstraat 15 2513 BH The Hague, Netherlands	+31.70.345.8570 - phone +31.70.346.7299 - fax service@emergogroup.com
CONFORMITY ASSESSMENT		
Device classification	Route to compliance	Standards applied
Class: Self-Certify		ISO 13485:2003

LumiQuick Diagnostics, Inc. declares that the above mentioned products meet the provision of the Council Directive 98/79/EC for In Vitro Diagnostic Medical Devices and Directive 98/79/EC as transposed in the national laws of the Member States.

COMPANY REPRESENTATIVE: Jeff C. Wang

TITLE: Quality Systems Manager

SIGNATURE:

Date:
2017.02.23
11:59:53-08'00'

DATE: 23/02/2017

Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 13485:2016

This is to certify that:

LumiQuick Diagnostics, Inc.
2946 Scott Blvd
Santa Clara
California
95054
USA

Holds Certificate No:

FM 574919

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 for the following scope:

The design, development, manufacture and distribution of in vitro diagnostics test kits and reagents used in the diagnosis and management of disease status, including Infectious Diseases tests, Drugs of Abuse tests, Cardiac Monitor tests, Cancer Marker tests, Fertility Hormone tests, ELISA tests & Urine Chemistry tests.

For and on behalf of BSI:

Gary E Slack, Senior Vice President - Medical Devices

Original Registration Date: 2011-10-20

Latest Revision Date: 2020-08-31

Effective Date: 2020-10-20

Expiry Date: 2023-10-19

Page: 1 of 1

...making excellence a habit.™

LumiQuick Diagnostics, Inc.

2946 Scott Blvd., Santa Clara, CA 95054, USA

Tel: 1-408-855-0061

Fax: 1-408-855-0063

E-mail: info@lumiquick.com

Website: www.lumiquick.com

LETTER OF AUTHORIZATION

We, LumiQuick Diagnostics, Inc., having a registered office at 2946 Scott Blvd, Santa Clara, CA 95054, USA assign SRL SANMEDICO, having a registered office at A. Corobceanu street 7A, apt. 9, Chişinău MD-2012, Moldova, as our authorized representative in correspondence with the conditions of directive 98/79/EEC.

We declare that the company mentioned above is authorized to register, notify, renew or modify the registration of medical devices on the territory of the Republic of Moldova.

This letter is valid through December 31, 2023 and will automatically renewed upon the agreement of both companies. Should you have questions, please contact us.

Best regards,

A handwritten signature in black ink, appearing to read 'Charles Yu', is written over a light blue horizontal line.

Charles Yu

President

Date: January 19, 2022

LumiQuick
DIAGNOSTICS, INC.

QuickProfile™ HELICOBACTER PYLORI ANTIGEN TEST

Movies available at YouTube  : www.youtube.com/lumiquickinc

Cat.# 71020

Immunochromatographic rapid assay for the Detection of Helicobacter pylori Antigens in Human Stool Specimens

INTENDED USE

QuickProfile™ H. pylori Antigen Test is an *in vitro* qualitative immunochromatographic assay for the rapid detection of *Helicobacter pylori* antigens in human stool specimen. This test is intended to aid in the diagnosis of *H. pylori* infection, to monitor the effectiveness of therapeutic treatment and to confirm the eradication of *H. pylori* in peptic ulcer patients.

INTRODUCTION

Helicobacter pylori is a corkscrew-shaped, gram-negative rod that lives in the mucous layer of the stomach. *H. pylori* infection is now accepted as the most common cause of gastritis, and is etiologically involved in gastric ulcer, duodenal ulcer, gastric adenocarcinoma and primary gastric B-cell lymphoma.^{1,2}

The organism is very common, infecting at least half of the world's population. *H. pylori* infection is typically acquired in childhood. Once acquired, infection persists chronically, probably continuing in the stomach throughout life. The damage to gastric structure and function of stomach is constant and direct. Approximately one in six of *H. pylori* infections develop peptic ulcer disease and a small portion of *H. pylori* infection leads to gastric cancer.³

The diagnostic tests for *H. pylori* can be classified into two categories: Invasive and Noninvasive tests. Direct detection by invasive test procedures requires an endoscopy and biopsy specimens from antrum and stomach body.⁴ The presence of *H. pylori* is then confirmed by direct culture, histological examination or rapid urease test. The endoscopy and biopsy specimens offer direct detection of active *H. pylori* infections. Although the procedure is highly specific and high positive predictive value, the cost and discomfort to the patients are very high.

The most widely available noninvasive test is probably the serological based test. The serology test detects *H. pylori* specific IgG antibody in patient serum with current or prior infection.^{5,6} Serology test is a simple, convenient test with relative high sensitivity. The main limitation of serology test is the inability to distinguish current and past infections. Antibody may be present in the patient's serum long after eradication of the organism.⁶ The urease breath test (UBT) with ¹⁴C or ¹³C labeled urea, is a noninvasive test based on the urease activity of the organism. UBT detects active *H. pylori* infection and is highly sensitive and specific. The UBT requires a high density and active bacteria and should not be performed until 4 weeks after therapy to allow residual bacterial to increase to a sufficient number for detection.⁷

QuickProfile™ H. pylori Antigen Test is an immunochromatographic assay that uses antibody-coated colloidal gold to detect the presence of *H. pylori* antigens in stool


specimens. The test detects directly antigens in specimens for an active infection. The test is simple and easy to perform and the test results can be visually interpreted within 10 minutes

PRINCIPLE OF THE TEST

QuickProfile™ H. pylori Antigen Test is a sandwich solid phase immunochromatographic assay. To perform the test, an aliquot of diluted stool sample is added to the sample well of the test cassette. The sample flows through a label pad containing H. pylori antibody coupled to red-colored colloidal gold. If the sample contains H. pylori antigens, the antigen will bind to the antibody coated on the colloidal gold particles to form antigen-antibody-gold complexes. These complexes move on the nitrocellulose membrane by capillary action toward the test line region on which H. pylori specific antibodies are immobilized. As the complexes reach the test line, they will bind to the antibody on the membrane in the form of a line. A second red control line will always appear in the result window to indicate that the test has been correctly performed and the test device functions properly. If H. pylori antigen is not present or lower than the detection limit of the test, only the control line will be visible. If the control line does not develop, the test is invalid.

MATERIALS PROVIDED

1. QuickProfile™ H. Pylori Antigen test card
Each cassette contains a test strip with H. pylori specific antibody on the test region of the membrane and colored H. pylori antibody-gold conjugate pad.

H. Pylori Antigen
Test Card

Sample bottle

2. Sample bottle
Each sample bottle contains 1.5 ml of stool specimen collection buffer. Store at 4-30°C

MATERIALS NOT PROVIDED

1. Specimen collection container
2. Timer.

WARNINGS AND PRECAUTIONS

1. For *in vitro* diagnostic use.
2. Wear protective glove while handling kit components and test specimens.
3. Patient specimens and inactivated Positive Control may contain infectious agents and should be handled and disposed of as potential biohazards.
4. Do not use kit components beyond expiration date.
5. Dispose all used materials in appropriate container. Treat as potential biohazard.

STORAGE INSTRUCTION

1. The expiration date is indicated on the package label.
2. Sample Collection Tubes without introducing the sample can be stored at 4-30°C.
3. Test device can be stored at 4-30 °C.

SPECIMEN COLLECTION AND STORAGE

Stool specimens should be collected in containers that do not contain media, preservatives, animal serum or detergents as any of these additives may interfere with the QuickProfile™ H. pylori Antigen Test.

Specimens may be stored at 2-8°C for 3 days without interfering with the assay performance. For long-term storage of specimens, -20°C or colder is recommended. Repeated freezing and thawing of specimens is not recommended and may cause erroneous results. Do not store specimens in self-defrosting freezers.

REAGENT PREPARATION

Bring all reagents, including test device, to room temperature before use.

SPECIMEN PREPARATION

1. Unscrew the sample bottle, use the attached applicator stick attached on the cap to transfer small piece of stool (5-6 mm in diameter; approximately 100 mg – 200 mg/0.1-0.2 g) into the sample bottle containing specimen preparation buffer.
2. Replace the stick in the bottle and tighten securely. Mix stool sample with the buffer thoroughly by shaking the bottle for a few seconds.

Note: Watery or diarrhea specimens are inappropriate for testing.

PROCEDURE

1. Bring all materials and specimens to room temperature.
2. Remove the test card from the sealed foil pouch.
3. Hold the sample bottle upright with the tip pointed away from the test performer, snap off the tip.
4. Hold the bottle in a vertical position over the sample well of the test card, deliver 3 drops (120 -150 µL) of diluted stool sample to the sample well.
5. Read the result within 10 to 15 minutes. A strong positive sample may show result earlier.


Test results after 15 minutes may not be accurate.

INTERPRETATION OF RESULTS

Positive result: A distinct pink colored band appears on test line regions, in addition to a pink line on the control line region.

Negative result: No line appears in the test line region. A distinct pink line shows on the control line region.

Invalid: The control line next to the test line does not become visible within 10 minutes after the addition of the sample.

QUALITY CONTROL

1. The control band is an internal reagent and procedural control. It will appear if the test has been performed correctly and the reagents are reactive.
2. Good Laboratory Practice recommends the daily use of control materials to validate the reliability of the device. Control materials are not provided with this test kit but may be commercially available.

LIMITATIONS

1. The test is for qualitative detection of H. pylori antigen in stool sample and does not indicate the quantity of the antigens.
2. The test is for *in vitro* diagnostic use only.
3. The test result should be used only to evaluate with patient with signs and symptoms of gastrointestinal disease. A definitive clinical diagnosis should only be made by the physician after all clinical and laboratory findings have been evaluated.

EXPECTED VALUES

Helicobacter pylori infects more than half the people in the world.⁹ The prevalence of the infection varies among countries and among different groups within the same country.¹⁰ The prevalence rate in the United States suggests an incidence of infection of 2%. The lifetime prevalence of peptic ulcer disease is about 12% in men and 9% in women.¹¹ Studies have found that more than 90% of patients with duodenal ulcer and 80% of patients with gastric ulcer are infected with *H. pylori*.^{12,13}

QuickProfile™ H. pylori Antigen Test detects the presence of *H. pylori* antigens in stool specimens. Expected values for any given population should be determined for each laboratory. The positivity rate of any given laboratory may vary depending on geographic location, ethnic group, and living environment.

PERFORMANCE CHARACTERISTICS

Accuracy

QuickProfile™ Helicobacter pylori Antigen Test was evaluated on 1049 stool samples. The test results were compared with an approved predicate kit.

		Predicate Kit		Total
		Pos.	Neg.	
QuickProfile™ H. Pylori Ag Test	Pos.	508	11	519
	Neg.	5	525	530
Total		513	536	1049

Out of five hundred and thirteen (513) samples that were tested positive by the predicate kit, five hundred and eight (508) were positive on QuickProfile™ H. Pylori Antigen Test. Out of five hundred and thirty six (536) samples that were tested negative by the predicate kit, five hundred and twenty five (525) were negative on QuickProfile™ H. Pylori Antigen Test. Sixteen (16) samples that had disparity in results were verified by ELISA. Seven (7) samples had results in agreement with QuickProfile™ H. Pylori Antigen Test while nine (9) samples agreed with the predicate kit. The agreement with the predicate kit is summarized as below.

Agreement of positive = 508/513 = 99.03%
 Agreement of Negative = 525/536 = 97.95%
 Total Agreement = 1033/1049 = 98.47%

Assay Specificity


Following bacterial and viral strains were used to test the specificity of QuickProfile™ H. pylori Antigen Test. Positive and negative controls spiked with the bacteria or virus at the indicated concentration showed no interference on the test results.

Adenovirus type 40	1x10 ⁶ TCID ₅₀
Adenovirus type 41	1x10 ⁶ TCID ₅₀
Rotavirus Wa	1x10 ⁶ TCID ₅₀
Campylobacter jejuni	7.63 x 10 ⁷ CFU/ml
Candida albicans	1x10 ⁸ CFU/ml
Clostridium perfringens A	1x10 ⁸ CFU/ml
Citrobacter freundii	1x10 ⁸ CFU/ml
Enterococcus faecalis	1x10 ⁸ CFU/ml
Escherichia coli	1x10 ⁸ CFU/ml
Klebsiella pneumonia	1x10 ⁸ CFU/ml
Listeria monocytogenes	1x10 ⁸ CFU/ml
Moraxella catarrhalis	9.9x10 ⁶ CFU/ml
Neisseria gonorrhoeae	1x10 ⁸ CFU/ml
Pseudomonas aeruginosa	1x10 ⁸ CFU/ml
Staphylococcus epidermidis	1x10 ⁸ CFU/ml
Staphylococcus aureus	1x10 ⁸ CFU/ml
Shigella flexneri	1x10 ⁸ CFU/ml
Shigella sonnei	1x10 ⁸ CFU/ml
Streptococcus dysgalactiae	1x10 ⁸ CFU/ml
Streptococcus agalactiae	1x10 ⁸ CFU/ml
Streptococcus pyogenes	1x10 ⁸ CFU/ml

REFERENCES

1. Marshall, B.J. and Warren, J.R. Unidentified curved bacilli in the stomach of patients with gastric and peptic ulceration. *Lancet I*:1984: 1311-1314.
2. Graham K.S and Graham D.Y. 1999. Contemporary Diagnosis and Management of *H. pylori*-Associated Gastrointestinal Diseases, Handbooks in Health Care Co., Newtown, PA., 1999: 39-67.
3. Howden C.W. Clinical expressions of *Helicobacter pylori* infection. *Am J Med*; 1996;100:27S-33S.
4. El-Zimaity HM, Al-Assi MT, Genta RM, Graham DY. Confirmation of successful therapy of *Helicobacter pylori* infection: number and site of biopsies or a rapid urease test. *Am J Gastroenterol*. 1995;90:1962-1964.

5. Talley NJ, Newell DG, Ormand JE, et al. Serodiagnosis of *Helicobacter pylori*: Comparison of enzyme-linked immunosorbent assays. *J. Clin Microbiol*. 1991;29:1635-1639.
6. Cutler AF. Testing for *Helicobacter pylori* in clinical practice. *Am J. Med*. 1996;100:35S-41S.
7. Klein PD, Malaty HM, Martin RF, et al. Noninvasive detection of *Helicobacter pylori* infection in clinical practice: the ¹³C urea breath test. *Am J. Gastroenterol*. 1996;91:690-694.
8. National Committee for Clinical Laboratory Standards. Internal quality control: Principles and definitions; Approval Guideline, NCCLS document C24-A (NCCLS, 771 East Lancaster Ave, Vällanova, PA 19085, 1991).
9. Marshall BJ. *JAMA*. 1995;274:1064-1066
10. Breuer T, Malaty HM, Graham DY. The epidemiology of *H. pylori*-associated gastroduodenal diseases. In: Ernst PB, Michetti P, Smith PD, eds. The Immunobiology of *H. pylori*: From Pathogenesis to Prevention. Philadelphia. Lippincott-Raven, 1997:1-14.
11. Graham DY, Malaty HM, Evans DG, Evans, Jr. DJ, Klein PD, and Adam E. Epidemiology of *Helicobacter pylori* in a asymptomatic population in the United States. Effect of age, race, and socioeconomic status. *Gastroenterology*, 1991;100:1495-1501.
12. Anand BS, Raed AK, Malaty HM, et al. Low point prevalence of peptic ulcer in normal individual with *Helicobacter pylori* infection. *Am J Gastroenterol*. 1996;91:1112-1115.
13. Tytgat GNJ, Noach LA, Rauws EAJ. *Helicobacter pylori* infection and duodenal ulcer disease. *Gastroenterol Clin North Am*. 1993;22:27-139.

LumiQuick Diagnostics, Inc.
 2946 Scott Blvd.
 Santa Clara, CA 95054 USA

Emergo Europe
 Molenstraat 15
 2513 BH The Hague
 The Netherlands

Tel : (408) 855.0061
 Fax: (408) 855.0063
 Email: info@lumiquick.com
www.lumiquick.com
 Movies available at YouTube

 :
www.youtube.com/lumiquickinc

LumiQuick
DIAGNOSTICS, INC.

QuickProfile™ HBsAg TEST

FOR THE QUALITATIVE ASSESSMENT OF HBsAg IN HUMAN SERUM, PLASMA OR WHOLE BLOOD

REF 71003 HBsAg Test Strip
REF 71004 HBsAg Test Card

For In Vitro Diagnostic Use Only

INTENDED USE

QuickProfile™ HBsAg Test is an immunochromatography assay for the qualitative detection of Hepatitis B virus surface antigen (HBsAg) in human serum, plasma or whole blood specimen.

TEST PRINCIPLE

QuickProfile™ HBsAg Test is a double antibody sandwich immunoassay. Colloidal gold conjugated anti-HBsAg antibody complexes are dry-immobilized in the test device. When the sample is added, it migrates by capillary diffusion through the strip re-hydrating the gold conjugate complexes. If present, HBsAg will react with the gold conjugate complexes forming particles. These particles will continue to migrate along the strip until the Test Zone (T) where they are captured by anti-HBsAg antibodies immobilized there and a visible red line appears. If there is no HBsAg in sample, no red line will appear in the Test Zone (T). The gold conjugate complexes will continue to migrate alone until they are captured in the Control Zone (C) by immobilized goat anti-mouse IgG antibody aggregating a red line, which indicates the validity of the test.

MATERIALS PROVIDED

1. QuickProfile™ HBsAg Test
2. Instructions for use
3. Disposable transfer pipet

MATERIALS REQUIRED BUT NOT SUPPLIED

1. Whole blood or plasma: Vacutainer tube, or other appropriate tube, containing heparin or EDTA as an anticoagulant
2. Serum: Vacutainer tube, or other appropriate tube, without anticoagulant
3. Timer or clock

STORAGE

The sealed pouches in the test kit may be stored between 4-30°C for the duration of the shelf life as indicated on the pouch. The test must be used immediately after being removed from the sealed pouch.

PRECAUTIONS

1. This kit is for **in vitro** diagnostic use only.
2. This kit is for **PROFESSIONAL** use only.
3. Read the instructions carefully before performing the test.
4. This product does not contain any human source materials.
5. Do not use kit contents after the expiration date.
6. Handle all specimens as potentially infectious.
7. Follow standard Lab procedure and biosafety guidelines for handling and disposal of potentially infective material. When the assay procedure is completed, dispose of specimens after autoclaving them at 121° C for at least 20 min. Alternatively, they can be treated with 0.5% Sodium Hypochlorite for 1-2 hours before disposal.
8. Do not pipette by mouth. Do not smoke, eat, or drink in areas where reagents or specimens are handled.

SPECIMEN COLLECTION AND PREPARATION

1. The serum, plasma or whole blood specimen should be collected under standard laboratory conditions.
2. Heat inactivation of specimens, which may cause hemolysis and protein denaturation, should be avoided.
3. Patient samples performed best when tested immediately after collection. If specimens are to be stored, the red blood cells should be removed to avoid hemolysis. If the sample cannot be tested within 24 hours, serum or plasma should be frozen until the test can be performed. Allow sample to reach room temperature before proceeding.
4. Sodium azide can be added as a preservative up to 0.1% without effecting the test results.

PROCEDURE

For HBsAg Test Strip (Catalog Number: 71003)

1. Bring all materials and specimens to room temperature.
2. Remove the test strip from the sealed foil pouch.
3. Label the test strip with specimen identity by writing the ID on the top label of the strip.
4. Place the test strip on a flat horizontal surface.
5. Use the transfer pipet to draw up the sample.
6. Hold the transfer pipet in a vertical position over the sample pad and dispense 2 drops (80-100 µl) of sample onto the sample pad.
7. Read the result at 20 minutes after adding the sample.


For HBsAg Test Card (Catalog Number: 71004)

1. Bring all materials and specimens to room temperature.
2. Remove the test card from the sealed foil pouch.
3. Label the test card with specimen identity on the "ID ____" area of the cassette.
4. Place the test card on a flat horizontal surface.
5. Use the transfer pipet to draw up the sample.
6. Hold the transfer pipet in a vertical position over the sample well and dispense 2 drops (80-100 µl) of sample into the sample well.
7. Read the result at 20 minutes after adding the sample.


Note: Some positive samples may show positive results before 20 minutes. Results after 30 minutes may not be accurate.

INTERPRETATION OF RESULTS

Test Strip:

Test Card:

Positive:

Two colored bands appear within 20 minutes. One colored band appears in the Control Zone (C) and another colored band appears in the Test Zone (T). The test result is positive and valid. No matter how faint the colored band appears in the Test Zone (T), the test result should be considered as positive result.

Negative:

One colored band appears in the Control Zone (C) within 20 minutes. No colored band appears in the Test Zone (T). The test result is negative and valid.

Invalid result:

No colored band appears in the Control Zone (C) within 20 minutes. The test result is invalid. Repeat the test with a new test device.

QUALITY CONTROL

1. The control band is an internal reagent for procedural control. It will appear if the test has been performed correctly and the reagents are reactive.
2. Good Laboratory Practice recommends the daily use of control materials to validate the reliability of the device. Control materials are not provided with this test kit but may be commercially available.

LIMITATIONS

1. The test is for in vitro diagnostic use only.
2. Negative results do not rule out the possibility of hepatitis B exposure or infection. Infection through recent exposure to HBV may not be detectable.
3. The presumptive positive result obtained with QuickProfile™ HBsAg Test alone cannot be the final diagnosis of hepatitis B infection. As in case of all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test but should rather be made after all the clinical findings have been evaluated.
4. This test is intended ONLY for testing of an individual serum, plasma or whole blood sample. DO NOT use it for testing of other body fluids or pooled blood samples.
5. The test is for qualitative detection of HBsAg in human serum, plasma or blood sample and does not indicate the quantity of the antibodies.

PERFORMANCE CHARACTERISTICS:

Sensitivity:

QuickProfile™ HBsAg Test can detect HBsAg with a concentration of 1.0 ng/ml.

Accuracy:

Twelve hundred and forty-nine (1249) ELISA confirmed samples, including four hundred and sixty-nine (469) positive samples and seven hundred and eighty (780) negative samples were used for the clinical evaluation of QuickProfile™ HBsAg Test. The results are summarized in the following table. The sensitivity, specificity and accuracy are all 99.6%.

QuickProfile™ HBsAg Test	ELISA HBsAg Test	
	Positive	Negative
Positive	467	3
Negative	2	777
Agreement	99.6%	99.6%

BIBLIOGRAPHY

1. Sehulster, L. et al. Immunological and biophysical alteration of Hepatitis B virus antigens by sodium hypochlorite disinfection, Appl. And Envir. Microbiol., 42:762-767, 1981.
2. U.S. Department of Health and Human Services. Biosafety in microbiological and biomedical laboratories. HHS Publication(NIH) 88-8395. Washington:U.S. Government Printing Office, May 1988.

LumiQuick Diagnostics, Inc.
2946 Scott Blvd.
Santa Clara, CA 95054 USA
Tel : (408) 855.0061
Fax: (408) 855.0063
Email: info@lumiquick.com
www.lumiquick.com

LumiQuick
DIAGNOSTICS, INC.

QUICK PROFILE™ HCV ANTIBODY TEST

FOR THE QUALITATIVE ASSESSMENT OF HCV ANTIBODY
IN HUMAN SERUM, PLASMA OR WHOLE BLOOD

REF 71027 HCV Test Strip

REF 71030 HCV Test Card

For In Vitro Diagnostic Use Only

INTENDED USE

QuickProfile™ HCV Ab Test is a chromatographic immunoassay for qualitative detection of the antibodies against hepatitis C virus (HCV Ab) in human serum, plasma or whole blood samples. It is intended for professional use as an aid for diagnosis and management of patients related to infection with hepatitis C as well for primary screening of blood from volunteer donors.

SUMMARY

Hepatitis C virus (HCV) is an envelope, single stranded positive sense RNA (9.5 kb) virus belonging to the family of Flaviviridae. Six major genotypes and series of subtypes of HCV have been identified. Isolated in 1989, HCV is now recognized as the major cause for transfusion associated non-A, non-B hepatitis. The disease is characterized with acute and chronic form. More than 50% of the infected individuals develop severe, life threatening chronic hepatitis with liver cirrhosis and hepatocellular carcinomas. Since the introduction in 1990 of anti-HCV screening of blood donations, the incidence of this infection in transfusion recipients has been significantly reduced. Clinical studies show that significant amount of HCV infected individuals develop antibodies to NS5 non-structural protein of the virus. For this, the third generation tests include antigens from the NS5 region of the viral genome in addition to NS3 (c200), NS4 (c200) and the Core (c22). Third generation tests have improved sensitivity and shorten the time between infection with HCV and the appearance of detectable antibodies (window period) to 60 days.

TEST PRINCIPLE

QuickProfile™ HCV Ab Test employs a chromatographic lateral flow device in a strip or cassette format. Recombinant HCV antigens are immobilized at the Test Zone (T) and goat anti mouse IgG antibodies are immobilized at the Control Zone (C) on the nitrocellulose membrane. When the sample is added, it migrates by capillary diffusion and rehydrating the colloidal gold conjugated recombinant HCV antigens (Au-Ag) dried onto the fiberglass strip. If present in sample, HCV antibodies will bind the gold conjugated antigens forming complexes. These complexes will continue to migrate along the strip until the Test Zone (T) zone where they are captured by the HCV antigens to form a visible red line. The colloidal gold-mouse IgG is used as the indicator for control line. A red line formed by gold-mouse IgG and goat anti-mouse IgG at the Control Zone (C) indicates the validity of the test.

MATERIAL PROVIDED

1. QuickProfile™ HCV Ab Test
2. Sample buffer
3. Instructions for use

MATERIALS REQUIRED BUT NOT PROVIDED

Clock or timer
Specimen collection container
Centrifuge
Biohazard waste container

STORAGE

The sealed pouches in the test kit may be stored between 4-30°C for the duration of the shelf life as indicated on the pouch. The test must be used immediately after being removed from the sealed pouch.

DCR 15-078 71027+71030
5089 E1R1 01/06/2016

PRECAUTIONS

1. This kit is for *in vitro* diagnostic use only.
2. This kit is for **PROFESSIONAL** use only.
3. Read the instructions carefully before performing the test.
4. This product does not contain any human source materials.
5. Do not use kit contents after the expiration date.
6. Handle all specimens as potentially infectious.
7. Follow standard Lab procedure and biosafety guidelines for handling and disposal of potentially infective material. When the assay procedure is completed, dispose of specimens after autoclaving them at 121° C for at least 20 min. Alternatively, they can be treated with 0.5% Sodium Hypochlorite for 1-2 hours before disposal.
8. Do not pipette reagent by mouth and do not no smoke, eat or drink while performing assays.
9. Wear gloves during the whole procedure.


SPECIMEN COLLECTION AND PREPARATION

Fresh serum, plasma or whole blood samples can be used for this assay. Blood collected by venipuncture should be allowed to clot naturally and completely – the serum/plasma must be separated from the clot as early as possible as to avoid hemolysis of the red blood cell. Care should be taken to ensure that the serum samples are clear and not contaminated by microorganisms. Any visible particulate matter in the sample should be removed by centrifugation at 3000 RPM for at least 20 minutes at room temperature, or by filtration with 0.22u filters. Plasma samples collected into EDTA, sodium citrate or heparin may be tested, but highly lipaemic, icteric, or hemolyzed samples should not be used as they could give erroneous results in the assay. Do not inactivate samples by heat. This can cause deterioration of the target proteins in the sample.

ASSAY PROCEDURE


For HCV Test Strip (Catalog Number: 71027)

1. Allow the test strip and sample to reach room temperature if necessary.
2. Open the pouch, Take out the test strip and transfer pipet.
3. Using the transfer pipet to draw up the sample, dispense one drop (approx 40µl) specimen to the sample pad as shown in the illustration, and wait for a few seconds until the sample is completely absorbed by sample pad.
4. Add one drop (approx 40µl) sample buffer to the sample pad as shown in the illustration.
5. Read the results at 20 minutes.

For HCV Test Card (Catalog Number: 71030)

1. Allow the test card and sample to reach room temperature if necessary.
2. Open the pouch, Take out the test card and transfer pipet.
3. Using the transfer pipet to draw up the sample, dispense one drop (approximately 40µl) of specimen to the sample well marked as "S" and wait for a few seconds until the sample is completely absorbed by sample pad.
4. Add one drop (approx. 40 µl) sample buffer into the sample well marked as "S".
5. Read the results at 20 minutes.

**Some positive samples may show positive results before 20 minutes.
Results after 30 minutes may not be accurate.**

INTERPRETATION OF RESULTS

TEST STRIP		<p style="text-align: center;">Positive</p> <p>If two color bands are visible within 20 minutes, the test is positive and valid. The test result can be read as soon as a distinct colored band appears in test zone and control zone respectively. It indicates that HCV Ab has been detected using QuickProfile™ HCV Ab Test.</p>
	<p style="text-align: center;">NEGATIVE</p> <p>Only one red line appears in the Control Zone (C), indicating that no HCV Ab have been detected with QuickProfile™ HCV Ab Test. However, this does not exclude the possibility for infection with HCV.</p>	
TEST CARD		<p style="text-align: center;">INVALID</p> <p>One red line should always appear at Control Zone (C) if no red line appears in the Control Zone (C), the test is invalid. The sample must be re-tested using a new device</p>

QUALITY CONTROL

1. The control band is an internal reagent for procedural control. It will appear if the test has been performed correctly and the reagents are reactive.
2. Good Laboratory Practice recommends the daily use of control materials to validate the reliability of the device. Control materials are not provided with this test kit but may be commercially available.

LIMITATIONS

1. The test is for in vitro diagnostic use only.
2. Negative results do not rule out the possibility of hepatitis C exposure or infection. Infection through recent exposure to HCV may not be detectable.
3. The positive result obtained with QuickProfile™ HCV Ab Test alone cannot be the final diagnosis of hepatitis C infection. As in the case of all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test but should rather be made after all the clinical findings have been evaluated.
4. This test is intended ONLY for testing of individual serum, plasma or whole blood samples. DO NOT use it for testing of other body fluids or pooled blood samples.
5. The test is for qualitative detection of anti-HCV antibody in human serum, plasma or blood sample and does not indicate the quantity of the antibodies.

PERFORMANCE CHARACTERISTICS:

1. Accuracy

In clinical evaluation of the QuickProfile™ HCV Ab Test, 727 confirmed negative and 327 positive samples were tested. A sensitivity of 99.08% (324/327) and a specificity of 99.17% (721/727) were obtained. Overall, agreement with the Predicate Test is 99.15%.


QuickProfile™ HCV Ab Test	Predicate HCV Ab Test	
	Positive	Negative
Positive	324	6
Negative	3	721
Agreement	99.08%	99.17%

2. Interference

No cross reactivity was observed with specimens from patients infected with HAV, HBV, HIV, HTLV, CMV, and TP.

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