

ALL TEST
TCA Rapid Test Cassette
(Whole Blood/Serum/Plasma)
Package Insert
REF DTC-402 English

A rapid test for the qualitative detection of Tricyclic Antidepressants in human whole blood or serum or plasma.

For medical and other professional *in vitro* diagnostic use only.

[INTENDED USE]

The TCA Rapid Test Cassette (Whole Blood/Serum/Plasma) is a lateral flow chromatographic immunoassay for the detection of TCA in whole blood or serum or plasma at a cut-off concentration of 300ng/mL. This test will detect other related compounds, please refer to the Analytical Specificity table in this package insert.

This assay provides only a qualitative, preliminary test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

[SUMMARY]

TCA (Tricyclic Antidepressants) are commonly used for the treatment of depressive disorders. TCA overdose can result in profound CNS depression, cardiotoxicity and anticholinergic effects. TCA overdose is the most common cause of death from prescription drugs. TCA are taken orally or sometimes by injection. TCA are metabolized in the liver. Both TCA and their metabolites are excreted in whole blood or serum or plasma mostly in the form of metabolites for up to ten days.¹

[PRINCIPLE]

The TCA Rapid Test Cassette (Whole Blood/Serum/Plasma) is an immunoassay based on the principle of competitive binding. Drugs that may be present in the whole blood/serum/plasma specimen compete against the drug conjugate for binding sites on the antibody.

During testing, a whole blood/serum/plasma specimen migrates upward by capillary action. TCA, if present in the whole blood/serum/plasma specimen below the cut-off level, will not saturate the binding sites of the antibody in the test. The antibody coated particles will then be captured by immobilized TCA-protein conjugate and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the TCA level exceeds the cut-off level because it will saturate all the binding sites of anti-TCA antibodies.

A drug-positive whole blood/serum/plasma specimen will not generate a colored line in the test line region because of drug competition, while a drug-negative whole blood/serum/plasma specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

[REAGENTS]

The test contains mouse monoclonal anti-TCA antibody coupled particles and TCA-protein conjugate. A goat antibody is employed in the control line system.

[PRECAUTIONS]

- For professional *in vitro* diagnostic use only. Do not use after the expiration date.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Do not use test if pouch is damaged.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing and follow the standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being tested.
- The used test should be discarded according to local regulations.
- Humidity and temperature can adversely affect results.

[STORAGE AND STABILITY]

Store as packaged in the sealed pouch at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

[SPECIMEN COLLECTION AND PREPARATION]

- The TCA Rapid Test Cassette can be performed using whole blood (from venipuncture or fingerstick) or serum or plasma.
- To collect **Fingerstick Whole Blood specimens:**
 - Wash the patient's hand with soap and warm water or clean with an alcohol swab. Allow to dry.
 - Massage the hand without touching the puncture site by rubbing down the hand towards the fingertip of the middle or ring finger.
 - Puncture the skin with a sterile lancet. Wipe away the first sign of blood.
 - Gently rub the hand from wrist to palm to form a rounded drop of blood over the puncture site.

- Add the Fingerstick Whole Blood specimen to the test by using a **capillary tube**:
 - Touch the end of the capillary tube to the blood until filled to approximately 40 µL. Avoid air bubbles.
 - Place the bulb onto the top end of the capillary tube, then squeeze the bulb to dispense the whole blood to the specimen well of the test cassette.
- Testing should be performed immediately after specimen collection. Do not leave the specimens at room temperature for prolonged periods. Serum and plasma specimens may be stored at 2-8°C for up to 3 days, for long-term storage, specimens should be kept below -20°C. Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 2 days of collection. Do not freeze whole blood specimens. Whole blood collected by fingerstick should be tested immediately.
- Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.
- If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiologic agents.

[MATERIALS]

- Materials Provided**
- Test cassettes
 - Droppers
 - Buffer
 - Package insert
- Materials Required But Not Provided**
- Specimen collection containers
 - Lancets (for fingerstick whole blood only)
 - Heparinized capillary tubes and dispensing bulb (for fingerstick whole blood only)
 - Centrifuge
 - Timer

[DIRECTIONS FOR USE]

- Allow the test, specimen, buffer and/or controls to reach room temperature (15-30°C) prior to testing.**
- Bring the pouch to room temperature before opening it. Remove the cassette from the sealed pouch and use it within one hour.
 - Place the cassette on a clean and level surface.

For serum or plasma specimen:



- Hold the dropper vertically and transfer **1 full drop of serum or plasma** (approximately 40µL), then add **2 drops of buffer** (approximately 80µL) to the specimen well(S) of the cassette, and then start the timer. Avoid trapping air bubbles in the specimen well. See illustration below.

For Venipuncture Whole blood specimen:

- Hold the dropper vertically and transfer **1 drop of whole blood** (approximately 40µL) to the specimen well(S), then add **2 drops of buffer** (approximately 80µL), and start the timer. See illustration below.

For Fingerstick Whole blood specimen:

- To use a capillary tube: Fill the capillary tube and transfer approximately 40µL of fingerstick whole blood specimen to the specimen well(S) of test cassette, then add **2 drops of buffer** (approximately 80µL) and start the timer. See illustration below.
- Wait for the colored line(s) to appear. **Read the result at 5 minutes.** Do not interpret the result after 10 minutes.



[INTERPRETATION OF RESULTS]

(Please refer to the illustration above)

NEGATIVE: Two colored lines appear. One colored line should be in the control line region (C) and another colored line should be in the test line region (T). This negative result indicates that the TCA concentration is below the detectable cut-off level.

Interfering Substances

The TCA Rapid Test Cassette (Whole Blood/Serum/Plasma) has been tested for possible interference from visibly hemolyzed and lipemic specimens. In addition, no interference was observed in specimens containing up to 100 mg/dL hemoglobin, up to 100 mg/dL bilirubin and up to 200 mg/dL human serum albumin.

[BIBLIOGRAPHY]

- Tietz NW. *Textbook of Clinical Chemistry*. W.B. Saunders Company, 1998; 1735
- Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*, 2nd Ed. Biomedical Publ., Davis, CA, 1982; 488

Index of Symbols					
	Consult Instructions For Use		Tests per kit		Authorized Representative
	For <i>in vitro</i> diagnostic use only		Use by		Do not reuse
	Store between 2-30 °C		Lot		Lot Number
	Do not use if package is damaged				Catalog #

Hangzhou AllTest Biotech Co., Ltd.
 6500 Yuhang Street
 Hangzhou Economic & Technological Development Area
 Hangzhou, 311019 P.R. China
 Web: www.alltest.com.cn Email: info@alltest.com.cn

CE **EC REP**
 Metshel EC-REP GmbH
 Borkstrasse 10,
 48163 Münster,
 Germany

Number: 145320702
 Revision date: 2023-12-20

blood/Serum/Plasma by the TCA Rapid Test Cassette (Whole Blood/Serum/Plasma) at 5 minutes.

Compound	Concentration (ng/mL)
Nortriptyline	300
Nordoxepine	150
Trimipramine	1,300
Amipriptyline	600
Promazine	1,300
Desipramine	80
Cyclobenzaprine	600
Imipramine	140
Clomipramine	18,000
Doxepine	600
Maprotiline	600
Promethazine	18,000
Perphenazine	18,000

Precision

A study was conducted at three hospitals using three different lots of product to demonstrate the within run, between run and between operator precision. An identical panel of coded specimens, containing no TCA and 50% TCA above and below the 300ng/mL cut-off was provided to each site. The following results were tabulated:

TCA Concentration (ng/mL)	n	Site A	Site B	Site C
0	10	10	0	10
150	10	8	2	9
450	10	1	9	2

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free Whole blood/Serum/Plasma or determine positive whole blood/serum/plasma. The following compounds show no cross-reactivity when tested with the TCA Rapid Test Cassette (Whole blood/Serum/Plasma) at a concentration of 100 µg/mL.

Non Cross-Reacting Compounds

Acetophenetidin	Dextromethorphan	Methadone	Phenylpropanolamine
N-Acetylcocaineamide	Diazepam	D-methamphetamine	Prednisolone
Acetylsalicylic acid	Diclofenac	(l)-methamphetamine	Prednisone
Amisopyrine	Difenisal	Methoxyphenamine	Procaine
Amobarbital	Digoxin	3,4-Methylenedioxymethylamphetamine	D-Proparanolol
Amoxicillin	Diphenhydramine	amphetamine	D-Propoxyphene
Ampicillin	Doxylamine	(±) 3,4-Methylenedioxy-methamphetamine	D-Pseudoephedrine
l-Ascorbic acid	Ecgonine hydrochloride	methamphetamine	Quinidine
Apomorphine	Ecgoninemethylester	Methylphenidate	Quinine
Aspartame	(R,2S)-(-)-Ephedrine	Morphine-3-β-D-glucuronide	Ranitidine
Atropine	l-Ephedrine	glucuronide	Salicylic acid
D,l-Amphetamine	Erythromycin	Nalidixic acid	Sesazolbin
l-Amphetamine	Ethyl-p-aminobenzoate	Naloxone	Serotonin
Benzoic acid	Fenfluramine	Naltrexone	(S)-Hydroxyfluramine
Benzoic acid	Fenpropofen	Naproxen	Sulfamethazine
Benzofegonine	Furosemide	Nacinamide	Salindac
Benzphetamine	Gentamic acid	Nifedipine	Temazepam
Bilirubin	Hemoglobin	Norcocaine	Tetracycline
(±)-Brompheniramine	Hydralazine	(-)-α-Ephedrine	Tetrahydrocortisone
Caffeine	Hydrochlorothiazide	Norethindrone	3 Acetate
Cannabidiol	Hydrocodone	D-Norpropoxyphene	Tetrahydrocortisone
Cannabinol	Hydrocortisone	Nascepine	3 (β-D glucuronide)
Chloralhydrate	p-Hydroxyamphetamine	D,l-Octopamine	Tetrahydrozoline
Chloramphenicol	O-Hydroxyhippuric acid	Oxiclic acid	Thebaine
Chloridazepoxide	3-Hydroxyfluramine	β-Estradiol	Thiamine
Chlorothalidate	p-Hydroxy-methamphetamine	Oxycodone	Thioridazine
(±) Chlorpheniramine	methamphetamine	Oxymetazoline	Tolbutamide
Chlorpromazine	Propafenone	Paracetamol	Triamterene
Chlorquine	(±)-Isoproterenol	Penicillin-G	Trifluoperazine
Cholesterol	Isosucrine	Pentazocine	Trimethoprim
Clonidine	Ketamine	Pentobarbital	D, l-Tryptophan
Cocaine hydrochloride	Ketoprofen	Phencyclidine	Tyramine
Codeine	labetalol	Phenelzine	D, l-Tyrosine
Cortisone	levorphanol	Phenobarbital	Uric acid
(-) Cotinine	loperamide	Phentermine	Verapamil
Creatinine	Meperidine	l-Phenylephrine	Oxazepam
Deoxycorticosterone	Maproamate	β-Phenylethylamine	Zomepirac

*NOTE: The shade of color in the test line region (T) may vary, but it should be considered negative whenever there is even a faint colored line.

POSITIVE: One colored line appears in the control line region (C). No line appears in the test line region (T). This positive result indicates that the TCA concentration exceeds the detectable cut-off level.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

[QUALITY CONTROL]

A procedural control is included in the test. A colored line appearing in the control region (C) is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

[LIMITATIONS]

- The TCA Rapid Test Cassette (Whole blood/Serum/Plasma) provides only a qualitative, preliminary result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.²
- It is possible that technical or procedural errors, as well as other interfering substances in the whole blood or serum or plasma specimen may cause erroneous results.
- A positive result indicates presence of the drug or its metabolites but does not indicate level of intoxication, administration route or concentration in whole blood or serum or plasma.
- A negative result may not necessarily indicate drug-free whole blood/serum/plasma. Negative results can be obtained when drug is present but below the cut-off level of the test.
- Test does not distinguish between drugs of abuse and certain medications.

[PERFORMANCE CHARACTERISTICS]

Accuracy

A side-by-side comparison was conducted using the TCA Rapid Test Cassette and GC/MS at the cut-off of 300ng/mL. Testing was performed on 90 clinical specimens previously collected from subjects present for Drug Screen Testing. The following results were tabulated:

Clinic Result of Whole Blood				
Method	Results	GC/MS		Total Results
		Positive	Negative	
TCA Rapid Test Cassette	Positive	23	2	25
	Negative	2	63	65
Total Results		25	65	90
% Agreement		92.0%	96.9%	96.6%

Clinic Result of Serum or Plasma				
Method	Results	GC/MS		Total Results
		Positive	Negative	
TCA Rapid Test Cassette	Positive	23	2	25
	Negative	2	63	65
Total Results		25	65	90
% Agreement		92.0%	96.9%	96.6%

Analytical Sensitivity

A drug-free whole blood/serum/plasma pool was spiked with TCA at the following concentrations of ±50% cut-off and 3x cut-off, the data are summarized below:

For whole blood:				
TCA Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0	30	30	0
150	-50%	30	30	0
300	Cut-off	30	15	15
450	+50%	30	0	30
600	3X	30	0	30

For serum or plasma:				
TCA Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0	30	30	0
150	-50%	30	30	0
300	Cut-off	30	15	15
450	+50%	30	0	30
600	3X	30	0	30

Analytical Specificity

The following table lists compounds that are positively detected in Whole