
VANCOGENX[®] - SPACE HIP


TECRES[®]
ADVANCING HIGH TECHNOLOGY

Vancogen[®]-Space Hip & Vancogen[®]-Space Hip XL

Ref	A (mm)	B (mm)	C (mm)	D (mm)	Gentamicin base	Vancomycin base
SPC0030	46	54,5	96	-	1,1 g	1,1 g
SPC0130	54	60	94	-	1,9 g	1,9 g
SPC0230	60	73	98	-	3,0 g	3,0 g
SPC0330	46	54,5	211	10	1,3 g	1,3 g
SPC0430	54	60	209	10,5	2,1 g	2,1 g
SPC0530	60	73	211	11	3,2 g	3,2 g

Vancogen[®]-Space Hip Flat Stem & Vancogen[®]-Space Hip Flat Stem XL

Ref	A (mm)	B (mm)	C (mm)	D (mm)	E (mm)	F (mm)	Gentamicin base	Vancomycin base
SPC0630	46	54,6	96	10	11	9	1,1 g	1,1 g
SPC0730	54	60	94,3	10,5	16	9	1,6 g	1,6 g
SPC0830	60	73	95,8	11	16	9	2,6 g	2,6 g
SPC0930	46	54,6	211	10	11	9	1,2 g	1,2 g
SPC1030	54	60	209,2	10,5	16	9	1,8 g	1,8 g
SPC1130	60	73	211	11	16	9	2,8 g	2,8 g

VANCOGENX[®]-SPACE HIP

EN

INSTRUCTIONS FOR USE

To the Medical Staff's attention

PRESENTATION

Vancogenx[®]-Space Hip is a temporary implantable spacer for the hip joint.

Vancogenx[®]-Space Hip includes and releases two medicinal substances: gentamicin sulfate and vancomycin hydrochloride.

Vancogenx[®]-Space Hip is sterile and single use. The device is sterilized with gamma irradiation.

Vancogenx[®]-Space Hip is available in the following variants:

	Stem profile	Stem length
Vancogenx [®] -Space Hip	Round	Standard
Vancogenx [®] -Space Hip XL	Round	Long
Vancogenx [®] -Space Hip Flat Stem	Rectangular	Standard
Vancogenx [®] -Space Hip Flat Stem XL	Rectangular	Long

DESCRIPTION

Vancogenx[®]-Space Hip is a temporary single use and sterile pre-formed implantable spacer made from fully formed polymethylmethacrylate (PMMA) radiopaque bone cement loaded with gentamicin and vancomycin coating an inner core of stainless steel.

The polymethylmethacrylate (PMMA) bone cement includes barium sulphate as radiopaque agent and gentamicin sulphate and vancomycin hydrochloride as antibiotics.

The formulation of the fully formed resin in the finished product is:

	Formulation [% w/w]
PMMA*	84,90
BaSO ₄	10,00
Gentamicin sulphate	3,20
Vancomycin hydrochloride	1,90

* Polymethyl Methacrylate, N,N-dimethyl-p-toluidine, Hydroquinone, Benzoyl peroxide, Methmethacrylate

The inner structure is in stainless steel AISI 316 ESR, in compliance with the standard ASTM F138.

The one-piece design mimics a "hemi-hip" prosthesis and includes head, neck and distal region.

The design allows motion of the leg in relation to the hip for basic mobility (sitting, standing, walking) under limited weight-bearing conditions (e.g. crutches, walkers, canes).

Dimensions, gentamicin content and vancomycin content of each size and variant of Vancogenx[®]-Space Hip are reported in the tables placed at the beginning of this leaflet.

INTENDED USE & INDICATIONS FOR USE

Vancogenx[®]-space Hip is indicated for temporary use (maximum 180 days) as a total hip replacement (THR) in skeletally mature patients undergoing a two-stage procedure due to a septic process.

The device is inserted into the femoral canal and acetabular cavity following removal of the existing femoral and acetabular implants and debridement. The device is intended for use in conjunction with systemic antimicrobial antibiotic therapy (standard treatment approach to an infection) as defined by physician and/or surgeon.

Vancogenx[®]-space Hip is not intended for use for more than 180 days, at which time it must be explanted and a permanent device implanted or another appropriate treatment performed (e.g. resection arthroplasty, fusion etc.)

It is particularly indicated for operations which present risks of or existing infections caused by organisms susceptible to Gentamicin and/or Vancomycin.

CONTRAINDICATIONS

- The patient's condition is such that a two-stage arthroplasty procedure is contraindicated due to decreased immune response or other relevant systemic clinical conditions.
- Lack of adequate bone structure precludes adequate support of the prosthesis in the proximal femur or acetabular region.
- The procedure is unjustified due to deficiencies in the patient's muscular, nervous or vascular systems.
- Poor bone quality (as in osteoporosis) could cause the prosthesis to migrate or to fracture host bone.
- Infection of the THR cannot be confirmed.
- The infected THR devices cannot be removed.
- Other osteosynthesis devices present in the joint interfere mechanically with the device.
- The infecting pathogens are resistant to gentamicin and vancomycin.
- The patient is sensitive (allergic) to gentamicin, aminoglycosides, vancomycin, glycopeptides or PMMA bone cement.
- A systemic or secondary remote infection is suspected or confirmed.
- The patient does not have a THR and the infection is secondary to trauma or other surgical procedures.
- The patient does not have sufficient bone stock to allow insertion and fixation of the prosthesis
- The patient has an anatomical condition which may favour the acetabular protrusion (interrupted, thin or weak acetabular roof)
- The patient has neuromuscular disorders that do not allow control of the hip joint.
- The patient's weight, age or activity level would cause the surgeon to expect early failure of the system.
- The patient is unwilling or unable to use protected weight bearing mobility throughout the implantation period (e.g. crutches, canes, walkers etc).
- Myasthenia gravis.
- Severe renal impairment.
- Pre-existing calcium metabolism disorder.

INTENDED USERS

The device is intended to be used by suitably qualified orthopaedic surgeons, experienced in two-stage arthroplasty procedures, in an operating theatre environment.

PATIENT TARGET GROUP(S)

The device is intended to be used in skeletally mature patients suffering from a peri-prosthetic joint infection (PJI) of the hip joint.

CLINICAL BENEFITS & PERFORMANCES

Clinical benefits

- The patient can maintain a certain joint mobility, walk with partial weight-bearing (according to surgeon indications based on device stability and clinical/anatomical condition) and perform some physical therapy (active and passive) and muscle strengthening exercises: this will make easier the second-stage and quicker the functional recovery after the second stage since muscle retraction and joint stiffness are avoided,
- The patient can perform basic movements, generally without the need of help: self-independence is a great psychological help; moreover, this may limit the need of an external help (caregiver),
- The first operation (1st stage) is shorter since the device is ready-made, and the second operation is generally facilitated and shorter as well (2nd stage) thanks to the muscle tonicity. This means that also the anaesthesiological risks and any surgery-related risk are reduced,
- Earlier hospital dismissal is generally allowed to the patient.
- The presence of the antibiotics protects from bacterial adhesion and the local release of the antibiotics participates to the treatment of the infection.

Clinical performance

Main

- Maintain the correct joint space between the resected femur and the pelvis following prosthesis removal,
- Allow partial weight-bearing in relation to device stability and anatomical condition,
- Restore a (partial) joint functionality which allows movement and possibility to perform physical therapy, providing a higher quality of life,
- Facilitate the second stage revision as muscle retraction and joint stiffness are avoided,
- Reduce functional recovery times after implantation of the revision prosthesis thanks to the maintenance of a certain muscle tonicity.

Ancillary

- Help to oppose the present infection with the controlled and effective release of antibiotics directly in the infection site all over the implantation period.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCES

The Summary of Safety and Clinical Performances (SSCP) is available on the European Database on medical devices (EUDAMED) at the link <https://ec.europa.eu/tools/eudamed>, within six months after the date of publication of the notice of its full functionality.

ADVERSE EVENTS

The following adverse events may arise when using a spacer. Some are typical of any surgery, others are specific of arthroplasty operations. Most are not directly attributable to the spacer as such; however, surgeons must be aware of these complications and be ready to treat them if they occur.

General Surgical Risks

- Arrhythmias
- Myocardial infarction
- Pulmonary embolism
- Sudden death
- Transitory hypotension
- Venous thrombosis

Hip replacement Surgery Risks

- Acetabular protrusion
- Arthrofibrosis
- Damage to blood vessels
- Damage to femur or acetabulum
- Extensive blood loss
- Hematoma
- Limb length discrepancy
- Nerve damage, bone bed damage
- Phlebitis, thrombophlebitis
- Wound healing problems

Vancogen[®]-Space Hip Devices Risks

- Adverse soft tissue reaction
- Allergic pyrexia
- Debris release
- Device breakage
- Device dislocation
- Device loosening
- Difficult device removal
- Gentamicin-associated toxicity (nephro- / ototoxicity)
- Recurrent infection
- Vancomycin-associated toxicity (nephro- / ototoxicity)

Other adverse reactions associated with systemic use of Gentamicin therapy include hypersensitivity, anaphylactic reactions, nausea, vomiting, urticaria, reversible granulocytopenia, anemia, blood dyscrasia, convulsions, central nervous toxicity, abnormal hepatic function, hypomagnesaemia, stomatitis, purpura, allergic contact sensitization and neuromuscular blockade.

Adverse reactions associated with systemic use of Vancomycin include hypersensitivity, anaphylactic reactions, rash, urticarial, pruritic, toxic epidermal necrolysis.

The postoperative follow-up should be closely monitored by the surgeon for patients reporting a previous experience of sensitivity or allergic reaction to other antibiotics, in particular of the aminoglycosides' family.

Attention: since the devices can/must be fixed with antibiotic loaded bone cement, the surgeon must be aware of its negative effects.

INFORMATION FOR USE

The present instructions for use leaflet must be read carefully.

Familiarity with and attention to Vancogen[®]-Space Hip surgical technique and familiarity with proper use of the Vancogen[®]-Space Hip is essential for successful use of the device.

Only surgeons who have read Vancogen[®]-Space Hip surgical technique and are aware of the limitations of its application should utilize the device.

WARNINGS - PRECAUTIONS FOR USE

Do not use after the expiration date since the effectiveness of the device may be compromised.

Altering the structural identity of the device can cause damages to the same (creation of cracks or fragments and risk of sterility loss leading to infection risks for the patient).

Surgeons must make sure that the device to be implanted is of the right size for the site of implantation.

Remove the device not later than six months after its implantation.

The device was tested to be safely used for not more than 6 months. If this period is extended for too long this can lead to wear, development of debris and eventually to breakage that can cause pain, inflammation and bone re-absorption.

Using the device under conditions other than the indicated use is unlikely to provide benefit to the patient and increases the risk of the development of drug-resistant bacteria.

Do not re-sterilize and/or re-use the device. It is designed for single-use on a single patient.

Re-sterilization should not be carried out since it can cause infection risks for the patient. Re-sterilization can also alter the device's morphology, the efficiency of the antibiotics and of the mechanical features, causing a malfunction of the same with serious risks for the patient's health.

Re-using the device after its extraction must by all means be avoided since it can cause contamination and worsening of the patient's infection. During the extraction, the spacer can also be damaged and residues of the cement can remain in the operatory site: in this case all residues of cement shall be fully removed.

Vancogenx®-Space Hip should be used with caution:

- in conjunction with other nephrotoxic or ototoxic drugs (e.g. other aminoglycosides, amphotericin B, cephalosporins, ciclosporins, cisplatin, clindamycin, colistin, polymyxin B, foscarnet, vancomycin, viomycin) since they can increase gentamicin toxicity (the low serum level of gentamicin makes this interaction very unlikely to occur);
- in patients who are predisposed to or who have pre-existing clinical conditions that would put them at risk for gentamicin and vancomycin toxicity (e.g. renal dysfunction, hearing difficulties, dehydration, advanced age, taking drugs which may affect kidneys, undergoing general anaesthesia, etc).

In such cases, patients should be monitored for toxic blood levels of gentamicin and vancomycin (mainly in the first days following implantation), as well as for kidney function.

Simultaneous use of gentamicin and/or vancomycin and strong diuretics, such as ethacrynic acid or furosemide, may intensify the ototoxic effect of gentamicin and vancomycin by altering the antibiotic concentrations in serum and tissue.

In cases of significant obesity gentamicin serum concentrations should be closely monitored and a reduction in dose should be considered.

Co-administration with the following agents should be avoided:

- Neuromuscular blocking agents such as succinylcholine, botulinum toxin risk of toxicity due to enhanced neuromuscular block
- Anticoagulants such as warfarin and phenindione as Gentamicin has been known to potentiate them
- Cholinergics: antagonism of effect of neostigmine and pyridostigmine.
- Piperacilline and tazobactam as they can potentiate vancomycin nephrotoxicity.
- Bisphosphonates: increased risk of hypocalcaemia.
- Sulfites in susceptible people, especially those with a history of asthma or allergy: risk of allergic-type reactions including anaphylactic symptoms and bronchospasm.

The use of Vancogenx®-Space Hip may result in overgrowth of non-susceptible organisms. If overgrowth of non-susceptible organism occurs, appropriate therapy should be initiated.

The use of antibiotic-loaded bone cement for fixation of the device may increase the potential of toxic drug reactions.

In vitro elution studies (microbiological method) has shown that the daily release of gentamicin and vancomycin never exceeds the recommended systemic adult dose for gentamicin (5-7 mg/kg/day, adult with normal renal function) and vancomycin (30 – 45 mg/Kg/day) according to Goodman and Gilman's recommendations¹. It is therefore unlikely that the amount of gentamicin and vancomycin absorbed locally from Vancogenx®-Space Hip will result in serum levels in the toxic range.

Owing to the presence of gentamicin and vancomycin, adverse events of these substances may occur, even if occurrence is unlikely because of the very low serum levels reached.

Gentamicin and vancomycin are potentially nephrotoxic and/or ototoxic.

Despite the risk of accumulation is low, care should be taken, and the level of gentamicin and vancomycin should be monitored in patients with impaired renal function, and in patients with previous neuromuscular disorders (e.g. myasthenia gravis, Parkinson's disease), in particular if they also have a history of renal insufficiency.

In case of administration of muscle relaxants and ether, the neuromuscular blocking properties of gentamicin may be intensified (using calcium salts it may be reversed). However, this is quite unlikely in patients with healthy kidneys, whereas in patients with impaired renal function, the probability of interaction increases in proportion to serum levels.

Plasma trough concentrations of gentamicin higher than 2µg/ml for periods longer than 10 days have been associated to toxicity.

Ototoxicity is associated to an excessively high concentration of vancomycin in plasma (between 60 and 100 µg/ml). Nephrotoxicity induced by vancomycin is unusual, when normal doses are used. It depends on kidney function and measurement of antibiotic concentration in blood.

USE DURING PREGNANCY, BREAST-FEEDING AND IN CHILDREN

There are no available tests which demonstrate the safe use of the device during pregnancy or breast-feeding, and in children.

During pregnancy and breast-feeding the surgeon should weigh the benefit for the mother against the potential risk to the child before using the device. It is advisable not to implant the device during pregnancy and breast-feeding. In these patients, the device should only be used in critical, life endangering situations.

¹ *The Pharmacological Basis of Therapeutics 13th Edition (Brunton L, Hilal-Dandan R, Knollmann B) - McGraw Hill, New York. Chapter 58: Aminoglycosides (McDougall C), Chapter 59: Protein Synthesis Inhibitors and Miscellaneous Antibacterial Agents (McDougall C)*

In the event that the device is used during pregnancy or if the patient becomes pregnant while the device is in use, she should be notified of the possible dangers to the foetus. Women of childbearing potential should consider the benefits and dangers associated with the use of the product. The device is indicated for use in children (skeletal young patients) only when the surgeon deems that saving the joint through other forms of intervention is not possible.

STERILE PACKAGING INFORMATION

The device is contained inside a double blister packaging configuration and it is sterilised by gamma irradiation.

Please discard the device if the blisters are damaged or unintentionally opened before use because the sterility of the device may be compromised.

Do not re-sterilize.

STORAGE & TRANSPORT CONDITIONS

- Store and transport at a temperature below 25 °C and in a cool and dry environment. However, if the temperature during the transport cannot be maintained under 25°C with a relative humidity < 90%, the medical device can tolerate a temperature between 0 °C and 40 °C with a relative humidity less than 90% for maximum 5 days.
- The device must be shipped and transported in a double wall corrugated box with ECT \geq 44 lbf/in (or ECT \geq 7.71 kN/m).

UTILIZATION, IMPLANTATION and EXPLANATION

It is important to maintain strictly aseptic surgical techniques.

The surgeon shall become thoroughly familiar with the technique of implantation of the prostheses by: (1) reading appropriate literature, (2) training in the operative skills and techniques required for hip arthroplasty revision surgery, and (3) reviewing information regarding use of instrumentation for sizing, implantation and explanation of the component.

In order to prevent dislocation, the same measures taken for a permanent total hip replacement (THR) are advised, plus other specifics such as:

- 1) choice of the correct length of the stem (refer to Section "Vancogen[®]-Space Hip size choice"),
- 2) proximal cementation of the stem (with Vancogen[®] bone cement) to avoid spacer rotation in case of lysis or fragmentation of the proximal femoral bone or stem instability,
- 3) insertion with appropriate joint tension of the soft tissues around the hip joint,
- 4) obtaining adequate head support in the case of severe acetabular bone loss,
- 5) application, in cases at risk, of an orthopaedic abduction brace (possibly articulated) to assist flexion without dislocation.

No specific instrumentation is required to successfully implant the device; however, it is recommended that the following instruments and devices be available for the procedure:

- trial and templates,
- Vancogen[®] bone cement (in case of proximal neck cementation),
- equipment for hip arthroplasty surgery (such as osteotomies, rasps, reamers, femoral head impactor without metallic tips, etc.).

For the application of the device all the routes of access to the hip may be used.

The infected prosthesis components and the residual bone cement (if any) must be removed from the femoral canal and acetabular cavity before the device implantation.

Clean the host site by aggressive debridement, removing all residues.

Reaming of the femoral canal and of the acetabulum should be performed for a better fit of the device on the basis of the patient anatomical condition.

Do not insert the device if the joint hosts other osteosynthesis devices that can interfere mechanically with the device.

Vancogen[®]-Space Hip size choice

Selection of the proper size of the device depends on the judgement of the surgeon with relationship to the requirements of the patient.

For the implant size selection, templates and trial devices are available.

Size appropriateness can be evaluated preoperatively with the superimposition of Templates on patient x-ray imaging.

Size appropriateness should be confirmed through trial reduction, achieving a compromise between stability and joint mobility.

The appropriate size for Vancogen[®]-Space Hip device is based on the fit of the head into the acetabular cavity and the condition of the femur.

Size selection for Vancogen[®]-Space Hip device can be selected by considering:

- Dimensions of the removed implant,
- Remaining bone stock,
- State of ligamentous apparatus,
- Flexion and extension spaces.

a) Head size choice

The head size which best articulates within the acetabulum should be chosen in order to avoid possible dislocation. The basis for measurement could be the removed acetabular cup. A reaming of the remaining acetabular dome, if the residual bone quality permits, may allow the application of a larger head diameter. This may be useful both for the removal of the infected tissue, and for deepening of the spacer head that may help prevent possible dislocation.

b) Stem length choice

A long-stem is advised when distal anchorage is required. This is important in the absence of proximal support, in the presence of large metaphyseal defects or after a trans-femoral approach for implant removal.

Trial placement

Insert the Trial into the femoral canal to verify the stem fit. Once properly seated, reduce the hip joint to determine the correct fit in the acetabular cavity. An adequate fit of the device head into the acetabular cavity will assist in reducing the incidence of dislocation. Once the trial has been used, remove it from the joint, and use it as a reference to select the implantable spacer device size.

Caution: Trials must not be implanted.

Implantation

Before inserting the device, the host site must be irrigated with Ringer or physiological solution and a thorough debridement must be performed. In particular, any residual cement from the previous implant must be removed.

Avoid washing with aqueous solutions the spacer before or after implantation in order to maintain optimal levels of antibiotics release.

Care should be taken to keep the wound dry before the spacer is inserted to avoid loss of antibiotics at the surface of the implant, as the spacer antibiotics are activated in an aqueous environment.

In case of instability, and in order to avoid device rotation and prevent head dislocation, the proximal cementation of the neck with Vancogen® bone cement is mandatory.

Note: Proximal neck cementation with Vancogen® bone cement is mandatory when Vancogen®-Space Hip XL, Vancogen®-Space Hip Flat Stem and Vancogen®-Space Hip Flat Stem XL are used.

Note: full cementation of the stem shall not be performed.

Insert the device completely in the femoral canal.

A head impactor can be used to properly seat the device. However, a metal head impactor may fracture the device. Do not subject the device to excessive forces.

Make sure it is in the correct position by ensuring that the lower protrusion of the device collar rests on the outer diaphysis cortex resected at the inter-trochanteric level.

Turn the device head towards the acetabulum cavity and proceed by reducing the head.

POST-OPERATIVE TREATMENT

As a general rule, post-operative treatment is superimposable with a primary hip prosthesis, with the difference that weight-bearing can be only partial (use of crutches and/or braces). In some clinical-anatomical conditions only touch-down or no weight-bearing can be prescribed by the surgeon.

Partial weight-bearing must be assessed on an individual basis in relation to the anatomic conditions of the femur and acetabulum, bone growth and the clinical conditions of the patient during rehabilitation stages. In particular, one should avoid the risk that excessive weight-bearing or forced mobilisation cause the structure of the spacer to damage the biological structure. If the surgeon deems it necessary, in cases at risk of dislocation, an orthopaedic abduction brace (possibly articulated) to assist flexion may be prescribed.

The physician shall make the patient aware of all the warnings, precautions and measures to be taken and shall provide patient with the filled implant card delivered with the device.

The follow-up may include x-rays and/or removal of fluid from the hip joint to see if the infection is healing. The follow-up may also include drug monitoring tests to make sure the level of drugs in the system is correct. After the infection is healed, a second surgery to remove the device and implant a permanent hip joint implant will be needed.

MRI SAFETY INFORMATION

Non-clinical testing has demonstrated that the Vancogen® -Space Hip and Vancogen®-Space Hip Flat Stem are "MR Conditional" in accordance with the ASTM F2503-20 standard definitions. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Horizontal bore MRI system with a static magnetic field of 1.5 Tesla or 3 Tesla.
- Gradient magnetic fields lower or equal to 19T/m.
- BO*IdBO/drl product lower or equal to 48T²/m.
- RF whole body transmit/receive coil use only.
- Whole body averaged SAR (Specific Absorption Rate) limited to Normal operating mode (WB-SAR ≤2W/kg).
- Vancogen®-Space Hip is expected to produce a maximum temperature rise of **2.1 ± 1.5 °C** at 1.5T and of **1.5 ± 1.1 °C** at 3T.
- Vancogen®-Space Hip Flat Stem is expected to produce a maximum temperature rise of **2.7 ± 2.2 °C** at 1.5T and of **1.5 ± 2.2 °C** at 3T.
- MR image quality may be compromised if the imaging area of interest is in the exact same area of the implant. In vitro assessment at 1.5T showed on spin echo sequences a maximal artifact length of 302 ± 15 mm (including the device) and a maximal artifact width of 18.00 ± 7.38 mm (surrounding the device) and on the corresponding gradient echo sequences, a maximal artifact length of 356 ± 15 mm (including the device) and a maximal artifact width 45.00 ± 7.38 mm (surrounding the device).

IMPLANT CARD

The device is provided with an implant card to be filled by the Healthcare Institution/provider with patient data, and to be supplied to the patient. This implant card reports the information allowing the identification of the implanted device. Further information to be supplied to the patient about warnings, precautions and measures to be taken is reported on the patient leaflet and at the following website: <https://ifu.tecres.com>. The patient leaflet included in the device package must be supplied to the patient.

Instruction for Implant Card completion

The following data shall be filled by the HealthCare Institution/provider on the implant card front side:

1. Name of the patient or patient ID.
2. Date of implantation (day/month/year).
3. Name and address of the Healthcare Institution/provider.

The image shows a rectangular implant card with a white background and a black border. On the left side, there is a vertical column of four icons: a person icon, a calendar icon, a person with a plus sign icon, and a laptop icon. To the right of these icons are four horizontal lines for text entry. Above the first line is a circled number '1', above the second is a circled number '2', and above the third is a circled number '3'. The top left corner features the Tecres logo (a cube) and the text 'TECRES' and 'www.tecres.com'. The top right corner says 'IMPLANT CARD'. On the right side, there is a large, stylized black and white geometric logo. At the bottom left, there is a laptop icon and the URL 'https://ifu.tecres.com'.

Once filled the Implant Card shall be detachment from its support following the dotted lines and provided to the patient.

Explanation of overall symbols present on the implant card (front and back sides) are reported under section "SYMBOLS on Implant Card" of the present leaflet.

EXPLANATION

The device is not intended for permanent implantation and must be removed within 180 days of implantation. Osteotomes, mallets and other revision instruments may be used to aid in the explanation procedure. Care should be taken to assure that the wound site is thoroughly cleaned of all bone cement debris prior to implantation of a definitive prosthesis or performing an alternative surgical procedure (e.g. resection arthroplasty, fusion etc.). Failure to remove cement and/or bone debris may shorten the survival of the revision implant.

NOTE: *It may happen that the device will be left indefinitely in the patient when the patient becomes unfit for further surgery or when the surgeon agrees with the patient under their responsibility that this is the best treatment option in a benefit/risk evaluation.*

DISPOSAL

Vancogenx®-Space Hip once explanted shall be disposed as biological waste according to local waste regulations, since the device could have been contaminated with potentially infectious substances of human origin. In case of expired shelf-life, damaged or unintentionally opened packaging, please dispose Vancogenx®-Space Hip according to local waste regulations for unused or expired drugs since the device contain antibiotics.

SERIOUS INCIDENT REPORTING

Any serious incident potentially attributable to Vancogenx®-Space Hip shall be reported to the manufacturer, Tecres Spa, and to the competent authority in which the patient is established.

The manufacturer can be reached at the following contacts:

- email: info@tecre.it
- tel.: +39 045 9217311

The contact information of the competent authority located in the European Member State in which the patient is established is reported on European Commission website, under section "medical devices" at the following link: https://health.ec.europa.eu/medical-devices-sector/new-regulations/contact_en

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The image shows a rectangular card titled "TECRES IMPLANT CARD". On the left side, there is a logo of a cube with the text "TECRES" and "ADVANCED MEDICAL TECHNOLOGIES" below it. To the right of the logo, the title "TECRES IMPLANT CARD" is printed. Below the title, there are three numbered fields for data entry: 1. A field with a person icon and a question mark, followed by a horizontal line and a circle containing the number 1. 2. A field with a calendar icon showing the number 31, followed by a horizontal line and a circle containing the number 2. 3. A field with a person icon and a plus sign, followed by a horizontal line and a circle containing the number 3. Below these fields, there is a URL "https://ifu.tecres.com" with a small computer icon to its left. On the right side of the card, there is a large, stylized QR code.

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The image shows a rectangular card with a white background and a black border. On the left side, there is a logo for 'TECRES' which consists of a stylized cube. To the right of the logo, the text 'TECRES' is written in a bold, sans-serif font, followed by 'IMPLANT CARD' in a smaller font. Below the logo and text, there are three rows of information, each starting with a small icon and followed by a line for text. The first row has a person icon and a circle containing the number '1'. The second row has a calendar icon and a circle containing the number '2'. The third row has a hospital icon and a circle containing the number '3'. At the bottom left, there is a small icon of a computer monitor and the URL 'https://ifu.tecres.com'. On the right side of the card, there is a large, stylized black graphic that looks like a triangle or a set of lines forming a shape.

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TECRES S.p.A.