<p>BIOGROUP MEDICAL SYSTEM S.r.l. Loc. Campiano 9/b Talamello (RN) P.Iva 00964170419 Divisione Prove Interlaboratorio</p> 	<p>Istruzioni d'uso</p> <p>EMATOLOGIA</p>	<p>I.O. 06 – IFU/E</p>
		<p>Rev. 01</p> <p>Data 07/2022</p>

REF: MSQUALITYE8 - MSQSE812

Contenuto:

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REF: MSEQSE8 – MSEQUALITYE8 - MSEQUALITYE12

Includes:

1 vial, 2.5 ml containing liquid stabilized material at different concentrations.

Intention:

The reagent is used for accuracy testing of haematology analyzers, and also for the testing of accuracy and reproducibility of their pipetting systems.

COMPOSITION:

Contents
<p>Water solution:</p> <p>Potassium bichromate</p> <p>Preservative (sodium azide)</p> <p>E121</p>
Concentration
<p>Vial and lot specific</p> <p>0.05 %</p> <p>0.03%</p>

Safety:

- For IN VITRO diagnostics only. Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.
- Contains sodium azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.
- Sodium azide reacts with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10% sodium hydroxide.

Storage and stability:

Store the kit at 2-8°C in vertical position. **Do not freeze.**

PROCEDURE:

- Leave the material to stand for at least 15 minutes at room temperature before measurement; Do not use mechanical stirrers.
- Roll the tube holding it horizontally between your hands. Carry out the operation with firm movements but without jerking.
- Every 20/30 seconds perform an inversion movement.
- Continue until all the corpuscular part appears homogeneously suspended.
- Before the test, perform 8-10 inversion movements.

Enter the values managed on the **qs-veq.it** website within the expected deadline.

LIMITATIONS

The stabilizing substances and the components used do not make differential analysis possible under a microscope. An incorrect sample preparation for the first reading consequently alters all subsequent ones.

N.B. ATTENTION THE BOTTLES CONTAIN MATERIAL OF HUMAN ORIGIN. IT IS RECOMMENDED TO IMPLEMENT ALL THE RECOMMENDED PROCEDURES FOR THE PURPOSE OF PREVENTION AND DISPOSAL

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