**Disposable Multistage  
Dilation Balloon Catheter**

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FR (Mode d'emploi) ..... 22

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ES (Instrucciones de uso) ..... 52

Instruction page for Electronic instructions for use website ..... 62

## WARNINGS

1. **For single use only.** Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which in turn may result in inpatient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or crossinfection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Micro-Tech assumes no liability with respect to instruments reused, reprocessed or resterilized.
2. This device shall not be used for any purpose other than its intended use.
3. Make sure the endoscope is clear before the instrument is inserted into the endoscope. Do not use the instrument if you cannot see the distal part of the instrument inserted in the endoscope. Unclear vision when the instrument is inserted into the endoscope may result in injury to the patient, such as perforation, hemorrhage or mucosal damage, as well as damage to the endoscope and/or instrument.
4. The device may cause injury, permanent disability or death to the patient if it is not used properly.
5. The operation should be performed under good surveillance and necessary rescue equipment should be prepared in advance to ensure the safety and success of the operation.
6. Contents supplied **STERILE**. Do not use if sterile barrier is damaged. If damage is found contact the manufacturer: Micro-Tech (Nanjing) Co., Ltd.

## DEVICE NAME

Disposable Multistage Dilation Balloon Catheter

# DEVICE DESCRIPTION

## 【CATALOGUE NUMBER】

Table 1 Catalogue number

|             |             |             |             |
|-------------|-------------|-------------|-------------|
| MBD-AM-BB-2 | MBD-BM-BB-2 | MBD-CM-BB-2 | MBD-DM-BB-2 |
| MBD-EM-BB-2 | MBD-FM-BB-2 | MBD-AM-BA-2 | MBD-BM-BA-2 |
| MBD-CM-BA-2 | MBD-DM-BA-2 | MBD-EM-BA-2 | MBD-FM-BA-2 |
| MBD-AL-BA-2 | MBD-BL-BA-2 | MBD-CL-BA-2 | MBD-DL-BA-2 |
| MBD-EL-BA-2 | MBD-FL-BA-2 | MBD-FL-BB-2 | MBD-AS-BB-2 |
| MBD-BS-BB-2 | MBD-CS-BB-2 | MBD-DS-BB-2 | MBD-AS-BA-2 |
| MBD-BS-BA-2 | MBD-CS-BA-2 | MBD-DS-BA-2 |             |

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The specifications of the product are represented by the Balloon OD, Balloon Length and Catheter Length, as indicated in Table 2.

Table 2 Specification and parameter

| Specifications | Balloon OD:<br>Φ (mm) | Dilation<br>pressure:<br>P1-P2-P3<br>(atm) | Balloon<br>length:<br>(mm) | Catheter<br>length:<br>(mm) | Compatible Guidewire | Catheter<br>OD: (Fr) | Maxium<br>Pressure<br>(atm) | Working<br>Channel<br>(mm) |
|----------------|-----------------------|--------------------------------------------|----------------------------|-----------------------------|----------------------|----------------------|-----------------------------|----------------------------|
| MBD-AM-BB-2    | 6-7-8                 | 3-6-10                                     | 55                         | 2300                        | ≤0.035inch(0.89mm)   | 7                    | 11                          | ≥2.8                       |
| MBD-BM-BB-2    | 8-9-10                | 3-5.5-9                                    | 55                         | 2300                        | ≤0.035inch(0.89mm)   | 7                    | 10                          | ≥2.8                       |
| MBD-CM-BB-2    | 10-11-12              | 3-5-8                                      | 55                         | 2300                        | ≤0.035inch(0.89mm)   | 7                    | 9                           | ≥2.8                       |
| MBD-DM-BB-2    | 12-13.5-15            | 3-4.5-8                                    | 55                         | 2300                        | ≤0.035inch(0.89mm)   | 7                    | 9                           | ≥2.8                       |
| MBD-EM-BB-2    | 15-16.5-18            | 3-4.5-7                                    | 55                         | 2300                        | ≤0.035inch(0.89mm)   | 7                    | 8                           | ≥2.8                       |
| MBD-FM-BB-2    | 18-19-20              | 3-4.5-6                                    | 55                         | 2300                        | ≤0.035inch(0.89mm)   | 7                    | 7                           | ≥2.8                       |
| MBD-AM-BA-2    | 6-7-8                 | 3-6-10                                     | 55                         | 1800                        | ≤0.035inch(0.89mm)   | 7                    | 11                          | ≥2.8                       |
| MBD-BM-BA-2    | 8-9-10                | 3-5.5-9                                    | 55                         | 1800                        | ≤0.035inch(0.89mm)   | 7                    | 10                          | ≥2.8                       |
| MBD-CM-BA-2    | 10-11-12              | 3-5-8                                      | 55                         | 1800                        | ≤0.035inch(0.89mm)   | 7                    | 9                           | ≥2.8                       |
| MBD-DM-BA-2    | 12-13.5-15            | 3-4.5-8                                    | 55                         | 1800                        | ≤0.035inch(0.89mm)   | 7                    | 9                           | ≥2.8                       |
| MBD-EM-BA-2    | 15-16.5-18            | 3-4.5-7                                    | 55                         | 1800                        | ≤0.035inch(0.89mm)   | 7                    | 8                           | ≥2.8                       |
| MBD-FM-BA-2    | 18-19-20              | 3-4.5-6                                    | 55                         | 1800                        | ≤0.035inch(0.89mm)   | 7                    | 7                           | ≥2.8                       |
| MBD-AL-BA-2    | 6-7-8                 | 3-6-10                                     | 55                         | 1800                        | ≤0.035inch(0.89mm)   | 7                    | 11                          | ≥2.8                       |
| MBD-BL-BA-2    | 8-9-10                | 3-5.5-9                                    | 55                         | 1800                        | ≤0.035inch(0.89mm)   | 7                    | 10                          | ≥2.8                       |
| MBD-CL-BA-2    | 10-11-12              | 3-5-8                                      | 80                         | 1800                        | ≤0.035inch(0.89mm)   | 7                    | 9                           | ≥2.8                       |
| MBD-DL-BA-2    | 12-13.5-15            | 3-4.5-8                                    | 80                         | 1800                        | ≤0.035inch(0.89mm)   | 7                    | 9                           | ≥2.8                       |
| MBD-EL-BA-2    | 15-16.5-18            | 3-4.5-7                                    | 80                         | 1800                        | ≤0.035inch(0.89mm)   | 7                    | 8                           | ≥2.8                       |

|             |            |         |    |      |                                        |   |    |            |
|-------------|------------|---------|----|------|----------------------------------------|---|----|------------|
| MBD-FL-BA-2 | 18-19-20   | 3-4.5-6 | 80 | 1800 | $\leq 0.035\text{inch}(0.89\text{mm})$ | 7 | 7  | $\geq 2.8$ |
| MBD-FL-BB-2 | 18-19-20   | 3-4.5-6 | 80 | 2300 | $\leq 0.035\text{inch}(0.89\text{mm})$ | 7 | 7  | $\geq 2.8$ |
| MBD-AS-BB-2 | 6-7-8      | 3-6-10  | 30 | 2300 | $\leq 0.035\text{inch}(0.89\text{mm})$ | 7 | 11 | $\geq 2.8$ |
| MBD-BS-BB-2 | 8-9-10     | 3-5.5-9 | 30 | 2300 | $\leq 0.035\text{inch}(0.89\text{mm})$ | 7 | 10 | $\geq 2.8$ |
| MBD-CS-BB-2 | 10-11-12   | 3-5-8   | 30 | 2300 | $\leq 0.035\text{inch}(0.89\text{mm})$ | 7 | 9  | $\geq 2.8$ |
| MBD-DS-BB-2 | 12-13.5-15 | 3-4.5-8 | 30 | 2300 | $\leq 0.035\text{inch}(0.89\text{mm})$ | 7 | 9  | $\geq 2.8$ |
| MBD-AS-BA-2 | 6-7-8      | 3-6-10  | 30 | 1800 | $\leq 0.035\text{inch}(0.89\text{mm})$ | 7 | 11 | $\geq 2.8$ |
| MBD-BS-BA-2 | 8-9-10     | 3-5.5-9 | 30 | 1800 | $\leq 0.035\text{inch}(0.89\text{mm})$ | 7 | 10 | $\geq 2.8$ |
| MBD-CS-BA-2 | 10-11-12   | 3-5-8   | 30 | 1800 | $\leq 0.035\text{inch}(0.89\text{mm})$ | 7 | 9  | $\geq 2.8$ |
| MBD-DS-BA-2 | 12-13.5-15 | 3-4.5-8 | 30 | 1800 | $\leq 0.035\text{inch}(0.89\text{mm})$ | 7 | 9  | $\geq 2.8$ |

### 【STRUCTURE】

The Disposable Multistage Dilation Balloon Catheter is composed of Handle Junction, Dilating Balloon, Guidewire, Guidewire Locker and other components (See Figure 1).


Figure 1 Product Drawing of Disposable Multistage Dilation Balloon Catheter  
 1.Flexible tip 2.Dilating Balloon 3.Inner Sheath 4.Marker Band  
 5.Handle Junction 6.Capture Tube 7.1-way Stopcock  
 8. Stainless steel guidewire 9. Guidewire Locker

Note: This is a sketch map of Disposable Multistage Dilation Balloon Catheter, it is for reference only, as the actual device may vary.

### 【USER INFORMATION/TRAINING/QUALIFICATIONS】

This device is intended for use by physicians or medical personnel who have received appropriate training, in accordance with the applicable laws. As a consumable of digestive endoscope, the product shall be used by professionals familiar with the operation technique of digestive endoscope.

**【RADIOPACITY INFORMATION】**

The Device contains Radiopaque Marker Band and is visible on X-ray.

**INTENDED USE/INDICATIONS**

The device is indicated for use to dilate strictures of the gastrointestinal tract.

**PATIENT TARGET GROUP**

This device is intended to be used for adults and adolescents.

**CONTRAINDICATIONS**

1. Inability to pass the endoscope or tolerate the endoscopy;
2. Severe cardiopulmonary disease;
3. Severe coagulopathy;
4. Known or suspected perforation;
5. Acute abdomen;
6. Deeply ulcerated stenosis;
7. Recent esophageal surgery;
8. Large aortic aneurysm;
9. Colonic strictures with acute diverticulitis;
10. Allergic reaction to contrast medium;
11. Any other conditions the physician judges unsuitable for use.

**COMPLICATIONS**

- ◆ Perforation
- ◆ Hemorrhage
- ◆ Hematoma
- ◆ Swelling/Edema
- ◆ Infection
- ◆ Pain
- ◆ Mucosal laceration
- ◆ Aspiration
- ◆ Inflammation
- ◆ Complications which are not currently known or observed may be present

## **CLINICAL BENEFIT**

---

In the endoscopic view, the Disposable Multistage Dilation Balloon Catheter reaches the stricture position of digestive tract. Then dilate the balloon under determined pressure, the expansion force drives the stricture position to recover to normal. After procedure, the patient's gastrointestinal function recovered, and the patient's quality of life improved. The operation is simple, and less invasive.

## **PERFORMANCE CHARACTERISTICS**

---

1. Three different diameters with three separate pressure levels.
2. 2 X-ray markers on opposite ends to indicate the exact working length of the balloon.
3. Rounded Balloon shoulders for clear endoscopic view of the balloon's distal end.
4. Conical, atraumatic catheter tip.
5. Preloaded stiffed wire with hydrophilic, radiopaque tip.

## **CAUTIONS**

---

1. A thorough understanding of the technical principles, clinical applications and associated risks are required prior to the use of the device.
2. The patient shall be informed and respond with his/her consent regarding the details of the operation as well as the potential risks and complications, which may lead to injury, illness or death.
3. Please read the instruction carefully before using it.
4. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## **NOTES**

---

A single copy of the instructions for use is provided in a case box. If further copies are requested, please contact Micro-Tech, we will provide free of charge.

## **HOW SUPPLIED**

---

1. The device packaging is intact and sterilized.
2. The device is intended for single use only.

## **STORAGE**

---

This device should be stored in a clean, well-ventilated, non-corrosive gas environment with a service life of 3 years.

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1. Refer to package label for the minimum channel size required, the recommended working pressure and other specifications.
2. Refer to package label for the expiration date in order to ensure that the device is not used after such date.
3. Do not pre-inflate the balloon prior to its insertion.
4. If the preload guidewire in dilation balloon catheter is bent or the guidewire tip is frayed, please do not use directly, it needs to be removed from dilation balloon catheter, and then the device matches a pre-positioned guidewire to use.
5. If the package is opened or damaged upon receipt, do not use.
6. Check for kinks (preload guidewire or shaft) and damages. If an abnormality is found, contact Micro-Tech for further assistance.
7. Inflate the balloon only with sterile water, sterile saline or contrast medium mixture.
8. If the deflated balloon or the preloaded guidewire cannot be withdrawn from the working channel after the procedure, removal of the endoscope along with the deflated balloon and the preloaded guidewire may be necessary in order to prevent damages to the endoscope. After removal, withdraw the balloon sheath by cutting the balloon at its proximal end, as shown in Figure 4.

### 【INSTRUCTIONS FOR USE】

#### Instructions For Use of Disposable Multistage Dilation Balloon Catheter

Note: Use guidewire with 0.035 inch (0.89mm).

1. Choose the appropriate endoscope, guidewire (see label for specifications) and the pressure gauge with the required pressure range.
2. Choose the dilation balloon of appropriate size. Remove the protective sheath from the balloon. The dilation balloon comes with a preloaded guidewire. The guidewire can be used as is to access the stricture or removed from the balloon catheter prior to use and back loaded over a pre-positioned guidewire. Prior to advancing the balloon catheter over a pre-positioned guidewire flush the guidewire lumen. (See Figure 2)


Figure 2

3. Apply a lubricating agent to the balloon surface to facilitate passage through the endoscopy working channel.

4. Attach the balloon inflation lumen to an inflation device with pressure gauge of appropriate range.
5. Insert the balloon into the endoscope working channel, and advance it until it gets through the stricture. Continue to advance the balloon until the desired position is achieved. The position can be verified through fluoroscopy with the help of the radiopaque markers (See Figure3).


Figure 3

**WARNING:**

In order to minimize balloon volume, keep the balloon under negative pressure throughout the insertion process. Do not forcefully advance the balloon catheter if any resistance is encountered. Evaluate the cause of resistance before continuing to advance the device.

6. Once the balloon's position is confirmed, inflate the balloon with the inflation device. Adjust the pressure of the balloon according to the recommended working pressure on the package label. Larger diameters can be achieved with higher inflation pressures based on the pressure/-diameter tag label near the hub on the proximal end of the catheter. The balloon may burst if excessive pressure is applied.

**WARNING:**

If the balloon bursts during the procedure, stop immediately and take remedial action.

7. After the operation completely deflate the balloon. Ensure that all liquid is withdrawn from the balloon by maintaining negative pressure until the balloon is fully withdrawn from the working channel.

**WARNING:**

Do not forcefully withdraw the balloon or the preloaded guidewire from the endoscope if significant resistance is encountered. Removal of the endoscope along with the compromised balloon and the preloaded guidewire may be required in order to prevent damages to the endoscope. Remove the balloon catheter from the endoscope by cutting the proximal end of the balloon (See Figure 4).


Figure 4

## PRODUCT DISPOSAL

After use, dispose of product and packaging in accordance with hospital, administrative, and/or local government policy.

## POST-PROCEDURE

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

## COMPATIBILITY

### 【Applicable Endoscopes】

Endoscope which is legal listing in European Union is recommended, working channel of Endoscope  $\geq 2.8\text{mm}$

### 【Compatible Guidewire Diameter】

Guidewire Diameter  $\leq 0.035''$

## SYMBOL INDICATION

|                                                                                   |                                                                         |                                                                                   |                                            |
|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------|-----------------------------------------------------------------------------------|--------------------------------------------|
|  | Do not re-use                                                           |  | Do not re-sterilize                        |
|  | Date of manufacture                                                     |  | Manufacturer                               |
|  | Use-by date                                                             |  | Caution                                    |
|  | Catalogue number                                                        |  | Batch code                                 |
|  | Consult instructions for use or consult electronic instructions for use |  | Sterilized using ethylene oxide            |
|  | Authorized representative in the European Community                     |  | Keep away from sunlight                    |
|  | Keep dry                                                                |  | Do not use if package is damaged           |
|  | Medical device                                                          |  | Sterile barrier system / Sterile packaging |
|  | Importer                                                                |  | Not made with natural rubber latex         |
|  | Unique device identifier                                                |  | Open here                                  |

**【Packaging】** Packed in flexible peel pouch

**【Production Date】** See label

**【Sterilization】** Sterilized by EO (ethylene oxide) gas

**【Period of Validity】** 3 Years

## WARRANTY

EN

Limited Warranty to Buyer. Micro-Tech warrants to Buyer that, for the earlier of one (1) year from the date of purchase, or until the product is used by Buyer, the products will be free from defects in materials and workmanship when stored and used in accordance with the instructions for storage and use provided by Micro-Tech and in accordance with applicable regulatory requirements. Descriptions or specifications appearing in Micro-Tech's literature are meant to generally describe the products and do not constitute any express warranties. Any technical advice with respect to the product and guarantee of specific properties of or in the products shall only be effective if and to the extent specifically confirmed by Micro-Tech in writing. These warranties shall not apply for product failure or deficiency due to improper storage, alteration, or the consequences of uses for which the products were not designed or that adversely affect the products' integrity, reliability, or performance.

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**MICRO-TECH Europe GmbH**
**Shanghai International Holding Corp. GmbH (Europe)**
**Lotus Global Co., Ltd.**