



[www.biocare.com.cn](http://www.biocare.com.cn)



# iM 12&iM 15

## Technical Specification

## Appendix B Product Specifications

### B.1 Safety Specifications

#### B.1.1 Product Classification

For classification of this series of monitors comply with IEC60601-1, please refer to Table B.1.

Table B.1 Module Classification

Components	Type of Protection Against Electric Shock	Degree of Protection Against Electric Shock	Degree of Protection Against harmful ingress of water	Degree of Protection Against hazards of Explosion	Mode of Operation
Main unit	I	Not marked	IPX2	Not suitable	Continuous Operation
ECG (Resp) Module	NA	CF(*)			
IBP Module (Optional)					
NIBP Module					
Temp Module					
SpO <sub>2</sub> Module					
CO <sub>2</sub> Module (Optional)		BF(*)			
AG Module (Optional)					

#### ATTENTIONS:

- I: Class I equipment
- BF: Type BF applied part (The symbol ‘\*’ indicates the availability of defibrillation-proof function).
- CF: Type CF applied part (The symbol ‘\*’ indicates the availability of defibrillation-proof function).
- NA: Not applicable.
- IPX2: degree of protection against harmful ingress of water or particulate matter.
- Not suitable: the equipment is not suitable for use in an environment with air, oxygen or nitrous oxide mixed with flammable anesthetic gas.

**B.1.2 Environment Specifications**

Equipment Environment			
Item	Temperature	Humidity	Atmospheric Pressure
Operating	0℃～40℃ (32°F～104°F) If the machine includes CO <sub>2</sub> module, the operating temperature is 5 ℃～ 40 ℃ (41 ℉～104 ℉)	15%～80%, Non-Condensing	442.5 mmHg～805.5 mmHg (59 kPa～107.4 kPa)
Storage&Transport	-20℃～+55℃ (-4°F～140°F)	≤93%, Non-Condensing	525 mmHg～795 mmHg (70 kPa～106 kPa)

**B.1.3 Power Specifications**

(AC) Input Voltage	100 V～240 V
Input Power	160 VA
Frequency	50 Hz/60 Hz (Allowable frequency error ±1Hz)
Fuse	3.15A/250V

**B.2 Physical Specifications**

Host	12 inches monitor	15 inches monitor
Weight	About 4.5 kg	About. 5.5 kg
Size (L×W×H)	310 mm×163 mm×285 mm	370 mm×187 mm×313 mm

**B.3 Hardware Specifications**

Display	
Type	TFT LCD Screen
Dimensions	12.1 inches (12 inches monitor), 15 inches (15 inches monitor)
Resolution	800×600 pixels (12 inches monitor), 1024×768 pixels (15 inches monitor)
Screen Brightness	10-level, adjustable
LCD View Angle	Horizontal / vertical view angle at least 150°/120°
Recorder	
Type	Thermal array recorder
Horizontal Resolution	16 dots/mm (Paper Speed: 25.0 mm/s)
Vertical Resolution	8 dots/mm
Printing Paper Size	50 mm×20 mm
Paper Speed	12.5 mm/s; 25.0 mm/s; 50.0 mm/s
Waveform	Max. 3 waveforms
Battery	
Dimensions	182 mm×71 mm×25.5 mm
Weight	0.3 kg
Type	Rechargeable lithium battery
Rated voltage	14.8 V
Battery Capacity	4.4 Ah
Length of Power Supply	In environment temperature 25 ℃ and in standard configuration (the SpO <sub>2</sub> sensor connects, the ECG cable and Temp cable disconnect, the “Measure Mode” of NIBP is

	“Auto” and the “Interval” is 15 minutes), the continuous working time of a single battery is not less than 5 hours.
Time for recharging battery to 90% from zero power state	The charging time is not more than 12 hours to charge the battery to 90%.
Shutdown Delay	0 s, 0.5 s, 1 s, 1.5 s, 2 s
<b>Host LED</b>	
Physiological Alarm Indicator Lamp	1 (Dual color yellow & red)
Technical Alarm Indicator Lamp	1 (Blue)
Power Switch Indicator Lamp	1 (Green)
AC Power Indicator Lamp	1 (Green)
Battery Power Indicator Lamp	1 (Green)
Battery Charging Indicator Lamp	1 (Green) (Only for 12 inches monitor)
Keypad Backlight	5 (White)
Alarm Pause Key Backlight	1 (Red)
Speaker	Give out alarm sound (45 dB~85 dB), keystroke sound and QRS sound. Alarm sound complies with IEC 60601-1-8
<b>Interface</b>	
Power	1 AC power port
Network	Standard RJ45 network port, which can network with the central monitoring system and transmit all the animal Monitored data to the central monitoring system.
USB	USB disk supported. For the manufacturer to upgrade and service the application software, and export data (Structurally 2 USB host interfaces supported)
VGA	Supported, for connection of external display
Analog Output Port	1 piece. It can be connected to oscilloscope for output of the analog signals.
Nurse Call System Interface and Defibrillation Synchronization Interface	1 piece. It can be connected to port of the nurse call system or the defibrillator.
Equipotential Terminal Port	1 piece
<b>ECG Analog Signal Output</b>	
Bandwidth (-3 dB, reference 10Hz)	Surgery mode: 1 Hz~15 Hz Monitor mode: 0.5 Hz~40 Hz Diagnose mode: 0.05 Hz~150 Hz
Max. Transmission Delay	25ms (Wave filter closed under diagnose mode)
Sensitivity	1 V/mV $\pm 5\%$
Accuracy of input signal reproduction	Using the method described in 4.2.7.1 of AAMI EC11 to test the overall system error, which is within $\pm 5\%$ ; Using method A and D described in 4.2.7.1 of AAMI EC11 to test frequency response. Because of sampling characteristics and the asynchronism between sample rate and signal rate of the ECG module, digital systems may produce a noticeable modulating effect from one cycle to the next. This phenomenon, which is not physiologic, shall be clearly described in the operator's and service manuals.
<b>IBP Analog Signal Output</b>	
Bandwidth (-3 dB, reference 10Hz)	0 Hz~50 Hz
Max. Transmission Delay	30 ms (Filter closed)
Sensitivity	0.01 V/mmHg $\pm 5\%$

## B.4 Data Storage

Trend Data	Short Trend (Trend Window Time 4 min, 40 min, 2 h) Resolution of Trend Chart 5 s, 30 s, 1 min, 10 min): Max. storage time: 72h. Long trend (Trend Window Time 4 h, 16 h, 32 h, 48 h) Resolution of Trend Chart 15 min, 30 min, 1 h, 2 h, 3 h): Max. storage time: 480h.
Parameter Alarm Event	700 parameter alarm events and manual events, as well as the parameter waveform related to the occurring time, wave length 10s
NIBP Measuring Result	Max. 1000 groups
Single-Channel ECG Waveform	Max. 2h
Holographic Waveform	Max. 2 min (Power cutoff storage not supported)

## B.5 Wireless Network

Applicable Standard	IEEE 802.11b/g/n (2.4G)	IEEE 802.11a/n (5G)
Frequency Range	2.412 GHz~2.472 GHz	4.9 GHz~5.975 GHz
Band Width	20~40MHz	20~40MHz
Radiated Power	+18dBm	+13.5dBm
Signal Path	1-13 (China)	
Type and Frequency Characteristics of the Modulation	CCK/DSSS/OFDM/MCS7/MCS0	

## B.6 Measuring Specifications

### B.6.1 ECG Monitoring

Input Mode	3-Lead ECG input (Optional) 5-Lead ECG input (Standard)
Lead Selection	I , II , III (Optional) I , II , III, aVR, aVL, aVF, V I , II , III, aVR, aVL, aVF, V1~V6 (Optional)
Lead Standard	AHA, IEC
Measuring Range of Heart Rate	>20 kg: 15 bpm~300 bpm 10~20 kg: 15 bpm~350 bpm <10 kg: 15 bpm~350 bpm
Heart Rate Display Tolerance	±1% or ±1 bpm, whichever is higher
Sensitivity	1.25 mm/mV (×1/8), 2.5 mm/mV (×1/4), 5.0 mm/mV (×1/2), 10.0 mm/mV (×1), 20.0 mm/mV (×2), 40.0 mm/mV (×4), Auto. Error: ±5%
Resolution Stability	The resolution change 1 minute after the instrument is powered on does not exceed 0.66% per minute. The total change within 1h does not exceed any available fixed gain setting by ±10%.
Sweep Speed	6.25 mm/s, 12.5 mm/s, 25.0 mm/s, 50.0 mm/s. Error: ±10%
Noise Level	≤30 μV <sub>p-p</sub>
Input Circuit Current	≤0.1 μA
Input Impedance	≥2.5 MΩ
Patient Leakage Current	< 10μA

ESU Proof	Cutting Mode: 300 W Coagulation Mode: 100 W Recovery Time: $\leq 10$ s
ESU Noise Inhibition	Tested acc. to 5.2.9.14 of ANSI/AAMI EC 13:2002: 1) The ECG signal track does not disappear; 2) Change in heart rate does not exceed 10% of the heart rate when the electrosurgical knife is not activated.
CMRR	Diagnose Mode: $\geq 89$ dB Surgery & Monitor Mode: $\geq 100$ dB
Time Constant	Monitor Mode: $\geq 0.3$ s Diagnose Mode: $\geq 3.2$ s
Frequency Response	Surgery Mode: 1 Hz-15 Hz; Monitor Mode: 0.5 Hz-40 Hz; Diagnose Mode: 0.05 Hz-150 Hz.
ECG Parameter Frequency Characteristics	Surgery Mode: Meet (+0.4 dB ~ (-3.0 dB)) requirements at 15 Hz. Monitor Mode: Meet (+0.4 dB ~ (-3.0 dB)) requirements at 0.5 Hz ~ 40 Hz. Diagnose Mode: Meet (+0.4 dB ~ (-1.0 dB)) requirements at 0.05 Hz ~ 60 Hz. Meet (+0.4 dB ~ (-3.0 dB)) requirements at 61 Hz ~ 150 Hz.
Notch	Monitor & Surgery Mode: notch filter automatically activated at 50 Hz/60 Hz Diagnose Mode: Notch filter manually activated or deactivated at 50 Hz/60 Hz
Range of Electrode Polarized Voltage	$\pm 300$ mV d.c.
Lead Fall Testing Current	Measuring Electrode: $< 0.1$ $\mu$ A Drive Electrode $< 1$ $\mu$ A
<b>Pacemaker Pulse</b>	
Pacemaker Pulse Display Capacity	Pace-making mark can be displayed for the following pacemaker pulses: Pulse Amplitude: $\pm 2$ mV ~ $\pm 100$ mV Pulse Width: 0.1 ms ~ 2 ms Pulse Rise Time: 10 $\mu$ s ~ 100 $\mu$ s Pacemaker pulse should be no overshoot
Pacemaker Pulse Suppression Capacity	The monitor can inhibit the pacemaker pulse that conforms to the following conditions: Pulse Amplitude: $\pm 2$ mV ~ $\pm 100$ mV Pulse Width: 0.1 ms ~ 2 ms Pulse Rise Time: 10 $\mu$ s ~ 100 $\mu$ s Pacemaker pulse should be no overshoot
<b>Alarm Limit Specifications</b>	<b>Range</b>
Upper Limit of ECG Heart Rate	Alarm upper limit for $> 20$ kg: (Lower limit+2) bpm ~ 300 bpm Alarm upper limit for 10 ~ 20 kg: (Lower limit+2) bpm ~ 350 bpm Alarm upper limit for $< 10$ kg: (Lower limit+2) bpm ~ 350 bpm
Lower Limit of ECG Heart Rate	Alarm lower limit for $> 20$ kg: 15 bpm ~ (Upper limit-2)bpm Alarm lower limit for 10 ~ 20 kg: 15 bpm ~ (Upper limit-2)bpm Alarm lower limit for $< 10$ kg: 15 bpm ~ (Upper limit-2)bpm
Resolution	$\pm 1$ bpm
Accuracy	The tolerance of alarm limit setting is $\pm 1$ bpm. In addition, the ECF signal alarm below the publicized lower limit of the alarm will not fail. If the alarm is not disabled, the alarm will not fail if you enter the ECG input signal higher than the upper limit of alarm up to 300 bpm (350 bpm for $< 10$ kg and 10 ~ 20 kg).
<b>HR</b>	
Heart Rate Testing Amplitude	$\pm 0.3$ mV ~ $\pm 5$ mV
Resolution	1 bpm
Alarm Time for Tachycardia	Acc. to ANSI/AAMI EC13:2002 Part 4.1.2.1 g).

	4ah-Range: 11 s 4a-Range: 11 s 4ad-Range: 11 s 4bh-Range: 11 s 4b-Range: 11 s 4bd-Range: 11 s
Heart Rate Average	Acc. to ANSI/AAMI EC13:2002 Part 4.1.2.1 g). The average heart rate is obtained by the method below: If the interval of the last continuous 3 RR is higher than 1200ms, the heart rate is averaged based on the most recent 4 RR intervals; otherwise, the heart rate is averaged based on the most recent 12 RR intervals. The heart rate displayed on the screen is refreshed every second.
Response to Irregular Rhythm of the heart	Acc. to ANSI/AAMI EC13:2002 Part 4.1.2.1 e). The heart rate displayed after 20s stabilizing period is: 3a (Ventricular bigeminy) $\sim 80\pm 1\text{bpm}$ 3b (Slow alternating ventricular bigeminy) $\sim 60\text{ bpm}\pm 1\text{ bpm}$ 3c (Rapid alternating ventricular bigeminy) $\sim 120\text{ bpm}\pm 1\text{ bpm}$ 3d (Bidirectional systoles) $\sim 90\text{ bpm}\pm 6\text{ bpm}$
Response Time to Heart Rate Change	Acc. to ANSI/AAMI EC13:2002 Part 4.1.2.1 f). Increase of heart rate: response time $\leq 11\text{ s}$ Decrease of heart rate: response time $\leq 11\text{ s}$
High T-wave Suppression Capacity	Acc. to ANSI/AAMI EC13:2002 Part 4.1.2.1 c). The heart rate monitor inhibits all T-waves with amplitude lower than 1.2 mV, 100msQRS wave groups, T-wave period 180 ms and QT period 350ms.
Arrhythmia Type	a) Monitoring type: Asystole, VFib/VTac, VTac, Ventricular bradycardia, Extreme-Tachy, Extreme-Brady, Non-Sustained VT, PVC, Tachycardia, Bradycardia, VR(ventricular rhythm), V-Bigeminy, V-Trigeminy, Irr.Rhythm, PVCs/min, Run PVCs>2, Couplet, R on T, Multiform, HeartBeat Pause, Missed Beats b) Pace-making: Pacemaker not captured (PNC), Pacemaker not paced (PNP).
<b>ST Interval Measuring</b>	
Range	(-2.0 mV)~(+2.0 mV)
Accuracy	Measuring Tolerance: measuring tolerance within (-0.8 mV)~(+0.8 mV) is $\pm 0.02\text{ mV}$ or $\pm 10\%$ , whichever is higher. It not defined for other ranges.
ST Interval Updating Interval	A single heart beat interval or 1s, whichever is higher.

### B.6.2 Respiration (Resp) Monitoring

Measuring Method	Chest Impedance Method
Measuring Lead	Lead I and II for selection. Lead I defaulted.
Respiration Exciting Waveform	$< 300\text{ }\mu\text{A}$ , Sine signal, 62.8 kHz ( $\pm 10\%$ )
Range of Respiration Impedance	$0.5\text{ }\Omega \sim 3\text{ }\Omega$
Range of Base Impedance	$250\text{ }\Omega \sim 2000\text{ }\Omega$ (Use of ECG cable with $1\text{ k}\Omega$ resistor)
Differential Input Impedance	$> 2.5\text{ M}\Omega$
Bandwidth	$0.2\text{ Hz} \sim 2\text{ Hz}$ (-3 dB)
Waveform Sensitivity	$\times 1/4$ , $\times 1/2$ , $\times 1$ , $\times 2$ , $\times 4$ , Auto
Sweep Speed	$6.25\text{ mm/s}$ ; $12.5\text{ mm/s}$ ; $25.0\text{ mm/s}$
Resolution	1 rpm
Accuracy	$\pm 2\text{ rpm}$
Asphyxia Alarm	Off, 10 s, 15 s, 20 s, 25 s, 30 s, 35 s, 40 s

RR	
Range	Monitoring Range for >20 kg: 0 rpm~120 rpm Monitoring Range for 10 ~20 kg: 0 rpm~150 rpm Monitoring Range for <10 kg: 0 rpm~150 rpm
Resolution	1 rpm
Respiration Monitoring Tolerance	Within 7 rpm~150 rpm, the measuring error is $\pm 2$ rpm or $\pm 2\%$ , whichever is higher. The tolerance is not defined for other ranges.
Asphyxia Alarm Tolerance	Within 10 s~40 s (Increase/decrease by 5s for each rotation of the knob), the asphyxia alarm tolerance is $\pm 5$ s.
Alarm Limit Specifications	Range
RR Upper Limit	Alarm upper limit for >20 kg: (Lower limit+2) rpm ~100 rpm Alarm upper limit for 10 ~20 kg: (Lower limit+2) rpm ~100 rpm Alarm upper limit for <10 kg: (Lower limit+2) rpm ~100 rpm
RR Lower Limit	Alarm lower limit for >20 kg: 0 rpm ~ (Upper limit-2) rpm Alarm lower limit for 10 ~20 kg: 0 rpm ~ (Upper limit-2) rpm Alarm lower limit for <10 kg: 0 rpm ~ (Upper limit-2) rpm

### B.6.3 SpO<sub>2</sub> Monitoring

Alarm Limit Specifications	Range
SpO <sub>2</sub> Upper Limit	(Lower limit+1)% ~100%
SpO <sub>2</sub> Lower Limit	80% ~ (Upper limit-1)%
Accuracy Tolerance	$\pm 1\%$ of the setting
Sensing element	Optical power <15 mW Red light wavelength: 658 nm~664 nm, infrared light: 897 nm~915 nm Information on the wavelength range is particularly useful for clinicians (e.g. in optical dynamic therapy)

### SpO<sub>2</sub> Module

Monitoring Parameters	SpO <sub>2</sub> and Pulse Rate (PR)
Range	0%~100%
Resolution	1%
Data update period	1 s
Accuracy	Within 70%~100%, the measuring tolerance is $\pm 2\%$ . Within 0%~69%, the measuring tolerance is not defined.

### Masimo Oximeter Module

Monitoring parameter	Pulse oximetry (SpO <sub>2</sub> ) and pulse rate (PR)
Range	1%~100%
Resolution	1%
Accuracy	> 20 kg and <10~20 kg: In the range of 70%~100%, the measurement error is $\pm 2$ ; <10 kg: In the range of 70%~100%, the measurement error is $\pm 3$ ; In the range of 0%~69%, the measurement error is not defined.
Average time	2 s-4 s, 4 s-6 s, 8 s, 10 s, 12 s, 14 s, 16 s
Data update peiriod	1 s
Weak perfusion condition	Pulse amplitude: >0.02%; Light transmittance: >5%.



## Appendix B Product Specifications

Weak perfusion SpO <sub>2</sub> accuracy	> 20 kg and <10~20 kg:±2% <10 kg:±3%.
--	--

### Nellcor Oximeter Module

Monitoring parameter	Pulse oximetry (SpO <sub>2</sub> ) and pulse rate (PR)
Range	1%~100%
Resolution	1%
Data update peiriod	1 s
Accuracy	> 20 kg: In the range of 70 %~100 %, the measurement error is ±2; <10 kg: In the range of 70 %~100 %, the measurement error is ±3; Insufficiency: In the range of 70 %~100 %, the measurement error is ±2; In the range of 0 %~69 %, the measurement error is not defined.

### B.6.4 PR Specifications

Alarm Limit Specifications	Range
PR Upper Limit	Alarm upper limit for > 20 kg: (Lower limit+2) bpm~250 bpm Alarm upper limit for <10~20 kg: (Lower limit+2) bpm~250 bpm Alarm upper limit for <10 kg: (Lower limit+2) bpm~250 bpm
PR Lower Limit	Alarm lower limit for > 20 kg: 25 bpm~ (Upper limit-2)bpm Alarm lower limit for <10~20 kg: 25 bpm~ (Upper limit-2)bpm Alarm lower limit for <10 kg: 25 bpm~ (Upper limit-2)bpm

### PR from SpO<sub>2</sub> Module

Range	30 bpm~250 bpm
Resolution	1 bpm
Measuring Tolerance	±2 bpm
Average Time	8 s

### PR from Masimo SpO<sub>2</sub> Module

Range	25 bpm~240 bpm
Resolution	1 bpm
Measuring Tolerance	The measuring tolerance is ±3 bpm or ±1%, whichever is higher.
Average Time	2 s-4 s, 4 s-6 s, 8 s, 10 s, 12 s, 14 s, 16 s

### PR from Nellcor SpO<sub>2</sub> Module

Range	20 bpm~300 bpm
Resolution	1 bpm
Measuring Tolerance	> 20 kg and <10 kg: 20 bpm~250 bpm: ±3 bpm Insufficiency: 251 bpm~300 bpm: not defined.

### PR from IBP

Range	30 bpm~350 bpm
Resolution	1 bpm
Measuring Tolerance	30 bpm~200 bpm: ±1 bpm or ±1%, whichever is higher; 201 bpm~350 bpm: ±2%.

**B.6.5 NIBP Monitoring**

Measuring Method	Automatic oscillometric method				
Safety Requirements	Acc. to ANSI/AAMI SP-10 Non-invasive Automated Blood Pressure Monitor, Part 4.4				
Work Mode	Manual, Auto, STAT Measuring				
Measuring Time under Continuous Mode	5 min				
Measuring Interval under Auto Mode	1 min, 2 min, 3 min, 4 min, 5 min, 10 min, 15 min, 30 min, 60 min, 90 min, 2 h, 4 h, 3 h, 8 h, Timer interval error: < 10 s				
Resolution	1 mmHg (0.133kPa)				
Nominal Range of Monitoring	<b>Blood Pressure (unit)</b>		<b>&gt;20 kg</b>	<b>10~20 kg</b>	<b>&lt;10 kg</b>
	Systolic Pressure	mmHg	40~270	40~200	40~135
		kPa	5.3~35.9	5.3~26.6	5.3~18.0
	Mean Pressure	mmHg	20~230	20~165	20~110
		kPa	2.7~30.6	2.7~22.0	2.7~14.7
	Diastolic Pressure	mmHg	10~210	10~150	10~100
		kPa	1.3~27.9	1.3~20.0	1.3~13.3
Range of Initial Inflation Pressure Setting	>20 kg: 80 mmHg~280 mmHg (10.7 kPa~37.3 kPa) 10~20 kg: 80 mmHg~210 mmHg (10.7 kPa~27.9 kPa) <10 kg: 60 mmHg~140 mmHg (8.0 kPa~18.6 kPa)				
Default of Initial Inflation Pressure	>20 kg: 160 mmHg (21.3 kPa) 10~20 kg: 140 mmHg (18.6 kPa) <10 kg: 90 mmHg (12.0 kPa)				
Measuring Tolerance of Pressure Source Testing	±3 mmHg (±0.4 kPa)				
Overpressure Protection	>20 kg state: When the pressure in cuff exceeds 297 mmHg (39.5 kPa)±3 mmHg (0.4 kPa), the control valve shall relieve the pressure. 10~20 kg state: When the pressure in cuff exceeds 240 mmHg (31.9 kPa)±3 mmHg (0.4 kPa), the control valve shall relieve the pressure. <10 kg state: When the pressure in cuff exceeds 147 mmHg (19.6 kPa)±3 mmHg (0.4 kPa), the control valve shall relieve the pressure.				
<b>Alarm Limit Specifications</b>	<b>Range</b>				
Upper Limit of Systolic Blood Pressure	>20 kg: (Lower limit+5) mmHg~270 mmHg ( (Lower limit+0.7) kPa~35.9 kPa) 10~20 kg: (Lower limit+5) mmHg~200 mmHg ( (Lower limit+0.7) kPa~26.6 kPa) <10 kg: (Lower limit+5) mmHg~135 mmHg ( (Lower limit+0.7) kPa~18.0 kPa)				
Lower Limit of Systolic Blood Pressure	>20 kg: 41 mmHg~ (Upper limit-5) mmHg (5.3 kPa~ (Upper limit -0.7) kPa) 10~20 kg: 40 mmHg~ (Upper limit-5) mmHg (5.3 kPa~ (Upper limit-0.7) kPa) <10 kg: 40 mmHg~ (Upper limit-5) mmHg (5.3 kPa~ (Upper limit-0.7) kPa)				
Upper Limit of Mean Blood Pressure	>20 kg: (Lower limit+5) mmHg~230 mmHg ( (Lower limit+0.7) kPa~30.6 kPa) 10~20 kg: (Lower limit+5) mmHg~165 mmHg ( (Lower limit+0.7) kPa~21.9.0 kPa) <10 kg: (Lower limit+5) mmHg~110 mmHg ( (Lower limit+0.7) kPa~14.6 kPa)				
Lower Limit of Mean Blood Pressure	>20 kg: 20 mmHg~ (Upper limit-5) mmHg (2.7 kPa~ (Upper limit-0.7) kPa) 10~20 kg: 20 mmHg~ (Upper limit-5) mmHg (2.7 kPa~ (Upper limit-0.7) kPa) <10 kg: 20 mmHg~ (Upper limit-5) mmHg (2.7 kPa~ (Upper limit-0.7) kPa)				
Upper Limit of Diastolic Blood	>20 kg: (Lower limit+5) mmHg~210 mmHg ( (Lower limit+0.7) kPa~27.9 kPa)				

## Appendix B Product Specifications

Pressure	10~20 kg: (Lower limit+5) mmHg~150 mmHg ( (Lower limit+0.7) kPa~20.0 kPa) <10 kgv: (Lower limit+5) mmHg~100 mmHg ( (Lower limit+0.7) kPa~13.3 kPa)
Lower Limit of Diastolic Blood Pressure	>20 kg: 11 mmHg~ (Upper limit-5) mmHg (1.4 kPa~ (Upper limit-0.7) kPa) 10~20 kg: 11 mmHg~ (Upper limit-5) mmHg (1.4 kPa~ (Upper limit-0.7) kPa) <10 kg: 10 mmHg~ (Upper limit-5) mmHg (1.3 kPa~ (Upper limit-0.7) kPa)

### B.6.6 Temperature (Temp) Monitoring

Range	0℃~50℃ (32°F~122°F)
Measuring Method	Thermal resistance method
Accuracy	The measuring tolerance is ±0.1℃ (exclusive of probe tolerance)
Updating Interval	1 s
Nominal Resistance of Temp. Sensor	2252 Ω (25℃)
Type of Temp. Sensor	YSI400 Sensor or its Compatible Sensor (Precision±0.1℃)
Channel Number	2 channels
Resolution	0.1℃
Alarm Indication	Audible & visual alarm, data and parameter blinking, alarm message displayed in the screen, 3 levels of alarm.
<b>Alarm Limit Specifications</b>	<b>Range (℃)</b>
Upper Limit	(Lower Limit +1)℃~50℃
Lower Limit	0℃~(Upper Limit -1)℃

### B.6.7 IBP Monitoring

Measuring Method		Invasive direct measuring
Volume displacement (Abbott)		<0.04 mm <sup>3</sup> /100mmHg
IBP		
Measuring Range		-50 mmHg~350 mmHg
Resolution		1 mmHg
Accuracy		±2% or ±1 mmHg, whichever is higher (exclusive of the sensor)
Updating Interval		1 s
<b>Alarm Limit Specifications</b>		<b>Range</b>
Art P1 P2	Upper Limit of Systolic Blood Pressure	(Lower limit+2) mmHg~350 mmHg ((Lower limit+0.3)kPa~46.7 kPa)
	Upper Limit of Mean Blood Pressure	
	Upper Limit of Diastolic Blood Pressure	
PA	Upper Limit of Systolic Blood Pressure	(Lower limit+2) mmHg~120 mmHg ( (Lower limit+0.3)kPa~16.0 kPa)
	Upper Limit of Mean Blood Pressure	
	Upper Limit of Diastolic Blood Pressure	
Art	Lower Limit of Systolic Blood Pressure	0 mmHg~(Upper limit-2)mmHg (0 kPa~(Upper limit-0.3)kPa)
	Lower Limit of Mean Blood Pressure	
	Lower Limit of Diastolic Blood Pressure	
P1 P2	Lower Limit of Systolic Blood Pressure	-50 mmHg~(Upper limit-2)mmHg (-6.7 kPa~(Upper limit -0.3)kPa)
	Lower Limit of Mean Blood Pressure	
	Lower Limit of Diastolic Blood Pressure	
PA	Lower Limit of Systolic Blood Pressure	-6 mmHg~(Upper limit-2)mmHg (-0.8 kPa~(Upper limit-0.3)kPa)
	Lower Limit of Mean Blood Pressure	
	Lower Limit of Diastolic Blood Pressure	
LAP RAP	Upper Limit of Mean Blood Pressure	(Lower limit+2)mmHg~40 mmHg ((Lower limit+0.3)kPa~5.3 kPa)

ICP CVP	Lower Limit of Mean Blood Pressure	-10 mmHg~(Upper limit-2)mmHg (-1.3 kPa~ (Upper limit-0.3)kPa)
------------	------------------------------------	--

### B.6.8 CO<sub>2</sub> Monitoring (Optional)

Measuring Mode	Sidestream type (support 50ml/min pumping rate), mainstream type
Measuring Method	Infrared radiation absorption technique

#### Phasein Sidestream ISA Module

Measuring Method	Infrared Spectrum Method
Measuring Mode	Sidestream
Range	0%~25%
Accuracy	At 0%~25%: $\pm (0.2\%+2\%$ of reading) At 15%~25%: undefined
Unit selection	%, mmHg, kPa
Preheating time	< 10 s (Report the concentration and reach the highest precision)
Total System Response Time	< 3 s (use of 2m sampling tube)
Primary agent threshold (ISA OR+/AX+)	0.15 vol%. When an agent is identified, concentrations will be reported even below 0.15 vol%
Secondary agent threshold (ISA OR+/AX+)	0.2 vol% + 10% of total agent concentration
Airway Leakage	$\leq 0.5$ ml/min
Range of Breathing Rate	0 rpm~150 rpm
Accuracy of Breathing Rate	$\pm 1$ rpm
Asphyxia Alarm Delay	20 s, 25 s, 30 s, 35 s, 40 s, 45 s, 50 s, 55 s, 60 s
Sampling Flow Rate	50 ml/min $\pm 10$ ml/min
Automatic Pressure Compensation	yes
<b>Alarm Limit Specifications</b>	<b>Range</b>
EtCO <sub>2</sub> Upper Limit	(Lower Limit +2)mmHg~99 mmHg
EtCO <sub>2</sub> Lower Limit	0 mmHg~(Upper Limit -2)mmHg
FiCO <sub>2</sub> Upper Limit	0 mmHg~99 mmHg
awRR Upper Limit	(Lower limit+2) rpm~100 rpm
awRR Lower Limit	0 rpm~ (Upper limit-2) rpm

#### Phasein Mainstream IRMA Module

Measuring Method	Infrared Spectrum Method
Measuring Mode	Mainstream
Range	0%~25%
Accuracy	Range:0 %~15 %, Default: $\pm(0.2\%+$ reading 2 %); Range:15%~25%, Default: Undefined.
Resolution	1 mmHg (0.133 kPa)
Unit selection	%, mmHg, kPa
Total System Response Time	< 1 s
Primary agent threshold	0.15 vol%. When an agent is identified, concentrations will be reported even below 0.15 vol% as long as apnea is not detected.
Secondary agent threshold	0.2 vol% + 10% of total agent concentration
Range of Breathing Rate	0 rpm~150 rpm
Accuracy of Breathing Rate	$\pm 1$ rpm

## Appendix B Product Specifications

Asphyxia Alarm Delay	20 s, 25 s, 30 s, 35 s, 40 s, 45 s, 50 s, 55 s, 60 s
Automatic Pressure Compensation	yes
<b>Alarm Limit Specifications</b>	<b>Range</b>
EtCO <sub>2</sub> Upper Limit	(Lower Limit +2)mmHg~99 mmHg
EtCO <sub>2</sub> Lower Limit	0 mmHg~(Upper Limit -2)mmHg
FiCO <sub>2</sub> Upper Limit	0 mmHg~99 mmHg
awRR Upper Limit	(Lower limit+2) rpm~100 rpm
awRR Lower Limit	0 rpm~ (Upper limit-2) rpm

### Respironics Sidestream LoFlo Module

Measuring Method	Infrared Spectrum Method
Measuring Mode	Sidestream
Preheating time	Max. length of waveform is 20s. Full accuracy requirements satisfied after 2min (environment temp.: 25℃)
Range	0%~19.7% (0 mmHg ~ 150 mmHg) (0 kPa~20 kPa)
Resolution	0.1 mmHg 0 mmHg~69 mmHg 0.25 mmHg 70 mmHg~150 mmHg
Stability	Short-term drift: ≤0.8 mmHg (0.1 kPa) within 4h Long-term drift: accuracy maintained within 120h.
Unit selection	%, mmHg, kPa
Accuracy (Gas Temp. at 25℃)	0 mmHg~40 mmHg (0 kPa~5.3 kPa), ±2 mmHg (0.27 kPa) 41 mmHg~70 mmHg (5.5 kPa~9.3 kPa), ±5% of the reading 71 mmHg~100 mmHg (9.4 kPa~13.3 kPa), ±8% of the reading 101 mmHg~150 mmHg (13.4 kPa~20 kPa), ±10% of the reading (When the breathing rate is > 80 rpm, all ranges are ±12% of the reading)
Total System Response Time	< 3 s
Range of Breathing Rate	2 rpm~150 rpm
Accuracy of Breathing Rate	±1 rpm
Asphyxia Alarm Delay	20 s, 25 s, 30 s, 35 s, 40 s, 45 s, 50 s, 55 s, 60 s
Sampling Flow Rate	≥50 ml/min(100Hz)
Automatic Pressure Compensation	no
<b>Alarm Limit Specifications</b>	<b>Range</b>
EtCO <sub>2</sub> Upper Limit	(Lower Limit +2) mmHg~99 mmHg
EtCO <sub>2</sub> Lower Limit	0 mmHg~(Upper Limit -2) mmHg
FiCO <sub>2</sub> Upper Limit	0 mmHg~99 mmHg
awRR Upper Limit	(Lower limit+2) rpm~100 rpm
awRR Lower Limit	0 rpm~ (Upper limit-2) rpm

**Respironics Mainstream CAPNOSTAT5 Module**

Measuring Method	Infrared Spectrum Method
Measuring Mode	Mainstream
Preheating time	Max. length of waveform is 15s. Full accuracy requirements satisfied after 2min (environment temp.: 25℃)
Range	0%~19.7% (0 mmHg~150 mmHg) (0 kPa~20 kPa)
Resolution	0.1 mmHg 0 mmHg~69 mmHg 0.25 mmHg 70 mmHg~150 mmHg
Stability	Short-term drift: ≤0.8 mmHg (0.1 kPa) within 4h Long-term drift: accuracy maintained within 120h.
Rise Time	< 60 ms
Unit selection	%, mmHg, kPa
Accuracy (Environment Temp. at 35℃)	0 mmHg~40 mmHg (0 kPa~5.3 kPa), ±2 mmHg (0.27 kPa) 41 mmHg~70 mmHg (5.5 kPa~9.3 kPa), ±5% of the reading 71 mmHg~100 mmHg (9.4 kPa~13.3 kPa), ±8% of the reading 101 mmHg~150 mmHg (13.4 kPa~20 kPa), ±10% of the reading
Range of Breathing Rate	0 rpm~150 rpm
Accuracy of Breathing Rate	±1 rpm
Asphyxia Alarm Delay	20 s, 25 s, 30 s, 35 s, 40 s, 45 s, 50 s, 55 s, 60 s
Sampling Flow Rate	100 Hz
Automatic Pressure Compensation	no
<b>Alarm Limit Specifications</b>	<b>Range</b>
EtCO <sub>2</sub> Upper Limit	(Lower Limit +2)mmHg~99 mmHg
EtCO <sub>2</sub> Lower Limit	0 mmHg~(Upper Limit -2)mmHg
FiCO <sub>2</sub> Upper Limit	0 mmHg~99 mmHg
awRR Upper Limit	(Lower limit+2) rpm~100 rpm
awRR Lower Limit	0 rpm~ (Upper limit-2) rpm

**Kingst KM7002-V33/KM7003-V40 Sidestream Module**

Measuring Method	Non-scattering Infrared Gas Analysis
Measuring Technology	Non-dispersive Infrared Gas Analysis (NIDR)
Range	0%~20% (0 mmHg~150 mmHg) (0 kPa~20 kPa)
Protection Level / Type	BF
Preheating time	2 min at 25 °C
Response Time	50 ml/min
Delay Time	50 ml/min
Fully-automatic Drift Calibration	Automated according to the time and temperature. Time 5 s~8 s
Airway Leakage	< 0.1% (within the flow range above)
Accuracy	When < 5.0%: ±0.3% (±2.0 mmHg) (0.27 kPa) When ≥5.0%: < 6% of the reading
Range of Breathing Rate	3 rpm~150 rpm
Accuracy of Breathing Rate	1% or ±1 rpm, whichever is higher.
Asphyxia Alarm Delay	30 s, 35 s, 40 s, 45 s, 50 s, 55 s, 60 s

## Appendix B Product Specifications

Automatic Pressure Compensation	yes
<b>Alarm Limit Specifications</b>	<b>Range</b>
EtCO <sub>2</sub> Upper Limit	(Lower Limit +2)mmHg~99 mmHg
EtCO <sub>2</sub> Lower Limit	0 mmHg~(Upper Limit -2)mmHg
FiCO <sub>2</sub> Upper Limit	0 mmHg~99 mmHg
awRR Upper Limit	(Lower limit+2) rpm~100 rpm
awRR Lower Limit	0 rpm~ (Upper limit-2) rpm

### B.6.9 C.O. Specifications(Optional)

Measurement method	Thermodilution method	
Measuring range	C.O.:	0.01~20L/min
	TB:	23~43℃
	TI:	0~27℃
Resolution	C.O.:	0.01L/min
	TB, TI:	0.1℃
Accuracy	C.O.:	±5% or ±0.1 L /min, whichever is greater
	TB, TI:	±0.1℃ (without sensor)

Alarm Limit Specifications	Range
TB Upper Limit	(Lower Limit + 1.1)~43℃ (Lower Limit + 2)~109.4°F
TB Lower Limit	23~(Upper Limit - 1.1)℃ 73.4~(Upper Limit - 2)°F

### B.6.10 AG Specifications (Optional)

Measurement method	Infrared radiation absorption characteristics	
Warm-up time	30 s	
Measuring range	CO <sub>2</sub> :	0%~25%
	O <sub>2</sub> :	0%~100%
	N <sub>2</sub> O:	0%~100%
	Des:	0%~25%
	Sev:	0%~25%
	Enf:	0%~25%
	Iso:	0%~25%
	Hal:	0%~25%
	awRR:	0 rpm~254 rpm
Resolution	CO <sub>2</sub> : 1 mmHg awRR: 1 rpm	
Measurement accuracy drift	Meet the accuracy requirements within 6 hours	
Suffocation alarm delay	20 s, 25 s, 30 s, 35 s, 40 s, 45 s, 50 s, 55 s, 60 s	
Update time	1 s	
IRMA AX+	Primary agent threshold	0.15 vol%. When an agent is identified, concentrations will be reported even below 0.15 vol% as long as apnea is not detected.

	Secondary agent threshold	0.2 vol% + 10% of total agent concentration			
ISA OR+/AX+	Primary agent threshold	0.15 vol%. When an agent is identified, concentrations will be reported even below 0.15 vol% as long as apnea is not detected.			
	Secondary agent threshold	0.2 vol% + 10% of total agent concentration			
Interfering gases and steam effect					
gases and steam	Gas concentration	Carbon dioxide			
		IRMA CO <sub>2</sub> 、OR	IRMA AX+/OR+	Anesthetic gas	Nitrous oxide
N <sub>2</sub> O <sup>4)</sup>	60 vol%	_1&2)	_1&2)	_1)	_1)
Hal <sup>4)</sup>	4	_1)	_1)	_1)	_1)
Enf, Iso, Sev <sup>4)</sup>	5	Reading of +8% <sup>3)</sup>	_1)	_1)	_1)
Des <sup>4)</sup>	15	Reading of +12% <sup>3)</sup>	_1)	_1)	_1)
Xe (Xenon) <sup>4)</sup>	80	Reading of -10% <sup>3)</sup>		_1)	_1)
He (Helium) <sup>4)</sup>	50	Reading of -6% <sup>3)</sup>		_1)	_1)
Quantitative spray <sup>4)</sup>	Not for quantitative spray				
Ethanol <sup>4)</sup>	0.3	_1)	_1)	_1)	_1)
Isopropano <sup>4)</sup>	0.5	_1)	_1)	_1)	_1)
Acetone <sup>4)</sup>	1	_1)	_1)	_1)	_1)
Methane <sup>4)</sup>	3	_1)	_1)	_1)	_1)
Carbon monoxide <sup>4)</sup>	1	_1)	_1)	_1)	_1)
Nitric oxide <sup>5)</sup>	0.02	_1)	_1)	_1)	_1)
Oxygen <sup>5)</sup>	100	_1&2)	_1&2)	_1)	_1)
1): “Accuracy _ All conditions” The specification contains negligible interference and influence. 2): for the probe which cannot be measured, nitrous oxide and / or the concentration of oxygen should be set. (IRMA CO <sub>2</sub> not measure Nitrous oxide or oxygen, IRMA AX+ not measure the oxygen) 3): the gas concentration interference indicated, such as 50vol% of the helium usually leads to a decrease of 6% carbon dioxide readings. That is, if the measurements contain 5.0% vol% of carbon dioxide and 50vol% of nitrogen mixed gas, the actual measured concentration of carbon dioxide is usually as follows: (1-0.06)*5.0vol%=4.7vol% Carbon dioxide. 4): meet the to EN ISO 21647:2004 standard. 5): supplement EN ISO 21647:2004 standard.					

Alarm Limit Specifications	Range
EtCO <sub>2</sub> Upper Limit	(Lower Limit +2)mmHg~99 mmHg
EtCO <sub>2</sub> lower limit	0 mmHg~(Upper Limit -2)mmHg
FiCO <sub>2</sub> Upper Limit	0 mmHg~99 mmHg
awRR Upper Limit	(lower limit+2) rpm~100 rpm
awRR lower limit	0 rpm~(upper limit-2) rpm
FiEnf Upper Limit	(lower limit+0.2)%~8%
FiEnf lower limit	0%~(upper limit-0.2)%
EtEnf Upper Limit	(lower limit+0.2)%~8%
EtEnf lower limit	0%~(upper limit-0.2)%
EtHal Upper Limit	(lower limit+0.2)%~8%
EtHal lower limit	0%~(upper limit-0.2)%
Filso Upper Limit	(lower limit+0.2)%~8%
Filso lower limit	0%~(upper limit-0.2)%
Etlso Upper Limit	(lower limit+0.2)%~8%
Etlso lower limit	0%~(upper limit-0.2)%
EtSev Upper Limit	(lower limit+0.2)%~10%
EtSev lower limit	0%~(upper limit-0.2)%
FiSev Upper Limit	(lower limit+0.2)%~10%
FiSev lower limit	0%~(upper limit-0.2)%



## Appendix B Product Specifications

EtDes Upper Limit	(lower limit+0.2)%~22%
EtDes lower limit	0%~(upper limit-0.2)%
FiDes Upper Limit	(lower limit+0.2)%~22%
FiDes lower limit	0%~(upper limit-0.2)%
FiO <sub>2</sub> Upper Limit	(lower limit+16) mmHg~760 mmHg ((lower limit+2.1) kPa~101.1 kPa)
FiO <sub>2</sub> lower limit	137 mmHg~(upper limit-16) mmHg (18.3 kPa~(upper limit-2.1) kPa)
EtO <sub>2</sub> Upper Limit	(lower limit+16) mmHg~760 mmHg ((lower limit+2.1) kPa~101.1 kPa)
EtO <sub>2</sub> lower limit	137 mmHg~(upper limit-16) mmHg (18.3 kPa~(upper limit-2.1) kPa)
FiN <sub>2</sub> O Upper Limit	(lower limit+2)%~82%
FiN <sub>2</sub> O lower limit	0%~(upper limit-2)%
EtN <sub>2</sub> O Upper Limit	(lower limit+2)%~100%
EtN <sub>2</sub> O lower limit	0%~(upper limit-2)%

### B.6.11 Recorder Specifications

Recorder	To record the animal information, the hospital information, waveform, parameters and others displayed in the screen
Method	Thermal array recorder
Printing Paper	Thermal paper
Print Resolution	8 dots/mm on Y-Axis
Delay Characteristics	≤0.5 mm
Amplitude-frequency Characteristics	Monitor Mode: 0.5 Hz~40 Hz; Diagnose Mode: 0.05 Hz~150 Hz.
Time Constant	≥0.3 s



Care from  
hearts

**Shenzhen Biocare Bio-Medical Equipment Co., Ltd.**

Add:#16-1, Jinhui Road, Jinsha Community, Kengzi Sub-District, Pingshan New District, 518122 Shenzhen, P.R.CHINA

Tel: +86-755-2796 0888/3300 5899

Fax: +86-755-2796 0643

E-mail: [sales@biocare.com.cn](mailto:sales@biocare.com.cn)

Website: [www.biocare.com.cn](http://www.biocare.com.cn)

© 2025 BIO CARE all rights reserved. Specifications subject to changes without prior notice.

ENG-CATA-V2.1-20250122