**LAMIDEY  
NOURY**  
M E D I C A L

# MODULE ARGON

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
**INDEX**

1.	Introduction .....	4
1.1	Applicability .....	4
1.2	Classification of the medical device .....	4
1.3	Meaning of the symbols used.....	5
1.4	Composition of the Argon electrosurgical system .....	6
1.4.1	List of system elements .....	6
1.4.2	The electrosurgical unit .....	6
1.4.3	Argon cylinder recommended .....	6
1.4.4	The pressure regulator .....	7
1.4.5	The accessories for use.....	7
1.5	Définitions.....	9
1.6	General description .....	9
1.7	Argon plasma coagulation .....	9
2.	Receiving the device .....	10
2.1	Reception checking.....	10
2.2	Operating environment and installation .....	10
2.2.1.	Operating environment .....	10
2.2.2.	Installation.....	10
2.2.3.	Installation and connections of pressure regulator .....	11
2.2.4.	Connections on the back of the Argon module.....	13
2.2.5.	Connections on the front of the Argon module .....	14
2.3.	Controls before use .....	15
2.3.1.	Control of accessories.....	15
2.3.2.	Control of Argon system.....	15
3.	Précautions for use .....	16
3.1.	Precautions related to the use of Argon Gas.....	16
3.2.	Precautions related to the use of electrosurgical unit .....	16
3.3.	Recommandations related to the use of the Argon module.....	16
3.4.	Risks of serious adverse events .....	17
4.	Use of the Argon module .....	18
4.1.	Putting into service the argon module.....	18
4.2.	Initialization of the Argon module.....	18
4.3.	Error Messages .....	20
4.4.	Settings .....	20
4.4.1.	Open surgery .....	20
4.4.2.	Digestive endoscopy .....	20
4.5.	After use .....	21
5.	Maintenance, cleaning, sterilization.....	21
5.1.	Argon Module.....	21
5.2.	Argon accessories .....	21
5.3.	Maintenance.....	21
5.4.	Environment protection (Directive 2002-96-CE).....	21
6.	Electromagnetic compatibility .....	22
6.1.	General instructions .....	22
6.2.	Cable length.....	22
6.3.	Electromagnetic emission.....	23
6.4.	Electromagnetic and magnetic immunity .....	24
7.	Technical description .....	26
7.1.	Identification of the medical device .....	26
7.2.	Connection to the mains power network.....	26
7.3.	Other features .....	26
7.4.	Maintenance.....	26
7.4.1	Preventive Maintenance .....	27

---

7.4.2 Replacement of spare parts .....28

## 1. INTRODUCTION

### **IMPORTANT**

The instructions below are an integral part of the device that you have just bought. They describe the operation and the use of the device and must be read carefully before using it.

All the safety instructions and warning notes must be observed with scrupulously. Make sure that they are conveyed whenever the device is used by new personnel. No portion of this document may be photocopied, reproduced, or translated without the written approval of **LAMIDEY-NOURY MEDICAL**. If necessary, contact our Technical Assistance.

THE SYMBOL OPPOSITE IS RESERVED TO DRAW ATTENTION TO WARNINGS AND INSTRUCTIONS THAT MUST BE FOLLOWED ABSOLUTELY

#### 1.1 Applicability

The use of Argon Modules is reserved for surgeons. They are designed to be used intermittently for surgical operations in the operating room. It is intended for the following specialties: Visceral and digestive surgery, Digestive endoscopy, Thoracic surgery, Urology, gynaecology; ...

The Argon Modules are complementary to the HF electrosurgical units, intended to apply Argon gas jet current. When a fulguration current is used, this voltage allows to ionize the gas in order to make a contactless surface coagulation, by argon plasma. The use of argon plasma allows to target the tissues precisely, because the electric arc is channeled by the gas flow. The use of the Argon also allows to improve the visibility of the treated zone by blowing blood and residuals to the surface.

#### 1.2 Classification of the medical device

Following the criterias of the classification of the medical devices mentioned by the norm CEI 60601-1, the table below points out classifications that can be used with Argon module:

Criteria of classification	Classification
Protection against electric shock	Classe I
Applied parts	Device without applied parts
Protection against harmful effects of the penetration of water or solid bodies.	Device non protected (See § 3.3 and 5.1)
Method(s) of sterilization	Device not sterilizable
Suitability for use in an oxygen rich environment	Not applicable : use forbidden in an oxygen rich environment

The technical features are detailed in the section 7

**1.3 Meaning of the symbols used**

**SYMBOLS IN THE FRONT PLATE OF THE ARGON**



READ CAREFULLY THE USE INSTRUCTIONS: non-respect of these instructions can lead the patient or the surgeon to run a risk.

	Identification of the switching-on warning light (in the front plate)		Control keys of the Argon gas output (+ for incrementation ; - for decrementation)
	Start key of the Argon gas purge cycle		Selection key in the configuration menu

**OTHER SYMBOLS ON THE DEVICE**

	"On" position of the switch sector located on the back (seeing alight)		"Off" Position of the mains switch on the rear panel (light off)
	Equipotentiality terminal		Protective ground terminal (on the AC power cord)
	Reference code		Serial code
	Name and address of the producer		Date of manufacture (XXXX)
	Weight of the empty device		Maximum safety load
	Single Phase		Compliance with Directive 93/42/CEE (0459 : organism notified LNE/G-Med issuer of the certificate)
	Separate collection for electrical and electronic equipment		

**SYMBOLS ON PACKAGING CARTONS**

	FRAGILE : handle with care		Keep dry
	Transport and storage position : This way up		Temperature limits for storage (-20°C – 70°C)
	Humidity limits for storage (500-1060 hPa)		Hygrometry limits for storage (0-95%)

## 1.4 Composition of the Argon electrosurgical system

### 1.4.1 List of system elements

The Argon Coagulation Module is used with the replaceable elements listed below:

1 OPTIMA electrosurgical unit, Ref V10GOPT4, and its monopolar-use accessories
1 cylinder of medical-surgical Argon (to be supplied by a medical gas supplier)
1 pressure regulator with low pressure outlet (Ref A301 9813), equipped with a pressure sensor
1 main power cable, length 5 m (Ref V11F3AX)
1 control cable (Ref A409 0312), length 42 cm
1 gas connection between pressure regulator and module (Ref A 409 0345C)
1 sterile disposable antibacterial filter, with male / female Luer Lock fittings (Ref V11A025)

Paragraph 7.4.2 details the precautions for replacing spare parts.

### 1.4.2 The electrosurgical unit

The Argon module has been designed to work in association with the LAMIDEY NOURY MEDICAL OPTIMA unit equipped with the Argon Option. This unit is used placed on the top of the Argon module.

The use of other electrosurgical units near the Argon module must be excluded as it may cause unexpected Argon module malfunctions (see § 6.2)

### 1.4.3 Argon cylinder recommended


The Argon electrosurgical System must be used with pure Argon gas, with medical grade surgical. The system has been validated with gas supplied by AIR LIQUID (Purity > 99.9996%):

Reference Air Liquide : M2001S05R2A001	Charge : 1 m <sup>3</sup>	Diameter : 14cm
Tradename : CH1	Operating pressure : 200 bar	Height with hat :62 cm
Water capacity : 5 liters	Connection : ISO 5145, N°4	Full mass: 11.5 Kg

- Insufficient purity gas can damage the pressure regulator prematurely and cause pressure regulation malfunction.
- Do not install cylinders with a pressure higher than 200 bar.
- Do not exceed the expiry date of the gas.
- Handle the bottle carefully (see the gas supplier's operating instructions)
- Store away from heat

### 1.4.4 The pressure regulator

The Argon module is delivered with a regulator equipped with a pressure sensor. It is intended to limit the pressure of use to less than 2 Bars (2000 Hectopascal). The sensor is connected to the Argon module to provide an indication of the fill level in the cylindre.

#### Features :

Manufacturer : CAHOUIET	Model NOVA 40	Connection type ISO 5145 N ° 4
Maximum inlet pressure 200 bar	Release pressure 1,7 Bars	Operating temperatures and storage -20°C à 60°C

Medical device marked CE 0459.

#### Safety instructions of the pressure reducer:

- Do not attempt to disassemble the pressure regulator or modify its settings
- The assembly on the cylinder must be done strictly without grease or oil
- Tighten by hand, without tools, to avoid damaging the O-ring
- After assembly, always stand on the side and not facing the pressure regulator and check the operation by opening the gas valve very gradually.
- Never lift the cylinder by the pressure regulator
- Preventive control every 5 years
- Disposal after 15 years of use
- Only replace the pressure regulator with a model designed for Argon gas

### 1.4.5 The accessories for use

Antibacterial filter (item code V11A025): The Argon module and its elements are not sterile. To avoid contaminating the patient by the gas jet, it is imperative to install a disposable antibacterial filter on the gas outlet on the front of the module.

Features : - filtration 0,22 µ - Sterile single use – Luer Lock male/female connection

*The Argon handle (item code V11A201)*

The Argon handle, with manual controls, with 3.5m cable is reusable and sterilizable with steam at 134 ° C / 20 min. It is equipped with a double connection: a tri-plug type for the electro-surgical unit, and a luer lock connection for argon module.

Argon electrodes: 3 models of electrodes, mounted on the argon handle, allow the application of argon plasma coagulation:

In laparoscopic surgery, there is a risk of embolism by intra-peritoneal overpressure. Do not use argon gas without overpressure safety system.

Argon probes for flexible endoscopy:

Digestive endoscopy, length 2.2m:

Reusable probe diam 2,3mm – code art. V11A401

Disposable probe diam 2,3mm – code art. VSUA222

Reusable probe diameter 3.2 mm – code art V11A402

Disposable probe diam. 3.2 mm – code art. VSUA223


In bronchial fibroscopy:

-Reusable probe, long 1,5m, diam. 1.5 mm - art. V11A403

-Disposable probe, long 1,5m, diam. 1.5 mm - art. VSUA150

-Disposable probe, long 1m, diam. 2.3 mm - art. VSUA102

Connection cable for Argon probes: 3.5 m long cable with tri-plug for the electro-surgical unit, and Luer lock fitting for the Argon module. Article number V11A400

## 1.5 Définitions

*HF Currents:* Electrosurgical units are high frequency alternating current electrosurgical generators (HF currents or radiofrequency currents). According to the norm EN 60601-2-2, the fundamental emission frequency of electrosurgical units must be greater than or equal to 200 kHz. The waveform is pure sinusoidal for a pure section current. It is sinusoidal modulated by bursts in the other cases.

*Active electrode:* An active electrode is an electrode in contact with which the HF energy is converted into thermal energy (Joule effect) to obtain cutting or coagulation effects. In Argon electrosurgery, the active electrodes are in the form of a needle, inserted into an insulating tube and mounted at the end of a manually operated electrode holder, or in the form of a flexible probe gastrointestinal or bronchial fibroscopy probe.

*Neutral electrode:* A neutral electrode is a large surface electrode placed in contact with the patient to close the HF circuit without producing any significant thermal effect. This electrode is commonly called plate, dispersive electrode or plate.

*Monopolar mode:* This is the mode of application of the HF current for which the patient is traversed by the HF current between an active electrode and a neutral electrode. Argon electrosurgery works in monopolar mode.

*Argon Plasma (or APC: Argon Plasma Coagulation):* argon plasma is obtained by ionization of the gas subjected to a very high electric voltage, such as that produced by the coagulation by fulguration. The flow of ionized gas then conducts the HF current on the targeted tissues, without contact with these tissues.

## 1.6 General description

The Argon LAMIDEY NOURY Medical module is a complementary module to electrosurgical units equipped with fulguration coagulation, designed to control the flow of Argon gas delivered by an active Argon monopolar electrosurgical electrode.

The gas flow is emitted during the entire duration of activation of the HF current, at the end of the electrode mounted on the Argon handle, or at the end of an Argon probe.

The use of the ARGON module requires the application of a neutral electrode in contact with the patient. Only the neutral electrodes of the double zone type, divided into two conductive parts, allow monitoring of the patient-plate contact, to avoid the risk of burns caused by a current density too high.

The mains power supply and the mains switch are located on the rear panel of the module, the controls for key settings and the LCD control panel are on the front panel.

Automatic control systems of the internal parameters, report the anomalies detected.

The filling level of the bottle is indicated on the screen of the module. When the gas load falls below <10%, at each activation an audible signal and a message warn the user to replace the cylinder. The module allows an adjustment of the gas flow from 0 to 10 Liters / minutes.

## 1.7 Argon plasma coagulation

Argon plasma coagulation is a form of monopolar coagulation without contact: the current is transmitted by a jet of ionized gas.

By selecting an argon plasma coagulation mode on the OPTIMA unit, it delivers a high-voltage coagulation of fulguration type, able to ionize the gas. The higher the power setting, the greater the distance needed to prime the arc. It is indeed necessary to approach the electrode of application sufficiently close to the tissues so that the ionization begins.

The intensity of the electric arc depends on the flow of gas. By finely adjusting this flow, the coagulation effect is adjusted precisely. The benefits of Argon plasma coagulation on contact coagulation are as follows:

- It allows to coagulate large areas with diffuse bleeding in a reduced time.
- Argon reduces the effects of carbonization.
- Argon plasma coagulation acts on a shallow depth.

## 2. RECEIVING THE DEVICE

### 2.1 Reception checking

During reception, inspect the unit, in order to note any damage during transport. Claims will only be accepted if they have been reported to the deliverer or directly to the sender. If the unit is returned, it is imperative to use the original packaging or any other that can ensure its integrity for his return. Remove the device from its packaging and read carefully the documentation and instructions provided. The local mains voltage must match the value indicated below the power plug (230 Vac).

### 2.2 Operating environment and installation

#### 2.2.1. Operating environment

The ARGON module is suitable for use in hospitals, in the operating room, with an electrosurgical unit, model OPTIMA, stacked on the module. This device is intended for use in medical premises of Group 2 from the standard IEC 60364-7-710, isolated from the public mains network according to the medical IT scheme.

The module must be connected to the mains power supply dedicated to the operating room. The use of the electrosurgical system with argon module, connected to the public mains network, is likely to cause electromagnetic disturbances on equipment not designed for use in the operating room.

**FIRE HAZARD:** The use of the system near a source of flammable gas or oxygen is prohibited due to the arcing at the tip of the active electrode. A safety distance greater than 25 cm from any source of oxygen must be respected.

Electromagnetic environment: see paragraph 6.2 detailing the electromagnetic environment and its exclusions.

#### 2.2.2. Installation.


**WARNING :** Castors fitted to the device are intended to facilitate its movement within the operating room, for convenience.

**THEY ARE NOT INTENDED FOR DISPLACEMENTS INCLUDING THRESHOLD OR ROLLING ON CABLE.** It is therefore necessary to take all the necessary precautions to avoid overbalancing of the device in case of obstacle on the path of the Argon module.

The Argon module is intended to carry only the electrosurgical units made by LAMIDEY NOURY MEDICAL, the defined Argon gas cylinder and the accessories necessary for the use of the Argon system (double pedal, connecting cables, application electrodes). The addition of any additional devices on the Argon module is not advised.

The unit must be installed on the trolley in which the ARGON MODULE is integrated. The Argon cylinder with its gaz regulator is installed on the lower platform of the trolley.

**2.2.3. Installation and connections of pressure regulator**

Before installing the pressure regulator on the cylinder for the first time, check that the cylinder is equipped with an intact tamper-proof guarantee system when using it for the first time. Remove this system to mount the regulator.


**DANGER, HIGH PRESSURE:** IT IS STRICTLY FORBIDDEN TO REMOVE THE HEAD OF TAP  
DO NOT POSITION AGAINST THE EXIT OF THE GAS WHEN OPENING THE TAP.


- a- Slightly open the valve before mounting the regulator to remove any dust.
- b- Remove the yellow cover of the regulator and present it against the threaded connection of the bottle.
- c- Tighten MANUALLY the high pressure connection clockwise. The pressure sensor should be held up. DO NOT TIGHTEN TO THE CLAMP (risk of damage to the seal)
- d- Connect the low pressure hose to the regulator
- e- Present the pressure sensor connector facing the sensor, aligning the coding pin with the corresponding groove in the socket. TIGHTEN CLOCKWISE, WITHOUT FORCING

**PRESSURE SENSOR CONNECTION CODING**

### 2.2.4. Connections on the back of the Argon module

a- Connect the low pressure gas hose to the quick coupling. Just push the tube into the fitting as far as it will go. (To remove the tube, simply push on the fitting ring, which releases the tube)

b- Connect the Argon control cable between the unit and the module. This cable is used to transmit the unit activation command to the Argon module. ATTENTION: the two connectors have a polarizer

c- Connect the power cable to the corresponding socket on the rear panel of the unit: the unit, normally supplied with 230 V 50 Hz, is supplied with a three-conductor detachable power cord terminated normalized 10/16 A & Earth. The envelope being metallic, it must be connected to a protective earth. The continuity of the earth connection, which guarantees the safety of the patient and the operators, must be regularly checked.

**IT IS FORBIDDEN TO OPERATE THE UNIT IF IT IS NOT EARTHED**

d- The equipotential connection terminal can be connected to another equipotential terminal, such as that of the HF electro-surgical generator, in order to reduce the effects of possible electromagnetic disturbances.

e- If necessary, connect the control pedal of the electro-surgical unit

### 2.2.5. Connections on the front of the Argon module

HF socket of the Argon handle

Luer lock connector type for Argon gas

M2

Plate connection

M1

HF connection

Screen of control

Antibacterial filter

Argon connection

- a- **PLATE CONNECTION:** apply a neutral plate to the patient and connect it to the electrosurgical unit (see section 3.4)
- b- Connect the sterile antibacterial filter to the Argon gas outlet of the module (Luer Lock connection)
- c- Then connect the ARGON handle Luer Lock fitting, or ARGON probe adapter cable, to the antibacterial filter.
- d- Lastly connect the 3-pin HF connection socket: if it is an Argon handle, it can be connected either in the M1 or M2 jacks of the OPTIMA unit. In the case of a gastrointestinal probe connection, only the M1 socket should be used (with red socket).

### 2.3. Controls before use

Daily control of the module before the first intervention of the day and a control of the accessories before each intervention make it possible to optimize the safety by preventing any failures.

#### 2.3.1. Control of accessories

THE USE OF DEFECTIVE ELECTROSURGICAL ACCESSORIES PRESENTS RISK OF THERMAL DAMAGE AND FIRE

Before each intervention, make sure that all the accessories that are to be used are available and, except pedals and plate cable, are in the sterile state.

By their use, the electrosurgical accessories undergo an inevitable wear due to mechanical aggression (pulls and torsions of cables, cuts, falls) and chemical. This requires careful visual examination of the items and a functional test before each use. Items with visible signs of deterioration should not be used, especially in case of defective insulation.

It is possible to test the operation of the hand controls or pedals by connecting them to an electrocautery, after setting the HF powers to 0. Activating the controls (yellow button in section or yellow lever for the pedal, or the blue button in coagulation or blue lever for the pedal), the unit indicates by an acoustic and visual signal (blue message for the coagulation or yellow

#### 2.3.2. Control of Argon system

Before the first intervention of the day, a visual check of the equipment should be carried out. Make sure it is well grounded. If the device shows signs of shock, or if a plug is damaged, it should not be used without a thorough technical inspection by a **LAMIDEY-NOURY MEDICAL** approved service.

When the Argon module is powered up, an automatic control of the module and the pressure of the connected bottles is performed. If a malfunction is detected, it is indicated on the display and by alarm tones, before the automatic shutdown of the device. Paragraph 4 describes how to start the module and how to remedy the faults reported.

**3. PRECAUTIONS FOR USE****3.1. Precautions related to the use of Argon Gas**

- a. VERIFY THE CONTENT OF BOTTLES by checking their labeling. The color of the bottle warheads is not an indication of their contents.
- b. NEVER USE THE ARGON WITHOUT PRESSURE REGULATOR to reduce the gas pressure.
- c. DO NOT LIFT THE BOTTLES OF GAS BY THEIR TAP OR BY THEIR PRESSURE REGULATOR.
- d. Surgical Medical Argon is only intended for use in argon cutting and coagulation. It should not be inhaled or used for other applications.
- e. The accessories of use (sleeves, electrodes, probes) must not contain residues of liquid after cleaning, decontamination and rinsing. The presence of liquid in the circuit compromises the efficiency of argon plasma coagulation. Prior to sterilization of the accessories, the liquid residues must be removed by means of a filtered air jet applied into the gas connection nozzle.

**3.2. Precautions related to the use of electrosurgical unit**

- a. Non observance of the IFU of electrosurgical unit can lead to damages or serious damages. Consult the IFU of the associate unit of the Argon module, relative to the use of the monopolar mode and neutral electrodes.
- b. ATTENTIVE EXAMINATION OF ACCESSORIES, particularly cables, electrodes and probes, is recommended before use, especially to ensure that they are properly insulated (see section 2.3.1)

**3.3. Recommendations related to the use of the Argon module**

- a. Although the enclosure of the apparatus is designed and intended to provide protection against the ingress of water or particles of material, the user must take all necessary measures to avoid any contact with water, which, as for any electrical equipment, could cause injury (electric shock).
- b. Improper storage and / or physical stress during transport can lead to failure or malfunction that may cause damage. Therefore, it is necessary that the user checks the correct operation of the device before using it and strictly follows the mandatory verification operations described in the instructions for use.
- c. Failure to observe general operating instructions, including warnings, such as the strict application of inspection and maintenance procedures, or the disposal of devices, may result in malfunctions that may induce adverse effects or adverse effects. Serious adverse effects (intended use not achieved), or environmentally hazardous situations.
- d. The device must be used strictly according to the instructions for use: by persons corresponding to the profile of the user, and for indications and intended use as described. Any other use is prohibited and contraindicated by LAMIDEY NOURY MEDICAL as it could lead to dangerous situations that could induce adverse effects (adverse effects or serious adverse effects) or ineffective (intended use not achieved) or dangerous situations for the patient environment.
- e. Like any device, a malfunction and / or failure may occur during the life of the device. Therefore, it is necessary to ensure the proper functioning of the device before use and strictly follow the periodic maintenance as described in the instructions for use. Failure to observe the operations mentioned in the instructions for use may lead to dangerous situations that could lead to a danger (undesirable effects or serious undesirable effects) or a non-effective device (intended use not achieved) or dangerous.

- f. As indicated by the labeling (crossed-out wheeled bin): Disposal must be done with care and must not be disposed of with consumer waste. Waste electrical and electronic equipment (DEEE) can pollute the environment and must be disposed of according to a specific circuit suitable for DEEE there.
- g. ELECTRICAL SAFETY is only guaranteed if the device is correctly connected to the mains and earthed in accordance with the CURRENT safety standards. It is imperative to check this basic condition. If in doubt, ask qualified personnel to carry out a meticulous inspection of all electrical equipment. LAMIDEY NOURY MEDICAL can not be held responsible for damages caused by the lack of an effective grounding of the installation. It is also forbidden to carry out an operation with a device that is not effectively grounded.
- h. If there is any incompatibility between the AC power cable and the power outlet, only replace with a cable that complies with current regulations. It is inadvisable to use adapters, multiple sockets or extenders. If you can not help it, make sure they comply with current safety standards.
- i. In use, the device must be positioned so that the M / A switch on the rear panel remains accessible at all times.
- j. Before using the Argon module, always check the gas flow setting, and adjust it as low as possible according to the desired effect (see § 4.4)
- k. The Argon module and its internal tubing are not sterile. Therefore, to avoid contaminating the patient by the flow of gas, it is imperative to connect a sterile disposable antibacterial filter on the Argon gas outlet on the front of the module, and to connect a user accessory sterile.
- l. Never touch the device with wet or damp hands or feet - Do not use the device barefoot.
- m. Do not leave the device exposed to ATMOSPHERIC AGENTS (rain, sun, etc.).
- n. The use of the device and its accessories is prohibited in EXPLOSIVE or FLAMMABLE ATMOSPHERE. As with any electrosurgical device, it must be ensured that the active electrode is not used in an environment with a concentration of endogenous digestive gases or Oxygen.
- o. Do not leave the device connected unnecessarily. Turn it off at the end of use.
- p. It is dangerous to modify or attempt to modify the characteristics of the device.
- q. Before any cleaning or maintenance operation, disconnect the Argon module either by removing the plug from the socket or by pressing the main switch not a word.
- r. In the event of damage or malfunction, turn off the power by flipping the M / A switch on the rear panel, or disconnecting the AC power cord. For a repair call a service center approved by the manufacturer by requiring original spare parts. Any other solution may compromise the safety of the device and especially that of its user.

### **3.4. Risks of serious adverse events**

The benefits of the device outweigh the occurrence of serious adverse events. However, serious adverse events are reported in the scientific and medical literature, and are reviewed periodically to ensure compliance with the state of the art. Serious adverse events reported and known are:

- a. ARGON CAN ACT AS ASPHYXIAN by displacing air in confined spaces and making the atmosphere dangerous. In the event of a major leak, an atmosphere containing less than 19% oxygen becomes dangerous in a confined space, and presents risks of loss of consciousness and suffocation. In case of discomfort, bring the person back to a normal atmosphere. If necessary, administer oxygen.
- b. During argon plasma coagulation, because of the use of a high voltage current, contact of the electrode with the tissue can cause deep necrosis. AVOID CONTACT ON FABRICS.
- c. If the tip of an Argon electrode or an Argon probe is too close to tissue (distance <1mm), there is a RISK OF DISTENSION of tissue. Maintain a distance between the tip and the targeted tissues during activation of coagulation.

- d. Do not place the tip of an Argon electrode or probe facing a section of a vessel to NOT INJECT GAS INTO THE VASCULAR SYSTEM.
- e. The use of Argon in LAPAROSCOPY presents RISKS OF GAS EMBOLISM because Argon gas is added to the insufflation of CO<sup>2</sup> and may create an overpressure. It is imperative to have a security against the risk of intra-peritoneal overpressure.
- f. Do not direct the argon jet at conductive metal parts: the passage of current through a clip, metal clip, stimulation electrode, or monitoring electrode may cause thermal injury to contact tissue.
- g. The use of argon plasma coagulation may present a HAZARDOUS SITUATION FOR SAFETY at low power settings: there is increased risk of air embolism if there is not enough HF power to quickly produce a waterproof eschar on the targeted tissue.

These serious adverse events should be immediately treated by medical personnel.

#### 4. USE OF THE ARGON MODULE

##### 4.1. Putting into service the argon module

- a. Connect accessories as specified in paragraph 2.2.3
- b. Open the tap of the Argon bottle
- c. Switch on the OPTIMA electrosurgical unit with the on / off switch on the rear panel.
- d. Set the OPTIMA electrosurgical unit: Select the "Argon Plasma" function in monopolar coagulation choices.
- e. Switch on the Argon module with the on / off switch on the rear panel.

##### 4.2. Initialization of the Argon module

When the power is turned on, the device displays a welcome message and then performs a self-test phase (Initialization). The display indicates the software version of the device.

At the end of the initialization, if no anomaly is detected, the device will then display the preset flow rate and the percentage of gas fill (100% = load at 200 bar)

Keys  and  allow to adjust the flow rate from 0 to 10 Liters / minute as follows:


- From 0.3 to 2 l / min, settings in steps of 0.05 l / min
- From 2 to 10 l / min, settings in steps of 0.1 l / min

The key with circular arrow  allows to start the purge of the gas circuit to expel the air from the tubings: a short press triggers a purge of 4 s.


At startup, simultaneously pressing the keys  and  until the end of initialization, the display goes into configuration mode:

Press a key  or  to display the available setting choices:

Press the key  to change the parameter or  to exit the configuration mode.

Press a key  or  to display the available setting choices:


Press the key  to register the setting or  to return to the configuration menu.

Open Surgery mode limits the maximum flow rate to 10 l / min.

Endoscopy mode limits the maximum flow rate to 2 l / min. This mode is intended for applications in digestive endoscopy and bronchoscopy, using a flexible fibroscope.

During activation, the display shows the flow setpoint for 750 ms.

Beyond 750ms, the display will show the flow measured at the output

NOTE: It is possible that the diameter of an argon probe connected does not reach the desired level.

### 4.3. Error Messages

Messages	Instructions
Warning Flow rate : <80%	The application electrode is closed, probably because of tissue contact. The gas flow is interrupted for safety reasons. Release the activation control and move the electrode away to resume normal use. If necessary, clean the tip of the electrode to restore the normal flow of gas.
Command default Please restart	This message means that an activation command occurred during the initialization phase. Check that no command is activated on the electrosurgical unit and check the condition of the Argon control cable on the rear panel.
Warning Empty tank	Check that the cylinder valve is in the open position. Otherwise, replace the empty bottle.
Pressure default Please restart	Check that the pressure sensor is connected properly. After reconnecting the sensor, restart the argon module. If the fault persists, the module must be entrusted to the after-sales service of LAMIDEY NOURY MEDICAL.
Com. default Please restart	Internal communication fault between cards. Restart the argon module. If the fault persists, return the Argon module to the LAMIDEY NOURY MEDICAL Service Department.
Valve default Please restart	Failure of the mass flow valve. Restart the argon module. If the fault persists, return the Argon module to the LAMIDEY NOURY MEDICAL Service Department.
Flow rate default Please restart	The measured flow has an excessive deviation from the expected flow. Restart the argon module. If the fault persists, return the Argon module to the LAMIDEY NOURY MEDICAL Service Department.

### 4.4. Settings

#### 4.4.1. Open surgery

Select the ARGON PLASMA or ARGON PULSE PLASMA mode on the OPTIMA unit.

The power settings for Argon plasma are between 60 and 100 W. The gas flow should be adjusted from 4 to 7 Liters / minutes

Some Argon electrodes allow cutting and coagulation under Argon, without plasma (the tip electrode or the knife protrudes from the tube). For these electrodes, use argon cutting or coagulation currents, and adjust the Argon flow rate to less than 2 L / min.

#### 4.4.2. Digestive endoscopy

Use low power for argon plasma, and FLOW <2L / MIN. Recommended settings:

Interventions	Power	Flow rate
Right colon; Duodenum; Radical Rectitis	20 à 30 W	0,6 L/minute
Esophagus; Stomach	40 W	0,8 L/minute
Left colon	30-40 W	1 L/minute
Tumor	40-50W	1,2 L/minute

#### **4.5. After use**

As a precaution, CLOSE SYSTEMATICALLY THE TAP OF THE ARGON CYLINDRE AFTER USE. This will prevent any leakage from emptying the gas cylinder. Connections on the back of the module can remain in place. Disconnect front-use accessories for cleaning and sterilization.

### **5. MAINTENANCE, CLEANING, STERILIZATION**

#### **5.1. Argon Module**

**CLEANING:** The Argon Module itself does not require any special maintenance other than cleaning the body and front panel, which can be done with a clean, soft cloth soaked in a non-flammable cleaning agent.

**DECONTAMINATION:** use a surface decontaminant such as CIDALKAN wipes (Alkapharm laboratory) or WIP ANIOS or LINGET ANIOS (Anios Laboratory).

The Argon module is not sterilizable.

#### **5.2. Argon accessories**

Clean and sterilize reusable accessories before first use. Before cleaning, separate the sleeve electrodes or the probe from the connection cable.

**CLEANING:** Argon accessories (electrodes, probes, handles and probe connection cables) can be cleaned in an instrument washing machine. See the instructions on the instruction manuals of the accessories concerned.

After cleaning, flush out any excess liquid that may remain in the tubing, using compressed air at a maximum pressure of 4 bar. For this, the compressed air nozzle is loosely fitted to the Luer-Lock connection piece of the connection cable (do not blow from the distal tip to the proximal). Argon sleeves, flexible argon probe connection cables and flexible probes are sterilized in an autoclave up to +134 °C for 20 minutes.

Accessories must remain separate during sterilization. Regular inspection of accessories is recommended (see § 2.3.1).

#### **5.3. Maintenance**

Maintenance at the expense of the user consists mainly in the cleaning and sterilization of accessories and in the control of the proper functioning of the device before use. The annual periodic checks, as well as the safety checks and the verification of the parameters, must be carried out by specialized technicians (see paragraph 7.4):

#### **5.4. Environment protection (Directive 2002-96-CE)**

At the end of their lifetime, the electronic and electrical equipments must be disposed off in accordance with the legislation in force: they must be collected separately and sent to an approved reprocessing or recycling center. LAMIDEY NOURY MEDICAL undertakes to take back at its expense the Argon modules of its manufacture, put on the market in France. Contact LAMIDEY NOURY MEDICAL at 01 69 20 69 69 for any removal request. For devices placed on the European market, contact the local distributor of LAMIDEY NOURY MEDICAL.

## 6. ELECTROMAGNETIC COMPATIBILITY

### 6.1. General instructions

All the information below is derived from the normative requirements to which the manufacturers of electromedical devices are subject, in the sense of IEC60601-1-2 Ed4.

The medical device complies with applicable electromagnetic compatibility standards, however, the user will ensure that any electromagnetic interference does not create an additional hazard, such as radio frequency transmitters or other electronic devices.

In this chapter you will find information necessary to ensure the installation and commissioning of your medical device under the best conditions in terms of electromagnetic compatibility.

The different cords of the medical device must be separated from each other. Certain types of mobile telecommunication devices such as mobile phones are likely to interfere with the medical device. The separation distances recommended in this chapter must therefore be strictly observed.

**WARNING:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

**WARNING:** Use of accessories, transducers and cables other than those specified or provided by the MANUFACTURER of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

The ARGON MODULE is compatible with the OPTIMA electrosurgical unit (Ref V10GOPT4). Its intended use includes this unit placed over the module, the two devices being used together.

In the event of a mains failure, the argon module restarts on the gas flow rate on which it was set before the shutdown.

The use of the ARGON MODULE in an electromagnetic environment that has not been mastered may alter the control of the gas flow and the display of the device.

### 6.2. Cable length

Cable and accessoriess	Maximum length	Type of test	In conformity with
Power cord Control cable	5 m	RF emission	CISPR 11, Classe A
		Harmonic current emissions	Non applicable
		Voltage changes, voltage fluctuations and flicker	Non applicable
		Electrostatic discharge immunity	IEC 61000-4-2
		Radiated, radio-frequency, electromagnetic field immunity	IEC 61000-4-3
	0.42 m	Electrical fast transient/burst immunity	IEC 61000-4-4
		Surge immunity	IEC 61000-4-5
		Immunity to conducted disturbances, induced by radio-frequency fields	IEC 61000-4-6
		Power frequency magnetic field immunity	IEC 61000-4-8
		Voltage dips, short interruptions and voltage variations immunity	IEC 61000-4-11

### 6.3. Electromagnetic emission

The medical device is intended for use in the electromagnetic environment described in the table below. The user and the installer must therefore ensure that the medical device is used in the environment described below.

Emission test	Conformity	Electromagnetic environment - remarks
Radio-frequency disturbance (Radiated emissions) (CISPR 11)	Group 1	The medical device uses RF energy for its internal function.
Mains terminal disturbance voltage (conducted EMISSIONS) (CISPR 11)	Class A	The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
Harmonic current EMISSIONS (IEC61000-3-2)	Non applicable	
Voltage changes, voltage fluctuations and flicker EMISSIONS (IEC61000-3-3)	Non applicable	

#### 6.4. Electromagnetic and magnetic immunity

The medical device is intended to be used in the magnetic and electromagnetic environment described in the table below. The user and installer must comply with the electromagnetic environment.

EXCLUSIONS - The environment describes the following exceptions:

- Installations close to HF ELECTROSURGICAL EQUIPMENT in operation, other than the OPTIMA unit.
- CONTROLLED ACCESS AREA of a MAGNETIC RESONANCE SYSTEM, where the intensity of Electromagnetic DISTURBANCES is high

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) (IEC 61000-4-2)	±8kV contact ±15kV air	±8kV contact ±15kV air	Environment of a professional health care facility
Electrical fast transient/burst (IEC 61000-4-4)	±2 kV for power supply lines	±2 kV for power supply lines	
Surge (IEC 61000-4-5)	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	
Power frequency magnetic field (IEC 61000-4-8)	30 A/m	30 A/m	
Voltage dips, short interruptions and voltage variations (IEC 61000-4-11)	0% UT for 0.5 cycles at 0°, 45°, 90°, 135°, 180°, 225°, 270° et 315°  0% UT for 1 cycle  70% UT for 25 cycles à 50 Hz For 30 cycles à 60 Hz Phase: to 0°	0% UT for 0.5 cycles at 0°, 45°, 90°, 135°, 180°, 225°, 270° et 315°  0% UT for 1 cycle  70% UT for 25 cycles à 50 Hz For 30 cycles à 60 Hz Phase: to 0°	Environment of a professional health care facility.  If the use of the system requires continued operation during mains power interruptions, it is recommended that the medical device be powered from a separate power source (UPS, etc.).
Voltage interruptions (IEC61000-4-11)	0 % UT; for 250 cycles to 50 Hz for 300 cycles to 60 Hz	0 % UT; for 250 cycles to 50 Hz for 300 cycles to 60 Hz	

### 6.5. Electromagnetic immunity, radiofrequency

The medical device is intended to be used in the magnetic and electromagnetic environment described in the table below. The user and the installer must ensure compliance of the electromagnetic environment.

Immunity test	Test level	Compliance level	Electromagnétique environment / remarks
	WARNING: PORTABLE RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ARGON MODULE, including cables specified by LAMIDEY NOURY MEDICAL. Otherwise, degradation of the performance of this equipment could result.		
Radiated radiofrequency electromagnetic fields (IEC61000-4-3)	3 V/m 80 MHz to 2.7 GHz 80 % MA to 1 kHz  10 V/m 80 MHz to 2.7 GHz 80 % MA to 1 kHz	3 V/m 80 MHz to 2.7 GHz 80 % MA to 1 kHz 10 V/m 80 MHz to 2.7 GHz 80 % MA to 1 kHz	Environment of a professional health care facility, excluding installation close to HF ELECTROSURGICAL EQUIPMENT in operation, other than the OPTIMA unit
Proximity fields from RF wireless communications equipment IEC 61000-4-3 (interim method)	9 V/m 710 MHz, 745 MHz, 780 MHz, 5240 MHz, 5550 MHz, 5785 MHz  27 V/m 385 MHz  28 V/m 450 MHz, 810 MHz, 870 MHz, 930 MHz, 1720 MHz, 1845 MHz, 1970 MHz, 2450 MHz	9 V/m 710 MHz, 745 MHz, 780 MHz, 5240 MHz, 5550 MHz, 5785 MHz  27 V/m 385 MHz  28 V/m 450 MHz, 810 MHz, 870 MHz, 930 MHz, 1720 MHz, 1845 MHz, 1970 MHz, 2450 MHz	
Conducted disturbances, induced by RF fields (IEC610004-6)	3 V 150KHz to 80MHz  6 V in ISM band and band between 0.15 MHz and 80 MHz, amateur radio band included 80% MA à 1 KHz	3 V 150KHz à 80MHz  6 V in ISM band and band between 0.15 MHz and 80 MHz, amateur radio band included 80% MA à 1 KHz	

## 7. TECHNICAL DESCRIPTION

### 7.1. Identification of the medical device

See § 1.3 for the meaning of the symbols used for the markings on the device.

Manufacturer	LAMIDEY NOURY MEDICAL 3, rue des petits ruisseaux - ZA des Godets F 91370 Verrières le Buisson (France) Tel +33 (0)1 69 20 69 69 – e-mail : <a href="mailto:info@lamidey-noury.fr">info@lamidey-noury.fr</a>
Device	Argon Module, reference V11A100
Serial number	See the markings on the nameplate on the rear panel
Date of manufacture	
Software	The software version is displayed during the initialization phase

### 7.2. Connection to the mains power network

Voltage, frequency	230 V~, 50 Hz
Absorbed power	60 VA
Earth connection	Via the mains cable
Power switch	Double light rocker switch located at the back of the module.

Note: The device does not deliver output voltage.

### 7.3. Other features

Protection against electric shocks :	Class 1
Applied parts :	Non applicable
Overcurrent protection :	External fuses 2x T 1 AL - 250 V
Physiological effect :	Not applicable
Cooling :	Not applicable
Environmental conditions : transportation and storage:	Temperatures : -20 à +70°C, Relative Humidity : 5 à 95%
Environmental conditions Operating :	Temperature : +10 à +40°C Relative humidity : 0 à 85%
External pressure source:	2 maximum bar (for Argon gas of the medical / surgical type)
Net weight :	19 Kg
Maximum safety load:	30 Kg
Pressure regulator :	Manufacturer: Sté CAHOUEUET - Model NOVA 40 Maximum inlet pressure: 200 bar max Output pressure: 1.7 bar
Argon cylinder recommended :	Manufacturer: Air Liquide Santé - Ref M2001S05R2A001 CH1 / 5 l / ISO connector 5145 N ° 4
Argon flow rate :	Adjustable from 0 to 10 l / min
Precision of the display :	- 0,01 l/min
Precision of measurement :	- 2% of full scale

### 7.4. Maintenance

**WARNING: NO MODIFICATION OF THIS EQUIPMENT IS AUTHORIZED.** The warranty will lapse if any modification is made to the Argon module without authorization. In case of need for repair or adjustment, it must be sent to the technical assistance center of the company LAMIDEY NOURY MEDICAL, with the description of the incident.

#### 7.4.1 Preventive Maintenance

LAMIDEY NOURY MEDICAL recommends AN ANNUAL PREVENTIVE CONTROL of the device, which includes an electrical safety check, and a functional check.

- Required equipment includes:
  - ° An electrical safety tester according to the requirements of IEC 60601-1 or IEC 62353
  - ° A gas flow meter, in the range of 0 to 10 l / min (accuracy <1%)
- The technical staff responsible for preventive maintenance operations must be qualified in biomedical technical maintenance and master the use of these devices.

NOTE: The regulator must be factory tested every 5 years (see next §)

Electrical safety checks: The argon module has no part applied to the patient. The control operations are as follows:

- Check the earth connection via the power cord, between the equipotential terminal and the earth terminal of the mains plug. Compliance: Impedance  $\leq 200 \text{ m}\Omega$ .
- Measurement of contact current, from accessible metal parts - Compliance:

Normal Conditions : $\leq 100 \mu\text{A}$	First default condition : $\leq 500 \mu\text{A}$
--	--

- Measurement of leakage current to earth - Compliance:

Normal Conditions : $\leq 5 \text{ mA}$	First default condition : $\leq 10 \text{ mA}$
---	--

#### Functional checks:

##### *Checks,*

- Operation of the screen (backlight and readability of information)
- Operation of keys  and 
- Operation of key purge 
- Indication of the filling level of the bottle (in%)
- Error message if empty bottle (= gas valve closed on power up)
- Air proof of the gas circuit (no leak noise)

##### *Gas flow control:*

- Connect a tubing between the gas outlet fitting and the control flow meter.
- Measure the flow rate for the displayed values by incrementing 0.5, from 0 to 2 l / min and incrementing by 1, from 2 to 10 l / min. Compliance range: difference between measurement and display  $\leq 0.2 \text{ l / min}$ .


### 7.4.2 Replacement of spare parts

NOTE: Electrical diagrams, component list, adjustment instructions and all necessary information to assist biomedical technical personnel in the troubleshooting of the ARGON module are available upon request from LAMIDEY NOURY MEDICAL.

The following spare parts are replaceable by hospital biomedical technical staff. They are available from the after-sales service of LAMIDEY NOURY MEDICAL:

- Mains supply fuses:

The 2 fuse holders are located on the back, in the block the M / A switch, the power base and the fuse holders.  
To access the fuses, use a flat screwdriver and make a quarter turn to the left

WARNING: the protection provided by the fuses is guaranteed only if they are replaced by fuses with the same features (2x T 1 AL - 250 V).

- Power cord :

. The connection to the ground must be regularly checked. In the event of a fault, the mains cable must only be replaced by a cord conforming to the IEC 60227 standard, equipped with 3 conductors with a cross-section of at least 1mm<sup>2</sup> each, and a maximum length of 5m.

- Argon gas cylindre :

FAILURE TO COMPLY WITH THE INSTRUCTIONS FOR USING THE ARGON BOTTLE MAY RESULT IN HAZARDOUS SITUATIONS. The gas cylinder must only be replaced by an equivalent bottle, containing pure Argon, of medical-surgical quality, under a maximum pressure of 200 bar, and an ISO 5145, N ° 4 connection. (See the installation of the bottle in 2.2.2).

- Gaz pressure reducer

The regulator must be returned to LAMIDEY NOURY MEDICAL every 5 years for preventive control. THIS DEVICE MUST BE REPLACED EVERY 15 YEARS, by a regulator meeting the specifications of LAMIDEY NOURY MEDICAL. If the pressure sensor is defective, return the regulator + sensor assembly to the LAMIDEY NOURY MEDICAL Service Department.

DO NOT USE PLIERS for mounting the regulator on the bottle. This could damage the seal. Tighten only by hand.

- Tubing (Ref A409 0345C) :

The gas connection between the regulator and the Argon module uses a special pressure-resistant tubing, which is compatible with the quick connectors.

- Control cable (Ref A409 0312) :

This shielded cable transmits the gas activation command between the Optima unit and the module. It must be replaced in case of failure.