# GIMA

PROFESSIONAL MEDICAL PRODUCTS

## **ASTE PORTAFLEBO I.V STANDS PIED À PERFUSION INFUSIONSSTÄNDER POSTE DE GOTEIO INTRAVENOSO HASTES SUPORTE DE SORO**

**27841/27843**

**27842**

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In accord with Regulation (EU) 2017/745, this product is part of the CLASS I MEDICAL DEVICES. You should carefully read this manual before using this medical device.

**DESCRIPTION**

Telescopic I.V. stand with 2/4 hooks for the support of bags and infusions. Composed of two coaxial tubes in chromed or aluminum steel. Tubular adjustable up to 210 by plastic regulation crank. Equipped with black plastic base with five wheels Ø 50 mm. The device is an aid for the care of the patient who is bedded or with limits functional autonomy.

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**ASSEMBLY INSTRUCTIONS:**

1. Take the wheels and insert them into the base.
2. Take the tubular with a larger section (25 x 1.2) equipped with a knob, insert it first into the black cap, then into the base, fix the parts, and screw it into the area below.
3. Take the tubular section with 20 x 1 and insert the tube inside the other until it stops using the plastic regulation crank.
4. For 27841/27843 codes, take the plastic support with two hooks and press to insert the cap at the end of the tubular. If necessary, use a soft hammer to insert it all completely.
5. For 27842 code, fix the support with 4 hooks by castor supplied.

**TECHNICAL DATA**

<b>Measurements</b>	27841 – 27842	Adjustable height 160-210 cm
	27843	Adjustable height 165-210 cm
<b>Weight</b>	27841	2,900 kg
	27842	3,150 kg
	27843	2,080 kg
<b>Max capacity</b>	27841 – 27843	Plastic hook 5 kg each hook
	27842	Plastic hook 2 kg each hook

**WARNINGS!**

- Fix well with the supplied screws.
- Check the locking nuts are tightened.
- Make sure that the parts are locked in their housing.
- Observe the maximum flow rates indicated in the data sheets.
- Avoid direct contact with open wounds.
- Not for use in open environments.

**MAINTENANCE AND CLEANING**

- Use only original spare parts.
- Periodically check all fixing points.
- To preserve at best, it is recommended to avoid prolonged exposure to sunlight, to contact and deposit in saline and very humid environments.
- Wash the stainless steel surfaces with soapy water and use appropriate disinfectants.
- Rinse with a damp cloth and dry thoroughly.
- Do not spill on the surface cleaners containing chlorine or ammonia, fatty and oily substances.
- Do not use cleaning utensils containing iron. If necessary disinfection only use suitable non-corrosive products.

	Keep away from sunlight		Medical Device compliant with Regulation (EU) 2017/745		Keep in a cool, dry place
	Product code		Consult instructions for use		Manufacturer
	Date of manufacture		Lot number		Medical Device
	Caution: read instructions (warnings) carefully				

**GIMA WARRANTY TERMS**

The Gima 12-month standard B2B warranty applies.

## DECLARATION OF CONFORMITY

We, undersigned GIMA S.p.A., with operational headquarters in Gessate (MI), Via Marconi 1, and registered office in Milano, Via Tommaso Grossi 2, acting as manufacturer of the medical device:

GIMA Single Registration Number (SRN):

<b>Medical Device ( Trade Name and description)</b>	<b>Code</b>	<b>Basic UDI-DI code</b>
I.V.STAND ON 5 WHEELS TROLLEY - 4 hooks	27842	80232790000V90994500000C3

Risk class I (Not sterile), according to the Rule 1 Annex VIII of Regulation (EU) 2017/745 (MDR), declares, under its own responsibility, that this medical device:

- comply with essential requirements and dispositions of Regulation (EU) 2017/745 (MDR), as from the Technical File filed at the company;
- common Specifications have not been used for the compliance of the above medical device;


Gessate, 5/28/2021

**GIMA S.p.A.**

The legal Representative  
(Nicola Manzoni)

A handwritten signature in black ink, appearing to read "N. Manzoni", is written over a horizontal line.

# I.V. STANDS

GIMA code	IV STANDS	Adjustable height	Weight kg	Base Material	Hook		Load kg		Castors with brake
					N	Material	Per hook	Total	
27839	Stainless steel AISI 304 IV stand	160-210 cm	3.3	stainless steel	4	stainless steel	1.2	5	1
27841	Chromed-steel IV stand	160-210 cm	2.9	plastic	2	plastic	5	10	-
27842	Chromed-steel IV stand	160-210 cm	3.2	plastic	4	metal	2	8	-
27843	Aluminium IV stand	165-210 cm	2.1	plastic	2	plastic	5	10	-
27840	Epoxy painted IV stand	160-220 cm	4.0	painted steel	2	plastic	5	10	-
27844	Stainless steel AISI 304 IV stand	135-215 cm	2.9	stainless steel	4	plastic	1.2	5	1
27820	Stainless steel AISI 202 IV stand	135-215 cm	3.5	aluminium	4	stainless steel	3	12	1

**ACCESSORIES FOR IV STANDS**

27845	Plastic support for bottle holder - spare	2	plastic	5	10	 <b>27849 - optional</b>
27846	Metal support for bottle holder - spare	4	metal	2	8	
27849	Stainless steel AISI 304 support for bottle holder - optional	2	stainless steel			
27847	Plastic bottle holder - box of 40 pcs.	-				
27848	Stainless steel AISI 304 bottle holder for bottle ≥ 1 l				2.5	

# INFUSION STANDS

- 27826 INFUSION STAND
  - 27829 INFUSION STAND with hydraulic pump
- Made of chrome plated steel. Four chrome plated steel hooks. Five spoke stainless steel base Ø 60 cm with castors Ø 50 mm, two with brake. In conformity with MDR.

TECHNICAL SPECIFICATIONS		
	27826	27829
Adjustable height (cm):	160-250	189-260
Load (kg)*:	30	30
Weight (kg):	10	18
Frame for infusion pump (cm):	— 39x75 —	— —
Packaging size (cm):	— 185x18x43 —	— —

\*Load 5 kg each hook, 10 kg frame for infusion pump

# HYPODERMOCLYSIS STAND

- 27857 HYPODERMOCLYSIS STAND
- Stand with hydraulic height adjustment for hypodermoclysis, for washing and surgical irrigation. Hydraulic pump with pedal for height adjustment and gripping handle for easy transport.

Stainless steel adjustable gripping handle for easy transport

**LOAD 28 KG**

Pedal driven shock absorbing pump. Height adjustable: 147-247 cm. Weight: 22 kg

Stainless steel base with 5 pivoting castors Ø 75 mm (2 antistatic with brake, 3 standard) and grey vipla bumpers.