Generic Name : Microalbumin Kit (Immunoturbidimetric Assay Method)

Abbreviated name: MALB

Order Information

Cat. No.	Package size
MAL0202	R1 2×18 mL + R2 2×5 mL
MAL0203	R1 3×36 mL + R2 1×20 mL
MAL0204	R1 3×50 mL + R2 3×11 mL
MAL0302	R1 2×18 mL + R2 2×5 mL + Cal 5×1 mL
MAL0303	R1 3×36 mL + R2 1×20 mL + Cal 5×1 mL
MAL0304	R1 3×50 mL + R2 3×11 mL + Cal 5×1 mL
MAL1302	R1 1×30 mL + R2 1×7 mL + Cal 5×1 mL
MAL1303	R1 1×30 mL + R2 1×7 mL + Cal 5×1 mL
MAL1304	R1 1×30 mL + R2 1×7 mL + Cal 5×1 mL

Intended use

In vitro test for the quantitative determination of MALB concentration in urine on photometric systems.

Summary

A healthy kidney infiltrated and excreted unnecessary albumin which sustain in blood as valuable protein. Microalbumin (MALB) is defined as a condition characterized by urinary albumin excretion above 20 mg/L in normal adults. Investigations have established a range of 20 mg/L-200 mg/L as diagnostic for MALB.

The presence of MALB in the urine is a reliable and sensitive evidence to diagnose earliest stage of nephropathy. MALB can be considered as an indicator of deteriorating renal function in diabetic subjects. It also may beneficial to control diabetes treatment program.

Method

Immunoturbidimetric Assay Method

Reaction Principle

When the sample is added to the R1 and R2, turbidity changes caused by antigen-antibody reaction. The turbidity is proportional to albumin quantity in the sample. The actual concentration is then determined by a calibration curve prepared from calibrators of known concentration.

Reagents

Components and concentration

R1	Tris buffer	25 mmol/L
	Preservative	0.5 g/L
R2	Tris buffer	20 mmol/L
	goat anti-human MALB antibody	500 mg/L
	Preservative	1 g/L

Warnings and precautions

1. For in vitro diagnostic use only.
2. All products derived from human source have been tested HBsAg, HIV antibody, and HCV antibody negative. However, as no testing method can rule out the potential risk of infection with absolute certainty, the material should be considered potentially infectious.
3. Take appropriate precautions for the use of laboratory reagents. Do not swallow and avoid contact with skin and mucous membranes.
4. When the reagent accidentally enters the eyes and mouth, or contact with the skin, immediately wash with plenty of water. If necessary, visit the doctor for further medical treatment.
5. Disposal of all waste material should be in accordance with local guidelines.
6. Do not allow to mix reagents with different lots.
7. Avoid direct exposure to sunlight and freezing. The results can not be assured when stored at inappropriate condition.
8. Do not use the product after its expiration date, as the reliability of the measured values cannot be assured after this date.
9. Material safety data sheet is available for professional user on request.

Reagent Preparation

R1: Ready to use

R2: Ready to use

Calibrator Preparation

Calibrators are ready to use.

Storage and stability

Up to expiration date indicated on the label, when stored unopened at 2-8°C and protected from light.

Once opened, the R1 and R2 reagent is stable for 30 days when refrigerated on the analyzer or refrigerator.

Once opened, the calibrator is stable for 14 days at 2-8°C and protected from light.

Contamination of the reagents must be avoided.

Do not freeze the reagents.

Reagent blank absorbency

The absorbance of MALB reagent blank at 340 nm should be <0.3 A.

Materials required but not provided

1. Calibrators: MALB Calibrator, Cat.No.OSA181 or within MALB kit.
2. Controls: MALB Control, Cat.No.OQA181.
3. NaCl solution 9 g/L.
4. General laboratory equipments.

Specimen collection and preparation

- a)Urine is suitable for samples.
 b)Samples Stability: 7 days at 20-25°C
 1 month at 2-8°C
 6 months at -20°C (only freeze once)

Assay procedure

	Blank	Sample
Reagent 1	250 µL	250 µL
Dist. water	16 µL	—
Sample	—	16 µL
Mix	and incubate at 37°C for 4 min, read the blank absorbance(A Blank), then add:	
Reagent 2	40 µL	40 µL
Mix	thoroughly at 37°C, and read the absorbance (A Sample) again 5 min. later.	
$\Delta A = [A \text{ Sample}] - [A \text{ Blank}]$		

Calibration

- 1.It is recommended to use the Mindray calibrators and 9 g/L NaCl for multipoint calibration.
- 2.Traceability of the calibrator: This method has been standardized against ERMDA470K.
- 3.Calibration frequency:30days
 It is recommended to calibrate when new working reagent is made.
 After reagent lot changed.
 As required following quality control procedures.

Quality control

One level control material should be analyzed with each batch of samples. In addition, the control should be run with each new calibration, with each new reagent cartridge, and after specific maintenance or troubleshooting procedures as detailed in the appropriate system manual.

We recommend using the supporting control to verify the performance of the measurement procedure; other suitable control material can be used in addition.

Each laboratory should establish its own internal quality control scheme and procedures for corrective action if control can not recover within the acceptable tolerances.

Calculation

The analyzer calculates the MALB concentration of each sample automatically after calibration.

Conversion factor: $\mu\text{g/L} \times 1000 = \text{mg/L}$

Reference Intervals

Each laboratory should establish its own reference intervals based upon its patient population. The reference intervals measured at 37°C listed below

were taken from literature:

Sample Type	S.I. Units
Urine	≤20 mg/L

Performance Characteristics

Representative performance data obtained from Mindray system (Mindray BS series analyzers / Mindray MALB Reagent) is given below. Results may vary if a different instrument, an individual laboratory or a manual procedure is used.

Limitations-interference

The following substances were tested for interference with this methodology. Criterion: Recovery within ±10 % of initial value.

Substance	Level Tested	Observed Effect
Acetone	350 mg/dL	NSI
Ascorbic acid	30 mg/dL	NSI*
Bilirubin	66 mg/dL	NSI*
Creatinine	500 mg/dL	NSI*
Glucose	2000 mg/dL	NSI*
Hemoglobin	300 mg/dL	NSI
Urea	4200 mg/dL	NSI*
Urobilinogen	20 mg/dL	NSI*

* NSI: No Significant Interference (within ±10%)

Linearity Range

The Mindray System provides the following analytical ranges:

Sample Type	S.I. Units
Urine	4~300 mg/L

If the value of sample exceeds 300 mg/L, the sample should be diluted with 9 g/L NaCl solution (e.g. 1+ 9) and the result should be multiplied by 10.

Analytic Sensitivity/Limit of Detection

The lowest measurable MALB concentration that can be distinguished from zero is 4 mg/L with 99.7% confidence.

Precision

Precision performance using the CLSI Approved Guideline EP5-A2 to assay control appears in the table below. U: mg/L

Type of Imprecision	Level 1			Level 2		
	Mean	SD	CV (%)	Mean	SD	CV (%)
Within-run		0.27	1.31		2.60	0.72
Between-run	20.86	0.40	1.90	361.81	1.343	0.37
Between-day		0.46	2.19		1.65	0.46
Within-device		0.66	3.18		3.36	0.93

Method Comparison

A comparison between Mindray System (Mindray BS series analyzers

/Mindray MALB Reagent) (y) and Hitachi System (Hitachi analyzer/ MALB Reagent) (x) using 40 samples gave following correlation (mg/L):
 $y=1.0224x-0.6087$, $R^2=0.9987$

Details of the comparison experiments are available on request.

References

- 1.Rothschild MA, Oratz M, Schreiber SS. Serum albumin. Hepatology 1988; 8:385-401.
- 2.Hofmann W, Guder WG. A diagnostic program for quantitative analysis of proteinurea. J Clin Chem Clin Biochem 1989; 27:589-600.
- 3.Greiling H, Gressner AM,eds. Lehrbuch der Klinischen Chemie und Pathobiochemie,3rd ed. Stuttgart/New York: Schattauerr, 1995:223-224, 749-750.
- 4.Guder WG, Narayanan S, Wisser H, Zawta B. The Quality of Diganostic Samples. Brochure in: Samples: From the patient to the Laboratory,2nd edition. Darmstadt: GIT Verlag,2001.

Graphical symbols

In Vitro Diagnostic
medical device

Batch Code

European
Conformity

Authorized representative in
the European Community

Use By

Consult
Instructions for use

Biological
Risks

Temperature
Limit

Manufacturer

Catalogue
number

© 2017 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building,Keji 12th Road South,Hi-tech industrial park,Nanshan,Shenzhen 518057,P.R.China

Website: www.mindray.com

E-mail Address: service@mindray.com

Tel:+86-755-81888998

Fax:+86-755-26582680

EC-Representative: Shanghai International Holding Corp.
GmbH(Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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© 2017 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Tüm hakları saklıdır.
Üretici: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
AT Temsilcisi: Shanghai International Holding Corp. GmbH(Europe)