1. DESCRIPTION

Embosoft is a medical device that is done to promote the vascular occlusion after selective delivery by microcatheter. The product is composed of microspheres made of polyvinyl alcohol (PVA) and bovine gelatin dispersed in saline, stored in a cyclic olefin copolymer (COP) syringe, packaged individually in a polycarbonate blister thermally sealed with Tyvek® and sterilized by steam. Each syringe is intended for single use only on one patient and should not be re-sterilized. Embosoft is available in 1 or 2 mL of microspheres, completed with 9 or 8 mL of saline and in five size range: 100-300 μm , 300-500 μm , 500-700 μm , 700-900 μm , 900-1100 μm .

2. INDICATIONS

Embosoft microspheres can be used for vascular occlusions being indicated for use in:

- Embolization of hypervascular tumors and process, including uterine fibroids, meningiomas, etc.
- Embolization of the prostate arteries for relief of symptoms related to Benign Prostatic Hyperplasia.
- Haemostatic embolization.

3. CONTRAINDICATIONS

- Patients who do not tolerate occlusion procedure;
- Positioning the microcatheter in the target lesion is not possible;
- The arteries leading to the injured part are small enough not to allow the passage of the microspheres;
- Presence of probable principle of vasospasm;
- Final orientation of arteries directed to the cranial nerves;
- Presence of extra anastomosis to the pronounced intracranial and deviations;
- High-flow arteriovenous shunts or with a diameter greater than the selected microspheres.
- Use in the pulmonary vasculature.
- Severe atherosclerosis.
- Patients with known allergy to gelatin.
- 100-300 μm microspheres are not recommended for use in the bronchial circulation

Precautions, restrictions, warnings, special care and clarifications on the use of the medical products as well as its storage and transport.

4. WARNINGS

The product contains material of animal origin.

These warnings apply to any treatment of embolism in blood vessels at high risk procedures:

1. The procedure must be performed by or under the supervision of a doctor who has been specifically trained

to perform these procedures and has knowledge of angiography and embolization techniques.

2. An evaluation of the patient's condition, along with the vascular assessment of the injury and the place must be made before the embolization procedure.
3. A preliminary assessment of angiography shall be done to determine the potential danger of collateral damage, visible in extra and intracranial areas as well as in area of brain nerves. The Embolization shall not be performed before the areas, cranial nerves or veins are protected.
4. Appropriate equipment shall be evaluated to avoid occurrence of complications during the procedure.
5. Select a microsphere in an appropriate size according to the size delineated by the angiography of the lesion to be treated. Improper choice of the size of the microspheres is a potential cause of its entry in an inappropriate place, preventing vascularization of other normal tissues, or other vessel out of the lesion.
6. Ensure the appropriate suspension of microspheres in contrast solution to prevent the microspheres to clog the infusion microcatheter. A higher concentration of microspheres may also clog the microcatheter during the infusion.
7. If the microcatheter clog during the procedure, do not try unblocking it with infusion injection or injection of pressure while inside the body. Remove the microcatheter and determine the extent of the obstruction. Clean the microcatheter before reintroducing

- it.
8. Special care shall be taken to avoid introduction of contaminants in the angiography room.
9. Incomplete occlusion of the vascular territory can increase the possibility of post-operative hemorrhage, ischemia, heart attack, developing alternative vascular disorders, or recurrence of symptoms.
10. Inadvertent flow for cranial arteries can create a potential ischemic stroke in normal vessels.
11. The patient monitoring to assess the ongoing level of vascular occlusion may be clinically indicated.
12. Do not re-use. For single use of the patient. Discard after the procedure.
13. Scitech may not be deemed responsible for any consequential, incidental or direct damages, resulting from the re-use of this product.
14. Do not connect a 20mL syringe directly to the microcatheter to release Embosoft because it can occur its occlusion.

5. PRECAUTIONS

- Inspect the packaging process before use to ensure its integrity.
- Do not use after expiration of the specified date on the package. Do not use if the package is damaged, opened or violated.
- Discard any unused material.

6. ADVERSE POTENTIAL COMPLICATIONS

Specific complications on embolization procedures may occur during or after the procedure and may include, however not limiting to the following:

- Ischemia in unwanted place
- Ischemic attack or ischemic stroke
- Undesired reflux or passage of microspheres in normal vessels adjacent to the lesion
- Pulmonary embolism
- Rupture or injury to vessels and bleeding
- Neurological deficit
- Paralysis of the cranial nerve
- Damage to vessels and tissues
- Infection
- Allergic reaction
- Bleeding
- Spasms in the vessels
- Migration of microspheres causing stroke or cerebral infarction
- Recanalization by collateral vascular pathways
- Death
- Migration of microspheres causing stroke or cerebral infarction

7. FORMS OF PRESENTATION

Each syringe contains individually 1 or 2ml of Embosoft microspheres with 8 or 9ml of saline solution for each total volume of 10ml.

Sterilization

Sterilized by steam. Do not use if sterile barrier is

breached. Do not re- use or re-sterilize.

8. PACKAGING AND STORAGE

Store in a cool and dry place, between temperatures of 15 to 25°C.

The product will be packaged in 20 mL syringes made of cyclic olefin copolymers (COP), each of which contains 1 or 2 ml of microspheres for embolization + 9 or 8ml of saline solution for each total volume of 10ml, individually packed in a polycarbonate blister and heat sealed with Tyvek®. It is delivered to the sterile, and nonapyrogenic consumption.

Packaging:

Primary package: Assembled syringe and blister heat sealed with Tyvek

Secondary package: Cardboard Box.

Size range of microspheres Embosoft and color code:

Size Range (µm)	Color Code
100-300	Yellow
300-500	Blue
500-700	Red
700-900	Green
900-1100	Purple

9. TRACEABILITY

Because it is an implantable material, as set out in art. 4th. IX c / c art. 18, Paragraph 1st, both from RDC 14/2011 (Brazil local laws), are available from the

manufacturer detachable labels sufficient for the traceability of the product with the identification of the material inside the package.

10. INSTRUCTIONS

- Perform an initial angiogram to determine the vascular supply of the injury. The angiogram must demonstrate the access route to the microcatheter and identify relevant side circulations.
- Position the microcatheter using standard technique. Enter as close as possible to the target area to prevent improper release of microspheres to normal veins.
- Select an appropriate size of microspheres, based on observation of the angiogram of the treatment area.
- Visually inspect the integrity of the packaging to ensure sterility. Do not use if packaging is damaged or broken. Inspect the microspheres visually to check size uniformity.
- Mix the EMBOSOFT microspheres with the appropriate contrast agent until the desired dilution is achieved. Use a 3-way tap to inject contrast in the syringe EMBOSOFT. Ensure that the microspheres are equally suspended and floating free.
- Remove all the air from a syringe 1 to 3cc and connect it to the three-way tap.
- Transfer the EMBOSOFT syringe solution partially to the syringe 1 to 3cc. Check that concentration and microsphere suspension prior to injection.
- Remove all air from the syringe (1 to 3 cc) solution

containing the microspheres Connect the syringe (1 or 3cc) microcatheter.

- Release the microspheres through the microcatheter infusion under direct fluoroscopic control. The injection speed must be controlled so as not to occur reflux.
- Continue releasing microspheres until achieved the desired devascularization.
- Remove the microcatheter out of the infusion, while applying gentle suction to prevent the microspheres to remain in the lumen of the microcatheter.
- Remove any remaining Microspheres EMBOSOFT and used syringes.

Do not use the syringe 20 ml pre-filled to inject the Embosoft directly through the microcatheter.

11.EXPIRY DATE

The expiry date of the product is 2 years from sterilization date.

12.LABEL SYMBOLS

The following symbols are used:

No.	Description	Symbol
1.	CATALOGUE NUMBER	REF
2.	BATCH CODE	LOT
3.	USE BY (EXPIRY DATE)	
4.	STERILIZED USING STEAM OR DRY HEAT	STERILE
5.	DO NOT REUSE	
6.	DO NOT RESTERILIZE	

7.	CAUTION	
8.	TEMPERATURE LIMITATION	15°C 25°C
9.	CONSULT INSTRUCTIONS FOR USE	
10.	DO NOT USE IF PACKAGE IS DAMAGED	
11.	KEEP AWAY FROM SUNLIGHT	
12.	KEEP DRY	
13.	NON-PYROGENIC	
14.	MANUFACTURER	
15.	CE MARKING WITH IDENTIFICATION NUMBER OF THE NOTIFIED BODY	CE 1434

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