ASO/CRP/RF Triple Control

Order Information

Cat. No.	Package size
105-004650-00	L: 3×1 mL + H: 3×1 mL

Intended Purpose

ASO/CRP/RF Triple Control is used for quality control of quantitative determination of rheumatism analytes by monitoring accuracy and precision of Mindray BS series chemistry analyzers and test ability of clinical laboratory.

Summary

ASO/CRP/RF Triple Control contains the following analytes:

Antistreptolysin "O" (ASO)	Rheumatoid Factor (RF)	C-reactive Protein (CRP)
-------------------------------	---------------------------	-----------------------------

Mindray BS measurement system is composed of Mindray BS series chemistry analyzers, Mindray reagent kits, calibrators and control materials.

Components

ASO/CRP/RF Triple Control is a liquid control based on ASO, RF and CRP solution which contains 2 levels: L, H with different and specific concentration, and 0.095% Sodium azide and 0.1% ProClin300 as preservative. The control values are lot-specific with the matched models listed in the value sheet.

Storage and stability

The control is stable up to expiration date indicated on the label when stored in unopened vial at 2-8°C and protected from light.

Once opened, it is stable for 5 days when capped tightly, avoiding microbial contamination and protected from light at 2-8°C.

Please store control tightly capped when not in use.

Preparation

1. Control is ready to use.
2. It is recommended to take the vial out of the refrigerator and allow reaching room temperature.
3. According to the quality control procedure described in the operation manual for the analyzers.
4. Dispense the required volume into a sample cup, and then run QC test you have requested or have it analyzed in the same way as patient samples.
5. Quality control is recommended to run daily in parallel with the patient samples, after new calibration or reagent cartridge change, after specific maintenance or trouble shooting procedures.

6. Please perform scheduled maintenance and standard operation including calibration and analysis to assure the performance of measurement system.

Materials required but not provided

1. General laboratory materials: NaCl solution 9 g/L (saline).
2. Mindray reagent Kit and Mindray calibrator. The catalogue number of reagent and calibrator used with ASO/CRP/RF Triple Control can be obtained from the corresponding Mindray reagent kits and calibrators instruction for use.
3. Mindray BS series chemistry analyzers, and General laboratory equipment.

Performance Characteristics

■ Assay Values

The control assay values (the target value and range) determined by Mindray standard transfer procedure and routine method are listed in the value sheet. The target value was obtained from Mindray measurement system and the range was calculated as the target values ± 3 standard deviations. The control values are specific for different lot numbers and different model of chemistry analyzers, so please check the lot number and model of analyzer before use.

Mindray systems (BS series chemistry analyzers/Mindray reagents) are calibrated by Mindray calibrators, then test the control on Mindray systems, the assayed values should be within the defined range.

■ Homogeneity

Analytes	Intra-vial homogeneity	Inter-vial homogeneity
ASO	SD $\leq$ 8.0 IU/mL or CV $\leq$ 6.0%	SD $\leq$ 11.97 IU/mL or CV $\leq$ 8.98%
FR-CRP	SD $\leq$ 0.5 mg/L or CV $\leq$ 5.0%	SD $\leq$ 0.75 mg/L or CV $\leq$ 7.48%
RF	SD $\leq$ 2.0 IU/mL or CV $\leq$ 8.0%	SD $\leq$ 2.99 IU/mL or CV $\leq$ 11.97%

Result interpretation

If the control result falls outside the range, the measurement system should be checked. For instance, position of reagent or sample placed on the analyzer; expiration date or storage condition of the calibrator, reagent and control; the parameter setting in software; performance of the analyzer; or the calibration process. Each laboratory should establish its own internal quality control scheme and procedure for corrective action.

Warnings and precaution

1. For in vitro diagnostic use only. For laboratory professional use.

2. Please Take the necessary precautions for handling all laboratory reagents.
3. Please confirm the integrity of the package before use. Do not use the controls with damaged packages. The controls avoid direct exposure to sunlight and freezing. The results can't be assured when stored at inappropriate condition.
4. If unintentionally opened before used, store the controls tightly capped at 2-8°C and protected from light, and the stability is equally to in-use stability.
5. Do not mix controls with different lots and bottles.
Do not use the control beyond the expiration date and the in-use date.
Do not mix new opened control with in-use control.
Avoid the formation of foam.
6. Instability or deterioration should be suspected if there are visible signs of leakage, precipitates or microbial growth, or if calibration or controls do not meet the insert and/or the Mindray System criteria.
7. Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
8. Preservative contained. Do not swallow. Avoid contact with skin and mucous membranes.
9. When the control accidentally enters the eyes and mouth, or contact with the skin, immediately wash with plenty of water. If necessary, visit the doctor for further medical treatment.
10. Safety data sheet is available for professional user on request.
11. Disposal of all waste material should be in accordance with local guidelines.
12. The control was tested with CE-marked methods and shown to be negative for antibodies to HIV and HCV and nonreactive for HBsAg. However, as no testing method can rule out the potential risk of infection with absolute certainty, so this material should be handled as a patient specimen¹.
13. All human material should be considered potentially infectious.
14. All identified risks have been reduced as far as possible without adversely affecting the benefit-risk ratio, and the overall residual risk is acceptable.
15. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.
16. This kit contains components classified as follows in accordance with the Regulation (EC) No 1272/2008:

**Warning**

H317	May cause an allergic skin reaction.
------	--------------------------------------

H412	Harmful to aquatic life with long lasting effects.
------	--

Prevention:

P280	Wear protective gloves and protective clothing.
------	---

P261	Avoid breathing mist/vapours/spray.
------	-------------------------------------

P273	Avoid release to the environment.
------	-----------------------------------

P272	Contaminated work clothing should not be allowed out of the workplace.
------	--

Response:

P302+P352	IF ON SKIN: Wash with plenty of water.
-----------	--

P333+P313	If skin irritation or rash occurs: Get medical advice/attention.
-----------	--

P362+P364	Take off contaminated clothing and wash it before reuse.
-----------	--

Disposal:

P501	Dispose of contents/container to authorised hazardous or special waste collection point in accordance with any local regulation.
------	--

References

1. Biosafety in Microbiological and Biomedical Laboratories; U.S. Department of Health and Human Services; US Government Printing Office; Washington:2007.

Graphical symbols

In Vitro Diagnostic
medical device

Unique device
identifier

European
Conformity

Consult Instructions
For use

Use-by
date

Authorized representative in
the European Community

Batch Code

Temperature
limit

Manufacturer

Catalogue
number

Control

Biological Risks

Keep away from sunlight
Indicates a medical device that needs
protection from light sources

© 2022-2023 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-81888998; **Fax:** +86-755-26582680

EC-Representative: Shanghai International Holding Corp. GmbH
(Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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Fabricant : Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Représentant en Europe : Shanghai International Holding Corp. GmbH (Europe)