[EN]PROBE CLEANSER	
[Packing Specifications]	
50mL×1	
[Intended Use]	
It is used to clean the instrument regularly.	
[Test Principle]	
Not applicable.	
[Major Components]	
Surfactant	≤0.2%
Sodium hypochlorite	≤12%
Sodium hydroxide	≤5%

- [Storage and Stability]
- It is important to store in dark place and keep container closed when not in use.
 - The product will be stable for 18 months when stored at 2°C to 30°C (35.6 °F to 86 °F) in a dark place, with relative humidity not above 90%.
 - The working temperature range of the product is consistent with that of its applicable instruments.
 - The open-vial validity is 60 days.

[Applicable Instruments]

The product applies to Hematology Analyzers, Automated Glycohemoglobin Analyzers, Auto Slide Maker & Stainer, Specific Protein Analyzers, Auto Specific Protein and ESR Analyzers, Chemistry Analyzers, Chemiluminescence Immunoassay Analyzers, Flow Cytometers and Auto Coagulation Analyzers manufactured by Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

[Sample Requirements]

Not applicable.

[Materials Required but Not Provided]

The following materials are required but not provided with the product: Mindray-manufactured measurement instruments as well as matched reagents.

[Test Procedure]

For detailed usage instruction, please refer to the Operator's Manual of the instruments.

[Cut-off Value/Reference Intervals]

Not applicable.

[Result Elaboration]

Not applicable.

[Limitations]

Not applicable.

[Performance Characteristics]

Not applicable.

[Precautions and Warnings]

- For *in-vitro* diagnostic use only. For laboratory professional use.
- Read the package insert carefully before using this product. It shall be used before the expiry date and disposed of properly when expired.
- Confirm the integrity of the package before use. Do not use the product if the package is damaged. Otherwise, the test results may not be accurate.
- Do not use the reagent if it is frozen.
- If the reagent is polluted or affected by other factors and becomes abnormal, stop using it and replace it with a normal one.
- Dispose of the waste, residual and contaminated packing based on local regulations.
- The following factors can affect the sample analysis: expired or ineffective reagent; reagent polluted by dust in the air; improper disposal of the sample; mixed with or used with reagents produced by other company; mixed use of residual from the old container and the newly-opened; used in condition other than specified.
- Take the necessary precautions for the use of the product. Do not swallow. Avoid contact with skin and mucous membranes. If you accidentally take the reagent into your mouth, or the reagents accidentally spill on your skin or into your eyes, wash them off with plenty of water and go seek medical treatment if necessary.
- Disposal of waste liquid and materials should be in accordance with local guidelines.
- The Material Safety Data Sheet (SDS) is available upon request.
- All identified risks have been reduced as far as possible by generally acknowledged state of art, and the overall residual risk is acceptable.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.
- This product contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

  Danger	
H411	Toxic to aquatic life with long lasting effects.
H400	Very toxic to aquatic life.
H290	May be corrosive to metals.
H314	Causes severe skin burns and eye damage.
EUH031	Contact with acids liberates toxic gas.
EUH206	Warning! Do not use together with other products. May release dangerous gases (chlorine).
Prevention	
P260	Do not breathe mist/vapours/spray.
P264	Wash all exposed external body areas thoroughly after handling.
P280	Wear protective gloves, protective clothing, eye protection and face protection.
P234	Keep only in original packaging.
P273	Avoid release to the environment.
Response	
P301+P330+P331	IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
P303+P361+P353	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower].
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310	Immediately call a POISON CENTER/doctor/physician/first aider.
P363	Wash contaminated clothing before reuse.
P390	Absorb spillage to prevent material damage.
P391	Collect spillage.
P304+P340	IF INHALED: Remove person to fresh air and keep comfortable for breathing.
Storage	
P405	Store locked up.
Disposal	
P501	Dispose of contents/container to authorised hazardous or special waste collection point in accordance with any local regulation.

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