Testosterone

Enzyme immunoassay for the quantitative determination of total Testosterone in human serum or plasma

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Product Number: DNOV002 (96 Determinations)

1. INTRODUCTION

Testosterone (17 β -Hydroxy-4-androstene-3-one) is a steroid hormone from the androgen group mainly produced by the Leydig cells located in the testes and, minimally, by the ovaries and the adrenal cortex. It is also present in women who, compared to men, have a greater tendency to convert into estrogen this hormone. .

In postpubertal males, testosterone is secreted primarily by the testes with only a small amount derived from peripheral conversion of androstenedione. In humans is deputy to the development of the sexual organs (differentiation of the testis and the whole genital apparatus) and of secondary sexual characteristics, such as beard, hair distribution, the timbre of the voice and muscles. The testosterone, during puberty, is also involved on skeletal development, limiting the elongation of the long bones and avoiding, in this way, a disproportionate growth of the limbs. In adult humans, the levels of testosterone have a very important role as regards the sexuality, the musculoskeletal system, the vitality and good health (mainly understood as protection from metabolic diseases such as hypertension and diabetes mellitus); helps to ensure fertility, as it stimulates the maturation of sperm in the testes. It also influences the quality and quantity of sperm produced, since the seminal work on the streets and on the prostate, deputies to the production of sperm. Daily production of testosterone in men varies from 5 to 7 milligrams, but exceeded 30 years, tends to decrease annually by 1%.

In adult women over 50% of serum testosterone is derived from peripheral conversion of androstenedione secreted by the adrenal and ovary, with the remainder from direct secretion of testosterone by these glands. The majority of circulating testosterone is bound by SHBG and a smaller portion is bound by albumin. Only a small percentage (<1%) exists in circulation as unbound or free testosterone.

Testosterone effects can be classified as virilizing and anabolic effects, although the distinction is somewhat artificial, as many of the effects can be considered both. Anabolic effects include growth of muscle mass and strength, increased bone density and strength, and stimulation of linear growth and bone maturation. Virilizing effects include maturation of the sex organs, and after birth (usually at puberty) a deepening of the voice, growth of the beard and axillary hair (male secondary sex characteristics).

Testosterone levels decline gradually with age in men (andropause). The signs and symptoms are non-specific, and are generally associated with aging such as loss of muscle mass and bone density, decreased physical endurance, decreased memory ability and loss of libido.

In females of all ages, elevated testosterone levels can be associated with a variety of virilizing conditions, including adrenal tumors and polycystic ovarian disease.

2. INTENDED USE

Competitive immunoenzymatic colorimetric method for quantitative determination of total Testosterone in human serum or plasma.

3. PRINCIPLE OF THE ASSAY

Testosterone in the blood is bound to SHBG (60%) and in lower quantity to other proteins (for example albumin); the unbound Testosterone (< 1% of total Testosterone) is known as "free Testosterone". The chemical formulation of this assay allows to release completely the Testosterone from bound proteins; thus the NovaTec Testosterone kit allows to measure the concentration of total Testosterone (bound + free) in the sample. For the measurement of free Testosterone only, the NovaTec "Free Testosterone" kit is available.

Microtiter strip wells are precoated with anti-Testosterone antibodies (solid-phase). Testosterone in the sample competes with added horseradish peroxidase labelled Testosterone (enzyme-labelled antigen) for antibody binding. After incubation a bound/free separation is performed by solid-phase washing. The immune complex formed by enzyme-labelled antigen is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product. The intensity of this product is **inversely** proportional to the amount of Testosterone in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorption at 450 nm is read using an ELISA microwell plate reader.

4. MATERIALS

4.1. Reagents supplied

- **Anti-Testosterone IgG Coated Wells:** 12 breakapart 8-well snap-off strips coated with anti-Testosterone antibody; in resealable aluminium foil.
- **Stop Solution:** 1 bottle containing 15 ml sulphuric acid, 0.15 mol/l (avoid any skin contact).
- **Testosterone conjugate:** 1 bottle containing 12 ml of horseradish peroxidase labelled Testosterone.
- **TMB Substrate Solution:** 1 bottle containing 15 ml 3, 3', 5, 5'-tetramethylbenzidine (H₂O₂-TMB 0.26 g/l) (avoid any skin contact).
- **Wash solution 10x conc.:** 1 bottle containing 50 ml of a 10x concentrated solution of phosphate buffer 0.2 M, Proclin < 0.0015%
- **Testosterone control A:** 1 bottle containing 1 ml of a lot-specific control solution. The concentration is indicated on the label of the bottle.

- **Testosterone control B:** 1 bottle containing 1 ml of a lot-specific control solution. The concentration is indicated on the label of the bottle.
- **Testosterone Standards:** 5 bottles, 1 ml each
 - Standard 0: 0.0 ng/ml
 - Standard 1: 0.2 ng/ml
 - Standard 2: 1.0 ng/ml
 - Standard 3: 4.0 ng/ml
 - Standard 4: 16.0 ng/ml

4.2. Materials supplied

- 1 Strip holder
- 1 Cover foils
- 1 Test protocol
- 1 Distribution and identification plan

4.3. Materials and Equipment needed

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450 nm, 620-630 nm)
- Incubator 37°C
- Manual or automatic equipment for rinsing wells
- Pipettes
- Rotating mixer
- Distilled water
- Timer
-

5. STABILITY AND STORAGE

The reagents are stable up to the expiry date stated on the label when stored at 2...8 °C .

6. REAGENT PREPARATION

It is very important to bring all reagents, samples and standards to room temperature (22...28°C) before starting the test run for at least 30 minutes! At the end of the assay store the reagents immediately at 2-8° C; avoid long exposure to room temperature.

6.1. Coated snap-off Strips

The ready to use break apart snap-off strips are coated with anti-Testosterone antibodies. Store at 2...8 °C in the dark. *Immediately after removal of strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C in the dark,; once opened, the microplate is stable until expiry date of the kit.*

6.2. Testosterone-HRP Conjugate

Ready to use. Mix gently for 5 minutes with a vortex mixer.
After first use the conjugate is still stable for another 6 months if stored at 2...8 °C.

6.3. Testosterone Standards

The standards are ready to use and have the following concentration of Testosterone:

Standard 0:	0.0 ng/ml
Standard 1:	0.2 ng/ml
Standard 2:	1.0 ng/ml
Standard 3:	4.0 ng/ml
Standard 4:	16.0 ng/ml

The solutions have to be stored at 2...8 °C. *After first use the standards are still stable for another 6 months if stored at +2...+8° C.*

6.4. TMB Substrate Solution

The bottle contains 15 ml of a tetramethylbenzidine/hydrogen peroxide system. The reagent is ready to use and has to be stored at 2...8°C in the dark. *The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away. After first use the TMB substrate solution is still stable for another 6 months if stored at 2...8 °C.*

6.5. Stop Solution

The bottle contains 15 ml 0.15 M sulphuric acid solution (R 36/38, S 26). This ready to use solution has to be stored at 2...8°C. *After first use stable until expiry date.*

6.6. Wash Solution

Dilute the content of each vial of the "10X Conc. Wash Solution" with distilled water to a final volume of 500 ml prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C. In concentrated wash solution it is possible to observe the presence of crystals; in this case mix at room temperature until the complete dissolution of crystals; for greater accuracy, dilute the whole bottle of concentrated wash solution to 500 ml, taking care to transfer completely the crystals, then mix until crystals are completely dissolved.

6.7. Testosterone Controls

Ready to use. The bottles contain 1 ml of a lot-specific control solution. The concentration is indicated on the label.

7. SPECIMEN COLLECTION AND PREPARATION

Use human serum or plasma samples with this assay. If the assay on the same day as sample collection, the specimen should be aliquoted and stored deep-frozen (-20°C). If samples are stored frozen, mix thawed samples well before testing. *Avoid repeated freezing and thawing.*

Treatment of the patient with cortisone, natural or synthetic steroids can impair Testosterone determination.

8. ASSAY PROCEDURE

8.1. Test Preparation

Please read the test protocol carefully **before** performing the assay. Result reliability depends on strict adherence to the test protocol as described. Prior to commencing the assay, the distribution and identification plan for all specimens and controls should be carefully established on the result sheet supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than one plate is used, it is recommended to repeat the dose response curve. Please allocate at least:

1 well	(e.g. A1)	for the substrate blank
2 wells	(e.g. B1+C1)	for standard 0
2 wells	(e.g. D1+E1)	for standard 1
2 wells	(e.g. F1+G1)	for standard 2
2 wells	(e.g. H1+A2)	for standard 3
2 wells	(e.g. B2+C2)	for standard 4
2 wells	(e.g. D2+E2)	for control A
2 wells	(e.g. F2+G2)	for control B

It is recommended to determine standards, controls and patient samples in duplicate.

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each standard and each patient sample.

Adjust the incubator to 37° ± 1°C.

1. Dispense 25 µl standards, controls and samples into their respective wells. Add 100 µl conjugate to each well. Leave well A1 for substrate blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour at 37°C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well 3 times with 300 ml of diluted wash solution. Avoid overflows from the reaction wells. During each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: in case you use an automatic washer, it is advised to do 6 washing steps.

Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.

5. Dispense 100 µl TMB Substrate Solution into all wells.
6. **Incubate for exactly 15 min at room temperature (22...28°C) in the dark.**
7. Dispense 100µl Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution. Shake the microplate gently. *Any blue colour developed during the incubation turns into yellow.*
8. Measure the absorbance (E) of the specimen at 450 nm against a reference wavelength of 620-630 nm or against blank within 5 minutes.

8.2. Measurement

Adjust the ELISA Microwell Plate Reader **to zero** using the **substrate blank in well A1**.

If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the substrate blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard and patient sample in the distribution and identification plan.

Where applicable calculate the **mean absorbance values** of all duplicates.

9. RESULTS

9.1. Calculation of results

Calculate the mean absorbance for each point of the standard curve and each sample. Plot the mean value of absorbance of the standards against concentration. Draw the best-fit curve through the plotted points. (Four Parameter Logistic).

Interpolate the values of the samples on the standard curve to obtain the corresponding values of the concentrations expressed in ng/ml.

9.2. Reference values

The serum-Testosterone reference values are:

	n° sera	Median pg/ml	Range pg/ml
Male (age)			
< 12	16	< 0,10	< 0.10 - 1.01
12 - 18	16	5,02	0.56 - 8.63
19 - 55	16	3,54	2.12 - 6.01
> 55	16	1,51	0.11 - 7.25
Female (age)			
< 12	16	< 0,10	< 0.10 - 0.16
12 - 18	15	0,23	< 0.10 - 0.63
19 - 55	16	0,20	< 0.10 - 0.63
> 55	16	0,13	< 0.10 - 0.32

Please pay attention to the fact that the determination of a range of expected values for a “normal” population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

10. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of Testosterone for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the standard curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

11. SPECIFIC PERFORMANCE CHARACTERISTICS

11.1. Precision

Intra Assay Variation

Within run variation was determined by replicate determination (16x) of three different sera in one assay. The within assay variability is $\leq 7.0\%$.

Inter Assay Variation

Between run variation was determined by replicate measurements (9x) of three different sera in different lots. The between assay variability is $\leq 8.3\%$.

11.2. Specificity

The cross reaction of the antibody calculated at 50% according to Abraham is:

Testosterone	100%
DHT	2.03%
Androstenedione	0.01%
Androsterone	0.05%
DHEA-S	0%
Cortisol	0.01%
Cortisone	0%
17 α Estradiol	0.16%
Estrone	0%
Prednisone	0.01%

11.3. Sensitivity

The lowest detectable concentration of testosterone that can be distinguished from the zero standard is 0.10 ng/ml at the 95 % confidence limit.

11.4. Accuracy

The recovery of 0.4 - 0.8 - 4.0 - 14.0 ng/ml of Testosterone added to the sample gave an average value of 98.9% with reference to the original concentrations.

The dilution test performed on 3 samples diluted up to 4 times gave an average value of 99.7%.

11.5. Correlation with RIA

The new Testosterone ELISA kit was compared to the old Testosterone ELISA kit. Serum samples from 24 females and 28 males were analysed.

The linear regression curve was calculated:

$$Y = 1.13 \cdot X + 0.06$$

$$r^2 = 0.92$$

The new Testosterone ELISA kit was compared to another commercially available Testosterone assay. Serum samples of 24 females and 29 males were analysed according in both test systems.

The linear regression curve was calculated:

$$Y = 0.89 \cdot X + 0.24$$

$$r^2 = 0.93$$

12. PRECAUTIONS AND WARNINGS

- In compliance with article 1 paragraph 2b European directive 98/79/EC the use of the in vitro diagnostic medical devices is intended by the manufacturer to secure suitability, performances and safety of the product. Therefore the test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use by professional persons. Not for internal or external use in humans or animals
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive. Nevertheless, all materials should still be regarded and handled as potentially infectious.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and control vials for microbial contamination prior to further use.

- When using automatic equipment, the user has the responsibility to make sure that the kit has been appropriately tested. To improve the performance of the kit on ELISA automatic systems, it is recommended to increase the number of washes.
- Avoid the exposure of TMB substrate to direct sunlight, metal or oxidants. Do not freeze the solution.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- Maximum precision is required for dispensation of the reagents.
- This method allows the determination of Testosterone from 0.2 to 16.0 ng/ml
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Microbiologically contaminated samples should not be used in the assay. Highly lipemic or haemolysed specimens should also not be used.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background.
- Some reagents contain small amounts of Sodium Azide (NaN₃) or Proclin 300[®] as preservative. Avoid the contact with skin or mucosa.
- Sodium Azide may be toxic if ingested or absorbed through the skin or eyes; moreover it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, allow scroll through large amounts of water to prevent azide build-up.
- The NovaTec ELISA is only designed for qualified personnel who are familiar with good laboratory practice.

12.1. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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

Testosterone

Test Preparation

Prepare reagents and samples as described.
Establish the distribution and identification plan for all specimens and controls on the result sheet supplied in the kit.
Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate blank	Standard 0 - 4	Controls	Sample
Standard 0 - 4	-	25 µl	-	-
Controls	-	-	25 µl	-
Sample	-	-	-	25 µl
Conjugate	-	100 µl	100 µl	100 µl
Cover wells with foil supplied in the kit Incubate for 1 h at 37°C Wash each well three times with 300 µl diluted wash solution In case you use an automatic washer, it is advised to do 6 washing steps.				
TMB Substrate	100 µl	100 µl	100 µl	100 µl
Incubate for exactly 15 min at room temperature (22...28°C) in the dark				
Stop Solution	100 µl	100 µl	100 µl	100 µl
Shake the microplate gently Photometric measurement at 450 nm against a reference wavelength of 620-630 nm or against blank within 5 minutes.				

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NovaTec Immundiagnostica GmbH Technologie & Waldpark