## BIOCHEMISTRY CONTROL SERUM (Human) I, II

BioSystems

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#### INTENDED USE

The Biochemistry Control Serum (human) is intended for quality control in the laboratory only and is supplied with suggested intervals of acceptable values. The control can be used manually or in BioSystems analyzers. For *in vitro* professional use only in the clinical laboratory.

#### COMPOSITION

**Biochemistry Control Serum.** For 5 mL. Lyophilized human serum containing component concentrations suitable for the quality control of the clinical laboratories, and without preservatives which might interfere with the tests.

Components from human origin have been tested and found to be negative for the presence of antibodies anti-HIV and anti-HCV, as well as for HBs antigen. However, they should be handled cautiously as potentially infectious.

#### PREPARATION AND USE

1. Open a vial very carefully, avoiding any loss of the lyophilized material.
2. Pipette exactly 5.00 mL of distilled water into the vial. The component values depend on the accuracy of this reconstitution step.
3. Close the vial with the stopper and let stand for 20 minutes at room temperature. For alkaline phosphatase measurement, it is recommended to let the reconstituted material stand for 1 hour at room temperature before use.
4. Swirl gently, avoiding the formation of foam, to ensure complete dissolution of contents.
5. If the material is to be used for analysis of trace elements, avoid contact of the reconstituted material with the stopper to prevent a possible contamination.
6. The reconstituted control serum is to be treated like the patient serum.

#### STORAGE AND STABILITY

Store at 2-8°C.

The lyophilized Biochemistry Control Serum is stable until the expiration date given in the label.

The components of the reconstituted material are stable for at least 7 days at 2-8°C and 30 days at -20°C (when frozen once), excepting:

- AST is stable 8 hours at 2-8°C and 30 days at -20°C.
- Alkaline phosphatase is stable 5 hours at 2-8°C and 30 days at -20°C.
- NEFAs are stable 2 days at 2-8°C and 30 days at -20°C.

CK and bilirubin are sensitive to light. Store the vials protected from light.

Discard the vial if there are signs of microbial contamination or excessive turbidity in the reconstituted product.

#### WARNING AND PRECAUTIONS

Exercise the normal precautions required for handling all laboratory reagents. Safety data sheet available for professional user on request. Disposal of all waste material should be in accordance with local guidelines. Any serious incident that might occur in relation to the device shall be reported to BioSystems S.A.

#### SPECIFIED VALUES

The assigned concentration values for components are batch-dependent and their traceability is shown in the enclosed value sheets. Traceability of the results can be assured only if the BioSystems reagents and recommended measurement procedures are used.

The suggested intervals of acceptable values have been calculated from previous experience in interlaboratory variability and are given for orientation only. Each laboratory should establish its own precision parameters.

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8435287140BIOCHCoSeHHG

M18042i-09

The lateral lines mark the modifications in the current version.

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BioSystems S.A. Costa Brava, 30. 08030 Barcelona (Spain)