


Document ID: 60886352

**DRYSTAR 5302, DRYSTAR 5301, DRY IMAGER,
CRYSTAL IMAGER, DRYSTAR 5302 XL,
DRYSTAR 5301 XL, DRY 2.4, DRY 1.4,
DRYSTAR AXYS, CRYSTAL-M IMAGER,
DRYSTAR AXYS XL, DRY 2.4M**

Type 5366/100/101/110/300/400/410/500/510 and 5367/100/300/400/500

Includes the Maintenance Checklist, attached in chapter 9.5.3.



DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Meding Library. Any printed copy of this document is uncontrolled.



Manufacturer

Agfa NV

Published by

Agfa-Gevaert HealthCare GmbH
Tegernseer Landstraße 161
D - 81539 München
Germany

Copyright © 2021 Agfa NV

All rights reserved.

Technical modifications reserved.

Agfa and the Agfa rhombus are trademarks of Agfa-Gevaert N.V., Belgium, or its affiliates.

All other trademarks mentioned in this document are held by Agfa NV or the respective owners and are used in an editorial fashion with no intention of infringement.

Nothing contained in this legal notice nor in any text in this document shall be construed as granting by implication, estoppel or otherwise, any license or right to use any of the trademarks, service marks, trade names or logos appearing in this document without the express prior written consent of their respective owner.



WARNING:

Improper operation or service activities may cause damage or injuries.

- (1) Read the "Generic Safety Directions" prior to attempting any operation, repair or maintenance task on the equipment.
Refer to Document ID [11849633](#).
- (2) Strictly observe all safety directions within the "Generic Safety Directions" and on the product.



IMPORTANT:

The installation and service of the product(s) described herein is to be performed by qualified personnel who are employed by Agfa NV or one of its affiliates or who are otherwise authorized by Agfa NV or one of its affiliates to provide such services.

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.

LIST OF CONTENTS

0	ABOUT THIS MANUAL	13
0.1	Purpose of this document	13
0.2	Changes compared to previous revision	13
0.3	Referenced documents	13
0.4	Explanation of notes	14
0.5	Conventions	14
1	PRODUCT DESCRIPTION	15
1.1	Intended use	15
1.2	Type overview	15
1.3	Main printer characteristics	16
1.4	Printer components	17
1.4.1	Components Location	17
1.4.2	List of Components	21
1.4.3	Tray detection sensors in Printers with SN > 151.001/751.001	24
1.5	Overview of control boards	25
1.6	Overview of major hardware changes	26
1.7	Film run	28
1.8	Backside Printing Recognition	34
1.9	Left / Right Shift Print Algorithm	35
1.10	Mammo Printing Functionality (type 5367 only)	37
1.11	Printer calibrations overview	39
1.12	Printer Modules	40
1.13	Hardware	42
1.14	Modules and E-Labels	44
1.15	Printer online diagnosis software	46
1.16	Thermal model	48
2	SAFETY DIRECTIONS	49
2.1	Safety Notes for Installation Planning and Installation	49

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.

2.2	Safety Notes for preventive and corrective Maintenance Work.....	49
2.3	General Safety Note for all Service Activities	50
2.3.1	Safety switch	50
2.3.2	Electrical Check after Power Supply Replacements.....	51
3	PRE-INSTALLATION.....	53
3.1	Items to be organized before installation	53
3.1.1	Power cables.....	53
3.1.2	Tools	53
3.1.3	Documents	53
3.2	Scope of delivery	55
3.3	Accessories.....	56
3.4	Remote service	56
3.5	Required data for printer setup on site	57
3.6	Specifications	57
3.6.1	Printer dimensions	57
3.6.2	Printer Interfaces.....	58
3.6.3	Electrical connection	59
3.6.4	Resolution	60
3.6.5	Pixel Matrices.....	60
3.6.6	Interface Data.....	61
3.6.7	User access to the printer	62
3.6.8	Other technical data.....	62
4	INSTALLATION AND CONFIGURATION	63
4.1	Installation overview.....	64
4.1.1	Installation steps	64
4.2	Installation wizard	66
4.3	Printer configuration.....	66
4.3.1	Setting film calibration parameters	67
4.3.2	Setting the Default Densitometer	68
4.3.3	Setting the Film View	69

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.

4.3.4	Setting Beeps on empty Tray.....	70
4.3.5	Activating the Mammo License (type 5367 only)	71
4.3.6	Entering Site data	72
4.4	Film calibration	73
4.5	Printer - Modality connection	73
4.5.1	Modality set-up.....	73
4.6	Logical printers Set-up	74
4.7	Printing of diagnostic and / or test images	78
4.8	Image adjustments.....	78
4.8.1	For type 5367 only: Mammography specific Image Adjustments (optional)	79
4.8.2	Image Adjustments as described in the Connectivity Release Document	81
4.8.3	Image Adjustments at the Modality.....	81
4.8.4	Basic Image Adjustment at the Printer via Key-Operator Menu	83
4.9	Performing the Acceptance Test.....	87
4.10	Completing installation.....	87
4.10.1	Customer training.....	87
4.10.2	Backup of Printer Settings	88
4.10.3	Preparing the printer for remote service	88
5	CALIBRATION.....	89
5.1	Film Calibration	90
5.2	CDM (Internal Densitometer) Calibration.....	92
5.3	Print Head (TH) Calibration.....	94
5.4	Smart Profile Calibration	96
5.5	Compliance Film Calibration as part of Afssaps / EUREF Calibration Procedure.....	101
5.5.1	Checking the Densitometer Calibration	104
5.5.2	Checking the Densitometer Configuration in the Printer	104
5.5.3	Resetting the Compliance Film Calibration, Internal Densitometer Calibration and Thermal Print Head Calibration	105
5.5.4	Performing Internal Densitometer Calibration, Film Calibrations and Thermal Print Head Calibration	106
5.5.5	Performing a Density Reproduction Check.....	110

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.

5.5.6	Performing the Afssaps / EUREF Calibration Check.....	111
5.5.7	Performing a Compliance Film Calibration	112
6	ACCEPTANCE TEST	116
7	TROUBLESHOOTING	117
7.1	Film density, sharpness or contrast not OK.....	117
7.1.1	Dmin (minimum Density) too high.....	117
7.1.2	Dmax too high or too low	118
7.1.3	Density Reproduction not okay	119
7.1.4	Optimal perceivable Levels not okay	120
7.1.5	Contouring visible; 'Landscaping'.....	121
7.1.6	Contrast and/or Sharpness not okay	122
7.2	Line Artifacts in Transport Direction.....	123
7.2.1	Dust Stripes on the Film.....	123
7.2.2	Dark Line(s) < 1mm in Transport direction	124
7.2.3	White or bright Line(s) < 1mm in Transport direction	125
7.2.4	Bright vertical Lines over the whole Film Length	126
7.2.5	Narrow (between 1 and 5 mm) to wide (> 5 mm) Lines in Transport direction.....	127
7.2.6	Persistent dark line in transport direction.....	128
7.2.7	Persistent blank line in transport direction.	129
7.2.8	Intermittent blank line in transport direction.	130
7.2.9	Wide Bands with lower Density at the left and/or right Hand Side of the Film	131
7.2.10	Regular Line Structure vertical to Transport direction: Banding	132
7.3	Point Artifacts	133
7.3.1	Point Artifacts - Pinholes.....	133
7.3.2	Point Artifacts - Bubbles.....	134
7.3.3	Dust Particles on the Print Drum.....	135
7.3.4	Dust Particles in the Film	136
7.4	Surface Artifacts.....	137
7.4.1	Blemish: The Artifact looks like Dirt spread out on the Surface of the Film.....	137
7.4.2	Scratches: One or more bright Line over the whole Length of the Film	138

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.

7.4.3	Film Curl if Film is hanging at the Lightbox.....	139
7.4.4	Ghost image visible on films of printers which mainly print small film formats	139
7.5	Display Problems	140
7.5.1	Blank LCD Display	140
7.5.2	Keys at LCD Display do not work	140
7.6	Error and Warning Messages	141
7.6.1	Display Message: MODXX: Please replace Module.....	141
7.6.2	Display message: New Error EXXXXX	142
7.6.3	Display Message: New Error E312000	143
7.6.4	Display Message: New Maintenance.....	143
7.6.5	Error message: "Film Jam: Open upper input tray and clear film path". Error HC314003 in the infocounter.	145
7.6.6	Message "Error System Failure"	146
7.6.7	Keypad shows "Please wait - Selftest Proceeding" for more than 5 Minutes.....	148
7.7	Image Transfer Problems	149
7.7.1	Images are not shown in the Printer Queue	149
7.8	Film Jam.....	150
7.8.1	Message: Remove Film Jam	150
7.9	Printer error codes	151
7.9.1	Error Codes of printers with serial number < 100.000	151
7.9.2	Error Codes of printers with serial number > 100.000	190
8	REPAIR	227
8.1	Required tools.....	227
8.2	Connection of the Service PC to the Printer	228
8.2.1	Connection of the Service PC to the Printer	229
8.2.2	Adaptation of the Service PC IP Address	229
8.3	Software tools reference	231
8.3.1	Key-Operator Menu on the local Keypad.....	232
8.3.2	Infocounter	235
8.3.3	Test images.....	236
8.3.4	Show Error Codes.....	244

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.

8.3.5	Printer Browser Pages	245
8.3.6	Key-Operator Tools.....	247
8.3.7	Service Engineer Tools	252
8.3.8	Printer Logging	277
8.3.9	Printer Bootup-Message	278
8.3.10	ADVP	280
8.4	Fuses, LEDs and jumpers.....	286
8.4.1	Fuses	286
8.4.2	LEDs	286
8.4.3	Hard disk jumper settings	293
8.5	Replacements	296
8.6	Adjustments	296
8.6.1	Horizontal Film Alignment Procedure	296
8.7	Checks	299
8.7.1	Density reproduction check.....	299
8.7.2	Check for dust stripes / TPH Cleaning Procedure	302
8.7.3	Resistor Check.....	305
9	MAINTENANCE.....	307
9.1	Introduction: Maintenance frequency, required tools and time.....	307
9.2	Diagnostics.....	309
9.2.1	Questioning of the Customer	309
9.2.2	Infocounter Analysis.....	309
9.2.3	Image Quality Check.....	313
9.3	Maintenance Points - Power off	314
9.3.1	Check of general Condition.....	314
9.3.2	Check of Condition of Belts and Gears.....	314
9.3.3	Visual Check of the Cables.....	315
9.3.4	Cleaning of the Suction Cups	315
9.3.5	Cleaning of Electronics Compartment	316
9.3.6	Cleaning of Power Supply Air Inlet	316

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.

9.3.7	Check and Cleaning of Drum.....	317
9.3.8	Cleaning of Thermal Print Head.....	318
9.3.9	Electrical Check	319
9.4	Maintenance Points - Power on	321
9.4.1	Smart Profile Calibration	321
9.4.2	Check of Safety Switch	323
9.4.3	SW Patch Installation.....	323
9.4.4	Film Calibration	323
9.4.5	Density Reproduction Check	324
9.4.6	Test Image Evaluation	326
9.4.7	Host System Image Quality Check	327
9.5	Completion of the maintenance	327
9.5.1	Backup printer settings	327
9.5.2	Reset of Maintenance Action counters and relative Counters in the Infocounter.....	328
9.5.3	Customer Conversation	329
10	RELEASE INFORMATION	330
10.1	Software release 1.6.0_C1 (DRYSTAR AXYS only)	330
10.1.1	Supported Features	330
10.1.2	Limitations.....	330
10.1.3	Known Software Issues	330
10.2	Software update 1.6.1_C1 (DRYSTAR AXYS only)	331
10.2.1	Solved Software Issues	331
10.3	Software update 1.6.1_C1 (DRYSTAR AXYS only)	332
10.3.1	Solved Software Issues	332
10.3.2	New Features.....	333
10.4	Software update 1.6.3_C1 (DRYSTAR AXYS only)	336
10.4.1	Solved Issue.....	336
10.5	Software Release 1.8.0_C1 (DRYSTAR 5302 only).....	337
10.5.1	Supported Features	337
10.5.2	Limitations.....	337

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.

10.5.3	Known Issues.....	337
10.6	Software Release 1.8.0_C2 (DRYSTAR 5302 only).....	338
10.6.1	Solved Software Issue	338
10.6.2	New Feature.....	338
10.7	Software Release 1.8.1_C1 (DRYSTAR 5302 only).....	339
10.7.1	Solved Software Issues	339
10.8	Software Release 1.9.0_C2 (DRYSTAR 5302 only).....	340
10.8.1	New Features.....	340
10.9	Software Release 1.9.0_C3 (DRYSTAR 5302 only).....	340
10.9.1	Solved Software Issue	340
10.10	Software Release 1.9.0_C4 (DRYSTAR 5302 only).....	341
10.10.1	Solved Software Issues	341
10.11	Software Release 1.9.1_C1 (DRYSTAR 5302 only).....	342
10.11.1	Solved Software Issues	342
10.12	Software Release 1.9.2_C1 (DRYSTAR 5302 only).....	344
10.12.1	Solved Software Issues	344
10.13	New Feature.....	345
10.14	Software update 2.0.0_C2	345
10.14.1	New Features of Software 2.0.0_C2.....	345
10.14.2	Solved Issue in Software 2.0.0_C2.....	346
10.14.3	Known Software Issue	347
10.15	Software update 2.0.0_C3	347
10.15.1	New Feature of Software 2.0.0_C3	347
10.16	Software update 2.1.0_C2 (DRYSTAR AXYS only)	347
10.16.1	New Feature of Software 2.1.0_C2	348
10.17	Software update 2.1.0_C3	348
10.17.1	Solved Issues in Software 2.1.0_C3.....	348
10.18	Software update 2.2.0_C1	349
10.18.1	Solved Issues in Software 2.2.0_C1.....	349
10.19	Software update 2.3.0_C2	350

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.

10.19.1	Solved Issues and Improvements in Software 2.3.0_C2.....	350
10.20	Software update 3.0.0_C1	351
10.20.1	New Features in Software 3.0.0_C1.....	351
10.20.2	Solved Issues and Improvements in Software 3.0.0_C1.....	351
10.21	Software update 3.1.0_C1	352
10.21.1	Solved Issues and Improvements in Software 3.1.0_C1.....	352
10.22	Software Patch 3.1.1_C1	353
10.23	Software Release 4.0.0_C1	354
10.24	Software Release 5.0.0_C1	355
10.25	Software Release 5.1.0_C1	356
10.26	Software Release 6.0.0_C1 (Printer SN > 100.000 only)	357
10.27	Software Release 6.1.0_C1 (Printer SN > 100.000 only)	357
10.28	Software Release 6.1.1_C1 (Printer SN > 100.000 only)	358
10.29	Software Release 6.1.2_C1 (Printer SN > 100.000 only)	359
10.30	Software Release 6.1.3_C1 (Printer SN > 100.000 only)	360
10.31	Software Release 6.1.4_C1 (Printer SN > 100.000 only)	361
10.32	Software Release 6.2.0_C1 (Printer SN > 100.000 only)	362
10.33	Software Release R6.2.1_C1 (Printer SN > 100.000 only)	363
10.34	Software Release R6.2.2_C1 (Printer SN > 100.000 only)	364
11	UPGRADE PROCEDURE	365
12	SPARE PARTS.....	365
13	WIRING DIAGRAM.....	366
13.1	DRYSTAR 5302 Type 5366/100/300/400/500 SN < 100.000 with SCSI hard disk.....	367
13.2	DRYSTAR 5302 Type 5366/100/300/400/500 SN < 100.000 with SATA hard disk.....	375
13.3	DRYSTAR 5302 Type 5366/100/300/400/500 SN < 100.000 with Micro SATA flash module	383
13.4	DRYSTAR Axys Type 5367/100/400/500 SN < 100.000 with SATA hard disk.....	391
13.5	DRYSTAR Axys Type 5367/100/400/500 SN < 100.000 with Micro SATA flash module	399
13.6	DRYSTAR 5302 Type 5366/100/300/400/400 SN > 100.000	407
13.7	DRYSTAR 5301 Type 5366/110/410/510 SN > 100.000	413
13.8	DRYSTAR Axys Type 5367/100/300/400/500 SN > 100.000	417

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.

13.9	DRYSTAR 5302 Type 5366/100/300/400/400 SN > 151.001/751.001	423
13.10	DRYSTAR 5301 Type 5366/110/410/510 SN > 151.001/751.001	427
13.11	DRYSTAR Axys Type 5367/100/300/400/500 SN > 151.001/751.001	430

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.

0 About this manual

0.1 Purpose of this document

This document contains all information that service engineers need for installation as well as for corrective and preventive maintenance of the DRYSTAR 5302, DRYSTAR 5301, DRY IMAGER, CRYSTAL IMAGER, DRYSTAR 5302 XL, DRYSTAR 5301 XL, DRY 2.4, DRY 1.4, DRYSTAR AXYS, CRYSTAL-M IMAGER, DRYSTAR AXYS XL, DRY 2.4M.

It does **not** contain:

- Specific service information for other target readers, for example Clinical Application Specialists or Service Managers
- Service Bulletins

0.2 Changes compared to previous revision

Revision 5:

- Updated chapter 1.4
- Added chapter 1.4.3: Tray detection sensors in printers with SN 151.001/751.001
- Adapted chapter 1.6: Overview of major hardware changes
- Added Release Information in chapters 10.32, 10.33, 10.34
- Added chapters 13.9, 13.10, 13.11: Wiring diagrams of printers with SN $\geq$ 151.001/751.001

0.3 Referenced documents

Document	Reference
List of Service Documents,	Document ID 62000370
Generic Safety Directions	Document ID 11849633
User Manuals	Section 3.1.3
Additional documents are referenced in the corresponding sections.	

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.

0.4 Explanation of notes

Safety-relevant Notes		
Icon	Signal Word	Situation
	CAUTION:	Hazardous situation which, if not avoided, can lead to a minor injury to a user, engineer, patient or any other person.
	WARNING:	Hazardous situation which, if not avoided, can lead to a potential serious injury to a user, engineer, patient or any other person.
	DANGER:	Direct, immediate danger: If not avoided, it can lead to a serious injury to a user, engineer, patient or any other person.
	-	Instruction to avoid damage to equipment and/or environmental pollution.
	-	Prohibition to avoid damage to equipment and/or environmental pollution.

None-Safety-relevant Notes		
Icon	Name	Type of Information
	IMPORTANT:	Highlights very important actions which have to be carried out to prevent malfunction.
	NOTE:	- Indicates advice to facilitate the following step or action without having a direct influence on the step or action. - Highlights unusual points. - Indicates background information. - Can be used to explain or highlight displays of the graphical user interface.

0.5 Conventions

Highlighting of tasks		
Task number	Task Description	Remark
(1)	Connect the cable.	Examples for working steps to be performed in the listed sequence.
(2)	Switch the machine on.	

Highlighting of buttons, functions and names within a task		
(1)	Press F9 or double-click the Refresh button.	The bold typeface is used for menu topics, keyboard keys, icons, device buttons, commands etc.
(2)	If the Install client or server? prompt appears, type server .	The Courier New bold typeface is used for system messages and prompts.

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.

1 Product description

1.1 Intended use

The intended use statement is listed in the User Manual of the printer.

1.2 Type overview

The printer is available in several variants.

Main difference is the print head resolution (320 or 508 dpi (dots per inch)) and number of film trays (1 or 2).

Printer name	Type	DPI	Number of film trays
DRYSTAR 5302	Type 5366/100	320	2
DRY IMAGER	Type 5366/101	320	2
DRYSTAR 5301	Type 5366/110	320	1
CRYSTAL IMAGER	Type 5366/300	320	2
DRYSTAR 5302 XL	Type 5366/400	320	2
DRYSTAR 5301 XL	Type 5366/410	320	1
DRY 2.4	Type 5366/500	320	2
DRY 1.4	Type 5366/510	320	1
DRYSTAR AXYS	Type 5367/100	508	2
CRYSTAL-M IMAGER	Type 5367/300	508	2
DRYSTAR AXYS XL	Type 5367/400	508	2
DRY 2.4M	Type 5367/500	508	2

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.

1.3 Main printer characteristics

Purpose of the printer type 5366

- Medium resolution (320 dpi) DICOM network dry printer for 14x17", 14x14", 14x11", 10x12", 8x10" DT2 blue base or clear base films

Purpose of the printer type 5367

- High resolution (508 dpi) DICOM network dry printer for:
 - 14x17", 14x14", 14x11", 10x12", 8x10" DT2 blue base or clear base films
 - 14x11", 10x12", 8x10" DT2 Mammo blue base films

Printing principle

- Direct thermal on thermo-sensitive, dry films
- 1 or 2 input trays and 1 output tray

Main application printer type 5365

- Well suited to all medical applications, especially CR and DR applications. Cannot be used for mammography.

Main application printer type 5367

- Well suited to all medical applications, especially Mammo, CR and DR applications.

Connections

- 1 Gb/s Network connection (Default: Automatic detection of 10Mb/s and 100Mb/s network is enabled, automatic detection of 1Gb/s network is disabled)

Service Concept

- Customer installable
- No service PC needed for basic installation, maintenance and service, but is convenient for entering detailed host profile settings.
- Strict modular design: Production assemblies = service assemblies
- All switches, sensors and motors are monitored via current measurement on main engine board (Devil / Devius and Devnix PCB*)
- If online diagnosis does not detect the exact reason for the error, an offline check of the modules can be performed via browser pages

*see overview of boards

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.