




USER MANUAL [PRODUCT DESCRIPTION SHEET]

Issue No. : 1 Date : 01-01-2017

Rev. No. : 1 Date : 01-09-2017

Product Name	:	Angeltouch Blood Administration Set [Vented]		
Product Code	:	TT-11		
CE Mark with No.	:			
Manufacturer's Name	:	Angioplast Private Limited		
Manufacturer's Address	:	Plot No. 4803, Phase IV, Vatva GIDC, Ahmedabad 382 445. [India]		
EC REP Address	:	mdi Europa GmbH, Langenhagener Str. 71, 30855, Langenhagen, Germany		
Intended Use	:	The sterile Blood Administration Set is a medical device used for transfusion of whole blood or blood derivatives from blood bag to the patients.		
Product Description	:	The sterile Set is made up different components made up from ABS, DEHP Free PVC, PP, HDPE etc. It is made up of Latex Tube injection site and the needle is made of Stainless steel tube.		
Product Specification	:	Spike Guard + Vented Spike + Drip Chamber + Fluid Filter + Tubing (155cm) + Flow Regulator + Tube Latex + Luer Tip + Needle with Protector (18G x 1.5").		
Material List& Make	:	Sr. No.	Name of Components	No.
		1	Spike Guard	1
		2	Spike (Closure Piercing Device)	1
		3	Housing for air-vent filter	1
		4	Blood Chamber	1
		5	Blood Filter	1
		6	Tubing (155cm)	1
		7	Flow Regulator (A) Clamp (B) Duddy	1
		8	Injection Site (Tube)	1
		9	Male Fitting Adaptor [Luer Tip]	1
		10	Needle with Protector (18G x 1.5")	1
				Material Used
				Polypropylene (PP)
				Acrylonitrile Butadiene Styrene (ABS)
				DEHP Free Poly-vinyl Chloride (PVC)
				DEHP Free Poly-vinyl Chloride (PVC)
				High Density Polyethylene (HDPE)
				DEHP Free Poly-vinyl Chloride (PVC)
				(A) Polypropylene (B)High Density Polyethylene (HDPE)
				Latex
				Poly-propylene (PP)
				Cannula : Stainless Steel + Poly-propylene (PP)



Prepared by		Approved by
	Page 1 of 2	
S. A Christian		M. M Patel

USER MANUAL [PRODUCT DESCRIPTION SHEET]

Issue No. : 1 Date : 01-01-2017

Rev. No. : 1 Date : 01-09-2017

Instruction For Use	:	<ol style="list-style-type: none">1. Disinfect inlet of container.2. Close the air-vent and roller clamp respectively.3. Insert the administration set holding the container in an upright position.4. Suspend container. Half fill drip chamber by squeezing and releasing. Do not fill completely.5. Open the snap vent and roller clamp to expel air from tubing. Close clamp.6. Perform vein-puncture and regulate flow by opening clamp.7. Designed to deliver approximately 20 drops/ml. <p>WARNING:</p> <ol style="list-style-type: none">1. Do not re-use. <p>CAUTION:</p> <p>* Do not use if package is open or damaged. * Do not re-sterilize. * Discard after single use. * Gravity Feed Only.</p>
Product Use	:	For Single Use Only
Is product Sterilized	:	Yes
Sterilization Method	:	Ethylene Oxide Gas Sterilization
Storage Conditions	:	Do not store in direct sunlight, at extreme temperature or in high humidity. Keep away from getting wet. Store in cool & dry place.
Packing Detail	:	
Unit Pack [Pouch]	:	HMHDPE, 215 mm x 110 mm x 225G
Gross Pack [Size & Qty.]	:	HMHDPE Pouch, 260 mm [10.25"] x 520 mm [16.5"] x 150G - 25 Nos.
Shipper Pack [Size & Qty.]	:	Corrugated 5 ply, 595 mm x 375 mm x 325 mm - 425 Nos.

Prepared by		Approved by
	Page 2 of 2	
S. A Christian		M. M Patel