Ref Electrode

Reference Electrode

(For use with direct ISE module)

INSTRUCTIONS FOR USE

mindray

Ref Electrode Reference Electrode (For use with direct ISE module)

Order Information

Cat. No.	Package size
115-084088-00	1 electrode

Intended Purpose

Reference electrode is used for in vitro test for forming electromotive force with the Sodium/Potassium/Chloride electrodes and providing the potential for measuring the concentration of $\text{Na}^+/\text{K}^+/\text{Cl}^-$ in human serum, plasma and urine samples on Mindray BS series chemistry analyzers.

Summary

It forms electromotive force with the $\text{Na}^+/\text{K}^+/\text{Cl}^-$ measuring electrodes, providing the potential for measuring the concentration of $\text{Na}^+/\text{K}^+/\text{Cl}^-$ in human serum, plasma and urine samples.

Assay Principle

The sample concentration of $\text{Na}^+/\text{K}^+/\text{Cl}^-$ in serum or urine on Mindray chemistry analyzer is measured by ion-selective electrode method (direct method). In the measurement of $\text{Na}^+/\text{K}^+/\text{Cl}^-$, when the sample flows through the ion measuring electrode and the reference electrode, the potential difference is formed with the

measuring potential to assist in the measurement. ¹

Key Components

The main components of the reference electrode are the sensitive membrane in contact with the sample flow path, the reference solution in the carrier, and the internal reference electrode Ag/AgCl.

Storage and stability

Storage and stability: Refer to expiration date indicated on the label when stored unopened for 12 months at 2-32°C and without condensation.

On-board stability: After the electrodes are installed, 10,000 samples can be tested for nine months. The expiration date is determined by the date when the test period or tested sample volume reaches the limit first.

For replacement refer to instructions on the Mindray Chemistry Analyzer Operator's Manual.

Preparation

ISE Reagent Pack and K⁺/Na⁺/Cl⁻ electrode are ready to use.

Please perform scheduled maintenance and standard operation including calibration and analysis to assure the performance of measurement system.

Materials required but not provided

Cat. No.	Product Description	Package size
115-084089-00	Potassium Electrode	1 electrode
115-084090-00	Sodium Electrode	1 electrode

Cat. No.	Product Description	Package size
115-084091-00	Chloride Electrode	1 electrode
105-025539-00	ISE Reagent Pack	Buffer A: 460 mL Buffer B: 360 mL
105-025540-00	Urine Diluent	1x50 mL
105-025543-00	ISE Cleaning Solution	1x50 mL
105-041321-A0	PROBE CLEANSER	1x50 mL
105-002225-A0	PROBE CLEANSER	1x50 mL
105-025028-A0	DETERGENT C	6x20 mL

For the PROBE CLEANSER and DETERGENT C, only the names and specifications are different. The formula and use are the same.

Calibrator and Control: Please check the section of reagent instruction of Calibration and Quality Control.

Please refer to the manual for information about storage and stability.

For information on materials required for operation of the instrument, and maintenance procedures refer to the Mindray Chemistry Analyzer Operator's Manual.

Assay procedure

Parameters Item		BS-600M chemistry analyzer
Sample types	Serum, plasma	Urine: Dilution first. Mix sample with urine Diluent with a proportion of 1:9.

Sample	Serum, plasma	70 μ L
	Urine	140 μ L

For Mindray BS series chemistry analyzers, the Parameters Item are the same. For information on viewing and editing assay parameters or for a detailed description of system procedures, refer to the Mindray Chemistry Analyzer Operator's Manual.

Calibration

1. For detailed calibration procedures, please refer to Mindray Chemistry Analyzer Operator's Manual.
2. Calibration frequency: 8h by default.

Recalibration may be necessary when the following occur:

- As ISE reagent lot changed;
- As the calibration factors are expired;
- As $\text{Na}^+/\text{K}^+/\text{Cl}^-$ /reference electrodes are replaced;
- As important components are replaced, such as sample tube, pump tube and tubing on the ISE module;
- As QC is out of control;
- After the ISE module and electrodes are maintained, including electrode cleaning and tube cleaning.

Each laboratory can set its own calibration procedures according to specific circumstances.

Quality control

Two levels of control material are recommended to analyze each batch of samples. It is recommended to use the Mindray ClinChem

Multi Control (Cat. No. 105-009119-00 and 105-009120-00) for internal serum/plasma quality control and MR Urine Quality Control (105-001691-00, Low: 1x100mL, High: 1x100mL) for internal urine quality control.

Controls provided by different manufacturers may be different. If a third-party QC is used, it is recommended to accumulate the target value and re-accumulate the target value when the lot of reagent and control is changed. The measured control value should be within the specified range. If the result is beyond the specified range, run the selected control at least once. If the result cannot be changed, it is recommended to replace the control or run it after recalibration. If necessary, take corresponding measures or contact the manufacturer.

Warnings and precautions

1. For in vitro diagnostic use only. For laboratory professional use. The installation and maintenance of electrodes must be operated by skilled/trained clinical professionals.
2. Please take the necessary precautions for handling all laboratory reagents.
3. Reliability of assay results cannot be guaranteed if the instructions in this package insert is not followed.
4. Please confirm the integrity of the package before use. Do not use the kits with damaged packages. The results can't be assured when stored at inappropriate condition.
5. Do not freeze the electrodes. Do not use electrodes and reagent

that have expired. For the expiration date, see the package label. For electrode replacement procedures, refer to the Mindray Chemistry Analyzer Operator's Manual. If the electrodes are stored refrigerated (2-8°C), they should be kept for a period of time and restored to room temperature before use.

6. After a large number of tests are finished, the proteins and lipids in the sample will adhere to the electrode surface, which may affect the measurement performance. It is recommended to clean the electrodes after finishing the ISE tests every day or maintain the electrodes before powering off the analyzer. For the specific cleaning steps, please refer to the Mindray Chemistry Analyzer Operator's Manual.
7. Biohazards will exist as electrodes may touch the patient's sample during the test, and the residual serum sample may contain virus and serum. Wear gloves when using the electrodes to avoid infection. If the sample touches the skin, wash the body area immediately with soap and water.
8. Before installing the electrodes, check visually whether the black sealing rubber ring of each electrode falls off, or it may cause poor sealing of the electrode tubes and affect the reliability of the results.
9. Before installing the electrodes, check visually if there are bubbles in the reference solution. Shake the electrodes gently to float up the bubbles when they exist.
10. Disposal of all samples and waste shall be in accordance with

national and local laws and regulations.

11. Safety data sheet is available for professional user on request.
12. Instability or deterioration should be suspected if there are visible signs of leakage, precipitates or microbial growth, or if calibration/controls do not meet the insert and/or the Mindray System criteria.
13. Install the electrodes according to the labels on the electrode compartment. Incorrect installation may result in calibration failure
14. All identified risks have been reduced as far as possible by generally acknowledged state of art, and the overall residual risk is acceptable.
15. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

References

1. Eisenman G. Glass Electrodes for Hydrogen and Other Cations, Principles and Practice. New York: Marcel Dekker Inc.; 1967:2.

Graphical symbols

In vitro
diagnostic
medical
device

Batch code

Unique
device
identifier

Authorized
representative in the
European
Community

Temperature limit

Consult
instructions for
use

Manufacturer

European
Conformity

Use-by date

© 2022-2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-81888998; **Fax:** +86-755-26582680

EC-Representative: Shanghai International Holding Corp.
GmbH(Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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