M708348B234 Rev D 2019-03

English

KYPHON® HV-R® (High Viscosity, Radiopaque) Bone Cement Instructions for Use

Carefully read all instructions prior to use.

INDICATIONS FOR USE

KYPHON® HV-R® Bone Cement is indicated for use as a fixation material in vertebral compression fractures undergoing Kyphoplasty treatment.

STERILE

KYPHON® HV-R® Bone Cement is supplied sterile. The powder and package are sterilized with gamma radiation. The liquid is sterilized using filtration and is contained in a glass vial. The outside of the glass vial is sterilized with ethylene oxide gas. This device is intended for single use only. Do not re-sterilize. Do not use if package is open or damaged.

DEVICE NAME

KYPHON® HV-R® Bone Cement

DESCRIPTION

KYPHON® HV-R® Bone Cement is a polymethylmethacrylate (PMMA) that contains approximately 30% barium sulfate. It is designed for delivery in a highly viscous state.

**Nominal Composition of KYPHON® HV-R® Bone Cement:
(Actual weight percentages of individual components will vary within accepted ranges.)**

POWDER (20 g of sterile powder in a packet)	LIQUID (9.0 g of sterile liquid in a vial)
Methylmethacrylate-styrene-copolymer 68.0% w/w	Methylmethacrylate (monomer) 99.1% w/w
Barium sulfate 30.0% w/w	N, N-dimethyl-p-toluidine 0.9% w/w
Benzoyl peroxide 2.0% w/w	Hydroquinone 75 ppm

CONTRAINDICATIONS

PMMA bone cement is contraindicated in the presence of active or incompletely treated infection at the site where the bone cement is to be applied.

WARNINGS

- The sterile KYPHON® HV-R® Bone Cement is designed and intended for one time use only. DO NOT re-sterilize and/or re-use. Sterility is assured only if the unit container is not damaged.
- Maintain aseptic surgical technique to prevent possible infection including treatment site infection.
- Monitor patients carefully for any change in blood pressure during and immediately following the application of bone cement.
- Adverse patient reactions affecting the cardiovascular system, including Bone Cement Implantation Syndrome (BCIS), have been associated with the use of bone cements. Hypotensive reactions have occurred between 10 and 165 seconds following application of bone cement; they have lasted from 30 seconds to 5 or more minutes. Some have progressed to cardiac arrest. Patients should be monitored carefully for any change in blood pressure during and immediately following the application of bone cement, especially those potentially at increased risk for peri-operative death, including elderly patients, patients with underlying cardiac or pulmonary compromise, and patients being treated for multiple vertebral body fractures in one procedure.
- Caution should be exercised during the mixing of the two components to prevent excessive exposure to the concentrated vapors of the monomer, which may produce irritation of the respiratory tract, eyes, and possibly the liver.
- Polymerization of the bone cement is an exothermic reaction, which occurs while the cement is hardening *in situ*. The released heat may damage bone or other tissues surrounding the implant.
- The liquid component is a powerful lipid solvent. It should not be allowed to come into contact with rubber or latex gloves. Should contact occur, the gloves may dissolve and tissue damage may occur. Wearing a second pair of gloves and adherence to the mixing instructions may diminish the possibility of hypersensitivity reactions. The mixed bone cement should not make contact with the gloved hand until the bone cement is being tested for the consistency of dough.
- Personnel wearing permeable contact lenses should NOT be near or involved in mixing the bone cement; soft contact lenses are permeable.
- The handling characteristics of bone cements are affected by operating room conditions, including the room temperature, temperature of the cement components prior to mixing, humidity, the geometry of the mixing apparatus, time spent mixing, and the geometry of the delivery device. Any change in one or more of these conditions can alter the handling characteristics of the bone cement, including the following:
 - Handling period - the time it takes for the bone cement to reach the doughy state (the cement has reached the doughy state when it no longer sticks to surgical gloves.)
 - Working period - the time the bone cement remains in the doughy state and can be delivered.
 - Hardening period - the time it takes for the bone cement to harden or until it can no longer be delivered.
- The user must be aware of these factors and adjust technique to account for variability in operating room conditions.
- Avoid over-pressurization of the bone cement and do not insert the bone cement into the cavity of the vertebral body until the bone cement has reached the doughy state. The cement has reached the doughy state when it no longer sticks to surgical gloves. Maintain patient positioning until the end of the polymerization process. Failure to do so may lead to extravasation of the bone cement beyond the site of its intended application and may or could damage the surrounding tissues.
- The safety of the bone cement in pregnant women or in children has not been established. Bone cement may adversely affect bone growth and fetal health.
- Difficulty in swallowing and blistering of the throat, alleged to be an allergic reaction to the use of bone cement, have also been reported.
- DO NOT mix more than one vial of liquid and one packet of powder together at any one time. Never modify the ratios between the liquid and solid components. Doing so can affect bone cement properties, including handling characteristics.
- Long-term follow-up is advised for all patients on a regularly scheduled basis. Inadequate filling or unanticipated postoperative events may affect the stability of the bone cement and/or the bone cement-bone interface. A fibrous tissue layer may develop between the cement and the bone, and loosening of the bone cement may occur, leading to failure.
- KYPHON® HV-R® Bone Cement is provided in finished form with all the necessary components for use. The addition of radiopacifier, e.g., barium; antibiotics; or other drugs or materials to KYPHON® HV-R® Bone Cement is not recommended. Never add other substances or foreign bodies to the acrylic resin. The safety and effectiveness of adding such drugs or materials has not been evaluated and may cause patient harm.

- Modifying the polymerization time by either warming or cooling the bone cement and/or associated delivery devices has not been tested and could affect bone cement properties, including handling characteristics.

PRECAUTIONS

- Store product below 25°C. Keep the product at a temperature of 23 ± 1°C for a period of 24 hours prior to use.
- Never deliver bone cement into the vertebral body without the use of high quality fluoroscopic guidance capable of visualizing movement of bone cement.
- Only physicians thoroughly trained in the surgical use of bone cement and balloon kyphoplasty should use KYPHON® HV-R® Bone Cement. The operator should have specific training and be familiar with the properties, handling characteristics, and application of the bone cement and adhere to the instructions for use.
- Do not use after the expiration date printed on the package. The device may not be safe or effective beyond its expiration date.
- Do not use if package is opened or damaged because product integrity including sterility may be compromised.
- Do not use damaged products. Prior to use, inspect the packaging and product to verify that no damage has occurred.
- The device is single use and for single patient: DO NOT RE-STERILIZE AND/OR RE-USE.
- Re-sterilization should not be performed because there is a possible risk of infection to the patient. Re-sterilization could modify the device's morphology and the mechanical features, leading to a malfunction of the device with serious risks for the patient's health.
- Never modify the ratios between the liquid and the solid components.
- Do not divide the product into two or more parts because the use of these parts in different times would be a re-use and may lead to an error in the ratio of the powder and liquid components. Moreover this could cause a sterility loss.
- Follow the mixing and handling instructions to avoid contact dermatitis. Strict adherence to the instructions for mixing the powder and liquid components may reduce the incidence of this complication.
- Adequately ventilate the operating room to eliminate as much monomer vapor as possible. The liquid monomer is highly volatile and flammable. Ignition of monomer fumes caused by use of electrocautery devices in surgical sites near freshly implanted bone cements has been reported.
- All leftover materials must be considered waste surgical material and consequently should be eliminated at the end of the surgical technique.**
- Dispose of the polymer component in an authorized waste facility. The liquid component can be evaporated under a well ventilated hood or absorbed by an inert material and transferred in a suitable container for disposal.

ADVERSE EVENTS

Serious adverse events, some with fatal outcome, associated with the use of acrylic bone cements include:

- cardiac arrest
- cerebrovascular accident
- myocardial infarction
- pulmonary embolism

Other reported adverse events relevant to the anatomy being treated with acrylic bone cements include:

- deep or superficial wound infection
- fistula
- hematoma
- hemorrhage
- heterotopic new bone formation
- extravasation of bone cement potentially resulting in but not limited to:
 - compression or irritation of nerve structures, such as the spinal cord or nerve roots, causing radiculopathy, parasthesia, paraplegia or paralysis and/or;
 - introduction into the vascular system resulting in embolism of the lung and/or heart or other clinical sequelae
- pyrexia due to allergy to bone cement
- short-term conduction irregularities
- thrombophlebitis
- transitory fall in blood pressure
- migration of hardened bone cement bolus

INSTRUCTIONS FOR THE PREPARATION AND CLINICAL USE OF KYPHON® HV-R® BONE CEMENT

Timing for the Preparation and Application of KYPHON® HV-R® Bone Cement

The handling characteristics of bone cements are affected by operating room conditions, including the room temperature, temperature of the cement components prior to mixing, humidity, the geometry of the mixing apparatus, time spent mixing, and the geometry of the delivery device. Any change in one or more of these conditions can alter the handling characteristics of the bone cement, including the following:

- Handling period - the time it takes for the bone cement to reach the doughy state (the cement has reached the doughy state when it no longer sticks to surgical gloves)
- Working period - the time the bone cement remains in the doughy state and can be delivered
- Hardening period - the time it takes for the bone cement to harden or until it can no longer be delivered

The user must be aware of these factors and adjust technique to account for variability in operating room conditions.

Under specific conditions in our laboratory, KYPHON® HV-R® Bone Cement had the following dough time, hardening time, and working period, using KYPHON® Bone Filler Devices, which allowed sufficient time for minimally invasive surgical introduction of the bone cement.

Average Room Temperature at Start of Mixing (±1°C)	Average Dough Time (minutes)	Average Hardening Time (minutes)	Average Working Period (minutes)
15	23.8	40.0	16.2
19	14.0	33.8	19.8
23	8.9	19.3	10.4
25	7.6	18.4	10.8
29	4.6	12.3	7.7

Lower temperatures or other changes in operating room conditions can increase the handling, doughy (working) and hardening periods. Conversely, higher temperatures or other changes in operating room conditions can decrease the handling, doughy (working) and hardening periods.

Handling Characteristics of KYPHON® HV-R® Bone Cement at 23 ± 1°C in our Laboratory

Period	Activity	Approximate Cumulative Time From Initiation of Mixing*
Mixing	Mix liquid and powder	0-2 minutes
Handling	Transfer into delivery system	2-8 minutes
Working (doughy state)	Fill cavity in vertebral body	8-16 minutes
Hardening	Wait before completing procedure	16-20 minutes

Note: These cumulative time periods will vary depending on temperature and other factors. For example, the colder the environment, the longer the time necessary for the cement to develop the required doughy consistency. Warmer temperatures require more rapid preparation and handling (see Figure 1). Ensure the cement's viscosity is high enough (doughy) before delivery begins. The cement has reached the doughy state when it no longer sticks to surgical gloves.

* These times are based on cement prepared by Medtronic. Times were obtained from bone cement mixed in a KYPHON® Mixer and delivered through a KYPHON® Bone Filler Device (inner diameter = 0.114"). Times may vary when other mixing methods, delivery devices and/or vacuum are used.

Figure 1


Figure 1. The effect of room temperature on the handling characteristics of KYPHON® HV-R® Bone Cement is shown in the above time vs. temperature graph. At room temperatures of 15°C, 19°C, 23°C, 25°C and 29°C, units of KYPHON® HV-R® Bone Cement were mixed in a KYPHON® Mixer and transferred to KYPHON® Bone Filler Devices (BFD, size 3). Dough time and hardening time were measured. The dashed line represents dough time (when dispensed cement no longer sticks to a surgical glove) and the solid line represents hardening time (when cement is too hard to dispense).

Preparation Procedures

- Per the note above, temperature can affect the handling of the cement. Prior to use, it is advised to keep the product at a temperature of 23 ± 1°C for a period of 24 hours.
- Prior to use, KYPHON® HV-R® Bone Cement packaging should be examined for damage and the presence of all required components.
- Maintain aseptic transfer surgical technique to prevent possible infection.
- Maintain strict adherence to the instructions for mixing the powder and liquid to prevent possible dermatitis from the liquid monomer during handling.
- Assure that the inner package is undamaged, that the powder is not discolored (yellow or brown) and the liquid is not syrupy. These conditions indicate that the product has not been stored correctly.
- Assure that the preparation accessories are specifically compatible with the bone cement product.
- Do not open the vial of liquid over the mixing bowl to avoid the risk of glass fragments entering the dough.
- Never add other substances or foreign bodies to the acrylic resin. Never modify the ratios between the liquid and solid components.
- Care should be taken in the mixing of the liquid and powder components such that the entire contents of the vial and packet are utilized. The liquid monomer and the powder component should be thoroughly mixed.

Mixing Procedure

1. Always use sterile technique when mixing the bone cement.
2. Open powder packet and break open the vial. Immediately combine the powder and liquid into a mixing bowl or KYPHON® Mixer.
3. When using the KYPHON® Mixer, follow the KYPHON® Mixer Instructions For Use.
4. When mixing by hand, mix with a spatula until the ingredients are smoothly and uniformly combined, approximately one minute or longer as necessary. The liquid must moisten all the powder; use the spatula to carefully work any lumps of unmoistened powder into the overall mass of moist dough.
5. After the powder and liquid have been thoroughly mixed, the transfer of the mixture into a cement delivery device(s) may begin. Allow for complete dissolution of the powder into the liquid prior to transferring the bone cement into the cement delivery device(s). This can take up to approximately eight minutes after the start of mixing.

Clinical Application Procedures

When the bone cement is in a doughy state, it has reached the appropriate level of viscosity and is ready to be placed into the pathological fracture. Do not insert the bone cement into the cavity of the vertebral body until the bone cement has reached the doughy state. The cement has reached the doughy state when it no longer sticks to surgical gloves.

1. Before delivery, dispense a small sample of the bone cement from the delivery device to ensure proper viscosity. The bone cement is not ready to deliver until it does not drip from the cement delivery device, does not stick to surgical gloves and has lost its sheen.
2. Always create a cavity within the vertebral body before delivering the cement.

3. Delivery of the bone cement into the vertebral body during balloon kyphoplasty should begin in the anterior or central region of the bone cavity. Continue by allowing the bone cement to flow into the posterior area without moving the delivery device posteriorly. Care should be taken to avoid placing the delivery device against the anterior cortex. Do not deliver the bone cement into the vertebral body before creating a cavity.
4. While the bone cement hardens, it is important to maintain patient positioning until the end of the polymerization or hardening process.

HOW SUPPLIED

KYPHON® HV-R® Bone Cement is supplied sterile. The powder and package are sterilized with gamma radiation. The liquid is sterilized using filtration and is contained in a glass vial. The outside of the glass vial is sterilized with ethylene oxide gas. This device is intended for single use only. Do not resterilize. Do not use if package is open or damaged.

STORAGE

KYPHON® HV-R® Bone Cement should be stored in its original shipping materials. Proper care should be taken to ensure that KYPHON® HV-R® Bone Cement will not be damaged. Store below 25°C and away from sunlight.

LIMITATION OF LIABILITY

MEDTRONIC WILL NOT BE RESPONSIBLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL OR EXEMPLARY DAMAGES RESULTING FROM REUSE OF KYPHON® HV-R® BONE CEMENT.

IN NO EVENT SHALL MEDTRONIC BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, CONSEQUENTIAL, OR EXEMPLARY DAMAGES ARISING OUT OF OR IN CONNECTION WITH KYPHON® HV-R® BONE CEMENT, BASED UPON BREACH OF CONTRACT (INCLUDING BREACH OF WARRANTY).

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© 2011 Medtronic Sofamor Danek USA, Inc. All rights reserved.
TECRES S.p.A.
Medtronic Sofamor Danek USA, Inc.
Tel: 800 933 2635 (In U.S.A.)
901 396 3133 (Outside of U.S.A.)