PTCA Dilatation Catheter

Cathy n°4

*Translumina -
Balance is the key*

Cathy n° 4 - Trackability and pushability in balance

The Translumina **Cathy n°4** represents the next generation catheter technology with an optimum in pushability and trackability. This Balloon Catheter has been developed in order to increase the success rate when crossing the lesions. Improvements with soft and hard hightech materials result in a perfect product balance for easy navigation.

Technical Data

Shaft material proximal	PTFE coated Hypo-Tube
Shaft material distal	Polyamide
Balloon material	Modified Polyamide
Balloon marker material	Platinum / Iridium
Proximal shaft diameter	1,9 F
Distal shaft diameter	2,7 F
Entry profile	0,41 mm
Balloon crossing profile	0,64 mm
Recommended guide wire	0,014"
Usable length	143 cm
Length of GW-Lumen	27 cm
Guiding Catheter	5 F

The next generation in balance

Translumina makes great demand on its product development. The new **Cathy n°4** stands for a well balanced concept. It combines outstanding design features to offer our customers a perfect product for the most challenging cases.

What's new?

The new distal shaft

Our high performance shaft consisting of a new material provides excellent pushability and kink resistance. This feature allows for high manoeuvrability.

The new designed luer

The transparent luer has a positive, tactile feel assisting in navigation of the system. It is designed with an integrated protection to minimize any kinking.

The new flexible tip

The soft tip material combined with an improved robust segment ensures perfect crossability and trackability. This feature allows easy access to all lesions.

Cathy n°4

Product matrix / Ordering information*

Small vessel sizes

Balloon Ø [mm]	Balloon length [mm] & Article number			
	12	15	20	30
Ø 1,50		T-CCC1515	T-CCC1520	
Ø 2,00	T-CCC2012	T-CCC2015	T-CCC2020	T-CCC2030
Ø 2,25	T-CCC22512	T-CCC22515	T-CCC22520	T-CCC22530
Ø 2,50	T-CCC2512	T-CCC2515	T-CCC2520	T-CCC2530

new sizes available

Medium vessel sizes

Balloon Ø [mm]	Balloon length [mm] & Article number			
	12	15	20	30
Ø 2,75	T-CCC27512	T-CCC27515	T-CCC27520	T-CCC27530
Ø 3,00	T-CCC3012	T-CCC3015	T-CCC3020	T-CCC3030
Ø 3,50	T-CCC3512	T-CCC3515	T-CCC3520	T-CCC3530
Ø 4,00	T-CCC4012	T-CCC4015	T-CCC4020	T-CCC4030

* Please contact our customer service for available sizes

Compliance chart

Ballon Ø [mm]	Inflation pressure [bar or 10 ⁻⁵ Pascal]														
	NP											RBP			
	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
Ø 1,50	1,39	1,43	1,46	1,50	1,53	1,57	1,60	1,64	1,67	1,71	1,74	1,78	1,81	1,85	1,88
Ø 2,00	1,90	1,94	1,97	2,00	2,03	2,07	2,10	2,13	2,17	2,20	2,23	2,27	2,35	2,38	2,41
Ø 2,25	2,15	2,18	2,21	2,25	2,28	2,32	2,35	2,38	2,42	2,45	2,49	2,52	2,55	2,59	2,62
Ø 2,50	2,40	2,43	2,47	2,50	2,54	2,57	2,60	2,64	2,67	2,71	2,74	2,77	2,81	2,84	2,88
Ø 2,75	2,65	2,68	2,72	2,75	2,78	2,82	2,85	2,88	2,92	2,95	2,98	3,01	3,28	3,31	3,34
Ø 3,00	2,89	2,93	2,97	3,00	3,04	3,08	3,12	3,15	3,19	3,23	3,26	3,30	3,34	3,37	3,41
Ø 3,50	3,37	3,42	3,46	3,50	3,54	3,58	3,63	3,67	3,71	3,75	3,79	3,84	3,88	3,92	3,96
Ø 4,00	3,85	3,90	3,95	4,00	4,05	4,10	4,16	4,21	4,26	4,31	4,36	4,41	4,46	4,51	4,56

CE 0124

© 2009 Translumina GmbH

All Rights Reserved

Please refer to the Instruction for Use supplied with these devices for indications, contraindications, adverse effects, suggested procedures, warnings and precautions.

Translumina GmbH | Neue Rottenburger Strasse 50
72379 Hechingen | Germany | info@translumina.de |
t + 49 7471 98 94 - 0 | f + 49 7471 98 94 - 360

EC Design- Examination Certificate

according to directive 93/42/EEC,
annex II (4)

As a notified body of the European Union, DEKRA Certification GmbH certifies for the manufacturer

translumina GmbH

Neue Rottenburger Straße 50, 72379 Hechingen, Germany

that the design dossier for the product(s) described in the annex complies with the requirements of the directive 93/42/EEC. This certificate is based on the result of the examination of the design dossier according to the directive 93/42/EEC annex II.4 as documented in the report mentioned in the Annex.

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Annex to the EC design examination certificate 50668-23-M2 dated 2016-02-23

Revision status: 0

Date: 2016-03-21

Page 1 of 1

Report number: 50668-P13-04

Product: Cathy No. 4 PTCA Catheter

Intended use:

The PTCA Balloon Catheter is indicated for dilatation in stenosed coronary artery or aorto-coronary bypass grafts as well as for dilatation of stenosed tibial anterior and posterior artery of the lower limb with a maximum lesion length of 45 mm.

Technical data:

Catheter lengths of use (mm)	1430
Balloon diameter (mm)	1,50 / 2,00 / 2,25 / 2,5 / 2,75 / 3,0 / 3,5 / 4,0
Balloon lengths (mm)	9,5 to 41,5
Article codes PTCA Catheter	T-CCC1509 T-CCC2012, T-CCC22512, T-CCC2512, T-CCC27512, T-CCC3012, T-CCC3512, T-CCC4012 T-CCC1515, T-CCC2015, T-CCC22515, T-CCC2515 T-CCC27515, T-CCC3015, T-CCC3515, T-CCC4015 T-CCC1520, T-CCC2020, T-CCC22520, T-CCC2520 T-CCC27520, T-CCC3020, T-CCC3520, T-CCC4020 T-CCC2030, T-CCC22530, T-CCC2530, T-CCC27530 T-CCC3030, T-CCC3530, T-CCC4030

EC Declaration of Conformity

By the responsible manufacturer:

translumina®

translumina GmbH
Neue Rottenburger Str. 50
72379 Hechingen
Germany

Zertifikat – Certificate – Certificat

Translumina GmbH declares under its sole responsibility that the below entitled Medical Device complies to the requirements of the Medical Device Directive 93/42/EEC of the European Union.

The Quality Management System is in conformity with the requirements and the Conformity Assessment Procedure was followed in accordance to MDD 93/42/EEC, Annex II.

The device **Cathy No.4 PTCA Catheter** in all variations is classified as Class III as per Annexure IX, Rule 6 of MDD 93/42/EEC.

Issue of this certificate:

Cathy No.4 PTCA Catheter

Ref. No.:

T-CCC1509
T-CCC1515
T-CCC1520


T-CCC2012	T-CCC22512	T-CCC2512	T-CCC27512	T-CCC3012	T-CCC3512	T-CCC4012
T-CCC2015	T-CCC22515	T-CCC2515	T-CCC27515	T-CCC3015	T-CCC3515	T-CCC4015
T-CCC2020	T-CCC22520	T-CCC2520	T-CCC27520	T-CCC3020	T-CCC3520	T-CCC4020
T-CCC2030	T-CCC22530	T-CCC2530	T-CCC27530	T-CCC3030	T-CCC3530	T-CCC4030

Certificate Registration Number 50668-23-M2


Notified Body's No. 0124
Notified Body DEKRA Certification
Handwerkstrasse 15
70565 Stuttgart
Germany

This certificate is valid from **2016-03-21** to **2021-03-20**

Hechingen, 2016-03-21

Juergen Ettenhuber
Manager QA/RA

Joerg Weber
Manager R&D

[Vers.2014/01]

This document has to be maintained at least five years after the last production of the concerned product was achieved.

Cathy N° 4

Angioplasty Balloon Dilation Catheter
Ballon-Dilatations-Katheter für Angioplastie
Sonde à ballonnet de dilatation pour angioplastie
Catetere a palloncino da dilatazione per PTCA
Catéter de dilatation con balón para angioplastia
Caterer Balão Dilatador para Angioplastia

INSTRUCTIONS FOR USE

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MODE D'EMPLOI

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
Table of Contents:

1. Description	3
2. Indications	3
4. Warnings and Precautions	4
5. Packaging and Storage	5
6. Specifications	5
7. Application	5
8. Potential Complications.....	8
9. Warranty and Limitation of Liability	8

These instructions are intended for all personnel who will handle the „Cathy N° 4“, a PTCA Catheter, and should be read carefully before the product is used!

Product Label – Symbol Description

	Manufacturer		Batch code
	Date of manufacture		Catalogue number
	Expiry Date		Sterilized using ethylene oxide
	Do not use if packaging is damaged		Consult instructions for use
	Do not re-use		Caution
	Keep away from sunlight		Keep dry
	Nominal diameter of balloon		Rated burst pressure
	Nominal length of balloon		Nominal pressure
	Do not re-sterilize		

Please Note:

These instructions for use relate to all available balloon lengths and diameters of the PTCA Catheter Cathy N° 4.

Information on Sterilization:

The product has been sterilized with ethylene oxide. It is non-pyrogenic. Do not use if the packaging has been opened, damaged or the expiration date has been exceeded!

1. Description

The Cathy N° 4 PTCA balloon catheter contains of a shaft system with an integrated balloon near the distal tip. The shaft has a combination of a single lumen hose, which divides into two lumina in the area of the distal tip. The longer lumen, which is accessed through the Luer connector, is used for inflation of the balloon with diluted contrast medium. The second, distal, shorter lumen which is accessed through the catheter tip, is used to insert guide wire, which permits the exact positioning of the catheter in the area of the stenosis. The balloon has been designed to expand to an exact defined diameter at a specific pressure (semi-compliant balloon material). The nominal diameter stated is reached at nominal pressure.

In order to facilitate the positioning of the balloon in the area of the stenosis, several markers have been integrated. Two platinum / iridium markers indicate the balloon length. Two additional markers on the proximal catheter shaft facilitate estimation of the positioning of the dilatation catheter in relation to the tip of the guiding catheter. The marker nearer to the Luer-Lock adapter of the dilatation catheter is used in femoral application, the other one in brachial application.

The design of this catheter does not incorporate a lumen for distal injections and distal measurements.

2. Indications

- The PTCA catheter Cathy No 4 is indicated for balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis or the tibial anterior and posterior artery of the lower limb with a maximum lesion length of 45 mm.
- The perfusion can be improved by the dilatation.
- Patients being considered for balloon angioplasty should be acceptable candidates for such treatment.

3. Contraindications**3.1 The coronary application of the Cathy N° 4 is contra-indicated for:**

- Patients with hemorrhagic diathesis or other disorders in which the application of anticoagulant or antiplatelet therapy is contraindicated (e.g. septic ulceration, cerebrovascular complications)
- Patients in whom aortocoronary bypass procedures are contraindicated
- Patients with imminent thrombus formation and alteration of flow after myocardial infarction
- Patients who have experienced a myocardial infarction less than one (1) week prior to the proposed procedure
- Patients with proximal atherosclerosis and/or significant vessel tortuosity in whom guide wire access and adequate guiding catheter support is prohibited
- Coronary artery spasm in the absence of a significant stenosis
- Patients with unprotected left main coronary vessel disease
- Lesions that cannot be successfully pre-dilated
- Patients with known allergies to the used materials
- Patients in whom the necessary pharmaceutical treatments are contraindicated
- Severe allergic reaction to the contrast medium, stent and/or delivery system material and/or the necessary pharmaceutical treatment (e.g. aspirin)
- Women during pregnancy

3.2 The peripheral application of the Cathy N° 4 is contra-indicated for:

- Limited Access to the target vessel
- Diffuse vessel disease
- extremely calcified lesions
- angiographic indication of a massive thrombus amount next to or inside the target lesion
- lesions located in or next to an aneurysm
- extremely tortuous vessels
- extremely angled lesion
- contra-indication against thrombocyte aggregation inhibitors and/or coagulation inhibitors
- Allergy against contrast medium, stent and/or delivery system material and/or the necessary pharmaceutical treatment

4. Warnings and Precautions

4.1 In General:

The medicinal techniques and procedures described here do not entail all the procedures of PCTA acknowledged in medicine.

- Use the catheter prior to the expiration date noted on the packaging.
- Check before use that the packaging, in particular the sterile packaging, has not been damaged or opened.
- The device is intended for single use only. It is designed for transitional use (<60 min.) Do not resterilize and/or reuse it, as this can potentially result in device malfunction and increased risk of a lack of sterility and cross contamination.
- The Cathy N° 4 PTCA catheter should only be used by physicians trained in PTCA procedures.
- Only use diluted contrast medium.
- Do not use air or any gaseous medium to inflate the balloon.
- The balloon pressure must not exceed the rated burst pressure (RBP) specified on the product label.
- Only use an inflation syringe with a manometer to inflate the balloon.
- Administer appropriate anticoagulation and vasodilatation therapy before inserting the dilatation catheter.

4.2 During Procedure

- PTCA should only be performed in hospitals where emergency coronary artery bypass graft surgery can be performed immediately in the event of injuries or life threatening complications. During the entire PTCA procedure, an emergency bypass surgery team should stand by.
- Only use fluoroscopic control to insert the Cathy N° 4 balloon catheter into the vessel. Only move the balloon catheter forward and backward if it has been completely deflated under vacuum.
- Do not overtighten the hemostatic valve, as lumen constriction may occur, impairing the inflation / deflation of the balloon.
- The inflation diameter of the balloon may not exceed the diameter of the artery proximal and distal to the stenosis.
- Do not attempt to straighten a catheter if the shaft is bent or kinked, as it breaks. Instead, prepare a new catheter.
- Avoid contact with organic solvents (such as alcohol)

4.3 Compatibility with MRI

The Cathy No. 4 PTCA catheter must not be used during an MRI (magnetic resonance imaging) examination.

4.4 Advice due to the single-use

The product „Cathy N° 4 PTCA catheter“ is designed for the single-use only. The reconditioning for a re-use is prohibited. – The product is exposed high mechanical loads during its use and possibly affected thereafter. Because of that only for the single-use a save application can be warranted.

5. Packaging and Storage

Store in a dry, dark and cool place. Use before the expiration date (indicated on the product label).

6. Specifications

6.1 Balloon Compliance

Tab.: Ratio pressure / balloon diameter.

- The grey areas mark the diameter at rated burst pressure.
- Nominal pressure = 9 [bar, atm or 10⁵ Pa]

Pressure [bar, atm or 10⁵ Pascal]

Balloon ø [mm]	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	20
ø 1.5	1.29	1.32	1.36	1.39	1.43	1.46	1.5	1.53	1.57	1.6	1.64	1.67	1.71	1.74	1.78	1.81	1.88
ø 2.0	1.8	1.84	1.87	1.9	1.94	1.97	2	2.03	2.07	2.1	2.13	2.17	2.2	2.23	2.27	2.35	2.41
ø 2.25	2.04	2.08	2.11	2.15	2.18	2.21	2.25	2.28	2.32	2.35	2.38	2.42	2.45	2.49	2.52	2.55	2.62
ø 2.5	2.3	2.33	2.37	2.4	2.43	2.47	2.5	2.54	2.57	2.6	2.64	2.67	2.71	2.74	2.77	2.81	2.88
ø 2.75	2.55	2.59	2.62	2.65	2.68	2.72	2.75	2.78	2.82	2.85	2.88	2.92	2.95	2.98	3.01	3.28	3.34
ø 3.0	2.78	2.82	2.86	2.89	2.93	2.97	3	3.04	3.08	3.12	3.15	3.19	3.23	3.26	3.3	3.34	3.41
ø 3.5	3.25	3.29	3.33	3.37	3.42	3.46	3.5	3.54	3.58	3.63	3.67	3.71	3.75	3.79	3.84	3.88	3.96
ø 4.0	3.7	3.75	3.8	3.85	3.9	3.95	4	4.05	4.1	4.16	4.21	4.26	4.31	4.36	4.41	4.46	4.56

6.2 Material

- Catheter shaft proximal: PTFE coated hypo-tube
- Catheter shaft distal: Polyamide
- Balloon: modified polyamide

7. Application

Please Note:

All personnel handling this device should read current literature about PTCA, as published, for example, by the American Heart Association (AHA) or the American College of Cardiology (ACC).

7.1 Accessories

- Introduction sheath, compatible with the guiding catheter
- Femoral or brachial guiding catheter (5 F min.)
- Hemostatic valve
- 60% contrast medium diluted 1:1 with normal saline
- Sterile, heparinized normal saline
- 20 cc Luer-Lock syringe (optional)
- Inflation syringe with manometer
- 3-way stopcock for a minimum of 25 bar

- Guide wire 0.014"
- Guide wire introducer
- Guide wire turning device

7.2 Prior to use

- Carefully check the sterile packaging of all equipment to be used for this procedure. Do not use if the packaging has been opened or damaged.
- Do not use damaged or kinked instruments.
- Please remember that the balloon diameter stated at nominal pressure should not be bigger than the diameter of the artery above and below the stenosis.
- Check the guide wire size according to the product label.

Please Note:

When using a hydrophilic coated guide wire, the coating could swell. In this case, the guide wire should be selected smaller.

7.3. Preparation

1. Prepare the inflation syringe according to the manufacturer's instructions.
2. Prepare diluted contrast medium and sterile saline solution.
3. Fill a 10 or 20 cc syringe with 4 cc of contrast medium.
4. Carefully remove the balloon catheter from the sterile packaging. Work in sterile technique. Pull the catheter out of the transport tubing.
5. Connect a pressure-resistant 3-way stopcock (minimum 25 bar) to the Luer connector of the balloon catheter.
6. Connect the syringe to the 3-way stopcock, with the result that the plunger points vertically downwards.
7. Generate a strong vacuum in the balloon by pulling the syringe plunger to the back. Maintain the vacuum for 15–20 seconds and make sure that no more bubbles are seen within the diluted contrast medium.
8. Release the vacuum carefully.
9. Repeat this procedure if necessary until no more bubbles can be seen in the system.

Warning:

If it is not possible to create or to maintain the vacuum, check the connections. If there are no errors here, do not use this catheter under any circumstances and send it back to translumina.

10. Create another vacuum in the catheter and close the 3-way stopcock.
11. Remove the air from the tubing and connector of the inflation syringe.
12. Connect the inflation to the 3-way stopcock of the balloon catheter.
13. Open the stopcock carefully. The inflation syringe should be in a vertical position in order to avoid getting air into the system.
14. Make sure that the transport protection has been removed.
15. For flushing the guide wire lumen, attach a needle to the syringe filled with saline solution and insert it into the distal tip of the catheter. Then apply slight pressure to the syringe plunger.
16. The balloon catheter is now prepared and ready for use.

7.4 Use of the Cathy N° 4 PTCA Catheter

1. Insert a suitable sheath into the femoral or brachial artery using conventional techniques.

2. Insert the guiding catheter into the vessel through the introduction sheath and fix the hemostatic valve to it.
3. Insert the guide wire into the guiding catheter and carefully push it forwards under fluoroscopic control. Please pay attention to the admissible guide wire size on the product label.
4. To prevent unrestricted guide wire movement, tighten the hemostatic valve so that it closes firmly around the guide wire after the guide wire has passed the stenosis.
5. Push the proximal end of the guide wire into the distal tip of the Cathy N° 4 catheter, ensuring that the guide wire exits through the exit of the guide wire lumen. The balloon catheter must move without resistance.
6. Open the valve and carefully push the Cathy N° 4 over the guide wire to the distal end of the guiding catheter.
7. Continue the procedure under fluoroscopic control using accepted angioplasty techniques. The radiopaque balloon markers aid the positioning of the balloon in the area of the stenosis.
8. When the balloon has been positioned correctly in the stenosis, inflate it with nominal pressure to start with and keep the pressure for 10 to 90 seconds. Due to the restricted blood flow, the patient can feel heart pain. Adapt the inflation time to the patient's condition.
9. The inflation procedure can be repeated two to three times, but do not exceed the rated burst pressure (RBP) under any circumstances.

Caution:

- The balloon catheter may only be inserted into the guiding catheter if it has been fully deflated and only over the guide wire.
- Do not overtighten the hemostatic valve since this may restrict the flow of contrast medium into the balloon and thus increase the inflation/deflation time.

Please Note:

- It is recommended that the guide wire and/or the balloon catheter remain in the area of the stenosis until the procedure is completed.
- Contrast media may influence the inflation/deflation time due to different viscosities.

7.5 Catheter Exchange Procedure

1. Open the hemostatic valve.
2. Hold the guide wire and the hemostatic valve firmly in one hand while grasping the balloon catheter shaft with the other hand. Maintain the guide wire's position in the coronary vessel by holding the wire stationary and begin pulling the dilatation catheter out of the guiding catheter carefully.
3. Pull the balloon catheter until the opening in the guide wire lumen has been reached. Then proceed in the same way as with an over-the-wire system until the Cathy N° 4 has been totally removed.
4. Close the hemostatic valve.
5. Prepare the next balloon catheter (cf. 7.3).
6. Push the balloon catheter over the guide wire.
7. Open the hemostatic valve and push the balloon catheter forward and leave the guide wire in the coronary artery.
8. Be careful not to twist or rotate the balloon catheter around the guide wire.

Caution:

Before removing the balloon catheter, make sure that the balloon is completely deflated.

8. Potential Complications

8.1 Potential complications of the coronary application:

The following complications due to PTCA and stent implantations have been reported and may occur. Potential undesired effects are possible, but not limited to the following:

- Dissection, injury, rupture, sudden and total occlusion or other damages to the coronary arteries
- Thrombosis of the coronary arteries
- Unstable angina
- Arterial spasm
- Myocardial infarction
- Cardiac dysrhythmia, including ventricular fibrillation
- Embolism
- Restenosis of the dilated artery
- Hemorrhage or hematoma
- Bleeding from anticoagulation/antiplatelet medication
- Side effects of the concomitant medication
- Allergic reaction to contrast medium
- Allergic reaction to the stent and/or delivery system material
- Infections
- Arteriovenous fistula
- Emergency surgery of coronary vessel bypass graft
- Deformation of stent symmetry on deployment or on re-crossing with another device (IVUS, balloon catheter, etc.)
- Stent migration
- renal failure
- Death

8.2 Potential complications of the peripheral application:

- Sepsis / infection
- stent migration
- embolism of atherosclerotic or thrombosis material
- acute or sub-acute stent thrombosis
- amputation
- arterio-venous fistula
- Emergency surgery to treat vascular complications
- Hemorrhage or hematoma
- Side effects of the concomitant medication
- Allergic reaction to the contrast medium
- Allergic reaction to the stent and/or delivery system material
- injury of the artery including perforation and dissection
- Restenosis of the dilated artery
- Thrombus
- tissue necrosis
- total occlusion
- renal failure

9. Warranty and Limitation of Liability

translumina GmbH guarantees that this product has been manufactured, sterilized and packaged with the due care prescribed in EU Directive for Medical Devices 93/42/EEC. Each product is individually tested before it is released for packaging. **translumina GmbH** will replace any device if defects which can be proven to have originated during manufacture or packaging are immediately brought to the attention of **translumina GmbH** or its distributors. This warranty applies in lieu of and ruling out any other verbal or written, express, implicit or other warranty.

As a result of biological differences in individuals, no product is 100% effective under all circumstances. Because of this fact and since we have no control over the conditions, the instruments used, the operation and the diagnosis, **translumina GmbH** and its distributors do not guarantee either good effectivity in application; unsatisfactory results also cannot be ruled out. **translumina GmbH** and its distributors shall not be liable for any damage or losses of any nature arising directly or indirectly from the use of the device. The product liability expires if damage caused by transport, storage or improper handling by the operator or damage caused by factors outside the control of the manufacturer occur.

translumina GmbH and its distributors shall not be liable for damages arising from resterilisation or reuse of the product.

Important Information:

If more than 2 years have elapsed between the revision date stated below and today's date, please contact the manufacturer in order to find out if any changes have occurred.

For further information please contact:

translumina GmbH

Neue Rottenburger Strasse 50
D-72379 Hechingen
Phone: ++ 49 74 71 98 94 - 0
Fax: ++ 49 74 71 98 94 - 3 80
E-mail: info@translumina.de
Internet: www.translumina.de

© 2015 translumina GmbH
All rights reserved.

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® YUKON® and translumina® are registered trademarks of translumina GmbH