

Foresight[®]

HBeAg EIA Test Kit

Package Insert

REF

I231-1061

English

An enzyme immunoassay (EIA) for the qualitative detection of Hepatitis B Envelope Antigen (HBeAg) in human serum or plasma.

For professional *in vitro* diagnostic use only.

| INTENDED USE |
|--|
| The HBeAg EIA Test Kit is a one step enzyme immunoassay for the qualitative detection of Hepatitis B Envelope Antigen (HBeAg) in human serum or plasma. It is intended as a screening tool for exposure to the Hepatitis B virus (HBV), and as an aid in the evaluation of infectivity, the disease state, and management of patients. |
| SUMMARY |

Hepatitis B virus is a spherical enveloped, partially double-stranded DNA virus of the Hepadnaviridae family.¹ The Hepatitis B infection of the liver is transmitted through sexual contact, blood borne exposure, transmission from mother to child during delivery, sharing of objects that pierce the skin, child-to-child and household contact.^{2,3,4,5} HBV infection has been linked to a variety of mild to chronic liver diseases, including cirrhosis, and hepatocellular carcinoma.⁴ In some cases, the virus may persist for a lifetime. Annually, 1 million people die from chronic active hepatitis, cirrhosis or primary liver cancer. Hepatitis B affects millions of people worldwide and is considered a global public health problem.^{4,5} Hepatitis B envelope antigen (HBeAg) is produced by hepatocytes infected with the Hepatitis B virus (HBV) and is one of the serological markers of Hepatitis B infection. It appears during the viral replication period, at the most infectious phase of the disease, and is regarded as an indicator of infectivity. HBeAg appears 3-6 weeks after exposure and remains detectable up to several weeks. The presence of HBeAg indicates that the patient is infectious and its long term presence is a sign of chronic liver disease.⁵ In the chronic carrier state, HBeAg may persist in the blood for years. Upon recovery, HBeAg declines and its specific antibody (HBeAb) appears. In cases of infection with mutant strains of HBV, pre-core mutants, HBeAg is not synthesized and viral replication may be ongoing in the absence of HBeAg.^{3,5} In these cases it is recommended that the test is combined with other HBV serological marker tests for the evaluation of the disease state. The HBeAg EIA Test Kit is an immunoassay for the qualitative detection of the presence of Hepatitis B Envelope Antigen in serum or plasma specimen. The test utilizes monoclonal antibodies to selectively detect various subtypes of HBeAg in serum or plasma.

The HBeAg EIA Test Kit is a solid phase qualitative enzyme immunoassay based on a sandwich principle for the detection of HBeAg in human serum or plasma. The microwell plate is coated with monoclonal antibodies specific to HBeAg. During testing, the specimen and the enzyme-conjugated HBeAg antibodies are added to the antibody coated microwell plate and then incubated. If the specimen contains HBeAg, it will bind to the antibodies coated on the microwell plate and simultaneously bind to the conjugate to form immobilized antibody-HBeAg-conjugate complexes. If the specimen does not contain HBeAg, the complexes will not be formed. After initial incubation, the microwell plate is washed to remove unbound materials. Substrate A and substrate B are added and then incubated to produce a blue color, indicating the amount of HBeAg present in the specimen. Sulfuric acid solution is added to the microwell plate to stop the reaction which produces a color change from blue to yellow. The color intensity, which corresponds to the amount of HBeAg present in the specimen, is measured with a microplate reader at 450/630-700 nm or 450 nm.

PRECAUTIONS

- For professional *in vitro* diagnostic use only. Do not use after expiration date.
- Do not mix reagents from other kits with different lot numbers.
- Avoid cross contamination between reagents to ensure valid test results.
- Follow the wash procedure to ensure optimum assay performance.
- Use plate sealer to cover microwell plate during incubation to minimize evaporation.
- Use a new pipet tip for each specimen assayed.
- Ensure that the bottom of the plate is clean and dry and that no bubbles are present on the surface of the liquid before reading the plate. Do not allow wells to dry out during the assay procedure.
- Do not touch the bottom of the wells with pipette tips. Do not touch the bottom of the microwell plate with fingertips.
- Do not allow sodium hypochlorite fumes from chlorine bleach or other sources to contact the microwell plate during the assay as the color reaction may be inhibited.
- All equipment should be used with care, calibrated regularly and maintained following the equipment manufacturer's instructions.

HEALTH AND SAFETY INFORMATION

- Some components of this kit contain human blood derivatives. No known test method can offer complete assurance that products derived from human blood will not transmit infectious agents. Therefore, all blood derivatives should be considered potentially infectious. It is recommended that these reagents and human specimens be handled using established good laboratory working

practices.

- Wear disposable gloves and other protective clothing such as laboratory coats and eye protection while handling kit reagents and specimens. Wash hands thoroughly when finished.
- ProClin™ 300 is included as a preservative in the Conjugate, Concentrated Wash Buffer, Substrate and Controls. Avoid any contact with skin or eyes.
- Do not eat, drink or smoke in the area where the specimens or kits are handled. Do not pipette by mouth.
- Avoid any contact of the Substrate A, Substrate B, and Stop Solution with skin or mucosa. The Stop Solution contains 0.5M sulfuric acid which is a strong acid. If spills occur, wipe immediately with large amounts of water. If the acid contacts the skin or eyes, flush with large amounts of water and seek medical attention.
- Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for one hour at 121°C. Disposables should be autoclaved or incinerated. Do not autoclave materials containing sodium hypochlorite.
- Handle and dispose all specimens and materials used to perform the test as if they contained infectious agents. Observe established precautions against microbiological hazards throughout all the procedures and follow the standard procedures for proper disposal of specimens.
- Observe Good Laboratory Practices when handling chemicals and potentially infectious material. Discard all contaminated material, specimens and reagents of human origin after proper decontamination and by following local, state and federal regulations.
- Neutralized acids and other liquids should be decontaminated by adding sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to a 1.0% sodium hypochlorite may be necessary to ensure effective decontamination.

STORAGE AND STABILITY

- Unopened test kits should be stored at 2-8°C upon receipt. All unopened reagents are stable through the expiration date printed on the box if stored between 2-8°C. Once opened, all reagents are stable for up to 3 months after the first opening date if stored between 2-8°C. Return reagents to 2-8°C immediately after use.
- Allow the sealed pouch to reach room temperature before opening the pouch and remove the required number of strips to prevent condensation of the microwell plate. The remaining unused strips should be stored in the original resealable pouch with desiccant supplied at 2-8°C and can be used within 3 months of the opening date. Return the remaining unused strips and supplied desiccant to the original resealable pouch, firmly press the seal closure to seal the pouch completely and immediately store at 2-8°C.
- Concentrated Wash Buffer may be stored at room temperature to avoid crystallization. If crystals are present, warm up the solution at 37°C. Working Wash Buffer is stable for 2 weeks at room temperature.
- Do not expose reagents especially the Substrate to strong light or hypochlorite fumes during storage or incubation steps.
- Do not store Stop Solution in a shallow dish or return it to the original bottle after use.

SPECIMEN COLLECTION AND PREPARATION

- The HBeAg EIA Test Kit can be performed using only human serum or plasma collected from venipuncture whole blood.
- EDTA, sodium heparin, and ACD collection tubes may be used to collect venipuncture whole blood and plasma specimens. The preservative sodium azide inactivates horseradish peroxide and may lead to erroneous results.
- Separate serum or plasma from blood as soon as possible to avoid hemolysis. Grossly hemolytic, lipidic or turbid samples should not be used. Specimen with extensive particulate should be clarified by centrifugation prior to use. Do not use specimens with fibrin particles or contaminated with microbial growth.
- Serum and plasma specimens may be stored at 2-8°C for up to 7 days prior to assaying. For long term storage, specimens should be kept frozen below -20°C.
- Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.
- If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiologic agents.

| REAGENTS AND COMPONENTS | | | | |
|-------------------------|--------------------------------|---|--------------------------|---------------------------|
| Materials Provided | | | | |
| No. | Reagent | Component Description | Quantity | |
| | | | 96 wells/kit | 480 wells/kit |
| | HBeAg Microwell Plate | Microwell plate coated with Anti-HBeAg | 1 plate (96 wells/plate) | 5 plates (96 wells/plate) |
| 1 | HBeAg Conjugate | Anti-HBeAg bound to peroxidase; Preservative: 0.1% ProClin™ 300 | 1 x 8 mL | 5 x 8 mL |
| 2 | Concentrated Wash Buffer (25x) | Tris-HCl buffer containing 0.1% Tween 20; Preservative: 0.1% ProClin™ 300 | 1 x 40 mL | 5 x 40 mL |

| | | | | |
|---|------------------------|--|----------|----------|
| 3 | Substrate A | Citrate-phosphate buffer containing hydrogen peroxide; Preservative: 0.1% ProClin™ 300 | 1 x 8 mL | 5 x 8 mL |
| 4 | Substrate B | Buffer containing tetramethylbenzidine (TMB); Preservative: 0.1% ProClin™ 300 | 1 x 8 mL | 5 x 8 mL |
| 5 | Stop Solution | 0.5M Sulfuric acid | 1 x 8 mL | 5 x 8 mL |
| 6 | HBeAg Negative Control | Normal serum non-reactive for HBeAg, HCV, HIV-1, and HIV-2; Preservative: 0.1% ProClin™ 300 | 1 x 1 mL | 5 x 1 mL |
| 7 | HBeAg Positive Control | Inactivated serum containing HBeAg and negative for HCV, HIV-1, and HIV-2; Preservative: 0.1% ProClin™ 300 | 1 x 1 mL | 5 x 1 mL |
| | Plate Sealers | | 2 | 10 |
| | Package Insert | | 1 | 1 |

Materials Required But Not Provided

- Freshly distilled or deionized water
- Sodium hypochlorite solution for decontamination
- Absorbent paper or paper towel
- Water bath or incubator capable of maintaining 37°C ± 2°C
- Calibrated automatic or manual microwell plate washer capable of aspirating and dispensing 350 µL/well
- Disposable gloves
- Calibrated micropipettes with disposable tips capable of dispensing 50 µL
- Graduated cylinders for wash buffer dilution
- Vortex mixer for specimen mixing (optional)
- Timer
- Disposable reagent reservoirs
- Calibrated microplate reader capable of reading at 450 nm with a 630-700 nm reference filter, or reading at 450 nm without a reference filter
- Automated processor (optional)

DIRECTIONS FOR USE

Allow reagents and specimens to reach room temperature (15-30°C) prior to testing. The procedure must be strictly followed. Assay must proceed to completion within time limits. Arrange the controls so that well A1 is the Blank well. From well A1, arrange the controls in a horizontal or vertical configuration. The procedure below assigns specific wells arranged in a vertical configuration. Configuration may depend upon software.

| Step | Detailed Procedure | Simplified Procedure |
|------|---|--|
| | <ul style="list-style-type: none">Prepare Working Wash Buffer by diluting the Concentrated Wash Buffer 1:25. Pour the contents of the bottle containing the concentrated wash buffer in a graduated cylinder and fill it with freshly distilled or deionized water to 1000 mL for 96 wells/plate testing. The Working Wash Buffer is stable for 2 weeks at 15-30°C.Note: If crystals are present in the Concentrated Wash Buffer, warm it up at 37°C until all crystals dissolve.Remove unused strips from the microwell plate, and store in the original resealable pouch at 2-8°C. | <ul style="list-style-type: none">Prepare Working Wash Buffer by diluting the Concentrated Wash Buffer 1:25Remove and store unused strips at 2-8°C |
| 0 | <ul style="list-style-type: none">Leave A1 as Blank well. | <ul style="list-style-type: none">Leave A1 as Blank well |
| 1 | <ul style="list-style-type: none">Add 50 µL of Negative Control in wells B1 and C1. (Blue Reagent)Add 50 µL of Positive Control in wells D1 and E1. (Red Reagent)Add 50 µL of specimen to assigned wells starting at F1. | <ul style="list-style-type: none">B1 and C1: Add 50 µL Negative ControlD1 and E1: Add 50 µL Positive ControlStarting F1: Add 50 µL specimen |
| 2 | <ul style="list-style-type: none">Add 50 µL of Conjugate to each well except for the Blank well. (Red Reagent) | <ul style="list-style-type: none">Add 50 µL of Conjugate to each well except for the Blank well |
| 3 | <ul style="list-style-type: none">Mix gently by swirling the microwell plate on a flat bench for 30 seconds.Cover the microwell plate with the Plate Sealer and incubate in a water bath or an incubator at 37°C ± 2°C for 30 minutes ± 2 minutes. | <ul style="list-style-type: none">Mix gentlyCover the microwell plate with the Plate Sealer and incubate at 37°C for 30 min |
| 4 | <ul style="list-style-type: none">Remove the Plate Sealer.Wash each well 5 times with 350 µL of Working Wash Buffer per well, then remove the liquid.Turn the microwell plate upside down on absorbent tissue for a few seconds. Ensure that all wells have been completely washed and dried.Note: Improper washing may cause false positive results | <ul style="list-style-type: none">Remove the Plate SealerWash each well 5 times with 350 µL of Working Wash BufferTurn the microwell plate upside down on absorbent tissue |

