

The image shows a laboratory setting with various pieces of equipment. In the foreground, a white control panel for a centrifuge is visible, featuring a green LCD screen and several buttons. The screen displays '08:15', 'Power: 100%', and 'Temp: 4°C'. Below the screen, there are buttons for '088203', a power button, a stop button, and a start button. The 'orto alresa' logo is printed on the bottom right of the panel. In the background, a stainless steel centrifuge rotor is mounted on a black base, and a dark-colored centrifuge rotor is also visible. The overall scene is brightly lit, suggesting a clean and professional laboratory environment.

orto alresa

About centrifugation

Other laboratory products | DISTILLERS

DISTILLERS



Ortoalresa's distillers allow obtaining distilled water from the water supply, with ideal characteristics for its use in other equipment and appropriate for laboratory use: reagents preparation, bacteriological cultures, final cleaning of glassware, etc.

Designed for applications according the standard: ASTM D 1193

Features

- Automatic distillers of continuous production and with water flow control.
- Stainless steel for the inner and the inner heater, stainless steel painted with epoxy for housing.
- Reduced size.
- Easily removable for its cleaning.
- High quality of distillation: obtains water of types III and IV (laboratory degree water)
- Conductivity at 20°C: 1,5 microsiemens/ cm.
- Resistivity at 20°C: 0,67 megaohms/ cm.
- PH outcoming water between 5-8 (DA 006) at 25° C.

Easy to use

- Control panel with general switch and light signals.
- Water in-put connection adjustable to the feed tube.
- Cooling water drain outlet connector adaptable to containers.

Safety

- Safety mechanism that disconnects the resistances in case of cooling water failure and connects them again when the cooling water recovers.
- Electric: earthed.
- Sealed with a silicone gasket.
- Safety system with hydraulic manager of temperature that protects the distiller in case of overheating.
- Optional accessory for limiting the flow of water optimizing the flow down to the necessary minimum.

EU Directives: 2011/65/EU, 2012/19/EU, 2014/30/EU, 2014/35/UE.
Standards: EN 61326-1, EN 61010-1, EN 61010-2-101.



Versions

	Dimensions (mm) (w x d x h)			Net Weight (Kg)	Voltage (V)	Frequency (Hz)	Power (W)	Capacity (l/h)	Refrigeration water (l/h)
DA 005	370	220	440	8	220	50-60	3.000	4	60
DA 006	370	260	640	14	220	50-60	6.000	8	84
DA 007	370	220	440	8	110	50-60	3.000	4	60

Accessories

	Description
PP 354	Plastic tank for 30 litres with tap.
PV 192	Flow water limiter DA 005/007
PV 193	Flow water limiter DA 006

Water purification

After the distillation process we obtain water with quality type III and IV, this depends directly on incoming water quality and other environmental factors.

Below we show the effectiveness of the distillation process against different compounds and organisms:

Distillation

Pyrogen and virus	3 water droplets
Bacteria	3 water droplets
Particles	3 water droplets
Dissolved inorganic solutes	2 water droplets
Dissolved inorganic gases	2 water droplets
Dissolved organic	2 water droplets






orto
alresa
About centrifugation

Álvarez Redondo S.A. 

Los frailes, 121
Pol. Ind. Los Frailes
28814 Daganzo | Madrid

www.ortoalresa.com 

(+34) 91 884 40 16 

info@ortoalresa.com to ask for information 
sat@ortoalresa.com to contact our technical service
sales@ortoalresa.com to contact our comercial department
marketing@ortoalresa.com to contact our marketing department