Instruction for use

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**«HSV 1,2 IgG-ИФА»**

**A SOLID-PHASE ENZYME IMMUNOASSAY  
FOR THE QUALITATIVE DETERMINATION  
OF IgG ANTIBODIES TO HSV 1, 2  
IN HUMAN SERUM OR PLASMA**

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
Authorized Representative in EU:  
Polmed.de  
Steinacker 20, D-73773  
Aichwald, Germany  
e-mail: info@polmed.de

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**CONTENT**

1. INTENDED USE	9
2. SUMMARY AND EXPLANATION	9
3. PRINCIPLE OF THE TEST	9
4. WARNINGS AND PRECAUTIONS	10
5. KIT COMPONENTS	11
6. SPECIMEN COLLECTION AND STORAGE	12
7. TEST PROCEDURE	12
8. QUALITY CONTROL	14
9. CALCULATION OF RESULTS	14
10. EXPECTED VALUES	15
11. LITERATURE	15

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*Instruction for use*

## **A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUALITATIVE DETERMINATION OF IgG ANTIBODIES TO HSV 1,2 IN HUMAN SERUM OR PLASMA**

### **1. INTENDED USE**

A solid-phase enzyme immunoassay for the qualitative determination of IgG antibodies to HSV 1,2 in blood serum or plasma.

This kit is designed for measurement of IgG antibodies to HSV 1,2 – in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 45 unknown samples in duplicates.

### **2. SUMMARY AND EXPLANATION**

Herpes simplex virus (HSV) is one of the most common pathogens in humans. HSV is transmitted by all secretions of infected body, especially via saliva, semen and cervical fluid. Acute HSV infections appear as a vesicular rash of labial or genital area. In immunocompromised hosts, HSV may cause life-threatening sequela in central nervous system. HSV is rarely fully eradicated after acute infections and persists in human organism lifelong, showing the periodic reactivation. In case of acute infection or reactivation during pregnancy, HSV may cross the placental barrier and cause severe fetal malformations.

Specific IgG-antibodies to HSV are not protective; their titer usually raise in response to the reactivation of virus and therefore may be used to monitor the actual status of HSV activity. In pregnant women, the absence of specific IgG-antibodies (seronegativity) requires tight restrictions of the lifestyle during the pregnancy minimizing contacts to seropositive humans. The seronegative individuals should not receive the blood transfusions and organ transplants from seropositive donors.

There are two very similar serotypes of HSV – HSV I and HSV II, showing different distribution in affected human tissues and organs. The present test system does not detect the differences between these two serotypes.

### **3. PRINCIPLE OF THE TEST**

This test is based on indirect enzyme immunoassay principle. Tested specimen is placed into the microwells coated by the antigen. Antibodies from the specimen bind coated antigen on the microwell surface. Unbound material is removed by washing procedure. Second antibodies directed towards species specific Ig, labelled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

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## 5.1. Contents of the Kit

## 5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	HSV 1,2 IgG EIA strips, 8x12 wells	1	pcs		until exp.date
2	Calibrator set, 1.5 ml	1	pcs	blue	2 months
3	Control sera (1.5 ml)	2	pcs	colourless and red	2 months
4	Conjugate, 14 ml	1	pcs	red	until exp.date
5	EIA buffer 50ml	1	pcs	blue	until exp.date
6	Substrate solution, 14 ml	1	pcs	colourless	until exp.date
7	Washing solution concentrate 26X, 22 ml	1	pcs	colourless	Concentrate – until exp.date Diluted washing solution – 45 days at 2-8 °C or 15 days at RT
8	Stop solution, 14 ml	1	pcs	colourless	until exp.date
9	Plate sealing tape	2	pcs		N/A
10	Instruction HSV 1,2 IgG EIA	1	pcs		N/A
11	QC data sheet HSV 1,2IgG EIA	1	pcs		N/A

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- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 5–250 µl;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0;
- Dry thermostat for +37 °C ±0.1 °C.

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#### 7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 6 wells for the calibrators CAL and control samples CONTROL - , CONTROL + and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Dilute all samples using buffer DIL SPE (EIA buffer) 101 fold. (5 µl of sample + 500 µl of diluent). Do not dilute control sample and calibrators.
3	Pipet 100 µl of calibrators, control samples CONTROL - , CONTROL + and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
4	Incubate 30 minutes at +37 °C.
5	Prepare washing solution by 26X dilution of washing solution concentrate BUF WASH 26X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
6	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
7	Incubate 30 minutes at +37 °C.
8	Wash the strips 5 times.
9	Dispense 100 µl of SUBS TMB into the wells.
10	Incubate 10-20 minutes at +18...+25 °C.
11	Dispense 100 µl of STOP into the wells.
12	Measure OD (optical density) at 450 nm.
13	Set photometer blank on air.
14	Apply coefficient method for data reduction. 1. Calculate a cutoff value by multiplying mean OD of Calibrator by the Q value (see Lot QC insert). Cutoff = OD(Calibrator) * Q value 2. Divide mean OD of each sample by Cutoff. $K = \text{OD}(\text{sample}) / \text{Cutoff}$ 3. If the K value is greater than 1.1, the result is POSITIVE. If the K value is less than 0.9, the result is NEGATIVE. If the K value is between 0.9 and 1.1, the result is EQUIVOCAL.

### 7.5. Sample processing

Material type	Notes on material collection, storage and handling	Sample dilution example	Sample into the well, $\mu$ l	Calculation factor
blood serum or plasma	grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing	5 $\mu$ l of sample + 500 $\mu$ l of diluent	100	1

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The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of *GLP* (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

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## 9. CALCULATION OF RESULTS

Some laboratories, based on their population studies, set up a second cutoff, which stands between anamnestic ('normal') IgG antibody level and 'high' IgG antibody level characteristic of reactivation or late period of primary infection. Recommended values for this second cut-off for two age groups (8 months – 3 year, > 3 years) are presented in the table below.

## 10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for HSV 1,2 IgG. Based on data obtained by XEMA, the following normal range is recommended (see below).

Sex, age	Units, K	
	Lower limit	Upper limit
Seronegative	<0.1	0.9
Seropositive > 3 years	1.1	4.9
newborn*	<0.1	1.5
under 8 months*	<0.1	3.5
8 months – 3 years	<0.1	4.8

\*antibodies of maternal origin

## 11. LITERATURE

1. Isakov V. A., Borisova V. V., Isakov D. V. Herpetic infection: pathogenesis and laboratory diagnostics. Manual for physicians. St.-Petersburg, Lanj publishing, 1999, p. 192.
2. Cowan F. M., Johnson A. M., Ashley R. et al. Relation ship between antibodies to herpes simplex virus (HSV) and symptoms of HSV infection // J. Infect. Dis. 1996- v.174 p.470-475
3. Munday P. E., Vuddamalay J., Slomka M. J. et al. Role of type specific herpes simplex virus serology in the diagnosis and management of genital herpes // Sex Transm. Inf 1998 v.74 p.175-178
4. Roest R. W., van der Meijden W. I., van Dijk G. et al. Prevalence and association between herpes simplex viruses types 1- and 2-specific antibodies in attendees at a sexually transmitted disease clinic. // Int Epidemiol. 2001 v.30 p.580-588
5. Whitley R/I/ Herpes simplex virusinfections. In: Infection Diseases of the Fetus and Newborn Infants (Remington JS, Klein JO eds.) Philadelphia: WB Saunders, 1990 - p.282-305

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Russian Diagnostic  
Manufacturers Association

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Russian Association  
of Medical Laboratory  
Diagnosticians

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## **«HSV 1,2 IgM-ИФА»**

**A SOLID-PHASE ENZYME IMMUNOASSAY  
FOR THE QUALITATIVE DETERMINATION  
OF IgM ANTIBODIES TO HSV 1/2  
IN HUMAN SERUM OR PLASMA**

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
**CONTENT**

1. INTENDED USE	11
2. SUMMARY AND EXPLANATION	11
3. PRINCIPLE OF THE TEST	11
4. WARNINGS AND PRECAUTIONS	12
5. KIT COMPONENTS	13
6. SPECIMEN COLLECTION AND STORAGE	14
7. TEST PROCEDURE	14
8. QUALITY CONTROL	15
9. CALCULATION OF RESULTS	15
10. EXPECTED VALUES	15
11. PERFORMANCE CHARACTERISTICS	16
12. LITERATURE	16

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*Instruction for use*

## **A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUALITATIVE DETERMINATION OF IGM ANTIBODIES TO HSV 1/2 IN HUMAN SERUM OR PLASMA**

### **1. INTENDED USE**

A solid-phase enzyme immunoassay for the qualitative determination of IgM antibodies to HSV 1/2 in serum or plasma.

This kit is designed for measurement of IgM antibodies to HSV 1/2 in serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 46 unknown samples in duplicates.

### **2. SUMMARY AND EXPLANATION**

Herpes simplex virus (HSV)infection is a chronic recurrent infectious disease caused by Herpes simplex viruses, types 1 and 2. The disease may be manifested as a local skin and/or mucosal vesicular rash or as a generalized form with multi-organ damage.

During HSV infection, a sequential synthesis of HSV-specific IgM, IgG and IgA occurs. IgM antibodies appear 10 or more days following infection, their level continuously rising during two weeks. HSV-specific IgG appear 7–10 days following appearance of IgM and persist lifelong.

A primary HSV infection is diagnosed by appearance of HSV-specific IgM and/or a 4-fold increase of HSV-specific IgG in sequential samples taken with 14-20 days interval. Appearance of HSV-specific IgM in an HSV-infected person indicates exacerbation of the disease.

### **3. PRINCIPLE OF THE TEST**

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by murine monoclonal antibodies to IgM. Antibodies from the specimen bind coated murine monoclonal antibodies to IgM on the microwell surface. Unbound material is removed by washing procedure. HSV 1,2, labelled with peroxidase enzyme, is then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Positivity index (PI, %) is calculated by the formula (see Calculations).Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

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## 5.1. Contents of the Kit

## 5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/ diluted components
1	HSV 1, 2 IgM EIA strips, 8x12 wells	1	pcs		until exp.date
2	Control sera (0.5 ml and 0.2 ml, resp.)	2	pcs	colourless; red	2 months
3	Conjugate concentrate, 1.25 ml	1	pcs	blue	Concentrate - until exp.date Diluted - 1 day at 2-8°C
4	Conjugate dilution buffer, 12 ml	1	pcs	bright yellow	until exp.date
5	EIA buffer 14 ml	1	pcs	blue	until exp.date
6	Substrate solution, 14 ml	1	pcs	colourless	until exp.date
7	Washing solution concentrate 26X, 22 ml	1	pcs	colourless	Concentrate - until exp.date Diluted washing solution - 45 days at 2-8 °C or 15 days at RT
8	Stop solution, 14 ml	1	pcs	colourless	until exp.date
9	Plate sealing tape	2	pcs		N/A
10	Instruction HSV 1, 2 IgM EIA	1	pcs		N/A
11	QC data sheet HSV 1, 2 IgM EIA	1	pcs		N/A

**5.2.** Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 90–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 10–250 µl;
- Dry thermostat or thermostat shaker for 37 °C ±2 °C;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

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**7. TEST PROCEDURE****7.1.** Reagent Preparation

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**7.4. Assay procedure**

1	Put the desired number of microstrips into the frame; allocate 4 wells for control samples CONTROL - CONTROL + (3 and 1 wells resp.) and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Pipet 90 µl of EIA buffer into each well.
3	Pipet 10 µl of control samples CONTROL - CONTROL + and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
4	Incubate 30 minutes at +37 °C and continuous shaking at 500-600 rpm.
5	Prepare washing solution by 26X dilution of washing solution concentrate BUF WASH 26X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
6	Dispense 100 µl of working conjugate solution into the wells. Cover the wells by plate adhesive tape.
7	Incubate 30 minutes at +37 °C and continuous shaking at 500-600 rpm.
8	Wash the strips 5 times.
9	Dispense 100 µl of SUBS TMB into the wells.
10	Incubate 10–20 minutes at +18...+25 °C.
11	Dispense 100 µl of STOP into the wells.
12	Measure OD (optical density) at 450 nm.
13	Set photometer blank on air.

Alternative incubation:

1	Put the desired number of microstrips into the frame; allocate 4 wells for control samples CONTROL - CONTROL + (3 and 1 wells resp.) and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Pipet 90 µl of EIA buffer into each well.
3	Pipet 10 µl of control samples CONTROL - CONTROL + and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
4	Incubate 60 minutes at +37 °C.
5	Prepare washing solution by 26X dilution of washing solution concentrate BUF WASH 26X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
6	Dispense 100 µl of working conjugate solution into the wells. Cover the wells by plate adhesive tape.
7	Incubate 30 minutes at +37 °C.

(Continuation see page 15)

8	Wash the strips 5 times.
9	Dispense 100 µl of SUBS TMB into the wells.
10	Incubate 10–20 minutes at +18...+25 °C.
11	Dispense 100 µl of STOP into the wells.
12	Measure OD (optical density) at 450 nm.
13	Set photometer blank on air.

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For the assay to be valid, the following requirements should be met:

1. OD450 for CONTROL+ should be nlt 1 AU.
2. OD450 for CONTROL- should not be more than 0.15 AU for all replicates.
3. OD450 for any CONTROL- replicate should be within 50%-150% of the mean OD450 value for CONTROL-. If any value lies outside this range (although meets requirement #2), it should be discarded and not used for calculation of the mean OD450 value for CONTROL-.

## 9. CALCULATION OF RESULTS

**9.1.** Calculate the mean absorbance values (OD450) for CONTROL- in triplicates and each pair of samples.

**9.2.** Calculate the cut-off value: (mean OD450 for CONTROL-) +0.2.

**9.3.** Calculate Positivity Index (PI) for each sample:

PI = mean OD450(sample) / Cut-off

## 10. EXPECTED VALUES AND ASSAY LIMITATIONS

If PI value is greater than 1.1, the result is POSITIVE.

If PI value is less than 0.9, the result is NEGATIVE.

If PI value is between 0.9 and 1.1, the result is EQUIVOCAL. Such samples should be retested. If the result is equivocal again, a new sample should be obtained 2-4 weeks later and tested again. If the result remains equivocal, the sample should be considered negative.

Therapeutic consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutic measures.

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## 11. PERFORMANCE CHARACTERISTICS

### 11.1. Analytical specificity

Specificity of the test was evaluated on 47 serum specimens found negative in DiaSorin (Italy) and NovaTec (Germany). All tested specimens were found negative. Based on these data, specificity of the test is 100%.

### 11.2. Analytical sensitivity

The kit HSV1,2 IgM EIA was evaluated using panel of naturally occurring serum and plasma specimens from Boston Biomedica, Inc ( Anti-HSV1/2 Mixed Titer Performance Panel PTH201).

### 11.3. Precision

Intra-assay precision is shown below:

Serum, no	duplicated	value, K	CV1, %	CV2, %
1	32	0.547	7.4	6.2
2	32	1.548	5.1	4.8

Inter-assay precision is shown below:

Serum, no	duplicated	value, K	CV, %
1	8	0.785	8.2
2	8	1.95	6.4

## 12. LITERATURE

1. Isakov V.A., Borisova V.V., Isakov D.V. Herpes infection: pathogenesis and laboratory diagnostics. A guideline for physicians – St.-Petersburg, Lanj - 1999 - p.192
2. Cowan F.M., Johnson A.M., Ashley R. et al. Relation ship between antibodies to herpes simplex virus (HSV) and symptoms of HSV infection // J. Infect. Dis. - 1996 - v.174-p.470-475
3. Munday P.E., Vuddamalay J., Slomka M.J. et al. Role of type specific herpes simplex virus serology in the diagnosis and management of genital herpes // Sex Transm. Inf - 1998-v.74 - p.175-178
4. Roest R.W., van der Meijden W.I., van Dijk G. et al. Prevalence and association between herpes simplex viruses types 1- and 2-specific antibodies in attendees at a sexually transmitted disease clinic. // Int Epidemiol. - 2001 - v.30 - p.580-588
5. Whitley R/I/ Herpes simplex virusinfections. In: Infection Diseases of the Fetus and Newborn Infants (Remington JS, Klein JO eds.) - Philadelphia: WB Saunders, 1990-p.282-305

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RUSSIAN ASSOCIATION  
OF MEDICAL LABORATORY  
DIAGNOSTICS

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*Instruction for use*

## **A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUALITATIVE DETERMINATION OF IgM ANTIBODIES TO *HELICOBACTER PYLORI* IN HUMAN SERUM OR PLASMA**

### **1. INTENDED USE**

A solid-phase enzyme immunoassay for the qualitative determination of IgM antibodies to *Helicobacter pylori* in serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 46 unknown samples in duplicates.

### **2. SUMMARY AND EXPLANATION**

*H. pylori* is detectable in nearly 100% of adult patients with duodenal ulcer, about 80% of patients with gastric ulcer, more than 90% of gastric cancer patients and more than 90% of patients with active chronic gastritis. An association between *H. pylori* and gastric cancer is confirmed. In developing countries, where prevalence of *H. pylori* infection is high (ca. 70-80% of population for Eastern Europe and Asia), most children become infected by the age of 10, and gastric cancer rates are very high. In the USA and other developed countries, standards of hygiene and the increasing socioeconomic status of the population have reduced the incidence of infection, and in parallel, the rates of peptic ulcers and gastric cancer have declined. There is an excellent correlation between the clinical presentation of gastritis, the presence of *H. pylori* in the stomach and elevated serum *H. pylori* IgG antibody. Moreover, eradication of *H. pylori* from a patient leads to a recovery from gastritis and a significant reduction in duodenal ulcer recurrence. *H. pylori*-specific IgG antibody also falls significantly after successful antibacterial therapy. *H. pylori* strains are classified into two broad groups - those that express both VacA (vacuolising toxin) and CagA (cytotoxin-associated gene) (type I) and those that produce neither (type II). Type I strains are predominate in patients with ulcers and cancer, while type II strains are mostly associated with asymptomatic infections. 80-100% of patients with duodenal ulcer disease produce CagA antibodies against a 128 kd antigen compared with 60-63% of *H. pylori*-infected persons with gastritis only, indicating that serologic responses to the 128 kd protein are more prevalent among *H. pylori*-infected persons with duodenal ulcers than infected persons without peptic ulceration. In *H. pylori*-infected patients who develop gastric cancer or duodenal ulcer, serum IgG against CagA are found in 94% and 80-100% of patients, resp., thus indicating that detection of antibodies to CagA is a useful marker for diagnosis of gastric cancer and duodenal ulcer. For gastritis, prevalence of IgG-antibodies to *H. pylori* is lower – ca. 60-63%.

The IgM test detects antibodies formed 2-4 weeks after infection with *H. pylori*. IgG and IgA antibodies require up to two months after infection before becoming detectable. If the immune system fails to eradicate the bacterium, it will begin to colonise the mucosal lining of the upper gastrointestinal tract. In the absence of *H. pylori* IgG, IgM positivity should be considered as a reason for treatment only in case of existing active symptoms of gastric ulcer. If treatment is not chosen in the presence of *H. pylori* IgM antibodies, it is advisable to carry out IgG testing after 2-3 months to confirm infection eradication by the immune system.

### 3. PRINCIPLE OF THE TEST

This test is based on indirect enzyme immunoassay principle. Tested specimen is placed into the microwells coated with CagA protein of *Helicobacter pylori*. Antibodies from the specimen bind coated CagA protein of *Helicobacter pylori* on the microwell surface. Unbound material is removed by washing procedure. Second antibodies directed towards species specific IgG, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Positivity index (PI, %) is calculated by the formula (see Calculations). Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

### 4. WARNINGS AND PRECAUTIONS

- 4.1.** For professional use only.
- 4.2.** This kit is intended for in vitro diagnostic use only.
- 4.3. INFECTION HAZARD:** There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.
- 4.4.** Avoid contact with stop solution containing 5.0% H<sub>2</sub>SO<sub>4</sub>. It may cause skin irritation and burns.
- 4.5.** Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.
- 4.6.** Do not use the kit beyond the expiration date.
- 4.7.** All indicated volumes have to be performed according to the protocol. Optimal test results are obtained only when using calibrated pipettes and microplate readers.
- 4.8.** Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
- 4.9.** Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.
- 4.10.** Do not mix reagents from different lots.
- 4.11.** Replace caps on reagents immediately. Do not swap caps.
- 4.12.** Do not pipette reagents by mouth.
- 4.13.** Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.
- 4.14.** Material Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.
- 4.15.** The Material Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

## 5. KIT COMPONENTS

### 5.1. Contents of the Kit

	Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	SORB MTP	<i>Helicobacter pylori</i> IgM EIA strips, 8x12 wells	1	pcs		until exp. date
2	CONTROL - CONTROL +	Control sera (0.5 ml and 0.2 ml, resp.)	2	pcs	colourless; red	2 months
3	CONJ HRP	Conjugate, 14 ml	1	pcs	green	until exp. date
4	DIL	EIA Sample buffer 14 ml	1	pcs	purple	until exp. date
5	SUBS TMB	Substrate solution, 14 ml	1	pcs	colourless	until exp. date
6	BUF WASH 26X	Washing solution concentrate 26x, 22 ml	1	pcs	colourless	Concentrate - until exp.date
7	STOP	Stop solution, 14 ml	1	pcs	Diluted washing solution - 1 month at 2-8 °C or 5 days at RT	until exp. date
8	N003	Plate sealing tape	2	pcs	colourless	until exp. date
9	K119MI	Instruction <i>Helicobacter pylori</i> IgM EIA	1	pcs	-	N/A
10	K119MQ	QC data sheet <i>Helicobacter pylori</i> IgM EIA	1	pcs	-	N/A

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- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 90-250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 10-250 µl;
- Dry thermostat for 37 °C ±2 °C;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0-3.0.

**5.3.** Storage and stability of the Kit

Store the whole kit at 2 to 8 °C upon receipt until the expiration date.

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Specimens may be stored for up to 48 hours at 2-8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

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- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18 to +25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the vials with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 26X by 26 dilution in distilled water.

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**7.4.** Alternative units:

**7.5.** Assay procedure:

1	Put the desired number of microstrips into the frame; allocate 4 wells for control samples CONTROL - CONTROL + (3 and 1 wells resp.) and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Pipet 90 µl of EIA Sample buffer into each well.
3	Pipet 10 µl of control samples CONTROL - CONTROL + and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
4	Incubate 30 minutes at 37 °C.
5	Prepare washing solution by 26x dilution of washing solution concentrate BUF WASH 26X with distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
6	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
7	Incubate 30 minutes at 37 °C.
8	Wash the strips 5 times.
9	Dispense 100 µl of SUBS TMB into the wells.
10	Incubate 10-20 minutes at 18-25 °C.
11	Dispense 100 µl of STOP into the wells.
12	Measure OD (optical density) at 450 nm.
13	Set photometer blank on air.

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1. OD450 for CONTROL+ should be nit 0.6 AU.
2. OD450 for CONTROL- should not be more than 0.15 AU for all replicates.

## 8. EXPECTED VALUES

### CALCULATION OF RESULTS.

1. Calculate the mean absorbance values (OD450) for CONTROL- in triplicates and each pair of samples.
2. Calculate the cut-off value: (mean OD450 for CONTROL-) + 0.25
3. Calculate Positivity Index (PI) for each sample:  

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If PI value is between 0.9 and 1.1, the result is EQUIVOCAL. Such samples should be retested. If the result is equivocal again, a new sample should be obtained 2-4 weeks later and tested again. If the result remains equivocal, the sample should be considered negative. Therapeutic consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutic measures. NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

## 11. PERFORMANCE CHARACTERISTICS

### 11.1 Analytical specificity.

Specificity of the test was evaluated on 93 negative serum specimens. Based on these data, specificity of the test is 98%.

### 11.2 Analytical sensitivity.

Specificity of kit was evaluated using internal panel "H.Pylori Positive Control".

### 11.3 Precision.

Intra-assay precision for two different lots (CV1, CV2) is shown below:

Serum, no	Replicates	Value, IP average	CV1, %	CV2, %
1	32	3.4	2.9	2.8
2	32	9.1	4.8	5

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*Instruction for use***A SOLID-PHASE ENZYME IMMUNOASSAY  
FOR THE QUANTITATIVE DETERMINATION OF HCG  
IN HUMAN SERUM OR PLASMA****1. INTENDED USE**

A solid-phase enzyme immunoassay for the quantitative determination of hCG in blood serum or plasma.

This kit is designed for measurement of hCG in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 41 unknown samples in duplicates.

**2. SUMMARY AND EXPLANATION**

Human chorionic gonadotropin (HCG) is a glycoprotein hormone secreted by trophoblastic cells of placenta. A molecule of HCG consists of two noncovalently bound subunits: alpha- and beta-HCG. Beta-subunit is specific for HCG hormone.

Determination of HCG is widely used for early diagnosis of pregnancy. Multiple pregnancy results in correspondent elevation of serum HCG; while ectopic pregnancy and placental insufficiency cause decreased serum HCG levels. Determination of HCG in serum during second trimester is used for pregnancy monitoring, especially in screening for Down syndrome, along with other laboratory tests (AFP and Estriol).

Serum HCG is also a laboratory marker of trophoblastic tumours – chorionepitheliomas, some seminomas and teratomas. Serial determination of serum HCG can be used for therapy monitoring in these cancers.

The present test system uses beta ( $\beta$ )-chain specific monoclonal antibody as the capture, and alpha ( $\alpha$ )-chain specific monoclonal antibody as the tracer; therefore only the whole intact HCG molecule is detected.

### 3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to  $\beta$  chain of hCG-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Second antibodies – murine monoclonal to  $\alpha$  chain of hCG, labelled with peroxidase enzyme, are then added into the microwells. After washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

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## 5.1. Contents of the Kit

## 5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	hCG EIA strips, 8x12 wells	1	pcs		until exp.date
2	Calibrator set, 0.8 ml each. The set contains 6 calibrators: 0; 15; 60; 125; 250, 500 IU/l	6	pcs	blue (C1 - colourless)	2 months
3	Control serum (0.8 ml)	1	pcs	colourless	2 months
4	Conjugate, 14 ml	1	pcs	red	until exp.date
5	EIA buffer 100 ml	1	pcs	blue	until exp.date
6	Substrate solution, 14 ml	1	pcs	colourless	until exp.date
7	Washing solution concentrate 26X, 22 ml	1	pcs	colourless	Concentrate - until exp.date Diluted washing solution - 45 days at 2-8 °C or 15 days at RT
8	Stop solution, 14 ml	1	pcs	colourless	until exp.date
9	Plate sealing tape	2	pcs		N/A
10	Instruction hCG EIA	1	pcs		N/A
11	QC data sheet hCG EIA	1	pcs		N/A

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- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 25–250 µl;
- Dry thermostat for +37 °C ±0.1 °C;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

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**7.2. Procedural Note:** It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

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**7.4. Alternative units:** 1 IU/l = 0.2 ng/ml

**7.5. Highdose hook effect.** A high-dose hook effect is possible in kits of this nature. To avoid falsely low results, samples obtained from pregnant women should be diluted accordingly before assaying:

Week 1-2 of gestation – no dilution

Week 3-4 – 1:20

Week 4-14 – 1:400

Week 14-21 – 1:400 and 1:1000 (both dilutions should be assayed)

Week 22, 3rd trimester – 1:400.

Dilution examples:

1:20 – 25 ul of sample + 475 ul of EIA buffer

1:400 – 25 ul of 1:20 dilution + 475 ul of EIA buffer

1:1000 – 200 ul of 1:400 dilution + 300 ul of EIA buffer.

If pregnancy status is not clear but suspected, the sample should be assayed both undiluted and diluted 1:20.

### 7.6. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 14 wells for the calibrators CAL 1–6 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using DIL SPE (EIA buffer) (see p 7.5). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
3	Dispense 100 µl of CONJ HRP into the wells.
4	Pipet 50 µl of calibrators CAL 1–6, control samples CONTROL and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
5	Incubate 60 minutes at 37 °C.
6	Prepare washing solution by 26X dilution of washing solution concentrate (BUF WASH 26X) by distilled water. Wash the strips 5 times.
7	Dispense 100 µl of SUBS TMB into the wells
8	Incubate 10–20 minutes at +18...+25 °C
9	Dispense 100 µl of STOP into the wells.
10	Measure OD (optical density) at 450 nm.
11	Set photometer blank on first calibrator
12	Apply point-by-point method for data reduction.

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
**9.2.** Plot a calibration curve on graph paper: OD versus hCG concentration.

**9.3.** Determine the corresponding concentration of hCG in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

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## 10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for hCG. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units, IU/l		Units alternative, ng/ml	
	Lower limit	Upper limit	Lower limit	Upper limit
Males	-	15	-	3.0
Females	-	15	-	3.0
Pregnancy week:				
1	-	50	-	10.0
2	20	500	4.0	100
3	500	5000	100	1000
4	3000	19000	600	3800
5-8	14000	169000	2800	33800
9-13	16000	180000	3200	36000
22	4500	70000	900	14000
23	3000	69500	600	13900
3rd trimester	2400	50000	480	10000

Pregnancy, week	Median, kIU/l	SKO
14	55	0.83
15	42.8	0.78
16	38.4	0.74
17	32.8	0.75
18	25.6	0.69
19	23	0.56
20	21	0.52
21	18	0.51

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### 11.2. Analytical sensitivity

Sensitivity of the assay was assessed as being 1.25 IU/l.

### 11.3. Linearity



Linearity was checked by assaying dilution series of 5 samples with different hCG concentrations. Linearity percentages obtained ranged within 90 to 110%.

### 11.4. Recovery

Recovery was estimated by assaying 5 mixed samples with known hCG concentrations. The recovery percentages ranged from 90 to 110%.

## 12. LITERATURE

1. Krieg, A.F. Pregnancy Tests and Evaluation of Placental Function. Clinical Diagnosis and Management by Laboratory Methods, 16th ed. Henry, J.B. editor, W.B. Saunders Co. Philadelphia. app. 680, 1979.
2. Brody, S., and Carlstrom, G., Immunoassay of Human Chorionic Gonadotropin in Normal and Pathologic Pregnancy. J. Clin. Endocrinol. Metab., 22: 564, 1962.
3. Braunstein, G.D., Rasor, J., Adleer, D., et al: Serum Human Chorionic Gonadotropin Levels Throughout Normal Pregnancy. Am. J. Obstet. Gynecol., 126: 678, 1976.
4. Skinner, M.S., Seckinger, D., Evaluation of Beta-Subunit Chorionic Gonadotropin as an Aid in Diagnosis of Trophoblastic Disease. Ann. Clin. Lab. Sci., 9 (4): 347-52, 1979.
5. Braunstein, G.D., Vaitukaitis, J.L., Carbone, P.P., and Ross, G.T., Ectopic Production of Human Chorionic Gonadotropin by Neoplasms. Ann. Intern. Med., 78: 39-45, 1973.

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1. Hall, P. F., Testicular Steroid Synthesis: Organization and Regulation. In: The Physiology of Reproduction, Ed.: Knobil, E., and Neill, J. et al., pp 975-98. Raven Press, New York (1988).
2. Tsang, B. K., Armstrong, D. T. and Whitfield, J. F., Steroid biosyntheses by isolated human ovarian follicular cells in vitro, J. Clin. Endocrinol. Metab. 51:1407 - 11 (1980).
3. Gore-Langton, R. E. and Armstrong, D. T., Follicular steroidogenesis and its control. In: The physiology of Reproduction, Ed.: Knobil, E., and Neill, J. et al., pp. 331-85. Raven Press, New York (1988).
4. Siiteri, P. K. Murai, J. T., Hammond, G. L., Nisker, J. A., Raymoure, W. J. and Kuhn, R. W., The serum transport of steroid hormones, Rec. Prog. Horm. Res. 38:457- 510 (1982).
5. Martin, B., Rotten, D., Jolivet, A. and Gautray, J-P-. Binding of steroids by proteins in follicular fluid of the human ovary, J.Clin. Endocrinol. Metab. 35: 443-47 (1981).
6. Baird, D. T., Ovarian steroid secretion and metabolism in women. In: The Endocrine Function of the Human Ovary. Eds.: James, V.H:T., Serio, M. and Giusti, G. pp. 125-33, Academic Press, New York (1976).
7. McNastty, K. P., Baird, D. T., Bolton, a., Chambers, P., Corker, C.S. and McLean, H., concentration of oestrogens and androgens in human ovarian venous plasma and follicular fluid throughout the menstrual cycle, J. Endocrinol. 71:77-85 (1976).

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*Instruction for use*

## **A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF ESTRADIOL IN HUMAN SERUM OR PLASMA**

### **1. INTENDED USE**

A solid-phase enzyme immunoassay for the quantitative determination of estradiol in blood serum or plasma.

This kit is designed for measurement of estradiol in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 41 unknown samples in duplicates.

### **2. SUMMARY AND EXPLANATION**

Estradiol (E2) is a steroid hormone with molecular weight 272.4 Da. In humans, E2 shows the highest physiological activity among the estrogens. In males, minor quantities of E2 are produced by adrenals and testicles. In females, E2 is produced by ovarian follicles. The physiological activity of E2 involves multiple functions in female reproductive system. Regulation of E2 production and secretion is complex and depends on pituitary and ovarian hormones.

Serum E2 level is low in follicular phase of the menstrual cycle; 3–5 days before serum LH peak, serum E2 level begins to rise and reaches a maximum ca. 12 hours before LH peak. After LH peak, E2 level drops dramatically to the minimal level and starts to rise again. The maximal E2 level in serum is observed in luteal phase, at day 9 after ovulation; then the decline of serum E2 reflects the degradation of corpus luteum.

During pregnancy, the determination of serum E2 reflects the status of fetoplacental system. In first trimester, serum E2 level is in the range corresponding to the ovulation levels. Sharp increase of serum E2 in pregnant women is observed between 9th and 10th week; then the increase continues less sharply by the end of pregnancy.

Increased levels of serum estradiol are characteristic for metrorrhagias in post-menopausal age; adrenal hyperplasia; estrogen-secreting tumours; liver cirrhosis; feminization in children and males; intake of gonadotropins and estrogens.

Decreased levels of serum estradiol are observed in Turner syndrome, primary or secondary hypogonadism; germaphroditism; post-climacteric syndrome; fetal dysfunctions; intake of oral contraceptives.

### 3. PRINCIPLE OF THE TEST

This test is based on competition enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific rabbit polyclonal to estradiol-antibodies simultaneously with conjugated Estradiol-peroxidase. Estradiol from the specimen competes with the conjugated Estradiol for coating antibodies. After washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is inversely related to the quantity of the measured analyte in the specimen.

### 4. WARNINGS AND PRECAUTIONS

**4.1.** For professional use only.

**4.2.** This kit is intended for in vitro diagnostic use only.

**4.3. INFECTION HAZARD:** There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

**4.4.** Avoid contact with stop solution containing 5.0% H<sub>2</sub>SO<sub>4</sub>. It may cause skin irritation and burns.

**4.5.** Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

**4.6.** Do not use the kit beyond the expiration date.

**4.7.** All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

**4.8.** Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

**4.9.** Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

**4.10.** Do not mix reagents from different lots.

**4.11.** Replace caps on reagents immediately. Do not swap caps.

**4.12.** Do not pipette reagents by mouth.

**4.13.** Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

**4.14.** Material Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

**4.15.** The Material Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

## 5.1. Contents of the Kit

## 5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1 SORB MTP	Estradiol EIA strips, 8x12 wells	1	pcs		until exp.date
2 CAL 1-6	Calibrator set, 0.8 ml each. The set contains 6 calibrators: 0: 0.1; 0.3; 1; 3; 20 nmol/l	6	pcs	bright red (C1 – colourless)	2 months
3 CONTROL	Control serum (0.8 ml)	1	pcs	colourless	2 months
4 CONJ HRP	Conjugate, 14 ml	1	pcs	red	until exp.date
5 SUBS TMB	Substrate solution, 14 ml	1	pcs	colourless	until exp.date
6 BUF WASH 26X	Washing solution concentrate 26X, 22 ml	1	pcs	colourless	Concentrate – until exp.date Diluted washing solution – 45 days at 2-8 °C or 15 days at RT
7 STOP	Stop solution, 14 ml	1	pcs	colourless	until exp.date
8 N003	Plate sealing tape	2	pcs		N/A
9 K2081	Instruction Estradiol EIA	1	pcs		N/A
10 K208Q	QC data sheet Estradiol EIA	1	pcs		N/A

**5.2. Equipment and material required but not provided**

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 25–250 µl;
- Dry thermostat for 37 °C ±0.1 °C (optional: with shaking 600 rpm);
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0.

**5.3. Storage and stability of the Kit**

Store the whole kit at +2...+8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED) to minimize exposure to moisture.

**6. SPECIMEN COLLECTION AND STORAGE**

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at +2...+8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

**7. TEST PROCEDURE****7.1. Reagent Preparation**

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 26X by 26 dilutions in distilled water.

**7.2. Procedural Note:**

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

See the example of calibration graphic in Quality Control data sheet.

**7.3. Assay flowchart**

1 nmol/l = 272 pg/ml

#### 7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 14 wells for the calibrators CAL 1–6 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Pipet 25 µl of calibrators CAL 1–6, control samples CONTROL and unknown samples into the wells.
3	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape (included into the kit).
4	Incubate 120 minutes at 37 °C. OR: Incubate at +37 °C with continuous shaking (ca 600 rpm) during 60 min.
5	Prepare washing solution by 26X dilution of washing solution concentrate (BUF WASH 26X) by distilled water. Wash the strips 5 times.
6	Dispense 100 µl of SUBS TMB into the wells.
7	Incubate 10–20 minutes at +18...+25 °C.
8	Dispense 100 µl of STOP into the wells.
9	Measure OD (optical density) at 450 nm.
10	Set photometer blank on air.
11	Apply lin-log method for data reduction.

#### 7.5. Handling notes

Calibrators and control sample(s) – only one freezing/thawing cycle is allowed

Incubation at +37 °C may be replaced by incubation at +37 °C with continuous shaking (ca 600 rpm) during 60 min

#### 8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

## 9. CALCULATION OF RESULTS


1. Calculate the mean absorbance values (OD450) for each pair of calibrators and samples.

2. Plot a calibration curve on graph paper: OD versus estradiol concentration.

3. Determine the corresponding concentration of estradiol in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

4. Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

Calibrators	Value	Absorbance Units (450 nm) x 1000
CAL 1	0 nmol/l	2579
CAL 2	0.1 nmol/l	2408
CAL 3	0.3 nmol/l	1766
CAL 4	1 nmol/l	1300
CAL 5	3 nmol/l	934
CAL 6	20 nmol/l	352

## 10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for Estradiol. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units, nmol/l		Units alternative, pg/ml	
	Lower limit	Upper limit	Lower limit	Upper limit
Children under 11 yrs	-	0.2	-	54.4
Males	0.029	0.3	7.9	81.6
Females				
Pregnancy week:				
1st trimester	0.1	10.5	27	2856
2nd trimester	3.0	21	816	5712
3rd trimester	6.0	80	1632	21760
Menstrual cycle:				
follicular phase	0.05	0.70	13.6	190.4
luteinic phase	0.1	1.10	27.2	299.2
ovulation	0.34	1.80	92.5	489.6
post menopausal	-	0.23	-	62.6

## 11. PERFORMANCE CHARACTERISTICS

### 11.1. Analytical specificity / Cross reactivity

Analyte	Cross-reactivity, % wt/wt
Estradiol	100
Estrone	0.2
Estriol	0.6
Cortisol	0.06
Prednisolone	0.09
Corticosterone	<0.01
Progesterone	<0.01
17-Hydroxyprogesterone	<0.05
Pregnenolone	<0.05
Testosterone	<0.01

### 11.2. Analytical sensitivity.

Sensitivity of the assay was assessed as being 0.025 nmol/l.

### 11.3. Linearity.

Linearity was checked by assaying dilution series of 5 samples with different estradiol concentrations. Linearity percentages obtained ranged within 90 to 110%.

### 11.4. Recovery.

Recovery was estimated by assaying 5 mixed samples with known estradiol concentrations. The recovery percentages ranged from 90 to 110%.

## 12. LITERATURE

- Hall, P. F., Testicular Steroid Synthesis: Organization and Regulation. In: The Physiology of Reproduction, Ed.: Knobil, E., and Neill, J. et al., pp 975-98. Raven Press, New York (1988).
- Tsang, B. K., Armstrong, D. T. and Whitfield, J. F., Steroid biosyntheses by isolated human ovarian follicular cells in vitro, J. Clin. Endocrinol. Metab. 51:1407 – 11 (1980).
- Gore-Langton, R. E. and Armstrong, D. T., Follicular steroidogenesis and its control. In: The physiology of Reproduction, Ed.: Knobil, E., and Neill, J. et al., pp. 331-85. Raven Press, New York (1988).
- Siiteri, P. K. Murai, J. T., Hammond, G. L., Nisker, J. A., Raymoure, W. J. and Kuhn, R. W., The serum transport of steroid hormones, Rec. Prog. Horm. Res. 38:457 - 510 (1982).
- Martin, B., Rotten, D., Jolivet, A. and Gautray, J-P-. Binding of steroids by proteins in follicular fluid of the human ovary, J.Clin. Endocrinol. Metab. 35: 443-47 (1981).
- Baird, D. T., Ovarian steroid secretion and metabolism in women. In: The Endocrine Function of the Human Ovary. Eds.: James, V.H:T., Serio, M. and Giusti, G. pp. 125-33, Academic Press, New York (1976).
- McNastty, K. P., Baird, D. T., Bolton, a., Chambers, P., Corker, C. S. and McLean, H., concentration of oestrogens and androgens in human ovarian venous plasma and follicular fluid throughout the menstrual cycle, J. Endocrinol. 71:77-85 (1976).

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*Instruction for use*

## **A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF THYROGLOBULIN IN HUMAN SERUM OR PLASMA**

### **1. INTENDED USE**

A solid-phase enzyme immunoassay for the quantitative determination of thyroglobulin in blood serum or plasma.

This kit is designed for measurement of thyroglobulin in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 42 unknown samples in duplicates.

### **2. SUMMARY AND EXPLANATION**

Thyroglobulin (TG) is a high MW (ca. 650-700 kDa) glycoprotein synthesized by the thyroid epithelial cells. In normal thyroid gland, TG is secreted to the follicular lumen and undergoes iodination of tyrosine residues leading to formation of thyroid hormones (T3 and T4). Minor quantities of TG penetrate to the circulation in normal donors. Synthesis of TG is regulated by hormones (TSH, TRH, exogenous thyroid therapy).

In differentiated thyroid carcinoma, serial determination of serum TG is used for post-treatment monitoring. An elevation of serum TG in such patients indicates a presence of residual thyroid tissue, relapse or metastatic growth of the tumour.

The elevated serums TG are also observed in benign thyroid diseases, e.g. thyroiditis, hyperthyroidism and non-toxic goiter. The monitoring of serum TG is also used for prognostic evaluation of thyrostatic treatment of Graves' disease.

### **3. PRINCIPLE OF THE TEST**

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to human TG-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Second antibodies – murine monoclonal to human TG, labelled with peroxidase enzyme, are then added into the microwells. After washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

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## 5.1. Contents of the Kit

## 5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1 SORB MTP	Thyroglobulin EIA strips, 8x12 wells	1	pcs		until exp.date
2 CAL 1-5	Calibrator set, 0.8 ml each. The set contains 5 calibrators: 0; 10; 25; 100; 400 ng/ml	5	pcs	bright blue (C1 – colourless)	2 months
3 CONTROL	Control serum (0.8 ml)	1	pcs	colourless	2 months
4 CONJ HRP	Conjugate, 7 ml	1	pcs	bright blue	until exp.date
5 SUBS TMB	Substrate solution, 14 ml	1	pcs	colourless	until exp.date
6 BUF WASH 26X	Washing solution concentrate 26X, 22 ml	1	pcs	colourless	Concentrate – until exp.date Diluted washing solution – 45 days at 2-8 °C or 15 days at RT
7 STOP	Stop solution, 14 ml	1	pcs	colourless	until exp.date
8 N003	Plate sealing tape	2	pcs		N/A
9 K232I	Instruction Thyroglobulin EIA	1	pcs		N/A
10 K232Q	QC data sheet Thyroglobulin EIA	1	pcs		N/A

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- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 25–250 µl;
- Dry thermostat for 37 °C ±0.1 °C;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0.

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**7.4. Assay procedure**

1	Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1–5 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Pipet 50 µl of calibrators CAL 1–5, control samples CONTROL and unknown samples into the wells.
3	Dispense 50 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape (included into the kit).
4	Incubate 60 minutes at 37 °C.
5	Prepare washing solution by 26X dilution of washing solution concentrate (BUF WASH 26X) by distilled water. Wash the strips 5 times.
6	Dispense 100 µl of SUBS TMB into the wells.
7	Incubate 10–20 minutes at +18...+25 °C.
8	Dispense 100 µl of STOP into the wells.
9	Measure OD (optical density) at 450 nm.
10	Set photometer blank on first calibrator.
11	Apply point-by-point method for data reduction.

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
**9.2.** Plot a calibration curve on graph paper: OD versus thyroglobulin concentration.

**9.3.** Determine the corresponding concentration of thyroglobulin in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

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## 10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for Thyroglobulin. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units, ng/ml	
	Lower limit	Upper limit
Healthy donors	–	50

## 11. PERFORMANCE CHARACTERISTICS

### 11.1. Analytical sensitivity.

Sensitivity of the assay was assessed as being 0.5 ng/ml.

### 11.2. Linearity.

Linearity was checked by assaying dilution series of 5 samples with different thyroglobulin concentrations. Linearity percentages obtained ranged within 90 to 110%.

### 11.3. Recovery.

Recovery was estimated by assaying 5 mixed samples with known thyroglobulin concentrations. The recovery percentages ranged from 90 to 110%.

## 12. LITERATURE

1. CA Spencer, M Takeuchi, and M Kazarosyan – Current status and performance goals for serum thyroglobulin assays. Clin. Chem., Jan 1996; 42: 164 – 173.
2. PW Ladenson – Optimal laboratory testing for diagnosis and monitoring of thyroid nodules, goiter, and thyroid cancer. Clin. Chem., Jan 1996; 42: 183 – 187.
3. Anthony P. Weetman – Graves' Disease. N. Engl. J. Med., Oct 2000; 343: 1236 – 1248.

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1. DJ Byrne, MC Browning, and A Cuschieri – CA72-4: a new tumour marker for gastric cancer. Br J Surg, Sep 1990; 77(9): 1010-3.
2. Ian J. Jacobs and Usha Menon – Progress and Challenges in Screening for Early Detection of Ovarian Cancer. Mol. Cell. Proteomics, Apr 2004; 3: 355 – 366.
3. R Hamazoe, M Maeta, T Matsui, S Shibata, S Shiota, and N Kaibara – CA72-4 compared with carcinoembryonic antigen as a tumour marker for gastric cancer. Eur J Cancer, Jan 1992; 28A(8-9): 1351-4.

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*Instruction for use*

## **A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF CA72-4 IN HUMAN SERUM OR PLASMA**

### **1. INTENDED USE**

A solid-phase enzyme immunoassay for the quantitative determination of CA72-4 in blood serum or plasma.

This kit is designed for measurement of CA72-4 in blood serum or plasma. The kit contains reagents sufficient for 96 determinations and allows to analyze 42 unknown samples in duplicates.

### **2. SUMMARY AND EXPLANATION**

CA 72-4, or a carbohydrate antigen 72-4, is a high MW (230-1000 kD) antigen (epitope) associated to gastric and ovarian cancer as well as some other malignancies and not expressed in noticeable quantities in tissues of healthy adult individuals. Quantitative determination of CA 72-4 in serum or plasma is helpful (particularly, in combination with CA 19-9 – see Xema, Cat.# K223) for monitoring of gastric cancer and its therapy, while combined determination of CA 72-4 and CA 125 (see Xema, Cat.# K222) is used for monitoring of ovarian cancer.

Elevated levels of CA 72-4 are often seen in adenocarcinomas of the gastro-intestinal tract, ovaries (mucinous type) and lungs. Besides, raised CA 72-4 is sometimes also seen in patients with benign pathology (chronic inflammation, cysts, fibrosis). That is why, results of CA 72-4 determination should always be interpreted in conjunction with other laboratory and clinical data.

NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy may develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies gives false results in the present assay. Sera from HAMA positive patients should be treated by depleting adsorbents before assaying.

### **3. PRINCIPLE OF THE TEST**

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal antibodies to human CA72-4. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Unbound material is removed by washing procedure. Second antibodies – murine monoclonal Ab to human CA72-4, labelled with peroxidase, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the concentration of the measured analyte in the specimen.

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### 5. KIT COMPONENTS

#### 5.1. Contents of the Kit

Symbol	Description	Qty	Units	Colour code	Stability of opened/ diluted components
1	SORB MTP CA72-4 EIA strips, 8x12 wells	polystyrene microwells coated with murine monoclonal Ab to human CA72-4	1	pcs	until exp.date
2	CAL 1-5 Calibrator set, 0.5 ml each. The set contains 5 calibrators: 0; 5; 15; 50; 200 U/ml	human CA72-4 diluted in a preselected human serum preservative - 0.01% Bronidox L, 0.01% 2-Methyl-4-isothiazolin-3-one-hydrochloride; also contains blue dye. Ready-to-use	5	pcs	blue (C1 - colourless) 2 months
3	CONTROL Control serum (0.5 ml)	dilution of preselected human serum, with high content of CA72-4 with preservative - 0.01% Bronidox L, 0.01% 2-Methyl-4-isothiazolin-3-one-hydrochloride; also contains purple dye. Ready-to-use	1	pcs	colourless 2 months
4	CONJ 11X Conjugate concentrate, 1.2 ml	Aqueous solution of murine monoclonal to human CA72-4 coupled with horseradish peroxidase diluted on phosphate buffered saline with casein from bovine milk and detergent (Tween-20), contains 0,1% phenol as preservative and green dye	1	pcs	green Concentrate - until exp.date Diluted - 1 day at 2-8 °C
5	DIL CONJ Conjugate dilution buffer, 12 ml	Aqueous Tris-buffered BSA solution, preservative - 0,01% Bronidox L, 0,01% 2-Methyl-4-isothiazolin-3-one-hydrochloride; contains light blue dye	1	pcs	light blue until exp.date
6	SUBS TMB Substrate solution, 14 ml	ready-to-use single-component tetramethylbenzidine (TMB) solution	1	pcs	colourless until exp.date
7	BUF WASH 26X Washing solution concentrate 26x, 22 ml	aqueous solution of sodium chloride and detergent (Tween 20), contains proClin300 as a preservative.	1	pcs	colourless Concentrate - until exp.date Diluted washing solution - 45 days at 2-8°C or 15 days at RT
8	STOP Stop solution, 14 ml	5,0% vol/vol solution of sulphuric acid. Ready-to-use	1	pcs	colourless until exp.date
9	N003 Plate sealing tape		2	pcs	N/A
10	K244I Instruction CA72-4 EIA		1	pcs	N/A
11	K244Q QC data sheet CA72-4 EIA		1	pcs	N/A

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- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 20–250 µl;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

**5.3. Storage and stability of the Kit**

Store the whole kit at 2 to 8 °C upon receipt until the expiration date. A short-term storage at higher temperatures (nmt 25 °C for nmt 5 days) is allowed (e.g., during transportation).

After opening the pouch keep unused microtiter wells **TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED)** to minimize exposure to moisture.

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Specimens should be tested during the day of sampling. If test run is scheduled for another day, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

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- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18 to +25 °C) before use.
- Prepare washing solution from the concentrate BUF WASH 26X by 26 dilutions in distilled water.

## 7.2. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1-5 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Prepare working conjugate solution by dilution of conjugate concentrate 11 fold by conjugate dilution buffer. ATTENTION: working conjugate solution is unstable and should not be stored! Prepare the volume required for actual assay run.
3	Pipet 20 µl of calibrators CAL 1-5, control samples CONTROL and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
4	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
5	Incubate 60 minutes at 37 °C.
6	Prepare washing solution by 26x dilution of washing solution concentrate BUF WASH 26X by distilled water. Wash the strips 5 times.
7	Dispense 100 µl of SUBS TMB into the wells. Dispensing of TMB should be completed within 3 minutes.
8	Incubate in a dark place for 10-20 minutes, depending on blue coloration of wells.
9	Dispense 100 µl of STOP into the wells in the same order as TMB.
10	Measure OD (optical density) at 450 nm.
11	Set photometer blank on first calibrator.
12	Apply point-by-point method for data reduction.

## 8. QUALITY CONTROL

It is recommended to use external control samples (not included) according to state and federal regulations. The use of external control samples is advised to assure the day to day validity of results.

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
**9.2.** Plot a calibration curve on graph paper: OD versus CA72-4 concentration.

**9.3.** Determine the corresponding concentration of CA72-4 in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

**9.4.** is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

Calibrators	Value	Absorbance Units (450 nm)
CAL 1	0 U/ml	0.08
CAL 2	5 U/ml	0.17
CAL 3	15 U/ml	0.32
CAL 4	50 U/ml	0.67
CAL 5	200 U/ml	2.13

## 10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutical measures. Each laboratory should establish its own normal range for CA72-4. Based on data obtained by XEMA, the following normal range is recommended (see below).

NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

NOTE2: if analyte concentration in a specimen is lower than the analytical sensitivity (2 U/ml) or higher than the highest calibrator (200 U/ml), it is recommended to express the results in the following manner: "CA 72-4 concentration is lower than 1.0 U/ml" or "CA 72-4 concentration is higher than 200 U/ml".

Patient group	Units, U/ml	
	Lower limit	Upper limit
Healthy individuals	-	6

## 11. PERFORMANCE CHARACTERISTICS

### 11.1. Analytical sensitivity

Sensitivity of the assay was assessed as being 1.0 U/ml.

### 11.2. Linearity

Linearity was checked by assaying dilution series of 5 samples with different CA72-4 concentrations. Linearity percentages obtained within the range of 5-200 U/ml fell within 90–110%.

### 11.3. Recovery

Recovery was estimated by assaying a mixed sample (Contol + calibrator 3, 1:1) with known CA72-4 concentration. The recovery percentages ranged from 90% to 110%.

11.4. High-dose hook effect was not seen up to 50 000 IU/ml.

## 12. LITERATURE

1. DJ Byrne, MC Browning, and A Cuschieri - CA72-4: a new tumour marker for gastric cancer. *Br J Surg*, Sep 1990; 77(9): 1010-3.
2. Ian J. Jacobs and Usha Menon - Progress and Challenges in Screening for Early Detection of Ovarian Cancer. *Mol. Cell. Proteomics*, Apr 2004; 3: 355 - 366.
3. R Hamazoe, M Maeta, T Matsui, S Shibata, S Shiota, and N Kaibara - CA72-4 compared with carcinoembryonic antigen as a tumour marker for gastric cancer. *Eur J Cancer*, Jan 1992; 28A(8-9): 1351-4.

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*Instruction for use*

## **A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF TOTAL PSA IN HUMAN SERUM OR PLASMA**

### **1. INTENDED USE**

A solid-phase enzyme immunoassay for the quantitative determination of C-peptide in blood serum or plasma.

This kit is designed for measurement of C-peptide in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 41 unknown samples in duplicates.

### **2. PRINCIPLE OF THE TEST**

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to human C-peptide-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Second antibodies - murine monoclonal to human C-peptide, labelled with peroxidase enzyme, are then added into the microwells. After washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

### 3. WARNINGS AND PRECAUTIONS

**3.1.** For professional use only.

**3.2.** This kit is intended for in vitro diagnostic use only.

**3.3. INFECTION HAZARD:** There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

**3.4.** Avoid contact with stop solution containing 5.0 % H<sub>2</sub>SO<sub>3</sub>. It may cause skin irritation and burns.

**3.5.** Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

**3.6.** Do not use the kit beyond the expiration date.

**3.7.** All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

**3.8.** Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

**3.9.** Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

**3.10.** Do not mix reagents from different lots.

**3.11.** Replace caps on reagents immediately. Do not swap caps.

**3.12.** Do not pipette reagents by mouth.

**3.13.** Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

**3.14.** Material Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

**3.15.** The Material Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

## 4.1. Contents of the Kit

## 4. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1 SORB MTP	C-peptide EIA strips, 8x12 wells	1	pcs		until exp.date
2 CAL 1 - 6	Calibrator set, 0.5 ml each. The set contains 6 calibrators: 0; 0.2; 0.5; 2; 5; 20 ng/ml	6	pcs	blue (C1 - colourless)	2 months
3 CON-TROL	Control serum (0.5 ml)	1	pcs	colourless	2 months
4 CONJ HRP	Conjugate, 14 ml	1	pcs	blue	until exp.date
5 DIL SPE	EIA buffer 22 ml	1	pcs	blue	until exp.date
6 SUBS TMB	Substrate solution, 14 ml	1	pcs	colourless	until exp.date
7 BUF WASH 26X	Washing solution concentrate 26x, 22 ml	1	pcs	colourless	Concentrate - until exp. date Diluted washing solution - 45 days at +2...+8 °C or 15 days at RT
8 STOP	Stop solution, 14 ml	1	pcs	colourless	until exp.date
9 N003	Plate sealing tape	2	pcs	-	N/A
10 K267CI	Instruction C-peptide EIA	1	pcs	-	N/A
11 K267CQ	QC data sheet C-peptide EIA			-	N/A

**4.2.** Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 25–250 µl;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

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**5. SPECIMEN COLLECTION AND STORAGE**

Specimens of serum/plasma may be stored for up to 72 hours at +2...+8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

Urine:

Centrifuge sample to clear. Urine samples may be stored for up to 36 hours at +2...+8 °C prior to assaying.

Dilute sample by DIL SPE in 21 times. Do not store diluted samples!

**6. TEST PROCEDURE****6.1.** Reagent Preparation

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18 to +25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 26X by 26 dilutions in distilled water.

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#### 6.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 14 wells for the calibrators CAL 1 - 6 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Before first use of the kit dissolve the calibrators: add 0.5 ml deionized water to each vial and mix thoroughly avoiding foam formation. Liquid calibrators should be assayed within an 72 hour, aliquoted and stored frozen below -15 °C. For long term storage it is recommended to freeze the calibrators in aliquoted portions and store below -15 °C. ATTENTION: AVOID MORE THAN 1 FREEZE-THAW CYCLES!
3	Pipet 25 µl of calibrators CAL 1 - 6, control samples CONTROL and unknown samples into the wells.
4	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape (included into the kit).
5	Incubate 60 minutes at 18-25°C.
6	Prepare washing solution by 26x dilution of washing solution concentrate (BUF WASH 26X) by distilled water. Wash the strips 5 times.
7	Dispense 100 µl of SUBS TMB into the wells.
8	Incubate 10-20 minutes at 18-25 °C.
9	Dispense 100 µl of STOP into the wells.
10	Measure OD (optical density) at 450 nm.
11	Set photometer blank on first calibrator.
12	Apply point-by-point method for data reduction.
13	Apply point-by-point method for data reduction.

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
**8.2.** Plot a calibration curve on graph paper: OD versus C-peptide concentration.

**8.3.** Determine the corresponding concentration of C-peptide in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

**8.4.** Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

Calibrators	Value	Absorbance Units (450 nm)
CAL 1	0 ng/ml	0.05
CAL 2	0.2 ng/ml	0.23
CAL 3	0.5 ng/ml	0.33
CAL 4	2 ng/ml	0.90
CAL 5	5 ng/ml	1.60
CAL 6	20 ng/ml	2.57

## 9. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for C-peptide. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units, ng/ml		Units alternative, pM	
	Lower limit	Upper limit	Lower limit	Upper limit
serum (plasma)	0.9	5.0	298	1655
urine	2 ug/day	260 ug/day	0	0

## 10. PERFORMANCE CHARACTERISTICS

### 10.1. Analytical specificity / Cross reactivity

Analyte	Cross-reactivity, % wt/wt
Human insulin	<0.01
Human proinsulin (intact)	2.7
Humalog	0.3

### 10.2. Precision

Intra-assay precision is shown below:

Serum, no	duplicated	value, ng/ml	CV1, %
1	10	0.6	6.7
2	10	5.8	7,7

Inter-assay precision is shown below:

Serum, no	duplicated	value, ng/ml	CV, %
1	10	0.85	4.5
2	10	4.3	5.6

### 10.3. Analytical sensitivity

Sensitivity of the assay was assessed as being 0.015 ng/ml.

### 10.4. Linearity

Linearity was checked by assaying dilution series of 5 samples with different C-peptide concentrations. Linearity percentages obtained ranged within 90 to 110%.

### 10.5. Recovery

Recovery was estimated by assaying 5 mixed samples with known C-peptide concentrations. The recovery percentages ranged from 90 to 110%.

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*Instruction for use*

## **A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF TOTAL IgG IN HUMAN BIOLOGICAL FLUIDS**

### **1. INTENDED USE**

A solid-phase enzyme immunoassay for the quantitative determination of total IgG in biological fluids.

This kit is designed for measurement of total IgG in biological fluids. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 42 unknown samples in duplicates.

### **2. SUMMARY AND EXPLANATION**

Immunoglobulin G (IgG) is the main part of serum  $\gamma$  - globulin fraction. IgG is secreted during secondary immune response and plays a key role in humoral immunity. Decrease of serum IgG concentration below 5 g/l is a marker of severe life-threatening immunodeficiency. Determination of serum IgG concentration and IgG/IgA/IgM ratios can be used for monitoring of humoral immune status. Marked elevation of serum IgG may be observed in chronic inflammation, autoimmune diseases and myeloma.

### **3. PRINCIPLE OF THE TEST**

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to human total IgG-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Unbound material is removed by washing procedure. Second antibodies - murine monoclonal to human total IgG, labelled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

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## 5.1. Contents of the Kit

## 5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	SORB MTP total IgG EIA strips, 8x12 wells	1	pcs		until exp.date
2	CAL 1-5 Calibrator set, 1 ml each. The set contains 5 calibrators: 0; 1; 5; 10; 25 g/l	5	pcs	blue (Cl – colourless)	2 months
3	CONTROL Control serum (1 ml)	1	pcs	colourless	2 months
4	CONJ HRP Conjugate, 14 ml	1	pcs	red	until exp.date
5	DIL SPE EIA sample buffer 100 ml	1	pcs	blue	until exp.date
6	SUBS TMB Substrate solution, 14 ml	1	pcs	colourless	until exp.date
7	BUF WASH 26X Washing solution concentrate 26X, 22 ml	1	pcs	colourless	Concentrate – until exp.date Diluted washing solution – 45 days at 2-8 °C or 15 days at RT
8	STOP Stop solution, 14 ml	1	pcs	colourless	until exp.date
9	N003 Plate sealing tape	2	pcs		N/A
10	K2711 Instruction total IgG EIA	1	pcs		N/A
11	K271Q QC data sheet total IgG EIA	1	pcs		N/A

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- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250  $\mu\text{l}$ , is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 10–250  $\mu\text{l}$ ;
- Dry thermostat for 37  $^{\circ}\text{C} \pm 0.1$   $^{\circ}\text{C}$
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

**5.3.** Storage and stability of the Kit

Store the whole kit at +2...+8  $^{\circ}\text{C}$  upon receipt until the expiration date.

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Specimens may be stored for up to 48 hours at +2...+8  $^{\circ}\text{C}$  before testing.

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- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25  $^{\circ}\text{C}$ ) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 26X by 26 dilutions in distilled water.

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#### 7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1-5 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Dilute samples using buffer DIL SPE (EIA sample buffer) 5000 fold. See table M for dilution modes and factors for different types of analyzed material. Do not dilute control sample and calibrators.
3	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using DIL SPE (EIA sample buffer). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
4	Pipet 100 µl of calibrators CAL 1-5 and control samples CONTROL into allocated wells. For testing of blood serum or plasma pipet 100 µl of the unknown diluted sample (DILUTION 2) into the allocated wells. See table M for the volumes of other materials. Pipetting should be made within 3 minutes, to ensure an uniform incubation time for all samples. Carefully mix the contents of the wells by short horizontal rotating of the plate for 5-7 seconds and cover the wells by plate adhesive tape (included into the kit).
5	Incubate 30 minutes at +37 °C.
6	Prepare washing solution by 26X dilution of washing solution concentrate BUF WASH 26X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
7	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
8	Incubate 30 minutes at +37 °C.
9	Wash the strips 5 times.
10	Dispense 100 µl of SUBS TMB into the wells.
11	Incubate 10-20 minutes at +18...+25 °C.
12	Dispense 100 µl of STOP into the wells.
13	Measure OD (optical density) at 450 nm.
14	Set photometer blank on first calibrator.
15	Apply point-by-point method for data reduction. Use Calculation factor listed in table M to calculate analyte concentration in different material types.

### 7.5. Sample processing

Material type	Notes on material collection, storage and handling	Sample dilution example	EIA sample buffer into the well, $\mu\text{l}$	Sample into the well, $\mu\text{l}$	Calculation factor
blood serum or plasma	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	10 $\mu\text{l}$ of sample + 990 $\mu\text{l}$ of diluent = DILUTION 1. 10 $\mu\text{l}$ of DILUTION 1 + 490 $\mu\text{l}$ of diluent = DILUTION 2	0	100	1
saliva	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		90	10	0.002
urine	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		50	50	0.0004
cerebrospinal fluid	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	10 $\mu\text{l}$ of sample + 500 $\mu\text{l}$ of diluent	0	100	0.01

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
**9.2.** Plot a calibration curve on graph paper: OD versus total IgG concentration.

**9.3.** Determine the corresponding concentration of total IgG in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is

recommended due to non-linear shape of curve.

**9.4.** Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

Calibrators	Value	Absorbance Units (450 nm)
CAL 1	0 g/l	0.08
CAL 2	1 g/l	0.38
CAL 3	5 g/l	1.15
CAL 4	10 g/l	1.78
CAL 5	25 g/l	2.30

## 10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for total IgG. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units, g/l	
	Lower limit	Upper limit
newborn	7.0	15
1-3 month	2.7	8.0
4-6 month	1.8	8.5
7-12 month	3.5	12
1-6 yrs	6.5	18
7-11 yrs	8.5	15
> 11 yrs	9.0	20

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**11.2.** Analytical sensitivity. Sensitivity of the assay was assessed as being 0.06 g/l.

**11.3.** Linearity. Linearity was checked by assaying dilution series of 5 samples with different total IgG concentrations. Linearity percentages obtained ranged within 90 to 110%.

**11.4.** Recovery. Recovery was estimated by assaying 5 mixed samples with known total IgG concentrations. The recovery percentages ranged from 90 to 110%.

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*Instruction for use*

## **A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF TOTAL IGA IN HUMAN BIOLOGICAL FLUIDS**

### **1. INTENDED USE**

A solid-phase enzyme immunoassay for the quantitative determination of total IgA in biological fluids.

This kit is designed for measurement of total IgA in biological fluids. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 42 unknown samples in duplicates.

### **2. SUMMARY AND EXPLANATION**

Immunoglobulin A (IgA) is a main factor of mucosal immune response to bacteria and viruses. Selective IgA deficiency is one of the most frequent hereditary disorders causing chronic infections inflammation in gastrointestinal, urinary and respiratory systems. Determination of IgA concentration in serum and other biological fluids can be used as screening for selective IgA deficiency and other immunodeficiency syndromes. Marked elevation of serum IgA is observed in some autoimmune diseases and IgA myeloma.

### **3. PRINCIPLE OF THE TEST**

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to human total IgA-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Unbound material is removed by washing procedure. Second antibodies – murine monoclonal to human total IgA, labelled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

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**4.4.** Avoid contact with stop solution containing 5.0%  $H_2SO_4$ . It may cause skin irritation and burns.

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## 5.1. Contents of the Kit

## 5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	SORB MTP total IgA EIA strips, 8x12 wells	1	pcs		until exp.date
2	CAL 1-5 Calibrator set, 1 ml each. The set contains 5 calibrators: 0; 0.1; 0.5; 2; 5 g/l	5	pcs	blue (C1 – colourless)	2 months
3	CONTROL Control serum (1 ml)	1	pcs	colourless	2 months
4	CONJ HRP Conjugate, 14 ml	1	pcs	blue	until exp.date
5	DIL SPE EIA buffer 100 ml	1	pcs	blue	until exp.date
6	SUBS TMB Substrate solution, 14 ml	1	pcs	colourless	until exp.date
7	BUF WASH 26X Washing solution concentrate 26X, 22 ml	1	pcs	colourless	Concentrate – until exp.date Diluted washing solution - 45 days at 2-8 °C or 15 days at RT
8	STOP Stop solution, 14 ml	1	pcs	colourless	until exp.date
9	N003 Plate sealing tape	2	pcs		N/A
10	K275I Instruction total IgA EIA	1	pcs		N/A
11	K275Q QC data sheet total IgA EIA	1	pcs		N/A

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- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250  $\mu\text{l}$ , is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 5–250  $\mu\text{l}$ ;
- Dry thermostat for +37 °C  $\pm$ 0.1 °C
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

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Specimens may be stored for up to 48 hours at +2...+8 °C before testing.

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#### 7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1-5 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Dilute samples using buffer DIL SPE (EIA buffer) 5000 fold. See table M for dilution modes and factors for different types of analyzed material. Do not dilute control sample and calibrators.
3	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using DIL SPE (EIA buffer). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
4	Pipet 100 µl of calibrators CAL 1-5 and control samples CONTROL into allocated wells. For testing of blood serum or plasma pipet 100 µl of the diluted unknown sample into the allocated wells. See table M for the volumes of other materials. Pipetting should be made within 3 minutes, to ensure an uniform incubation time for all samples. Carefully mix the contents of the wells by short horizontal rotating of the plate for 5-7 seconds and cover the wells by plate adhesive tape (included into the kit).
5	Incubate 30 minutes at +37 °C.
6	Prepare washing solution by 26X dilution of washing solution concentrate BUF WASH 26X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
7	Dispense 100 µl of CONJ HRP into the wells.
8	Incubate 30 minutes at +37 °C.
9	Wash the strips 5 times.
10	Dispense 100 µl of SUBS TMB into the wells.
11	Incubate 10-20 minutes at +18...+25 °C.
12	Dispense 100 µl of STOP into the wells.
13	Measure OD (optical density) at 450 nm.
14	Set photometer blank on first calibrator.
15	Apply point-by-point method for data reduction. Use Calculation factor listed in table M to calculate analyte concentration in different material types.

### 7.5. Sample processing

Material type	Notes on material collection, storage and handling	Sample dilution example	EIA buffer into the well, µl	Sample into the well, µl	Calculation factor
blood serum or plasma	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	10 µl of sample + 990 µl of diluent = DILUTION 1. 10 µl of DILUTION 1 + 490 µl of diluent = DILUTION 2	0	100	1
saliva	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	5 µl of sample + 500 µl of diluent	90	10	0.2
urine	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		80	20	0.001
cerebrospinal fluid	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		50	50	0.0004

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
**9.2.** Plot a calibration curve on graph paper: OD versus total IgA concentration.

**9.3.** Determine the corresponding concentration of total IgA in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

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## 10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for total IgA. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units, g/l	
	Lower limit	Upper limit
Healthy donors	0.9	5.0
>61 yr	1.0	6.5
newborn	-	0.05
1-3 month	0.06	0.6
4-6 month	0.1	1.0
7-12 month	0.35	1.7
1-6 yrs	0.8	2.2
7-11 yrs	0.9	2.6

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**11.3.** Linearity. Linearity was checked by assaying dilution series of 5 samples with different total IgA concentrations. Linearity percentages obtained ranged within 90 to 110%.

**11.4.** Recovery. Recovery was estimated by assaying 5 mixed samples with known total IgA concentrations. The recovery percentages ranged from 90 to 110%.

## 12. LITERATURE

1. Heiddis B. Valdimarsdottir and Arthur A. Stone – Psychosocial Factors and Secretory Immunoglobulin A. Critical Reviews in Oral Biology & Medicine, Jan 1997; 8: 461 – 474.

2. Amir H Abdul Latiff and Michael A Kerr – The clinical significance of immunoglobulin A deficiency. Ann Clin Biochem, Mar 2007; 44: 131 – 139.

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*Instruction for use*

## **A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF TOTAL IgM IN HUMAN BIOLOGICAL FLUIDS**

### **1. INTENDED USE**

A solid-phase enzyme immunoassay for the quantitative determination of total IgM in biological fluids.

This kit is designed for measurement of total IgM in biological fluids. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 42 unknown samples in duplicates.

### **2. SUMMARY AND EXPLANATION**

Immunoglobulin M (IgM) is secreted during primary immune response and exists in monomeric and pentameric forms. Elevated serum IgM is observed in chronic inflammation, macroglobulinemia and IgM myeloma. Decreased IgM level may occur in some immunodeficiency syndromes.

### **3. PRINCIPLE OF THE TEST**

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to human total IgM-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Unbound material is removed by washing procedure. Second antibodies – murine monoclonal to human total IgM, labelled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

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## 5.1. Contents of the Kit

## 5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	total IgM EIA strips, 8x12 wells	1	pcs		until exp.date
2	total IgM EIA strips, 8x12 wells Calibrator set, 1 ml each. The set contains 5 calibrators: 0; 0.5; 2; 5; 10 g/l	5	pcs	purple (C1 – colourless)	2 months
3	Control serum (1 ml)	1	pcs	colourless	2 months
4	Conjugate, 14 ml	1	pcs	purple	until exp. date
5	EIA buffer 100 ml	1	pcs	blue	until exp. date
6	Substrate solution, 14 ml	1	pcs	colourless	until exp.date
7	Washing solution concentrate 26X, 22 ml	1	pcs	colourless	Concentrate – until exp.date Diluted washing solution – 45 days at 2-8 °C or 15 days at RT
8	Stop solution, 14 ml	1	pcs	colourless	until exp. date
9	Plate sealing tape	2	pcs		N/A
10	Instruction total IgM EIA	1	pcs		N/A
11	QC data sheet total IgM EIA	1	pcs		N/A

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- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 5–250 µl;
- Dry thermostat for 37 °C ±0.1 °C
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

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#### 7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1–5 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Dilute samples using buffer DIL SPE (EIA sample buffer) 5000 fold. See table M for dilution modes and factors for different types of analyzed material. Do not dilute control sample and calibrators.
3	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using DIL SPE (EIA sample buffer). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
4	Pipet 100 µl of calibrators CAL 1–5 and control samples CONTROL into allocated wells. For testing of blood serum or plasma pipet 100 µl of the unknown diluted sample (DILUTION 2) into the allocated wells. See table M for the volumes of other materials. Pipetting should be made within 3 minutes, to ensure an uniform incubation time for all samples. Carefully mix the contents of the wells by short horizontal rotating of the plate for 5–7 seconds and cover the wells by plate adhesive tape (included into the kit).
5	Incubate 30 minutes at 37 °C.
6	Prepare washing solution by 26X dilution of washing solution concentrate BUF WASH 26X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
7	Dispense 100 µl of CONJ HRP into the wells.
8	Incubate 30 minutes at 37 °C.
9	Wash the strips 5 times.
10	Dispense 100 µl of SUBS TMB into the wells.
11	Incubate 10–20 minutes at +18...+25 °C.
12	Dispense 100 µl of STOP into the wells.
13	Measure OD (optical density) at 450 nm.
14	Set photometer blank on first calibrator.
15	Apply point-by-point method for data reduction. Use Calculation factor listed in table M to calculate analyte concentration in different material types.

## 7.5. Sample processing

Material type	Notes on material collection, storage and handling	Sample dilution example	EIA sample buffer into the well, $\mu\text{l}$	Sample into the well, $\mu\text{l}$	Calculation factor
blood serum or plasma	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	10 $\mu\text{l}$ of sample + 990 $\mu\text{l}$ of diluent = DILUTION 1. 10 $\mu\text{l}$ of DILUTION 1 + 490 $\mu\text{l}$ of diluent = DILUTION 2	0	100	1
saliva	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		90	10	0.002
urine	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		50	50	0.0004
cerebrospinal fluid	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		80	20	0.001

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
**9.2.** Plot a calibration curve on graph paper: OD versus total IgM concentration.

**9.3.** Determine the corresponding concentration of total IgM in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

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## 10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for total IgM. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units, g/l	
	Lower limit	Upper limit
Healthy donors	0.7	3.7
newborn	0.1	0.35
1-3 month	0.12	0.9
4-6 month	0.25	1.2
7-12 month	0.35	1.0
1-6 yrs	0.55	2.2
7-11 yrs	0.65	1.7

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**11.3.** Linearity. Linearity was checked by assaying dilution series of 5 samples with different total IgM concentrations. Linearity percentages obtained ranged within 90 to 110%.

**11.4.** Recovery. Recovery was estimated by assaying 5 mixed samples with known total IgM concentrations. The recovery percentages ranged from 90 to 110%.

## 12. LITERATURE

Erik J. Wiersma, Cathy Collins, Shafie Fazel, and Marc J. Shulman Structural and Functional Analysis of J Chain-Deficient IgM J. Immunol., Jun 1998; 160: 5979 – 5989.

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