



SYPHILIS SEROLOGY KIT
DIRECTIONS FOR USE

RPR CARBON KIT: For Detection Of Syphilis.

SUMMARY

At one time, syphilis was a major medical disease with a host of different manifestations transmitted primarily through sexual contact. The advent of penicillin in 1943 changed this. The etiologic agent of syphilis is *Treponema pallidum*, a spiral bacterium (spirochete). The spirochete causes some damage to the heart and the liver, releasing some tissue fragments. The patient's immune system produces antibodies, called reagins, against these fragments. There are two different techniques for the detection of syphilis. TPHA tests, which detect antibodies to *Treponema pallidum*, and non-treponemal serologic tests, which detect Reagin in infected people.

PRINCIPLE

When used by the recommended techniques, the reagent will agglutinate (clump) in the presence of reagin. No agglutination usually indicates the absence of reagin (see **Limitations**).

KIT DESCRIPTION

Lorne RPR Carbon Kit is a non-treponemal serologic test for the detection of syphilis. The RPR Carbon Antigen contains micro particulate carbon, which aids in the microscopic reading of results. All the reagents are supplied at optimum dilution for use with all recommended techniques without the need for further dilution or addition. For lot reference number and expiry date see **Vial Labels**.

STORAGE

Do not freeze. Reagent vials should be stored at 2 - 8°C on receipt. Prolonged storage at temperatures outside this range may result in accelerated loss of reagent reactivity.

SPECIMEN COLLECTION

Specimens should be drawn with or without anticoagulant using an aseptic phlebotomy technique. If testing is delayed specimens can be stored at 2-8°C for 7 days or for up to 3 months at or below -20°C. Specimens must be free from bacterial contamination, fibrin, haemolysis and lipaemia.

PRECAUTIONS

1. The kit is for *in vitro* diagnostic use only.
2. Do not use kit past expiration date (see **Vial and Box Labels**).
3. Protective clothing should be worn when handling the reagents, such as disposable gloves and a laboratory coat.
4. No known tests can guarantee products derived from human or animal sources are free from infectious agents. Care must be taken in the use and disposal of each vial and its contents.
5. RPR Positive Control: H319 - Causes serious eye irritation. Follow the precautionary statement given in the SDS.

DISPOSAL OF KIT REAGENT AND DEALING WITH SPILLAGES

For information on disposal of kit reagent and decontamination of a spillage site see **Material Safety Data Sheets**, available on request.

CONTROLS AND ADVICE

1. It is recommended the RPR Positive and Negative Controls be tested in parallel with each batch of tests. Tests must be considered invalid if controls do not show expected results.
2. Shake all the reagents well before use to ensure homogeneity.
3. Do not interchange components between different kits.
4. The circles on the agglutination cards should never be touched with fingers, as this may invalidate the test results.
5. Use of kit and interpretation of results must be carried out by properly trained and qualified personnel in accordance with the requirements of country where reagents are in use.
6. The user must determine suitability of the kit for use in other techniques.

KIT COMPONENTS PROVIDED

- 1) RPR Carbon Antigen (Red Label): Carbon particles coated with a lipid complex (cardiolipin, lecithin and cholesterol) in phosphate buffer 20 mmol/L, pH 7.0 containing a preservative.
- 2) RPR Positive Control (Red cap): Artificial serum with reagin titer $\geq 1/4$.
- 3) RPR Negative Control (Blue cap): Animal serum containing a preservative
- 4) Dispensing bottle (1 x 2 ml).
- 5) Dispensing Needle (x1).
- 6) Disposable agglutination slides.
- 7) Plastic stirrers.

MATERIALS AND EQUIPMENT NOT SUPPLIED

- a) Pipette capable of accurately delivering 50 μ l
- b) Mechanical rotating table capable of rotating at 80-100 rpm.
- c) 9 g/L saline solution.

QUALITATIVE TECHNIQUE

1. Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
2. Place 50 μ L of the sample and one drop of each Positive and Negative controls into separate circles on the slide test.
3. Swirl the RPR-carbon reagent gently before using. Invert the dropper assembly and press gently to remove air bubbles from the micropipette.
4. Place the micropipette in a vertical position and perpendicular to the slide, and add one drop (20 μ L) of this reagent next to the samples to be tested.
5. Mix the drops with a stirrer, spreading them over the entire surface of the circle. Use different stirrers for each sample
6. Place the slide on a mechanical rotating table at 80-100 r.p.m. for 8 min. False positive results could appear if the test is read after more than 8 minutes.

INTERPRETATION OF QUALITATIVE RESULTS

1. **Reactive:** Visible agglutination (medium to large clumps) constitutes a positive result and within the accepted limitations of the test procedure, indicates the presence of reagin.
2. **Weak-Reactive:** Weak agglutination (small clumps) around the periphery of the test area constitutes a weak positive result and within the accepted limitations of the test procedure, indicates the presence of reagin.
3. **Negative:** No agglutination constitutes a negative result and within the accepted limitations of the test procedure, indicates the absence of reagin.

SEMI QUANTITATIVE TECHNIQUE

1. The semi-quantitative test can be performed in the same way as the quantitative technique using dilutions of the serum in 9 g/L saline solution.
2. Make doubling dilutions of specimen as follows:

Dilution	Serum	Saline
1/2	100 μ l undiluted serum	100 μ l
1/4	100 μ l 1/2 diluted serum	100 μ l
1/8	100 μ l 1/4 diluted serum	100 μ l
1/16	100 μ l 1/8 diluted serum	100 μ l

3. Test the specimen dilutions in the same way as for the quantitative technique above.
4. Read the test and note the last positive dilution series.

STABILITY OF THE REACTIONS

Slide tests should be interpreted straight after the 8-minute rotating period to avoid the possibility that a negative result may be incorrectly interpreted as positive due to drying of the reagent.

LIMITATIONS

1. RPR carbon test is non-specific for syphilis. All Reactive samples should be retested with treponemic methods such as TPHA and FTA-Abs to confirm the results.
2. A Non Reactive result by itself does not exclude a diagnosis of syphilis. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.
3. False positive results have been reported in diseases such as infectious mononucleosis, viral pneumonia, toxoplasmosis, pregnancy and autoimmune diseases.
4. Bilirubin (≤ 20 mg/dL), hemoglobin (≤ 10 g/L) and lipids (≤ 10 g/L), do not interfere. Rheumatoid factors (≥ 300 IU/mL), interfere. Other substances may interfere⁵.
5. False positive or negative results may also occur due to:
 - a) Not expelling air from end of needle
 - b) Not maintaining dispensing bottle and needle in a vertical position when dispensing the antigen.
 - c) When transferring the specimen from the collecting tube some of the specimen being drawn up in to the teat
 - d) Contamination of test materials
 - e) Improper storage of test materials or omission of reagents
 - f) Deviation from the recommended techniques

SPECIFIC PERFORMANCE CHARACTERISTICS

1. The kit has been characterised by all the procedures mentioned in the **Recommended Techniques**.
2. Prior to release, each lot of Lorne RPR Syphilis Kit is tested by the **Recommended Techniques** to ensure suitable reactivity.
3. The reagent sensitivity is calibrated against the "Human Reactive Serum" from the CDC (Centres for Disease Control) and comparable to the RPR reagent from Becton Dickinson.
4. **Prozone effect:** No prozone effect was detected up to titers $\geq 1/128$.
5. **Diagnostic sensitivity:** 100%
6. **Diagnostic specificity:** 100 %.

DISCLAIMER

1. The user is responsible for the performance of the kit by any method other than those mentioned in the **Recommended Techniques**.
2. Any deviations should be validated prior to use using established laboratory procedures.

BIBLIOGRAPHY

1. George P. Schmid. Current Opinion in Infectious Diseases 1994; 7: 34-40.
2. Sandra A Larsen et al. Clinical Microbiology Reviews 1995; 8 (1): 1-21.
3. Sandra Larsen et al. A manual of Test for Syphilis American Public Health Association 1990: 1-192.
4. Joseph Earle Moore et al. Gastrointestinal Haemorrhage 1952; 150(5): 467-473.
5. Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.

AVAILABLE KIT SIZES








Kit Size	Catalogue Number
150 Tests Per Kit	044150A
500 Tests Per Kit	044500A

For the availability of other sizes, please contact:

Lorne Laboratories Limited

Unit 1 Cutbush Park Industrial Estate
Danehill
Lower Earley
Berkshire, RG6 4UT
England
Tel: +44 (0) 118 921 2264
Fax: +44 (0) 118 986 4518
E-mail: info@lornelabs.com

TABLE OF SYMBOLS

	Batch Number		<i>in-vitro</i> Diagnostic
	Catalogue Reference		Store At
	Expiry Date		Manufacturer
	Read Pack Insert		

DECLARATION OF CONFORMITY

PRODUCT IDENTIFICATION

Product name	Catalogue number
RPR Carbon kit	044150A 044500A

MANUFACTURER

Name	Lorne Laboratories
Address	Unit 1 Cutbush Park Industrial Estate Danehill Lower Earley Berks, RG6 4UT
Country	United Kingdom

MEANS OF CONFORMITY

I hereby declare that the products listed above comply with the essential requirements and provisions of Directive 98/79/EC of the European Parliament and of the Council (also SI 2002 No.618 which transposes the requirements of Directive 98/79/EC).

This declaration is valid from 17 May 2015.



Eddy Velthuis
Technical Director

Certificate of Approval

This is to certify that the Management System of:

Inzek International Trading B.V.

Vissenstraat 32, 7324 AL Apeldoorn, Netherlands

has been approved by LRQA to the following standards:

ISO 13485:2016



P.G. Cornelissen - Area Manager North Europe

Issued by: Lloyd's Register Nederland B.V.

for and on behalf of: Lloyd's Register Quality Assurance Limited

Current issue date: 9 March 2019
Expiry date: 8 March 2022
Certificate identity number: 10177988

Original approval(s):
ISO 13485 – 9 March 2019

Approval number(s): ISO 13485 – 00019238

The scope of this approval is applicable to:

Design, development, production and distribution of In Vitro Diagnostics Medical device - reagents and instrument for point of care testing.



001

Certificate



Quality Management System EN ISO 13485:2016

Registration No.: SX 1034230-1

Organization: nal von minden GmbH
Carl-Zeiss-Str. 12
47445 Moers
Germany

Scope: Design and development, manufacture and distribution of in vitro diagnostic test kits and reagents for the detection or determination of cardiac markers, tumor markers, infections, inflammation, allergies, endocrine disorders, diabetes, hormones, vitamins, special proteins, metabolic disorders, drug misuse, immune status, vaginal pH levels, pregnancy, kidney function, for urine analysis, sperm testing, coagulation management systems, for use in clinical laboratories, as near-patient tests and for self-testing as well as associated in vitro diagnostic devices for sampling and analysis systems for rapid tests.

Distribution of control materials and tests for blood group determination, medical masks, medical gloves, blood pressure monitoring devices, medical thermometers, swabs and lancets.

The Certification Body of TÜV Rheinland LGA Products GmbH certifies that the organization has established and applies a quality management system for medical devices. Proof has been furnished that the requirements specified in the abovementioned standard are fulfilled. The quality management system is subject to yearly surveillance.

Report No.: 1089325-40
Effective date: 2021-12-02
Expiry date: 2024-12-01
Issue date: 2021-11-29



Dipl.-Ing. Sven Hoffmann
TÜV Rheinland LGA Products GmbH
Tillystraße 2 · 90431 Nürnberg · Germany



Certificate



Quality Management System EN ISO 13485:2016

Registration No.: SX 1034230-1



Organization: nal von minden GmbH
Carl-Zeiss-Str. 12
47445 Moers
Germany

The scope of certification also covers the following:

No.	Facility	Scope
/01	c/o nal von minden GmbH Carl-Zeiss-Str. 12 47445 Moers Germany	Manufacture and distribution
/02	c/o nal von minden GmbH Friedenstr. 32 93053 Regensburg Germany	Design and development and distribution
/03	c/o nal von minden GmbH Robert-Bosch-Breite 34 37079 Göttingen Germany	Design and development and manufacture
/04	c/o nal von minden GmbH Raseweg 4 37124 Rosdorf Germany	Administration and distribution

Report No.: 1089325-40
Effective date: 2021-12-02
Expiry date: 2024-12-01
Issue date: 2021-11-29





Dipl.-Ing. Sven Hoffmann
TÜV Rheinland LGA Products GmbH
Tillystraße 2 · 90431 Nürnberg · Germany

EC Certificate

Directive 98/79/EC Annex IV, excluding Sections 4 and 6
Full Quality Assurance System
In Vitro Diagnostic Medical Devices

Registration No.: HL 60131398 0001

Report No.: 21200072 015

Manufacturer: nal von minden GmbH
Carl-Zeiss-Str. 12
47445 Moers
Deutschland

Products:

- IVDs for the detection of infectious disease markers
- IVDs for the detection of the tumor marker PSA
- Urine tests for self-testing

(see attachment for products and sites included)

Replaces Certificate, Registration No.: HL 60114562 0001

Expiry Date: 2023-11-27

The Notified Body hereby declares that the requirements of Annex IV, excluding section 4 and 6 of the directive 98/79/EC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex IV, section 5 of the aforementioned directive. For placing on the market of List A devices covered by this certificate an EC design-examination certificate according to Annex IV, section 4 and a verification of manufactured products according to section 6 is required.

Effective Date: 2018-11-28

Date: 2018-11-27

Notified Body



Dipl.-Ing. Sven Hoffmann

TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg
TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 98/79/EC concerning in vitro diagnostic medical devices with the identification number 0197.

TÜV Rheinland
LGA Products GmbH
Tillystraße 2, 90431 Nürnberg

**Attachment to
Certificate**

Registration No.: HL 60131398 0001
Report No.: 21200072 015

Manufacturer: nal von minden GmbH
Carl-Zeiss-Str. 12
47445 Moers
Deutschland

Products included:

In vitro diagnostica for self-testing:

- HCG pregnancy tests
- LH ovulation tests
- Single- and multi-constituent test strips for urinalysis

In vitro diagnostica rapid tests:

- Chlamydia trachomatis Rapid Tests
- PSA Rapid Tests

Site included:

nal von minden GmbH
Friedenstr. 32
93053 Regensburg
Germany

Activities: Design and development

Date: 2018-11-27

Notified Body




Dipl.-Ing. Sven Hoffmann

Certificate

Standard **ISO 9001:2015**

Certificate Registr. No. **01 100 1810016**

Certificate Holder: **nal von minden GmbH**
Carl-Zeiss-Str. 12
47445 Moers
Germany

including the locations according to annex

Scope: Design and development, manufacture and distribution of in-vitro diagnostic test kits and reagents for the detection of cardiac markers, tumour markers, infections, inflammation, allergies, endocrine disorders, diabetes, hormones, vitamins, special proteins, metabolic disorders, drug misuse, immune status, vaginal pH levels, pregnancy, kidney function, for urine analysis, performance diagnostics, sperm testing, coagulation management systems, for use in clinical laboratories, as near-patient tests and for self-testing and associated in-vitro diagnostic devices for sampling and analysis systems, as well as veterinary diagnostic devices and testing for narcotics and substances.
Distribution of control materials and tests for blood group determination, medical masks, medical gloves, blood pressure monitoring devices, medical thermometers, swabs, lancets and laboratory analyses.

Proof has been furnished by means of an audit that the requirements of ISO 9001:2015 are met.

Validity: The certificate is valid from 2021-09-10 until 2024-09-09.
First certification 2018

2021-09-10



TÜV Rheinland Cert GmbH
Am Grauen Stein · 51105 Köln

Annex to certificate

Standard **ISO 9001:2015**

Certificate Registr. No. **01 100 1810016**

No.	Location	Scope
/01	c/o nal von minden GmbH Carl-Zeiss-Str. 12 47445 Moers Germany	Design and development, manufacture and distribution of in-vitro diagnostic test kits and reagents for the detection of cardiac markers, tumour markers, infections, inflammation, allergies, endocrine disorders, diabetes, hormones, vitamins, special proteins, metabolic disorders, drug misuse, immune status, vaginal pH levels, pregnancy, kidney function, for urine analysis, performance diagnostics, sperm testing, coagulation management systems, for use in clinical laboratories, as near-patient tests and for self-testing and associated in-vitro diagnostic devices for sampling and analysis systems, as well as veterinary diagnostic devices and testing for narcotics and substances. Distribution of control materials and tests for blood group determination, medical masks, medical gloves, blood pressure monitoring devices, medical thermometers, swabs, lancets and laboratory analyses.

Annex to certificate

Standard **ISO 9001:2015**

Certificate Registr. No. **01 100 1810016**

/02 c/o nal von minden GmbH
Friedenstr. 32
93053 Regensburg
Germany

Design and development, distribution of in-vitro diagnostic test kits and reagents for the detection of cardiac markers, tumour markers, infections, inflammation, allergies, endocrine disorders, diabetes, hormones, vitamins, special proteins, metabolic disorders, drug misuse, immune status, vaginal pH levels, pregnancy, kidney function, for urine analysis, performance diagnostics, sperm testing, coagulation management systems, for use in clinical laboratories, as near-patient tests and for self-testing and associated in-vitro diagnostic devices for sampling and analysis systems, as well as veterinary diagnostic devices and testing for narcotics and substances. Distribution of control materials and tests for blood group determination, medical masks, medical gloves, blood pressure monitoring devices, medical thermometers, swabs, lancets and laboratory analyses.

Annex to certificate

Standard **ISO 9001:2015**

Certificate Registr. No. **01 100 1810016**

/03 c/o nal von minden GmbH
Robert-Bosch-Breite 34
37079 Göttingen
Germany

Design and development, manufacture and distribution of in-vitro diagnostic test kits and reagents for the detection of cardiac markers, tumour markers, infections, inflammation, allergies, endocrine disorders, diabetes, hormones, vitamins, special proteins, metabolic disorders, drug misuse, immune status, vaginal pH levels, pregnancy, kidney function, for urine analysis, performance diagnostics, sperm testing, coagulation management systems, for use in clinical laboratories, as near-patient tests and for self-testing and associated in-vitro diagnostic devices for sampling and analysis systems, as well as veterinary diagnostic devices and testing for narcotics and substances. Distribution of control materials and tests for blood group determination, medical masks, medical gloves, blood pressure monitoring devices, medical thermometers, swabs, lancets and laboratory analyses.

2021-09-10


TÜV Rheinland Cert GmbH
Am Grauen Stein · 51105 Köln

EC DECLARATION OF CONFORMITY

ZAO "Vector-Best" hereby ensures under own responsibility and declares that the products listed on pages 2,4 are in conformity with applicable provisions and fulfill the essential requirements of Annex I Directive 98/79/EC of 27 October 1998 regarding in vitro diagnostic medical devices.


Classification of products: Other devices (all devices except Annex II and self-testing devices)

Conformity assessment procedure: Annex III (not including section G).

Manufacturer:
ZAO "Vector-Best"
Address: AHC, Koksovo,
Novosibirsk Region, 630559, Russia,
Tel: +7 (383) 363 20 60,
Fax: +7 (383) 363 35 55

European authorized representative:
Biorion GmbH,
Rheinhorststr. 18, D-67071
Ludwigshafen, Germany.
Tel.: +49 (0) 621 5720 915,
fax: +49 (0) 621 5720 916

Date: 2013/04/12



Murat Khuzainov
General Director ZAO «Vector-Best»

No.	Product name	Identification data	REF
1.	Vectohep A-IgM	ELISA kit for determination of IgM to hepatitis A virus	D-0352
2.	Vectohep A-IgG	ELISA kit for quantitative and qualitative determination of IgG to hepatitis A virus	D-0362
3.	Vectohep TTV-IgG	ELISA kit for determination of IgG to TT virus	D-0802
4.	Vectohep E-IgG	ELISA kit for determination of IgG to hepatitis E virus	D-1056
5.	Vectohep E-IgM	ELISA kit for determination of IgM to hepatitis E virus	D-1058
6.	Vectohep G-IgG	ELISA kit for determination of IgG to hepatitis G virus	D-1252
7.	LymeBest-IgG	ELISA kit for determination of IgG to infectious borreliosis agents	D-1452
8.	LymeBest-IgM	ELISA kit for determination of IgM to infectious borreliosis agents	D-1454
9.	RecombBest antipallidum-IgG	ELISA kit for determination of IgG to Treponema pallidum	D-1852
10.	RecombBest antipallidum-total antibodies	ELISA kit for determination of total antibodies to Treponema pallidum	D-1856
11.	RecombBest antipallidum-IgM	ELISA kit for determination of IgM to Treponema pallidum	D-1858
12.	RecombBest antipallidum-total antibodies	ELISA kit for determination of total antibodies to Treponema pallidum	D-1857
13.	VectohSV-1,2 - IgG	ELISA kit for determination of IgG to herpes simplex virus types 1 and 2	D-2152
14.	VectohSV - IgM	ELISA kit for determination of IgM to herpes simplex virus types 1 and 2	D-2154
15.	VectohHV-8 - IgG	ELISA kit for determination of IgG to human herpes virus type 8	D-2160
16.	VectohHV-8 - IgG	ELISA kit for determination of IgG to human herpes virus type 8	D-2166
17.	Ureaplasma urealyticum - IgG-EIA-BEST	ELISA kit for determination of IgG to Ureaplasma urealyticum antigens	D-2254
18.	Ureaplasma urealyticum - IgA-EIA-BEST	ELISA kit for determination of IgA to Ureaplasma urealyticum antigens	D-2258
19.	VectoParotitis-IgG	ELISA kit for determination of IgG to parotitis virus	D-2602
20.	VectoParotitis-IgM	ELISA kit for determination of IgM to parotitis virus	D-2604
21.	Toxocara-IgG-EIA-BEST	ELISA kit for determination of IgG to toxocara antigens	D-2752
22.	Opisthorchiasis - IgG-EIA-BEST	ELISA kit for determination of IgG to opisthorchiasis antigens	D-2952
23.	Echinococcus-IgG-EIA-BEST	ELISA kit for determination of IgG to Echinococcus	D-3356

24.	Ascend-IgG-EIA-BEST	antigens	ELISA kit for determination of IgG to Ascans lumbricoides	D-3452
25.	Lambli-a-antibodies-EIA-BEST		ELISA kit for determination of IgG, IgM and IgA to Lambli-a antibodies	D-3552
26.	Lambli-a-IgM-EIA-BEST		ELISA kit for determination of IgM to Lambli-a antibodies	D-3554
27.	Lambli-a-antigen-EIA-BEST		ELISA kit for determination of Lambli-a antigen	D-3556
28.	Helicobacter pylori-Caga-antigen-EIA-BEST		ELISA kit for determination of total antibodies to Caga Helicobacter pylori	D-3752
29.	TSHEIA-BEST		ELISA kit for determination of concentration of thyroid-stimulating hormone	X-3952
30.	T3 total-EIA-BEST		ELISA kit for determination of concentration of total triiodothyronine	X-3954
31.	T4 total-EIA-BEST		ELISA kit for determination of concentration of total thyroxine	X-3956
32.	Anti-TPO-EIA-BEST		ELISA kit for determination of antibody concentration to thyroperoxidase	X-3968
33.	PAPP-A-EIA-BEST		ELISA kit for determination of concentration of pregnancy-associated plasma protein A	D-4160
34.	Mycoplasma hominis-IgG-EIA-BEST		ELISA kit for determination of IgG to Mycoplasma hominis	D-4352
35.	Mycoplasma hominis-IgA-EIA-BEST		ELISA kit for determination of IgA to Mycoplasma hominis	D-4356
36.	Mycoplasma pneumoniae-IgG-EIA-BEST		ELISA kit for determination of IgG to Mycoplasma pneumoniae	D-4362
37.	Mycoplasma pneumoniae-IgM-EIA-BEST		ELISA kit for determination of IgM to Mycoplasma pneumoniae	D-4366
38.	Vedochmean - CHF - IgG		ELISA kit for determination of IgG to Crimsean-Congo hemorrhagic fever virus	D-5052
39.	Vedochmean - CHF - IgM		ELISA kit for determination of IgM to Crimsean-Congo hemorrhagic fever virus	D-5054
40.	CEA-EIA-BEST		ELISA kit for determination of concentration of carcinoembryonic antigen	T-8454
41.	AFP-EIA-BEST		ELISA kit for determination of concentration of Alpha-Fetal Protein	T-8456
42.	CA-125-EIA-BEST		ELISA kit for determination of concentration of oncomarker CA-125	T-8466
43.	CA 19-9-EIA-BEST		ELISA kit for determination of concentration of CA 19-9	T-8470
44.	CA 15-3-EIA-BEST		ELISA kit for determination of concentration of oncomarker CA 15-3	T-8472
45.	NSE-EIA-BEST		ELISA kit for determination of concentration of neuron specific enolase	T-8476

46.	Ferritin-EIA-BEST		ELISA kit for determination of concentration of ferritin	T-8552
47.	IgE total-EIA-BEST		ELISA kit for determination of concentration of total IgE	A-8660
48.	IgG total-EIA-BEST		ELISA kit for determination of concentration of total IgG	A-8662
49.	IgM total-EIA-BEST		ELISA kit for determination of concentration of total IgM	A-8664
50.	IgA total-EIA-BEST		ELISA kit for determination of concentration of total IgA	A-8666
51.	Gamma-Interferon-EIA-BEST		ELISA kit for determination of concentration of gamma-interferon	A-8752
52.	Interleukine-4-EIA-BEST		ELISA kit for determination of concentration of interleukine-4	A-8754
53.	Alpha-TNF-EIA-BEST		ELISA kit for determination of concentration of alpha-tumor necrosis factor	A-8756
54.	Alpha-Interferon-EIA-BEST		ELISA kit for determination of concentration of alpha-interferon	A-8758
55.	Interleukine-6-EIA-BEST		ELISA kit for determination of concentration of interleukine-6	A-8768
56.	Interleukine-2-EIA-BEST		ELISA kit for determination of concentration of interleukine-2	A-8772
57.	Procalcitonin-EIA-BEST		ELISA kit for determination of concentration of procalcitonin	A-9004
58.	NTproBNP-EIA-BEST		ELISA kit for determination of concentration of N-terminal prohomone of brain natriuretic peptide	A-9102
59.	Troponin I-EIA-BEST		ELISA kit for determination of concentration of troponin I	A-9106