Helios

DRF

Helios DRF

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*Manufactured by Assing,
created for excellence*

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**Helios DRF: lightness, compactness
and reduced dimensions.**

Each movement has been optimized to facilitate the use by the operator and to ensure maximum patient safety, preventing his movement during the examination

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The minimum height from ground of the patient tabletop is just 45 cm, the lowest of the category.

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X-ray tube head motorized rotation $\pm 180^\circ$


It allows the directing of the beam towards any external receptor device. As an example, when the table is tilted at $+90^\circ$, the X-ray tube can be directed towards a stretcher or a wall stand. It is activated from remote or from the two control boards placed on DRF (B) and tube front (A).

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Focal distance up to 200cm

Adjustable focus to film distance from 105 to 200cm. Helios DRF allows chest examinations directly on the table, be it in horizontal or vertical position.

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Multifocal grid

A multifocal grid is provided, which can be extracted and stored in the housing located beneath the flat panel holder.

The operator can select among the following modes:

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Elevated Coverage

The large dimensions of the tabletop 250x72cm (optional 250x80cm) linked to the longitudinal excursion of the X-ray tube column and detector allows a patient total scan of 210cm. In this case there is no need to reposition the patient. **It is a perfect solution when the patient is not collaborating.**

Tomography

The tube motorized rotation is suitable to the execution of tomography and oblique projections without the mechanical bar support. The tomography angles (for a total of 5) can be chosen by the customer at the order.

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
Table Remote Control Console.

Remote console through a 10,4" touch screen monitor (19" as option). It allows all table movements and exam settings. The console includes an additional set of 4 joysticks to control the main movements as tilting, lateral run, longitudinal movement, focal distance, etc.

Best ergonomics:

all table controls are further duplicated on the two control boards placed on flat panel and X-ray tube front.

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TLK Movement

When the TLK option is chosen, the tilting movement of the patient tabletop clears the exam area from itself so allowing contact radiographies as well as exams on patients on stretchers (radiography and fluoroscopy).

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Abdominal compressor

The abdominal compressor with automatic parking is optionally available. It is controlled from console with display of the compression force and from the table control boards. Retrofit available.

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

Helios DRF can be combined to the use of a tube stand, a vertical stand and a 35x43 or 43x43 cm cassette-sized wireless flat panel.

With the use of a wireless portable detector and a ceiling tube stand the system achieves maximum flexibility making imaging possible virtually everywhere in the room, as conducting examinations on patients in wheelchairs or stretchers. Adding a vertical bucky stand would complete the spectrum of upright examinations such as lateral projections.

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The Helios DRF remote-controlled system is fully integrated and validated with the ERACLE software platform, developed by our partner Digitec, for the acquisition, processing, visualization, and management of digital radiological images in radiography and fluoroscopy modes (continuous and pulsed).

The ERACLE platform supports imaging with both static and dynamic detectors.

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*Opzionale

Dynamic Imaging for fluoroscopic and angiographic applications:

- Gastrointestinal examinations
- Urogenital tract examinations
- Endoscopy
- Arthrography
- Venography
- Lymphography
- Myelography
- Step Angio
- Digital angiography*
- Digital subtracted angiography*

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
*Opzionale

Static Imaging for radiographic examinations

- Skeletal radiography
 - Skull
 - Thorax
 - Abdomen
 - Vertebral column
 - Pelvis
 - Upper and lower extremities
 - Tomography
 - Stitching *
 - Tomosynthesis*

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The stitching function consists of an automatic reconstruction of an X-ray image starting from a series of images acquired during the scanning of the patient. The image is reconstructed, keeping all original pixels, and can be viewed on the monitor, processed, printed or sent to the network. This function is usually used for the scan of spine and lower arts.

Using the "stitching" command, the system automatically generates a sequence of acquisitions according to the selected number of frames (2 for 70 cm, 3 for 100 cm, 4 for 130 cm, 5 for 160 cm length). Afterwards, the software generates one "linked" image in about 5 seconds.

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DSA is a multiple shot radiographic technique used visualize vessels. Radiopaque structures, such as bones, are digitally eliminated – subtracted- from the image, thus enabling an accurate depiction of the blood vessels

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Advanced applications for dose reduction and image quality improvement

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Identifies the presence of the typical disturbances caused by an anti-scatter grid and eliminates the artefact.

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Through the Software Grid, it is possible to take examinations without a grid so reducing patient dose. An algorithm simulates the physical grid resulting into high resolution images with no artifacts. It is mainly recommended for thorax or abdomen exams where the x-ray diffusion is generally very high.

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
CHEST X-RAY IMAGE
WITHOUT BOOST LINES


CHEST X-RAY IMAGE
WITH BOOST LINES

Boost Lines is an image-processing algorithm used to enhance the visibility of tubes and catheters in chest X-ray images, giving a better information about their position.

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SOFT IMAGE
soft tissue component


BONE IMAGE
hard tissue component

Enables to separate soft tissue and bones in different images obtained from double chest X-ray exposure. It is an ideal technique for chest X-rays to detect pulmonary nodules reducing dosage compared to a dedicated CT scan.

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CHEST X-RAY IMAGE
WITHOUT BONE SUPPRESSION



CHEST X-RAY IMAGE
WITH BONE SUPPRESSION

Artificial Intelligence algorithm that automatically removes bones from single shot chest X-ray images.

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The Tomosynthesis is a technique which allows the reconstruction of volumetric images starting from a finite number of bi-dimensional projections taken at different tube angles. The iterative approach is a method that gradually refines the result, repeating the process multiple times until a clearer and noise-free image is obtained.

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
The position of the collimator and AEC chambers is identified directly on the patient overview display.


- Collimated area
- Active AEC chambers
- Inactive AEC chamber

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
The Tracking and Autocollimation function assists the system by calculating the necessary movement after the APR has been selected to proceed with the exposure.

CHECK THE POSITIONING
OF THE PATIENT

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
Patient's body thickness is measured and resulting values are transmitted to the acquisition software for a possible size suggestion.

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Assing has 3 branch offices in Italy and a network of authorized dealers in Italy and abroad.

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A good product is nothing without the right after-sales service.

Assing dedicated after sales service is our flagship. It is operating from our branch offices all over Italy and can constantly assist our customers.

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The after sales service can be contacted h24.