# INSTRUCTIONS FOR USE

Apollo™  
Onyx Delivery Micro Catheter

## TABLE OF CONTENTS

# *Apollo™ Onyx Delivery Micro Catheter*

---

English . . . . .	3
Français . . . . .	4
Deutsch . . . . .	5
Italiano . . . . .	6
Español . . . . .	8
Svenska . . . . .	9
Nederlands . . . . .	10
Português . . . . .	12
Suomi . . . . .	13
Dansk . . . . .	14
Ελληνικά . . . . .	15
Čeština . . . . .	17
Magyar . . . . .	18
Русский . . . . .	19
Polski . . . . .	21
Türkçe . . . . .	22
Norsk . . . . .	23
Slovenčina . . . . .	24
Română . . . . .	25
Български . . . . .	27
한국어 . . . . .	28
عربي . . . . .	30
<b>Symbol Glossary . . . . .</b>	<b>31</b>

# English

## Instructions for Use

EN

### Apollo™ Onyx Delivery Micro Catheter

#### CAUTION

Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician.

This device should be used only by physicians with a thorough understanding of angiography and/or percutaneous neurointerventional procedures.

#### CONTENTS

One Apollo Onyx Delivery Micro Catheter with Syringe Adapter (50414).


#### DESCRIPTION

The Apollo Onyx Delivery Catheter is a single-lumen, endhole catheter designed for the superselective infusion of physician-specified therapeutic agents such as embolization materials and diagnostic materials such as contrast media in tortuous, distal vessels. The catheter has a semi-rigid proximal shaft and a highly flexible distal shaft to facilitate the advancement of the catheter in the anatomy. The proximal end of the catheter incorporates a standard luer adapter to facilitate the attachment of accessories. The outer surfaces of the catheter are coated to increase lubricity.

The Apollo catheter is designed to facilitate catheter retrieval in the event the catheter becomes entrapped within the vasculature. The distal section of the catheter incorporates a detachment zone that allows detachment of the distal tip when the force required to extract the catheter exceeds the force to detach the tip. The catheter has 2 radiopaque marker bands to visualize the position of the catheter and the detach zone area:

- Proximal to the detach zone
- At the distal end of the catheter

REF	A	B - Detachable Tip Length	C	D
105-5095-000	0.5 mm	1.5 cm / 15 mm	1.25 mm	1.9 F/0.63 mm
105-5096-000	0.5 mm	3 cm / 30 mm	1.25 mm	1.9 F/0.63 mm
105-5097-000	0.5 mm	5 cm / 50 mm	1.25 mm	1.9 F/0.63 mm

The Apollo is packaged with a Syringe Adapter (50414). This device, attached to a 1 mL syringe filled with Onyx™, will reduce the dead space within the micro catheter luer hub. Reducing the dead space within the hub is intended to minimize the potential mixing of Onyx® and DMSO in the hub of the catheter during connection and injection.

#### INDICATIONS FOR USE

The Apollo Onyx Delivery Micro Catheter is intended to access the neuro vasculature for the controlled selective infusion of physician-specified therapeutic agents such as embolization materials and of diagnostic materials such as contrast media.

#### CONTRAINDICATIONS

- The Apollo Onyx Delivery Micro Catheter is contraindicated when, in the medical judgment of the physician, use of such product may compromise the patient's condition.
- Not intended for use in the coronary vasculature.
- The Apollo Onyx Delivery Micro Catheter is contraindicated for neonatal and pediatric use.

#### POTENTIAL COMPLICATIONS

Potential complications include, but are not limited to:

- Puncture site hematoma
- Vessel perforation
- Vessel spasm
- Hemorrhage
- Pain and tenderness
- Thrombotic episodes
- Neurological deficits including stroke and death
- Vascular thrombosis

#### WARNINGS

- Not intended for use with embolization particles, detachable coils or Onyx HD500.
- Do not steam shape the tip of the microcatheter. Steam shaping the catheter tip may cause damage to the detachment zone and unintended detachment.

#### WARNINGS

- Do not use a cannula or needle to introduce the end of the catheter into a connector. Insertion of a cannula or needle may cause damage to the detachment zone and result in unintended detachment.
- Always handle the distal end of the catheter with care to avoid damage to the detachment zone and unintended detachment.
- Infusion pressure with this device should not exceed 690 kPa/100 psi. Pressure in excess of 690 kPa/100 psi may result in catheter rupture, possibly resulting in patient injury.
- If flow through the catheter becomes restricted, do not attempt to clear the device by high pressure infusion. Remove the catheter and replace it with a new catheter. Excessive pressure may cause catheter rupture, possibly resulting in patient injury.
- Never advance or withdraw an intraluminal device against resistance. Excessive force against resistance may result in damage to the device or vessel perforation.
- When injecting contrast for angiography, ensure that the catheter is not kinked, prolapsed or occluded. Remove excess slack in the catheter to reduce the potential for catheter kink or prolapse.
- Verify catheter integrity prior to re-inserting guidewire or injecting embolic material to prevent vascular damage or unintended embolization. Catheter integrity is verified by angiographically confirming that contrast agent is exiting only from the catheter tip while viewing the entire distal section of the catheter.
- Regardless of the liquid embolic being used, leave a gap between the reflux and the proximal marker band. Excessive reflux may result in difficult catheter removal.
- This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilization increase the risks of patient infection and compromised device performance.

#### PRECAUTIONS

- Select tip size based on angioarchitecture. The detachment zone should never be distal to the last tortuous curve of the vessel. Refluxing over the detachment zone distal to the last tortuous curve may result in catheter entrapment. Do not place catheter such that the detached tip could interfere with patent vessels.
- Prior to use, carefully examine the Apollo and its packaging to verify that it has not been damaged during shipment. Do not touch or manipulate the catheter tip prior to use.
- Prior to use, all accessory devices and agents should be fully prepared according to the manufacturer's instructions.
- During navigation, check that the distal tip of the catheter is not kinked before passing the guidewire through it. Kinking or prolapsing of the catheter may result in unintended rupture of the catheter.
- Always monitor infusion rates when using the catheter.
- The Apollo has a hydrophilic coating on the outside of the catheter which must be kept hydrated.
- This catheter is not intended for use with chemotherapy agents.
- When the infusion catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing the resultant tip response.
- Navigating or repositioning the catheter while it is in a wedged position or with vessels that are in vasospasm may cause premature tip detachment.
- Do not reposition catheter after start of Onyx injections.
- When performing angiography, it is recommended to use a 3cc syringe rather than a 1cc syringe to reduce the risk of catheter over-pressurization.
- The Apollo is a flow directed micro catheter that can optionally be used with hydrophilic, 0.010" or less sized guidewires. The Apollo is not compatible with non-hydrophilically coated guidewires or guidewires greater than 0.010" in diameter.
- It is recommended that the Apollo be used with an appropriately sized guiding catheter which allows adequate clearance (minimum internal diameter of 0.053" or 1.35mm).
- When withdrawing the catheter, monitor the distal tip under angiography. Pulling the catheter against significant resistance can cause patient injury. If significant catheter resistance is felt, refer to the precaution in the procedure section below for guidance.
- If catheter entrapment is suspected (with any embolic agent), fast catheter retrieval technique may result in catheter shaft separation and potential vascular damage. Follow catheter retrieval instructions at the end of instructions for use.

#### STORAGE

Store the Apollo Onyx Delivery Micro Catheter in a dry place between 50°F (10°C) and 90°F (32°C).

#### DIRECTIONS FOR USE

##### Procedure

1. Place the appropriate guiding catheter following recommended procedures. Connect a one-way stop-cock to the guiding catheter to prevent backflow of blood during insertion of catheter. Connect a hemostatic side arm adapter and a one-way stop-cock to allow for insertion of catheter and to facilitate continuous flush of the guiding catheter with saline.  
**NOTE:** It is recommended that a guiding catheter with a minimum internal diameter of 0.053" (1.35 mm) be used with the Apollo Onyx Delivery Micro Catheter.
2. Gently remove the packaging coil for the Apollo Onyx Delivery Micro Catheter from the package. Flush the coil with heparinized saline through the female luer attached to the coil.

3. Flush the coil with saline through the female luer attached to the coil.
4. Remove the catheter by removing the hub from the clip and gently pulling on the hub.
5. Prior to use, flush the catheter lumen with heparinized saline by attaching a saline filled syringe to the catheter hub.
6. Remove the appropriate steerable guidewire from its packaging (the Silverspeed .010 and the Mirage .008 have been qualified for use with the Apollo) and check for damage. Follow the manufacturer's instructions for preparing and using the guidewire.
7. Inspect the catheter, taking care to not touch or manipulate the tip prior to use to verify that it is undamaged. Retain the packaging coil for storage of the catheter when it is not in use during the procedure.
8. Carefully insert the guidewire into the hub of the micro catheter and advance the guidewire into the catheter lumen.
9. A split introducer is loaded on the proximal end of the catheter to aid in insertion into the hemostatic side-arm adapter. To use, slide the Split Introducer from the proximal end until it covers the distal end of the catheter.
10. Close the one-way stop-cock.
11. Loosen the hemostatic valve.
12. Introduce the guidewire, micro catheter as a unit through the hemostatic valve. Advance the guidewire/catheter assembly to the distal tip of the guiding catheter. Slide the introducer back to the proximal end of the catheter next to the hub.
13. Tighten the valve around the catheter to prevent backflow, but allowing some movement through the valve by the catheter.
14. Open the one-way stop-cock.
15. The catheter may be advanced through the guidecatheter and the vasculature by gently pushing the proximal shaft. It is recommended to use a guidewire during navigation to reduce the risk of catheter kink or prolapse.

#### WARNING

- Verify catheter integrity prior to re-inserting guidewire or injecting embolic material to prevent vascular damage or unintended embolization. Catheter integrity is verified by angiographically confirming that contrast agent is exiting only from the catheter tip while viewing the entire distal section of the catheter.

16. To infuse, connect a syringe with infusate to the catheter luer, and infuse as required.

#### SYRINGE ADAPTER (50414) DIRECTIONS FOR USE:

Usable Length	Apollo Minimum Dead Space Volume	Catheter Dead Space with Syringe Adapter (50414)	Approximate Infusion Rate at 100 psi (690 kPa)	
			Water	Contrast (76% Renografin)
165 cm	>.23 ml	>.20ml	12 ml/min	3 ml/min

1. Prepare the Onyx according to the Onyx Instructions for Use.
2. Following aspiration of Onyx into syringe, detach the needle and attach the Syringe Adapter to the syringe. Purge air from the Syringe Adapter.
3. Immediately connect the Onyx syringe with the Syringe Adapter to the hub of the previously DMSO primed micro catheter according to the Onyx Instructions for Use. Make sure the connection is tight and there is no air in the hub following connection.
4. Proceed with the Onyx injection according to the Onyx Instructions for Use.
5. If a second Onyx syringe is required, **do not remove the Syringe Adapter from the catheter.** When the new Onyx syringe is ready, simply remove the empty syringe from the interface device's proximal end and connect the next syringe to the Syringe Adapter being careful to avoid introducing air.

#### CATHETER RETRIEVAL INSTRUCTIONS:

**PRECAUTION:** If catheter entrapment is suspected (with any embolic agent), fast catheter retrieval technique may result in catheter shaft separation and potential vascular damage.

1. Slowly remove any 'slack' in the distal catheter shaft.
2. Gently and slowly apply 3-5 cm traction to the catheter to begin catheter retrieval.
3. Detachment of the catheter may be observed by visualizing the separation of the catheter between the distal and proximal marker bands.
4. In the event that the catheter does not detach, assess the following parameters by observing the distal shaft of the catheter:

- Vessel straightening
- Catheter tip releasing from embolic cast
- Traction on embolic cast

**NOTE:** Do not apply more than 20 cm of traction to catheter to minimize risk of catheter separation, proximal to the detachment zone.

5. Under some difficult clinical situations, it may be safer to leave a flow-directed catheter in the vascular system rather than risk rupturing the malformation and, consequently a hemorrhage, by exercising too much traction on an entrapped catheter.

This is accomplished by stretching the catheter and cutting the shaft near the entry point of vascular access allowing the catheter to remain in the artery.

If catheter breaks during removal, distal migration or coiling of the catheter may occur. Some day surgical resection should be considered to minimize risk of thrombosis.

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Micro Therapeutics, Inc.