# COAGULATION CONTROL I, II

## CONTROL DE COAGULACIÓN I, II

NIVEL I		NIVEL II	
COD 61007	COD 61012	COD 61008	COD 61013

### USO PREVISTO

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### INTENDED USE

The Coagulation Control is intended for quality control in the laboratory only and is supplied with suggested intervals of acceptable values. For *in vitro* professional use only in the clinical laboratory.

### CONTENTS

	COD 61007	COD 61012	COD 61008	COD 61013
Cl. Control	4 x 1 mL	2 x 1 mL	-	-
Cli. Control	-	-	4 x 1 mL	2 x 1 mL

### COMPOSITION

**Coagulation Control.** 4 x 1.0 mL. Lyophilized human plasma.

*Components from human origin have been tested and found to be negative for the presence of antibodies anti-HIV and anti-HCV, as well as for HBs antigen. However, they should be handled cautiously as potentially infectious.*

### PREPARATION AND USE

1. Open a vial very carefully, avoiding any loss of the lyophilized material.
2. Pipette exactly 1.00 mL of distilled water into the vial. The component values depend on the accuracy of this reconstitution step.
3. Close the vial with the stopper and let stand for 30 minutes at room temperature.
4. Swirl gently, avoiding the formation of foam, to ensure complete dissolution of contents.
5. The reconstituted control plasma is to be treated like the patient plasma.

### STORAGE AND STABILITY

Store at 2-8°C.

The lyophilized Coagulation Control is stable until the expiration date given in the label.

The components of the reconstituted Coagulation Control are stable for 4 hours at 22°C.

Reconstituted Coagulation Control is stable for at least 8 days at -20°C. Freeze only once.

### WARNINGS AND PRECAUTIONS

Exercise the normal precautions required for handling all laboratory reagents. Safety data sheet available for professional user on request. Disposal of all waste material should be in accordance with local guidelines. Any serious incident that might occur in relation to the device shall be reported to BioSystems S.A.

### SPECIFIED VALUES

The assigned concentration values for components are batch-dependent and their traceability is shown in the enclosed value sheets. Traceability of the results can be assured only if the BioSystems reagents and recommended measurement procedures are used.

The intervals of suggested acceptable values have been calculated from previous experience in interlaboratory variability and are given for orientation only; each laboratory should establish its own precision parameters.

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