NanoCross™ Elite

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Instructions for use

Device description

Dilatation balloon catheters are used to exert radial force to dilate narrow vessel segments. The NanoCross™ Elite 0.36 mm (0.014 in) over-the-wire PTA balloon dilatation catheter is an over-the-wire (OTW) coaxial lumen catheter (1). A semi-compliant inflatable balloon (2) is mounted at the distal end. The distal atraumatic tip is tapered (3). The distal portion of the catheter has a lubricious coating (4).

The luer marked THRU (6) on the manifold (5) is the proximal opening for the central lumen. The central lumen ends at the distal tip of the catheter. This lumen is used to pass the catheter over a guidewire with a maximum diameter of 0.36 mm (0.014 in). The luer marked BALLOON (7) is used to inflate and deflate the dilatation balloon, using a mixture of contrast medium and saline solution. The balloon has 2 radiopaque marker bands (8) for positioning the balloon relative to the stenosis. The radiopaque marker bands indicate the dilating or working section of the balloon (9). On the 150 mm and longer balloons, 2 additional radiopaque marker bands (11) denote the middle of the balloon body.


Figure 1: 0.36 mm (0.014 in) NanoCross Elite OTW PTA balloon dilatation catheter

The NanoCross Elite catheter is available in multiple balloon sizes. Nominal balloon diameter and length are printed on the strain relief (10). Refer to the package label for additional information about catheter length and sheath compatibility.

All balloons distend to sizes above the nominal size at pressures greater than the nominal pressure. Consult the balloon compliance card that is packaged with the device for the diameters of the balloons at given pressures.

Indications for use

The NanoCross™ Elite 0.36 mm (0.014 in) OTW PTA balloon dilatation catheter is intended to dilate stenoses in the iliac, femoral, ilio-femoral, popliteal, infra-popliteal, and renal arteries, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. This device is also indicated for stent post-dilatation in the peripheral vasculature.

Contraindications

The NanoCross Elite catheter is contraindicated for use in the coronary arteries or in the neurovasculature, or when unable to cross the target lesion with a guidewire.

Warnings

- This device is supplied **sterile**, for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risks of patient infection and compromised device performance.
- This device should only be used by physicians who have experience with and who understand the clinical and technical aspects of percutaneous transluminal angioplasty.
- To reduce the potential for vessel damage, select a balloon size with an inflation diameter that approximates the diameter of the lumen at the intended inflation site.
- Before beginning a percutaneous transluminal angioplasty procedure, identify allergic reactions to contrast media.
- Use the catheter with caution for procedures involving calcified lesions or synthetic vascular grafts.
- Use direct fluoroscopic observation while manipulating the catheter in the vascular system.
- If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the device or the vessel lumen. If resistance is felt, carefully withdraw the dilatation catheter.
- Never use air or any gas to inflate the balloon.
- Do not manipulate the balloon when it is inflated.
- Consider precautions to prevent or reduce clotting. Physician experience and discretion determines the appropriate anticoagulation therapy for each patient.

Precautions

- Inspect the packaging and the device before use. Do not use the device if the packaging or the device is damaged.
- Do not use contrast media with this device that are contraindicated for intravascular use.
- Do not wipe the catheter surface with dry gauze.
- Do not attempt to pass the catheter through an introducer sheath that is smaller than the size indicated on the product label. Refer to the product labeling for sizing information.

- To avoid kinking, advance the dilatation catheter slowly, in small increments, until the proximal end of the guidewire emerges from the catheter.
- Ensure that the guidewire is in place before changing the position of the balloon catheter.
- Balloons should not be inflated in excess of the rated burst pressure.
- Larger and longer models of the catheter may exhibit longer deflation times, particularly on long catheter shafts.
- Do not retract the balloon catheter after stent postdilatation unless the balloon is free from the stent and fully deflated under vacuum.
- This device may be a biohazard after use. Handle and dispose of the device in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

Potential adverse events

The following complications may result from balloon dilatation procedures. Complications or adverse events may not be limited to this list:

- Allergic reaction to device materials or procedure medications
- Aneurysm
- Arrhythmia
- Arteriovenous fistula
- Artery injury (such as dissection, perforation, or rupture)
- Bleeding requiring transfusion
- Contrast medium reaction or renal failure
- Death
- Embolism
- Endocarditis
- Fever
- Hematoma
- Hemorrhage
- Hypertension or hypotension
- Infection
- Inflammation
- Intraluminal thrombus
- Myocardial infarction
- Pseudoaneurysm
- Sepsis
- Shock
- Stroke
- Transient ischemic attack
- Venous thromboembolism

How supplied

Warning: This device is supplied **sterile**, for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risks of patient infection and compromised device performance.

- The catheter is packaged with a protective sheath positioned over the balloon.
- The catheter is packaged with a protective wire that is positioned in the central lumen.

Storage

Store the catheter away from sunlight at room temperature in a dry place. Do not store catheters where they are directly exposed to organic solvents, ionizing radiation, or ultraviolet light. Rotate inventory so that the device is used before the Use-by date on the package label.

Device selection

Consider the following details when selecting a catheter:

- The diameter of the expanded balloon should not exceed the diameter of the artery immediately distal to or proximal to the stenosis.

Warning: To reduce the potential for vessel damage, select a balloon size with an inflation diameter that approximates the diameter of the lumen at the intended inflation site.

- Verify that the selected accessories accommodate the balloon catheter based on sizing information on the product label. Refer to the product labeling for size information.

Caution: Inspect the packaging and the device before use. Do not use the device if the packaging or the device is damaged.

Recommended items

Prepare the following items using sterile technique:

- 10 mL syringe filled with sterile heparinized saline
- Three-way stopcock
- Contrast medium (The standard inflation medium is a 1:1 mixture of contrast medium and sterile saline.)

Caution: Do not use contrast media with this device that are contraindicated for intravascular use.

- Appropriately sized exchange guidewire (refer to product labeling)
- Appropriately sized hemostatic introducer sheath (refer to product labeling)
- Inflation device with manometer

Preparation

1. Carefully remove the catheter from the inner package.

2. Remove and discard the protective sheath from the balloon and the protective wire from the guidewire lumen.

Caution: Do not wipe the catheter surface with dry gauze.

To displace air:

Note: As supplied, the balloon and the guidewire lumen of the catheter contain air. This air must be displaced to make certain that only liquid fills the balloon while the catheter is in the bloodstream.

1. Flush out the lumen for the guidewire with heparinized saline solution through the luer wire, marked THRU.
2. Mix equal volumes of contrast medium and normal saline.
3. Fill a 10 mL syringe with approximately 4 mL of the mixed contrast medium solution.

Warning: Before beginning a percutaneous transluminal angioplasty procedure, identify allergic reactions to contrast media.

Caution: Do not use contrast media with this device that are contraindicated for intravascular use.

4. Evacuate air from the balloon segment:
 - a. Attach the syringe to the luer marked BALLOON. Hold the syringe pointing down with the plunger pointed upward.
 - b. Apply negative pressure and aspirate the lumen for 15 seconds. Slowly release the pressure to neutral, so the contrast medium fills the shaft of the dilatation catheter.
 - c. Disconnect the syringe from the BALLOON luer of the dilatation catheter.
 - d. Remove all air from the syringe barrel. Reconnect the syringe to the BALLOON luer of the dilatation catheter. Maintain negative pressure on the balloon until air no longer returns to the device.
 - e. Slowly release the device pressure to neutral.

Note: On first use, avoid any positive pressure to avoid unfolding and partially inflating the balloon outside of the vessel.
5. Repeat steps 4b to 4e as necessary to remove all air from the system.
6. Replace the syringe with a manometer controlled-dilatation system. Do not introduce air into the system.

Insertion and dilatation

1. If applicable, use an appropriately sized introducer sheath to introduce the catheter percutaneously. Refer to the product labeling for sizing information.

Caution: Do not attempt to pass the catheter through an introducer sheath that is smaller than the size indicated on the product label. Refer to the product labeling for sizing information.

Warning: Use the catheter with caution for procedures involving calcified lesions or synthetic vascular grafts.

2. After the guidewire is in place through the lesion, place the prepared catheter over the guidewire and advance the catheter tip to the introduction site. Use a guidewire that is of suitable length to maintain control and position of the guidewire at all times.

Warning: Use direct fluoroscopic observation while manipulating the catheter in the vascular system.

Warning: If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the device or the vessel lumen. If resistance is felt, carefully withdraw the dilatation catheter.

Caution: To avoid kinking, advance the dilatation catheter slowly, in small increments, until the proximal end of the guidewire emerges from the catheter.

Note: To improve catheter movement, thoroughly wipe the guidewire with saline soaked gauze prior to using the guidewire in the catheter.

3. Position the catheter with the center of the balloon in the stenotic area of the vessel. The radiopaque marker bands indicate the working length of the balloon. On the 150 mm and longer balloons, 2 additional marker bands denote the middle of the balloon body.

Caution: Ensure that the guidewire is in place before changing the position of the balloon catheter.

4. When the radiopaque marker bands are appropriately positioned, inflate the balloon to dilate the target area.

Warning: Never use air or any gas to inflate the balloon.

Warning: Do not manipulate the balloon when it is inflated.

Caution: Balloons should not be inflated in excess of the rated burst pressure.

5. Using the inflation device, deflate the balloon by aspirating the catheter.

Caution: Larger and longer models of the catheter may exhibit longer deflation times, particularly on long catheter shafts.

6. Withdraw the catheter while maintaining a vacuum in the balloon.

Caution: Do not retract the balloon catheter after stent post dilatation unless the balloon is free from the stent and fully deflated under vacuum.

Note: Multiple balloon inflations and deflations can cause resistance during device withdrawal. If resistance is encountered, use a gentle clockwise twisting motion to ease withdrawal through the sheath. If the balloon cannot be withdrawn through the sheath, withdraw the catheter and sheath together, while maintaining vessel access with the guidewire.

7. Use angiography to verify dilatation results.

Disposal

Caution: This device may be a biohazard after use. Handle and dispose of the device in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

Warranty

This product is subject to standard Medtronic warranty terms. Please go to www.medtronic.com for standard terms and conditions of supply for Australia and New Zealand.

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©2017 Medtronic. All rights reserved. Medtronic and Medtronic logo are trademarks of Medtronic. ™* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
ev3 Inc.
Medtronic Australasia Pty Ltd