# Infusomat® Space® Line

## Infusomat® Space® Line SafeSet



### Infusomat® plus Line

### Infusomat® plus Line SafeSet



<b>en</b>	Instructions for use	<b>6</b>
<b>de</b>	Gebrauchsanweisung	<b>7</b>
<b>bg</b>	Инструкция за употреба	<b>8</b>
<b>cs</b>	Návod k použití	<b>9</b>
<b>da</b>	Brugsanvisning	<b>10</b>
<b>el</b>	Οδηγίες χρήσης	<b>11</b>
<b>es</b>	Instrucciones de uso	<b>12</b>
<b>et</b>	Kasutusjuhised	<b>13</b>
<b>fi</b>	Käyttöohje	<b>14</b>
<b>fr</b>	Mode d'emploi	<b>15</b>
<b>hr</b>	Upute za uporabu	<b>16</b>
<b>hu</b>	Használati utasítás	<b>17</b>
<b>id</b>	Petunjuk penggunaan	<b>18</b>
<b>it</b>	Istruzioni per l'uso	<b>19</b>
<b>kk</b>	Қолдану жөніндегі нұсқаулық	<b>20</b>
<b>ko</b>	사용 지침	<b>21</b>
<b>lt</b>	Naudojimo instrukcija	<b>22</b>
<b>lv</b>	Lietošanas instrukcija	<b>23</b>
<b>nl</b>	Gebruikersinformatie	<b>24</b>
<b>no</b>	Bruksanvisning	<b>25</b>
<b>pl</b>	Instrukcja użytkowania	<b>26</b>
<b>pt</b>	Instruções de utilização	<b>27</b>
<b>ro</b>	Instrucțiuni de utilizare	<b>28</b>
<b>ru</b>	Инструкция по применению	<b>29</b>
<b>sk</b>	Návod na použitie	<b>30</b>
<b>sl</b>	Navodila za uporabo	<b>31</b>
<b>sr</b>	Uputstvo za upotrebu	<b>32</b>
<b>sv</b>	Bruksanvisning	<b>33</b>
<b>th</b>	คำแนะนำสำหรับการใช้	<b>34</b>
<b>tr</b>	Kullanım Kilavuzu	<b>35</b>
<b>uk</b>	Інструкція з використання	<b>36</b>
<b>vi</b>	Hướng dẫn sử dụng	<b>37</b>
<b>zh</b>	使用说明	<b>38</b>

**B | BRAUN**

**Infusomat<sup>®</sup> Space<sup>®</sup> Line**  
**Infusomat<sup>®</sup> Space<sup>®</sup> Line SafeSet**

**Infusomat<sup>®</sup> plus Line**  
**Infusomat<sup>®</sup> plus Line SafeSet**

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

# Infusomat<sup>®</sup> plus Line

## en Instructions for Use

Infusomat® Lines can be used with devices which are compatible or in conformance with Luer connectors (e.g., ISO 80369-7) and compatible B. Braun infusion pumps.

### Sterility

Sterilized by EO. Please refer to the information on primary packaging.

Do not use if package is damaged. Visually inspect the sterile barrier system for breaches and integrity prior to use.

### Intended purpose

I.V. administration set. For infusion by gravity and compatible pumps. Single-use, do not re-use.

### Indication

Type IV-Standard

Indications: For standard infusion solutions and parenteral fluids (e.g., parenteral infusion solutions and lipids), and for low soluble IV-solutions and for cytostatic solutions.

Type Light Protect

Indications: For light sensitive drugs

### Patient population

The device can be used for all patients for which infusion therapy is prescribed. Infusomat® Lines can be used for adults, pediatric and neonates. Selection of pump and dedicated sets according to national standards and/or hospital protocols.

### Intended User

Infusomat® Lines should only be used by health-care professionals (professional nurses, licensed and/or certified physicians, healthcare workers, midwives, paramedics, pharmacists) who have been adequately trained within the educational background in this technique.

Upon medical assessment and adequate instruction, patients and/or caregivers may also be allowed to take over defined handling steps according to national guidelines. Healthcare workers are responsible for ensuring that patients and/or caregivers are instructed in the correct use of Infusomat® Lines and that the instruction of caregivers and/or patients is documented by the healthcare workers according to local documentation guidelines.

### Contraindications

IV administration sets shall not be used in case of proven incompatibility of the drug/IV solution to one of the used materials. SPCs of the drugs/solutions for contraindications should be considered. No administration of blood products, blood components and enteral nutrition solutions.

### Residual Risks/Side effects

General risks occurring during infusion therapy including air embolism, local- to systemic infection, particulate contamination, medication error including over and under supply. Kinking of the line may result in decrease or stop of infusion flow. Leakage of drugs/fluids can be an additional hazard to patients, visitors or staff.

### Precautions

- Visual check of the IV administration sets for damages before use.
- Do not re-attach the protective cap of the patient connector.
- Always refer to the manufacturer's instruction for use of infusion pumps.

- If pressure infusion without infusion pump is used, remove air from container beforehand. Observe general safety precautions.
- Mind the general guidelines for parallel infusions.
- Please also consider the corresponding SPCs of the drugs/ solutions regarding possible incompatibilities of drugs/ solutions and/or material of the used Infusomat® Line
- Monitor all connections for tightness over the whole infusion procedure
- Due to the structure of the fluid filter of Infusomat® Line SafeSet which achieves a particle retention with a filter efficiency > 80% for a particle size of 3 µm, the infusion of highly viscous solutions, lipid emulsions > 10% and certain suspensions can affect the flow rate. With regard to this, please also follow the information in the relevant summary of product characteristics.
- Use aseptic techniques in accordance to national guidelines and/or hospital protocols.

### Warning

Do not use if packaging or product is damaged or contaminated or if protective caps are loose or missing. Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

Administration of light sensitive drugs requires UV-protected IV administration sets. Always check, if incompatibilities with other drugs or with respective material of different types of Infusomat® Lines is stated in the SPC (Summary of product characteristics) of the drug. Always check, if the use of a 0.2 µm filter is stated in the SPC (Summary of product characteristics) of the drug.

### Operating Instructions

- ☉ Illustrations on page 4–5.

### Duration of use

Infusomat® Line Type IV should be changed acc. to national standards and / or hospital protocols, but latest after 96 h (24 h for sets with 0.2 µm Infusion filter).

### Disposal

Disposal acc. to local guidelines and or clinical protocols.

Use integrated spike protection at flow regulator housing for disposal.

### Storage and handling conditions

Storage conditions:

Keep dry.

Keep away from sunlight.

### Notice to the user

If, during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

If further Instructions for Use are required, they can be requested from the manufacturer or can be obtained from the B. Braun homepage: <https://eifu.bbraun.com/>.

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**AU** B. Braun Australia Pty. Ltd.  
Level 5, 7 – 9 Irvine Place  
Bella Vista NSW 2153  
Australia  
For sterilization type see primary  
packaging

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**ID** Distributed by:  
PT. B. Braun Medical Indonesia  
Jakarta, Indonesia

**HN** AGENCIA MATAMOROS S. de R.L  
For manufacturing site refer  
to packaging

**MY** Authorised representative:  
B. Braun Medical Industries  
Sdn. Bhd.  
Bayan Lepas Free Industrial Zone,  
11900 Penang, Malaysia.  
www.bbraun.com

**RS** Ovlašćeni predstavnik proizvođača:  
B. Braun Adria RSRB d.o.o. Beograd  
Milutina Milankovića 11g, Beograd

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B. Braun Melsungen AG