EIAgen

HBsAg Kit

REF	071011	 96
REF	071012	 192
REF	071015	 480

This package insert must be read carefully before product use.

Package insert instructions must be carefully followed.

Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package.

Manufacturer:

Adaltis S.r.l

Via Durini, 27

20122 Milano (Italy)

Tel. +39-0774-5791 - Fax +39-0774-353085

www.adaltis.net

en

SYMBOLS USED ON LABELS

English EN							
	In Vitro Diagnostic Medical Device	Catalogue Number	Lot Number	Attention, See Instructions For Use	Temperature Limitation	Use By	Number of Test
							
	Manufacturer	Keep away from Sunlight	Date of Manufacture	Biological Risk	Microplate	Positive Control	Negative Control
							
	Conjugate	Assay Diluent	Substrate A (Urea Peroxide)	Substrate B (TMB)	Stop Solution (0,5 M H ₂ SO ₄)	Wash Buffer Conc. (20x)	Danger
							
	Warning						

Attention:

Negative Control, Positive Control, Conjugate, Assay Diluent and Substrate Solution A classified as: Skin Sens. 1 and Aquatic Chronic 3

- Signal word:**
Warning
- Hazard-determining components of labelling:**
Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1)
- Hazard statements:**
H317 May cause an allergic skin reaction.
H412 Harmful to aquatic life with long lasting effects.
- Precautionary statements:**
P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
P273 Avoid release to the environment.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P362+P364 Take off contaminated clothing and wash it before reuse.
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
P321 Specific treatment (see medical advice on this label).
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Attention:

Substrate Solution B classified as: Carc. 1B and Repr. 1B

- Signal word:**
Danger
- Hazard-determining components of labelling:**
N,N-dimethylformamide
- Hazard statements:**
H350 May cause cancer.
H360D May damage the unborn child.
- Precautionary statements:**
P201 Obtain special instructions before use.
P202 Do not handle until all safety precautions have been read and understood.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P308+P313 IF exposed or concerned: Get medical advice/attention.
P405 Store locked up.
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Refer to www.adaltis.net for the Safety Data Sheets.

A. INTENDED USE

EIAgen HBsAg Kit is an enzyme-linked immunosorbent assay (ELISA) for qualitative detection of HBsAg in human serum or plasma (EDTA, sodium citrate or heparin). It is intended for screening of blood donors and for diagnosing of patients related to infection with hepatitis B virus.

For "in vitro" diagnostic use only.

B. INTRODUCTION

Hepatitis B virus (HBV) is an enveloped, double-stranded DNA virus belonging to the Hepadnaviridae family and is recognized as the major cause of blood transmitted hepatitis together with hepatitis C virus (HCV). Infection with HBV induces a spectrum of clinical manifestations ranging from mild, inapparent disease to fulminant hepatitis, severe chronic liver diseases, which in some cases can lead to cirrhosis and carcinoma of the liver. Classification of a hepatitis B infection requires the identification of several serological markers expressed during three phases (incubation, acute and convalescent) of the infection. Now several diagnostic test are used for screening, clinical diagnosis and management of the disease. Hepatitis B surface antigen or HBsAg, previously described as Australia antigen, is the most important protein of the envelope of Hepatitis B Virus. The surface antigen contains the determinant "a", common to all known viral subtypes and immunologically distinguished in two distinct subgroups (ay and ad). HBV has 10 major serotypes and four HBsAg subtypes have been recognized (adw, ady, ayw, and ayr). HBsAg can be detected 2 to 4 weeks before the ALT levels become abnormal and 3 to 5 weeks before symptoms develop. The serological detection of HBsAg is a powerful method for the diagnosis and prevention of HBV infection and ELISA has become an extensively used analytical system for screening of blood donors and clinical diagnosis of HBV in infected individuals.

C. PRINCIPLE OF THE TEST

For detection of HBsAg, EIAgen HBsAg Kit uses antibody "sandwich" ELISA method in which, polystyrene microwell strips are pre-coated with monoclonal antibodies specific to HBsAg. Patient's serum or plasma sample is added to the microwells. During incubation, the specific immunocomplex formed in case of presence of HBsAg in the sample, is captured on the solid phase. Then the second antibody conjugated the enzyme horseradish peroxidase (the HRP-Conjugate) directed against a different epitope of HBsAg is added into the wells. During the second incubation step, these HRP-conjugated antibodies will be bound to any anti-HBs-HBsAg complexes previously formed during the first incubation, and the unbound HRP-conjugate is then removed by washing. Chromogen solutions containing tetramethyl-benzidine (TMB) and urea peroxide are added to the wells. In presence of the antibody-antigen-antibody(HRP) "sandwich" immunocomplex, the colorless Chromogens are hydrolyzed by the bound HRP-conjugate to a blue-colored product. The blue color turns yellow after stopping the reaction with sulfuric acid. The amount of color intensity can be measured and it is proportional to the amount of antigen captured in the

wells, and to its amount in the sample respectively. Wells containing samples negative for HBsAg remain colorless.

D. COMPONENTS

The kit contains reagents for 96 tests (code 071011), 192 tests (code 071012), or 480 tests (code 071015).

Microplate	1
Negative Control	1 x 1.5 mL/vial
Positive Control	1 x 1.5 mL/vial
Conjugate	1 x 6 mL/vial
Assay Diluent	1 x 5 mL/vial
Substrate Solution A (Urea Peroxide)	1 x 6 mL/vial
Substrate Solution B (TMB)	1 x 6 mL/vial
Stop Solution (H₂SO₄ - 0.5M)	1 x 6 mL/vial
Wash Buffer Concentrate 20X	1 x 30 mL/vial
Plate Sealing Foils	2
Plastic Sealable Bag	1
Number of tests	96
Code	071011

Microplate	2
Negative Control	1 x 1.5 mL/vial
Positive Control	1 x 1.5 mL/vial
Conjugate	2 x 6 mL/vial
Assay Diluent	2 x 5 mL/vial
Substrate Solution A (Urea Peroxide)	2 x 6 mL/vial
Substrate Solution B (TMB)	2 x 6 mL/vial
Stop Solution (H₂SO₄ - 0.5M)	2 x 6 mL/vial
Wash Buffer Concentrate 20X	2 x 30 mL/vial
Plate Sealing Foils	3
Plastic Sealable Bag	2
Number of tests	192
Code	071012

Microplate	5
Negative Control	4 x 1.5 mL/vial
Positive Control	4 x 1.5 mL/vial
Conjugate	5 x 6 mL/vial
Assay Diluent	4 x 5 mL/vial
Substrate Solution A (Urea Peroxide)	5 x 6 mL/vial
Substrate Solution B (TMB)	5 x 6 mL/vial
Stop Solution (H₂SO₄ - 0.5M)	1 x 30 mL/vial
Wash Buffer Concentrate 20X	3 x 50 mL/vial
Plate Sealing Foils	6
Plastic Sealable Bag	5
Number of tests	480
Code	071015

1. Microplate

Blank microwell strips fixed on white strip holder.

12 strips of 8 microwells coated with monoclonal antibodies reactive to HBsAg (anti-HBs).

Plates are sealed into a aluminium pouch with desiccant.

The microwell strips can be broken to be used separately. Place unused wells or strips in the provided plastic sealable storage bag together with the desiccant and return to 2...8°C. Once open, stable for one month at 2...8°C.

2. Negative Control

Yellowish liquid filled in a vial with natural screw cap.

Protein-stabilized buffer tested non-reactive for HBsAg.

Ready to use as supplied. It contains 0,1% Proclin™ 300 as preservative.

Once open, stable for one month at 2...8°C.

3. Positive Control

Red-colored liquid filled in a vial with red screw cap.
HBsAg diluted in protein-stabilized buffer.
Ready to use as supplied. It contains 0,1% Proclin™ 300 as preservative.

Once open, stable for one month at 2...8°C.

Important Note: *The absence of viable pathogens in the Positive Control can not be fully ensured, and therefore, the reagent should be handled as potentially biohazardous, in accordance with good laboratory practices.*

4. Conjugate

Red-colored liquid in a white vial with green screw cap.
Horseradish peroxidase-conjugated anti-HBs antibodies.
Ready to use as supplied. It contains 0,1% Proclin™ 300 as preservative.

Once open, stable for one month at 2...8°C.

5. Assay Diluent

Green-colored in a vial with pink screw cap.
Buffer solution containing protein.
Ready to use as supplied. It contains 0,1% Proclin™ 300 as preservative.

Once open, stable for one month at 2...8°C.

6. Substrate Solution "A"

Colorless liquid filled in a white vial with grey screw cap.
Urea peroxide solution.
Ready to use as supplied.

Once open, stable for one month at 2...8°C.

7. Substrate Solution "B"

Colorless liquid filled in a black vial with black screw cap.
TMB (Tetramethyl benzidine solution) and N,N-dimethylformamide.

Ready to use as supplied.

Once open, stable for one month at 2...8°C.

Note: *To be stored protected from light as sensitive to strong illumination.*

8. Stop Solution

Colorless liquid in a white vial with red screw cap.
Diluted sulfuric acid solution (0.5M H₂SO₄).
Ready to use as supplied.

Once open, stable for one month at 2...8°C.

9. Wash Buffer Concentrate 20x

Colorless liquid filled in a clear bottle with natural screw cap. Buffer solution containing surfactant Tween-20.
The concentrate must be diluted **1 to 20** with distilled/deionized water before use.

Once diluted, stable for one week at room temperature, or for two weeks when stored at 2...8°C.

E. MATERIALS REQUIRED BUT NOT PROVIDED

1. Calibrated Micropipettes (20, 50, 100 µL) and disposable plastic tips.
2. EIA grade water (bidistilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.

5. Calibrated ELISA microplate thermostatic incubator capable to provide a temperature of 37 ± 0.5°C.
6. Calibrated ELISA microwell reader with 450nm (reading) and possibly with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
This package insert must be read carefully before product use.
2. Read carefully the Safety Data Sheet (SDS) before product use.
3. When the kit is used for the screening of blood units and blood components, it has to be used in a laboratory certified and qualified by the national authority in that field (Ministry of Health or similar entity) to carry out this type of analysis.
4. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
5. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
6. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-born microbial agents, when opening kit vials and microplates and when performing the test. Protect the Substrate Solution "B" (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.
7. Upon receipt, store the kit at 2...8°C into a temperature controlled refrigerator or cold room.
8. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
9. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.
10. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
11. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.
12. Do not use the kit after the expiration date stated on the external container and internal (vials) labels.
13. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical

Laboratories", ed. 1984. Materials from human origin may have been used in the preparation of the Negative Control of the kit. These materials have been tested with tests kits with accepted performance and found negative for antibodies to HIV 1/2, HCV, TP and HBsAg. However, there is no analytical method that can assure that infectious agents in the specimens or reagents are completely absent. Therefore, handle reagents and specimens with extreme caution as if capable of transmitting infectious diseases. Bovine derived sera have been used for stabilizing of the positive and negative controls. Bovine serum albumin (BSA) and fetal calf sera (FCS) are derived from animals from BSE/TSE free-geographical areas.

14. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.
15. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated before waste. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min.
16. Accidental spills from samples and operations have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
17. The Stop Solution contains 0.5M sulphuric acid. Avoid contact with skin and eyes. In the event of contact, rinse immediately with plenty of water.
18. ProClin™ 300, 0.1% used as preservative, can cause sensation of the skin. Wipe up spills immediately or wash with water if come into contact with the skin or eyes.
19. Do not smoke, eat, drink or apply cosmetics in areas in which specimens or kit reagents are handled
20. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.
21. Do not pipette by mouth.

G. SPECIMEN: PREPARATION AND RECOMMENDATIONS

1. Specimen Collection:

No special patient's preparation required. Collect the specimen in accordance with the normal laboratory practice. Either fresh serum or plasma specimens can be used with this assay. Blood collected by venipuncture should be allowed to clot naturally and completely - the serum/plasma must be separated from the clot as early as possible as to avoid haemolysis of the RBC. Care should be taken to ensure that the serum specimens are clear and not contaminated by microorganisms. Any visible particular matters in the specimen should be

removed by centrifugation at 3000 RPM (round per minutes) for 20 minutes at room temperature or by filtration.

2. Plasma specimens collected into EDTA, sodium citrate or heparin may be tested, but **highly lipaemic, icteric, or hemolytic specimens should not be used** as they can give false results in the assay. **Do not heat inactivate specimens.** This can cause deterioration of the target analyte. Samples with visible microbial contamination should never be used.
3. EIAgen HBsAg Kit is intended ONLY for testing of individual serum or plasma samples. Do not use the assay for testing of cadaver samples, saliva, urine or other body fluids, or pooled (mixed) blood.
4. **Transportation and Storage:** Store specimens at 2...8°C. Specimens not required for assaying within 7 days should be stored frozen (-20°C or lower). Multiple freeze-thaw cycles should be avoided. For shipment, samples should be packaged and labeled in accordance with the existing local and international regulations for transportation of clinical samples and ethological agents.

H. PREPARATION OF COMPONENTS AND WARNINGS

The components of the kit will remain stable through the expiration date indicated on the label and package when stored between 2...8°C, do not freeze. To assure maximum performance of EIAgen HBsAg Kit, during storage, protect the reagents from contamination with microorganism or chemicals.

1. Microplates:

Allow the microplate to reach room temperature (18...30°C), about 1 hr, before opening the pouch. Unused strips have to be placed inside the plastic sealable bag, with the desiccant supplied and stored at 2...8°C. After first opening, remaining strips are stable one month at 2...8°C.

2. Negative Control:

Ready to use. Mix well on vortex before use.

3. Positive Control:

Ready to use. Mix well on vortex before use. Handle this component as potentially infectious.

4. Conjugate:

Ready to use. Mix well on vortex before use. Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes. If this component has to be transferred use only plastic, possibly sterile disposable containers.

5. Assay Diluent:

Ready to use. Mix well on vortex before use.

6. Substrate Solution "A" and "B":

Ready to use. Mix well on vortex before use. Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes. Do not expose to strong illumination, oxidizing agents and metallic surfaces. If this component has to be transferred use only plastic, and if possible, sterile disposable container.

7. Stop Solution:

Ready to use. Mix well on vortex before use.

8. Wash Buffer Concentrate 20x (vial of 30 mL):

The whole content of the 20x concentrated solution has to be diluted with distilled/deionized water up to 600 mL (up to 1000 mL for the vial of 50 mL), the volume is reported on the label, and mixed gently end-over-end before use. As some salt crystals may be present into the vial, take care to dissolve all the content when preparing the solution. During preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles.

Note: Once diluted, the wash solution is stable for 1 week at room temperature and two weeks at 2...8°C.

I. INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT

1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample or the components of the kit. They should also be regularly maintained in order to show a precision of 1% and a trueness of $\pm 2\%$. Decontamination of spills or residues of kit components should also be carried out regularly.
2. The ELISA incubator has to be set at $+37^{\circ}\text{C}$ (tolerance of $\pm 0.5^{\circ}\text{C}$) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.

3. INSTRUCTIONS FOR WASHING

The ELISA washer is extremely important to the overall performances of the assay:

- A good washing procedure is essential in order to obtain correct and precise analytical data.
- It is therefore, recommended to use a good quality ELISA microplate washer, maintained at the best level of washing performances. In general, no less than 5 automatic washing cycles of 350-400 μL /well are sufficient to avoid false positive reactions and high background.
- To avoid cross-contaminations of the plate with specimen or HRP-conjugate, after incubation, do not discard the content of the wells but allow the plate washer to aspirate it automatically.
- Assure that the microplate washer liquid dispensing channels are not blocked or contaminated and sufficient volume of Wash buffer is dispensed each time into the wells.
- In case of manual washing, we suggest to carry out **5 washing cycles**, dispensing 350-400 μL /well and aspirating the liquid for **5** times. If poor results (high background) are observed, increase the washing cycles or soaking time per well.
- In any case, the liquid aspirated out the strips should be treated with a sodium hypochlorite solution at a final concentration of 2.5% for 24 hours, before they are wasted in an appropriate way.
- The concentrated Wash buffer should be diluted

1:20 before use. If less than a whole plate is used, prepare the proportional volume of solution.

4. The ELISA microplate reader has to be equipped with a reading filter of 450nm and ideally with a second filter (630nm) for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0 ; (c) linearity to ≥ 2.0 ; repeatability $\geq 1\%$. Blanking is carried out on the well identified in section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer's instructions.
5. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, shaking, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in section O "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing samples and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells due to strongly reactive samples, leading to false positive results. The use of ELISA automated work stations is recommended for blood screening and when the number of samples to be tested exceed 20-30 units per run. If using fully automated equipment, during incubation, do not cover the plates with the plate cover. The tapping out of the remainders inside the plate after washing, can also be omitted.

L. PRE ASSAY CONTROLS AND OPERATIONS

1. Check the expiration date of the kit printed on the external label of the kit box. Do not use if expired.
2. Check that the liquid components are not contaminated by naked-eye visible particles or aggregates. Check that the Substrate Solutions "A" and "B" are colorless or pale blue by aspirating a small volume of it with a sterile transparent plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box. Check that the aluminum pouch, containing the microplate, is not punctured or damaged.
3. Dilute all the content of the 20x concentrated Wash Buffer as described above.
4. Allow all the other components to reach room temperature ($18...30^{\circ}\text{C}$), about 1 hr and then mix as described.
5. Set the ELISA incubator at $+37^{\circ}\text{C}$ and prepare the ELISA washer by priming with the diluted wash solution, according to the manufacturers instructions. Set the right number of washing cycles as found in section I.3.
6. Check that the ELISA reader has been turned on at least 20 minutes before reading.

- If using an automated workstation, turn it on, check settings and be sure to use the right assay protocol.
- Check that the micropipettes are set to the required volume.
- Check that all the other equipment is available and ready to use.
- In case of problems, do not proceed further with the test and advise the supervisor.

M. ASSAY PROCEDURE

The assay has to be carried out according to what reported below, taking care to maintain the same incubation time for all the samples in testing.

MANUAL ASSAY:

- Preparation:** Mark three wells as Negative control (e.g. B1, C1, D1), two wells as Positive control (e.g. E1, F1) and one Blank (e.g. A1, neither samples nor Conjugate should be added into the Blank well). If the results will be determined by using dual wavelength plate reader, the requirement for use of Blank well could be omitted. Use only number of strips required for the test.
- Adding Diluent:** Add 20µL of Assay Diluent into each well except the Blank.
- Adding Sample:** Add 100µL of Positive control, Negative control, and Specimen into their respective wells except the Blank. **Note: Use a separate disposal pipette tip for each specimen, Negative Control, Positive Control to avoid cross-contamination. Mix by tapping the plate gently.**
- Incubating:** Cover the plate with the plate cover and incubate for **60 minutes at 37°C**.
- Adding Conjugate:** At the end of the incubation, remove and discard the plate cover. Add 50µL Conjugate into each well except the Blank, and mix by tapping the plate gently.
- Incubating:** Cover the plate with the plate cover and incubate for **30 minutes at 37°C**.
- Washing:** At the end of the incubation, remove and discard the plate cover. Wash each well **5 times** with diluted Washing buffer. Each time allow the microwells to soak for **30-60 seconds**. After the final washing cycle, turn down the plate onto blotting paper or clean towel and tap it to remove any remainders (see section I.3).
- Coloring:** Add 50µL of Substrate Solution "A" and 50µL Substrate Solution "B" into each well including the Blank. Incubate the plate at **37°C for 30 minutes avoiding light**. The enzymatic reaction between the Substrate Solutions and the Conjugate produces blue color in Positive control and HBsAg positive sample wells.
- Stopping Reaction:** Using a multichannel pipette or manually, add 50µL Stop solution into each well and mix gently. Intensive yellow color develops in Positive control and HBsAg positive sample wells.
- Measuring the Absorbance:** Calibrate the plate reader with the Blank well and read the absorbance at **450nm**. If a dual filter instrument is used, set the reference wavelength at **630nm**. Calculate the Cut-off value and evaluate the results. (**Note:** read the absorbance within **10 minutes** after stopping the reaction).

Important notes:

- If the second filter is not available, ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
- Reading has should ideally be performed immediately after the addition of the Stop Solution but definitely no longer than 10 minutes afterwards. Some self oxidation of the substrate can occur leading to a higher background.

N. ASSAY SCHEME

Steps	Operations
Assay Diluent	20 µL
Controls and Samples	100 µL
1st incubation	60 min
Temperature	+37°C
Conjugate	50 µL
2nd incubation	30 min
Temperature	+37°C
Wash Step	5 cycles (see section I.3)
Substrate "A"	50 µL
Substrate "B"	50 µL
3rd incubation	30 min (avoiding light)
Temperature	+37°C
Stop Solution	50 µL
Reading OD	450/630nm

An example of dispensation scheme is reported below (valid for both incubation time procedures):

Microplate												
	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	S3										
B	NC	S4										
C	NC	S5										
D	NC	S6										
E	PC	S7										
F	PC	S8										
G	S1	S9										
H	S2	S10										

Legenda: BLK = Blank NC = Negative Control
PC = Positive Control S = Sample

O. INTERNAL QUALITY CONTROL

The test results are valid if the Quality Control criteria are fulfilled. It is recommended that each laboratory must establish appropriate quality control system with quality control material similar to or identical with the patient sample being analyzed.

- The A value of the Blank well, which contains only Chromogen and Stop Solution, is < 0.080 at 450 nm.
- The A values of the Positive Control must be ≥ 0.800 at 450/630nm or at 450nm after blanking.
- The A values of the Negative Control must be < 0.100 at 450/630nm or at 450nm after blanking.

If one of the Negative control A values does not meet the Quality Control criteria, it should be discarded and the mean value calculated again using the remaining two values. If more than one Negative control A values do

not meet the Quality Control Range specifications, the test is invalid and must be repeated.

If the results of the test match the requirements stated above, proceed to the next section.

If they do not, do not proceed any further and perform the following checks:

Problems	Check
Blank well ≥ 0.080 OD at 450nm	1. that the Substrate Solutions have not become contaminated during the assay
Negative Control (NC) ≥ 0.100 OD at 450/630nm or at 450nm after blanking	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of the negative one); 4. that no contamination of the negative control or of the wells where the control was dispensed has occurred due to spills of positive samples or of the conjugate; 5. that micropipettes have not become contaminated with positive samples or with the conjugate 6. that the washer needles are not blocked or partially obstructed.
Positive Control < 0.800 OD at 450/630nm or at 450nm after blanking	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the control (dispensation of negative control instead of positive control); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

If any of the above problems have occurred, report the problem to the supervisor for further actions.

P. RESULTS

Each microplate should be considered separately when calculating and interpreting the results of the assay, regardless of the number of plates concurrently processed. The results are calculated by relating each specimen absorbance (A) value to the Cut-off value (C.O.) of the plate. If the Cut-off reading is based on single filter plate reader, the results should be calculated by subtracting the Blank well A value from the print report values of specimens and controls. In case the reading is based on dual filter plate reader, do not subtract the Blank well A value from the print report values of specimens and controls.

$$\text{Cut-Off (C.O.)} = \text{NC mean} + 0.06$$

The value found for the test is used for the interpretation of results as described in the next paragraph.

An example of calculation is reported below:

Example:

1. Quality Control

Blank well A value: A1= 0.025 at 450nm (Note: blanking is required only when reading with single filter at 450nm)

Well No.: B1 C1 D1
Negative control A values after blanking: 0.020 0.012 0.016

Well No.: E1 F1
Positive control A values after blanking: 2.421 2.369
All control values are within the stated quality control range

$$2. \text{ Calculation of Nc: } = \frac{(0.020+0.012+0.016)}{3} = 0.016$$

$$3. \text{ Calculation of the Cut-off: (C.O.)} = 0.016 + 0.06 = 0.076$$

Q. INTERPRETATION OF RESULTS

Negative Results (A / C.O. < 0.9):

Specimens giving absorbance less than the Cut-off value are negative for this assay, which indicates that no hepatitis B virus surface antigen has been detected with EIAGEN HBsAg Kit, therefore the patient is probably not infected with HBV and the blood unit do not contain hepatitis B virus surface antigen and could be transfused in case that other infectious diseases markers are also absent.


Positive Results (A / C.O. > 1.1):

Specimens giving an absorbance equal to or greater than the Cut-off value are considered initially reactive, which indicates that hepatitis B virus surface antigen has probably been detected using EIAGEN HBsAg Kit. All initially reactive specimens should be retested in duplicates using EIAGEN HBsAg Kit before the final assay results interpretation. Repeatedly reactive specimens can be considered positive for hepatitis B virus surface antigen with EIAGEN HBsAg Kit.

Borderline (A / C.O. = 0.9-1.1):

Specimens with absorbance to Cut-off ratio between 0.9 and 1.1 are considered borderline and retesting of these specimens in duplicates is required to confirm the initial results.

Follow-up, confirmation and supplementary testing of any positive specimen with other analytical system (e.g. PCR) is required. Clinical diagnosis should not be established based on a single test result. It should integrate clinical and other laboratory data and findings.

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If, after retesting of the initially reactive samples, both wells are negative results (A/C.O.<0.9), these samples should be considered as non-repeatable positive (or false positive) and recorded as negative. As with many very sensitive ELISA assays, false positive results can occur due to the several reasons, most of which are connected with, but not limited to, inadequate washing step. For more information regarding "Troubleshooting Guide" (section S).

If after retesting in duplicates, one or both wells are positive results, the final result from this ELISA test should be recorded as repeatedly reactive. Repeatedly reactive specimens could be considered positive for hepatitis B virus surface antigen and therefore the patient is probably infected with HBV and the blood unit must be discarded.

After retesting in duplicates, samples with values close to the Cut-off value should be interpreted with caution and considered as "borderline" zone sample, or uninterpretable for the time of testing.

Important notes:

1. Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgment errors and misinterpretations.
2. When test results are transmitted from the laboratory to another department, attention must be paid to avoid erroneous data transfer.
3. Diagnosis of viral hepatitis infection has to be taken and released to the patient by a suitably qualified medical doctor.

R. PERFORMANCES

Evaluation of Performances has been conducted in accordance to what reported in the Common Technical Specifications or CTS (art. 5, Chapter 3 of IVD Directive 98/79/EC).

Evaluation studies carried out in Paul-Ehrlich-Institut (PEI), German Red Cross Institute Baden-Württemberg – Hessen, and three blood banks, demonstrated the following performance characteristics of EIAGEN HBsAg Kit.

1. SPECIFICITY

When evaluated on European blood donors (n=5038), the overall diagnostic specificity of the kit was 99.78%.

During multi-center evaluation (Site A, B and C), EIAGEN HBsAg Kit demonstrated specificity of 99.92%.

Laboratory	Number	EIAGEN HBsAg Kit		
		-	+	Specificity
"A" blood bank	1958	1955	3	99.85%
"B" blood bank	2518	2516	2	99.92%
"C" blood bank	6344	6340	4	99.94%
Total	10820	10811	9	99.92%

2. SENSITIVITY

EIAGEN HBsAg Kit was evaluated for sensitivity on 22 HBV commercial available HBV seroconversion panels, and on total 403 HBsAg positive including 146 HBsAg HBV genotyped and HBsAg subtyped plasma samples available at the Paul-Ehrlich-Institut. With respect to seroconversion sensitivity, the results for EIAGEN HBsAg Kit on the 22 HBV seroconversion panels showed a sensitivity level at least equivalent with the range of current CE marked HBsAg screening assays for which PEI holds data. 10 additional seroconversion panels were tested in-house. The seroconversion sensitivity was comparable to other CE-marked HBsAg screening test. With respect to diagnostic sensitivity EIAGEN HBsAg Kit detected all positive samples as positive, including the HBV genotypes A-F or HBsAg subtypes examined.

In conclusion, the overall score of EIAGEN HBsAg Kit for the seroconversion sensitivity was comparable with other CE marked HBsAg test kits for which PEI holds data and all 403 HBsAg positive samples were reactive giving an overall sensitivity of 100%.

3. ANALYTICAL SENSITIVITY

0.067 IU/mL (NIBSC 00/588)

4. ANALYTICAL SPECIFICITY

No interference was observed with samples from patients with high-level of rheumatoid factor, and pregnant woman. Same day and frozen specimens have been tested to check for interferences due to collection and storage. Total of 100 samples reactive for anti-HBc, anti-HCV and anti-HIV-1 were screened for HBsAg with EIAGEN HBsAg Kit. 98 out of 100 samples were negative for HBsAg. 200 blood samples from patients were also tested with EIAGEN HBsAg Kit. 191 out of 200 samples had negative screening results for HBsAg. 8 out of 9 samples with initial reactive screening results had repeat reactive test results with EIAGEN HBsAg Kit but hepatitis B virus was not confirmed in all cases.

5. DETECTION OF MUTATIONS

Panel of 108 samples sequenced by PCR were tested to demonstrate the performance of EIAGEN HBsAg Kit in detection of HBsAg mutations. The results are given in the table below.

Background	Number	EIAGEN HBsAg Kit
adr (+) wild type	35	33
adr (+) 4 mutations	5	4
adw (+) wild type	37	34
adw (+) 16 mutations	25	24
ayw (+) wild type	2	2
ayw (+) 2 mutations	2	2
ayr (+) 2 mutations	2	2
Total	108	101

S. SUGGESTIONS FOR TROUBLESHOOTING

Adherence to assay procedure and specifications, as well as a correct use of reagents and proper pipetting, may help to avoid the following kinds of errors:

ERROR	POSSIBLE CAUSES / SUGGESTIONS
OD very different ($\pm 50\%$) from OD reported on QC	- incorrect dispensing volume of reagents (suggestion: check the correspondence between the volume dispensed by the pipette and the one required by the assay; re-calibrate again pipettes) - incorrect temperature or incorrect incubation time (suggestion: more care in the incubator maintenance; note down the beginning of the incubation) - error in washing or in photometer reading (suggestion: check operating or settings of respective instruments) - contamination of Substrate Solutions or Conjugate (suggestion: use only disposable and clean plastic containers)
Low reproducible results	- not constant dispensing volume of samples or reagents (suggestion: check the pipettes precision and the correspondence between the volume dispensed by the pipette and the one required by the assay; re-calibrate again pipettes) - error in washing or in reading (suggestion: check operating or settings of respective instruments) - contamination of Substrate Solutions (suggestion: use only disposable and clean plastic containers) - pollution or degradation of reagents (suggestion: use appropriate tips, disposable and clean plastic containers for reagents and high quality distilled or equivalent water)
no colorimetric reaction after addition of substrate solutions	- some reagent not pipetted - strong contamination of Conjugate or Substrate Solutions - errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)
too low reaction (too low ODs)	- incubation time too short, incubation temperature too low - incorrect conjugate dilution
too high reaction (too high ODs)	- incorrect conjugate dilution - incubation time too long, incubation temperature too high - water quality for wash buffer insufficient (low grade of deionization) - insufficient washing (conjugates not properly removed)
unexplainable outliers	- contamination of pipettes, tips or containers - inconstant and insufficient washing (conjugates not properly removed)
too high within-run CV%	- reagents and/or strips not pre-warmed to Room Temperature prior to use - plate washer is not washing correctly (suggestion: clean washer head)
too high between-run CV%	- incubation conditions not constant (time, temperature) - controls and samples not dispensed at the same time (with the same intervals) (check pipetting order) - person-related variation

T. AUTOMATION

The procedures identified in this Instruction for Use are for manual testing only. When using automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow their approved validation procedures to demonstrate compatibility of this product on automated systems.

U. LIMITATIONS

1. Positive results must be confirmed with another available method and interpreted in conjunction with the patient clinical information.
2. Antigens may be undetectable during the early stage of the disease. Therefore, negative results obtained with EIAgen HBsAg Kit are only indication that the sample does not contain detectable level of hepatitis B virus surface antigen and any negative result should not be considered as conclusive evidence that the individual is not infected with HBV or the blood unit is not infected with HBV.
3. If, after retesting of the initially reactive samples, the assay results are negative, these samples should be considered as non-repeatable (false positive) and interpreted as negative. As with many very sensitive ELISA assays, false positive results can occur due to the several reasons, most of which are related but not limited to inadequate washing step. For more information please refer to "Troubleshooting Guide", or contact Adaltis technical support for further assistance.
4. The most common assay mistakes are: using kits beyond the expiry date, bad washing procedures, contaminated reagents, incorrect assay procedure steps, insufficient aspiration during washing, failure to add specimens or reagents, improper operation with the laboratory equipment, timing errors, the use of highly hemolyzed specimens or specimens containing fibrin, incompletely clotted serum specimens.
5. The prevalence of the marker will affect the assay's predictive values.
6. This assay cannot be utilized to test pooled (mixed) plasma. EIAgen HBsAg Kit has been evaluated only with individual serum or plasma specimens.
7. EIAgen HBsAg Kit is a qualitative assay and the results cannot be used to measure antigen concentration.
8. INDICATIONS OF INSTABILITY DETERIORATION OF THE REAGENT:

Values of the Positive or Negative controls, which are out of the indicated quality control range, are indicators of possible deterioration of the reagents and/or operator or equipment errors. In such case, the results should be considered as invalid and the samples must be retested. In case of constant erroneous results and proven deterioration or instability of the reagents, immediately substitute the reagents with new one or contact Adaltis technical support for further assistance.

BIBLIOGRAPHY

1. Stevens, C. E., P. E. Taylor, and M. J. Tong. 1988. Viral hepatitis and liver disease. Alan R. Riss, New York, N.Y. 142. Stevens, C. E., P. E. Taylor, M. J. Tong, P. T. Toy, G. N. Vyas, P. V. Nair.
2. J. Y. Weissman, and S. Krugman. 1987. Yeast-recombinant hepatitis B vaccine. Efficacy with hepatitis B immune globulin in prevention of perinatal hepatitis B virus transmission. JAMA 257:2612–2616. 143. Stevens, C. E., P. T. Toy, P. E. Taylor, T. Lee, and H. Y. Yip. 1992. Prospects for control of hepatitis B virus infection: implications of childhood vaccination and long term protection. Pediatrics 90(Suppl.):170–173.
3. Hurie, M. B., E. E. Mast, and J. P. Davis. 1992. Horizontal transmission of hepatitis B virus infection to U.S. born children of Hmong refugees. Pediatrics 89:269–273.
4. Szmuness, W., C. E. Stevens, E. J. Harley, E. A. Zang, W. R. Olesko, D. C. Williams, R. Sadovsky, J. M. Morrison, and A. Kellner. 1980. Hepatitis B vaccine: demonstration of efficacy in a controlled trial in a high risk population in the U.S. N. Engl. J. Med. 303:833–841.
5. Bhatnagar, P. K., E. Papas, H. E. Blum, D. R. Milich, D. Nitecki, M. J. Karels, and G. N. Vyas. 1982. Immune response to synthetic peptide analogues of hepatitis B surface antigen specific for the a determinant. Proc. Natl. Acad. Sci. USA 79:4400–4404.

Claims regarding the quality of the kit should be addressed to:

Adaltis Srl
Via Durini, 27
20122 Milano (Italy)
Tel. +39-0774-5791 - Fax +39-0774-353085
www.adaltis.net
e-mail: info@adaltis.net

For further information and help, ask us at info@adaltis.net

ElAgen

HCV Ab (v.4) Kit

REF 071067


Σ 96

REF 071064

Σ 192

REF 071068

Σ 480

IVD

CE 0459

This package insert must be read carefully before product use.

Package insert instructions must be carefully followed.

Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package.

Manufacturer:

Adaltis S.r.l

Via Durini, 27

20122 Milano (Italy)

Tel. +39-0774-5791 - Fax +39-0774-353085

www.adaltis.net

en

SYMBOLS USED ON LABELS

English EN							
	In Vitro Diagnostic Medical Device	Catalogue Number	Lot Number	Attention, See Instructions For Use	Temperature Limitation	Use By	Number of Test
				MICROPLATE	CONTROL+	CONTROL-	CAL
	Manufacturer	Keep away from Sunlight	Date of Manufacture	Microplate	Positive Control	Negative Control	Calibrator
	CONJ	DILSPE	SUBSTMB	SOLNSTOP	WASHBUF20X	DILAS	RCNSx mL
	Conjugate	Sample Diluent	Substrate TMB	Stop Solution (H ₂ SO ₄ 0,3 M)	Wash Buffer Conc. (20x)	Assay Diluent	Reconstitute with x mL
	Biological Risk	Danger	Warning				

Attention:

Stop Solution classified as: Skin Corr. 1A

- Signal word:**
Danger
- Hazard-determining components of labelling:**
Sulphuric acid
- Hazard statements:**
H314 Causes severe skin burns and eye damage.
- Precautionary statements:**
P260 Do not breathe dust/fume/gas/mist/vapours/spray.
P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310 Immediately call a POISON CENTER/doctor.
P405 Store locked up.
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Attention:

Negative Control, Positive Control, Calibrator, Conjugate, Sample Diluent, Assay Diluent and Wash Buffer Concentrate 20X classified as: Skin Sens. 1

- Signal word:**
Warning
- Hazard-determining components of labelling:**
Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1)
- Hazard statements:**
H317 May cause an allergic skin reaction.
- Precautionary statements:**
P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P321 Specific treatment (see on this label).
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
P302+P352 IF ON SKIN: Wash with plenty of water.
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Refer to www.adaltis.net for the Safety Data Sheets.

A. INTENDED USE

Fourth generation Enzyme ImmunoAssay (ELISA) for the determination of antibodies to Hepatitis C Virus in human plasma (EDTA, Heparin and Citrate) and sera. The kit may be used for the screening of blood units of HCV-infected patients.

For "in vitro" diagnostic use only.

B. INTRODUCTION

The World Health Organization (WHO) define Hepatitis C infection as follows:

"Hepatitis C is a viral infection of the liver which had been referred to as parenterally transmitted "non A, non B hepatitis" until identification of the causative agent in 1989. The discovery and characterization of the hepatitis C virus (HCV) led to the understanding of its primary role in post transfusion hepatitis and its tendency to induce persistent infection.

HCV is a major cause of acute hepatitis and chronic liver disease, including cirrhosis and liver cancer. Globally, an estimated 170 million persons are chronically infected with HCV and 3 to 4 million persons are newly infected each year. HCV is spread primarily by direct contact with human blood. The major causes of HCV infection worldwide are use of unscreened blood transfusions, and re-use of needles and syringes that have not been adequately sterilized. No vaccine is currently available to prevent hepatitis C and treatment for chronic hepatitis C is too costly for most persons in developing countries to afford. Thus, from a global perspective, the greatest impact on hepatitis C disease burden will likely be achieved by focusing efforts on reducing the risk of HCV transmission from nosocomial exposures (e.g. blood transfusions, unsafe injection practices) and high-risk behaviours (e.g. injection drug use).

Hepatitis C virus (HCV) is one of the viruses (A, B, C, D, and E), which together account for the vast majority of cases of viral hepatitis. It is an enveloped RNA virus in the flaviviridae family which appears to have a narrow host range. Humans and chimpanzees are the only known species susceptible to infection, with both species developing similar disease. An important feature of the virus is the relative mutability of its genome, which in turn is probably related to the high propensity (80%) of inducing chronic infection. HCV is clustered into several distinct genotypes which may be important in determining the severity of the disease and the response to treatment.

The incubation period of HCV infection before the onset of clinical symptoms ranges from 15 to 150 days. In acute infections, the most common symptoms are fatigue and jaundice; however, the majority of cases (between 60% and 70%), even those that develop chronic infection, are asymptomatic. About 80% of newly infected patients progress to develop chronic infection. Cirrhosis develops in about 10% to 20% of persons with chronic infection, and liver cancer develops in 1% to 5% of persons with chronic infection over a period of 20 to 30 years. Most patients suffering from liver cancer who do

not have hepatitis B virus infection have evidence of HCV infection. The mechanisms by which HCV infection leads to liver cancer are still unclear. Hepatitis C also exacerbates the severity of underlying liver disease when it coexists with other hepatic conditions. In particular, liver disease progresses more rapidly among persons with alcoholic liver disease and HCV infection. HCV is spread primarily by direct contact with human blood. Transmission through blood transfusions that are not screened for HCV infection, through the reuse of inadequately sterilized needles, syringes or other medical equipment, or through needle-sharing among drug-users, is well documented. Sexual and perinatal transmission may also occur, although less frequently. Other modes of transmission such as social, cultural, and behavioural practices using percutaneous procedures (e.g. ear and body piercing, circumcision, tattooing) can occur if inadequately sterilized equipment is used. HCV is not spread by sneezing, hugging, coughing, food or water, sharing eating utensils, or casual contact.

In both developed and developing countries, high risk groups include injecting drug users, recipients of unscreened blood, haemophiliacs, dialysis patients and persons with multiple sex partners who engage in unprotected sex. In developed countries, it is estimated that 90% of persons with chronic HCV infection are current and former injecting drug users and those with a history of transfusion of unscreened blood or blood products. In many developing countries, where unscreened blood and blood products are still being used, the major means of transmission are unsterilized injection equipment and unscreened blood transfusions. In addition, people who use traditional scarification and circumcision practices are at risk if they use or re-use unsterilized tools.

WHO estimates that about 170 million people, 3% of the world's population, are infected with HCV and are at risk of developing liver cirrhosis and/or liver cancer. The prevalence of HCV infection in some countries in Africa, the Eastern Mediterranean, South-East Asia and the Western Pacific (when prevalence data are available) is high compared to some countries in North America and Europe.

Diagnostic tests for HCV are used to prevent infection through screening of donor blood and plasma, to establish the clinical diagnosis and to make better decisions regarding medical management of a patient. Diagnostic tests commercially available today are based on Enzyme immunoassay (EIA) for the detection of HCV specific antibodies. EIAs can detect more than 95% of chronically infected patients but can detect only 50% to 70% of acute infections. A recombinant immunoblot assay (RIBA) that identifies antibodies which react with individual HCV antigens is often used as a supplemental test for confirmation of a positive EIA result. Testing for HCV circulating by amplification tests RNA (e.g. polymerase chain reaction or PCR, branched DNA assay) is also being utilized for confirmation of serological results as well as for assessing the effectiveness of antiviral therapy. A positive result indicates the presence of active infection and a potential for spread of the infection and or/the development of chronic liver disease.

Antiviral drugs such as interferon taken alone or in combination with ribavirin, can be used for the treatment of persons with chronic hepatitis C, but the cost of treatment is very high. Treatment with interferon alone is effective in about 10% to 20% of patients. Interferon combined with ribavirin is effective in about 30% to 50% of patients. Ribavirin does not appear to be effective when used alone.

There is no vaccine against HCV. Research is in progress but the high mutability of the HCV genome complicates vaccine development. Lack of knowledge of any protective immune response following HCV infection also impedes vaccine research. It is not known whether the immune system is able to eliminate the virus.

Some studies, however, have shown the presence of virus neutralizing antibodies in patients with HCV infection. In the absence of a vaccine, all precautions to prevent infection must be taken including (a) screening and testing of blood and organ donors; (b) Virus inactivation of plasma derived products; (c) implementation and maintenance of infection control practices in health care settings, including appropriate sterilization of medical and dental equipment; (d) promotion of behaviour change among the general public and health care workers to reduce overuse of injections and to use safe injection practices; and (e) Risk reduction counselling for persons with high-risk drug and sexual practices. “

The genome encodes for structural components, a nucleocapsid protein and two envelope glycoproteins, and functional constituents involved in the virus replication and protein processing. The nucleocapsid-encoding region seems to be the most conservative among the isolates obtained all over the world.

C. PRINCIPLE OF THE TEST

Microplates are coated with HCV-specific antigens derived from “core” and “ns” regions encoding for conservative and immunodominant antigenic determinants (Core peptide, recombinant NS3, NS4 and NS5 peptides).

The solid phase is first treated with the diluted sample and HCV Ab are captured, if present, by the antigens. After washing out all the other components of the sample, in the 2nd incubation bound HCV antibodies, IgG and IgM as well, are detected by the addition of polyclonal specific anti hlgG&M antibodies, labelled with peroxidase (HRP).

The enzyme captured on the solid phase, acting on the substrate TMB mixture, generates an optical signal that is proportional to the amount of anti HCV antibodies present in the sample. A cut-off value let optical densities be interpreted into HCV antibody negative and positive results.

D. COMPONENTS

The kit contains reagents for 96 tests (code 071067), 192 tests (code 071064), or 480 tests (code 071068).

Microplate	1
Negative Control	1x4 mL/vial
Positive Control	1x2 mL/vial
Calibrator	2 vials
Wash Buffer Concentrate 20x	1x50 mL/vial
Conjugate	1x16 mL/vial
Sample Diluent	1x50 mL/vial
Substrate TMB	1x16 mL/vial
Stop Solution	1x15 mL/vial
Assay Diluent	1x8 mL/vial
Plate sealing foils	2
Number of tests	96
Code	071067

Microplate	2
Negative Control	2x4 mL/vial
Positive Control	1x4 mL/vial
Calibrator	3 vials
Wash Buffer Concentrate 20x	2x50 mL/vials
Conjugate	2x16 mL/vials
Sample Diluent	2x50 mL/vials
Substrate TMB	2x16 mL/vials
Stop Solution	2x15 mL/vial
Assay Diluent	2x8 mL/vial
Plate sealing foils	4
Number of tests	192
Code	071064

Microplate	5
Negative Control	1x20 mL/vial
Positive Control	1x10 mL/vial
Calibrator	7 vials
Wash Buffer Concentrate 20x	5x50 mL/vials
Conjugate	2x40 mL/vial
Sample Diluent	5x50 mL/vials
Substrate TMB	2x40 mL/vials
Stop Solution	2x40 mL/vial
Assay Diluent	1x40 mL/vial
Plate sealing foils	10
Number of tests	480
Code	071068

1. Microplate

12 strips of 8 microwells coated with Core peptide, recombinant NS3, NS4 and NS5 peptides. Plates are sealed into a aluminium pouch with desiccant. Bring the microplate to room temperature (18...24°C) before opening the bag. Unused strips have to be returned into the pouch and the pouch has to be sealed and stored back to 2...8°C, in presence of the desiccant.

2. Negative Control

Ready to use control. It contains 10 mM Na-citrate buffer pH 6.0 ± 0.1, 2% casein as protein base and 0.1% Proclin 150 as preservative. The negative control is olive green colour coded.

3. Positive Control

Ready to use control. It contains 1% goat serum proteins, human antibodies positive to HCV, 10 mM Na-citrate buffer pH 6.0 ± 0.1, 0.5% Tween 20, 0.09% Na-azide and 0.1% Proclin 150 as preservatives. The Positive Control is dark green colour coded.

Important Note: The absence of viable pathogens in the Positive Control can not be fully ensured, and therefore,

the reagent should be handled as potentially biohazardous, in accordance with good laboratory practices.

4. Calibrator

Lyophilized calibrator. To be dissolved with the volume of EIA grade water reported on the label. It contains foetal bovine serum proteins, human antibodies to HCV whose content is calibrated on the NIBSC Working Standard code 06/188-006, 10 mM Na-citrate buffer pH 6.0 ± 0.1 , 0.3 mg/mL gentamicin sulfate and 0.1% Proclin 150 as preservatives.

Important Note: The absence of viable pathogens in the Calibrator can not be fully ensured, and therefore, the reagent should be handled as potentially biohazardous, in accordance with good laboratory practices.

Note: The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label.

5. Wash Buffer Concentrate 20x

20x concentrated solution. Once diluted, the wash solution (wash buffer diluted) contains 10 mM phosphate buffer pH 7.0 ± 0.2 , 0.05% Tween 20 and 0.05% Proclin 150.

Once diluted, the wash solution is stable for 1 week at $2...8^{\circ}\text{C}$.

6. Conjugate

Ready to use and red colour coded reagent. It contains Horseradish Peroxidase conjugated goat polyclonal antibodies to human IgG and IgM, 5% BSA, 10 mM Citrate buffer pH 6.4 ± 0.1 , 0.1% Proclin and 0.05% Tween 20 as preservatives.

7. Substrate TMB

Ready-to-use component. It contains 50 mM citrate-phosphate buffer pH 3.5-3.8, 4% dimethylsulphoxide, 0.03% tetra-methyl-benzidine or TMB and 0.02% hydrogen peroxide or H_2O_2 . Mix gently before use.

Note: To be stored protected from light as sensitive to strong illumination.

8. Assay diluent

Ready-to-use component. It contains goat serum, 10 mM tris buffered solution pH 8.0 ± 0.1 containing 0.1% Proclin 150 and 0.09% Na azide for the pre-treatment of samples and controls in the plate, blocking interference.

9. Stop Solution

Ready-to-use component.

It contains 0.3 M H_2SO_4 solution. Mix gently before use.

10. Sample Diluent

Ready-to-use component and dark green colour coded. It contains 1% casein, 10 mM Na-citrate buffer pH 6.0 ± 0.1 and 0.1% Proclin 150 as preservative.

To be used to dilute the sample.

Note: The diluent changes colour from olive green to dark bluish green in the presence of sample.

E. MATERIALS REQUIRED BUT NOT PROVIDED

1. Calibrated Micropipettes (200 μL and 10 μL) and disposable plastic tips.
2. EIA grade water (bidistilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator capable to provide a temperature of $+37^{\circ}\text{C}$.
6. Calibrated ELISA microwell reader with 450nm (reading) and possibly with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory. This package insert must be read carefully before product use.
2. Read carefully the Safety Data Sheet (SDS) before product use.
3. When the kit is used for the screening of blood units and blood components, it has to be used in a laboratory certified and qualified by the national authority in that field (Ministry of Health or similar entity) to carry out this type of analysis.
4. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
5. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
6. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-born microbial agents, when opening kit vials and microplates and when performing the test. Protect the Substrate (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.
7. Upon receipt, store the kit at $2...8^{\circ}\text{C}$ into a temperature controlled refrigerator or cold room.
8. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
9. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.
10. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
11. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.

12. Do not use the kit after the expiration date stated on the external container and internal (vials) labels.
13. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
14. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.
15. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated before waste. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min.
16. Accidental spills from samples and operations have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
17. The Stop Solution contains 0.3M sulphuric acid. Avoid contact with skin and eyes. In the event of contact, rinse immediately with plenty of water.
18. Dispose of reagent solutions containing sodium azide as preservatives according to all local, state and national regulations. To dispose of reagents containing azide, flush away using copious amounts of water. Dispose with caution as sodium azide can form explosive compounds on prolonged contact with lead or copper piping.
19. Do not smoke, eat, drink or apply cosmetics in areas in which specimens or kit reagents are handled
20. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.
21. Do not pipette by mouth.

G. SPECIMEN: PREPARATION AND RECCOMANDATIONS

1. Blood is drawn aseptically by venepuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.
2. Avoid any addition of preservatives to samples; especially sodium azide as this chemical would affect the enzymatic activity of the conjugate, generating false negative results.
3. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. When the kit is used for the screening of blood units, bar code labeling and electronic reading is strongly recommended.

4. Haemolysed (red) and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.
5. Sera and plasma can be stored at +2...8°C for up to five days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.
6. If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8µ filters to clean up the sample for testing.

H. PREPARATION OF COMPONENTS AND WARNINGS

A study conducted on an opened kit has not pointed out any relevant loss of activity up to 1 re-use of the device and up to 6 months.

1. Microplates:

Allow the microplate to reach room temperature (about 1 hr) before opening the pouch. Check that the desiccant has not turned to dark green, indicating a defect in manufacturing. In this case, call Adaltis's customer service.

Unused strips have to be placed back inside the aluminium pouch, with the desiccant supplied, firmly zipped and stored at +2...8°C. After first opening, remaining strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

2. Negative Control:

Ready to use. Mix well on vortex before use.

3. Positive Control:

Ready to use. Mix well on vortex before use. Handle this component as potentially infectious.

4. Calibrator:

Dissolve carefully the content of the lyophilised vial with the volume of EIA grade water reported on its label. Mix well on vortex before use.

Handle this component as potentially infectious.

Note: When dissolved the Calibrator is not stable. Store in aliquots at -20°C.

5. Wash Buffer Concentrate 20x (vial of 50 mL):

The whole content of the 20x concentrated solution has to be diluted with EIA grade water up to 1000 mL (the volume is reported on the label) and mixed gently end-over-end before use. As some salt crystals may be present into the vial, take care to dissolve all the content when preparing the solution. During preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles.

Note: Once diluted, the wash solution is stable for 1 week at +2...8°C.

6. Conjugate:

Ready to use. Mix well on vortex before use.

Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes. If this component

has to be transferred use only plastic, possibly sterile disposable containers.

7. Substrate TMB:

Ready to use. Mix well on vortex before use.

Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes. Do not expose to strong illumination, oxidizing agents and metallic surfaces. If this component has to be transferred use only plastic, and if possible, sterile disposable container.

8. Assay Diluent:

Ready to use. Mix well on vortex before use.

9. Stop Solution:

Ready to use. Mix well on vortex before use.

10. Sample Diluent:

Ready to use. Mix well on vortex before use.

I. INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT

1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample or the components of the kit. They should also be regularly maintained in order to show a precision of 1% and a trueness of $\pm 2\%$. Decontamination of spills or residues of kit components should also be carried out regularly.
2. The ELISA incubator has to be set at $+37^{\circ}\text{C}$ (tolerance of $\pm 0.5^{\circ}\text{C}$) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. The ELISA washer is extremely important to the overall performances of the assay. The washer must be carefully validated and correctly optimized. Usually 4-5 washing cycles (aspiration + dispensation of $350\ \mu\text{L}/\text{well}$ of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to set correctly their number, it is recommended to run an assay with the kit controls/calibrator and well-characterized negative and positive reference samples, and check to match the values reported below in sections O "Internal Quality Control". Regular calibration of the volumes delivered and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.
4. Incubation times have a tolerance of $\pm 5\%$.
 - ✓ Short Incubation Method (for 1st / 2nd incubation tolerance between 43 min to 47 min; for 3rd incubation tolerance between 14 and 16 min).
 - ✓ Standard Incubation Method (for 1st incubation tolerance between 57 min to 63 min; for 2nd and 3rd incubation tolerance between 29 and 31 min).
5. The ELISA microplate reader has to be equipped with a reading filter of 450nm and ideally with a second

filter (620-630nm) for blanking purposes. Its standard performances should be (a) bandwidth $\leq 10\ \text{nm}$; (b) absorbance range from 0 to ≥ 2.0 ; (c) linearity to ≥ 2.0 ; repeatability $\geq 1\%$. Blanking is carried out on the well identified in section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer's instructions.

6. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, shaking, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in sections O "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing samples and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells due to strongly reactive samples, leading to false positive results. The use of ELISA automated work stations is recommended for blood screening and when the number of samples to be tested exceed 20-30 units per run.

L. PRE ASSAY CONTROLS AND OPERATIONS

1. Check the expiration date of the kit printed on the external label of the kit box. Do not use if expired.
2. Check that the liquid components are not contaminated by naked-eye visible particles or aggregates. Check that the Substrate TMB is colorless or pale blue by aspirating a small volume of it with a sterile transparent plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box. Check that the aluminum pouch, containing the microplate, is not punctured or damaged.
3. Dilute all the content of the 20x concentrated Wash Buffer as described above.
4. Dissolve the Calibrator as described above.
5. Allow all the other components to reach room temperature (about 1 hr) and then mix as described.
6. Set the ELISA incubator at $+37^{\circ}\text{C}$ and prepare the ELISA washer by priming with the diluted wash solution, according to the manufacturers instructions. Set the right number of washing cycles as found in section I.3.
7. Check that the ELISA reader has been turned on at least 20 minutes before reading.
8. If using an automated workstation, turn it on, check settings and be sure to use the right assay protocol.
9. Check that the micropipettes are set to the required volume.
10. Check that all the other equipment is available and ready to use.
11. In case of problems, do not proceed further with the test and advise the supervisor.

M. ASSAY PROCEDURE

The assay has to be carried out according to what reported below, taking care to maintain the same incubation time for all the samples in testing.

The Assay Procedure can be done with two times incubation procedures. Choose the one that is required by your regulation:

1. Standard Incubation (1st incubation 60 minutes, 2nd and 3rd incubation 30 minutes)
2. Short Incubation (1st and 2nd incubation 45 minutes, 3rd incubation 15 minutes)

1. Standard Incubation - Manual assay:

1. Place the required number of microwells in the microwell holder. Leave the 1st well empty for the operation of blanking.
2. Dispense 200 µL of Negative Control in triplicate, 200 µL Calibrator in duplicate and 200 µL Positive Control in single in proper wells. Do not dilute Controls and Calibrator as they are pre-diluted, ready to use!
3. Add 200 µL of Sample Diluent to all the sample wells; then dispense 10 µL sample in each properly identified well. Mix gently the plate, avoiding overflowing and contaminating adjacent wells, in order to fully disperse the sample into its diluent.

Important note: Check that the colour of the Sample Diluent, upon addition of the sample, changes from light green to dark bluish green, monitoring that the sample has been really added.

4. Dispense 50 µL Assay Diluent into all the controls/calibrator and sample wells. Check that the color of samples has turned to dark blue.
5. Incubate the microplate for **60 min at +37°C**.
Important note: Strips have to be sealed with the adhesive sealing foil, supplied, only when the test is carried out manually. Do not cover strips when using ELISA automatic instruments.
6. Wash the microplate with an automatic washer by delivering and aspirating 350 µL/well of diluted wash solution as reported in section I.3.
7. Pipette 100 µL Enzyme Conjugate into each well, except the 1st blanking well, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1.

Important note: Be careful not to touch the plastic inner surface of the well with the tip filled with the Enzyme Conjugate. Contamination might occur.

8. Incubate the microplate for **30 min at +37°C**.
9. Wash microwells as reported in section I.3.
10. Pipette 100 µL Substrate TMB mixture into each well, the blank well included. Then incubate the microplate at room temperature (18...24°C) for **30 minutes**.

Important note: Do not expose to strong direct illumination. High background might be generated.

11. Pipette 100 µL Stop Solution into all the wells using the same pipetting sequence as in step 10 to stop the enzymatic reaction. Addition of stop solution will

turn the positive control and positive samples from blue to yellow.

12. Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and possibly at 620-630nm (background subtraction), blanking the instrument on A1.

Important notes:

1. If the second filter is not available, ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
2. Reading has should ideally be performed immediately after the addition of the Stop Solution but definitely no longer than 20 minutes afterwards. Some self oxidation of the substrate can occur leading to a higher background.
3. Shaking at 350 ±150 rpm during incubation has been proved to increase the sensitivity of the assay of about 20%.

2. Short Incubation - Manual assay:

1. Place the required number of Microwells in the microwell holder. Leave the 1st well empty for the operation of blanking.
2. Dispense 200 µL of Negative Control in triplicate, 200 µL Calibrator in duplicate and 200 µL Positive Control in single in proper wells. Do not dilute Controls and Calibrator as they are pre-diluted, ready to use!
3. Add 200 µL of Sample Diluent to all the sample wells; then dispense 10 µL sample in each properly identified well. Mix gently the plate, avoiding overflowing and contaminating adjacent wells, in order to fully disperse the sample into its diluent.

Important note: Check that the colour of the Sample Diluent, upon addition of the sample, changes from light green to dark bluish green, monitoring that the sample has been really added.

4. Dispense 50 µL Assay Diluent into all the controls/calibrator and sample wells. Check that the color of samples has turned to dark blue.
5. Incubate the microplate for **45 min at +37°C**.
Important note: Strips have to be sealed with the adhesive sealing foil, supplied, only when the test is carried out manually. Do not cover strips when using ELISA automatic instruments.
6. Wash the microplate with an automatic washer by delivering and aspirating 350 µL/well of diluted washing solution as reported in section I.3.
7. Pipette 100 µL Enzyme Conjugate into each well, except the 1st blanking well, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1.

Important note: Be careful not to touch the plastic inner surface of the well with the tip filled with the Enzyme Conjugate. Contamination might occur.

8. Incubate the microplate for **45 min at +37°C**.
9. Wash microwells as reported in section I.3.
10. Pipette 100 µL Substrate TMB mixture into each well, the blank well included. Then incubate the

microplate at room temperature (18...24°C) for **15 minutes**.

Important note: Do not expose to strong direct illumination. High background might be generated.

- Pipette 100 µL Stop Solution into all the wells using the same pipetting sequence as in step 10 to stop the enzymatic reaction. Addition of stop solution will turn the positive control and positive samples from blue to yellow.
- Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and possibly at 620-630nm (background subtraction), blanking the instrument on A1.

Important notes:

- If the second filter is not available, ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
- Reading has should ideally be performed immediately after the addition of the Stop Solution but definitely no longer than 20 minutes afterwards. Some self oxidation of the substrate can occur leading to a higher background.
- Shaking at 350 ±150 rpm during incubation has been proved to increase the sensitivity of the assay of about 20%.

N. ASSAY SCHEME

Method	Operations (Standard Incubation)	Operations (Short Incubation)
Controls & Calibrator	200 µL	200 µL
Sample Diluent & Sample	200 µL diluent.+ 10 µL sample	200 µL diluent.+ 10 µL sample
Assay Diluent	50 µL	50 µL
1st incubation	60 min (± 3)	45 min (± 2)
Temperature	+37°C	+37°C
Wash step	4-5 cycles	4-5 cycles
Enzyme Conjugate	100 µL	100 µL
2nd incubation	30 min (± 1)	45 min (± 2)
Temperature	+37°C	+37°C
Wash step	4-5 cycles	4-5 cycles
Substrate TMB	100 µL	100 µL
3rd incubation	30 min (± 1)	15 min (± 1)
Temperature	Room temperature (18...24°C)	Room temperature (18...24°C)
Stop Solution	100 µL	100 µL
Reading OD	450/620nm	450/620nm

An example of dispensation scheme is reported below (valid for both incubation time procedures):

Microplate

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	S2										
B	NC	S3										
C	NC	S4										
D	NC	S5										
E	CAL	S6										
F	CAL	S7										
G	PC	S8										
H	S1	S9										

Legenda: BLK = Blank NC = Negative Control
CAL = Calibrator PC = Positive Control S = Sample

O. INTERNAL QUALITY CONTROL

A check is carried out on the controls and the calibrator any time the kit is used in order to verify whether their OD450/620nm values are as expected and reported in the table below.

Check	Requirements
Blank well	< 0.100 OD 450/620nm value
Negative Control (NC)	< 0.050 mean OD450/620nm after blanking
Calibrator	S/Co >1.1
Positive Control	>1.000 OD450/620nm value

If the results of the test match the requirements stated above, proceed to the next section.

If they do not, do not proceed any further and perform the following checks:

Problems	Check
Blank well > 0.100 OD450nm	1. that the Substrate solution has not become contaminated during the assay
Negative Control (NC) > 0.050 OD450nm after blanking	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of the negative one); 4. that no contamination of the negative control or of the wells where the control was dispensed has occurred due to spills of positive samples or of the enzyme conjugate; 5. that micropipettes have not become contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
Calibrator S/Co<1.1	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution (ex.: dispensation of negative control instead of calibrator) 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
Positive Control <1.000 OD450nm	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the control (dispensation of negative control instead of positive control. In this case, the negative control will have an OD450nm value > 0.150,too. 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

If any of the above problems have occurred, report the problem to the supervisor for further actions.

P. RESULTS

The tests results are calculated by means of a cut-off value determined with the following formula:

$$\text{Cut-Off (Co)} = \text{NC mean} + 0.350$$

The value found for the test is used for the interpretation of results as described in the next paragraph.

Q. INTERPRETATION OF RESULTS

Test results are interpreted as ratio of the sample OD450nm and the Cut-Off value (or S/Co) according to the following table:

S/Co	Interpretation
< 0.9	Negative
0.9 - 1.1	Equivocal
> 1.1	Positive

A negative result indicates that the patient has not been infected by HCV or that the blood unit may be transfused. Any patient showing an equivocal result should be tested again on a second sample taken 1-2 weeks later from the patient and examined. The blood unit should not be transfused.

A positive result is indicative of HCV infection and therefore the patient should be treated accordingly or the blood unit should be discarded.

Important notes:

1. Interpretation of results should be done under the supervision of the responsible of the laboratory to reduce the risk of judgement errors and misinterpretations.
2. Any positive result should be confirmed by an alternative method capable to detect IgG and IgM antibodies (confirmation test) before a diagnosis of viral hepatitis is formulated.
3. As proved in the Performance Evaluation of the product, the assay is able to detect seroconversion to anti HCV core antibodies **earlier** than some other commercial kits. Therefore a positive result, not confirmed with these commercial kits, does not have to be ruled out as a false positive result ! The sample has to be anyway submitted to a confirmation test.
4. As long as the assay is able to detect also IgM antibodies some discrepant results with other commercial products for the detection of anti HCV antibodies - lacking anti hIgM conjugate in the formulation of the enzyme tracer and therefore missing IgM reactivity - may be present. The real positivity of the sample for antibodies to HCV should be then confirmed by examining also IgM reactivity, important for the diagnosis of HCV infection.
5. When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.
6. Diagnosis of viral hepatitis infection has to be done and released to the patient only by a qualified medical doctor.

An example of calculation is reported below:

The following data must not be used instead of real figures obtained by the user.

Negative Control: 0.019 – 0.020 – 0.021 OD450nm

Mean Value: 0.020 OD450nm Lower than 0.050 – Accepted

Positive Control: 2.189 OD450nm

Higher than 1.000 – Accepted

Cut-Off = 0.020+0.350 = 0.370

Calibrator: 0.550 - 0.530 OD450nm

Mean value: 0.540 OD450nm S/Co = 1.4

S/Co higher than 1.1 – Accepted

Sample 1: 0.070 OD450nm

Sample 2: 1.690 OD450nm

Sample 1 S/Co < 0.9 = negative

Sample 2 S/Co > 1.1 = positive

R. PERFORMANCES

Evaluation of Performances has been conducted in accordance to what reported in the Common Technical Specifications or CTS (art. 5, Chapter 3 of IVD Directive 98/79/EC) and with both incubation procedures (standard and short).

1. LIMIT OF DETECTION

The limit of detection of the assay has been calculated using the short incubation assay by means of the British Working Standard for anti-HCV, NIBSC code 06/188-006. The table below reports the mean OD450nm values of this standard when diluted in negative plasma and then examined.

Dilution	Lot#1	Lot#2	Lot#3
Factor	S/Co	S/Co	S/Co
1 X	3,50	4,00	4,30
2 X	2,10	2,60	2,60
4 X	1,3	1,40	1,30
Negative Plasma	0,25	0,20	0,20

In addition the sample coded Accurun 1 - series 3000 - supplied by Boston Biomedica Inc., USA, has been evaluated "in toto" showing the results below:

Accurun 1 series	Lot#1	Lot#2	Lot#3
Factor	S/Co	S/Co	S/Co
1 X	2,90	3,04	3,40

2. DIAGNOSTIC SPECIFICITY AND SENSITIVITY

The Performance Evaluation of the device was carried out in a trial conducted on more than total 5000 samples.

2.1 Diagnostic specificity

It is defined as the probability of the assay of scoring negative in the absence of specific analyte. A total of more 5000 unselected donors, including 1st time donors, were examined.

The diagnostic specificity was assessed against a kit US FDA approved.

5043 blood donors were tested providing a specificity of 99.5%. 210 hospitalized patients were tested for HCV Ab; a diagnostic specificity of 99.5% was found.

Moreover, diagnostic specificity was assessed by testing 162 potentially interfering specimens (other infectious diseases, E.coli antibody positive, patients affected by non viral hepatic diseases, dialysis patients, pregnant women, hemolized, lipemic, etc.). A value of specificity of 100% was assessed.

No false reactivity due to the method of specimen preparation has been observed. Both plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and sera have been used to determine the value of specificity. Frozen specimens have been tested, as well, to check for interferences due to collection and storage. No interference was observed.

2.2 Diagnostic Sensitivity

It defined as the probability of the assay of scoring positive in the presence of specific analyte. The diagnostic sensitivity has been assessed externally on a total number of 348 specimens; a diagnostic sensitivity of 100% was found. Internally more than other 50 positive samples were tested, providing a value of diagnostic sensitivity of again 100%. Positive samples from infections carried out by different genotypes of HCV were tested as well.

Furthermore, most of seroconversion panels available from Boston Biomedica Inc., USA, (PHV) and Zeptometrix, USA, (HCV) have been studied.

Results are reported below for some of them.

Panel	N° samples	Adaltis ¹	Ortho ^{1,2}
PHV 901	11	9	9
PHV 904	7	2	4
PHV 905	9	3	4
PHV 906	7	7	7
PHV 907	7	3	2
PHV 908	13	10	8
PHV 909	3	2	2
PHV 910	5	3	3
PHV 911	5	3	3
PHV 912	3	1	1
PHV 913	4	2	2
PHV 914	9	5	5
PHV 915	4	3	0
PHV 916	8	4	3
PHV 917	10	6	6
PHV 918	8	2	0
PHV 919	7	3	3
PHV 920	10	6	6
HCV 10039	5	2	0
HCV 6212	9	6	7
HCV 10165	9	5	4

Note:

1. Positive samples detected
2. HCV v.3.0

Finally the Product has been tested on the panel EFS Ac HCV, lot n° 06.140817, supplied by the Etablissement Francais Du Sang (EFS), France, with the following results:

EFS Panel Ac HCV

sample	Lot#1 S/Co	Lot#2 S/Co	Lot#3 S/Co	Average value
HCV 1	0,53	0,52	0,55	Negative
HCV 2	3,28	5,91	3,04	Positive
HCV 3	2,17	3,18	2,56	Positive
HCV 4	2,26	2,23	2,35	Positive
HCV 5	6,10	7,06	6,90	Positive
HCV 6	1,66	1,77	1,67	Positive

3. PRECISION

It has been calculated on five samples, one negative and four positives examined in 4 replicates each in six separate runs. The precision, obtained for both protocols, is equivalent.

Results are reported as follows:

Intra-lot results: ElAgen HCV Ab (v.4) Kit - 1st Lot (short incubation assay)

		Precision - %CV		
Sample	S/Co Mean	Within Run	Between Run	Total
Negative	0.03	6.66	10.56	12.48
Positives	1.20	8.52	8.49	12.03
	1.51	7.69	12.22	14.44
	3.57	7.43	11.82	13.97
	11.87	3.42	9.32	9.92

Intra-lot results: ElAgen HCV Ab (v.4) Kit - 1st Lot (long incubation assay)

		Precision - %CV		
Sample	S/Co Mean	Within Run	Between Run	Total
Negative	0.04	4.67	12.34	13.19
Positives	1.47	9.62	11.40	14.92
	1.82	8.92	12.77	15.58
	4.31	4.59	12.88	13.67
	13.78	2.42	8.96	9.26

Inter-lot results: ElAgen HCV Ab (v.4) Kit - 1st, 2nd and 3rd Lot (short incubation assay)

		Precision - %CV		
Sample	Lot 1	Lot 2	Lot 3	
Negative	8,65	8,29	6,13	
Calibrator	4,98	4,44	5,38	
Positive	4,11	3,11	1,37	

The variability shown in the tables above did not result in sample misclassification.

S. SUGGESTIONS FOR TROUBLESHOOTING

Adherence to assay procedure and specifications, as well as a correct use of reagents and proper pipetting, may help to avoid the following kinds of errors

ERROR	POSSIBLE CAUSES / SUGGESTIONS
OD very different ($\pm 50\%$) from OD reported on QC	- incorrect dispensing volume of reagents (suggestion: check the correspondence between the volume dispensed by the pipette and the one required by the assay; re-calibrate again pipettes) - incorrect temperature or incorrect incubation time (suggestion: more care in the incubator maintenance; note down the beginning of the incubation) - error in washing or in photometer reading (suggestion: check operating or settings of respective instruments) - contamination of Substrate or Conjugate (suggestion: use only disposable and clean plastic containers)
Low reproducible results	- not constant dispensing volume of samples or reagents (suggestion: check the pipettes precision and the correspondence between the volume dispensed by the pipette and the one required by the assay; re-calibrate again pipettes) - error in washing or in reading (suggestion: check operating or settings of respective instruments) - contamination of Substrate (suggestion: use only disposable and clean plastic containers) - pollution or degradation of reagents (suggestion: use appropriate tips, disposable and clean plastic containers for reagents and high quality distilled or equivalent water)
no colorimetric reaction after addition of substrate	- some reagent not pipetted - strong contamination of Conjugate or Substrate - errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)
too low reaction (too low ODs)	- incubation time too short, incubation temperature too low - incorrect conjugate dilution
too high reaction (too high ODs)	- incorrect conjugate dilution - incubation time too long, incubation temperature too high - water quality for wash buffer insufficient (low grade of deionization) - insufficient washing (conjugates not properly removed)
unexplainable outliers	- contamination of pipettes, tips or containers - inconstant and insufficient washing (conjugates not properly removed)
too high within-run CV%	- reagents and/or strips not pre-warmed to Room Temperature prior to use - plate washer is not washing correctly (suggestion: clean washer head)
too high between-run CV%	- incubation conditions not constant (time, temperature) - controls and samples not dispensed at the same time (with the same intervals) (check pipetting order) - person-related variation

T. AUTOMATION

The procedures identified in this Instruction for Use are for manual testing only. When using automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow their approved validation procedures to demonstrate compatibility of this product on automated systems.

U. LIMITATIONS

Repeatable false positive results, not confirmed by RIBA or similar confirmation techniques, were assessed as less than 0.1% of the normal population.

Frozen samples containing fibrin particles or aggregates after thawing have been observed to generate some false results.

BIBLIOGRAPHY

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhoo MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human Tlymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon-α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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
1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhoo MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human T-lymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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BIBLIOGRAFIA

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhoo MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human T-lymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JL, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB. Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore—an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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
REFERENCES

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhoo MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human T-lymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JL, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998 ;44 :98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993 ;167 :66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993 ;269 :392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994 ;169 :90-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995 ;171 :768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993 ;153 :1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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
1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhoo MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human T-lymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. A cikk megtalálható: <http://www.hepnet.com/nih/gretch.html>. A Nemzeti Egészségügyi Hivatal által szervezett, C hepatitisz konferenciára készült, amit 1997 március 24-26.-án tartottak Bethesda Maryland-ban
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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
BIBLIOGRAFIA

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhoo MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human Tlymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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LITERATURA

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, wyd. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turyn, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, wyd. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Waszyngton, DC: US Government Printing Office, 1994; Nr publikacji NIH 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN i in. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH i in. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ i in. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW i in. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB i in. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ i in. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A i in. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK i in. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhoo MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human Tlymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK i in. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL i in. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW i in. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ i in. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland i A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre i P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. Artykuł dostępny pod adresem: <http://www.hepnet.com/nih/gretch.html>. Napisany w ramach konferencji Narodowego Instytutu Zdrowia dot wirusowego zapalenia wątroby typu C w dniach 24-26 marca 1997 roku w Bethesda w stanie Maryland (USA)
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli i E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams i J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada i T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata i T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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BIBLIOGRAFIE


1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhoo MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human Tlymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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
1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhoo MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human T-lymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon-α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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BIBLIOGRAFIJA

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhoo MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human T-lymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JL, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB. Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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BIBLIYOGRAPHI

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. MMWR 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. Hepatology 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. Viral Hepatitis and Liver Disease, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. Digestive diseases in the United States: epidemiology and impact. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. JAMA 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. Ann Intern Med 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. Lancet 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. Am J Med Sci 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. N Engl J Med 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. Science 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. Science 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. N Engl J Med 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. N Engl J Med 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. N Engl J Med 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. Semin Liver Dis 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. Am J Epidemiol 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. West J Med 1992;156:30-5.
18. Fingerhood MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. Arch Intern Med 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human T-lymphotropic viruses. Am J Pub Health 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. Blood 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. Blood 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. J Med Virology 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. ASAIO Journal 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. J Infect Dis 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. JAMA 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. J Infect Dis 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. J Infect 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. J Infect Dis 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. Arch Intern Med 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. Infect Control Hosp Epidemiol 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. J. Clin. Microbiol. 40:1656-1659.
32. Dusheiko, G., H. Schmilovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. Hepatology 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. J. Clin. Microbiol. 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. J. Clin. Microbiol. 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. J. Hepatol. 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. Hepatology 16:293-299.