CMV IgM

**For the Qualitative Detection of
IgM Antibodies to Cytomegalovirus (CMV)
in Human Serum and Plasma**

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 CMV IgM

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the qualitative detection of IgM antibodies to cytomegalovirus (CMV) in human serum or plasma (EDTA or heparinized), as an aid in the determination of an acute CMV infection.

This product has been licensed in Canada for diagnostic purposes and has not been licensed for blood donor screening.

Catalog Numbers: **L2KCM2** (200 tests)

Test Code: **CMM** Color: **Dark Green**

Summary and Explanation

Cytomegalovirus (CMV), a member of the herpes virus family, is found throughout the world. Humans of all ages are susceptible and infection is spread through sexual contact, direct exposure to infected body fluids, blood transfusions and organ transplants.^{1,2,5,6} The majority of infections are asymptomatic; however, CMV infections can be severe in neonates and immunocompromised individuals.^{6,7} Infection can also be severe in patients with congenital or acquired cellular immune defects, including cancer patients, organ recipients and AIDS patients.^{5,6,8}

CMV is the most common congenital infection, infecting between 0.5 and 2.5 percent of newborn infants. Five percent of these will develop classic cytomegalic inclusion disease with jaundice, pneumonia and central nervous system disorder. Infected infants may be asymptomatic at birth, but can develop neurological problems later in life.^{3,4}

Between 40 and 100 percent of people have detectable antibody,⁹ with the prevalence highest in developing countries. The presence of IgM antibody is a reliable indicator of an active infection.

Principle of the Procedure

IMMULITE 2000 CMV IgM is a solid-phase, enzyme-labeled chemiluminescent three-step immunoassay. The solid phase (bead) is coated with inactivated, purified CMV antigen (strain AD-169 from infected cell lysates). The liquid phase consists of three reagents: 1) polyclonal goat anti-human IgG antibody in buffer, 2) polyclonal goat anti-human IgG antibody in buffer, and 3) alkaline phosphatase (bovine calf intestine) conjugated to polyclonal goat anti-human IgM antibody in buffer.

In the first cycle, the patient sample and polyclonal goat anti-human IgG antibody are incubated together without bead for 30 minutes. During this time, anti-IgG antibodies block IgG from the patient sample.

In the second cycle, the pretreated sample and polyclonal goat anti-human IgG antibody are transferred to the second reaction tube. Anti-IgG antibodies block remaining IgG from patient sample from binding to the CMV antigen on the bead. During this time, IgM specific for CMV in the sample binds to CMV antigen on the bead. Unbound sample and reagent are then removed by centrifugal washes.

In the third cycle, the enzyme conjugated polyclonal goat anti-human IgM antibody is added to the second reaction tube. The enzyme conjugate binds to immobilized IgM to form the antibody sandwich complex. The unbound enzyme conjugate is removed by centrifugal washes. Finally, chemiluminescent substrate is added to the reaction tube and the signal is generated in proportion to the bound enzyme.

Incubation Cycles: 3 × 30 minutes

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Samples that are cloudy or have particulate material should be clarified by low-speed centrifugation.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 CMV IgM has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 10 µL serum or plasma (EDTA or heparinized)

Storage: 3 days at 2–8°C or 6 months at –20°C.¹²

No predilution is required for patient samples.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.¹³⁻¹⁵

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

**H319,
H302 + H312,
H412**

**P280, P264,
P273, P305 +
P351 + P338,
P301 + P312,
P302 + P312,
P501**

Warning! Causes serious eye irritation. Harmful if swallowed or in contact with skin. Harmful to aquatic life with long lasting effects. Wear protective gloves/protective clothing/eye protection/face protection. Wash hands thoroughly after handling. Avoid release to the environment. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. IF ON SKIN: Call a POISON CENTER or doctor/physician if you feel unwell. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: tetrasodium ethylene diamine tetraacetate, aminocaproic acid; CMV IgM Reagent Wedge
Contains: sodium azide; CMV IgM Adjustor, CMV IgM Controls

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

The CMV IgM Adjustor and Controls supplied with the kit should be fully dissolved after reconstitution. Failure to ensure the homogeneity of the solution may result in poor reproducibility of results.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

The bead is coated with *inactivated* CMV antigen. However, caution is advised because of the possible presence of residual virus when working with, or disposing of, the materials supplied.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

CMV IgM Bead Pack (L2CM12)

With barcode. 200 beads, coated with inactivated, purified CMV antigen (strain AD-169 from infected cell lysates). Stable at 2–8°C until expiration date.

L2KCM2: 1 pack

CMV IgM Reagent Wedge (L2CMA2)

With barcode. Reagent A containing 17.5 mL of a buffer solution with polyclonal goat anti-human IgG, in buffer. Reagent B containing 11.5 mL of a buffer solution with polyclonal goat anti-human IgG, in buffer. Reagent C containing 11.5 mL of alkaline phosphatase (bovine calf intestine) conjugated to polyclonal goat anti-human IgM antibody, in buffer. Stable at 2–8°C until expiration date.

L2KCM2: 1 wedge

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

CMV IgM Adjustor (L2CMJ3)

One vial of lyophilized human serum with IgM reactive to CMV, with preservative. Reconstitute with **2.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. (No further dilution is required.) Stable at 2–8°C for 14 days after reconstitution, or for 6 months (aliquotted) at –20°C.

L2KCM2: 1 vial

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

CMV IgM Controls (LCMC1, LCMC2)

LCMC1: (Negative Control): One vial of lyophilized human serum with IgM nonreactive with CMV, with preservative.

LCMC2: (Positive Control): One vial of lyophilized human serum with IgM reactive with CMV, with preservative. Reconstitute each vial with **2.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 14 days after reconstitution, or for 6 months (aliquotted) at –20°C.

L2KCM2: 1 set

Refer to the control insert for the ratio ranges.

No predilution is required for CMV IgM Adjustor and CMV IgM Controls.

Before running adjustors or controls, place the appropriate Aliquot Labels (supplied with the kit) on test tubes, so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

Also Required

Distilled or deionized water; test tubes; controls

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
2 weeks

Quality Control Samples: The controls supplied with the kit should be used as quality control material to monitor assay performance.

For the current control ratio ranges, please refer to the Control insert.

Calculation of Cutoff and S/CO Ratio:

The Master Cutoff of the assay was determined from representative samples to achieve optimal sensitivity and specificity for the assay.

The cutoff is set equal to the average counts per second (mean cps) of the Adjustor (from the most recent adjustment) multiplied by Curve Parameter 1. (See the “Low Adjustor CPS” and “Curve Parameter 1” fields in the IMMULITE 2000 Kit Information screen, which can be accessed from the menu via Data Entry: Kit Entry.)

Calculation of a signal/cutoff (s/co) ratio is done by using the following formula:

$$\text{S/CO Ratio} = \frac{\text{Sample or Control cps}}{\text{Mean Adjustor cps} \times \text{P1}}$$

Calculation and reporting of qualitative results (reactive / nonreactive / indeterminate) and signal/cutoff (s/co) ratio are handled automatically by the IMMULITE 2000.

The result for a sample is reported as “Indeterminate” if the counts per second for that sample fall within $\pm 10\%$ of the cutoff. The result is reported as “Reactive” if the sample’s counts are *above* the indeterminate range, and “Nonreactive” if *below* this range.

Interpretation of Results

The cutoff of the IMMULITE 2000 CMV IgM assay was determined with reactive and nonreactive patient samples by a ROC analysis, with a balanced consideration of sensitivity and specificity.

A result of “**Reactive**” (ratio of ≥ 1.1) indicates that the patient sample is reactive and that CMV IgM antibodies were detected in the sample.

A result of “**Nonreactive**” (ratio of < 0.9) indicates that the patient sample is nonreactive and that CMV IgM antibodies were not detected in the sample.

Any result of “**Indeterminate**” (ratio between 0.9 and < 1.1) should be retested. Samples which still test as “Indeterminate” should be tested by an alternate method, or a second sample should be taken — if possible — within a reasonable period of time (e.g., one week).

The presence of IgM antibodies to CMV virus is an indication of recent exposure to the virus.

The magnitude of the measured results (cps) above the Cutoff is not indicative of the total amount of antibodies detected.

The results reported by the laboratory to the physician should include: “The following results were obtained with the IMMULITE 2000 CMV IgM EIA. Results obtained from other manufacturers’ assay methods may not be used interchangeably.”

Expected Values

Individuals with active CMV infection will exhibit detectable levels of IgM antibody at the time of the onset of symptoms. Levels of IgM antibodies to CMV will subsequently decline and are undetectable after several months.

Limitations

The results of the test must be taken within the context of the patient’s clinical history, symptomology and other laboratory findings.

The presence of a primary infection, or reactivation of a past infection, can also be determined by isolation of the virus in tissue culture.^{3,9}

For the determination of seroconversion from *nonreactive* to *reactive*, a second serum sample should be drawn three to four weeks after the acute, nonreactive sample, during the convalescent stage of the infection. The acute phase sample should be stored and tested in parallel with the convalescent sample.

Samples containing antinuclear antibodies or other anticell antibodies may give false reactive results.

A rise in CMV antibody level may occur in patients with measles, varicella-zoster virus (VZV) or herpes simplex virus (HSV) due to antigenic crossreactivity within the herpes virus family. Patients with acute EBV (Epstein-Barr Virus) infection may show increased reactivity to other viruses, including CMV.

The results in HIV patients, in patients undergoing immunosuppressive therapy, or in patients with other disorders leading to immunosuppression, should be interpreted with caution.

The performance characteristics of this assay have not been established for use with specimens from neonates, cord blood, or pretransplant patients.

Antibodies to human IgG have been added to the reagents to remove specific IgG and rheumatoid factors which may cause false reactive responses.

Heterophilic antibodies in human serum/plasma can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed as a signal-to-cutoff ratio. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Specificity: A study was conducted to evaluate whether the measurement of cytomegalovirus IgM antibody is affected by closely related microorganisms. Sera containing antibodies to rheumatoid factor ($n = 10$), syphilis ($n = 10$), *Toxoplasma* ($n = 20$), Epstein-Barr virus ($n = 10$), varicella zoster virus ($n = 10$) and rubella ($n = 10$) were tested by IMMULITE 2000 CMV IgM. Sixty-six out of 70 samples yielded nonreactive results. Two out of 70 samples (1 syphilis and 1 VZV) yielded indeterminate results, and 2 out of 70 samples (1 *Toxoplasma* and 1 EBV) yielded positive results. These 4 samples were further tested by a commercially available CMV IgM assay from VIDAS and yielded 2 equivocal results (1 syphilis and 1 *Toxoplasma*), and 2 confirmed positive results (1 EBV and 1 VZV).

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 522 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 40 volunteers into plain, heparinized, EDTA and Becton Dickinson SST[®] plastic vacutainer tubes. 10 of the matched samples were spiked with various concentrations of CMV IgM, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 CMV IgM procedure.

(Heparin) = 0.98 (Serum) – 0.024
 $r = 0.997$

(EDTA) = 1.08 (Serum) – 0.055
 $r = 0.988$

(SST) = 0.99 (Plain Tubes) – 0.020
 $r = 0.996$

Means:
0.75 (Serum)
0.79 (Heparin)
0.77 (SST)
0.74 (EDTA)

Clinical Performance: Clinical investigations evaluated the performance of the IMMULITE 2000 CMV IgM assay in comparison to VIDAS CMV IgM (bioMerieux, Marcy d'Etoile, France). The study included a total of 523 samples from pregnant women, transplant patients, and patients with potential crossreacting diseases or conditions; samples submitted for routine CMV serology testing; and samples from seroconversion panels.

Further analysis including CMV IgM/IgG immunoblot and Epstein-Barr (EBV) virus testing was performed on samples initially (retrospectively) identified as CMV IgM reactive. Of the 116 samples tested, 69 were identified as primary CMV infection, 13 as latent CMV infection, 24 as primary, latent or reactivated EBV infection, and 10 were resolved as CMV IgM negative.

The overall resolved agreement between IMMULITE 2000 and VIDAS is presented in Table 1. More detailed comparison of IMMULITE 2000 and VIDAS results for the samples with potential crossreacting agents is presented in Table 2.

Table 1

Sample Type	VIDAS									Total
	R			I			NR			
	IMMULITE 2000									
	R	I	NR	R	I	NR	R	I	NR	
Pregnant	2	1	0	0	0	2	0	1	152	158
Transplant patients	29	0	1	1	0	1	0	0	5	37
Total CMV IgM negative	0	1	5	1	0	3	0	0	130	140
Primary CMV infection	65	1	1	0	0	1	0	0	0	68
Latent CMV infection	2	1	8	0	0	2	0	0	0	13
BBI sero-conversion panels	5	0	0	0	0	0	0	0	10	15
Total	103	4	15	2	0	9	0	1	297	431
Crossreactives (see Table 2)	13	1	7	0	0	8	0	0	63	92
Overall Total	116	5	22	2	0	17	0	1	360	523

R = Reactive
I = Indeterminate
NR = Nonreactive

All samples excluding crossreactive cases:
Overall agreement: 92.8% (400/431); Positive agreement: 84.4% (103/122); Negative agreement: 99.7% (297/298)

Pregnant cases:

Overall agreement: 97.5% (154/158);
Positive agreement: insufficient number of samples;
Negative agreement: 99.3% (152/153)

Transplant cases:

Overall agreement: 91.9% (34/37);
Positive agreement: 96.7% (29/30);
Negative agreement: 100% (5/5)

Primary CMV cases:

Overall agreement: 95.6% (65/68);
Positive agreement: 97.0% (65/67);
Negative agreement: insufficient number of samples

Table 2: Crossreactives

Disease state	Total n	IMMULITE 2000:	VIDAS
Rheumatoid Factor (+)	36	34/36 nonreactive	34/36 nonreactive
<i>Toxoplasma</i> IgM* (+)	8	7/8 nonreactive	6/8 nonreactive
Primary, latent or reactivated EBV*	47	37/49 nonreactive	23/49 nonreactive

* One patient was positive for both *Toxoplasma* IgM and EBV

References

1) Drew WL. Herpesviridae: cytomegalovirus. In: Lennette EH, Halonen P, Murphy FA, editors. Laboratory diagnosis of infectious diseases: principles and practice, volume II, viral, rickettsial and chlamydial diseases. New York: Springer-Verlag, 1988: 247-60. 2) Ho M. Characteristics of cytomegalovirus. In: Greenough WB, Merigan TC, editors. Cytomegalovirus biology and infection: current topics in infectious disease. New York: Plenum, 1982: 9-32. 3) Hodinka RL, Friedman HM. Human cytomegalovirus. In: Balows A, Hausler, Jr WJ, Hermann KL, Isenberg HD, Shadomy HJ. Manual of clinical microbiology, 5th ed. Washington, DC: American Society for Microbiology, 1991: 829-37. 4) Reynolds DW, et al. Laboratory diagnosis of cytomegalovirus infections. In: Lennette EH, Schmidt NJ, editors. Diagnostic procedures for viral, rickettsial and chlamydial infections, 5th ed. Washington: American Public Health Association, 1979: 399-439. 5) Van der Bij W, et al. Antigen test for early diagnosis of active cytomegalovirus infection in heart transplant recipients. J Heart Transplant 1988;7:106-9. 6) Zaia JA. Epidemiology and pathogenesis of cytomegalovirus disease. Seminars in Hematology 1990;27:5-10. 7) Yow MD. Congenital cytomegalovirus disease: a now problem. J Infect Dis 1989;159:163-7. 8) Miller H, et al. Prospective study of cytomegalovirus

antigenemia in allograft recipients. J Clin Microbiol 1991;29:1054-5. 9) Wiedbrauk DL, Johnston SLG. Manual of clinical virology. New York: Raven Press, 1993: 82-91. 10) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 11) Baltz ML, Searcy RL. Clinical significance and advanced serologic diagnosis of ToRCH infections. Am Clin Lab 1994;March/April:18-23. 12) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: WB Saunders, 1995:358. 13) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377-82, 387-8. 14) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 15) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

Technical Assistance

Available outside the United States only. For technical assistance, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2019 Siemens Healthcare Diagnostics. All rights reserved.
CONTROL AB En: Control Antibody