

## NEUTRAL DISINFECTANT DETERGENT FOR THE PRE-DISINFECTION OF MEDICAL DEVICES

### INSTRUCTIONS FOR USE

#### Preparation of the solution at 0.5%:

- Wear gloves, goggles and apron
- Dilute 25 ml of DDN 9 in 5L of water and homogenize the solution

#### Use:

- 1- Soiled instruments have to be fully immersed
  - 2- Scrub medical devices and use a swab if necessary
  - 3- Contact time : 15 minutes
  - 4- Rinse thoroughly
- Renew the solution after use*

### DISINFECTING PROPERTIES

Activity	Norms	Concentration	Contact time
<b>Bactericidal</b>	EN 13727	0.5%	5 minutes
	EN 14561	0.5%	15 minutes
	EN 13697	0.5%	5 minutes
	EN 1276	0.25%	5 minutes
<b>Yeasticidal</b>	EN 13624	0.5%	5 minutes
	EN 14562	0.5%	15 minutes
	EN 1650	0.25%	5 minutes
	EN 13697	0.5%	5 minutes
<b>Virucidal</b>	EN 14476 Virus PRV (hepatitis B model), EN 14476 + A1	0.5%	5 minutes
	Virus DBV (hepatitis C model), Virus VIH	0.5%	15 minutes
	EN 14476 + A2	0.5%	15 minutes
	Herpes Virus	0.5%	5 minutes
	EN 17111	0.5%	15 minutes
	Vaccinia Virus	0.5%	15 minutes

### PRECAUTIONS OF USE

See (SDS 010449)

### FEATURES

- Form: **Liquid**
- Aspect : **Opalescent**
- pH at 0.5%: 7
- Expiration date: **2 years**

### COMPOSITION

- Non-ionic surfactants
- Sequestrants
- Quaternary ammonium

### PACKAGING

- ① Box of 250 doses of 25 ml.....ref. **1044947**
- ② Box of 12 x1 L.....ref. **1044916**
- ③ Box of 4 x5 L.....ref. **1044920**



Hazardous product - Follow precautions of use. Always read label and product information before use. Medical device Class IIb.

# DDN 9

NEUTRAL DISINFECTANT DETERGENT FOR THE  
PRE-DISINFECTION OF MEDICAL DEVICES



## ■ AREAS OF APPLICATION

- Operating room, endoscopy, ophthalmology and contactology

## ■ PRODUCT ADVANTAGES

- Proven efficacy on *Pseudomonas aeruginosa* biofilm
- Instrumentation integrity guaranteed for 72 hours in immersion-bath
- Compatible with a wide range of materials
- Can be used in ultrasounds vats
- Effective on enveloped virus (including corona virus)

**Fb Franklab®**  
notre expertise l'Ultra-Propreté



CE  
0459

