



Sterilization item Guide

2020.12.30

CS

Contents Table

- ◆ Material to sterilization
- ◆ Lumen Guideline
- ◆ Cautions
- ◆ Device that can be sterilized
- ◆ Substances that cannot be sterilized
- ◆ Non-woven packaging method
- ◆ Pouch packing method

Material to sterilization

Sterilizable material

- ▶ Metal products (stainless steel, aluminum, titanium alloy, etc.)
- ▶ Plastic products (polyethylene, polycarbonate, polyvinyl chloride (PVC), polypropylene (PP), Viton, etc.))
- ▶ Glassware, Rubber products (silicon)
- ▶ Ceramic products

Prohibited material

- ▶ Hard-Lumen¹: units having a pipe with a diameter of $<1\varnothing$ and a length of more than 1,000mm
- ▶ Soft-Lumen¹: units having a tube with a diameter of $<1\varnothing$ and a length of more than 12,000mm
- ▶ Soft-Lumen¹: Units having a tube with a diameter of $<2\varnothing$ and a length of more than 1,500mm, blocked on one side
- ▶ Substances that absorb moisture¹ (paper, Cellulose, Cotton, Gauze, Wood, etc.)
- ▶ Liquid (oil, water, foam, etc.), equipment with water
- ▶ Powder, copper and zinc, Brass
- ▶ Acetal, Buna N(Nitrile), Hypalon (glove material)

Material to sterilization

Materials that need attention when sterilizing

- ▶ Disposable products that are not recommended (disposable tarpaulin, gown, syringe, etc.)
- ▶ Items damaged by vacuum such as marker pens
- ▶ Polyurethane, nylon and poor quality plastics¹
- ▶ Sponges (styrofoam, sponge for microscope protection, etc.)
- ▶ In the case of some painted products such as aluminum surface painting, it occurs when discoloration or color change
- ▶ Some cable products¹
- ▶ Some cables may stick to transparent vinyl of sterilized roll and pouch. Wrap it in sterilized wrap (non-woven) and wrap it.



<Cable item>



< Cable wrapped in non-woven sheet >

※ ¹ Materials with limited sterilization: Porous materials, absorbent materials and lumens interfere with vacuum and absorb hydrogen peroxide, a sterilant, to reduce sterilization effect..

Lumen Guideline

RENO-S20

Sterilization cycle	Lumen	S20	Max quantity
Eco Cycle	SUS Ø1	150mm	2ea
	Teflon Ø1	400mm	2ea

Sterilization cycle	Lumen	S20	Max quantity
Advanced Cycle	SUS Ø1	400mm	5ea
	Teflon Ø1	1,000mm	5ea
	Teflon Ø2	1,500mm (One-side closed)	1ea


Maximum of 2 lumens per load in ECO CYCLE

Maximum of 5 lumens per load in ADVANCED CYCLE

Maximum of 10 containers per load In ECO / ADVANCED CYCLE

Maximum of 2 silicon pads per load in ECO / ADVANCED CYCLE

※ Silicon pad size criteria: 250mm x 150mm



Validation load weight (With tray)

- Advanced cycle: 2.6kg
- Eco cycle: 1.5kg

Lumen Guideline

RENO-S30

Sterilization cycle	Lumen	S30	Max quantity
Eco Cycle	SUS Ø1	200mm	10ea
	Teflon Ø1	400mm	10ea
	Flexible Endoscope Teflon Ø1	900mm 1ea with camera and light cable	1set

Sterilization cycle	Lumen	S30	Max quantity
Advanced Cycle	SUS Ø 0.7	500mm	10ea
	Teflon Ø1	1,000mm	10ea
	Teflon Ø2	1,500mm (One-side closed)	2ea

Maximum of 10 lumens per load in ECO/ ADVANCED CYCLE

Maximum of 10 containers per load In ECO / ADVANCED CYCLE

Maximum of 2 silicon pads per load in ECO / ADVANCED CYCLE

※ Silicon pad size criteria: 250mm x 150mm



Validation load weight (With tray)

- Advanced cycle: 3.5kg
- Eco cycle: 2.1kg

Lumen Guideline

RENO-D50

Sterilization cycle	Lumen	S30	Max quantity
Eco Cycle	SUS Ø1	200mm	10ea
	Teflon Ø1	400mm	10ea
	Flexible Endoscope Teflon Ø1	900mm 1ea with camera and light cable	1set

Sterilization cycle	Lumen	S30	Max quantity
Advanced Cycle	SUS Ø 0.7	500mm	10ea
	Teflon Ø1	1,000mm	10ea
	Teflon Ø2	1,500mm (One-side closed)	2ea

Maximum of 10 lumens per load in ECO/ ADVANCED CYCLE

Maximum of 10 containers per load In ECO / ADVANCED CYCLE

Maximum of 2 silicon pads per load in ECO / ADVANCED CYCLE

※ Silicon pad size criteria: 250mm x 150mm



Validation load weight (With tray)

- Advanced cycle: 3.5kg
- Eco cycle: 2.1kg

Lumen Guideline

RENO-S90

Sterilization cycle	Lumen	S90	Max quantity
Eco Cycle	SUS Ø1	200mm	10ea
	Teflon Ø1	400mm	10ea
	Flexible Endoscope Teflon Ø1	900mm	2ea
	Flexible Endoscope Teflon Ø1	900mm 1ea with camera and light cable	1set

Sterilization cycle	Lumen	S90	Max quantity
Advanced Cycle	SUS Ø1	500mm	10ea
	Teflon Ø1	1,000mm	10ea
	Teflon Ø2	1,500mm (One-side closed)	2ea

Maximum of 10 lumens per load in ECO/ ADVANCED CYCLE

Maximum of 10 containers per load In ECO / ADVANCED CYCLE
 Maximum of 2 silicon pads per load in ECO / ADVANCED CYCLE
 ※ Silicon pad size criteria: 250mm x 150mm

Validation load weight (With tray)

- Advanced cycle: 8kg
- Eco cycle: 5.5kg
- Non lumen cycle: 3.5kg

Lumen Guideline

RENO-S130

Sterilization cycle	Lumen	S130	Max quantity
Eco Cycle	SUS Ø1	200mm	10ea
	Teflon Ø1	400mm	10ea
	Flexible Endoscope Teflon Ø1	900mm	2ea
	Flexible Endoscope Teflon Ø1	900mm 1ea with camera and light cable	1set

Sterilization cycle	Lumen	S130	Max quantity
Advanced Cycle	SUS Ø0.7	500mm	10ea
	Teflon Ø1	1,000mm	10ea
	Teflon Ø2	1,500mm (One-side closed)	2ea
	Teflon Ø1	1,000mm (Without load)	20ea

Maximum of 10 containers per load In ECO / ADVANCED CYCLE
 Maximum of 2 silicon pads per load in ECO / ADVANCED CYCLE
 ※ Silicon pad size criteria: 250mm x 150mm



Validation load weight (With tray)

- Advanced cycle : 12kg
- Eco cycle : 8kg
- Non lumen cycle : 5kg

Lumen Guideline

RENO-S130D

Sterilization cycle	Lumen	S130D	Max quantity
Eco Cycle	SUS Ø1	200mm	10ea
	Teflon Ø1	400mm	10ea
	Flexible Endoscope Teflon Ø1	900mm	2ea
	Flexible Endoscope Teflon Ø1	900mm 1ea with camera and light cable	1set

Sterilization cycle	Lumen	S130D	Max quantity
Advanced Cycle	SUS Ø1	500mm	10ea
	Teflon Ø1	1,000mm	10ea
	Teflon Ø2	1,500mm (One-side closed)	2ea
	Teflon Ø1	1,000mm (Without load)	20ea

Maximum of 10 containers per load In ECO / ADVANCED CYCLE
 Maximum of 2 silicon pads per load in ECO / ADVANCED CYCLE
 ※ Silicon pad size criteria: 250mm x 150mm



Validation load weight (With tray)

- Advanced cycle : 12kg
- Eco cycle : 8kg
- Non lumen cycle : 5kg

Cautions


Before loading the instruments

- ▶ Check if instruments includes prohibited materials.
- ▶ Check if there are prohibited instruments loaded.
- ▶ Follow the lumen guideline.
- ▶ The recommendation amount is less than 80% loading.
However it can be different according to the instrument kinds.
- ▶ Check the instrument compatibility.
Some instrument are not compatible with Plasma sterilizer such as container only for steam sterilizer with out holes.
- ▶ The quantity of the Silicon pad must not be over 2ea.
※ Silicon pad size criteria: 250mm x 150mm

Cautions

Rigid Endoscopes (Arthroscope, Laparoscope, Telescope, etc.)

**Moisture removal
caution items**



Laparoscopic Forceps

**Moisture removal
caution items**



Cautions

Light Cables

Moisture removal
caution items



Camera Head & Camera Cables

Moisture removal
caution items



* In the case of cables, scratches or holes on the surface may have an effect on the vacuum process.

Cautions

Bipolar forceps & Bipolar codes



Trocar (Reusable)



Cautions

Micro Instruments set



* In the case of silicon material, the more the amount is used, the more it needs to be managed due to the surface change.

Tubes, lines

Moisture removal caution items



Cautions

Drill & Battery

Moisture removal caution items



Arthroscopy Set



Cautions

Laparoscopy Set



Endo Urology Set



Moisture removal caution items

Device that can be sterilized



Cable Code



Electrosurgical Cable
(Bipolar code)



Pincete



Wire Guide



Twister Needle Holder
TC



PLIERS

Device that can be sterilized



Electric Coagulation Forceps



Electro surgical Cable



Stainless steel device



Biolar Forceps



Forcep



GS Micro

Device that can be sterilized



Skin Stapler



Endo Clip



Cryo-probes



Power Drills



Micro eye surgery
instrument



TUR Set

Device that can be sterilized



Micro Forceps



Forward/ reverse scanlan



Light Cable



Air drill code



Chest Tube



Sucker

Device that can be sterilized



Lonstar ret & Hook



Co2 line



Silicon Tube



Arthro suture set



Drill



Face Cover

Device that can be sterilized



Air drill



Remover



Trocar



Rectal Tube



Laparoscopic Trocar



Biopsy forceps

Device that can be sterilized



Optical Sugical Instrument



Cataract Instrument



Optical Sugical Instrument



Handle of surgical light



Arthroscope



Laparoscope

Substances that cannot be sterilized



wood, pulp material



Gauze, cotton, linen



Sponge, styrofoam,
foam material



Patient clothes,
gowns, bed sheets



Liquid, powder,



Paper



Gel, oil, clay



Disposable

* 일회용품 같은 경우 재사용이 제한되어 있어 멸균 불가능함.

Non-woven packaging method

- 1 After preparing a tray suitable for medical tools as shown in the picture, Put the sterilized product and CI strip in the tray.




- 2 Wrap the tray as shown in the picture using non-woven sheet. (After packing, fix with CI tape)




Non-woven packaging method

② Tray packing flow chart



Pouch packing method

1 Pouch packing method



Thank you

Expert in Infection Control
